

ISMATEC® Laboratory Pumps

2009 Catalog

For Dispensing and Filling Applications • 0.4 µl/min to 13 liters/min

Featuring Tubing, Gear and
Piston Pumps for:

Analytics
Sampling
Biotechnology
Process Control



ISMATEC, established in 1964

ISMATEC®
IDEX Health & Science

Precision engineered pumps

The high accuracy and long service-life of ISMATEC® pumps are well known worldwide. For more than 40 years, ISMATEC has worked closely with laboratory and production customers to develop state-of-the-art pumping solutions for applications from highly sensitive cellular analyses to industrial dispensing and filling. You can count on ISMATEC to provide individual consultation to match the right pump to your application, and to provide the service you need after the sale.

As part of Micropump and the IDEX Health and Science, ISMATEC now offers an even broader array of industry-leading products and engineering expertise to meet your fluidic and gas management needs. To learn more, please contact us at the locations indicated below, or visit our website at www.idex-hs.com

Objectives for quality

By providing competent service and products of high quality we contribute to our customers' success.

ISO 9001:2000 and 14001 certified



United States

ISMATEC - Micropump Headquarters
IDEX Health & Science
1402 NE 136th Ave.
Vancouver, WA 98684
Tel. +1 360 253 2008
Fax +1 360 253 8294

Asia

IDEX Health & Science
Room 1701A, CBD International Plaza
No.16 Young An Dong Li
Chao Yang District
Beijing, China 1000022
Tel. +86 10 6566 9090
Fax +86 10 6567 5152



Germany

ISMATEC European Headquarters
IDEX Health & Science
ISMATEC Laboratoriumstechnik GmbH
Futtererstrasse 16
D-97877 Wertheim-Mondfeld
Tel. +49 93 77 / 92 03 0
Fax +49 93 77 / 13 88
office.ismatec@idexcorp.com
www.ismatec.de

Switzerland

ISMATEC Labortechnik SA
IDEX Health & Science
Feldeggstrasse 6
CH-8152 Glattbrugg-Zürich
Tel. +41 44 874 94 94
Fax +41 44 810 52 92
www.ismatec.ch

United Kingdom

Micropump Europe
IDEX Health & Science
Howard Road, Eaton Socon
St. Neots, Cambridgeshire
England
Tel. +44 1480 356600
Fax +44 1480 356920

www.ismatec.com
www.idex-hs.com/ismatec



Please ask also for our OEM pump brochure or contact us.

Our engineers welcome your ideas.



Tubing Pumps for Multi-Purpose Use

Overview and Table of Contents
 IPC / IP / IPC-N / IP-N, REGLO, REGLO *Quick*, and Ecoline Pumps
 MCP / BVP Drives and Pump-heads
 MS-CA Panel Mounted Pumps
 Flowmaster Pump
 Tubing

Page 4
 Page 6
 Page 16
 Page 26
 Page 28
 Page 30



Gear Pumps for Pulsefree Pumping

Overview and Table of Contents
 Pump-heads
 REGLO / MCP / BVP Series

Page 40
 Page 41
 Page 44



Rotary Piston Pumps for Corrosive Media

Overview and Table of Contents
 REGLO Series
 MCP Series

Page 50
 Page 52
 Page 56



Miscellaneous

Cassettes, Foot Switch, Software, Fittings
 OEM Pumps

Page 60
 Page 64



Select the pump type that meets your requirements

Selection criteria	Tubing Pumps	Gear Pumps	Rotary Piston Pumps
Flow rate min. to max.	<0.001 ml/min to 13 l/min	1 ml/min to 12 l/min	0.025 ml/min – 2.3 l/min
Number of channels	1 – 24	1	1
Differential pressure	max. 2.5 bar ³	max. 5.6 bar	max. 6.9 bar
Suction lift (water)	7 to 8 m	<1 m	~5 m
Dead volume	practically none	5 to 45 ml	very small
Chemical resistance	depends on tubing material	high	very high
Accuracy and Repeatability	high	high ¹	very high
Self-priming	✓	possible ⁴	possible
Sensitive to dry-running	no	yes	yes
Syphoning effect	no	yes	no
Pumping gently = low shearing forces	✓	no	no
under sterile conditions	✓	no	no
in both directions	✓	no	✓
pulse-free	²	✓	²
contamination-free	✓	no	no
Media containing particles	✓ very good	no	✓ max. 0.8 mm Ø
viscous	✓ very good	possible	✓ good
containing living cells	✓ very good	no	no
foaming	✓	no	no
corrosive / aggressive	³	✓ good	✓ very good
gas	³	no	⁴

¹ Requires non-return valve

³ Depends on the tubing material

² Pumping with low pulsation possible; depends on the pump-head


⁴ Depends on the pump-head

ISMATEC® pumps are well known for:

- Excellent quality standard
- A maximum of performance and reliability over many years of service
- Extremely low amortization costs

The multi-purpose tubing pump

Safe and easy to use



Safe and easy to use

- Developed for continuous duty, 24 hrs./day, 7 days/week
- **No contamination of the fluid or the pump**
- **Immune to dry-running**
- No valves which clog, hang up or need servicing
- No seals which can leak
- No syphoning effect when the pump stops
- **No flushing and cleaning required**
- Immediately adaptable to new application
- New tubing = pump is ready for use
- **Fast filling and emptying thanks to a MAX key and reversible flow direction**

Application range of tubing pumps

Industries	Applications	Special media
<ul style="list-style-type: none"> - Beverage - Biotechnology - Chemical - Environmental - Food - Ink - Paint - Pharmaceutical - Plating - Printing 	<ul style="list-style-type: none"> - Analytical - Fermentation - Filtration and separation - Galvanic process - Research and development - Water and waste water treatment 	<ul style="list-style-type: none"> - Abrasives - Cell cultures - Corrosives - Paint and pigments - Printing inks - Solid content up to 60%

For a wide range of applications

The tubing pump is capable of pumping almost any media, e.g. highly viscous media or media containing a solid content of up to 60%, as well as corrosive or sterile media.

- More than 10 different tubing materials are available.
- Multi-channel pumps (max. 24 channels) **delivers independently yet simultaneously** various media at **different flow rates**
- Self-priming from a depth of up to 8 m
- Reversible rotation direction
- **Almost pulsefree pumping with 12 roller-pumpheads**
- Extremely accurate dispensing with calibrate-able, microprocessor controlled drives

Inexpensive to maintain

- **Practically no maintenance costs**
- The pumps are almost maintenance-free
- The tubing is the only part to wear
- No expensive seals, valves or diaphragms
- Easy to use
- Simply change the tubing and the pump is ready again

High quality and precision guarantee an optimum performance even after many years of intensive use




The tubing cassettes, an ISMATEC invention, have been continuously refined.

The result:
Click'n'go – the cassette with the automatic pressure setting.

▶ When it comes to versatility, there are no other pump systems which can compete with ISMATEC® tubing pumps.

Unique!

ISMATEC® pumps feature:



Convex rollers

- Treat the liquid gently (e.g. living cells)
- Improve the delivery stability
- Increase the repeatability
- Guarantee optimum tube centering

The tube is progressively closed, starting from the center outwards.



Changing the tube within 5 seconds

- The Flowmaster pump-head (see Page 28), with flow rates up to 13 l/min
- Disassembling the pump-head in less than 1 minute




Carrying out programs independently of a PC!

- Create the application profile in the PC (with ProgEdit software, Page 61)
- Download the file data into the pump memory
- Disconnect the pump from the PC
- Carry out your application on the spot, using the pump as a stand-alone unit



Unique coupling system

- The MCP and BVP drives enable the user to mount a wide range of different pump-heads
- The pump-heads are interchangeable within seconds without using a tool
 - More than 14 peristaltic pump-heads are available (1 with PTFE tubing)



UNIQUE!

- Convex rollers and concave tub-bed.
- Pump-heads with this sign are ideal for cell and media sensitive pumping.



Note

All microprocessor controlled drives are LabVIEW compatible and can easily be integrated into process control systems.

Selection by flow rate and model

Flow rate ml/min		Channel	Model	Page
min.	max.			
0.0004	44	4–24	IPC, IP, IPC-N, IP-N	6
0.001	68	2–4	REGLO Analog/Digital	8
2.1	230	1	REGLO Quick	10
0.003	5400	1–8	Ecoline	12

Versatile and flexible, the MCP/BVP system with instantly interchangeable pump-heads

Drives	Model	Page
	BVP Standard, MCP Standard	14
	BVP Process, MCP Process	16
		18

Pump-heads

Flow rate ml/min		Channel	Model	Page
min.	max.			
0.45	3700	1	Pro-280/281, Pro-380/381	20
0.072	3600	1	360, 380, 380AD	22
0.07	1100	1	Easy-Load®, Easy-Load II	23
0.07	45	1	PTFE-Tubing	23
0.002	230	4–12	CA-4, CA-8, CA-12	24
0.09	1100	2–3	SB 2V, SB 3V	24
0.001	100	4–24	MS/CA4-12, MS/CA8-6	25

Flow rate ml/min

Flow rate ml/min		Channel	Model	Page
min.	max.			
0.021	26	2–4	MS-CA	26
37	13'000	1	Flowmaster	28

Tubing

Rating Comparison	30
Material Properties	32
Special tubing for aggressive media and high pressure	34
Ordering Information	36
Tube connectors	38



This mark indicates dispensing functions

(You'll find tubing dispensing pumps on Pages 6, 8, 17, 19)

- Pumping by speed or flow rate
- Dispensing by volume or time
- Interval dispensing with a pause
- Dispensing a volume within a pre-set time
- Interval dispensing with a pre-set number of dispensing cycles
- **Calibrating** the flow rate and dispensing volume
- Roller back-steps for **drip-free dispensing**
- Factory-set tube sizes

1

IPC / IP and IPC-N / IP-N

Lowest pulsation, highest accuracy
Ideal for analytical applications

- Very high repeatability on all channels
- Planetary drive



Tablet dissolution

CE



Dispensing and calibrating function see Page 5

IPC, IPC-N with dispensing functions

- Microprocessor controlled
- Very accurate pumping
- With stand-by mode (prevents the tubing system from drying out, e.g. over night)

IPC (and IP)

0.002–44 ml/min (per channel)

IPC-N (and IP-N)

0.4 µl/min–11 ml/min (per channel)

- 4, 8, 12, 16 or 24 channels
- 8 actively driven stainless steel rollers
- Laquered stainless steel housing
- Membrane key-pad and LED display
- Click'n'go cassettes with automatic pressure setting
- Each channel can take various tubing sizes
- 2-stop tubing
- Differential pressure 1.0 bar depending on tubing material. Tubing with small i.d.'s and/or cassettes with pressure lever (Page 60) may enable higher pressures.



IP and IP-N

without dispensing functions

- Microprocessor controlled
- Very accurate pumping
- Flow rates as IPC and IPC-N

Ordering information

Model	Order No.	Model	Order No.	Flow rates	Channels	Speed
IPC / IPC-N	IPC / IPC-N	IP/IP-N	IP / IP-N	ml/min per channel		rpm
IPC 4	ISM 930	IP 4	ISM 940	0.002 – 44	4	0.4 – 45
IPC 8	ISM 931	IP 8	ISM 941	0.002 – 44	8	0.4 – 45
IPC 12	ISM 932	IP 12	ISM 942	0.002 – 44	12	0.4 – 45
IPC 16	ISM 933	IP 16	ISM 943	0.002 – 44	16	0.4 – 45
IPC 24	ISM 934	IP 24	ISM 944	0.002 – 44	24	0.4 – 45
IPC-N 4	ISM 935	IP-N 4	ISM 945	0.0004 – 11	4	0.11 – 11.25
IPC-N 8	ISM 936	IP-N 8	ISM 946	0.0004 – 11	8	0.11 – 11.25
IPC-N 12	ISM 937	IP-N 12	ISM 947	0.0004 – 11	12	0.11 – 11.25
IPC-N 16	ISM 938	IP-N 16	ISM 948	0.0004 – 11	16	0.11 – 11.25
IPC-N 24	ISM 939	IP-N 24	ISM 949	0.0004 – 11	24	0.11 – 11.25

Foot switch for IPC / IPC-N and IP / IP-N Order No. IS 10039 see Page 61

CA spare cassettes and adaptors Order No. see Page 60

LabVIEW driver download for free www.ismatec.com



Studies and simulation of the environmental influences of rain water in free-land plants.

- ### Applications
- Toxicological in-vitro use
 - Perfusion of animal tissue slices
 - Sampling from tablet dissolution systems

Interfaces



IPC, IPC-N
PC-controllable
Analog: same as IP, IP-N




IP, IP-N

- Speed control (0–5 or 0–10 V, 0–20 or 4–20 mA)
- Speed output (0–10 V or 0–11 kHz)
- Start/Stop
- Rotation direction
- Autostart

Flow rates and tubing

Model Channels Rollers	Tubing i.d. (mm)	IPC / IP		IPC-N / IP-N	
		4 / 8 / 12 / 16 / 24	8	4 / 8 / 12 / 16 / 24	8
Speed rpm		0.4	44.0	0.11	11.25
Tygon®ST R-3603/R-3607 Order No.		ml/min per channel min.	ml/min per channel max.	ml/min per channel min.	ml/min per channel max.
SC0188	0.13	0.002	0.15	0.0004	0.039
SC0002	0.25	0.005	0.41	0.001	0.10
SC0005	0.51	0.015	1.5	0.004	0.38
SC0008	0.76	0.032	3.2	0.009	0.81
SC0011	1.02	0.057	5.7	0.014	1.4
SC0014	1.22	0.079	7.9	0.020	2.0
SC0017	1.52	0.12	12	0.030	3.0
SC0020	1.85	0.17	17	0.043	4.3
SC0023	2.54	0.30	30	0.075	7.5
SC0222	3.17	0.44	44	0.11	11

Approx. values: determined with water, at 22°C, no differential pressure, Tygon tubing.



With the planetary drive system each roller is directly driven by the sun wheel. This prevents axial push-pull friction on the tubing. Result: increased service-life of the tubing, lower pulsation, high repeatability.

Ask for our Pump Tubing Selection Guide.
(see also Pages 30 to 39)



Specifications IPC and IPC-N

Motor type	DC motor
Speed IPC	0.4–44 rpm
Speed IPC-N	0.11–11.25 rpm
Speed setting	1–100 %, resolution 0.1%
Flow rate setting	µl/min or ml/min
Power consumption	30 W
Mains connection	230V _{AC} /50Hz, 115V _{AC} /60Hz adjustable
Protection rating	IP 30

Specifications IP and IP-N

Motor type	DC motor
Speed IP	0.4–44 rpm
Speed IP-N	0.11–11.25 rpm
Speed setting	1–100 %, resolution 0.1%
Speed setting IP	rpm, resolution 0.1 rpm
Speed setting IP-N	rpm, resolution 0.03 rpm
Power consumption	30 W
Mains connection	230V _{AC} /50Hz, 115V _{AC} /60Hz adjustable
Protection rating	IP 30

Dimensions / Weight

4 channels	
Depth/Width/Height (mm)	180x175x130
Weight (kg)	4.6
8 channels	
Depth/Width/Height (mm)	220x175x130
Weight (kg)	5.1
12 channels	
Depth/Width/Height (mm)	260x175x130
Weight (kg)	5.8
16 channels	
Depth/Width/Height (mm)	300x175x130
Weight (kg)	6.5
24 channels	
Depth/Width/Height (mm)	380x175x130
Weight (kg)	7.9

1

REGLO *Analog / Digital*

The smallest calibrateable dispensing pump
Footprint only 178x100 mm!

- Very low pulsation (model with 12 rollers)
- High repeatability
- 10 cm wide, 13.5 cm high



REGLO *Digital*
with dispensing functions
0.001–68 ml/min (per channel)
– Microprocessor controlled

REGLO *Analog / Digital*

- 2 or 4 channels
- Choice of 6, 8 or 12 rollers
- Click'n'go cassettes with automatic pressure mechanism (each channel can take various tubing sizes)
- 3-stop tubing
- Differential pressure 1.0 bar depending on tubing material. Tubing with small i.d.'s and/or cassettes with pressure lever (Page 60) may enable higher pressures.



REGLO *Analog*
without dispensing functions
0.002–68 ml/min (per channel)
– Variable speed drive

Specifications REGLO *Analog*

Motor type	DC motor
Speed	2-channel 3.2 – 160 rpm 4-channel 2.0 – 100 rpm
Speed setting	2–99%, resolution 1% 2-digit potentiometer
Power consumption	20 W
Mains connection	230V _{AC} /50Hz, 115V _{AC} /60Hz adjustable
Protection rating	IP 30
Depth/Width/Height	2-channel 178x100x143 mm 4-channel 190x100x143 mm
Weight	2-channel 2.0 kg 4-channel 2.1 kg

Specifications REGLO *Digital*

Motor type	DC motor
Speed	2-channel 1.6 – 160 rpm 4-channel 1.0 – 100 rpm
Speed setting	rpm, resolution 0.1 rpm
Flow rate setting	µl/min or ml/min
Power consumption	75 W
Mains connection	100–230V _{AC} /50–60Hz
Protection rating	IP 30
Depth/Width/Height	2-channel 178x100x135 mm 4-channel 190x100x135 mm
Weight	2-channel 2.0 kg 4-channel 2.1 kg

Interfaces



REGLO Analog
2-digit potentiometer
2–99%, resolution 1% (for speed setting)



REGLO Digital
6-button membrane key-pad, LED-display
Flow rate setting in $\mu\text{l}/\text{min}$ and ml/min

**REGLO Analog**

- Speed control (0–5 or 0–10 V, 0–20 or 4–20 mA)
- Speed output
2-channel: 0–8 kHz
4-channel: 0–5 kHz
- Start/Stop
- Rotation direction

**REGLO Digital**

- PC-controllable Analog: only speed output (see Reglo *Analog*), start/stop and autostart

Flow rates and tubing



Model	Channels	Rollers	REGLO Analog+Digital 2/6		REGLO Analog+Digital 2/8		REGLO Analog+Digital 2/12		REGLO Analog+Digital 4/6		REGLO Analog+Digital 4/8		REGLO Analog+Digital 4/12	
			Speed rpm	1.6 ¹	160	1.6 ¹	160	1.6 ¹	160	1.0 ¹	100	1.0 ¹	100	1.0 ¹
Tygon®ST R-3603/R-3607	Order No.	Tubing i.d. (mm)	ml/min per channel		ml/min per channel		ml/min per channel		ml/min per channel		ml/min per channel		ml/min per channel	
			min. ¹	max.	min. ¹	max.	min. ¹	max.	min. ¹	max.	min. ¹	max.	min. ¹	max.
SC0189		0.13	0.003	0.22	0.002	0.17	0.002	0.15	0.002	0.14	0.002	0.11	0.001	0.093
SC0050		0.25	0.008	0.76	0.007	0.65	0.007	0.61	0.005	0.48	0.005	0.41	0.004	0.38
SC0053		0.51	0.031	3.1	0.027	2.7	0.025	2.5	0.019	1.9	0.017	1.7	0.016	1.6
SC0056		0.76	0.067	6.7	0.058	5.8	0.053	5.3	0.042	4.2	0.036	3.6	0.033	3.3
SC0059		1.02	0.012	12	0.10	10	0.090	9.0	0.073	7.3	0.063	6.3	0.056	5.6
SC0062		1.22	0.16	16	0.14	14	0.12	12	0.10	10	0.088	8.8	0.075	7.5
SC0065		1.52	0.24	24	0.20	20	0.17	17	0.15	15	0.13	13	0.10	10
SC0068		1.85	0.34	34	0.28	28	0.21	21	0.21	21	0.17	17	0.13	13
SC0071		2.54	0.53	53	0.44	44	0.31	31	0.33	33	0.27	27	0.19	19
SC0224		3.17	0.68	68	0.57	57	0.38	38	0.43	43	0.35	35	0.24	24

Approx. values: determined with water, at 22°C, no differential pressure, Tygon tubing.

¹ Min. flow rate for REGLO Analog = 2 % of max. flow rate

Ask for our Pump Tubing Selection Guide (see also Pages 30 to 39).



Spare cassettes MS/CA
(see Page 60)

Ordering information (values in brackets are for REGLO Digital)

Model	Order No.	Order No.	Flow rates	Channels	Rollers	Speed rpm
	REGLO Analog	REGLO Digital	ml/min per channel			
MS-2/06	ISM 830	ISM 831	0.005 (0.003) – 68	2	6	1.6 (3.2)–160
MS-2/08	ISM 829	ISM 832	0.004 (0.002) – 57	2	8	1.6 (3.2)–160
MS-2/12	ISM 795	ISM 596	0.003 (0.002) – 38	2	12	1.6 (3.2)–160
MS-4/06	ISM 828	ISM 833	0.003 (0.002) – 43	4	6	1.0 (2.0)–100
MS-4/08	ISM 827	ISM 834	0.003 (0.002) – 35	4	8	1.0 (2.0)–100
MS-4/12	ISM 796	ISM 597	0.002 (0.001) – 24	4	12	1.0 (2.0)–100
Foot switch	ISM 891	ISM 894	see Page 61			

Spare cassettes MS/CA . . . see Page 60

LabVIEW driver for Reglo Digital download for free: www.ismatec.com

Applications

- Addition of a reagent to a reactor and simultaneous removal of the reaction product from the upper fraction. Ramp control combined with a thermostat to maintain the ΔT during the reaction.
- Simultaneous addition of both components of a 2-component adhesive in ratio 1:10 with two different tubing sizes.

1

REGLO *Quick*™

Very fast tubing change-over

Unique, very user-friendly pump-head with fast tubing change-over



REGLO *Quick*
2.1–230 ml/min

CE



Easily accessible tube-bed thanks to wide opening angle.

- Variable speed
- Standard tubing
- 10 cm wide, 14.3 cm high

REGLO *Quick*

- 1 channel / 4 rollers
- Easily accessible tube-bed, with wide opening angle
- Transparent protection cover
- Standard tubing
1.6 mm wall thickness
- Differential pressure 1.5 bar
Depends on tubing material.
Tubing with small i.d.'s may enable higher pressures.
- Higher max. flow rate than
REGLO *Analog/Digital*
- "Start/stop", "speed" and
"direction" functions

Specifications

Motor type	DC motor
Speed	3.2 – 160 rpm
Speed setting	1–99 %, resolution 1% 2-digit potentiometer
Power consumption	30 W
Mains connection	230V _{AC} /50Hz, 115V _{AC} /60Hz adjustable
Protection rating	IP 30

Dimensions

Model	Depth/Width/Height
REGLO <i>Quick</i>	178 x 100 x 143 mm (pump-head closed)
Weight	
REGLO <i>Quick</i>	2.2 kg



2-digit potentiometer
1–99%, resolution 1%(for speed setting)






Interfaces



- Speed control (0–5 or 0–10 V, 0–20 or 4–20 mA)
- Speed output (0–8 kHz)
- Start/Stop
- Rotation direction

Insert the tubing appropriate to your application.



Flow rates and tubing

Model/ Type Channels Rollers			REGLO Quick 1 4	
Speed rpm			3.2	160
Tygon®ST R-3603/R-3607 Order No.	Wall (mm)	Tubing i.d. (mm)	ml/min	
			min.	max.
MF0030	1.6	3.2	2.1	103
SC0379	1.6	4.8	4.6	230
Approx. values: determined with water, at 22°C, no differential pressure, Tygon® tubing.				

Recommended standard tubing

- Tygon® LFL
- Tygon® ST
- Tygon® HC
- Pharmed®
- Silicone peroxide

Do not use other tubing formulations.

Applications

Single-channel delivery processes with variable flow rates where frequent tubing change-over is required e.g.:

- Addition of dye stuffs with tubing exchange after each dispensing process
- Flushing cylinder heads of HPLC pumps

Ordering information

Model	Order No.	Flow rates ml/min per channel	Channels max.	Rollers	Speed rpm
REGLO Quick	ISM 897	2.1 – 230	1	4	3.2 – 160
Foot switch	ISM 891	see Page 61			

1

Ecoline

Economical and powerful

Stackable pumps for dispensing and filling applications with variable flow rates

- Robust stainless steel housing (stackable)
- 2-digit potentiometer adjustable in 1% steps



Ecoline VC-MS/CA8-6
 0.005–150 ml/min
 8 channels
 6 rollers

Ecoline VC-MS/CA4-12
 0.003–83 ml/min
 4 channels
 12 rollers (low pulsation)

- Click'n'go cassettes with automatic pressure setting
- 3-stop tubing
- Differential pressure 1.0 bar ¹

¹ Possible with appropriate tubing material; tubing with small i.d.'s and/or cassettes with the pressure lever (see Page 60) may enable higher pressures.



Ecoline VC-Easy-Load®
 0.23–1600 ml/min
 1 channel
 3 rollers

- Easily accessible pump-head
- Allows rapid tube change-over
- Pump-head PSF housing (Polysulfone)
- Standard tubing 1.6 mm WT
- Differential pressure 1.5 bar ¹

¹ Differential pressure depends on tubing material; tubing with small i.d.'s may enable higher pressures.



Ecoline VC-360
 0.25–1300 ml/min
 1 channel
 3 convex rollers treat the liquid and tubing gently

- Hinged tube-bed for easy and rapid tube change-over
- Standard tubing 1.6 mm WT
- Differential pressure 1.5 bar ¹

¹ Differential pressure depends on tubing material; tubing with small i.d.'s may enable higher pressures.

Specifications

Motor type	DC motor
Speed	3.5 to 350 rpm
Speed setting	1–99%, resolution 1% 2-digit potentiometer
Power consumption	100 W
Mains connection	230V _{AC} /50Hz, 115V _{AC} /60Hz adjustable
Protection rating	IP 30

Size

Model	Depth/Width/Height
Ecoline VC-280	256 x 169 x 138 mm
Ecoline VC-380	256 x 169 x 138 mm
Ecoline VC-360	238 x 169 x 138 mm
Ecoline VC-Easy-Load®	285 x 169 x 138 mm
Ecoline VC-MS/CA8-6	313 x 169 x 138 mm
Ecoline VC-MS/CA4-12	281 x 169 x 138 mm

Weight

Model	Weight
Ecoline VC-280	5.2 kg
Ecoline VC-380	5.3 kg
Ecoline VC-360	4.9 kg
Ecoline VC-Easy-Load®	5.2 kg
Ecoline VC-MS/CA8-6	5.5 kg
Ecoline VC-MS/CA4-12	5.4 kg



1



2



3

- 1 Available without pump-head: ISM 1077
- 2 Easy stackable
- 3 These pump-heads are also available as OEM versions. Ask for the detailed data sheet.

**Easily exchangeable rotor**

WT 1.6 mm, 2 or 3 rollers
VC-280 or VC-380

WT 2.4 mm, 2 or 3 rollers
VC-281 or VC-381

Ecoline VC-280

1.7–5400 ml/min

Ecoline VC-380

1.6–5000 ml/min

1 channel

2 or 3 convex rollers treat the liquid and tubing gently

- With exchangeable rotor e.g. for lower pulsation, higher flow rates, or elevated differential pressures
- Standard tubing 1.6 mm wall thickness (WT)
- Differential pressure 1.5 bar ¹

¹ Differential pressure depends on tubing material; tubing with small i.d.'s may enable higher pressures.

Flow rates / Tubing

Model Type Channels Rollers			Ecoline VC-280 1 2		Ecoline VC-380 1 3		Ecoline VC-360 1 3		Ecoline VC-Easy-Load® 1 3	
Speed rpm			3.5	350	3.5	350	3.5	350	3.5	350
Tygon®ST R-3603/R-3607 Order No.	WT (mm)	Tubing i.d.(mm)	ml/min per channel		ml/min per channel		ml/min per channel		ml/min per channel	
			min.	max.	min.	max.	min.	max.	min.	max.
MF0001	1.6	0.8					0.25	25	0.23	23
MF0028	1.6	1.6	1.7	170	1.6	160	0.9	90	0.86	86
MF0030	1.6	3.2	6.6	660	5.9	590	3.5	350	3.2	320
SC0379	1.6	4.8	51	1500	13	1300	7.7	770	6.5	650
MF0031	1.6	6.4	25	2500	23	2300	13	1300	11	1060
MF0032	1.6	8.0	37	3700	34	3400			16	1600
SC0383	1.6	9.5	48	4800	44	4400				
SC0384	1.6	11.1	54	5400	50	5000				

Approx. values: determined with water at 22°C, no differential pressure, Tygon®ST tubing.

Flow rates / Tubing

Model Type Channels Rollers	Speed rpm	Ecoline VC-MS/CA8-6 8 6		Ecoline VC-MS/CA4-12 4 12	
		3.5	350	3.5	350
Tygon ST R-3603/R-3607 Order No.	Tubing i.d.(mm)	ml/min per channel		ml/min per channel	
		min.	max.	min.	max.
SC0189	0.13	0.005	0.49	0.003	0.32
SC0050	0.25	0.017	1.7	0.013	1.3
SC0053	0.51	0.067	6.7	0.055	5.5
SC0056	0.76	0.15	15	0.12	12
SC0059	1.02	0.26	26	0.20	20
SC0062	1.22	0.36	36	0.26	26
SC0065	1.52	0.53	53	0.36	36
SC0068	1.85	0.73	73	0.47	47
SC0071	2.54	1.2	120	0.68	68
SC0224	3.17	1.5	150	0.83	83

Approx. values: determined with water at 22°C,

no differential pressure, Tygon®ST tubing.

Ordering information

Model	Order No.	Flow rates ml/min per channel	Channels max.	Rollers
Ecoline VC-280 WT 1.6	ISM 1078	1.7 to 5400	1	2
Ecoline VC-380 WT 1.6	ISM 1079	1.6 to 5000	1	3
Ecoline VC-281 WT 2.4*	ISM 1085	13 to 4600	1	2
Ecoline VC-381 WT 2.4*	ISM 1080	12 to 4200	1	3
Rotor with 2 rollers WT 1.6	IS 3762	1.7 to 5400	1	2
Rotor with 3 rollers WT 1.6	IS 3763	1.6 to 5000	1	3
Rotor with 2 rollers, WT 2.4	IS 3757	13 to 4600	1	2
Rotor with 3 rollers, WT 2.4	IS 3754	12 to 4200	1	3
Ecoline VC-360	ISM 1076	0.25 to 1300	1	3
Ecoline VC-Easy-Load®	ISM 1091	0.23 to 1600	1	3
Ecoline VC-MS/CA8-6	ISM 1089	0.005 to 150	8	6
Ecoline VC-MS/CA4-12	ISM 1090	0.003 to 83	4	12
Ecoline without pump-head	ISM 1077			
Foot switch	IS 3572	see Page 61		
Spare cassettes MS/CA	IS 3510	see Page 60		

*For standard tubing 2.4 mm wall thickness 4.8 – 9.5 mm (3/16 – 3/8") inner diameter

Interfaces

- Speed control (0–5 or 0–10 V, 0–20 or 4–20 mA)
- Start/Stop, rotation direction

1

1



Easy interchangeable pump-heads

Mount the pump-head without using a tool

The MCP and BVP drives enable the user to choose individually from a large variety of different pump-heads. These heads are interchangeable and can be mounted or exchanged within seconds.



BVP Standard

see page 16

Variable speed drive with closed loop control analog interface

BVP Standard, Order No. ISM 444



MCP Standard

see page 17

Microprocessor controlled dispensing drive
4 program memories for saving 4 different sets of operating parameters, protection rating IP 30, RS232 and analog interface
MCP Standard, Order No. ISM 404

Interchangeable pump-head with bayonet coupling. Just insert and rotate.



Single-channel

Single-channel

Single-channel



0.072–530 ml/min
Type 360



0.49–3700 ml/min
Type Pro-280
For 1.6 mm WT*
3.6–3100 ml/min
Type Pro-281
For 2.4 mm WT*



0.07–1100 ml/min
Type MF Easy-Load®



0.44–2800 ml/min
Type 380



0.45–3400 ml/min
Type Pro-380
For 1.6 mm WT*
3.3–2900 ml/min
Type Pro-381
For 2.4 mm WT*



0.24–1000 ml/min
Type MF Easy-Load II
(with adjustable pressure setting)



0.41–3600 ml/min
Type 380 AD

4 drives and more than 14 pump-heads
An unique pump system

An investment for the future

For many years ISMATEC® pumps have been well known for being very robust, accurate and reliable, and consequently for offering a very long service life. These qualities make this unique type of pump system with 4 drives and more than 20 interchangeable pump-heads especially interesting. Without high investments you can adapt the chosen pump system at a later date for new applications or different media.



BVP Process

IP 65
dust-tight and protected against water jets

see page 18

Microprocessor controlled drive with digital LED display
Protection rating **IP 65**, analog interface
BVP Process, Order No. ISM 920


MCP Process

IP 65
dust-tight and protected against water jets

Stand-alone application

see page 19

Microprocessor controlled dispensing drive
4 program memories for PC programmed application profiles, protection rating **IP 65** RS232 and analog interface
MCP Process, Order No. ISM 915



1

BVP Standard

Economical

- Robust, powerful drive
- Variable speed



BVP Standard drive
(pump-heads on Pages 20 to 25)

Interfaces



- Speed control
(0–5 or 0–10V, 0–20 or 4–20mA)
- Speed output
(0–10 V_{DC} or 0–12 kHz)
- Start/Stop, rotation direction

BVP Standard

without dispensing functions

- **3-digit potentiometer for speed setting**
- More than 14 pump-heads available
- Bayonet coupling system enables a system change without tools
Flow rates, channels, rollers and differential pressure depend on the mounted pump-head (see Pages 20 to 25)

Specifications

Motor type	DC motor
Speed	2.4 – 240 rpm
Speed setting	1–99.9%, resolution 0.1% 3-digit potentiometer
Power consumption	100 W
Mains connection	230V _{AC} /50Hz, 115V _{AC} /60Hz adjustable
Protection rating	IP 30
Depth/Width/Height	220 x 155 x 260 mm (without pump-head)
Weight	5.7 kg (without pump-head)

Ordering information

The complete pump system BVP Standard consists of:

Drive	ISM 444
Pump-head	Pages 20 to 25
Tubing	Pages 30 to 39
Accessories	Page 61
Foot switch	ISM 891

Applications

With pump-head CA12

Addition of tablet coating solution during pharmaceutical production (pump controlled by external flow-through measuring system)

With pump-head 380

Cross-flow filtration of dyes and suspensions over a number of days for subsequent use of the filtrate in the dyeing industry

MCP *Standard*

Multi-purpose

Saves individual application parameters!

- Robust, powerful drive
- Ideal for dispensing and filling



MCP *Standard* drive
(pump-heads on Pages 20–25)

Interfaces



PC-controllable:
– RS232



– Speed control
(0–5 or 0–10V, 0–20 or 4–20mA)
– Speed output
(0–10 V_{DC} or 0–7.2 kHz)
– Start/Stop, Rotation direction,
Autostart

MCP *Standard*

- Pre-programmed tube sizes and pump-heads allow you to work with flow rates
- Membrane key-pad, LED display
- **4 program memories for saving individual application parameters**
- More than 20 pump-heads available
- Bayonet coupling system enables a system change without tools

Flow rates, channels, rollers and differential pressure depend on the pump-head mounted (see Pages 20 to 25)



Dispensing and calibrating
function see Page 5

CE

Specifications

Motor type	DC motor
Speed	1 – 240 rpm
Speed setting	rpm, resolution 0.1 rpm
Flow rate setting	µl/min, ml/min, l/min
Power consumption	100 W
Mains connections	230V _{AC} /50Hz, 115V _{AC} /60Hz adjustable
Protection rating	IP 30
Depth/Width/Height	220 x 155 x 260 mm (without pump-head)
Weight	6.4 kg (without pump-head)

Ordering information

The complete pump system MCP *Standard* consists of:

Drive	ISM 404
Pump-head	Pages 20 to 25
Tubing	Pages 30 to 39
Accessories	Page 61
Foot switch	IS 10039

Applications

With pump-head 360

Filling of pharmaceutical solutions at small production scale (500 to 1000 units) in syrup bottles. Controlled via foot switch.

With PTFE-Tubing pump-head

Used as “trickle funnel”, time-controlled addition of organic solvents during synthesis (small scale pilot plant)

1

BVP Process

Washdown

Protection rating of IP 65

- Extremely robust drive
- Microprocessor controlled
- Ideal for applications in a dusty, humid or corrosive environment and in clean room areas (IP 65, dust-tight and protected against water jets)



CE

BVP Process drive
(pump-heads on Pages 20 to 25)

Specifications

Motor type	DC motor
Speed	1 – 240 rpm
Speed setting	rpm, resolution 0.1 rpm
Power consumption	120 W
Mains connection	85–264 V _{AC} / 47–60 Hz
Protection rating	IP 65
Depth/Width/Height	220 x 155 x 260 mm (without pump-head)
Weight	6.9 kg (without pump-head)

Ordering information

The complete pump system BVP Process consists of:

Drive	ISM 920
Pump-head	Pages 20 to 25
Tubing	Pages 30 to 39
Accessories	Page 61
Foot switch	IS10039

Interfaces



- Speed control
(0–5 or 0–10V, 0–20 or 4–20mA)
- Speed output
(0–10 V_{DC} or 0–7.2 kHz)
- Start/Stop, rotation direction, autostart

BVP Process

without dispensing functions

Flow rates, channels, rollers and differential pressure depend on the mounted pump-head (see Pages 20 to 25)

- **Membrane key-pad for speed setting, LED display**
- **Stainless steel housing**
- More than 20 pump-heads available
- Bayonet coupling system enables a system change without tools

Applications

With Easyload® II pump-head
Continuous addition of ink for roller
based printing machine

MCP Process

Programmable

Programs can be carried out on the spot independently of a PC!

Protection rating of IP 65

- Extremely robust drive, suitable for industries
- Ideal for dispensing and filling applications in a dusty, humid or corrosive environment, and in clean room areas (IP 65, dust-tight and protected against water jets)

Interfaces



PC controllable:
– RS232



- Speed control (0–5 or 0–10V, 0–20 or 4–20mA)
- Speed output (0–10 V_{DC} or 0–7.2 kHz)
- Start/Stop, Rotation direction, Autostart
- 2 universal inputs
- 2 universal outputs



MCP Process drive
(pump-heads on Pages 20 to 25)

MCP Process

- Pre-programmed tube sizes and pump-heads allow you to work with flow rates
 - Stainless steel housing, membrane key-pad, LED display
 - **4 program memories for saving individual application parameters or PC programmed command sequences**
 - More than 20 pump-heads available
 - Bayonet coupling system enables a system change without tools
- Flow rates, channels, rollers and differential pressure depend on the pump-head mounted (see Pages 20 to 25)



Get your free plugin for ProgEdit
(German/English switchable)
see Page 61



Dispensing and calibrating
function see Page 5

Specifications

Motor type	DC motor
Speed	1 – 240 rpm
Speed setting	rpm, resolution 0.1 rpm
Flow rate setting	µl/min, ml/min, l/min
Power consumption	100 W
Mains connection	100–230 V _{AC} / 50 – 60 Hz
Protection rating	IP 65
Depth/Width/Height	220 x 155 x 260 mm (without pump-head)
Weight	6.9 kg (without pump-head)

Ordering information

The complete pump system MCP Process consists of:	
Drive	ISM 915
Pump-head	Pages 20 to 25
Tubing	Pages 30 to 39
Accessories	Page 61
Software ProgEdit	free
Foot switch	IS 10039

LabVIEW driver
download for free: www.ismatec.com

Applications

With pump-head Pro 381 and pressure sensor

Pressure controlled sterile filtration of human plasma through a filter plate 0.2 µm. Constant monitoring is no longer necessary thanks to an optical alarm if pre-set pressure is exceeded.

With pump-head MS/CA 4-12 and switching valves

Stabilizing of an electrophysiological test system with Ringer solution and switching to the test solutions without delay at very low flow rates.

1

Gentle pumping

This pump-head is gentle enough for pumping highly concentrated viable cells



Pro-280 0.49–3700 ml/min

Order No. ISM 785

- Coated aluminum pump-head
- Can be dismantled for cleaning
- Self-centering tube-track thanks to concave tube-bed and convex rollers, which lengthens the tube-life
- 2 stainless steel rollers (higher max. flow rate but more pulsation than with 3 rollers)
- For tubing with 1.6 mm wall thickness
- 1.5 bar differential pressure¹

Pro-281 3.6–3100 ml/min

Order No. ISM 793

Same pump-head as Pro-280, but

- For tubing with **2.4 mm wall thickness**
- 2.5 bar differential pressure¹

Flow rates/ Tubing		Model Pro-280			
Tygon®ST R-3603/R-3607 Order No.	Wall	Tubing	ml/min		
	(mm)	i.d. (mm)	min.	max.	
MF0028	1.6	1.6	0.49	120	
MF0030	1.6	3.2	1.9	450	
SC0379	1.6	4.8	4.2	1000	
MF0031	1.6	6.4	7.2	1700	
MF0032	1.6	8.0	11	2600	
SC0383	1.6	9.5	14	3300	
SC0384	1.6	11.1	16	3700	

Flow rates/ Tubing		Model Pro-281			
Tygon ST R-3603/R-3607 Order No.	Wall	Tubing	ml/min		
	(mm)	i.d. (mm)	min.	max.	
MF0029	2.4	4.8	3.6	870	
MF0033	2.4	6.4	6.5	1600	
SC0502	2.4	8.0	9.9	2400	
SC0503	2.4	9.5	13	3100	

Instantly interchangeable pump system



Pro-380 0.45–3400 ml/min

Order No. ISM 791

Same pump-head as Pro-280, but

- 3 stainless steel rollers (less pulsation but lower max. flow rate than with 2 rollers)

Pro-381 3.3–2900 ml/min

Order No. ISM 797

Same pump-head design as Pro-280, but

- 3 stainless steel rollers (less pulsation but lower max. flow rate than with 2 rollers)
- For tubing with **2.4 mm wall thickness**
- 2.5 bar differential pressure¹

Flow rates/ Tubing		Model Pro-380			
Tygon ST R-3603/R-3607 Order No.	Wall	Tubing	ml/min		
	(mm)	i.d. (mm)	min.	max.	
MF0028	1.6	1.6	0.45	110	
MF0030	1.6	3.2	1.7	400	
SC0379	1.6	4.8	3.7	890	
MF0031	1.6	6.4	6.5	1600	
MF0032	1.6	8.0	9.7	2300	
SC0383	1.6	9.5	13	3000	
SC0384	1.6	11.1	14	3400	

Flow rates/ Tubing		Model Pro-381			
Tygon ST R-3603/R-3607 Order No.	Wall	Tubing	ml/min		
	(mm)	i.d. (mm)	min	max	
MF0029	2.4	4.8	3.3	800	
MF0033	2.4	6.4	5.8	1400	
SC0502	2.4	8.0	8.8	2100	
SC0503	2.4	9.5	12	2900	

Especially suitable for:

- Chemical, biotechnological and pharmaceutical applications
- Food industry
- Elevated differential pressures (Pro-281 and Pro-381)
- Viscous fluids
- Fluids containing a high content of sensitive solids
- Applications requiring hygienic conditions, durability and reliability

Comparisons to gear, piston and centrifugal pumps proved that peristaltic pumps are the only suitable and sterilisable pump system for gently pumping media containing living cells.

¹ Differential pressure depends on tubing material; tubing with small i.d.'s enable higher pressures.

The flow rates are based on a drive speed of 1 (or 2.4) to 240 rpm. For the BVP Standard drive the indicated min. flow rates must be multiplied by factor 2.4.

Approx. values: determined with water, at 22° C, no differential pressure, Tygon tubing.



Application

Pumping animal cell cultures

A comparison between different peristaltic pump systems was made by Ms Fan Guo and Prof. Dr. U. Graf-Hausner of Zurich Technical University Winterthur, and Dr. Joanne Laukart of Ismatec.


Experimental target: Gentle cell delivery

- Cell viability remains high
- Cell damage remains low


Pumping application

Production monitoring by analyzing by-pass cell streams and reactor inoculation.

Test results



With the pump-head Pro-280 cell death was observed only after 56 hours whilst with the Easy-Load[®] pump-head death occurred already after 29 hours.




With the pump-head Pro-280 a viability of 87% was observed after **56 hours**. With the Easy-Load head a viability of only 72% was observed already after 29 hours.

Method

Cultivation of hybridoma cells in spinner bottles.

The cells are continuously pumped through an external loop at a constant flow rate ($\frac{1}{2}$ reactor volume / hour = 250 ml/h) until a identifiable cell death is observed.



The secret of the Pro-280

Before the roller totally closes the squeezed tubing, the cells can escape through a gap towards the tubing wall and, hence, are neither «squashed» nor damaged.

The pump-head Pro-280 has 2 convex rollers and the tube-bed is concave.

Conclusion

The Pro-280 (mounted either on an MCP or BVP Process drive) is a uniquely designed pump system especially suited for use in biotechnology. It has been tested for use in laboratory scale and small production scale up to 30 liters.

Pumps tested

(Colors according to diagram)

MCP Process Pro-280 2 convex, non-spring-loaded rollers	MCP Process WM5 2 straight, spring-loaded rollers
MCP Process 380AD 3 convex, non-spring-loaded rollers	MCP Process Easy Load 3 straight, non-spring-loaded rollers

Comparative parameters

1. Number of pump rollers: 2 or 3
2. Rollers: spring-loaded vs. non-spring-loaded
3. Tube-bed / roller geometry: convex vs. straight rollers



The winner: Pro-280


Thanks to convex rollers and a concave tube-bed, the MCP Process Pro-280 proved to be the gentlest pump system with respect to cell viability and cell concentration.

1

Multi-purpose single-channel pump-heads

For MCP and BVP drives

Change the pump-heads within seconds without the use of tools



360¹ 0.072–530 ml/min

Order No. ISM 719

- Easily accessible flip-up tube-bed guarantees easy and rapid tube change-over
- Transparent protection cover allows monitoring the tube and the revolving rotor
- Self-centering tube-track design thanks to the concave tube-bed and convex rollers (lengthens tube-life)
- Rotor for tubing i.d. from 0.8 to 6.4 mm with 1.6 mm wall thickness
- 3 stainless steel rollers
- 1.5 bar differential pressure²

Flow rates / Tubing Model 360

Tygon®ST R-3603/R-3607 Order No.	Wall	Tubing	ml/min	
	(mm)	i.d. (mm)	min.	max.
MF0001	1.6	0.8	0.072	17
MF0028	1.6	1.6	0.26	62
MF0030	1.6	3.2	1.0	240
SC0379	1.6	4.8	2.2	530



¹ An OEM version of this pump-head is also available. Ask for the detailed data sheet.



380¹ 0.44–2800 ml/min

Order No. ISM 718

Same design as pump-head 360, but larger size

- For tubing i.d. from 1.6 to 9.5 mm with 1.6 mm wall thickness
- 1.5 bar differential pressure²

Ideal for sterile media

Flow rates / Tubing Model 380

Tygon ST R-3603/R-3607 Order No.	Wall	Tubing	ml/min	
	(mm)	i.d. (mm)	min.	max.
MF0028	1.6	1.6	0.44	100
MF0030	1.6	3.2	1.7	400
SC0379	1.6	4.8	3.6	860
MF0031	1.6	6.4	6.0	1400
MF0032	1.6	8.0	8.8	2100
SC0383	1.6	9.5	12	2800

² Differential pressure depends on tubing material; tubing with small i.d.'s may enable higher pressures.



380AD 0.41–3600 ml/min

Order No. ISM 725

- Pressure on tubing adjustable via rollers
- Rotor accepts tubing with 1.6 and 2.4 mm wall thickness and 1.6 to 11.1 mm i.d.
- 3 stainless steel rollers
- 1.5 bar differential pressure² (with tubing wall thickness 1.6 mm)
- 2.5 bar differential pressure² (with tubing wall thickness 2.4 mm)

Ideal for media with high viscosity or a certain level of solid content

Flow rates / Tubing Model 380AD

Tygon ST R-3603/R-3607 Order No.	WS	Tubing	ml/min	
	(mm)	iØ (mm)	min.	max.
MF0028	1.6	1.6	0.4	99
MF0030	1.6	3.2	1.5	370
SC0379	1.6	4.8	3.4	830
MF0031	1.6	6.4	6.2	1500
MF0032	1.6	8.0	9.5	2300
SC0383	1.6	9.5	13.0	3000
SC0384	1.6	11.1	15.0	3600
MF0029	2.4	4.8	3.4	830
MF0033	2.4	6.4	6.2	1500

▶ The flow rates are based on a drive speed of 1 (or 2.4) to 240 rpm. For the BVP Standard drive the indicated min. flow rates must be multiplied by factor 2.4. Approx. values: determined with water, at 22° C, no differential pressure, and Tygon tubing.

The pump-head for PTFE tubing



Easy-Load[®] 3 0.07–1100 ml/min

Order No. MF 0313/ISM738

- Easily accessible pump-head
- Allows rapid tube change-over
- PSF housing (Polysulfone)
- Rotor for tubing with 1.6 mm wall thickness
- Rotor with 3 stainless steel rollers
- 0.7 bar differential pressure²

Flow rates / Tubing Model Easy-Load

Tygon [®] ST R-3603/R-3607 Order No.	Wall	Tubing	ml/min	
	(mm)	i.d. (mm)	min.	max.
MF0001	1.6	0.8	0.066	16
MF0028	1.6	1.6	0.25	59
MF0030	1.6	3.2	0.91	220
SC0379	1.6	4.8	1.9	450
MF0031	1.6	6.4	3.1	730
MF0032	1.6	8.0	4.7	1100

² Differential pressure depends on tubing material; tubing with small i.d.'s may enable higher pressures.



Easy-Load II³ 0.24–1000 ml/min

Order No. MF 0446/ISM738

- Same specifications as Easy-Load, but
- Adjustable pressure setting
 - Improved, automatic tubing retention
 - PPS housing (Polyphenylene sulfide)
 - Rotor with 4 stainless steel rollers
 - 0.7 bar differential pressure²

Flow rates / Tubing Model Easy-Load II

Tygon ST R-3603/R-3607 Order No.	Wall	Tubing	ml/min	
	(mm)	i.d. (mm)	min.	max.
MF0028	1.6	1.6	0.24	58
MF0030	1.6	3.2	0.92	220
SC0379	1.6	4.8	1.9	460
MF0031	1.6	6.4	3.0	730
MF0032	1.6	8.0	4.2	1000

³ 2 pump-heads can be mounted on one drive. (Special mounting sets must be ordered separately).



PTFE-Tube 0.07–45 ml/min

Order No. MF 0330/ISM727

Pump-head for PTFE tubing

- 6 stainless steel rollers
- Stainless steel rotor
- Anodized aluminum body
- Adjustable tube-bed pressure with locking knob
- Up to 6.9 bar differential pressure

Ideal for dispensing and pumping aggressive chemicals and for the filtration of organic solvents.

Flow rates for PTFE-Tubing

Tubing i.d.	Flow rates ml/min
2 mm	0.07 to 15
4 mm	0.19 to 45



See ordering information and tubing below.

Ordering information for pump-head PTFE tubing

Article	Specifications	Order No.	Pack size
PTFE pump tubing, 38 cm long	2 mm i.d. / 4 mm o.d., for 0.07 – 15 ml/min	MF 0331	2 pieces
PTFE pump tubing, 38 cm long	4 mm i.d. / 6 mm o.d., for 0.19 – 45 ml/min	MF 0332	2 pieces
Tube connectors (straight) (2 connectors are needed for one tube)	for tubing with 2 mm i.d.	MF 0333	3 pieces
Tube connectors (straight) (2 connectors are needed for one tube)	for tubing with 4 mm i.d.	MF 0334	1 pieces
PTFE extension tubing, 3.65 m long	for tubing with 2 mm i.d.	SC 1017BO	1 x 3.65 m
PTFE extension tubing, 3.65 m long	for tubing with 4 mm i.d.	SC 1016BO	1 x 3.65 m
Tubing grooving tool (Use only tubing with 4 mm i.d.)	Important for connections which must withstand 2.8 bar (40 psi) or greater.	MF 0337	1 pieces

1

Pump-heads up to 24 channels

For MCP and BVP drives

Change the pump-heads in seconds without tools



CA-4, CA-8, CA-12 **0.002–230 ml/min**

Order No. ISM 721 CA 4 (with cassettes for 4 channels)
 Order No. ISM 732 CA 8 (with cassettes for 8 channels)
 Order No. ISM 733 CA 12 (with cassettes for 12 channels)

- CA Click'n'go cassettes (spare cassettes, see Page 60)
- Automatic pressure setting
- Easy and rapid tube change-over; each channel separately, even while pump is running
- 8 rollers
- **4 to 12 channels**, each channel can take different tube sizes
- 2-stop tubing, 0.13 to 3.17 mm i.d.
- 1.0 bar differential pressure¹

Ideal pump-head for multi-channel applications at higher flow rates

Flow rates / Tubing				
Tygon®ST R-3603/R-3607 Order No.	Tubing i.d. (mm)	ml/min per channel		
		min.	max.	
SC0188	0.13	0.002	0.31	
SC0002	0.25	0.008	1.8	
SC0005	0.51	0.034	8.2	
SC0008	0.76	0.074	18	
SC0011	1.02	0.13	31	
SC0014	1.22	0.18	42	
SC0017	1.52	0.26	62	
SC0020	1.85	0.36	86	
SC0023	2.54	0.62	150	
SC0222	3.17	0.94	230	



OEM versions of these pump-heads are also available. Ask for the detailed data sheet.



SB2V
ISM 734 + ISM 010
Standard tubing
3.2 to 8.0 mm i.d.
Wall thickness 1.6 mm



SB3V
ISM 734 + ISM 011
Standard tubing
0.8 to 4.8 mm i.d.
Wall thickness 1.6 mm

SB 2V / SB 3V **0.09–1100 ml/min**

Order No. ISM 734 (depends on tube-bed)

+ Tube-bed set:
 2-channel, SB 2V ISM 010 for Standard tubing 1.1 to 1100 ml/min
 3-channel, SB 3V ISM 011 for Standard tubing 0.09 to 530 ml/min

- Spring-loaded, exchangeable tube-beds
- Individual and continuously adjustable pressure setting per channel
- 6 rollers, **2 to 3 channels**, for Standard tubing
- Differential pressure 1.5 bar¹

Ideal for sensitive substances requiring a smooth pressure adjustment

Flow rates / Tubing SB2V				
Tygon ST R-3603/R-3607 Order No.	WT (mm)	Tubing i.d. (mm)	ml/min per channel	
			min.	max.
MF0030	1.6	3.2	1.1	260
SC0379	1.6	4.8	2.3	550
MF0031	1.6	6.4	3.7	890
MF0032	1.6	8.0	4.6	1100

Flow rates / Tubing SB3V				
Tygon ST R-3603/R-3607 Order No.	WT (mm)	Tubing i.d. (mm)	ml/min per channel	
			min.	max.
MF0001	1.6	0.8	0.09	22
MF0028	1.6	1.6	0.26	63
MF0030	1.6	3.2	0.99	240
SC0379	1.6	4.8	2.2	530



Instantly interchangeable pump system

Choose from 4 drives
Memory for up to 4 user-defined programs. Pump-heads and tubing are pre-programmed, which enables you to work with varying flow rates. See Pages 14 to 19.

Single-channel pump-heads for MCP and BVP drives:
See Pages 20 to 23.

▶ Tubing for aggressive media, see Page 35.



MS/CA 4-12 0.001–57 ml/min

- Order No. ISM 735
+ max. 3 extension blocks ISM 737
- Click'n'go cassettes (spare cassettes, Page 60)
 - Each channel can take different tube sizes. Tube change-over is possible even while pump is running.
 - **12 rollers for extremely low pulsation**
 - 4 channels, extendable **up to 16 channels**
 - 3 extension blocks with 4 channels each
 - 3-stop tubing, 0.13 to 3.17 mm i.d.
 - 1.0 bar differential pressure¹

Flow rates / Tubing			
Tygon®ST R-3603/R-3607 Order No.	Tubing i.d. (mm)	ml/min per channel	
		min.	max.
SC0189	0.13	0.001	0.22
SC0050	0.25	0.004	0.91
SC0053	0.51	0.016	3.8
SC0056	0.76	0.033	8.0
SC0059	1.02	0.056	13
SC0062	1.22	0.075	18
SC0065	1.52	0.10	25
SC0068	1.85	0.13	32
SC0071	2.54	0.19	46
SC0224	3.17	0.24	57

MS/CA 8-6 0.002–100 ml/min

- Order No. ISM 724
+ max. 2 extension blocks ISM 185
- Click'n'go cassettes (spare cassettes, Page 60)
 - Each channel can take different tube sizes. Tube change-over is possible even while pump is running.
 - 6 rollers
 - 8 channels, extendable **up to 24 channels**
 - 2 extension blocks with 8 channels each
 - 3-stop tubing, 0.13 to 3.17 mm i.d.
 - 1.0 bar differential pressure¹

Flow rates / Tubing			
Tygon ST R-3603/R-3607 Order No.	Tubing i.d. (mm)	ml/min per channel	
		min.	max.
SC0189	0.13	0.002	0.33
SC0050	0.25	0.005	1.1
SC0053	0.51	0.019	4.6
SC0056	0.76	0.042	10
SC0059	1.02	0.073	18
SC0062	1.22	0.10	24
SC0065	1.52	0.15	36
SC0068	1.85	0.21	50
SC0071	2.54	0.33	79
SC0224	3.17	0.43	100

¹ Differential pressure depends on tubing material; tubing with small i.d.'s and/or cassettes with the pressure lever (see Page 60) may enable higher pressures.

▶ The flow rates are based on a drive speed of 1 (or 2.4) to 240 rpm. For the BVP Standard drive the indicated min. flow rates must be multiplied by factor 2.4. The indicated flow rates are approx. values: determined with water, at 22° C, no differential pressure, and Tygon tubing.

1

MS-CA Stand-mounted pumps

Save space with stand-mounted pumps

- **Fixed speed**
- Excellent start/stop characteristics
- Highly reproducible
- Ideal for limited space situations e.g. cabinets



MS-CA

0.021–26 ml/min (per channel)

- 2 or 4 channels
- 6 or 8 rollers
- Click'n'go cassettes
- Each channel can take different tube sizes
- 3-stop tubing
- Differential pressure 1.0 bar
Possible with appropriate tubing material;
tubing with small i.d.'s and/or pressure lever
cassettes (see Page 60) may enable higher pressures.



Spare cassettes MS/CA
(see Page 60)

Specifications MS-CA

Motor type	Synchronous motor
Speed	20, 40 or 60 rpm
Power consumption	8 W
Mains connection ¹	230V _{ac} (50/60 Hz) or 115V _{ac} (50/60 Hz)
Protection rating	IP 30
Depth/Width/Height	2-channel 125x88x135 mm 4-channel 145x88x135 mm
Weight	2-channel 1.2 kg 4-channel 1.3 kg

¹ Please state required mains voltage and frequency on your order.

These pumps are also available as OEM versions.

Our engineers welcome your ideas.
Please contact us.



Contents of tubing selection guide:

(is supplied with each pump)

- Material and tubing properties
- Chemical resistance chart
- Ordering information (see also Pages 30 to 39 and Internet www.ismatec.com)

Flow rates and tubing

Model Channels Rollers	MS-CA _620 2 / 4 6	MS-CA 4/640 4 6	MS-CA 2/660 2 6	MS-CA _820 2 / 4 8	MS-CA 4/840 4 8	MS-CA 2/860 2 8	
	Speed rpm	20	40	60	20	40	60
Tygon® ST R-3603/R-3607 Order No.	Tubing i.d. (mm)	ml/min per channel	ml/min per channel	ml/min per channel	ml/min per channel	ml/min per channel	ml/min per channel
SC0189	0.13	0.028	0.055	0.083	0.021	0.043	0.064
SC0050	0.25	0.10	0.19	0.29	0.08	0.16	0.24
SC0053	0.51	0.38	0.77	1.2	0.34	0.67	1.0
SC0056	0.76	0.84	1.7	2.5	0.73	1.5	2.2
SC0059	1.02	1.5	2.9	4.4	1.3	2.5	3.8
SC0062	1.22	2.0	4.1	6.1	1.8	3.5	5.3
SC0065	1.52	3.0	6.0	9.1	2.6	5.1	7.7
SC0068	1.85	4.2	8.4	13	3.5	7.0	10
SC0071	2.54	6.6	13	20	5.5	11	16
SC0224	3.17	8.5	17	26	7.1	14	21

Approx. values: determined with water, at 22°C, no differential pressure, Tygon tubing.

Interfaces

No interfaces for external control

1



Long-term determination of environmental influences on asbestos-cement

Ordering information

Model	Order No.	Flow rates* ml/min per channel	Channels max.	Rollers	Speed rpm
MS-CA 2/620	ISM 844	0.028 – 8.5	2	6	20
MS-CA 2/640	ISM 845	0.055 – 17.0	2	6	40
MS-CA 2/660	ISM 846	0.083 – 26.0	2	6	60
MS-CA 2/820	ISM 847	0.021 – 7.1	2	8	20
MS-CA 2/840	ISM 848	0.043 – 14.0	2	8	40
MS-CA 2/860	ISM 849	0.064 – 21.0	2	8	60
MS-CA 4/620	ISM 850	0.028 – 8.5	4	6	20
MS-CA 4/640	ISM 851	0.055 – 17.0	4	6	40
MS-CA 4/820	ISM 852	0.021 – 7.1	4	8	20
MS-CA 4/840	ISM 853	0.043 – 14.0	4	8	40

* 26 different tubing i.d.'s give 26 flow rates!

Please state required mains connection 230V_{AC} (50/60 Hz) or 115V_{AC} (50/60 Hz) on your order.

Spare cassettes MS/CA . . . (see Page 60)

Applications

- Multi-channel delivery processes with constant flow rate, e.g.
- Sipper pump for flow-through cuvettes
- Feeding of overflow level control systems

1

Flowmaster®

Ideal for heavy-duty processes

The most convincing solution for fluid transfer in biotechnology, as well as pharmaceutical, food and beverage industries

- Ideal for dispensing and filling applications in a dusty, humid or corrosive environment and in clean room areas
- Convex rollers treat the liquid and tubing gently
- Protection rating of IP 65



CE



Optimized for increased hygienic requirements

- Stainless steel housing
- Tube-loading under sterile conditions without aspirating air
- **Easy disassembly of the pump-head**
Thorough cleaning thanks to easy and rapid disassembly and reassembly in less than a minute.

Flowmaster FMT300

- 37 ml/min – 13 l/min**
- 1 channel
 - 3 convex stainless steel rollers
 - Automatic tube retention
 - Standard tubing
6.4 to 15.9 mm i.d.,
wall thickness 3.2 mm,
differential pressure max. 2 bar
(depends on tubing material used)

Safety

- Pump stops when opening the tube-bed
- Multiple overload protection

Specifications

Motor type	DC motor
Speed	5 to 500 rpm
Speed setting	resolution 0.1 rpm membrane key-pad, LED display
Power consumption	500 W
Mains connection	230V _{AC} /50Hz, 115V _{AC} /60Hz adjustable
Protection rating	IP 65
Depth/Width/Height	500x220x262 mm
Weight	26 kg

flowmaster®

Tube exchange
in **5** seconds



- Insert the tube (easily and fast)
- Press down the lever (automatically correct pressure setting of the tube)
- Start the pump!

Flow rates / Tubing

Tubing Information				Flow rates in liters/min							
Tygon® LFL	PharMed®	Wall (mm)	Tubing i.d. (mm)	rpm 5	rpm 10	rpm 50	rpm 100	rpm 200	rpm 300	rpm 400	rpm 500
Order No.	Order No.										
SC0531	MF0015	3.2	6.4	0.037	0.074	0.37	0.74	1.5	2.2	3.1	3.7
SC0395	MF0016	3.2	9.5	0.08	0.16	0.80	1.6	3.2	4.8	6.4	8.0
SC0396	MF0034	3.2	12.7	0.10	0.20	1.0	2.0	4.0	6.0	8.0	10.0
	SC0696	3.2	15.9	0.13	0.26	1.3	2.6	5.2	7.8	10.4	13.0

Approx. values: determined with water, at 22°C, no differential pressure, PharMed tubing

Ordering information

Model	Order No.	Flow rates ml/min	Channels max.	Rollers	Speed rpm
Flowmaster FMT300 230V _{AC} / 50 Hz	ISM 1020	37 – 13'000	1	3	5–500
Flowmaster FMT300 115V _{AC} / 50 Hz	ISM 1022	37 – 13'000	1	3	5–500

Accessories

Foot switch (Page 61) IS 10279

▶ Tubing for aggressive media, see Page 35.

PLC compatible interface with status information for process control systems (the level of the inputs can be configured: 5, 12 or 24 V)



- Speed control (0–5 or 0–10V, 0–20 or 4–20mA)
- Start/Stop, Rotation direction
- Autostart
- Speed output
- Digital output (potential free) (Error, okay, busy)

Settings menu

- Configuration of analog interface
- Entry of basic settings, e.g.
 - Foot switch control
 - Rotation speed (% or rpm)
 - Service life of tubing
 - **Timer function, etc.**



Suited for:
1-channel delivery processes with large flow rates in clean room areas or a dusty, humid or wet environment.

Applications

Filling of large volumes of pharmaceutical solution from a seed reactor to bulk vessels under sterile conditions.

Rating comparison

Rating: + meets the stated property 1 not recommended
 ± meets the stated property to limited extent 10 excellent
 – does not meet the stated property

Properties	Tygon® LFL	Tygon® ST R-3603/R-3607	PharMed® Ismaprene	Tygon® MHSL 2001	Tygon® MHLL
FDA	+	+	+	+	+
US Pharmacopoeia Class VI	+	–	+	–	+
Transparency	+	+	–	+	–
Long Life	7	1	10	3	10
Gas Permeability					
CO ₂	8	7	5	5	5
O ₂	9	9	8	9	8
N ₂	10	9	8	6	8
Temperature, above 0°C	2	2	7	1	7
Temperature, below 0°C	3	4	8	7	8
Pressure	9	5	1	1	1
Absorption / Adsorption	6	6	9	10	9
Chemical Resistance,					
Acids (H ₂ SO ₄)					
10%	10	10	10	10	10
30%	10	10	10	10	10
95–98%	1	1	1	7	1
Bases (NaOH)					
10–15%	10	10	10	10	10
30–40%	4	4	10	10	10
Hydrocarbons					
(aliphatic)	1	1	1	1	1
Mineral Salts	10	10	10	10	10
Alcohols	1	1	10	10	10
Ketones (Acetone)	1	1	1	7	1

Maximum recommended operating pressure

WT* (mm)	i.d. (mm)	bar	bar	bar	bar	bar
1.6	0.8	8.7	8.7	3.7	N.A.	N.A.
1.6	1.6	4.8	4.8	2.1	3.1	N.A.
1.6	2.4	3.8	3.8	1.6	N.A.	N.A.
1.6	3.2	3.0	3.0	1.3	2.0	N.A.
1.6	4.8	2.2	2.2	0.9	1.5	N.A.
1.6	6.4	1.8	1.8	0.8	1.1	N.A.
1.6	8.0	1.5	1.5	0.6	0.9	N.A.
1.6	9.5	1.3	1.3	0.5	0.8	N.A.
1.6	11.1	1.2	1.2	0.5	N.A.	N.A.
1.6	12.7	1.1	1.1	0.5	N.A.	N.A.
1.6	15.9	1.0	1.0	0.4	N.A.	N.A.
2.4	4.8	3.0	3.0	1.3	N.A.	N.A.
2.4	6.4	2.4	2.4	1.0	N.A.	N.A.
2.4	8.0	2.0	2.0	0.8	N.A.	N.A.
2.4	9.5	1.8	1.8	0.8	N.A.	N.A.
2.4	11.1	1.5	1.5	0.6	N.A.	N.A.
2.4	12.7	1.3	1.3	0.6	N.A.	N.A.
2.4	15.9	1.2	1.2	0.5	N.A.	N.A.
3.2	6.4	3.0	3.0	1.3	N.A.	N.A.
3.2	9.6	2.2	2.2	0.9	N.A.	N.A.
3.2	12.7	1.8	1.8	0.8	1.1	N.A.
3.2	15.9	1.5	1.5	0.6	0.9	N.A.

*WT: Wall thickness

We recommend you to also observe the chart 'Tubing Properties' on the following pages. All information has been supplied to ISMATEC by the tubing manufacturers. It is for your guidance only. We recommend you to test the tubing before use.

Properties	Tygon® HC F-4040-A	Tygon® SI Silicone 3350 (Platinum)	Silicone Peroxid	Norprene® A-60-G	Viton® Fluran® HCA F-5500-A
FDA	-	+	+	-	-
US Pharmacopoeia Class VI	-	+	+	-	-
Transparency	±	±	±	-	-
Long Life	2	4	4	10	3
Gas Permeability					
CO ₂	9	1	1	5	10
O ₂	10	1	1	8	10
N ₂	10	1	1	8	10
Temperature, above 0°C	2	10	10	7	9
Temperature, below 0°C	1	10	10	8	4
Pressure	7	1	1	1	1
Absorption / Adsorption	6	1	1	9	7
Chemical Resistance,					
Acids (H ₂ SO ₄)					
10%	10	10	10	10	10
30%	7	7	8	10	10
95-98%	1	1	1	1	10
Bases (NaOH)					
10-15%	1	10	10	10	10
30-40%	1	10	10	10	10
Hydrocarbons					
(aliphatic)	7	1	1	1	7
Mineral Salts	10	7	7	10	10
Alcohols	7	7	10	10	1
Ketones (Acetone)	1	4	1	1	1

Maximum recommended operating pressure

WT* (mm)	i.d. (mm)	bar	bar	bar	bar	bar
1.6	0.8	10.9	1.9	1.9	3.7	3.7
1.6	1.6	6.1	1.0	1.0	2.1	2.1
1.6	2.4	4.8	0.8	0.8	1.6	1.6
1.6	3.2	3.8	0.6	0.6	1.3	1.3
1.6	4.8	2.7	0.5	0.5	0.9	0.9
1.6	6.4	2.2	0.4	0.4	0.8	0.8
1.6	8.0	1.8	0.3	0.3	0.6	0.6
1.6	9.5	1.6	0.3	0.3	0.5	0.5
1.6	11.1	1.5	0.3	0.3	0.5	0.5
1.6	12.7	1.4	0.2	0.2	0.5	0.5
1.6	15.9	1.2	0.2	0.2	0.4	0.4
2.4	4.8	3.8	0.6	0.6	1.3	1.3
2.4	6.4	3.0	0.5	0.5	1.0	1.0
2.4	8.0	2.5	0.4	0.4	0.8	0.8
2.4	9.5	2.2	0.4	0.4	0.8	0.8
2.4	11.1	1.8	0.3	0.3	0.6	0.6
2.4	12.7	1.7	0.3	0.3	0.6	0.6
2.4	15.9	1.5	0.3	0.3	0.5	0.5
3.2	6.4	3.8	0.6	0.6	1.3	1.3
3.2	9.6	2.7	0.5	0.5	0.9	0.9
3.2	12.7	2.2	0.4	0.4	0.8	0.8
3.2	15.9	1.8	0.3	0.3	0.6	0.6

*WT: Wall thickness

Properties of ISMATEC® tubing

Proven quality for a wide range of laboratory applications. Be sure to choose the tubing most suitable for your application.

Type	Tygon® LFL	Tygon® ST R-3603/R-3607	PharMed® Ismaprene	Tygon® MHSL 2001	Tygon® MHLL																																																																						
Special Properties	The tubing with the longest service-life of any clear Tygon® tubing	The inexpensive all-round tubing for general laboratory applications	The ideal tubing for pharmaceutical and medical applications, and for foodstuffs.	The transparent, plasticiser-free tubing with superior pump-life. Especially designed for MEK and other aggressive solvents.	Chemically resistant to acetone, MEK and other aggressive solvents. Long life tubing.																																																																						
Advantages	<ul style="list-style-type: none"> • Transparent • Broad chemical resistance • Tasteless • Extremely low particulate spallation • Meets USP Class VI and FDA criteria • Non-aging • High dielectric constant 	<ul style="list-style-type: none"> • Transparent • Resistant to almost all inorganic chemicals • Tasteless • Smooth polished inner wall • Low gas permeability • Non-aging and non-oxidizing • High dielectric constant 	<ul style="list-style-type: none"> • Recommendable for cell and tissue cultures • Ideal for production filtration, fermentation and bioreactor process lines • Very long service-life • Non-toxic and non-hemolytic • Impermeable to normal light and UV-radiation • Appropriate for medical products and foodstuffs • Low particulate spallation • Can be autoclaved repeatedly • Withstands repeated CIP and SIP cleaning and sterilization • Meets USP class VI, FDA and NSF criteria 	<ul style="list-style-type: none"> • Plasticizer and oil-free • Smooth inner-bore • Low sorption maintains fluid and tube integrity • Does not impart anything into the pumping medium • No release of hazardous materials when properly incinerated 	<ul style="list-style-type: none"> • Plasticizer-free • Smooth innerbore • Low sorption maintains fluid integrity • Minimal adhesion and diffusion • Suitable for MEK, Acetone and other corrosive solvents • Long life tubing 																																																																						
Limitations	<ul style="list-style-type: none"> • Potential leaching of plasticizers • Not recommended for human blood and tissue 	<ul style="list-style-type: none"> • Potential leaching of plasticizers • Short service-life 	<ul style="list-style-type: none"> • Potential leaching of additives (lubricants) 	None	<ul style="list-style-type: none"> • Cannot be repeatedly sterilized • Only available as stopper tubing 																																																																						
Physical Properties	<ul style="list-style-type: none"> • Thermoplastic • PVC-based material with plasticizer • Flexible, firm, transparent 	<ul style="list-style-type: none"> • Thermoplastic • PVC-based material with plasticizer • Flexible, firm, transparent 	<ul style="list-style-type: none"> • Thermoplastic Elastomer based on polypropylene. • Flexible, firm, opaque 	• Polyolfin	<ul style="list-style-type: none"> • Special thermoplastic of high purity • Without additives • Without plasticizer • Environmental-friendly disposal • Flexible, firm, opaque 																																																																						
Service Temperature Range	-50°C to +74°C (-58°F to +165°F)	-50°C to +74°C (-58°F to +165°F)	-60°C to +135°C (-75°F to +275°F)	-73°C to +57°C (-100°F to +135°F)	-70°C to +74°C (-94°F to +165°F)																																																																						
Applications	<table border="1"> <tr><td>Acids</td><td>good</td></tr> <tr><td>Alkaline solutions</td><td>good</td></tr> <tr><td>Solvents</td><td>not recommended</td></tr> <tr><td>Pressure</td><td>good</td></tr> <tr><td>Vacuum</td><td>good</td></tr> <tr><td>Viscous media</td><td>excellent</td></tr> <tr><td>Sterile Media</td><td>limited</td></tr> </table>	Acids	good	Alkaline solutions	good	Solvents	not recommended	Pressure	good	Vacuum	good	Viscous media	excellent	Sterile Media	limited	<table border="1"> <tr><td>Acids</td><td>good</td></tr> <tr><td>Alkaline solutions</td><td>good</td></tr> <tr><td>Solvents</td><td>not recommended</td></tr> <tr><td>Pressure</td><td>fair</td></tr> <tr><td>Vacuum</td><td>good</td></tr> <tr><td>Viscous media</td><td>excellent</td></tr> <tr><td>Sterile Media</td><td>limited</td></tr> </table>	Acids	good	Alkaline solutions	good	Solvents	not recommended	Pressure	fair	Vacuum	good	Viscous media	excellent	Sterile Media	limited	<table border="1"> <tr><td>Acids</td><td>good</td></tr> <tr><td>Alkaline solutions</td><td>good</td></tr> <tr><td>Solvents</td><td>not recommended</td></tr> <tr><td>Pressure</td><td>not recommended</td></tr> <tr><td>Vacuum</td><td>excellent</td></tr> <tr><td>Viscous media</td><td>good</td></tr> <tr><td>Sterile Media</td><td>excellent</td></tr> </table>	Acids	good	Alkaline solutions	good	Solvents	not recommended	Pressure	not recommended	Vacuum	excellent	Viscous media	good	Sterile Media	excellent	<table border="1"> <tr><td>Acids</td><td>excellent</td></tr> <tr><td>Alkaline solutions</td><td>excellent</td></tr> <tr><td>Solvents</td><td>good/excellent</td></tr> <tr><td>Pressure</td><td>–</td></tr> <tr><td>Vacuum</td><td>–</td></tr> <tr><td>Viscous media</td><td>–</td></tr> <tr><td>Sterile Media</td><td>–</td></tr> </table>	Acids	excellent	Alkaline solutions	excellent	Solvents	good/excellent	Pressure	–	Vacuum	–	Viscous media	–	Sterile Media	–	<table border="1"> <tr><td>Acids</td><td>excellent</td></tr> <tr><td>Alkaline solutions</td><td>excellent</td></tr> <tr><td>Solvents</td><td>excellent</td></tr> <tr><td>Pressure</td><td>not recommended</td></tr> <tr><td>Vacuum</td><td>good</td></tr> <tr><td>Viscous media</td><td>good</td></tr> <tr><td>Sterile Media</td><td>good</td></tr> </table>	Acids	excellent	Alkaline solutions	excellent	Solvents	excellent	Pressure	not recommended	Vacuum	good	Viscous media	good	Sterile Media	good
Acids	good																																																																										
Alkaline solutions	good																																																																										
Solvents	not recommended																																																																										
Pressure	good																																																																										
Vacuum	good																																																																										
Viscous media	excellent																																																																										
Sterile Media	limited																																																																										
Acids	good																																																																										
Alkaline solutions	good																																																																										
Solvents	not recommended																																																																										
Pressure	fair																																																																										
Vacuum	good																																																																										
Viscous media	excellent																																																																										
Sterile Media	limited																																																																										
Acids	good																																																																										
Alkaline solutions	good																																																																										
Solvents	not recommended																																																																										
Pressure	not recommended																																																																										
Vacuum	excellent																																																																										
Viscous media	good																																																																										
Sterile Media	excellent																																																																										
Acids	excellent																																																																										
Alkaline solutions	excellent																																																																										
Solvents	good/excellent																																																																										
Pressure	–																																																																										
Vacuum	–																																																																										
Viscous media	–																																																																										
Sterile Media	–																																																																										
Acids	excellent																																																																										
Alkaline solutions	excellent																																																																										
Solvents	excellent																																																																										
Pressure	not recommended																																																																										
Vacuum	good																																																																										
Viscous media	good																																																																										
Sterile Media	good																																																																										
Complies with the standards	FDA 21CFR175.300 US Pharmacopoeia Class VI	FDA 21CFR175.300	FDA 21CFR177.2600 US Pharmacopoeia Class VI NSF listed (Standard 51)	FDA	FDA 21CFR177.2600 USP Pharmacopoeia Class VI																																																																						
Sterilization	Autoclavable with steam 30 minutes at 1 bar (15 psi) and 121°C (250°F); tubing will appear milky Gas sterilization with Ethylene oxide Not recommended for sterilization with radiation.	Autoclavable with steam 30 minutes at 1 bar (15 psi) and 121°C (250°F); tubing will appear milky Gas sterilization with Ethylene oxide Not recommended for sterilization with radiation.	Autoclavable with steam 30 minutes at 1 bar (15 psi) and 121°C (250°F) Gas sterilization with Ethylene oxide / Sterilization with radiation up to 2.5 mrad / Caution: Use special tubing version (welded stoppers) when autoclaving 2 or 3-stop colour-coded tubing.	Autoclavable with steam 30 minutes at 1 bar (15 psi) and 121°C (250°F) Gas sterilization with Ethylene oxide Sterilization with radiation up to 2.5 mrad	Autoclavable with steam 30 minutes at 1 bar (15 psi) and 121°C (250°F) Gas sterilization with Ethylene oxide Sterilization with radiation up to 2.5 mrad Caution: Can not be repeatedly sterilized																																																																						
Gas permeability ¹	<table border="1"> <tr><td>CO₂</td><td>563</td></tr> <tr><td>O₂</td><td>124</td></tr> <tr><td>N₂</td><td>67</td></tr> </table>	CO ₂	563	O ₂	124	N ₂	67	<table border="1"> <tr><td>CO₂</td><td>360</td></tr> <tr><td>O₂</td><td>80</td></tr> <tr><td>N₂</td><td>40</td></tr> </table>	CO ₂	360	O ₂	80	N ₂	40	<table border="1"> <tr><td>CO₂</td><td>1200</td></tr> <tr><td>O₂</td><td>200</td></tr> <tr><td>N₂</td><td>80</td></tr> </table>	CO ₂	1200	O ₂	200	N ₂	80	<table border="1"> <tr><td>CO₂</td><td>1140</td></tr> <tr><td>O₂</td><td>76</td></tr> <tr><td>N₂</td><td>190</td></tr> </table>	CO ₂	1140	O ₂	76	N ₂	190	<table border="1"> <tr><td>CO₂</td><td>–</td></tr> <tr><td>O₂</td><td>–</td></tr> <tr><td>N₂</td><td>–</td></tr> </table>	CO ₂	–	O ₂	–	N ₂	–																																								
CO ₂	563																																																																										
O ₂	124																																																																										
N ₂	67																																																																										
CO ₂	360																																																																										
O ₂	80																																																																										
N ₂	40																																																																										
CO ₂	1200																																																																										
O ₂	200																																																																										
N ₂	80																																																																										
CO ₂	1140																																																																										
O ₂	76																																																																										
N ₂	190																																																																										
CO ₂	–																																																																										
O ₂	–																																																																										
N ₂	–																																																																										
Odour and Taste	none	none	low	–	–																																																																						
Toxicity	non-toxic	non-toxic	non-toxic and non-hemolytic	–	–																																																																						
Tubing life ³	<table border="1"> <tr><td>at 0 bar</td><td>800 hrs</td></tr> <tr><td>at 0.7 bar</td><td>700 hrs</td></tr> </table>	at 0 bar	800 hrs	at 0.7 bar	700 hrs	<table border="1"> <tr><td>at 0 bar</td><td>35 hrs</td></tr> <tr><td>at 0.7 bar</td><td>30 hrs</td></tr> </table>	at 0 bar	35 hrs	at 0.7 bar	30 hrs	<table border="1"> <tr><td>at 0 bar</td><td>1000 + hrs</td></tr> <tr><td>at 0.7 bar</td><td>1000 hrs</td></tr> </table>	at 0 bar	1000 + hrs	at 0.7 bar	1000 hrs	<table border="1"> <tr><td>at 0 bar</td><td>75 hrs</td></tr> <tr><td>at 0.7 bar</td><td>–</td></tr> </table>	at 0 bar	75 hrs	at 0.7 bar	–	<table border="1"> <tr><td>at 0 bar</td><td>800 + hrs</td></tr> <tr><td>at 0.7 bar</td><td>800 + hrs</td></tr> </table>	at 0 bar	800 + hrs	at 0.7 bar	800 + hrs																																																		
at 0 bar	800 hrs																																																																										
at 0.7 bar	700 hrs																																																																										
at 0 bar	35 hrs																																																																										
at 0.7 bar	30 hrs																																																																										
at 0 bar	1000 + hrs																																																																										
at 0.7 bar	1000 hrs																																																																										
at 0 bar	75 hrs																																																																										
at 0.7 bar	–																																																																										
at 0 bar	800 + hrs																																																																										
at 0.7 bar	800 + hrs																																																																										

¹ Permeability Coefficient = $\frac{\text{Amount of gas (cm}^3\text{)} \times \text{tubing wall thickness (cm)}}{\text{Surface area of tubing ID (cm}^2\text{)} \times \text{time (seconds)} \times \text{pressure drop across tubing wall (cmHg)}} \times 10^{-11}$

For more tubings for aggressive media, see page 35.

This information has been supplied to Ismatec by the tubing manufacturers. It is for your guidance only. We recommend you to test the tubing before use.

					
Tygon® HC F-4040-A	Tygon® SI Silicone 3350 (Platinum)	Silicone Peroxide	Norprene® A-60-G	Norprene® Chemical	Viton® Fluran® HCA F-5500-A
The special tubing for hydrocarbons, petroleum products and distillates.	The platinum-cured silicone tubing with an ultra-smooth inner surface for sanitary transfer of sensitive fluids.	Silicone tubing blended with organic peroxide for biological applications	The high performance tubing for industrial use.	Chemically resistant to acetone, MEK and other aggressive solvents. Long life tubing.	The special tubing for concentrated acids and corrosive solvents.
<ul style="list-style-type: none"> Specially formulated to transport hydrocarbons, petroleum products and distillates Ideal for gasoline, kerosene, heating oils, cutting liquids and coolants based on glycols High dielectric constant Low gas permeability 	<ul style="list-style-type: none"> Steam autoclaveability Excellent biological compatibility Ultra-smooth inner bore reduces potential for particle entrapment Lower level of protein binding. Entirely non-toxic, non-hemolytic and non-pyrogenic. Weather, ozone, sunlight and radiation resistant Resistant to fungus Odourless 	<ul style="list-style-type: none"> Steam autoclaveability Excellent biological compatibility Greater physical compression capability Not prone to mould Non-toxic Waterproof and resistant to ozone, radiation and sunlight Resistant to fungus Odourless 	<ul style="list-style-type: none"> Offers longest service-life with good flow consistency Good resistance to acids and alkaline chemicals Superior weathering Abrasion resistant Non-ageing and non-oxidizing Outstanding flexural fatigue resistance Low gas permeability versus rubber tubing Ozone (300 pphm) and UV light resistant Ideal for use in vacuum system 	<ul style="list-style-type: none"> Plasticizer-free Smooth innerbore Low sorption maintains fluid integrity Minimal adhesion and diffusion Suitable for MEK, Acetone and other corrosive solvents Long life tubing 	<ul style="list-style-type: none"> High chemical resistance Low gas permeability Wide temperature range
<ul style="list-style-type: none"> Not recommended for strong acids and alkalies, foodstuffs, beverages and medicines Potential leaching of plasticizers 	<ul style="list-style-type: none"> Not suitable for concentrated solvents, oils, acids or diluted sodium hydroxide Relatively high gas permeability 	<ul style="list-style-type: none"> Not recommended for concentrated solvent, oils, acids or diluted sodium hydroxide Relatively high gas permeability 	<ul style="list-style-type: none"> Potential leaching of blend material 	<ul style="list-style-type: none"> Cannot be repeatedly sterilized Only available as stopper tubing 	<ul style="list-style-type: none"> Limited service-life
<ul style="list-style-type: none"> Thermoplastic PVC-based material with plasticizer Flexible, firm, translucent, yellow 	<ul style="list-style-type: none"> Thermal set rubber Siloxane polymers and amorphous silica Soft, translucent, clear to light amber Excellent compression strength 	<ul style="list-style-type: none"> Polydimethylsiloxane with silica filler and silicone oil Excellent resistance to compression Soft, translucent, clear to light amber 	<ul style="list-style-type: none"> Thermoplastic elastomer based on polypropylene Excellent tensile strength Firm, opaque, black 	<ul style="list-style-type: none"> Special thermoplastic of high purity Without additives Without plasticizer Environmental-friendly disposal Flexible, firm, opaque 	<ul style="list-style-type: none"> Fluoropolymerelastomer Firm, opaque, black
-37°C to +74°C (-35°F to +165°F)	-60°C to +200°C (-75°F to +392°F)	-51°C to +238°C (-60°F to +460°F)	-60°C to +135°C (-75°F to +275°F)	-70°C to +74°C (-94°F to +165°F)	-31°C to +204°C (-25°F to +400°F)
limited	limited	limited	excellent	excellent	excellent
not recommended	limited	good	excellent	excellent	excellent
not recommended	limited	not recommended	not recommended	excellent	limited
good	not recommended	not recommended	not recommended	not recommended	not recommended
good	good	good	good	good	good
excellent	fair	fair	excellent	good	good
limited	excellent	excellent	not recommended	good	fair
None	US Pharmacopoeia XXIII Class VI FDA 21CFR177.2600 also exceeds 3A sanitary standards	US Pharmacopoeia Class VI FDA 21CFR177.2600	None	FDA 21CFR177.2600 USP Pharmacopoeia Class VI	None
Not recommended	Autoclaveable with steam 30 minutes at 1 bar (15 psi) and 121°C (250°F) Gas sterilization with Ethylene oxide Sterilization with radiation up to 2.5 mrad	Autoclaveable for 30 min at 1 bar (15 psi), 121°C (250°F) Radiation: Irradiate at up to 2.5 mrad Gas: Not recommended to sterilize with Ethylene oxide	Not recommended	Autoclaveable with steam 30 minutes at 1 bar (15 psi) and 121°C (250°F) Gas sterilization with Ethylene oxide Sterilization with radiation up to 2.5 mrad Caution: Can not be repeatedly sterilized	Not recommended
100	25147	25147	1200	–	38
22	4715	4715	200	–	14
12	2284	2284	80	–	5
²	none	–	²	–	²
²	non-toxic	–	²	–	²
60 hrs	200 hrs	–	1000 + hrs	800 + hrs	150 hrs
60 hrs	100 hrs	–	1000 hrs	800 + hrs	90 hrs

² Must not be used for foodstuffs, beverages, and drugs

³ Tubing 6.4 mm i.d., 1.6 mm wall, 3-roller pump head, 600 rpm, 23°C, service-life ending with rupture (ISMATEC® tubing pumps run at approx. 11 rpm – 500 rpm)

1

Tubing for special applications

From GORE™ for single channel tubing pumps



e.g. for the Flowmaster®
dispensing pump up to 13 liters/min

For high-pressure application
Gore Style 100

- Extremely stable flow rates
- Variability in flow rate within 1% during total life time
- Virtually eliminates spallation (ensures continuously high purity in fluid transfer)
- For differential pressures up to 4 bar (60 psi)
- High burst strength (up to 25 bar / 360 psi)
- Exhibits 18 times the life of silicone rubber tubing

Application
in pharmaceutical, food and biotech processes

- Tangential flow filtration and other high-pressure applications
- Addition of anti-foam
- Long-term fermentation: continuous media recirculation over 75 days
- Transfer of live-cells from one container into another featuring excellent service life at low temperatures
- Ultra-filtration: high pressure stability allows higher system pressure and flow rate, which results in longer service life and fewer down-times due to tube exchanges

Specifications

Unique pressed composite material, not extruded, produced in clean room

- Platinum cured Silicone and expanded PTFE
- Available in bore sizes up to 50 mm ID
- USP Class VI approved and classified nontoxic
- Cited in FDA Type II Material Master File (MMF)
- Operates at pressures up to 4 bar (60 psi)
- In-line steam sterilizable



e.g. for the tubing pump
REGLO Quick™ up to 230 ml/min

For aggressive media
Gore Style 100CR

- Extremely long life perfluoroelastomer tubing
- Stable flow rates, variation less than 1% over tubing life
- Low solvent swell
- Extreme long service life
- Suitable for almost all aggressive chemicals, including organic solvents, such as:
 - Methylenechloride
 - Toluene and Acetone

Application
in electronic, medical, textile, industry

- Solvent-based ink for gravure printing
- Coating of glass bottles
- Chemical coating of plastic plates and film
- Chemical-based flow in waste water treatment
- Solvent-based coating of tablets
- Synthesis with high through-put
- Laboratory analysis or dispensing

Specifications

Fluoroelastomer tubing with expanded PTFE, not extruded, produced in clean room

- 1,6 to 16 mm ID
- Permanently stable flow rates
- USP Class VI approved
- FDA for food contact
- Operates at pressures up to 4 bar (60 psi)

Special tubing for aggressive media

From Tygon® for single and multi-channel tubing pumps

1



Stopper tubing for aggressive media Tygon® MHLL

- Resistant to highly aggressive chemicals
- Meets USP Class VI criteria
- Low sorption maintains the fluid integrity
- Unequaled combination of chemical resistance, clarity and flexibility
- Color opaque, beige
- Extreme long service life (+800h)

Applications

- Battery acid filling
- Addition of anti-foam
- Hazardous material handling
- Applications with acids, bases, ketones, salts and alcohols

Its exceptionally smooth inner surface inhibits particulate buildup and reduces the potential for contamination.

Specifications

- Combination of MH- and Pharmed-tubing, opaque
- Available as stopper-tubing up to 2.79 mm i.d.
- USPXXIII Class VI
- FDA for food contact
- Autoclaveable, gas sterilisation, radiation
- Absolutely no odor or taste



Standard- and stopper tubing for chemical applications Tygon® MHSL

- Chemically resistant to a wide range of fluids
- Meets FDA criteria for food contact
- Plastizier and oil free
- No contamination of the fluids
- Color clear, transparent
- Long service life in peristaltic pumps (75h)

Applications

- Transparent for visible flow monitoring
- Coating of tablets
- Laboratory analysis and dispensing
- Chemical-based flow in waste water treatment

Smooth inner surface, low sorption maintains fluid and tube integrity

Specifications

- Ultra-pure tubing for peristaltic pumps
- Stopper tubing up to 2.79 mm i.d.
- Standard tubing up to 15.9 mm i.d.
- FDA for food contact
- Autoclaveable, gas sterilisation, radiation
- Absolutely no odor or taste

Pump tubing ordering information

1 Standard tubing for single channel pumps

i.d. mm	WT ¹ mm	Tygon® LFL	Tygon® ST R-3603/R-3607	PharMed® Ismaprene	Tygon® HC F-4040-A	Tygon® MHSL 2001	Tygon® SI Silicone 3350	Silicone Peroxide	Norprene® A-60-G	Norprene® Chemical	Viton® Fluran® HCA F-5500-A
0.8	1.6		MF0001	MF0009			MF0291	MF0044	MF0017		MF0048
1.6	1.6	SC0389	MF0028	MF0010	MF0002	SC0830	SC0580	MF0035	SC0357		MF0049
2.4	1.6		SC0691	SC1006			SC0590				SC0590
3.2	1.6	SC0390	MF0030	MF0012	MF0004	SC0831	SC0581	MF0037	SC0358	SC1022	MF0051
4.8	1.6	SC0391	SC0379	MF0011	MF0003	SC0832	SC0582	MF0045	SC0359	SC1023	MF0322
6.4	1.6	SC0392	MF0031	MF0013	MF0005	SC0833	SC0584	MF0046	SC0360	SC1024	MF0052
8.0	1.6	SC0394	MF0032	MF0014	MF0006	SC0834	SC0587	MF0047	SC0361		MF0053
9.5	1.6		SC0383			SC0835	SC0387		SC0385	SC1025	
11.1	1.6		SC0384				SC0697		SC0386		
4.8	2.4		MF0029	MF0448	MF0476		SC0583	MF0288	SC0362		MF0050
6.4	2.4		MF0033		MF0007		SC0585	MF0040	SC0363		MF0054
8.0	2.4		SC0502				SC0515		SC0511		
9.5	2.4		SC0503				SC0516		SC0512		
11.1	2.4		SC0504				SC0517				
12.7	2.4		SC0505				SC0518				
6.4	3.2	SC0393	SC0380	MF0015			SC0586	MF0314	SC0364		MF0323
9.5	3.2	SC0395	SC0381	MF0016	MF0008		SC0588	MF0041	SC0365		MF0055
12.7	3.2	SC0396	SC0382	MF0034	SC0725	SC0845	SC0589	MF0315	SC0366	SC1026	
15.9	3.2		SC0695	SC0696		SC0846	SC0532		SC0698		
9.5	6.0			MF0351			MF0359*	MF0355*			
12.7	5.0			MF0352			MF0360*	MF0356*			
19.0	6.0			MF0353			MF0361*	MF0357*			
25.4	4.0			MF0354			MF0362*	MF0358*			
Pack size		7.5 m	15 m	7.5 m	15 m	15 m	15 m	7.5 m	15 m	15 m	7.5 m

¹ WT = Wall thickness

(Other tubing sizes and wall thicknesses available on request.)

Tygon MHSL and Tygon MHLL special tubing for aggressive media

To improve the gliding properties and service-life of the tubing, we recommend you to periodically lubricate the pump rollers with silicone spray.

Order No. SC 0179




Extension tubing for 2 and 3-stop tubing

i.d. mm	WT* mm	o.d. mm	Tygon® ST R-3603/R-3607	PharMed® Ismaprene	Tygon® HC F-4040-A	Silicone Peroxide	Silicone Platin	Viton® Fluran® HCA F-5500-A	Tygon® MHSL 2001
0.13	0.91	1.95	SC0226						
0.19	0.91	2.01	SC0025						
0.25	0.91	2.07	SC0026	SC0337	SC0173				
0.38	0.91	2.20	SC0027	SC0338	SC0174				SC0854
0.44	0.91	2.26	SC0028						
0.51	0.91	2.33	SC0029	SC0339	SC0175		SC06040	SC0550	
0.57	0.91	2.39	SC0030						
0.64	0.91	2.46	SC0031	SC0340	SC0176	SC0448	SC06041	SC0551	SC0856
0.76	0.86	2.48	SC0032	SC0341	SC0177	SC0449	SC06042	SC0552	
0.89	0.86	2.61	SC0033	SC0342	SC0120	SC0450	SC06043	SC0553	
0.95	0.86	2.67	SC0034						
1.02	0.86	2.74	SC0035	SC0343	SC0121	SC0451	SC06044	SC0554	SC0858
1.09	0.86	2.81	SC0036						
1.14	0.86	2.86	SC0037	SC0344	SC0122	SC0452	SC06045	SC0555	
1.22	0.86	2.94	SC0038						
1.30	0.86	3.02	SC0039	SC0345	SC0123	SC0453	SC06046	SC0556	
1.42	0.86	3.14	SC0040	SC0346	SC0124	SC0454	SC06047	SC0557	
1.52	0.86	3.24	SC0041	SC0347	SC0125	SC0455	SC06048	SC0558	SC0860
1.65	0.86	3.37	SC0042	SC0348	SC0126	SC0456	SC06049	SC0559	
1.75	0.86	3.47	SC0043						
1.85	0.86	3.57	SC0044	SC0349	SC0127	SC0457	SC06050	SC0560	
2.06	0.86	3.78	SC0045	SC0350	SC0128	SC0458	SC06051	SC0561	SC0862
2.29	0.86	4.01	SC0046	SC0351	SC0129	SC0459	SC06052	SC0562	
2.54	0.86	4.26	SC0047	SC0352	SC0130	SC0460	SC06053	SC0563	
2.79	0.86	4.51	SC0048	SC0353	SC0131	SC0461	SC06054	SC0564	SC0864
3.17	0.86	4.89	SC0223						
Tube Length			10 m	3 m	3 m	15 m	10 m	10 m	10 m

* WT = Wall thickness

3-Stopper tubing for MS/CA cassettes





i.d. mm	Color-Code	Tygon® LFL iØ mm	Tygon® ST R-3603/R-3607	PharMed® Ismaprene	Tygon® HC F-4040-A	Tygon® MHSL 2001	Tygon® MHLL	Tygon® SI Silicone 3350	Silicone Peroxide	Viton® Fluran® HCA F-5500-A
0.13 ¹	orange-black		SC0189							
0.19 ¹	orange-red		SC0049							
0.25 ¹	orange-blue	SC0397 0.27	SC0050	SC0303	SC0730*	SC0286				
0.38 ¹	orange-green	SC0398 0.38	SC0051	SC0304		SC0287	SC0802**	SC0710**		
0.44 ¹	green-yellow		SC0052							
0.51 ¹	orange-yellow	SC0399 0.48	SC0053	SC0305	SC0731	SC0288		SC0600		SC0255
0.57 ¹	white-yellow		SC0054							
0.64 ¹	orange-white	SC0400 0.64	SC0055	SC0306		SC0289	SC0804**	SC0601	SC0106	SC0256
0.76 ²	black-black	SC0401 0.76	SC0056	SC0307		SC0290		SC0602	SC0107	SC0257
0.89 ²	orange-orange	SC0402 0.89	SC0057	SC0308	SC0732	SC0291		SC0603	SC0108	SC0258
0.95 ²	white-black		SC0058							
1.02 ²	white-white	SC0403 1.02	SC0059	SC0309	SC0737	SC0292	SC0806**	SC0604	SC0109	SC0259
1.09 ²	white-red		SC0060							
1.14 ²	red-red	SC0404 1.14	SC0061	SC0310		SC0293		SC0605	SC0110	SC0260
1.22 ²	red-grey		SC0062							
1.30 ²	grey-grey	SC0405 1.25	SC0063	SC0311	SC0733*	SC0294		SC0606	SC0111	SC0261
1.42 ²	yellow-yellow	SC0406 1.37	SC0064	SC0312		SC0295		SC0607	SC0112	SC0262
1.52 ²	yellow-blue	SC0407 1.52	SC0065	SC0313	SC0734	SC0296	SC0808**	SC0608	SC0113	SC0263
1.65 ²	blue-blue	SC0408 1.60	SC0066	SC0314		SC0297		SC0609	SC0114	SC0264
1.75 ²	blue-green		SC0067							
1.85 ²	green-green	SC0409 1.85	SC0068	SC0315		SC0298		SC0610	SC0115	SC0265
2.06 ²	purple-purple	SC0410 2.06	SC0069	SC0316	SC0735	SC0299	SC0810**	SC0611	SC0116	SC0266
2.29 ²	purple-black	SC0411 2.20	SC0070	SC0317		SC0300		SC0612	SC0117	SC0267
2.54 ²	purple-orange	SC0412 2.62	SC0071	SC0318		SC0301		SC0613	SC0118	SC0268
2.79 ²	purple-white	SC0413 2.79	SC0072	SC0319	SC0736*	SC0302	SC0812**	SC0614	SC0119	SC0269
3.17 ²	black-white		SC0224							
Pack size (pcs.)		12	12	6	12	6	6	6	6	12
Tube length (mm)		400	400	400	400	300	300	400	400	400

¹ Wall thickness 0.91 mm ² Wall thickness 0.86 mm

* Welded stoppers for autoclaving

**These tubes are only equipped with 2 stoppers

2-Stopper tubing for CA cassettes




i.d. mm	Color-Code	Tygon® LFL iØ mm	Tygon® ST R-3603/R-3607	PharMed® Ismaprene	Tygon® HC F-4040-A	Tygon® MHSL 2001	Tygon® MHLL	Tygon® SI Silicone 3350	Silicone Peroxide	Viton® Fluran® HCA F-5500-A
0.13 ¹	orange-black		SC0188							
0.19 ¹	orange-red		SC0001							
0.25 ¹	orange-blue	SC0414 0.27	SC0002	SC0320	SC0740*	SC0156				
0.38 ¹	orange-green	SC0415 0.38	SC0003	SC0321		SC0157	SC0814**	SC0716**		
0.44 ¹	green-yellow		SC0004							
0.51 ¹	orange-yellow	SC0416 0.48	SC0005	SC0322	SC0741	SC0158		SC0620		SC0132
0.57 ¹	white-yellow		SC0006							
0.64 ¹	orange-white	SC0417 0.64	SC0007	SC0323		SC0159	SC0816**	SC0621	SC0092	SC0133
0.76 ²	black-black	SC0418 0.76	SC0008	SC0324		SC0160		SC0622	SC0093	SC0134
0.89 ²	orange-orange	SC0419 0.89	SC0009	SC0325	SC0742	SC0161		SC0623	SC0094	SC0135
0.95 ²	white-black		SC0010							
1.02 ²	white-white	SC0420 1.02	SC0011	SC0326	SC0747	SC0162	SC0818**	SC0624	SC0095	SC0136
1.09 ²	white-red		SC0012							
1.14 ²	red-red	SC0421 1.14	SC0013	SC0327		SC0163		SC0625	SC0096	SC0137
1.22 ²	red-grey		SC0014							
1.30 ²	grey-grey	SC0422 1.25	SC0015	SC0328	SC0743*	SC0164		SC0626	SC0097	SC0138
1.42 ²	yellow-yellow	SC0423 1.37	SC0016	SC0329		SC0165		SC0627	SC0098	SC0139
1.52 ²	yellow-blue	SC0424 1.52	SC0017	SC0330	SC0744	SC0166	SC0820**	SC0628	SC0099	SC0140
1.65 ²	blue-blue	SC0425 1.60	SC0018	SC0331		SC0167		SC0629	SC0100	SC0141
1.75 ²	blue-green		SC0019							
1.85 ²	green-green	SC0426 1.85	SC0020	SC0332		SC0168		SC0630	SC0101	SC0142
2.06 ²	purple-purple	SC0427 2.06	SC0021	SC0333	SC0745	SC0169	SC0822**	SC0631	SC0102	SC0143
2.29 ²	purple-black	SC0428 2.20	SC0022	SC0334		SC0170		SC0632	SC0103	SC0144
2.54 ²	purple-orange	SC0429 2.62	SC0023	SC0335		SC0171		SC0633	SC0104	SC0145
2.79 ²	purple-white	SC0430 2.79	SC0024	SC0336	SC0746*	SC0172	SC0824**	SC0634	SC0105	SC0146
3.17 ²	black-white		SC0222							
Pack size (pcs.)		12	12	6	12	6	6	6	6	12
Tube length (mm)		400	400	400	400	381	381	400	400	180

¹ Wall thickness 0.91 mm ² Wall thickness 0.86 mm






* Welded stoppers for autoclaving

**These tubes are only equipped with 2 stoppers

Tube connectors


for ISMATEC® pump tubing

Standard tube connectors in plastic



Tubing i.d. mm	Order No.	Order No.	Order No.	Order No.	Order No.
1.5	ISM 566 ¹	ISM 558 ¹		ISM 694 ²	
2.5	ISM 557 ¹	ISM 693 ¹		ISM 510 ¹	
3.0	ISM 558	ISM 510	ISM 524	ISM 493	
4.0	ISM 559	ISM 511	ISM 525	ISM 494	ISM 535
5.0	ISM 560	ISM 512	ISM 526	ISM 495	ISM 536
6.0	ISM 561	ISM 513	ISM 527	ISM 496	ISM 537
7.0		ISM 514			
8.0	ISM 562	ISM 515	ISM 528	ISM 497	
10.0	ISM 563	ISM 516	ISM 529	ISM 498	
12.0	ISM 564	ISM 517	ISM 530	ISM 499	
13.0	ISM 565	ISM 518	ISM 531	ISM 500	
14.0	ISM 566	ISM 519	ISM 532	ISM 501	
15.0		ISM 520		ISM 502	
16.0	ISM 567	ISM 521	ISM 533	ISM 503	
19.0	ISM 568	ISM 522	ISM 534	ISM 523	
Pack size	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces

Material: Polypropylene except ¹ = Acetal, ² = White Nylon






How to make connections with tube sleeves

- The tubing material of the tube sleeves should be equal to the tubes to be connected.
- For the tube sleeve choose an inner diameter that equals the outer diameter of the tubes to be connected.
- If possible, clean the tube bore thoroughly with compressed air.

Adhesive for connections with tube sleeves

ergo® Adhesive for all tubing materials Complete Set including:	Content	Order number
- ergo® Adhesive	20 g	SC 1010
- ergo® Primer	10 ml	
ergo® Adhesive (separate) type 5861 Porous	20 g	SC 1011
ergo® Primer (separate) type 5150 Porous important for pretreatment of the tubing	10 ml	SC 1012


Reducer tube connectors in plastic




Tubing i.d. mm			Tubing i.d. mm			Tubing i.d. mm		
A	B	Order No	A	B	Order No	A	B	Order No
1.5	2.5	ISM 569 ¹	3	4	ISM 538	4	6	ISM 553
1.5	3.0	ISM 570 ¹	4	6	ISM 539	6	8	ISM 554
2.5	3.0	ISM 571 ¹	6	4	ISM 540			
3.0	4.0	ISM 572	8	4	ISM 541			
4.0	6.0	ISM 573	8	6	ISM 542			
4.0	8.0	ISM 574	10	6	ISM 544			
6.0	8.0	ISM 575	10	8	ISM 545			
6.0	10.0	ISM 576	10	13	ISM 546			
8.0	10.0	ISM 577	12	8	ISM 547			
8.0	12.0	ISM 578	12	10	ISM 548			
10.0	12.0	ISM 579	15	6	ISM 549			
Pack size			10 pieces			10 pieces		

Material: Polypropylene except ¹ = Acetal

Steel 18/8 Standard tube connectors



Tubing i.d. mm	Tubing o.d. mm	Connector length	Order No.
0.30	.63	x15 mm	ISM 580
0.58	.90	x15 mm	ISM 581
0.58	.90	x11 mm	ISM 582
0.84	1.27	x11 mm	ISM 583
0.84	1.27	x16 mm	ISM 584
0.30	.63	x25 mm	ISM 585
0.58	.90	x25 mm	ISM 586
0.58	.90	x19 mm	ISM 587
Pack size			6 pieces



Contents of tubing selection guide:


- Material and tubing properties
- Chemical resistance chart
- Ordering information

Gear Pumps

Pulsefree pumping

Gear pumps allow differential pressures up to maximum 5.6 bar

2



ISMATEC® gear pumps run only in the clockwise direction (Exception REGLO-Z *Digital*)



Safe and easy to operate

- Developed for continuous duty, 24 hours a day, 7 days a week
- Compact drives with hermetically sealed and magnetically coupled pump-heads
- Safe overload protection – magnetically driven pump-heads decouple when load exceeds the coupling torque
- Internal bypass valve limits the differential pressure
- Pump-heads are interchangeable within seconds
- MAX key enables rapid filling of the system (BVP-Z and MCP-Z pumps)

Application range of gear pumps

Industries	Applications	Special media
<ul style="list-style-type: none"> - Biotechnological - Chemical - Food - Mining - Power - Pulp and Paper - Semiconductor - Textile 	<ul style="list-style-type: none"> - Sampling - Refrigeration technology - Water treatment - Liquid chromatography - Surface treatment - Distillation systems 	<ul style="list-style-type: none"> - Biozides - Dye stuffs - Thixotropic products - Liquid waxes - Hydrogen peroxide - Flux <p>Not suited for media containing particulates</p>

Multifunctional

- Interchangeable pump-heads for different flow rates available in specific, media-resistant materials
- Virtually no pulsation
- Very accurate dispensing pumps due to calibrateable drives
- REGLO-Z *Digital* with reversible rotation direction
- Specially designed pump-heads (Suction Shoe Design) for elevated differential pressures
- Excellent media compatibility stainless steel housing gears available in PTFE, Graphite, PPS or PEEK™
- Pump-heads for media with elevated viscosities

Low operation costs

- Interchangeable, magnetically coupled pump-heads
- Maintenance-free drives
- Only few wearing parts (gears, seals)
- Service kits allow the user to exchange worn parts
- High quality and precision for an optimum performance even after many years of intensive use

MICROPUMP®

ISMATEC gear pumps are

- Easy to service
- Almost maintenance-free
- Leak-free
- Differential pressure up to 5.6 bar

Note

All microprocessor controlled drives are LabVIEW compatible and can easily be integrated into process control systems.

Overview of gear pumps

Flow rates and models

Flow rates	ml/min		Bar	Model	Page
	min.	max.			
1	3290		5.2	REGLO-Z <i>Analog, Digital</i>	44
1	3290		5.2	REGLO-ZS <i>Analog, Digital</i>	44
1	7241		5.2	BVP-Z <i>Standard</i>	46
1	7241		5.2	MCP-Z <i>Standard</i>	47
1	7241		5.2	MCP-Z <i>Process</i>	48

Unique!

Only the ISMATEC® gear pump

MCP-Z *Process features:*



Carrying out programs independently of a PC

- Create the application profile in the PC (with ProgEdit software, Page 61)
- Download the file data into the pump memory
- Disconnect the pump from the PC
- Carry out your application on the spot, using the pump as a stand-alone unit



This mark indicates dispensing function

(You'll find dispensing gear pumps on pages 44, 47, 48)

- Pumping by speed or flow rate
- Dispensing by volume or time
- Interval dispensing by volumes with a pause
- Interval dispensing by time with a pause
- Programming a number of dispensing cycles
- **Calibrating** the flow rate and dispensing volume
- Factory set gear pump-heads

NEW pump-head reference chart

The Micropump pump-head numbers have recently been updated to new series product code to simplify the global understanding of the products we provide. Below, you'll find an easy reference chart updating the old Pump-head numbers to the new series product code. Please note that ISMATEC order numbers for Micropump pump-heads have now been updated to the Micropump order number.

OLD Order No.	OLD Pump-head No.	NEW Pump-head No.	Order No.
Series GA			
MI0006	Z-186	GA-X21.CFS.B	82092
MI0007	Z-181	GA-V21.CFS.B	82114
MI0008	Z-183	GA-V23.CFS.B	82115
MI0131	Z-1830	GA-T23.PFS.B	81473
MI0312	Z-186 P	GA-X21.JFS.B	L20820
MI0280	Z-1830 P	GA-T23.JFS.B	L18489
MI0309	Z-186 HC	GA-X21.CFC.B	L17164
MI0310	Z-183 HC	GA-V23.CFC.B	L2383
Series GB			
MI0022	Z-200	GB-P25.PVS.A.B1	81281
MI0023	Z-201	GB-P35.PVS.A.B1	81281
MI0306	Z-200 P	GB-P25.JVS.B	220004
MI0378	Z-201 PKC	GB-P35.JKS.B	L22609
Series GJ			
MI0013	Z-120	GJ-N23.FF1S.B.B1	82004
MI0015	Z-122	GJ-N25.FF1S.B.B1	82006
MI0016	Z-140	GJ-N23.FF1S.B	82001
MI0018	Z-142	GJ-N25.FF1S.B	82003
MI0019	Z-130	GJ-N23.PF1S.B.B1	81529
MI0020	Z-150	GJ-N23.PF1S.B	81531
MI0313	Z-140 P	GJ-N23.JF1S.B	L197735
MI0284	Z-140 HC	GJ-N23.FF1C.B	L20284
MI0311	Z-142 HC	GJ-N25.FF1C.B	L21812

Selection criteria

Find the optimum pump-head design

	 Cavity Style	 Suction Shoe
Flow performance and pressure		
Only flow	✓	✓
Pre-pressure necessary	—	✓
Back-pressure		
Flow rate stable	—	✓
Back-pressure high		
BVP-Z and MCP-Z drives	✓	✓
Forward and reverse delivery	✓	—
Bypass depending on pump-head/Series	GJ	GB
Range of flow rates (ml/min)		
Series GA	1–560	—
Series GJ	33–3950 (Reglo-Z)	✓
Series GJ	55–5480 (MCP-Z)	✓
Series GB	35–7241	—
Max. operating temperature		
54°C (129°F), 77°C (170°F), 99°C (210°F)	✓	—
(depending on the seals)		
Up to 177°C (350°F)	—	✓
Max. suction height varies (depends on pump-head, speed and tubing)		
for water		
wetted gears		
1 m	✓	—
30 cm	✓	✓
flooded		
8 m	✓	—
3 m	✓	✓
Pumping out of vacuum		
recommended up to 200 mbar		
absolutely	✓	—
not suitable	—	✓
Viscosity		
0.2 to 1500 cp	✓	✓
max 2000 cp, depending on pump-head	✓	—
Particles up to 5 µm	Z-150WI	—
Gear material		
PTFE	✓	—
Graphite	—	✓
PEEK™	✓	✓
PPS	✓	✓
NiC	Z-150WI	—

Cavity style:

Series GJ



- Max. suction height with water and flooded pump-head: 8 m, depending on pump-head and tubing
- Pumping out of a vacuum of 200 mbar
- Based on the traditional gear pump technology
- For application with moderate differential pressure

In comparison to the Suction Shoe pump-heads, the Cavity style pump-heads can be used for viscous media and applications with a certain suction height

Advantages:

- Excellent chemical resistance
- Smooth operation at a low noise level
- Low internal friction



Suction shoe style:

Series GA and GB



- An exclusive Micropump product featuring a patented technology
- Modified pump chamber compared to the conventional gear pump technique

This type of pump-head design has a seal plate mounted with a deliberate play in the suction part of the pump chamber (hence the expression Suction Shoe). Discharge pressure keeps the Suction Shoe seated tightly on top of the gears which prevents flow from decreasing in high-pressure applications.

Advantages:

- Temperature range from –46 – 177 °C (–51 – 350°F)
- The Suction Shoe acts as a dynamic seal element which results in a temperature- and pressure-independent pump chamber.
- Ease of servicing due to fewer parts. The service kit, including the Suction Shoes, enables an extended pump life; conventional pumps require more frequent replacement.




Main features of pump-head designs

	Reverse delivery		Main sealing zone		Max. suction height - dry		Max. suction height flooded		Recommended for		Temperature range		Differential pressure		Suction from vacuum		Wearing material		Pulseless helical gearing		Pulseless spur gearing	
Cavity Style 	Yes	Tooth edges	1 m	8 m	High flow rate	54°C, 77°C, 99°C depending on seals	only low	200 mbar	– Service kit – Cavity plate	Head: GJ												
Suction Shoe 	No	Front surfaces of gears	30 cm	3 m	Preferred for diff. pressure	from –46 to 177°C (–29°C 200 series)	5.6 bar 8.7 bar*	not suitable	Service kit incl. suction shoe	Head: GB	Head: GA											

*Pump-head for Industrial drive provided by customer

The magnetically coupled drive principle

Consists of two magnets, a driving magnet that attaches to the motor shaft and a driven magnet that is completely sealed within the pump-head and is connected to the driving gear. The driven magnet is a wetted component and is totally encapsulated.



The two magnets couple automatically such that the driving magnet turns the driven magnet and gears without physical contact.

Decoupling occurs when the pump load exceeds the coupling torque between the two magnets. This feature can act as a safety device to prevent damage to the pump and motor as well as associated piping. The magnets can be recoupled by bringing the motor to a complete stop, then eliminating the cause of the decoupling and restarting.

Pump-head material options

Enhance the chemical compatibility and application potential

- Base material
Standard: Stainless steel 316
Options: e.g. Hastelloy® B2, Hastelloy C-276, Alloy 20 and Titan
- Gears
Standard: PPS, Graphite, PTFE (depends on pump-head)
Options: e.g. PEEK™, PPSKV
- Static seals
Standard: Viton®, PTFE (depends on pump-head)
Options: EP, Buna N, Kalrez®
- Magnets
Standard: Ferrite
Options: e.g. SmCo, NdFeB

PTFE = Polytetrafluoroethylene
PPS = Polyphenylenesulphide
PEEK = Polyetheretherketone

Internal bypass

- An adjustable fluid bypass valve helps protect against decoupling and system damage from high-pressure build-up
- It allows for adjustment of a max. differential pressure (from 0.7 bar up to the max. differential pressure, depending on the individual pump-head)
- Should only be used for safety purposes and not for pressure controlling (bypass conditions may create a sufficient temperature rise to cause significant swelling in PTFE-gear pumps)



MICROPUMP®

Further pump-head options

- Integral drive
- High system pressure
- Deck ports
- 1/4-18 NPT ports
- Tri-clamp fittings

2

REGLO-Z, REGLO-ZS

Compact and powerful
Footprint only 10 by 18 cm!

2



REGLO-Z Analog

REGLO-ZS Analog

REGLO-Z Analog
1–3290 ml/min

- Variable speed
- Differential pressure max. 5.2 bar

REGLO-ZS
Drive and pump-head are separated by a 2 m long cable.

- Pulseless fluid delivery
- 10 cm wide, 13.5 cm high
- Reversible rotation
(with Cavity Style Pump-heads)
- Interchangeable Micropump® pump-heads
- Excellent repeatability

Repetitive error (rel.)	5 ml	0.5%	(REGLO-Z Digital)
	20 ml	<0.2%	
	100 ml	<0.1%	



Interchangeable pump-heads



with Cavity Style pump-head, rotation direction is reversible



with Suction shoe pump-head, run only in the clockwise direction

Dispensing and calibrating function see Page 41



REGLO-Z Digital



REGLO-ZS Digital

REGLO-Z Digital
1–3290 ml/min

- with dispensing functions
- Membrane key-pad
 - LED display with setting menu
 - Differential pressure max. 5.2 bar



Specifications REGLO-Z/ZS Analog

Motor type	DC motor
Speed	50 – 5000 rpm
Speed setting	1–99%, resolution 1% 2-digit potentiometer
Power consumption	50 W
Mains connection	230V _{AC} /50Hz, 115V _{AC} /60Hz adjustable
Protection rating	IP 30
Depth/Width/Height	
Drive REGLO-Z	178 x 100 x 143 mm
Drive REGLO-ZS	175 x 65 x 80 mm
External control unit	178 x 100 x 143 mm
Weight	
Drive REGLO-Z	2.1 kg (without pumphead)
Drive REGLO-ZS	0.7 kg (without pumphead)
External control unit	1.7 kg

Specifications REGLO-Z/ZS Digital

Motor type	DC-Motor
Speed range	50 – 5000 rpm
Speed setting	rpm, resolution 1 rpm
Flow rate setting	ml/min, liters/min
Power Consumption	75 W
Mains connection	100 – 230V _{AC} /50 – 60Hz
Protection rating	IP 30
Depth/Width/Height	
Drive REGLO-Z	178 x 100 x 135 mm
Drive REGLO-ZS	175 x 65 x 80 mm
External control unit	178 x 100 x 135 mm
Weight	
Drive REGLO-Z	1.7 kg (without pumphead)
Drive REGLO-ZS	0.7 kg (without pumphead)
External control unit	1.2 kg

Ordering Information

Model	Order No.
REGLO-Z Analog	ISM 895
REGLO-ZS Analog	ISM 896
REGLO-Z Digital	ISM 901
REGLO-ZS Digital	ISM 1143
Foot switch	ISM 891
Pump-head	Page 45
2 Nozzles	Page 45

LabVIEW drivers for Reglo-Z / -ZS Digital download for free: www.ismatec.com

Never use a gear pump for media containing particulates.

Interfaces

MICROPUMP®

Ordering Information pump-heads for REGLO-Z / -ZS

Pump-heads »Suction shoe«

- Enhanced pumping performance at elevated differential pressures
- Suited for elevated temperature ranges
- Not recommended for applications requiring a suction lift



REGLO Analog

- Speed control (0–5 or 0–10 V, 0–20 or 4–20 mA)
- Speed output
2-channel: 0–8 kHz
4-channel: 0–5 kHz
- Start/Stop
- Rotation direction



REGLO Digital

- PC-controllable
- Analog: only speed output (see Reglo *Analog*), start/stop and autostart

Suction Shoe	Pump-head No.	Order No.	Flow rate (ml/min)		Differential pressure max. bar	Gear material	Seals	Stainless steel housing	System pressure, max. (bar)	Temperature range °C	Internal Bypass
			min.	max.							
	GA-X21.CFS.B	MI0006	0.85	85	1.4	Graphite	PTFE	SS316	21	-46 – +177	–
	GA-X21.JFS.B	MI0312	0.85	85	2.3	PEEK*	PTFE	SS316	21	-46 – +177	–
	GA-V21.CFS.B	MI0007	2.1	210	2.8	Graphite	PTFE	SS316	21	-46 – +177	–
	GA-V23.CFS.B	MI0008	4.2	420	2.8	Graphite	PTFE	SS316	21	-46 – +177	–
	GA-V23.PFS.B	MI0131	4.6	460	5.2	PPS	PTFE	SS316	22	-46 – +177	–
	GA-V23.JFS.B	MI0280	4.6	460	5.2	PEEK	PTFE	SS316	22	-46 – +177	–
For corrosive media	GA-X21.CFC.B	MI0309	0.85	85	1.4	Graphite	PTFE	Hastelloy®-C276	21	-46 – +177	–
	GA-V23.CFC.B	MI0310	4.2	420	2.8	Graphite	PTFE	Hastelloy-C276	21	-46 – +177	–

- Ports (internal thread) 1/8"-27NPT
- Flow rates without differential pressure
- Operating temperature: with other seals up to 99°C possible

Pump-heads »Cavity style«

- Can be used for viscous media and applications requiring a certain suction lift
- Excellent chemical resistance
- Smooth and precise flow

Cavity Style	Pump-head No.	Order No.	Flow rate (ml/min)		Differential pressure max. bar	Gear material	Seals	Stainless steel housing	System pressure, max. (bar)	Temperature range °C	Internal Bypass
			min.	max.							
	GJ-N23.FFS.B.B1	MI0013	32	3200	1*	PTFE	PTFE	SS316	21	-46 – +54	✓
	GJ-N23.FFS.B	MI0016	32	3200	1*	PTFE	PTFE	SS316	21	-46 – +54	–
	GJ-N23.JFS.B	MI0313	32	3200	1*	PEEK	PTFE	SS316	21	-46 – +54	–
	GJ-N23.JFS.B.B1	MI0019	32	3200	1*	PPS	PTFE	SS316	21	-46 – +54	✓
	GJ-N23.JFS.B	MI0020	32	3200	1*	PPS	PTFE	SS316	22	-46 – +54	–
For corrosive media	GJ-N23.FFC.B	MI0284	32	3200	1*	PTFE	PTFE	Hastelloy-C276	21	-46 – +54	–
For abrasive media	GJ-N23.9FD.B	MI0265	32	3200	1*	NiC	PTFE	Surface hardened	21	-46 – +54	–

- Ports (internal thread) 1/8"-27NPT
- Flow rates without differential pressure
- Operating temperature: with other seals up to 99°C possible

* For applications with differential pressures exceeding 1 bar we recommend using the MCP-Z drive.



Service Kits contain the wearing parts (bushings, seals, gears)



These pump-heads are also available as OEM versions. Ask for the detailed data sheet.

Tubing adaptors for gear pump-heads

Threaded stainless steel connectors

Order No.	External Tubing thread adaptor	Tubing i.d. mm
AR0001	1/8" NPT Tube nozzle	6
AR0002	1/8" NPT Tube nozzle	3
AR0004	3/8" NPT Tube nozzle	12
AR0008	1/8" NPT Tube nozzle	8
AR0009	1/8" NPT Tube nozzle	9.5
AR0024	1/8" NPT Pipe connection	6 (outside)

Threaded connectors in Hastelloy-C

AR0001-HC	1/8" NPT Tube nozzle	6
-----------	----------------------	---

BVP-Z Standard

Economical

Robust, powerful gear pump drive

2

- Variable speed (no dispensing functions)
- Pulseless pumping
- Up to 5.2 bar differential pressure

Interfaces



- Speed control (0–5 or 0–10V, 0–20 or 4–20mA)
- Speed output (0–10V_{DC} or 0–12 kHz)
- Start/Stop



BVP-Z Standard

without dispensing functions

- **3-digit potentiometer** (for speed setting)
- 12 interchangeable Micropump® pump-heads

Flow rates and differential pressure depend on the pump-head mounted

CE

BVP-Z Standard
with interchangeable gear pump-heads
(see Page 49)
(material options, see Page 43)

Specifications

Motor type	DC motor
Speed	60 – 6000 rpm
Speed setting	1–99.9%, resolution 0.1% 3-digit potentiometer
Power consumption	150 W
Mains connection	230V _{AC} /50Hz, 115V _{AC} /60Hz adjustable
Protection rating	IP 30
Depth/Width/Height	220 x 155 x 260 mm (without pump-head)
Weight	5.7 kg (without pump-head)

Ordering information

The complete pump system BVP-Z Standard consists of:

Drive (magnet included)	ISM 446
Pump-head	Page 49
2 Nozzles	Page 49
Accessories	Page 61
– Foot switch	ISM 891
– Valve	on request

Applications

- Single-channel delivery processes under pressure for particulate-free fluids, e.g.: addition of reagents/solvents in organic synthesis at laboratory scale.
- Pumping propylene oxide into a laboratory reactor with a dispensing precision of +/-1% and a differential pressure of up to max. 3 bar.

MCP-Z *Standard*

Multi-purpose

Saves individual application parameters!

- Microprocessor controlled
- Ideal for dispensing and filling
- Pulseless pumping
- Robust, powerful gear pump drive
- Up to 5.2 bar differential pressure



MCP-Z *Standard*
with interchangeable gear pump-head
(see Page 49)
(material options, see Page 43)

MCP-Z *Standard*

- with dispensing functions
- Membrane key-pad, LED display
 - **4 program memories for saving individual application parameters**
 - 12 interchangeable Micropump® pump-heads (pre-programmed)
- Flow rates and differential pressure depend on the pump-head mounted



Interfaces



- PC-controllable:
- RS232



- Speed control (0–5 or 0–10V, 0–20 or 4–20mA)
- Speed output (0–10V_{DC} or 0–12 kHz)
- Start/Stop
- Autostart



Dispensing and calibrating
function see Page 41

Specifications

Motor type	DC motor
Speed	60 – 6000 rpm
Speed setting	rpm, resolution 1 rpm
Flow rate setting	µl/min, ml/min, liters/min
Power consumption	150 W
Mains connection	230V _{AC} /50Hz, 115V _{AC} /60Hz adjustable
Protection rating	IP 30
Depth/Width/Height	220 x 155 x 260 mm (without pump-head)
Weight	6.4 kg (without pump-head)

Ordering information

The complete pump system MCP-Z *Standard* consists of:

Drive (magnet included)	ISM 405
Pump-head	Page 49
2 Nozzles	Page 49
Accessories	Page 61
– Foot switch	IS 10039

Applications

Single-channel delivery and dispensing processes of particulate-free fluids under pressure.

With pump-heads GJ-N23 and GA-X21

Pulseless dispensing under pressure of different reagents with 2 pumps in different quantity ratios via a mixing valve into a reactor.

MCP-Z Process

Programmable

Programs can be carried out on the spot independently of a PC! Protection rating of IP 65

2

- Suitable for industries, extremely robust gear pump drive
- For pulseless pumping (up to 5.2 bar)
- Ideal for dispensing and filling applications in a dusty, humid or corrosive environment, and in clean room areas (IP 65, dust-tight and protected against water jets)



CE

MCP-Z Process
with interchangeable gear pump-heads
(material options, see Page 43)

MCP-Z Process

- Stainless steel housing
 - Membrane key-pad with LED display
 - **4 program memories for saving individual application parameters or PC programmed command sequences**
 - pre-programmed pump-heads
 - 21 interchangeable Micropump® pump-heads
- Flow rates and differential pressure depend on the pump-head mounted



Interfaces



PC-controllable:
– RS232



- Speed control (0–5 or 0–10V, 0–20 or 4–20mA)
- Speed output (0–10V_{DC} or 0–12 kHz)
- Start/Stop
- Autostart
- 2 universal inputs
- 2 universal outputs



Software ProgEdit
LabVIEW drivers
Free download on www.ismatec.com



Dispensing and calibrating
function see page 41

Specifications

Motor type	DC motor
Speed	60–6000 rpm
Speed setting	rpm, resolution 1 rpm
Flow rate setting	µl/min, ml/min, liters/min
Power consumption	200 W
Mains connection	100 – 230 V _{AC} / 50 – 60 Hz
Protection rating	IP 65
Depth/Width/Height	260 x 160 x 262 mm (without pump-head)
Weight	6.9 kg (without pump-head)

Ordering information


The complete pump system MCP-Z Process consists of:	
Drive (magnet included)	ISM 918
Pump-head	Page 49
2 Nozzles	Page 49
Accessories	Page 61
– Software ProgEdit	free download
– Foot switch	IS 10039

LabVIEW driver
download for free: www.ismatec.com

Applications

- Single-channel delivery and dispensing processes under pressure, for particulate-free solutions
- Addition of various reagents in different quantity ratios via mixing valve into reactor

Ordering Information pump-heads for BVP-Z / MCP-Z



ISMATEC gear pumps run only in the clockwise direction.
Never use a gear pump for media containing particulates.

Pump-heads »Suction shoe«

- Enhanced pumping performance at elevated differential pressures
- Suited for elevated temperature ranges
- Not recommended for applications requiring a suction lift

Suction Shoe	Pump-head No.	Order No.	Flow rate (ml/min)		Differential pressure max. bar	Gear material	Seals	Stainless steel housing	System pressure, max. (bar)	Temperature range °C	Internal Bypass
			min.	max.							
	GA-X21.CFS.B	MI0006	1	99	1.4	Graphite	PTFE	SS316	21	-46 – +177	–
	GA-X21.JFS.B	MI0312	1	99	2.3	PEEK™	PTFE	SS316	21	-46 – +177	–
	GA-V21.CFS.B	MI0007	3	252	2.8	Graphite	PTFE	SS316	21	-46 – +177	–
	GA-V23.CFS.B	MI0008	5	504	2.8	Graphite	PTFE	SS316	21	-46 – +177	–
	GA-V23.PFS.B	MI0131	6	560	5.2	PPS	PTFE	SS316	21	-46 – +177	–
	GA-V23.JFS.B	MI0280	6	560	5.2	PEEK	PTFE	SS316	21	-46 – +177	–
	GB-P25.PVS.B	MI0022	35	3509	3.5	PPS	Viton	SS316	21	-29 – +177	✓
	GB-P25.JVS.B	MI0306	35	3480	3.5	PEEK	Viton	SS316	21	-29 – +177	–
	GB-P35.PVS.B	MI0023	70	7020	3.5	PPS	Viton	SS316	21	-29 – +177	✓
Organic solvents	GB-P35.JKS.B	MI0378	73	7241	3.5	PEEK	Kalrez®	SS316	21	-29 – +177	✓
For corrosive media	GA-X21.CFC.B	MI0309	1	99	1.4	Graphite	PTFE	Hastelloy®-C276	21	-46 – +177	–
	GA-V23.CFC.B	MI0310	5	504	2.8	Graphite	PTFE	Hastelloy-C276	21	-46 – +177	–

- Ports (internal thread) 1/8" -27NPT
- Flow rates without differential pressure
- Operating temperature: with other seals up to 99°C possible



Service Kits

contain the wearing parts (bushings, seals, gears)

Pump-heads »Cavity style«

- Can be used for viscous media and applications requiring a certain suction lift
- Excellent chemical resistance
- Smooth and precise flow

Cavity Style	Pump-head No.	Order No.	Flow rate (ml/min)		Differential pressure max. bar	Gear material	Seals	Stainless steel housing	System pressure, max. (bar)	Temperature range °C	Internal Bypass
			min.	max.							
	GJ-N23.FFS.B.B1	MI0013	40	3950	3.5	PTFE	PTFE	SS316	21	-46 – +54	✓
	GJ-N23.FFS.B	MI0016	40	3950	3.5	PTFE	PTFE	SS316	21	-46 – +54	–
	GJ-N23.JFS.B	MI0313	40	3950	5.6	PEEK	PTFE	SS316	21	-46 – +54	–
	GJ-N25.FFS.B	MI0018	55	5460	3.5	PTFE	PTFE	SS316	21	-46 – +54	–
	GJ-N23.JFS.B.B1	MI0019	40	3950	5.2	PPS	PTFE	SS316	21	-46 – +54	✓
	GJ-N23.JFS.B	MI0020	40	3950	5.2	PPS	PTFE	SS316	22	-46 – +54	–
For corrosive media	GJ-N23.FFC.B	MI0284	40	3950	3.5	PTFE	PTFE	Hastelloy-C276	21	-46 – +54	–
	GJ-N25.FFC.B	MI0311	55	5480	3.5	PTFE	PTFE	Hastelloy-C276	21	-46 – +54	–
For abrasive media	GJ-N23.9FD.B	MI0265	40	3950	5.2	NIC	PTFE	Surface hardened	21	-46 – +54	–

- Ports (internal thread) 1/8" -27NPT
- Flow rates without differential pressure
- Operating temperature: with other seals up to 99°C possible



Delivery pump
BVP-Z Standard
ISM 446



Dosing pump
MCP-Z Standard
ISM 405



Programmed dosing
MCP-Z Process IP65
ISM 918

Tubing adaptors for gear pump-heads

Threaded stainless steel connectors

Order No.	External Tubing thread	Tubing adaptor	Tubing i.d. mm
AR0001	1/8" NPT	Tube nozzle	6
AR0002	1/8" NPT	Tube nozzle	3
AR0004	3/8" NPT	Tube nozzle	12
AR0008	1/8" NPT	Tube nozzle	8
AR0009	1/8" NPT	Tube nozzle	9.5
AR0024	1/8" NPT	Pipe connection	6 (outside)

Threaded connectors in Hastelloy-C

AR0001-HC	1/8" NPT	Tube nozzle	6
-----------	----------	-------------	---

Rotary Piston Pumps

For corrosive media and very accurate dispensing

The pump-heads are available with ceramic pistons and ceramic cylinder liners, which makes these components very resistant even to highly aggressive chemicals.



Rotation direction reversible



Safe and easy to use

- Developed for continuous duty, 24 hrs./day, 7 days/week
- Positive displacement
- No valves to clog
- **Precision better than $\pm 1\%$**
- High repeatability of dispensing volume
- Calibrateable drives
- MAX key enables rapid filling of the system (BVP-Z and MCP-Z pumps)

Inexpensive to maintain

- Interchangeable pump-heads
- No valves
- Only one moving part – the piston
- High quality and precision guarantee an optimum performance even after many years of intensive use

Application range of piston pumps

Industries	Applications	Special media
<ul style="list-style-type: none"> - Biotechnology - Chemistry - Industry - Electronic - Food and Dairy - Perfume/Cosmetics - Rubber/Plastics - Glass / Ceramic - Pulp and Paper - Medical 	<ul style="list-style-type: none"> - Accurate dispensing e.g. into bioreactors - Emulsion and slurry dosing - Medical diagnostics production - Milk and beverage enrichment - Plating bath replenishment - Titration equipment 	<ul style="list-style-type: none"> - Biozides - Dyes - Flux compound - Hydrogen peroxide - Liquid wax - Thixotropic products
		Not suited for media containing particles larger than 0.8 mm

For a wide range of applications

- Interchangeable pump-heads
- Adjustable stroke volume
- Very accurate dispensing pumps due to calibrateable drives
- Rotation direction reversible
- Chemically inert to a great extent
- Ideal for corrosive media, suited for viscous media
- Differential pressure up to 6.9 bar

Advantages of the valveless piston pump

- No valves which clog or hang up
- Only one moving part – the piston
- Drift free precision, better than $\pm 1\%$
- Variation coefficient smaller than 0.17%
- Medium-contacted parts available in ceramic and fluorocarbon
- Viscosity independent
- Positive displacement up to 6.9 bar
- Self-priming to 4.5 Meter

Typical Applications for Rotary Piston Pumps

Medical

For precise dispensing, aspirating, rinsing and mixing systems and for syringe pump replacement in diagnostic, clinical chemistry, dialysis and medical equipment manufacturing. Also for dispensing adhesives and lubricants used in assembly of disposable medical components.

Industrial


For accurate metering and mixing of paint and pigment additives, catalyst for foundry resins, plating bath regeneration, petroleum additives, photo chemicals, inks, monomers and adhesives.

Electronics Manufacturing

For dispensing of ceramic slurries in the manufacture of capacitors and diodes. Also for dispensing of insulating and encapsulating materials used in electric motor manufacture, addition of flux for wave soldering equipment, dispensing of mercury for switch manufacturing and metering of semiconductor wash and etch solutions.

Food and Dairy

For candy coating and polishing, vitamin fortification for milk; addition of flavors, colors and preservatives, hops for brewing and sanitizing agents for aseptic packaging. Also used for sample and reagent fluid control in milk analyzers and other food quality control instrumentation.



Valveless pumping

The valveless pumping function is accomplished by the synchronous rotation and reciprocation of the ceramic piston in the precisely mated ceramic cylinder liner. One complete piston revolution is required for each suction/discharge cycle.

The piston always bottoms for maximum fluid and bubble clearing. Together with the drive speed the stroke volume, which can be pre-set by the adjustment of the pump-head angle, determines the actual flow rate.

Note

All microprocessor controlled drives are LabVIEW compatible and can easily be integrated into process control systems.

Rotary piston pumps

Selection by flow rate and model

Flow rates ml/min		Bar	Model	Page
min.	max.			
0.045	180	6.9	REGLO-CPF <i>Analog</i> RH-type pump-heads	52 54
0.1	180	6.9	REGLO-CPF <i>Digital</i>	53
0.025	2300	6.9	MCP-CPF <i>Process</i> RH-type pump-heads Q-type pump-heads	56 54 58

Unique

Only the ISMATEC® Rotary Piston Pump

MCP-CPF *Process* features:



Carrying out programs independently of a PC

- Create the application profile in the PC (with ProgEdit software, Page 61)
- Download the file data into the pump memory
- Disconnect the pump from the PC
- Carry out your application on the spot, using the pump as a stand-alone unit



This mark indicates dispensing function

(You'll find dispensing rotary piston pumps on pages 53, 56)

- Pumping by speed or flow rate
- Dispensing by volume or time
- Interval dispensing with a pause
- Dispensing a volume within a pre-set time
- Interval dispensing with a pre-set number of dispensing cycles
- **Calibrating** the flow rate and dispensing volume
- Piston stroke back-steps **for drip-free dispensing**
- Factory-set piston pump-heads

REGLO-CPF

Calibrateable dispensing pumps
Ideal for dispensing corrosive media

- High repeatability
- Differential pressure up to 6.9 bar
- 10 cm wide, 13.5 cm high
- Wide selection of ceramic piston pumps



Rotation direction reversible



REGLO-CPF Analog
2-digit potentiometer
1–99%, resolution 1% (for speed)

3



REGLO-CPF Analog
with piston pump-head RH 00.CKC-LF

REGLO-CPF Analog
without dispensing functions
0.045–180 ml/min
Variable speed



Overview of piston pump-heads
on Pages 54 to 55

Specifications REGLO-CPF Analog

Motor type	DC-Motor
Speed	18 to 1800 rpm
Speed setting	1–99%, resolution 1% 2-digit potentiometer
Power consumption	50 W
Mains connection	230V _{AC} /50Hz, 115V _{AC} /60Hz adjustable
Protection rating	IP 30
Depth/Width/Height	250x100x143 mm
Weight	2.5 kg

Specifications REGLO-CPF Digital

Motor type	DC-Motor
Speed	40 to 1800 rpm
Speed setting	rpm, resolution 0.1rpm
Flow rate setting	µl/min and ml/min
Power consumption	75 W
Mains connection	100 – 230V _{AC} / 50 – 60Hz adjustable
Protection rating	IP 30
Depth/Width/Height	250x100x135 mm
Weight	2.1 kg



Rotation direction reversible



REGLO-CPF *Digital*
6-button membrane key-pad, LED display
Flow rate setting in µl/min and ml/min

Interfaces



REGLO-CPF *Analog*

- Speed control (0–5 or 0–10 V, 0–20 or 4–20 mA)
- Speed output 0–9 kHz
- Start/Stop
- Rotation direction



REGLO-CPF *Digital*

- RS232
- Analog
- Speed output 0–9 kHz, Start/Stop, Autostart



Dispensing and calibrating function see Page 51

REGLO-CPF *Digital*
with dispensing functions
0.1–180 ml/min
Microprocessor controlled



Application

Highly reproducible, single-channel dispensing processes of organic solvents or acids/bases, e.g.:
Dispensing of hydrogen fluoride and other highly corrosive acids with an X-Y-Z dispenser. Remote controlled pump.

Ordering information

Model (Drive only)	Order No. (drive only)	Flow rates	Channels Channels	Speed
REGLO-CPF <i>Analog</i>	ISM 1014	0.045 – 180	1	18 to 1800
REGLO-CPF <i>Digital</i>	ISM 321	0.1 – 180	1	40 to 1800
Foot switch				
REGLO-CPF <i>Analog</i>	ISM 891	see Page 61		
REGLO-CPF <i>Digital</i>	ISM 894	see Page 61		

The complete pump system REGLO-CPF consists of:

- 1 Drive
- 1 Piston pump-head see on Pages 54 to 55

RH pump-heads

For REGLO-CPF drives (Pages 52 to 53)
For MCP-CPF Process drive (Pages 56 to 57)



3

Pump-head RH 00
Stroke volumes 2.5 – 25 µl

Drives and flow rates:
REGLO-CPF Analog
0.045 – 45 ml/min
REGLO-CPF Digital
0.1 – 45 ml/min
MCP-CPF Process
0.025 – 45 ml/min



Type	RH00.CKC-LF	RH00.SKY-LF	RH00.STY-LF	RH00.CTC-LF
Order No.	FMI 009	FMI 010	FMI 011	FMI 012
Piston	Ceramic	316 Stainless Steel	316 SS	Ceramic
Cylinder case	Kynar® (Fluorocarbon {PVDF})	Kynar (Fluorocarbon {PVDF})	Tefzel®	Tefzel
Cylinder liner	Ceramic	Carbon	Carbon	Ceramic
Lip seals	Rulon® AR	Rulon J	Rulon J	Rulon AR
Gland washers	PTFE	PTFE	PTFE	PTFE
Max. temperature	100°C	60°C	60°C	100°C
Max. differential pressure	6.9 bar	6.9 bar	6.9 bar	6.9 bar
Flow ports	Kynar UNF 1/4"–28 (female)	Kynar UNF 1/4"–28 (female)	UNF 1/4"–28 (female)	UNF 1/4"–28 (female)

PTFE tubing for pump-heads mentioned above (must be ordered separately)

1.6 mm i.d., 3.2 mm o.d. with 2 fittings UNF 1/4"–28 male

Length	Order No.						
0.25 m	IC 0053	0.75 m	IC 0061	0.50 m	IC 0057	1.00 m	IC 0065



Tubing adaptors for the following pump-heads
RH00.CKC
RH00.SKY
RH0.CKC
RH1.CKC

These adaptors enable the use of other tubing. The integrally molded port fittings on the standard FMI Type K pump-heads accept all tubing with 6.4 mm o.d. For other tubing arrangements, these special port adaptors are required.

Description	Order No.
1 R412-0K for tubing with 3.2 mm i.d.	FMI 050
2 R412-1K for tubing with 6.4 mm i.d.	FMI 051
3 R412-2K for tubing with 9.5 mm i.d.	FMI 052
4 R412-5K for tubing with 1/4"–28 ferrule fittings	FMI 053
5 H476K for tubing with 3.2 mm o.d.	FMI 054

Pump-head RH 0
Stroke volumes 5 – 50 µl

Drives and flow rates:
REGLO-CPF *Analog*
0.09 – 90 ml/min
REGLO-CPF *Digital*
0.2 – 90 ml/min
MCP-CPF *Process*
0.050 – 90 ml/min



(LF = Low Flow
for flow rates below 50 ml/min)

Type and Order No.	RH0.CKC / FMI 005	RH0.CKC-LF / FMI 013	RH0.CTC / FMI 006
Piston	Ceramic	Ceramic	Ceramic
Cylinder case	Kynar® (Fluorocarbon {PVDF})	Kynar (Fluorocarbon {PVDF})	Tefzel®
Cylinder liner	Ceramic	Ceramic	Ceramic
Lip seals	Rulon® AR	Rulon AR	Rulon AR
Gland washers	PTFE	PTFE	PTFE
Max. temperature	100°C	100°C	100°C
Max. differential pressure	6.9 bar	6.9 bar	6.9 bar
Flow ports	2 fixed tube fittings for PTFE tubing 6 mm o.d.	Kynar UNF 1/4"–28 (female)	2 fixed tube fittings for PTFE tubing 6 mm o.d.
Tubing (must be ordered separately)	PTFE tubing 4 mm i.d., 6 mm o.d. Order No. MF 0336 (For other tubing material; use tubing adaptors, see Page 54)	PTFE tubing 1.6 mm i.d., 3.2 mm o.d. with 2 fittings UNF 1/4"–28 male Length Order No. 0.25 m IC 0053 0.50 m IC 0057 0.75 m IC 0061 1.00 m IC 0065	PTFE tubing 4 mm i.d., 6 mm o.d. Order No. MF 0336 (For other tubing material; use tubing adaptors, see Page 54) This pump-head is also available as LF version.

Pump-head RH 1
Stroke volumes 10 – 100 µl

Drives and flow rates:
REGLO-CPF *Analog*
0.18 – 180 ml/min
REGLO-CPF *Digital*
0.4 – 180 ml/min
MCP-CPF *Process*
0.1 – 180 ml/min



(LF = Low Flow
for flow rates below 50 ml/min)

Type and Order No.	RH1.CKC / FMI 007	RH1.CKC-LF / FMI 015	RH1.CTC / FMI 008
Piston	Ceramic	Ceramic	Ceramic
Cylinder case	Kynar (Fluorocarbon {PVDF})	Kynar (Fluorocarbon {PVDF})	Tefzel
Cylinder liner	Ceramic	Ceramic	Ceramic
Lip seals	Rulon AR	Rulon AR	Rulon AR
Gland washers	PTFE	PTFE	PTFE
Max. temperature	100°C	100°C	100°C
Max. differential pressure	6.9 bar	6.9 bar	6.9 bar
Flow ports	2 fixed tube fittings for PTFE tubing 6 mm o.d.	Kynar UNF 1/4"–28 (female)	2 fixed tube fittings for PTFE tubing 6 mm o.d.
Tubing (must be ordered separately)	PTFE tubing 4 mm i.d., 6 mm o.d. Order No. MF 0336 (For other tubing material, use tubing adaptors; see Page 54)	PTFE tubing 1.6 mm i.d., 3.2 mm o.d. with 2 fittings UNF 1/4"–28 male Length Order No. 0.25 m IC 0053 0.50 m IC 0057 0.75 m IC 0061 1.00 m IC 0065	PTFE tubing 4 mm i.d., 6 mm o.d. Order No. MF 0336 (For other tubing material; use tubing adaptors, see Page 54) This pump-head is also available as LF version.

MCP-CPF *Process*

Programmable

Programmable without a PC!

Protection rating of IP 65

3

- Ideal for aggressive media
- High repeatability
- Differential pressure 6.9 bar
- Ideal for dispensing and filling applications in a dusty, humid or corrosive environment and in clean room areas (IP 65, dust-tight and protected against water jets)



Dispensing and calibrating function see Page 51



Rotation direction reversible

MCP-CPF *Process*

- **Pre-programmed pump-heads allow you to work with flow rates**
 - Stainless steel housing, membrane key-pad, LED display
 - **4 program memories for saving individual application parameters or PC programmed command sequences**
 - **Programming similar to PLC**
 - Wide selection of different, interchangeable pump-heads
- Flow rates and differential pressure depend on the pump-head mounted, see Pages 54, 55 and 57–59

MCP-CPF *Process*
with rotary piston pump-head QP Q0.SSY-LF

Specifications


Motor type	DC motor
Speed	10.0 to 1800 rpm
Speed setting	rpm, resolution 0.1 rpm
Flow rate setting	µl/min, ml/min, liters/min
Power consumption	100 W
Mains connection	100 – 230 V _{AC} / 50 – 60 Hz
Protection rating	IP 65
Depth/Width/Height	220 x 155 x 260 mm (without pump-head)
Weight	6.9 kg (without pump-head)

Ordering information

The complete pump system MCP-CPF *Process* consists of:

Drive	ISM 919
Pump-head and tubing	see Pages 54, 55 and 57–59 / 62–63
Accessories	
– Software ProgEdit (Page 61)	SOF 104
– Foot switch (Page 61)	IS 10039

LabVIEW driver
download for free: www.ismatec.com




Software ProgEdit
LabVIEW drivers
Free download on
www.ismatec.com



'RH' pump-heads (description see Pages 54 to 55)

Type	Flow rates ml/min	Stroke volumes µl
RH 00	0.025 – 45	2.5 – 25
RH 0	0.050 – 90	5.0 – 50
RH 1	0.10 – 180	10.0 – 100



'Q' pump-heads (description see Pages 58 to 59)

Type	Flow rates ml/min	Stroke volumes µl
QP Q0	0.04 – 144	3.2 – 80
QP Q1	0.13 – 576	12.8 – 320
QP Q2	0.29 – 1300	28.8 – 720
QP Q3	0.51 – 2300	51.2 – 1280

Interfaces



PC-controllable:
– RS232



- Speed control (0–5 or 0–10 V, 0–20 or 4–20 mA)
- Speed output (0–10 V_{DC} or 0–7.2 kHz)
- Start/Stop
- Rotation direction
- Autostart
- 2 universal inputs
- 2 universal outputs



Application

- Single-channel sterile delivery and dispensing processes under pressure for particulate-free solvents
- Addition of various reagents in different volume ratios through mixing valve into reactor

Q-type pump-heads

For MCP-CPF Process drive (Page 56)



MCP-CPF Process with Q pump-heads and Low Flow Kit R479

3

Pump-heads Q0 and Q3

Q0 = Stroke vol. 3.2 – 80 µl
Q3 = Stroke vol. 51.2 – 1280 µl

Q0 = Flow rate 0.04 – 144 ml/min
Q3 = Flow rate 0.51 – 2300 ml/min



Type and Order No.	QP Q0.SS / FMI 202	QP Q0.SKY / FMI 316	QP Q3.CKC / FMI 217
Piston	316 Stainless Steel	316 Stainless Steel	Ceramic
Cylinder case	316 Stainless Steel	Kynar® (Fluorocarbon {PVDF})	Kynar (Fluorocarbon {PVDF})
Cylinder liner	Carbon	Carbon	Ceramic
Lip seals	Rulon® J	Rulon J	Rulon AR
Gland washers	PTFE	PTFE	PTFE
Cylinder head seal	PTFE	none	none
Max. temperature	60°C	60°C	100°C
Max. differential pressure	6.9 bar	4.1 bar	1.7 bar (to 1600rpm) 0.5 bar (from 1600 rpm)
Flow ports	1/4 NPT (female) <u>Included:</u> 2 stainless steel adaptors with thread 1/4 NPT (male) and fitting for tubing with 6.4 mm i.d.	for tubing up to 12.7 mm i.d. <u>Included:</u> 2 Kynar (PVDF) adaptors for tubing with 6 mm o.d.	for tubing up to 12.7 mm i.d. or PTFE tubing 6 mm o.d. <u>Included:</u> 2 Kynar (PVDF) adaptors for tubing with 6 mm o.d.
Tubing (must be ordered separately)	Tubing Tygon® ST R-3603 6.4 mm i.d. Order No. MF 0031 Accessories Low Flow Kit R 479 (see below) Order No. FMI 056	Tubing Tygon ST R-3603 12.7 mm i.d. Order No. SC 0382 PTFE Tubing 4 mm i.d. / 6 mm o.d., 3.6 m long Order No. MF 0336	Tubing Tygon ST R-3603 12.7 mm i.d. Order No. SC 0382 PTFE Tubing 4 mm i.d. / 6 mm o.d., 3.6 m long Order No. MF 0336



Low Flow Kit R 479 Order No. FMI 056
suitable for the following pump-heads:
QP Q0.SSY QP Q1.SSY QP Q2.CSY
QP Q1.CSC QP Q2.CSC QP Q2.SSY
QP Q1.CSY

This Low Flow adaptor Kit enables the use of the above mentioned pump-heads for flow rates below 50 ml/min or in case that a minimum dead volume or a maximum of chemical compatibility are required. The adaptor features a 1/4-28 inner thread. These threads are used with low flow tube fittings for small bore tubing of 3.2 mm o.d. or less. Hence, this »Low Flow Kit« is also very interesting for chromatography applications.

PTFE tubing for Low Flow Kit R 479

1.6 mm i.d. / 3.2 mm o.d., with 2 fittings 1/4-28 (male)
0.25 m long Order No. IC 0053 0.75 m long Order No. IC 0061
0.50 m long Order No. IC 0057 1.00 m long Order No. IC 0065



Tubing adaptors for pump-heads with a Kynar cylinder case:
Q0.SKY Q2.CKC
Q1.CKC Q2.CKY
Q1.CKY Q2.SKY
Q1.SKY Q3.CKC

In addition to the tubing mentioned above, these adaptors enable the use of other tubing.

Description	Order No.
1 R412-0K for tubing with 3.2 mm i.d.	FMI 050
2 R412-1K for tubing with 6.4 mm i.d.	FMI 051
3 R412-2K for tubing with 9.5 mm i.d.	FMI 052
4 R412-5K for tubing with 1/4-28 ferrule fittings	FMI 053
5 H476K for tubing with 3.2 mm o.d.	FMI 054

Pump-heads Q1 and Q2

Q1 = stroke vol. 12.8 – 320 µl

Q2 = stroke vol. 28.8 – 720 µl

Q1 = flow rates 0.13 – 576 ml/min

Q2 = flow rates 0.29 – 1300 ml/min



Type Order No.	QP Q1.CSC FMI 205	QP Q2.CSC FMI 212	QP Q1.CSC-W FMI 320	QP Q2.CSC-W FMI 321	QP Q1.CSC-WT FMI 219	QP Q2.CSC-WT FMI 218
Piston	Ceramic		Material and design like QP Q1.CSC and CP Q2.CSC but with isolation gland (2 extra ports 10–32 – female)		Material and design like QP Q1.CSC and CP Q2.CSC but with isolation gland (2 extra ports 1/8" NPT – female) and heating mantel	
Cylinder case	316 Stainless Steel		Thanks to a barrier gland of fluid, gas, steam or whatever is needed, the pumped fluid can be isolated from the seal area and atmosphere. Slurries, particulates, crystal formers and anaerobics are easily handled.		Same barrier gland as described under CP Q1./Q2.CSC-W	
Cylinder liner	Ceramic		Included for barrier gland ports: 2 Polypropylene adaptors, thread 10–32 UNF and fitting for tubing with 3.2 mm i.d.		In addition, 2 cartridge heaters (1/4" diam. x 1 1/2" long) and 1 thermo-couple (1/8" diam. x 1" long) can be used for heating the pump-head.	
Lip seals	Rulon® AR				Not included: Tubing adaptors for: – barrier gland ports 1/8" NPT (female) – main flow ports 1/4" NPT (female)	
Gland washers	PTFE					
Cylinder head seal	PTFE					
Max. temperature	177°C					
Max. differential pressure	6.9 bar					
Main flow ports	1/4" NPT (female) Included: 2 stainless steel adaptors with thread 1/4" NPT (male) and fitting for tubing with 9.5 mm i.d.					
Tubing (must be ordered separately)	Tubing Tygon® ST R-3603 9.5 mm i.d. Order No. SC 0383 Accessories Low Flow Kit R 479 Order No. FMI 056					

Other materials for wetted parts for:

Pump-heads Q1 and Q2 (see table below)

Q1 = stroke vol. 12.8 – 320 µl

Q2 = stroke vol. 28.8 – 720 µl

Q1 = flow rates 0.13 – 576 ml/min

Q2 = flow rates 0.29 – 1300 ml/min

Tubing and connections for pump-heads with the suffix -W or -WT (must be ordered separately)	
Description	Order No.
–2 stainless steel fittings for inlet/outlet, thread 1/4" NPT male, with fittings for tubing with 6.4 mm i.d.	FMI 060
–Tubing for inlet/outlet (Tygon ST R-3603) 6.4 mm i.d., 15 m long	MF 0031

Type / Order No.	QP Q1.CKC/FMI 352	QP Q1.CKC-W/FMI 356	QP Q1.CKY/FMI 358	QP Q1.CSY/FMI 359	QP Q1.SKY/FMI 361	QP Q1.SSY/FMI 363	QP Q1.SAN'/FMI 365
Type/ Order No.	QP Q2.CKC/FMI 355	QP Q2.CKC-W/FMI 357	QP Q2.CKY/FMI 353	QP Q2.CSY/FMI 360	QP Q2.SKY/FMI 362	QP Q2.SSY/FMI 364	QP Q2.SAN'/FMI 366
Piston	Ceramic	Ceramic	Ceramic	Ceramic	316 Stainless Steel	316 Stainless Steel	Ceramic
Cylinder Case	Kynar® ²	Kynar ²	Kynar ²	316 Stainless Steel	Kynar ²	316 Stainless Steel	316 Stainless Steel
Cylinder liner	Ceramic	Ceramic	Carbon	Carbon	Carbon	Carbon	316 Stainless Steel
Lip seals	Rulon AR	Rulon AR	Rulon AR	Rulon AR	Rulon J	Rulon J	PTFE
Gland washers	PTFE	PTFE	PTFE	PTFE	PTFE	PTFE	PTFE
Cylinder head seal	none	none	none	PTFE	none	PTFE	PTFE
Max. temperature	100°C	100°C	100°C	177°C	60°C	60°C	177°C
Max. diff. pressure	4.1 bar	4.1 bar	4.1 bar	6.9 bar	4.1 bar	6.9 bar	6.9 bar
Main flow ports	For tubing up to 9.5 mm i.d.	For tubing up to 9.5 mm i.d. With isolation gland Fittings for tubing with 3.2 mm i.d.	For tubing up to 9.5 mm i.d.	1/4" NPT (female)	For tubing up to 9.5 mm i.d.	1/4" NPT (female)	PTFE tubing adaptor

¹ designed for sanitary applications

² Kynar = Fluorocarbon (PVDF)



Tubing cassettes

Developed and consistently improved by ISMATEC

Click'n'go cassette with automatic pressure mechanism



Click 'n' go Cassettes (The new standard with all ISMATEC® cassette pumps)

The new, innovative tubing cassette from ISMATEC, a consequence of the constant improvement of the cassettes with pressure lever, offers the following advantages:

- Automatic tubing pressure
- Calibrated, fatigue-free spring guarantees optimal, reproducible tubing pressure independent of diameter, material and state of the tubing
- Pumping conditions are defined and repeatable at a later date
- Incorrect adjustment of the pressure lever is no longer possible
- **Readjustment** of tubing pressure is **not necessary**
- Long-time channel to channel conformity
- **Ideal for non-monitored long-time use**

Click 'n' go cassettes are not suitable for:

- Hard tubing materials
- Differential pressure greater than 1 bar

For these conditions you should choose the pressure lever cassettes.

Ordering information for spare Click'n'go cassettes (included with the corresponding pump model).

Model	Material	Order No.
MS/CA Click'n'go	POM-C*	IS 3510
CA ¹ Click'n'go	POM-C*	IS 3710



¹ Replacement adaptors for CA cassettes

Adaptor for CA cassettes	POM-C	IS 0123
Adaptor for CA cassettes	PVDF	IS 3861

Please order 2 adaptors each.

Option
Cassettes with adjustable pressure lever



ISMATEC pumps equipped with Click'n'go cassettes

Pressure Lever Cassettes (Optional)

The successful, tried and tested tubing cassettes for many years. The pressure lever allows you to set a different tubing pressure for each channel. Depending on the application, tubing material and diameter, an optimally adjusted tubing pressure can be set, i.e. just as high as necessary (can be favourable to the tubing life). For constant flow rates it may be necessary to periodically adjust the tubing pressure.

The pressure lever should not be used to adjust the flow rate! Such a 'calibration' can be unstable and is detrimental for tubing life.

These cassettes are available as an option. Under certain circumstances we recommend the use of these cassettes rather than the automatic click 'n' go cassettes:

- Under particular system conditions e.g. varying or high differential pressure
- When using very hard tubing material

Ordering information for cassettes with pressure lever
Must be ordered separately.

Model	Material	Order No.
MS/CA pressure lever	POM-C*	IS 0649
MS/CA pressure lever	PVDF*	IS 3629
CA ¹ pressure lever	POM-C*	IS 0122
CA ¹ pressure lever	PVDF*	IS 3820

PVDF = higher chemical resistance

*Material properties

POM-C	Polyoxymethylen Copolymer
PVDF	Polyvinylidene Fluoride

Foot switches, software



Foot switch for start/stop

Very practical for use with pumps as dispensing systems e.g. for filling tubes, bottles etc. Footswitch provides the start/stop signal required.

You have both hands free for handling the filling system.

Protection rating

IP21 (higher safety categories e.g. IP67 on request)

Ordering information

Foot switch suitable for pump models:	Order No.
IPC and IPC-N (firmware version older than 4.00)	ISM 016
IP and IP-N, IPC and IPC-N (from firmware version 4.00)	IS 10039
REGLO Analog, REGLO Quick™	
REGLO-Z , REGLO-ZS, REGLO-CPF Analog	ISM 891
REGLO Digital, REGLO-CPF Digital, REGLO-Z Digital	ISM 894
ecoline Series	IS 3572
MCP Standard (up to firmware version 7.00)	
MCP-Z (firmware version older than 4.00)	ISM 016
MCP Standard (from firmware version 7.00)	
MCP Process, MCP-CPF Process, MCP-Z Process, and MCP-Z Standard (from firmware version 4.00)	IS 10039
BVP Process	IS 10039
BVP Standard, BVP-Z	ISM 891
Flowmaster®	IS 10279



Get your personal ProgEdit Software for free from www.ismatec.com

Create your application program in the PC and then work with your pump as a stand-alone system (for MCP-Process, MCP-Z Process, MCP-CPF Process)

With the ProgEdit Software you can individually define and program your dosing and pumping sequences.

The application programs written in the PC can be downloaded into the program memory of the pump drive. Up to 4 individual application programs can be saved. The pump can use PC downloaded programs as a stand-alone system independent of a PC, e.g.:

- Filling and dispensing (e.g. dispensing of a perfume proportionally to the speed of a conveyor belt).
- Sampling (e.g. volumetrically proportional sampling from a waste water channel).
- Pressure control (e.g. for safe filtering)
- Fully automatic 1- and 2-point niveau regulation (e.g. pump keeps liquid level constant between given boundaries)



Software English/German switchable

- Freely defined, time dependent rotation rate profile (e.g. ramps)
- Monitoring inputs (2) and setting outputs (2) for applications similar to SPS
- 'for/while' loops and 'if-then-else' commands
- Compare constants and variables with sensor signals

Download from our web-site (www.ismatec.com):

- Free ProgEdit software
- Operating manual (PDF-file) (1.0 MB)
- LabView Software Driver

IC-Tubing and fittings

Metered IC-tubing

The IC-tubing is made of high density PTFE. It is chemically inert and is available as a pack of 3 m length.

You can make the appropriate tubing connections for your system by flanging the IC-tubing and using combined with IC fittings.

For an alternative to flanging tubing, we highly recommend the Flangeless fittings on Page 63.

Ordering information

Dimensions			Order No.
i.d.	o.d.	Length	
0.3	1.5 mm (1/16")	3 m	IC 0046
0.5	1.5 mm (1/16")	3 m	IC 0047
0.8	1.5 mm (1/16")	3 m	IC 0048
1.5	3.0 mm (1/8")	3 m	IC 0049

Flanged IC-tubing with fittings

The pre-flanged IC-tubing is made of the same high density PTFE tubing and is available in 4 different lengths and diameters. Comes complete with 2 fittings.



IC-tubing with 2 fittings (UNF 1/4-28)

Ordering information

Tubing i.d./o.d. mm	Length 0.25 m Order No.	Length 0.50 m Order No.	Length 0.75 m Order No.	Length 1.00 m Order No.
0.3/1.5 (1/16")	IC 0050	IC 0054	IC 0058	IC 0062
0.5/1.5 (1/16")	IC 0051	IC 0055	IC 0059	IC 0063
0.8/1.5 (1/16")	IC 0052	IC 0056	IC 0060	IC 0064
1.5/3.0 (1/16")	IC 0053	IC 0057	IC 0061	IC 0065



1 IC-fitting for metered tubing

Each fitting comes with a washer.

Material	Polyacetate
Color	white, blue, red, green, yellow, black
Thread	UNF 1/4-28 for both tube sizes

IC-fitting for tubing Color	1.5 mm (1/16") o.d. Order No.	3 mm (1/8") o.d. Order No.
White	IC 0001	IC 0007
Blue	IC 0002	IC 0008
Red	IC 0003	IC 0009
Green	IC 0004	IC 0010
Yellow	IC 0005	IC 0011
Black	IC 0006	IC 0012
Spare U-washers	WA 0002	WA 0001

2 Coupling-sleeve

To join 2 fittings
 Thread UNF 1/4-28
 Material Delrin
 Order No. IC 0013

3 End-plug

To close a PTFE tubing (used with coupling-sleeve)
 Material PTFE
 Order No. IC 0021

Pack size: 5 pieces

Flangeless fittings



Upchurch Scientific® Flangeless fittings eliminate the need to flange tubing. This removable and reusable system provides several benefits:

Convenience: Flangeless fittings are easy to replace. Just slip the nut and ferrule over the tubing and finger tighten the nut.

Minimal Down-Time: Component replacement is quick, taking only a few seconds — unlike the significant time required to flange tubing.

Cost-Effectiveness: Repairing a flanged tubing assembly requires a costly flanging tool or the purchase of a complete replacement assembly, including a new length of tubing and a set of fittings. The Flangeless fittings system typically requires only one new ferrule at minimal cost.

The Flangeless fittings are available for 1.5 mm (1/16") and 3.0 mm (1/8") o.d. tubing. Nuts and ferrules are sold separately. Use these fittings with the metered IC-tubing on page 62.




4

The convenience of Flangeless fittings

The Flangeless nut provides fingertight convenience – no wrenches required.

The Flangeless ferrule provides a leak-proof seal. There is no need to spend time flanging tubing.



Material

Nut: Delrin® (acetal resin)
 Ferrule: Tefzel® (ETFE)

Color

Nut: black, red, white, cream, green, blue
 Ferrule: blue 1.5 mm (1/16")
 yellow 3.0 mm (1/8")

Thread

UNF 1/4-28 for both tubing sizes

Ordering information

Flangeless fittings	1.5 mm (1/16") o.d.	3 mm (1/8") o.d.
Color	Order No.	Order No.
Nuts		
Black	P-201x	P-301x
Red	P-202x	P-302x
White	P-203x	P-303x
Cream	P-204x	P-304x
Green	P-205x	P-305x
Blue	P-206x	P-306x
Ferrules		
Blue	P-200x	—
Yellow	—	P-300x

Pack size: 10 pieces


OEM tubing pumps

from 1 μ l/min up to 7240 ml/min


ISMATEC®




1 to 5 channels
Synchro-, DC- or Steppermotor



6 rollers
0.002–68.8 ml/min



8 rollers
0.002–56 ml/min



12 rollers
0.001–38.4 ml/min
low pulsation

IDEX® Health & Science Complete Fluidic Solutions

- Custom fluidic system and component development
- Integration of fluid handling technologies
- Micro and nanofluidic technologies
- Ultra-hard materials processing for durability and compatability
- Process management for sterility and traceability compliance
- Patented, proprietary processes and long-term strategic partnerships with materials suppliers

OEM gear pumps

up to 5005 ml/min

MICROPUMP®



Small gear pump
DC motor with optical encoder



Micropump I-Drive electromagnetic drive pump for Series GA, GJ PPS head
Brushless DC motor

Applications

Integrated into systems for uses such as:

- UV/Vis spectrophotometric on-line analysis of water from a sewage treatment plant
- Continuous photometric analysis of process water
- Liquid analysis on-line from various processes
- Automated sample preparation of gas-form hydrides for IPC analysis
- Fully automated analysis of process sample streams by fluorescence spectroscopy
- Automated sampler for sample dialysis and ion-chromatography
- Continuous BOD, COD, toxicity, TOC, P and N analyzers

Visit our home page:

www.idex-hs.com

or ask for the separate OEM pump brochure

A-Z

2-Stopper tubing 37
 3-Stopper tubing 37
 360 tubing pump-head 22
 380 tubing pump-head 22
 380AD tubing pump-head 22

A

Adaptor for CA cassettes 60
 Adhesive for tubing sleeves 38

B

BVP Process tubing pump drive 18
 BVP Standard tubing pump drive 16
 BVP-Z Standard gear pump drive 46
 Bypass for gear pumps 43

C

CA Click'n'go 60
 CA-12 tubing pump-head 24
 CA-4 tubing pump-head 24
 CA-8 tubing pump-head 24
 Cassettes 60
 Cassettes with pressure level 60
 Cavity Style gear pump 42
 Click'n'go cassettes 60
 Connector for tubing 38, 62, 63
 Coupling-sleeve 62

D

Drive principle gear pump 43

E

Easy-Load® II tubing pump-head 23
 Easy-Load tubing pump-head 23
 Ecoline tubing pump 12
 End-plug 62
 Extension tubing 36

F

Flangeless fittings 63
 Fittings 62, 63
 Flowmaster® tubing pump 28
 Fluran® tubing 33
 Foot switch 61

G

Gentle cell delivery 21
 GORE™ tubing 34

H

HC Tygon® tubing 33

I

ic-fittings for tubing 62
 ic-tubing 62
 IP multi-channel tubing pumps 6
 IPC multi-channel tubing pumps 6
 IPC-N multi-channel tubing pumps 6
 IP-N multi-channel tubing pumps 6
 Ismaprene tubing 32

L

LFL Tygon tubing 32

M

MCP Process tubing pump drive 19
 MCP Standard tubing pump drive 17
 MCP-CPF Process piston pump drive 56
 MCP-Z Process gear pump drive 48
 MCP-Z Standard gear pump drive 47

M

MHLL Tygon tubing 32, 35
 MHSL Tygon tubing 32, 35
 MS/CA 4-12 tubing pump-head 25
 MS/CA 8-6 tubing pump-head 25
 MS/CA Click'n'go 60
 MS/CA panel mounted pump 26

N

Norprene® tubing 33
 Nozzles 45, 49

O

OEM Pumps 64
 Operating pressure, tubing 30
 Ordering information, tubing 36
 Overview, gear pumps 40
 Overview, MCP/BVP tubing pump system 14
 Overview, piston pumps 50
 Overview, rotary piston pumps 50
 Overview, tubing pumps 4

P

Panel mounted pumps 26
 Peroxide silicone tubing 33
 PharMed® tubing 32
 Planetary drive 6, 7
 Platinum silicon tubing 33
 Pressure lever for cassettes 60
 Pro-280 tubing pump-head 20, 21
 Pro-281 tubing pump-head 20
 Pro-380 tubing pump-head 20
 Pro-381 tubing pump-head 20
 ProgEdit software 61
 Properties of gear pump-heads 43
 Properties of Ismatec® tubing 30, 32
 PTFE-Tube tubing pump-head 23

Q

Q-piston pump-heads 58

R

R-3603/R-3607 Tygon-tubing 32
 Rating comparison, tubing 30
 REGLO Analog multi-channel tubing pump 8
 REGLO Digital multi-channel tubing pump 8
 REGLO Quick™ tubing pump 10
 REGLO-CPF Analog piston pump 52
 REGLO-CPF Digital piston pump 52
 REGLO-Z Analog gear pump 44
 REGLO-Z Digital gear pump 44
 REGLO-ZS Analog gear pump 44
 REGLO-ZS Digital gear pump 44
 RH-piston pump-heads 54

S

SB 2V tubing pump-head 24
 SB 3V tubing pump-head 24
 Selection criteria gear pump 42
 Selection criteria pumps 3
 Selection gear pumps 41
 Service life tubing 32
 SI Tygon tubing 33
 Silicone 3350 (Platinum) Tygon tubing 33
 Silicone spray 36
 Silicone tubing (Peroxide) 33
 Software 61
 ST Tygon tubing 32
 Standard tubing 36
 Sterilization of tubing 32
 Suction shoe style gear pumps 42

T

Tubing connector 38, 62, 63
 Tubing for aggressive media 32, 34, 35
 Tubing for special application 32, 34, 35
 Tubing maximum operating pressure 34
 Tubing order information 36
 Tygon tubing 32, 33, 35

V

Viton® tubing 33

The following trademarks and registered trademarks are used in this publication and are owned by the companies as indicated:

Owned by IDEX® Corporation or Units of IDEX Corporation

Trademark	Organization Owning the Trademark
Flowmaster®	Ismatec SA
IDEX®	IDEX Corporation
ISMATEC®	Ismatec SA
ISMATEC logo (®)	Ismatec SA
Micropump®	Micropump, Inc.
Micropump logo (®)	Micropump, Inc.
Precision Engineered Fluidics™	IDEX Health & Science
REGLO Quick™	Ismatec SA
Upchurch Scientific®	Upchurch Scientific, Inc.
Upchurch Scientific logo (®)	Upchurch Scientific, Inc.

Owned by other companies

Trademark	Organization Owning the Trademark
Delrin®	E.I. du Pont de Nemours and Company
CHEM-SURE™	W.L. Gore & Associates, Inc.
Easy-Load®	Cole-Parmer Instrument Company
ergo®	Kisling AG
Fluran®	Saint-Gobain Performance Plastics
Gore™	W.L. Gore & Associates, Inc.
Hastelloy®	Haynes International, Inc.
Kalrez®	DuPont Performance Elastomers
Kynar®	Arkema, Inc.
Norprene®	Saint-Gobain Performance Plastics
PEEK™	Victrex plc
PharMed®	Saint-Gobain Performance Plastics
Rulon®	Saint-Gobain Performance Plastics
Tefzel®	E.I. du Pont de Nemours and Company
Tygon®	Saint-Gobain Performance Plastics
STA-PURE®	W.L. Gore & Associates, Inc.
Viton®	E.I. du Pont de Nemours and Company

INDEX OF ORDER NUMBERS
Excluding Tubing and Connectors (36–39)

Order No.	Page	Order No.	Page	Order No.	Page	Order No.	Page	Order No.	Page
AR0001	45, 49	IC0051	62	ISM719	22	ISM937	6	P-205x	63
AR0001-HC	45, 49	IC0052	62	ISM721	24	ISM938	6	P-206x	63
AR0002	45, 49	IC0053	54, 55, 58, 62	ISM724	25	ISM939	6	P-300x	63
AR0004	45, 49	IC0054	62	ISM725	22	ISM940	6	P-301x	63
AR0008	45, 49	IC0055	62	ISM727	23	ISM941	6	P-302x	63
AR0009	45, 49	IC0056	62	ISM732	24	ISM942	6	P-303x	63
AR0024	45, 49	IC0057	54, 55, 58, 62	ISM733	24	ISM943	6	P-304x	63
FMI005	55	IC0058	62	ISM734	24	ISM944	6	P-305x	63
FMI006	55	IC0059	62	ISM735	25	ISM945	6	P-306x	63
FMI007	55	IC0060	62	ISM737	25	ISM946	6	SC0179	36
FMI008	55	IC0061	54, 55, 58, 62	ISM738	23	ISM947	6	SC0382	58
FMI009	54	IC0062	62	ISM785	20	ISM948	6	SC0383	59
FMI010	54	IC0063	62	ISM791	20	ISM949	6	SC1010	38
FMI011	54	IC0064	62	ISM793	20	MF0031	58, 59	SC1011	38
FMI012	54	IC0065	54, 55, 58, 62	ISM795	9	MF0313	23	SC1012	38
FMI013	55	IS0122	60	ISM796	9	MF0330	23	WA0001	62
FMI015	55	IS0123	60	ISM797	20	MF0331	23	WA0002	62
FMI050	54, 58	IS0649	60	ISM827	9	MF0332	23		
FMI051	54, 58	IS10039	61	ISM828	9	MF0333	23		
FMI052	54, 58	IS10279	61	ISM829	9	MF0334	23		
FMI053	54, 58	IS1080	13	ISM830	9	MF0335	23		
FMI054	54, 58	IS1085	13	ISM831	9	MF0336	23, 55, 58		
FMI056	58, 59	IS3510	60	ISM832	9	MF0337	23		
FMI060	59	IS3572	61	ISM833	9	MF0446	23		
FMI202	58	IS3629	60	ISM834	9	MI0006	45, 49		
FMI205	59	IS3710	60	ISM844	27	MI0007	45, 49		
FMI212	59	IS3762	13	ISM845	27	MI0008	45, 49		
FMI217	58	IS3763	13	ISM846	27	MI0013	45, 49		
FMI218	59	IS3820	60	ISM847	27	MI0016	45, 49		
FMI219	59	IS3861	60	ISM848	27	MI0018	45, 49		
FMI316	58	ISM010	24	ISM849	27	MI0019	45, 49		
FMI320	59	ISM011	24	ISM850	27	MI0020	45, 49		
FMI321	59	ISM016	61	ISM851	27	MI0022	49		
IC0001	62	ISM1014	53	ISM852	27	MI0023	49		
IC0002	62	ISM1020	29	ISM853	27	MI0131	45, 49		
IC0003	62	ISM1022	29	ISM891	61	MI0132	45		
IC0004	62	ISM1076	13	ISM894	61	MI0265	45, 49		
IC0005	62	ISM1078	13	ISM895	44	MI0280	45, 49		
IC0006	62	ISM1079	13	ISM896	44	MI0284	45, 49		
IC0007	62	ISM1089	13	ISM897	11	MI0306	45, 49		
IC0008	62	ISM1090	13	ISM901	44	MI0308	45		
IC0009	62	ISM1091	13	ISM915	19	MI0309	45, 49		
IC0010	62	ISM1143	44	ISM918	48	MI0310	45, 49		
IC0011	62	ISM185	25	ISM919	56	MI0311	49		
IC0012	62	ISM321	53	ISM920	18	MI0312	45, 49		
IC0013	62	ISM404	17	ISM930	6	MI0313	45, 49		
IC0021	62	ISM405	47	ISM931	6	MI0378	49		
IC0046	62	ISM444	16	ISM932	6	P-200x	63		
IC0047	62	ISM446	46	ISM933	6	P-201x	63		
IC0048	62	ISM596	9	ISM934	6	P-202x	63		
IC0049	62	ISM597	9	ISM935	6	P-203x	63		
IC0050	62	ISM718	22	ISM936	6	P-204x	63		

Your ISMATEC® representative

Revised 12/10/08, Number Printed 7,500

sales.ismatec@idexcorp.com
www.ismatec.com
www.idex-hs.com/ismatec

United States

ISMATEC - Micropump Headquarters
IDEX Health & Science
1402 NE 136th Ave.
Vancouver, WA 98684
Tel. +1 360 253 2008
Fax +1 360 253 8294

Germany

ISMATEC European Headquarters
IDEX Health & Science
ISMATEC Laboratoriumstechnik GmbH
Futtererstrasse 16
D-97877 Wertheim-Mondfeld
Tel. +49 93 77 / 92 03 0
Fax +49 93 77 / 13 88
office.ismatec@idexcorp.com
www.ismatec.de

Switzerland

ISMATEC Labortechnik SA
IDEX Health & Science
Feldeggstrasse 6
CH-8152 Glattbrugg-Zürich
Tel. +41 44 874 94 94
Fax +41 44 810 52 92
www.ismatec.ch

United Kingdom

Micropump Europe
IDEX Health & Science
Howard Road, Eaton Socon
St. Neots, Cambridgeshire
England
Tel. +44 1480 356600
Fax +44 1480 356920

Asia

IDEX Health & Science
Room 1701A, CBD International Plaza
No.16 Young An Dong Li
Chao Yang District
Beijing, China 100002
Tel. +86 10 6566 9090
Fax +86 10 6567 5152



ISMATEC®

IDEX
HEALTH & SCIENCE

©2009 IDEX Health & Science