General purpose stirred thermostatic baths/circulators

GD100, GD120

Operating Manual

EN

FR

DE

IT

ES

Contents

1	Safety.....	2
2	Assembly.....	4
2.1	Unpacking.....	4
2.2	Assembling control unit for use as bath controller or immersion thermostat.....	4
2.3	Use at lower temperatures.....	8
3	Operation.....	9
3.1	Controls and indicator lamps.....	9
3.2	Filling.....	12
3.3	Switching on and setting up.....	12
3.4	GD100 and GD120 Preset temperatures.....	13
3.5	GD120 Alarm (ALAR), Clock (CLOC) and Liquid selection (H2O, OIL) and adjustment Fig 7.....	13
3.6	Calibration.....	15
		17
4	Fault diagnosis.....	19
5	Technical specification.....	20
5.1	Environmental.....	20
5.2	Heating.....	20
6	Maintenance.....	22
6.1	Routine maintenance.....	22
6.2	Cleaning.....	22
6.3	Replacement of Fuses.....	23
7	Guarantee & service.....	24
7.1	Guarantee.....	24
7.2	Service.....	24
8	Compliance	25
8.1	Disposal and WEEE.....	25
8.2	RoHS.....	25
8.3	Electrical Safety and Electromagnetic Compatibility.....	25
9	Declaration of Conformity.....	Inside rear cover

1. Safety

The following symbols mean:-

Caution: Read these operating instructions fully before use and pay particular attention to sections containing this symbol

Caution: Surfaces can become hot during use.

Always observe the following safety precautions

Use only as specified by the operating instructions, or the intrinsic protection may be impaired.

After transport or storage in humid conditions, dry out the unit before connecting it to the supply voltage.

Connect only to a power supply with a voltage corresponding to that on the serial number label.

Ensure that the mains switch and isolating device (power supply connector) are easily accessible during use.

Connect only to a power supply which provides a safety earth (ground) terminal.

Before moving, disconnect at the power supply socket.

Do not check the temperature by touch, use the temperature display (where applicable) or a thermometer.

To reduce the risk of eye injury during high temperature operation, use safety goggles or spectacles.

Do not touch surfaces which become hot during high temperature operation.

Ensure that the operating temperature is less than the maximum operating temperature of your sample material.

Do not block or restrict ventilation slots.

If liquid is spilt inside the unit, disconnect it from the power supply and have it checked by a competent person.

It is the user's responsibility to carry out appropriate decontamination if hazardous material is spilt on or inside the equipment.

EN

-  Do not connect to a power supply or switch on before filling the tank.
-  Take care when topping up or draining, as the liquid in the tank may be very hot or cold.
-  If the alarm lamp is illuminated do not touch the liquid or the tank base, they may be very hot. Refill carefully, a hot heater can cause a spattering of very hot water droplets and scalding steam.
-  Always use a lid or polypropylene spheres when operating above 60°C. Take care when raising and removing the lid, it may be hot. Steam and hot vapours can cause scalding.
-  Drain before moving the bath. Before draining allow the liquid to cool to below 50°C .
-  Use only liquids specified in the operating instructions, within the specified temperature range. If the alarm lamp is illuminated the liquid temperature may be above its recommended maximum: do not inhale the vapours given off as they may be toxic; liquids should be safely discarded and replaced.
-  For indoor use only. Do not use outdoors.
-  Do not place removable control units on a flammable surface.
-  Set the adjustable overtemperature cut-out after setting or changing the working temperature and reset it at monthly intervals to check that the cut-out is operating correctly (GD120 only).
-  When pumping to external apparatus always use connecting pipe or hose suitable for the operating temperature and liquid used, and check the pipe connections are secure (GD120 only).
-  Do not disconnect any pipes or hoses while they contain very hot or very cold liquid (GD120 only).

For safe operation is it critical that the **GD100** is used with **WATER ONLY**

2. Assembly

2.1 Unpacking

Remove packing materials carefully, and retain for future shipment or storage of the unit.
Packs should contain:

2.1.1 **Controller**

Pump outlet plastic moulding, Operating instructions
Mains lead
Clamp (Where applicable)

2.1.2 **Tanks (Stainless steel)**

Tank
Bridge plate Circulation tray (except S5and S12)
Three screws M3
Screw M5

2.1.3 **Tanks (Plastic)**

Tank
Bridge plate
Clamp

2.2 Assembling control unit for use as bath controller or immersion thermostat

2.2.1 Bath controller-stainless steel tanks, S12, S18, S26 and S38.

The controller can be fitted to the bridge plate in either of two positions. These are, on the left hand side facing into the working area of the bath or on the left hand side of the bath facing out. Fit the controller through the hole in the bridge plate and align the holes with the three screw fixings in the base of the controller. Then using the three M3 screws supplied fix the controller to the plate. Note; the controller position cannot be changed after placing through the hole, it must be lifted out and then rotated.


Fit the lip of the bridge plate into the gap between the cut out in the plastic case and the tank in the top left hand corner of the case such that the lip of the bridge plate is lodged under the lip of the tank. Fix the bridge plate to the bath using the M5 screw, this goes through hole in the rear of the bridge plate and screws into the case.

Fit the circulation tray in the base of the tank with the large cut out in the tray underneath the control unit (models 12L and above).

EN

2.2.1.1 Accessory cooling

Where Grant dip coolers are to be used, fit them to the bridge plate as shown:

2.2.2 Bath controller - stainless steel tank S5

Fit the controller through the smaller of the two cut outs in the bridge plate with the front of the controller facing over the larger cut out. Using the three M3 screws supplied, fit the controller to the bridge plate. With the plate correctly aligned with the controller, three holes in the plate line up with three screw fixings in the base of the controller. Fit the lip of the bridge plate into the gap between the cut out in the plastic case and the tank in the front of the bath such that the lip of the bridge plate is lodged under the lip of the tank. Fix the bridge plate to the bath using the M5 screw, this goes through hole in the rear of the bridge plate and screws into the screw fixing in the case.

2.2.3 Bath controller - plastic tanks P5, P12, P18.

The controller is fitted to the bridge plate so that the controller faces into the working area of the bath and the fixing clamp is to the rear of the controller. Fit the controller through the hole in the bridge plate and using the three M3 screws supplied, fit the controller to the plate. Loosen the knurled clamp fixing screw and push the bridge plate onto the end of the plastic tank, ensuring that the clamp fixing plate goes into the recess in the tank on the outside. Tighten the knurled screw on the clamp making sure that the tank is securely clamped by the two plates on the bridging plate.

2.2.4 Immersion thermostats (for use with non Grant baths).

WARNING: *The liquid container on which the unit is mounted must have the necessary mechanical, chemical and heat resistance.*

Use of control unit with clamp

The clamp is fitted to the controller with the clamp at the rear of the controller and the clamp fixing knob facing outwards.

Using the three M3 screws supplied, fit the clamp to the controller. With the clamp correctly aligned with the controller, three holes in the clamp line up with three screw fixings in the base of the controller.

Clamp the controller to the side of your vessel. Vessels can be used with a maximum wall width of 30mm.


Take care not to over tighten the clamp if you are attaching it to a glass vessel.

CAUTION: *When using the immersion thermostat above 60° C, use a tightly fitting lid around the controller. The purpose of the lid is to prevent vapours from the hot liquid escaping to the atmosphere or being drawn into the controller with the cooling air. Hot vapours should not be inhaled. The lid also improves the heat up rate and reduces liquid evaporation.*

2.2.5 Lids


FG5 for S5

Select the combination of holes on the lid that is suitable for the application.

EN

PF5 for P5 plastic tank

Plastic lids

PL12 for P12 plastic tank

PL18 for P18 plastic tank

Hinged lids (These lids are supplied with hinges for fitting lid to bath)

LG12 for S18 stainless steel tank

LG26 for S26 stainless steel tank

LG38 for S38 stainless steel tank

The lid pack contains:

Lid (with upper part of hinge fitted)

2 hinges, (lower half)

4 screws

Assembly Instructions:

Remove the four bungs from the side of the bath. Attach the support brackets and lower half of the hinges to the bath with the screws provided, but do not tighten screws. Position the lid on the bath and check that lid opens smoothly and is correctly lined up with the side of the bath before tightening screws.

Lid operation:

The lid can be opened by lifting the stainless steel handle at the END of the lid and rotating it around the hinge, or removed by lifting the lid upwards slightly using the handle on the TOP of the lid and moving it away from the hinge.

Polypropylene spheres:

A pack of spheres (PS20) is sufficient to give a single layer for a surface of 1000 cm².

2.2.6 Use of pump on GD120.

To use the pump, remove the pump bypass moulding and replace with the pump output moulding.

Remove the two screws which hold the pump bypass moulding to the control unit, and using the screws fit the pump output moulding. Ensure that the two sealing 'O' rings are correctly fitted into their grooves before fitting. Retain the pump bypass with its longer screws for refitting when the pump is no longer required.

Connect the outlet pipe, marked with an upward arrow on the left hand pipe, to the external circuit. Connect the return pipe to the inlet pipe, marked with a downward arrow. Secure the pipes with clips.

The pump will pump liquid around a closed external circuit, which is one not open to the atmosphere. It cannot be used for circulation through an external open tank.

WARNING: If the operating temperature is greater than 60° C, ensure that the pipes are securely clamped to pump outlet and to the external circuit, and that the tubing is suitable for both the liquid and the temperature.

2.3 Use at lower temperatures.

Due to heat input from the pump and stirrer motor, the minimum operating temperature without accessory cooling depends on the size of the bath. The smallest baths, P5 and S5, have an operating temperature of approximately 10° C above ambient without a lid (15° with a lid).

The other baths can be used at a temperature of 5° C above ambient without accessory cooling.

For lower temperature operation, accessory cooling is required.

EN

3. Operation

3.1 Controls and indicator lamps (see fig. 1)

3.1.1 Power

The power switch is located in the lower left hand corner of the front panel.

If the unit does not come on, check whether the "alarm" lamp is illuminated. (See 3.1.3).

3.1.2 Heater lamp

The heater lamp (orange) on the front panel indicates when the heater is on. While the bath is warming up, the lamp is on continuously, and it starts to flash as the temperature reaches set point. When the unit is controlling at set temperature it flashes intermittently.

3.1.3 Alarm lamp

The alarm lamp (red) on the front panel illuminates when either the liquid level is low, or on the GD120 the over temperature cut-out has operated.

The low liquid level switch consists of a stainless steel float, and switches the heating and stirring off if the liquid level falls.

When the alarm lamp is illuminated, first check the liquid level and refill if necessary. If, on the GD120 it is still illuminated, press the overtemperature cut-out reset, if the alarm lamp does not go out, set overtemperature as described in 3.3.1.

3.1.4 GD100/GD120 Identification of units and controls.

(See the following pages)


Fig. 1 Front Panel

3.1.4.1 Temperature display

The temperature display normally shows the temperature in °C of liquid being controlled. It displays set temperature when the set button "S" is pressed. The display flashes to indicate that it is in "Set temperature" mode and will return to actual temperature after a few seconds.

3.1.4.2 Set temperature controls

Press "S" to change display from liquid temperature to set temperature, rotate the "Navigator" knob to set the required operating temperature. Increasing the speed of rotation of the knob increases the step size of the set temperature increment. The display flashes when in the set temperature mode. When the required temperature is indicated, press "S" again. This value is then stored in memory and the unit will retain this value after switch off.

3.1.4.3 Set over temperature control (GD120)

The cut-out operating temperature is adjustable. If the liquid temperature rises above the cut-out set temperature the heater and stirrer are switched off. To reset the cut-out push the outer ring and release once the bath has cooled down.

EN

Setting the over temperature control.

There are two methods of setting the over temperature cut out:

Option one. (Quick method)

Set the temperature to the required value and leave the bath to stabilise for at least 5 minutes after the set point has been reached. Turn the control slowly anticlockwise, using a screwdriver, until the alarm lamp comes on. Press the reset and at the same time turn the control slowly clockwise until the alarm lamp goes out. This gives an over temperature trip point of approximately 10° to 30° C above set temperature.

Option two. (Precision method)

Set the temperature to the cut out level required and leave the bath to stabilise for a least 5 minutes after the set point has been reached. Turn the control slowly anticlockwise, using a screwdriver, until the alarm lamp comes on. Allow the bath to cool, then press the reset button and the unit will start. (Bath will have to cool by 15° to 30° before the reset will work). This gives a over temperature trip point of the set value. Now decrease the set temperature to the working set point.

3.1.4.4 Function control F

Press F to enter function select mode, for pre-set temperature, high alarm, timer and liquid type selection.

3.2 Filling

Before filling the bath, disconnect the mains power supply. Fill the tank to a maximum level of 10mm from the top with a liquid suitable for your working temperature.

Specified Liquids:-

-30° to 30° C 50% water, 50% antifreeze (inhibited ethylene glycol) **WARNING - TOXIC:** Follow manufacturer's instructions
For safe disposal consult your Local Environmental Health Office.

0° to 30° C 80% water, 20% antifreeze (inhibited ethylene glycol) **WARNING - TOXIC:** Follow manufacturer's instructions
For safe disposal consult your local Environmental Health Office.
Use a lid to reduce the dilution of the mixture caused by condensing water vapour from the air, and to maintain the cool down rate.

+5 to 80° C Water, recommended temperature range.
Use a lid or polypropylene spheres above 60° C to reduce evaporation and heat loss.

80° to 99.9° C Water can be used but care should be taken with its use as hot water vapour can be dangerous. A lid or polypropylene spheres must be used. At these temperatures the temperature control will be reduced due to localised boiling, and there may be a large loss of water due to evaporation. The control units should not be used to boil water.

Note: For the long term reliability of water baths it is important to use oxygenated water that is free from ions and minerals that can cause corrosion. Grant recommend the use of distilled water and de-ionised water from modern ion exchange systems with replaceable cartridges. Do not use water from a system that incorporates a salt back flush system as this can leave sodium ions that are corrosive to stainless steel.

70° to 120° C Silicone fluid with the following characteristics;
Viscosity 20 centistokes
Flash point ≥230°C
Fire Point ≥280
E.g. Dow Corning silicone fluid DC200/20 Follow manufacture's instruction.
For safe disposal consult your Local Environmental Health Office.

3.3 Switching on and setting up.

Connect the control unit to a grounded (earthed) electrical power supply with voltage and frequency within the range specified on the serial number plate.
DO NOT SWITCH ON BEFORE FILLING THE TANK.

Note: GD100/GD120 When changing settings the display will revert to home display if no changes are made within 30 seconds.

3.3.1 GD120.

Connect the pump to your external circuit if required Before switching on for the first time, use a screwdriver to turn the set over temperature control fully clockwise. Press the set over temperature reset to ensure that it is reset. Switch the power on Press the "S" button and rotate the "Navigator" knob to set the required operating temperature. Press "S" button again to fix this value. Wait until the liquid temperature has stabilised at the operating temperature. The display flashes when in the set temperature mode.

If required the over temperature should now be set using the instructions in 3.1.4.3.

EN

3.3.2 GD100.

Switch the power on Press the S button and rotate the “Navigator” knob to set the required operating temperature. Press S button again to fix this value.


3.4 GD100 and GD120 (Fig. 2) Preset temperatures

Setting preset temperature values

Press F, Rotate the “Navigator” knob until required pre set adjust number [E^1 , E^2 or E^3] is indicated on the display. (The display will alternate between the ‘t’ and the preset temperature). Press S Rotate the **Navigator** knob until desired temperature value is indicated on the display. Press S. Controller accepts this as current set temperature value and also stores it as a preset value. Running at preset temperature value Press F Rotate the set temperature control until required pre set number [E^1 , E^2 or E^3] is indicated on the display. Press S Press S Controller accepts this as current set temperature.

Preset temperatures can be adjusted over full range of unit (-20° to 120°C), but will only be accepted when over 100°C if set to [H_1L] and below 100°C when set to low temperature liquid [L_1L] When set to water [H_2O], preset values greater than 100°C will revert to 100°C when S is pressed to leave preset menu.

Fig. 2 Set point and Pre-Set Temperature selection and adjustment (GD100, GD120)

3.5 **GD120 Alarm (ALR), Clock (CLK) and Liquid selection ($\text{H2O}, \text{OL}$) and adjustment (Fig 3)**

Set High Level Alarm

Press **F** three times then **S**, (ALR) displayed rotate navigator knob control to select alarm (on or OFF), when (OFF) is displayed, pressing **S** or **F** returns to bath temperature display with no alarm set. When (ON) is displayed pressing **S** causes the alarm to sound and the display to flash with alarm level indicated. This can now be adjusted using the navigator knob, pressing **S** returns with alarm ON. When the bath temperature exceeds set level, the display shows ($-\text{AL}-$). The alarm does not affect the heating control. Alarms are cancelled after either button is pressed during alarm or if alarm condition has been removed.


Set clock countdown

Press **F** twice then **S**, (CLK) displayed use navigator knob to select clock (On or OFF), when (OFF) is displayed, pressing **S** or **F** returns to bath temperature display with no timer set. When (On) is displayed pressing **S** causes the display to flash with time (in minutes) indicated. Rotate navigator control to set time required. Pressing **F** returns with timer (OFF), pressing **S** returns with timer (On). With clock (On) display alternates bath temperature and time remaining. At end of the time period, display shows (End). The clock acts as a display timer and has no effect on the control.

Set Liquid

Press **F** four times, use navigator knob to select Water (H2O), Oil (OL) or low temperature liquid (LET), pressing **S** returns to bath temperature with liquid selected. When (H2O) is selected temperature range is 0° to 100°C, for oil the range is 0° to 120°C, low temperature liquid (for example 50% water / 50% glycol) can be set in the range -20° to 50°C.

Fig. 3 Alarm, time and liquid selection and adjustment (GD120)

3.6. Calibration

Calibration should be carried out using a traceable reference thermometer with an accuracy of at least 0.1°C.

Calibration procedure for GD120.

The calibration of the GD120 is in two parts, low temperature and high temperature.

Part one.

Liquid temperature of the bath must be below 20°C, recommended temperature of 5°C.

- 1** Place the unit on a bath, turn on and set unit to 5°C. (This is the recommended temperature anything from 0° to 20°C can be used.)
- 2** Allow the unit to stabilise for at least 5 minutes after a stable temperature condition has been achieved.
- 3** Measure the liquid temperature by placing a traceable reference thermometer with an accuracy of at least 0.1°C into the bath.
- 4** Press **F** and **S** simultaneously for approximately 5 seconds or until **LCAL** is displayed.
- 5** Press **S** and temperature will be displayed.
- 6** Turn 'Navigator' knob so that the average temperature (over a period of 15 seconds) of the display equals the temperature on the external thermometer.
- 7** Press **S** and display shows **HCAL**
- 8** Press **F** and display shows temperature


Part two

Liquid temperature of the bath must be above 70°C, recommended temperature of 70°C.

- 1** Place the unit on a bath, turn on and set unit to 70°C. (This is the recommended temperature anything from 70° to 90°C can be used.)
- 2** Allow the unit to stabilise for at least 5 minutes after a stable temperature condition has been achieved.
- 3** Measure the liquid temperature by placing a traceable reference thermometer with an accuracy of at least 0.1°C into the bath.
- 4** Press **F** and **S** simultaneously for approximately 5 seconds or until **LCAL** is displayed.
- 5** The display should read **LCAL**
- 6** Rotate 'Navigator' knob until display shows **HCAL**
- 7** Press **S** and temperature will be displayed.
- 8** Turn 'Navigator' knob so that the average temperature (over a period of 15 seconds) of the display equals the temperature on the external thermometer.
- 9** Press **S** and display shows temperature

That completes the calibration procedure. Fig 4 shows the menu structure.

Fig. 4

4. Fault diagnosis

SYMPTOM	POSSIBLE CAUSE	ACTION REQUIRED
Unit does not operate	Unit not switched on	Switch on
No Display	Unit not plugged into power supply Fuse blown in the unit or for UK units, in the plug Electrical power supply failure	Plug in, switch on Check and replace Check that other electrical appliances on the same circuit are working
Display shows PF 0	Probe open circuit	Have unit checked by a competent person
Display shows PF 5	Probe short circuit	Have unit checked by a competent person
Alarm light on	Overtemperature cut-out has operated (GD120) Low liquid level, float switch has operated	Reset cut-out Top up liquid
Temperature does not rise when expected	Set temperature is lower than liquid temperature Set temperature is too close to ambient temperature Temperature control circuit fault	Check set temperature Raise the set temperature or fit accessory cooling Have unit checked by a competent person
SYMPTOM	POSSIBLE CAUSE	ACTION REQUIRED
Temperature continues to rise when not expected	Set temperature is higher than liquid temperature Temperature control circuit fault	Check set temperature Have unit checked by a competent person
Stirred motor does not rotate	Stirrer or propeller obstructed Pump obstructed If stirrer motor shaft rotates freely then motor thermal fuse may have blown	Clear obstruction Clear obstruction Have unit checked by a competent person

5. Specification

5.1 Environmental

This equipment is for indoor use and will meet its performance figures within an ambient temperature range 10° to 35°C with maximum relative humidity of 80% and is safe for use in an ambient temperature up to 40°C.

Operation at Installation category II (transient voltages), pollution degree 2, in accordance with IEC 664 for altitudes of up to 2000 metres.

5.2 Heating

Temperature setting range:

GD100	0°C to 100°C
GD120	-20°C to 120°C

Minimum operating temperature:

GD100	0°C to 100°C
GD120	-20°C to 120°C

Stability (DIN58966) at 37°C:

GD100	+/- 0.05
GD120	+/- 0.05

Uniformity at 37°C:

+/- 0.1

Supply voltage range:

230V ± 10% 50Hz

120V ± 10% 60Hz

Power rating:

230V 1.6kW

120V 1.8kW

Heater power:

230V 1.4kW

120V 1.4kW

Pump performance:

GD120

Max. Head 310mbar
at zero flow.

Max. Flow 17 lit/min
At zero head.

EN

Over temperature protection is by variable resettable over temperature cut-out.

Low liquid level is sensed by a float switch.

	S5	S12	S18	S26	S38
Tank capacity (Litres)	5	12	18	26	38
Tank dimensions l/w/d (mm)	300/150/150	325/300/150	505/300/150	505/300/200	600/300/200
Liquid depth min/max (mm) *	80/140	70/130	70/130	120/180	120/180

*Maximum depth can be increased by 10mm, by removing the circulation tray in 12, 18,

	P5	P12	P18
Tank capacity (Litres)	5	12	18
Tank dimensions l/w/d (mm)	300/150/150	325/300/150	505/300/150
Liquid depth min/max (mm)	80/140	70/130	70/130

Heat up time in minutes from 25° C to 70° C with water and lid on:-

	230V	120V
S5	10	16
S12	25	35
S18	38	53
S26	50	70
S38	70	100

6. Maintenance

6.1 Routine Maintenance

All Grant laboratory products are designed to comply with IEC61010-1 and can be flash tested. Some are fitted with radio frequency interference suppressors. Therefore it is recommended that only a dc. test is performed.

The over temperature cut-out on the GD120 should be checked periodically by turning the set over temperature control with a screwdriver anticlockwise until the alarm lamp comes on. The cut-out should then reset and set-up (see 3.3.1). If the alarm lamp fails to light with the knob fully turned anticlockwise the unit should be checked by a competent person.

The liquid level protection should also be checked periodically by lowering the level of liquid in the bath and noting that the unit cuts out with the top turn of the heater still immersed in the liquid.

When hoses are fitted to the pump they should be inspected periodically and replaced as necessary to avoid hose failure.

No other routine maintenance is required.

6.2 Cleaning

The cases can be cleaned with a damp cloth after disconnection. Do not use solvents.

The immersed parts can be cleaned using proprietary heating element cleaners.

CAUTION: these may be toxic - follow the cleaner manufacturer's instructions.

Before using any decontamination or cleaning method except that recommended, check with our Service Department, or in other countries with our distributor, that the proposed method will not damage the equipment

EN

6.3 Replacement of Fuses

230v

Disconnect the unit from the power supply socket. Remove the IEC power plug from the rear of the unit. Press down the drawer catch. Pull out the fuse drawer, check and replace with the correct fuses if necessary, as follows:

Two off, 1.25 x 0.25 inch ceramic quick acting rated:-

250V 10AF

Push back the drawer and replace the IEC plug.

120v

Disconnect the unit from power supply socket.

Locate fuses on rear panel next to power cable inlet.

Use a screwdriver of the correct size to turn the fuse holder a quarter turn anti-clockwise to release the fuse holder.

Check and replace with the correct fuses if necessary. The fuses should be:

Two off, 1.25 x 0.25 inch ceramic quick acting rated:-

250V 20AF

Replace fuse holder.

7. Guarantee and Service

7.3 **Guarantee**

When used in laboratory conditions and according to these working instructions, this bath is guaranteed for **THREE YEARS** against faulty materials or workmanship.

Service repairs outside of the guarantee period, carried out by Grant Instruments, carry a further one year guarantee.

7.4 **Service**

For service, the control unit can be removed and can be returned for repair to our service department .Outside the UK, return to your local distributor.

DO NOT RETURN THE TANK

Return address

Grant Instruments (Cambridge) Ltd.
Service Department
Station Road
Shepreth
Cambridgeshire
SG8 6GB
England

EN

8. Compliance

8.0 Disposal and WEEE

In Europe, at the end of its life the GR150, GP200 must be disposed of in accordance with the WEEE directive. For information regarding WEEE collections in the UK please contact our B2B Compliance Scheme directly on 01691 676 124.

Grant Instruments complies fully with the Waste Electrical & Electronic Equipment (WEEE) regulations 2006. We are a member of the B2B compliance scheme (Scheme Approval Number WEE/MP3338PT/SCH), which handle our WEEE obligations on our behalf. Grant Instruments have been issued with a unique registration number by the Environmental Agency, this reference number is WEE/GA0048TZ.

For other countries please contact your equipment supplier.

For General WEEE information please visit: www.b2bcompliance.org.uk

8.1 RoHS

All the products covered by this manual comply with the requirements of the EU RoHS Directive (Directive 2002/95/EC). This means the products are free of Lead and other hazardous substances covered by the directive.

8.2 Electrical safety and electromagnetic compatibility (EMC)

All the products covered by this manual comply the requirements of the Low Voltage Directive (2006/95/EC) for Electrical safety and the EMC directive (2004/108/EC) for Electromagnetic compatibility. See the Declaration of Conformity on the inside back page

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GRANT INSTRUMENTS (CAMBRIDGE) LTD.
Service Department
Station Road
Shepreth
Cambridgeshire
SG8 6GB
Inglaterra

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Declaration of Conformity

Manufacturer	GRANT INSTRUMENTS (CAMBRIDGE) LTD, Shepreth, Cambridgeshire SG8 6GB
Equipment Name	GD100, GD100L, GD100IN, GD100INL, GD120, GD120L GR150, GR150L, GP200, GP200L
Equipment Description	Optima series of general purpose stirrers / circulators
Directives	EMC Directive 2004/108/EC LVD Directive 2006/95/EC
Including Accessories	FF17 and LL17 external temperature probes
CE mark first applied	2000

Applied Harmonised Standards	BS EN 61326-1:2006 Electrical equipment for measurement, control and laboratory use. EMC requirements. General requirements
	BS EN 61010-1:2001 Safety requirements for electrical equipment for measurement, control and laboratory use. General requirements
	BS EN 61010-2-010:2003 Safety requirements for electrical equipment for measurement, control and laboratory use. Particular requirements for laboratory equipment for the heating of material

This product complies with the requirements of the above Directives when used with leads up to 3m long, compliance may be affected by using longer leads.

I confirm that this apparatus conforms to the requirements of the above Directives

L.J.C.....

Dated... 20/07/09

Ludo Chapman
Managing Director
Grant Instruments (Cambridge) Ltd.

USA

This device complies with Part 15 of the FCC rules. Operation is subject to the following two conditions: (1) this device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

Grant

**Grant Instruments
(Cambridge) Ltd**

Shepreth,
Cambridgeshire
SG8 6GB

Tel: +44 (0)1763 260811
Fax: +44 (0)1763 262410
Email: labsales@grant.co.uk
www.grant.co.uk