VWR

Microscope  **VisiScope®** series 400

INSTRUCTION MANUAL

Model	European Catalogue Number
IT405	630-2080

Version: 6
Issued: 20, 10, 2014

Legal Address of Manufacturer

Europe

VWR International bvba
Researchpark Haasrode 2020
Geldenaaksebaan 464
B-3001 Leuven
+ 32 16 385011
<http://be.vwr.com>

Country of origin: **ITALY**

Table of Contents

Warning
Safety Information
Package Contents
Intended use
Symbols and conventions
Unpacking
Installation steps
Product Specifications
Overview
Operation
Troubleshooting
Repair and maintenance
User replaceable accessories and spare parts
Technical service
Warranty
Compliance with local laws and regulations
Disposal

Warning

This microscope is a scientific precision instrument designed to last for many years with a minimum of maintenance. It is built to high optical and mechanical standards and to withstand daily use.

We remind you that this manual contains important information on safety and maintenance, and that it must therefore be made accessible to the instrument users.

We decline any responsibility deriving from incorrect instrument use that does not comply with this manual.

Safety Information

Avoiding Electrical Shock

Before plugging in the power supply, make sure that the supplying voltage of your region matches with the operation voltage of the equipment and that the lamp switch is in off position.

Users should observe all safety regulations of the region. The equipment has acquired the CE safety label. However, users do have full responsibility to use this equipment safely.

Please follow the guidelines below, and read this manual in its entirety to ensure safe operation of the unit.

Package Contents

DESCRIPTION	QUANTITY
Microscope stand with head, stage, nosepiece, focus	1
WF 10X/22 eyepieces	2
IOS 10x objective	1
IOS 20x objective	1
IOS 40x objective	1
Illuminator	1
Sample holder (square)	3
Sample holder (round)	2
LBD and IF550 filter	2
Phase contrast slide	1
Mechanical stage	1
Dust cover	1
Power cable	1

Intended use

For research and teaching and teaching use only. Not intended for any animal or human therapeutic or diagnostic use.

Symbols and conventions

The following chart is an illustrated glossary of the symbols that are used in this manual.


CAUTION

This symbol indicates a potential risk and alerts you to proceed with caution

Unpacking

The following figure shows the installation sequence of the components. The number in the figure shows the installation steps.

- Before installing, be sure every component is clean. Take care not to scratch parts or surfaces.
- Store the supplied hexagon wrench in a safe place. When changing the components, you will need it again.

Installation steps



Installing the condenser illumination unit and the LED housing

1. Insert the condenser illumination unit (1) into the bracket, according to figure.
2. Turn the condenser illumination unit clockwise about 90°, with the “AS” mark of filter holder (3) facing forwards. Align the screw of the condenser illumination unit and the hole of the holder, then screw the bolt in the hole with the supplied hexagon spanner.
3. Insert the connector plug (2) into the connector jack (4).
4. Push the LED housing gently into the holes of the illumination unit.

Installing the objectives

1. Turning the coarse focusing knob (1) till the nosepiece reaches to its lowest position.
 - **For a safe transport, the nosepiece is placed in the lowest position and the tension adjustment collar (2) is adjusted to the appropriate tension when the microscope leaves the factory.**
2. Screw the lowest magnification objective on to the turret from the right side, then turn the turret clockwise. Mount the other objectives in the same way, following the sequence from low to high.
 - **Note: the objectives can also be installed through the stage opening.**
 - Clean the objectives regularly. In inverted microscopes, the objectives are very sensitive to dust.
 - To prevent dust and contamination from entering the microscope, cover all the unused holes with dust caps (3).
 - When operating, use the low magnification objective (10X) to search and focus the specimen, then switch to higher magnifications.
 - When switching between objectives, slowly turn the nosepiece until it clicks. The click means that the objective is in the right position, in the center of the light path.


Installing the stage extension and the mechanical stage

The stage extension can be installed on either side of the stage to enlarge the working surface. The mechanical stage must be installed on the side opposite the extension. The stage extension can be installed on either side of the stage to enlarge the working surface. The mechanical stage must be installed on the side opposite the extension.

For right-handed operators, the mechanical stage is normally installed on the right side.


1. Installing the stage extension: Screw the bolts (1) on to the extension, then mount the extension from below the stage. Screw it to a firm fit.
2. Installing the mechanical stage: As for the extension, the mechanical stage is fixed by two bolts under the stage.

Installing the stage insert


1. When using the glass stage (1), make sure that the stage is horizontal.
2. Install the stage insert in the stage opening.

Turn the disk until the V-groove faces the user. This simplifies objective identification.

Installing the eyepieces

1. Remove the cap of the eyepiece tubes (1).
2. Insert the eyepieces (2) into the tubes.

Installing the color filters

- **Be sure the color filter has cooled down completely before you change them. Remove the filter holder (1), then install the color filters (2) you need.**

Mount the color filter flat as shown in (3), verifying that they are not tilted.

- **If the color filter is tilted or otherwise out of place (4), it may fall.**

The color filters can be stacked in the holder. This allows to install as many filters as needed, as long as the whole thickness is less than 11 mm.

Connecting the power cord

- **Do not tension or otherwise stress the power cord. If bent, the wires may easily break.**



1. Turn the main switch (1) to "O" (off) before connecting the power cord.
2. Insert the plug into the power jack (2) of the microscope.
3. Plug the power cord (4) into the mains socket. Check for a safe connection.

- **Please use the supplied power cord. If lost or damaged, please refer to qualified service.**
- **Connect the power cord to a grounded (earthed) power supply only.**

Replacing the fuse

Before replacing the fuse, turn the main switch (1) to "O" (off) and unplug the power cord. Rotate the fuse support (3) out of the holder using a straight screwdriver. Insert a new fuse in the support, then rotate the support back into the holder.


- **Fuse rating: see back of the microscope.**


Product Specifications

Illumination	Light source type P-LED [®] with white 8W LED; light intensity control using a knob on left side of the frame. Color temperature: 6300K LED average life time approx. 50.000h Voltage: 110/240Vac, 50/60Hz, 1A ; Fuse: T500mA 250V Max power required: 13W
Observation Modes	Brightfield, phase contrast
Focusing	Coaxial coarse and fine focusing mechanism (graduated, 0.002mm) with upper stop, to prevent the contact between objective and specimen. Adjustable tension of coarse focusing knob.
Stage	Fixed stage, dimensions 250x230 mm. Mechanical stage mountable on the right side of the stage, X-Y translation range 120x80 mm, with metallic interchangeable inserts for slides, Petri dishes, Terasaki, multi-Well plates, etc. Pair of side extensions to expand the surface of the stage. Glass stage insert with hole for small dimension specimens.
Nosepiece	Quintuple revolving nosepiece, rotation on ball bearings.
Head	Trinocular observation head, inclined 45°. Diopter adjustment on left eyepiece. Interpupillary adjustment 55-75 mm. Splitting ratios (eyepieces-photo tube): 100-0 / 20-80
Eyepieces	Wide field eyepieces WF10X/22 with field number 22.
Objectives	Infinity corrected optical system IOS (Infinity Optical System). Plan-achromatic LWD objectives infinity corrected, for thickness 1.2 mm, made by following objectives: -) Plan-achromatic IOS LWD 10XPh, N.A. 0.25, W.D. 7.94 mm -) Plan-achromatic IOS LWD 20XPh, N.A. 0.40, W.D. 7.66 mm -) Plan-achromatic IOS LWD 40XPh, N.A. 0.60, W.D. 3.71 mm All objectives are treated with an anti-fungus treatment.
Condenser	LWD condenser, N.A. 0.30, working distance 72 mm. The condenser can be removed to extend the working distance up to 150 mm. Precentered slider with phase rings 4x/10x and 20x/40x.
Dimensions	HEIGHT: 485 mm WIDTH: 300 mm DEPTH: 550 mm WEIGHT: 9 kg
Accessories	Hard case, instruction manual and dust cover included.

OVERVIEW

Operation

MICROSCOPE BASE

Turning on the LED

Connect the power, turn on the main switch (1).



Adjusting the brightness

Turn the brightness adjustment knob (2) to increase and decrease the brightness.

Adjusting the tension

► **The coarse focusing knob (1) is pre adjusted to a tight tension upon leaving the factory.**

If the nosepiece drops down by itself, or the specimen de focuses while adjusting the fine focus knob (3), the coarse focus knob is too loose. Turning the tension adjustment collar (2) in the indicated direction tightens the coarse focus tension (1). Rotate in the opposite direction to decrease the tension.

STAGE

Setting the specimen



► **For the best image quality, use flasks, Petri dishes and slides with a 1.2 mm thickness.**

Using $\Phi 35$ mm culture dishes:

You can place a $\Phi 35$ mm culture dish on the stage by using the standard center board (1) of the stage.

Using the mechanical stage:

1. When observing 96-well or 24-well micro-titration plates, please fasten them directly with the stage clip (2).
2. When fastening other plates, please use the adaptor plates supplied with the mechanical stage:
 - Terasaki bracket (3) for Terasaki boards
 - Culture dish bracket (4) for $\Phi 35$ mm culture dishes
 - Object slide bracket (5) for object slides and $\Phi 54$ mm culture dishes
3. Turning the X and Y knobs (6,7), move the specimen to the required position. (Movement Range: 120 (width) $\times$ 78 (length) mm).

Moving the specimen

Move the specimen to the desired position by freehand or by turning the knobs of the mechanical stage.

- **When switching objectives, take care not to touch the adaptor plates with the objectives, as their weight may damage the front lens.**

VIEWING TUBE

Dioptric adjustment

1. Look into the right eyepiece with your right eye only, and focus on the specimen.
 2. Look into the left eyepiece with your left eye only. If the image is not sharp, use the dioptric adjustment ring (1) to compensate.
- **The adjustment range is ± 5 diopter. The number indicated on the adjustment ring graduation should correspond to the operator's dioptric correction.**

Adjusting the interpupillary distance

Observing with both eyes, hold the two eyepiece prism assemblies. Rotate them around their common axis until the fields of view coincide.

- **The graduation on the interpupillary distance indicator (3), pointed by the spot “.” (2) on the eyepiece holder, shows the distance between the operator's eyes.**

The range of the interpupillary distance is 48-75mm.

Selecting the light path

Pull the light path selector lever (1) sideways using your thumb, selecting the light path you need.

LIGHT PATH SELECTOR LEVER	BRIGHTNESS	APPLICATION
In	100% used for binocular observation	Binocular observation
Out	20% used for binocular observation, and 80% used for video or photography	Binocular observation, television, and micrography or video can be operated simultaneously


ILLUMINATION UNIT

Using color filters

Selecting the appropriate color filters according your need.

You can stack a group of color filters in the filter holder, if you ensure that they are level and that the whole thickness is less than 11mm.

COLOR FILTER	USE
In	Single contrast color filter (green) used for phase contrast microscopy
Out	Color temperature compensation color filter blue used for bright field observation and microphotography

Using the aperture diaphragm

When in bright field observation, the aperture diaphragm controls the numerical aperture of the illumination system. When the numerical aperture of the objective and the aperture of the illumination system match, the highest resolution is achieved.

To identify the aperture diaphragm, the eyepiece can be removed or the centering telescope can be used. The aperture can be changed by moving the aperture adjustment lever. ((1)is the image of the aperture diaphragm, (2) is the edge of the objective).

Generally, when observing a fully chromatic specimen, you need to set the size of the condenser to 70-80% of the aperture of the objective. When observing unstained samples (e.g. bacteria), start from 70% and slowly turn the aperture diaphragm lever clockwise.

IDENTIFYING THE COMPONENTS

Phase contrast objectives

The magnifications of the optional phase contrast objectives are: 10X and 20X. These objectives are marked “PH”. Mount the objectives on the turret using the same procedures as the standard objectives

Phase contrast slider

Adjustable phase slider.

- The light ring is pre-centered when the microscope leaves the factory. It should therefore need no further adjustment. If a re-centering is needed, it can be performed via the two side bolts.
- The 4X/10X light ring (1) must be used with 4X and 10X phase contrast objectives, the 20x/40x light ring (2) with the 20x and 40x and the opening (3) is used for bright field.

SETUP AND USE


Installing the phase contrast slider

1. Insert the slider into the illumination system, printed face up.
2. Pull the slider into the desired position, to the click stop.
3. When in phase contrast observation, keep the aperture diaphragm adjustment lever on the "O" (open) position.

Centering the ring

► **Usually this operation is not needed. If necessary, please proceed with the following steps:**

1. Place a specimen on the stage and focus it.
 2. Take out the eyepiece from the tube without the dioptic adjustment, and replace it with the centering telescope (CT).
 3. Check that the phase ring and the objective correspond, and that both are steadily set on a click stop.
 4. Use the CT to focus on the light ring's image (5) and the phase contrast ring's image (2). If the light ring's image is not sharp, adjust the CT's eyepiece until you can see a clear image of the light ring (2).
 5. Adjust the bolts of the two centering holes in the phase contrast slider using a screwdriver until the light ring center and the phase contrast ring center coincide.
 6. The 10X and the 20X phase contrast objectives use the same ring on the phase contrast slider. The coincidence of the light ring center and the phase contrast center must be verified with both objectives.
- **If the light ring is centered incorrectly, the contrast will be severely impaired.**
- **The phase ring may need recentering during and after observation of very thick specimens.**
- **The phase ring may show an apparent misalignment if the cover glass is not flat.**

VIDEO PORT (VIDEO)

Selecting the light path

1. To activate the video port, pull out the light path selector lever.
- ▶ **For observation of dark specimens, first focus, then pull out the lever.**

Installing the video adapter

- ▶ **Usually this operation is not needed. If necessary, please proceed with the following steps:**
1. Loosen the locking bolt (1) on the trinocular viewing tube, and take out the dust cap (2).
 2. Set up the video adapter on the camera according to its instructions.
 3. Install the adapter into the tri-through port, and screw down the bolt (1).

Focus

During a binocular observation at 20% brightness, look at the image on the video imaging system, re-focusing the adapter if necessary.

MICROPHOTOGRAPHY

Installing the photography adapter

1. Loosen the locking bolt (1) on the trinocular viewing tube, and take out the dust cap (2).
 2. Install the photography adapter into the tri-through port according to its instructions, and screw down the locking bolt (1).
 3. Attach the camera ring (if any) to the adapter.
 4. Attach the camera to the ring.
- Warning: for some cameras (mainly reflex) the ring is not included with the microscope, and it should be supplied by the user.
 - For the photography of dark specimens, obscure the eyepieces and the viewfinder with a dark cloth in order to reduce stray light.
 - The camera magnification can be calculated as objective magnification × camera + lens magnification.
- ▶ **When shooting with a SLR, the mirror movement may cause camera movement. Please lift the mirror, use long exposure times and use an extension cord.**

Troubleshooting

Review the information in the table below to troubleshoot operating problems.

PROBLEM	CAUSE	SOLUTION
I. Optical Section:		
1. The illumination is open, but the field of view is dark.	The plug of the LED holder is not connected to the illumination set	Connect them
	The brightness is too low	Adjust to a proper setting
	Too many colour filters have been stacked	Minimize the number of the filters
2. The edge of the field of view is vignetted or the brightness is asymmetric.	The nosepiece is not in the correct position	Turn the nosepiece to a click stop
	The color filter is partially inserted	Insert the filter to full depth
	The phase contrast slider is not in the proper position	Move the slider to a click stop
3. Dust and stains can be seen in the field of view.	There are stains and dust on the specimen	Clean the specimen
	There are stains and dust on the eyepiece	Clean the eyepiece
4. There is an apparent double image.	The size of the aperture diaphragm is too small	Open the aperture diaphragm
5. Poor image quality: The image is not sharp The contrast is not high The details are not clear The phase contrast is low.	The nosepiece is not in the center of the light path	Turn the nosepiece to a click stop
	The aperture diaphragm in the view of field is opened too much or too little	Adjust the aperture diaphragm
	The lenses (condenser, objective, eyepieces are culture dish) is dirty	Thoroughly clean all the optical system
	In phase contrast observation, the bottom thickness of the sample is more than 1.2mm	Use a sample holder whose bottom thickness is less than 1.2mm
	A bright field objective is used for phase contrast observation	Switch to a phase contrast objective
	The condenser ring is not aligned with the objective phase ring	Adjust the condenser ring to match the objective phase ring
	The light ring and/or the phase contrast ring is not centered	Adjust the bolts to center them
	The objective used is not compatible with the phase ring	Please use a compatible objective
	The phase contrast depends on the sample position	The sample holder is not flat. Move the sample around until a compatible area is found.

6. One side of the image is out of focus.	The nosepiece is not in the center of the light path	Turn the nosepiece to a click stop
	The specimen is out of place (tilted)	Place the specimen flat on the stage.
	The optical performance of the sample cover glass is poor	Use a cover glass of better quality
II. Mechanical Section:		
1. The coarse focus knob is hard to turn.	The tension adjustment collar is too tight	Loosen the tension adjustment collar
2. The focus is unstable.	The tension adjustment collar is too loose	Tighten the tension adjustment collar
III. Electric section		
1. The LED doesn't turn on.	No power supply	Check the power cord connection
2. The brightness is not enough	The brightness adjustment is low	Adjust the brightness
3. The light blinks	The power cord is poorly connected	Check the power cord
IV. Viewing tube assembly		
3. The field of view of the two eyes is different	The interpupillar distance is not correct	Adjust the interpupillar distance
	The dioptic correction is not right	Adjust the dioptic correction
	The viewing technique is not correct, and the operator is straining the eyesight	When look into the objective, do not stare at the specimen but look at the whole field of view. Periodically, move the eyes away to look at a distant object, then back into the objective
V. Microphotography and video		
1. The image is unfocused	Incorrect focussing	Adjusting the focus system as in the present manual
2. The edge of the image is unfocussed	To some degree, it is inherent to the nature of achromatic objectives	The problem can be minimized by a correct setting of the aperture diaphragm
3. Bright patches appear on the image	Stray light is entering the microscope through the eyepieces and through the camera viewfinder	Cover the eyepieces and the viewfinder with a dark cloth

Repair and maintenance

Microscopy environment

This microscope is recommended for use in a clean, dry and shock free environment with a temperature of 0-40°C and a maximum relative humidity of 85 % (non condensing). Use a dehumidifier if needed.

To think about when and after using the microscope

The microscope should always be kept vertical when moving it so that no moving parts, such as the eyepieces, fall out. Never mishandle or impose unnecessary force on the microscope.

Never attempt to service the microscope yourself.

After use, turn down the illumination intensity control and turn the light off. Cover the microscope with included the dust cover, and keep it in a dry and clean place.

Electrical safety precautions

Before plugging in the power supply, make sure that the supplying voltage of your region matches with the operation voltage of the equipment and that the lamp switch is in off-position.

Users should observe all safety regulations of the region. The equipment has acquired the CE safety label. However, users do have full responsibility to use this equipment safely.

Cleaning the optics

If the optical parts require cleaning first use compressed air.

If that is not sufficient use a soft lint-free cloth with water and a mild detergent.

And as a final option use the piece of cloth moistened with a 3:7 mixture of ethanol and ether.

Note: ethanol and ether are highly flammable liquids. Do not use them near a heat source, near sparks or near electric equipment. Use these chemicals in a well ventilated room.

Remember to never wipe the surface of any optical items with your hands. Fingerprints can damage the optics. Do not disassemble objectives or eyepieces in an attempt to clean them.

For the best results, use the VWR cleaning kit (see catalog number below).

If you need to send the microscope to manufacturer for maintenance, please use the original packaging.

User replaceable accessories and spare parts

DESCRIPTION	QUANTITY	CAT. NO.
Eyepiece EWF10X/22MM	2	630-2163
WF10x/22mm micrometer eyepiece (10mm, 0.1mm div.)	1	630-1968
Slide micrometric 26x76mm 1mm	1	630-1650
Objective IOS LWD PLAN achromatic 4x/0,10 (w.d. 22mm)	1	630-2164
Objective IOS LWD PLAN achromatic for phase contrast 10x/0,25 (w.d. 7,94mm)	1	630-2165
Objective IOS LWD PLAN achromatic for phase contrast 20x/0,40 (w.d. 7,66mm)	1	630-2166
Objective IOS LWD PLAN achromatic for phase contrast 40x/0,60 (w.d.3,71mm)	1	630-2167
Objective IOS LWD PLAN achromatic 60x/0,70 (w.d. 2,50mm)	1	630-2168
Cut-off filter (infrared)	1	630-2169
Photo adapter for REFLEX camera with FULL FRAME sensor	1	630-2170
1/3" CCD camera adapter	1	630-1891
1/2" CCD camera adapter	1	630-1892
Universal adapter for cameras (only models with eyepiece adapter, 23mm)	1	630-2172
Dust cover type 7	1	630-1769
Phototube adapter for APS-C sensor	1	630-1645
Set for cleaning and maintenance of microscopes	1	630-2081
Immersion oil 8ml	1	630-1808

Technical service

Web Resources

Visit the VWR's website at **www.vwr.com** for:

- Complete technical service contact information
- Access to VWR's Online Catalogue, and information about accessories and related products
- Additional product information and special offers

Contact us For information or technical assistance contact your local VWR representative or visit. **www.vwr.com**.

Warranty

VWR International warrants that this product will be free from defects in material and workmanship for a period of five (5) years from date of delivery. If a defect is present, VWR will, at its option and cost, repair, replace, or refund the purchase price of this product to the customer, provided it is returned during the warranty period. This warranty does not apply if the product has been damaged by accident, abuse, misuse, or misapplication, or from ordinary wear and tear. If the required maintenance and inspection services are not performed according to the manuals and any local regulations, such warranty turns invalid, except to the extent, the defect of the product is not due to such non-performance.

Items being returned must be insured by the customer against possible damage or loss. This warranty shall be limited to the aforementioned remedies. IT IS EXPRESSLY AGREED THAT THIS WARRANTY WILL BE IN LIEU OF ALL WARRANTIES OF FITNESS AND IN LIEU OF THE WARRANTY OF MERCHANTABILITY.

Compliance with local laws and regulations

The customer is responsible for applying for and obtaining the necessary regulatory approvals or other authorizations necessary to run or use the Product in its local environment. VWR will not be held liable for any related omission or for not obtaining the required approval or authorization, unless any refusal is due to a defect of the product.

Disposal

This equipment is marked with the crossed out wheeled bin symbol to indicate that this equipment must not be disposed of with unsorted waste.

Instead it is your responsibility to correctly dispose of your equipment the end of its life cycle by handling it over to an authorized facility for separate collection and recycling. It is also your responsibility to decontaminate the equipment in case of biological, chemical and/or radiological contamination, so as to protect from health hazards the persons involved in the disposal and recycling of the equipment.

For more information about where you can drop off your waste equipment, please contact your local dealer from whom you originally purchased this equipment.

By doing so, you will help to conserve natural and environmental resources and you will ensure that your equipment is recycled in a manner that protects human health.

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Ireland / Northern Ireland

VWR International Ltd /
VWR International (Northern Ireland) Ltd
Orion Business Campus
Northwest Business Park
Ballycoolin
Dublin 15
Tel.: 01 88 22 222
Fax: 01 88 22 333
E-mail: sales@ie.vwr.com

Italy

VWR International PBI S.r.l.
Via San Giusto 85
20153 Milano (MI)
Tel.: 02-3320311/02-487791
Fax: 800 152999/02-40090010
E-mail: info@it.vwr.com

The Netherlands

VWR International B.V.
Postbus 8198
1005 AD Amsterdam
Tel.: 020 4808 400
Fax: 020 4808 480
E-mail: info@nl.vwr.com

New Zealand

Global Science - A VWR Company
241 Bush Road
Albany 0632, Auckland
Tel.: 0800 734 100
Fax: 0800 999 002
E-mail: sales@globalscience.co.nz

Norway

VWR International AS
Haavard Martinsens vei 30
0978 Oslo
Tel.: 02290
Fax: 815 00 940
E-mail: info@no.vwr.com

Poland

VWR International Sp. z o.o.
Limbowa 5
80-175 Gdansk
Tel.: 058 32 38 200 do 204
Fax: 058 32 38 205
E-mail: info@pl.vwr.com

Portugal

VWR International -
Material de Laboratório, Lda
Edifício Neopark
Av. Tomás Ribeiro, 43- 3 D
2790-221 Carnaxide
Tel.: 21 3600 770
Fax: 21 3600 798/9
E-mail: info@pt.vwr.com

Singapore

VWR Singapore Pte Ltd
18 Gul Drive
Singapore 629468
Tel.: +65 6505 0760
Fax: +65 6264 3780
E-mail: sales@sg.vwr.com

Spain

VWR International Eurolab S.L.
C/ Tecnología 5-17
A-7 Llinars Park
08450 - Llinars del Vallès
Barcelona
Tel.: 902 222 897
Fax: 902 430 657
E-mail: info@es.vwr.com

Sweden

VWR International AB
Fagerstagatan 18a
163 94 Stockholm
Tel.: 08 621 34 00
Fax: 08 621 34 66
E-mail: kundservice@se.vwr.com

Switzerland

VWR International GmbH
Lerzenstrasse 16/18
8953 Dietikon
Tel.: 044 745 13 13
Fax: 044 745 13 10
E-mail: info@ch.vwr.com

Turkey

Pro Lab Laboratuvar Teknolojileri Ltd.Şti.
a VWR International Company
Orta Mah. Cemal Gürsel Caddesi
Ördekcioglu İşmerkezi No.32/1
34896 Pendik - İstanbul
Tel.: +90216 598 2900
Fax: +90216 598 2907
Email: info@pro-lab.com.tr

UK

VWR International Ltd
Customer Service Centre
Hunter Boulevard - Magna Park
Lutterworth
Leicestershire
LE17 4XN
Tel.: 0800 22 33 44
Fax: 01455 55 85 86
E-mail: uksales@uk.vwr.com

**GO TO VWR.COM FOR THE
LATEST NEWS, SPECIAL OFFERS
AND DETAILS OF YOUR LOCAL
VWR DISTRIBUTOR**

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