

NanoVue™ Plus Spectrophotometer

NanoVue Plus Spectrophotometer (Fig 1) is a stand-alone instrument with a novel drop-and-measure sample plate that can be used with a wide range of chemicals. NanoVue Plus can be used for the accurate and reliable quantitation of nucleic acids and protein samples. The ability to pipette samples of 0.5 to 5 μl directly onto the sample plate eliminates the need for cuvettes or other sample devices. After measurement, the sample can either be easily recovered with a pipette or quickly discarded by wiping the sample plate clean for the next measurement.

NanoVue Plus is pre-programmed with a range of preset methods for the quantitation of nucleic acids and proteins, using UV or dye-intercalation methods (Lowry, Bradford, BCA, and Biuret). In addition, wavelength scanning gives you the flexibility to design your own methods, which can be stored in a personal folder for easy access. The instrument does not require a computer because it has a direct user interface. Fast instrument start-up with automatic self-calibration means NanoVue Plus Spectrophotometer is always ready to use.

NanoVue Plus Spectrophotometer offers:

- **Novel gold-colored hydrophobic sample plate:** Save time and effort by directly pipetting your sample onto the sample plate for measurement. The hydrophobic coating facilitates sample recovery and minimizes the potential for contamination of the sample mechanism thus improving data precision (Fig 2). The plate is compatible with a wide range of chemicals including DMSO and DMF
- **Low sample volumes:** Reduce sample loss and eliminate the need for dilution by using low volumes of 0.5 to 5 μl for sample measurement. The appropriate pathlength for your sample volume is selected automatically unless you opt for manual selection
- **Fast measurements:** Read-time is typically < 5 s per sample and you save more time by not having to wash cuvettes or dilute samples



Fig 1. NanoVue Plus Spectrophotometer uses a novel drop-and-measure sample plate which eliminates the need for cuvettes.

- **Convenience and ease of use:** Large, high-resolution graphical display enables a quick read of relevant results. Calibration curves, kinetics or ratio measurements are all displayed at the touch of a button (Fig 3)
- **Flexible analytical performance:** Full wavelength scan in less than 5 s from 200 to 900 nm with zoom facility, peak identification, and on-peak confirmation. Visualization of nucleic acid scans allows impurities to be detected and this is especially useful with RNA samples
- **Choice of data output:** Print to an integrated printer (optional) or to any suitable PC via a USB, wireless (Bluetooth) connection or SD card
- **Reliability and robust instrumentation:** Press-to-read feature reduces the amount of time the lamp stays switched on. Optics with no moving parts minimizes the incidence of optical misalignment






Fig 2. Sample carry over study using 25 mg/ml BSA with successive additions of a water-protein-water cycle of samples, which was repeated 25 times. The data showed no significant change in absorbance readings during the study for water sample readings or BSA samples readings thus indicating that there was no sample carry over.




Fig 3. NanoVue Plus spectrophotometer screenshots showing: (A) Simple interface for easy selection of applications; (B) Protein methods interface; (C) DNA quantitation results and (D) An example of a wavelength scan.

The key features of NanoVue Plus Spectrophotometer are summarized in Table 1 below.

Table 1. Features of NanoVue Plus Spectrophotometer

FeatureSpecification

Wavelength range	200 to 1100 nm
Wavelength calibration	Automatic upon switch on
Pathlength selection	Automatic with the option for manual selection
Stored methods	90
Spectral bandwidth	5 nm
Wavelength accuracy	± 2 nm
Wavelength reproducibility	± 0.5 nm
Absorbance accuracy	Maximum ± 1% (at 259 nm) at 0.7 to 0.8 A using uracil
Light source	Long-life stabilized Xenon
Detector	Twin CCD Array
Dimensions	260 × 390 × 100 mm
Weight	4.5 kg
Power requirements	100 to 240 VAC ± 10%, 50/60 Hz, 50 VA

Note: Technical specifications were measured after the instrument had warmed up at a constant ambient temperature.

Assays

Nucleic acids

NanoVue Plus Spectrophotometer allows you to measure the concentration and purity of nucleic acids in a variety of units ($\mu\text{g}/\text{ml}$, $\text{ng}/\mu\text{l}$, $\mu\text{g}/\mu\text{l}$, $\text{pmol}/\mu\text{l}$, and pmol) and also, to correct for dilution factors where necessary.

NanoVue Plus displays both individual absorbance values and absorbance ratios (260/280 and 260/230) on the screen, along with the sample concentration value. A graphical display of the result is also available at the touch of a button. The results can be printed immediately with the on-board printer option.

You can examine the UV spectrum (220 to 330 nm) for hybridization, PCR, and sequencing studies or for quantitation of minipreps after isolation by chromatography.

The instrument enables you to characterize oligonucleotide primers by keying in the base sequence (up to 66-mer) to obtain a conversion factor ($\mu\text{g}/\text{ml}$), molecular weight, theoretical absorbance ($\text{AU}/\mu\text{mol}$), and theoretical T_m .

DNA linearity

NanoVue Plus shows excellent linearity with double-stranded DNA (dsDNA) from 5 to 6000 $\mu\text{g}/\text{ml}$ (Fig 4). This large linear dynamic range and high reproducibility means that you can proceed to the next stage of your workflow with a high degree of confidence in the data produced by NanoVue Plus Spectrophotometer.




Fig 4. DNA concentration curve (0.5 to 6000 $\mu\text{g}/\text{ml}$) generated from salmon sperm dsDNA samples measured at 260 nm using NanoVue Plus [$n = 20$ replicates; mean ± 1 standard deviation (SD)].

Proteins

NanoVue Plus Spectrophotometer can be used to determine the concentration of protein samples by a variety of methods including Bradford, BCA, Lowry, Biuret, and direct UV methods with a choice of line fit and the ability to run up to 27 standards (including replicates). These calibration curves can be viewed on the graphical display, printed or stored as a method for future use.

Protein concentration can be determined in the near UV at 280 nm due to absorption by tyrosine, tryptophan, and phenylalanine amino acids. The absorbance at 280 nm varies greatly for different proteins because of their amino acid content therefore, the specific absorption value for a particular protein must be determined.

NanoVue Plus has in-built extinction coefficients for BSA, IgG, and lysozyme thus allowing you to select the most appropriate for your protein sample or you can input your own extinction coefficients. Typical results for BSA and IgG values using the absorbance at 280 nm program are shown in Figures 5 and 6, respectively.




Fig 5. Protein concentration curve (0.01 to 100 mg/ml) generated from BSA protein samples measured at 280 nm using NanoVue Plus (n = 20 replicates; mean \pm 1 SD).




Fig 6. Protein concentration curve (0.025 to 50 mg/ml) generated from bovine IgG samples measured at 280 nm using NanoVue Plus protein absorbance 280 program with IgG preset option (n = 20 replicates; mean \pm 1 SD).

Background correction

NanoVue Plus provides increased confidence in your results by allowing you to use the background correction feature at 320 nm to compensate for the effects of background absorbance caused by turbidity, stray particulates, and high-absorbance buffer solutions. This feature of the instrument is critical because very small samples are susceptible to stray particulates. The results are displayed on the screen with other absorbance results for quick evaluation. You have the flexibility to turn the background correction function on or off.

CyDye™

You can also measure purity, yields, and brightness of fluorescently labeled in-situ hybridization probes. CyDye DNA function holds spectral data for 19 commonly used dyes, including Cy™2, Cy3, Cy3B, Cy5, and Amersham™ Hyper™5 and both sample concentration and frequency of incorporation measurements are presented.

Calibration

An accurate pathlength distance is important for the correct functioning of all low volume spectrophotometers. We recommend a regular pathlength check as a matter of good laboratory practice. The pathlength calibration fluid supplied with the instrument allows for simple checking of the instrument pathlength whenever convenient for added confidence in performance. Should re-calibration be necessary, you can easily and quickly recalibrate NanoVue Plus in your laboratory without the inconvenience and cost of sending the instrument away.

Drop. Measure. Done.

NanoVue Plus provides intelligent performance across all spectrophotometry applications (Fig 7).




Fig 7. Operation of NanoVue Plus Spectrophotometer.

Summary

NanoVue Plus Spectrophotometer is an easy-to-use and reliable instrument for measuring nucleic acid and protein samples. Samples of 0.5 to 5 µl can be pipetted directly onto a novel sample plate for measurement, and then simply recovered using a pipette. In certain cases where sample recovery is not necessary, the sample plate can be quickly and easily wiped clean. This novel feature eliminates the need for cuvettes, capillaries or other sample devices—just drop, measure, and you're done.

For local office contact information, visit
www.gelifesciences.com/spectros

GE Healthcare UK Limited
Amersham Place
Little Chalfont, Buckinghamshire
HP7 9NA, UK

www.gelifesciences.com



Ordering information

Product

Code No.

NanoVue Plus with integrated printer (US & CA only)28-9569-63

NanoVue Plus with integrated printer
(Worldwide except US & CA)28-9569-66

NanoVue Plus, print via computer (US & CA only)28-9569-62

NanoVue Plus, print via computer
(Worldwide except US & CA)28-9569-65

NanoVue Plus with wireless connectivity (US & CA only)28-9569-64

NanoVue Plus with wireless connectivity
(Worldwide except US & CA)28-9569-67

NanoVue Plus with SD card (US & CA only) 28-9569-61

NanoVue Plus with SD card (Worldwide except US & CA) 28-9569-60

Pathlength calibration kit (Worldwide)28-9244-05

Sample Plate Replacement Kit28-9569-58

Lower Sample Plate Replacement Kit28-9569-59

Related products

Ultrospec 10 Cell density meter	80-2116-30
illustra™ plasmidPrep Mini Spin Kit	28-9042-69
illustra GenomiPhi™ V2 DNA Amplification Kit	25-6600-30
illustra RNAspin Mini Kit	25-0500-71
illustra HotStart PCR Master Mix	25-1500-01
illustra GFX™ PCR DNA and Gel Band Purification Kit	28-9034-70

GE, imagination at work, and GE monogram are trademarks of General Electric Company.

Amersham, Cy, CyDye, GenomiPhi, GFX, HyPer, illustra, and NanoVue are trademarks of GE Healthcare companies.

CyDye: This product or portions thereof is manufactured under an exclusive license from Carnegie Mellon University under US patent number 5,268,486 and equivalent patents in the US and other countries. Cy3-UTP or Cy5-UTP, Cy3.5-dCTP or Cy5.5-dCTP, Cy3-CTP or Cy5-CTP: These products are manufactured for GE Healthcare UK Limited by Perkin Elmer Life Sciences under US patent numbers 5047519 and 5151507. The cyanine dyes in the product are manufactured under an exclusive license from Carnegie Mellon University under US patent numbers 5,268,486 and equivalent patents in the US and other countries. The purchase of CyDye products includes a limited license to use the CyDye products for internal research and development but not for any commercial purposes. A license to use the CyDye products for commercial purposes is subject to a separate license agreement with GE Healthcare. Commercial use shall include:

1. Sale, lease, license or other transfer of the material or any material derived or produced from it.
2. Sale, lease, license or other grant of rights to use this material or any material derived or produced from it.
3. Use of this material to perform services for a fee for third parties, including contract research and drug screening.

If you require a commercial license to use this material and do not have one, return this material unopened to GE Healthcare Bio-Sciences AB, Björkgatan 30, SE-751 84 Uppsala, Sweden and any money paid for the material will be refunded.

All third party trademarks are the property of their respective owners.

© 2008 - 2010 General Electric Company—All rights reserved
First published Jan. 2008

All goods and services are sold subject to the terms and conditions of sale of the company within GE Healthcare which supplies them. A copy of these terms and conditions is available on request. Contact your local GE Healthcare representative for the most current information.

GE Healthcare Europe GmbH
Munzinger Strasse 5, D-79111 Freiburg,
Germany

GE Healthcare Bio-Sciences Corp
800 Centennial Avenue, P.O. Box 1327, Piscataway, NJ 08855-1327
USA

GE Healthcare, Japan Corporation
Sanken Bldg., 3-25-1, Hyakunincho, Shinjuku-ku, Tokyo 169-0073
Japan

GE Healthcare Bio-Sciences AB
Björkgatan 30
751 84 Uppsala
Sweden



Want to know more?

If you have questions about this item, please give us a call on 01257 270 433.

Richmond Scientific are a family run business in Lancashire. We're a small team, and we're always just a phone call away. With over 30 years in the business there's not much we haven't seen before, and we're always happy to chat.

Have kit to sell?

Call us on 01257 270 433

Email: harry@richmondscientific.com

Visit: www.richmondscientific.com

[Visit Website](#)

[Sell used kit](#)

[Our newsletter](#)