

Operating Manual

- BD (E3.1) | Incubators Avantgarde.Line with natural convection
- BF (E3.1) | Incubators Avantgarde.Line with forced convection
- ED (E3.1) | Drying and heating ovens Avantgarde.Line with natural convection
- FD (E3.1) | Drying and heating ovens Avantgarde.Line with forced convection
- FED (E3.1) | Drying and heating ovens Avantgarde.Line with forced convection and enhanced timer functions

with microprocessor temperature controller

Model	Model version	Art. No.	Model	Model version	Art. No.
BD 56	BD056-230V	9010/ 9110-0323	ED 260	ED260-230V	9010/ 9110-0339
50 00	BD056UL-120V	9010/ 9110-0324	ED 200	ED260UL-240V	9010/ 9110-0340
BD 115	BD115-230V	9010/ 9110-0325	ED 720	ED720-400V	9010/ 9110-0341
115	BD115UL-120V	9010/ 9110-0326	FD 56	FD056-230V	9010/ 9110-0303
BD 260	BD260-230V	9010/ 9110-0329	FD 50	FD056UL-120V	9010/ 9110-0304
BD 200	BD260UL-120V	9010/ 9110-0330	FD 115	FD115-230V	9010/ 9110-0305
BD 720	BD720-230V	9010/ 9110-0331	FDTIS	FD115UL-120V	9010/ 9110-0306
BD 720	BD720UL-240V	9010/ 9110-0332	FD 260	FD260-230V	9010/ 9110-0309
BF 56	BF056-230V	9010/ 9110-0313	FD 200	FD260UL-240V	9010/ 9110-0310
DF 30	BF056UL-120V	9010/ 9110-0314	FD 720	FD720-400V	9010/ 9110-0311
	BF115-230V	9010/ 9110-0315		FED056-230V	9010/ 9110-0295
BF 115	BF115UL-120V	9010/ 9110-0316	FED 56	FED056UL-120V	9010/ 9110-0296
BF 260	BF260-230V	9010/ 9110-0319		FED115-230V	9010/ 9110-0293
DF 200	BF260UL-120V	9010/ 9110-0320	FED 115	FED115UL-120V	9010/ 9110-0294
BF 720	BF720-230V	9010/ 9110-0321	FED 260	FED260-230V	9010/ 9110-0299
DF 720	BF720UL-240V	9010/ 9110-0322	FED 200	FED260UL-240V	9010/ 9110-0300
	ED056-230V	9010/ 9110-0333		FED720-400V	9010/ 9110-0301
ED 56	ED056UL-120V	9010/ 9110-0334	FED 720	FED720UL-208V	9010/ 9110-0302
ED 115	ED115-230V	9010/ 9110-0335			
	ED115UL-120V	9010/ 9110-0336			

BINDER GmbH

- ► Address: Post office box 102, 78502 Tuttlingen, Germany ► Phone: +49 7462 2005 0
- ► Fax: +49 7462 2005 100 ► Internet: http://www.binder-world.com
- ► E-mail: info@binder-world.com ► Service Hotline: +49 7462 2005 555
- ► Service Fax: +49 7462 2005 93 555 ► Service E-Mail: customerservice@binder-world.com
- ► Service Hotline USA: +1 866 885 9794 or +1 631 224 4340 x3
- ► Service Hotline Asia Pacific: +852 390 705 04 or +852 390 705 03
- Service Hotline Russia and CIS: +7 495 988 15 16

Content

1.	SAFETY	5
1.1	Personnel Qualification	5
1.2	Operating manual	
1.3	Legal considerations	
1.3 1.4	3.1 Intellectual property Structure of the safety instructions	
1.4	•	
	1.2 Safety alert symbol	
	 4.3 Pictograms 4.4 Word message panel structure 	
1.5	Localization / position of safety labels on the chamber	
1.6	Type plate	9
1.7	General safety instructions on installing and operating the chambers	
1.8 1.9	Intended use Foreseeable Misuse	
1.10	Residual Risks	
2.	CHAMBER DESCRIPTION	
2.1	Chamber overview	
2.1	Triangular instrument panel	
2.3	Main power switch ED, FD, FED 720	
3.	COMPLETENESS OF DELIVERY, TRANSPORTATION, STORAGE, AND INSTALLATION	18
0.4		
3.1 3.2	Unpacking, and checking equipment and completeness of delivery Guidelines for safe lifting and transportation	
3.3	Storage	19
3.4	Location of installation and ambient conditions	19
4.	INSTALLATION	21
4.1	Mounting the tilt protection holders (chambers with window)	
4.2	Electrical connection	
4.3 4.4	Connection to an exhaust/ventilation system (optional) Inserting the racks	
5.	START UP	
5.1	Behavior when opening the door	-
5.2	Performance during and after power failures	
5.3	Loading	
6.	OVERVIEW AND GENERAL SETTINGS ON THE R4 CONTROLLER	26
6.1	Controller overview	
6.2	Normal display	
6.3 6.4	Setting the menu language	
6.5	Selecting the temperature unit	30
6.6	Set-point entry for temperature and fan speed	31
6.6	S.1 Set-point entry for temperature in two-door chambers (ED, FD, FED 720) Adjusting the air flap position	
6.7 6.8	Changing the passwords for user level and general controller functions	33
7.	OVERTEMPERATURE PROTECTION	
7.1	Overtemperature protective device (class 1)	
7.1	Safety controller	

BINDER

7.3 7.4 7.5	Setting the safety controller set-point Alarm message and proceeding in case of an alarm Function check	.39 .39
7.6 8.	Disconnectable audible over-temperature alarm (option)	
o. 8.1 8.2	Selecting the timer function	.41
8.2	2.1 Entry and activation of the timer run-time and fan setting	.42
8.2 8.3	2.2 Turning off the timer function or changing the settings Timer function "Temperature dependent Delayed Off" (BF, FED)	
8.3 8.3		
8.4	Timer function "Delayed On" (BF, FED)	.49
8.4 8.4	5	
8.5	Temperature programming example (BF, FED)	
9.	RAMP FUNCTION	52
9.1	General information	
9.2 9.3	Setting and displaying the ramp function Displaying the effective ramp set-point and changing the target ramp set-point	
9.4	Turning off the ramp function	
10.	DATA RECORDING VIA USB INTERFACE	56
10.1 10.2	Starting data recording Terminating data recording	
11.	NETWORK CONFIGURATION FOR CHAMBERS WITH ETHERNET	
• • •	INTERFACE	58
12.	OPTIONS	61
12.1	APT-COM™ 4 Multi Management Software (option)	.61
12.2	Object temperature display with additional Pt100 temperature sensor (option)	
12.3 12.4	Analog output for temperature (option) HEPA fresh air filter (option for FD, FED)	
12.5	Water protected disconnectable internal socket (option BD, BF)	.62
12.6 12.7	Mostly gas-tight version (option for BF, FD, FED) Inert gas connection with mostly gas-tight version (option for BF, FD, FED)	
13.	CLEANING AND DECONTAMINATION	
13.1	Cleaning	
13.1	Decontamination / chemical disinfection	
14.	MAINTENANCE AND SERVICE, TROUBLESHOOTING, REPAIR, TESTING	68
14.1	General information, personnel qualification	
14.2 14.3	Maintenance intervals, service	
14.3 14.4	Simple troubleshooting Sending the chamber back to BINDER GmbH	
15.	DISPOSAL	71
15.1	Disposal of the transport packing	.71
15.2	Decommissioning	
15.3 15.4	Disposal of the chamber in the Federal Republic of Germany Disposal of the chamber in the member states of the EU except for the Federal Republic of Germa	
15.5	Disposal of the chamber in non-member states of the EU	
16.	TECHNICAL DESCRIPTION	74
16.1 16.2	Factory calibration and adjustment Definition of usable volume	

BINDER

16.9	Over current protection	75 76 78 80 81 83
16.11 16.12	Dimensions size 56 Dimensions size 115	86 88
	Dimensions size 260 Dimensions size 720	
17.	CERTIFICATES AND DECLARATIONS OF CONFORMITY	. 95
17.1 17.2 17.3 17.4 17.5 17.6	EU Declaration of Conformity for BD EU Declaration of Conformity for BF. EU Declaration of Conformity for ED EU Declaration of Conformity for FD EU Declaration of Conformity for FED Certificate for the GS mark of conformity of the "VDE Prüf- und Zertifizierungsinstitut" (Testing a Certification Institute of the Association for Electrical, Electronic and Information Technologies). Certificate of Compliance for the UL Certification Mark from Underwriters Laboratories	97 99 101 103 and 105
18.	CONTAMINATION CLEARANCE CERTIFICATE1	111
18.1 18.2	For chambers located outside the USA and Canada For chambers located in the USA and Canada	

Dear customer,

For the correct operation of the chambers, it is important that you read this operating manual completely and carefully and observe all instructions as indicated. Failure to read, understand and follow the instructions may result in personal injury. It can also lead to damage to the chamber and/or poor equipment performance

1. Safety

1.1 Personnel Qualification


The chamber must only be installed, tested, and started up by personnel qualified for assembly, startup, and operation of the chamber. Qualified personnel are persons whose professional education, knowledge, experience and knowledge of relevant standards allow them to assess, carry out, and identify any potential hazards in the work assigned to them. They must have been trained and instructed, and be authorized, to work on the chamber.

The chamber should only be operated by laboratory personnel especially trained for this purpose and familiar with all precautionary measures required for working in a laboratory. Observe the national regulations on minimum age of laboratory personnel.

1.2 Operating manual

This operating manual is part of the components of delivery. Always keep it handy for reference in the vicinity of the chamber. If selling the unit, hand over the operating manual to the purchaser.

To avoid injuries and damage observe the safety instructions of the operating manual. Failure to follow instructions and safety precautions can lead to significant risks.



read and understood the Operating Manual.

This Operating Manual is supplemented and updated as needed. Always use the most recent version of the Operating Manual. When in doubt, call the BINDER Service Hotline for information on the up-to-date-ness and validity of this Operating Manual.

1.3 Legal considerations

This operating manual is for informational purposes only. It contains information for correct and safe installing, start-up, operation, decommissioning, cleaning and maintenance of the product. Note: the contents and the product described are subject to change without notice.

Understanding and observing the instructions in this operating manual are prerequisites for hazard-free use and safety during operation and maintenance. Images are to provide basic understanding. They may deviate from the actual version of the chamber. The actual scope of delivery can, due to optional or special design, or due to recent technical changes, deviate from the information and illustrations in these instructions this operating manual. In no event shall BINDER be held liable for any damages, direct or incidental arising out of or related to the use of this manual.



This operating manual cannot cover all conceivable applications. If you would like additional information, or if special problems arise that are not sufficiently addressed in this manual, please ask your dealer or contact us directly, e.g. by phone at the number located on page one of this manual

Furthermore, we emphasize that the contents of this operating manual are not part of an earlier or existing agreement, description, or legal relationship, nor do they modify such a relationship. All obligations on the part of BINDER derive from the respective purchase contract, which also contains the entire and exclusively valid statement of warranty administration and the general terms and conditions, as well as the legal regulations valid at the time the contract is concluded. The statements in this manual neither augment nor restrict the contractual warranty provisions.

1.3.1 Intellectual property

Trademark Information: All BINDER trademarks relating to products or service, as well as trade names, logos and product names used on the website, products and documents of BINDER company are trademarks or registered trademarks of BINDER company (including BINDER GmbH, BINDER Inc.) in the U.S. and other countries and communities of states. This includes word marks, position marks, word/figurative marks, design configurations, figurative marks, and design patents.

Patent Information: BINDER products, categories of products, and accessories may be covered by one or more patents and/or utility models in the U.S. and other countries and communities of states. Additional patent applications may also be pending in the U.S. and other countries and communities of states.

Please visit www.binder-world.com for more information.

1.4 Structure of the safety instructions

In this operating manual, the following safety definitions and symbols indicate dangerous situations following the harmonization of ISO 3864-2 and ANSI Z535.6.

1.4.1 Signal word panel

Depending on the probability of serious consequences, potential dangers are identified with a signal word, the corresponding safety color, and if appropriate, the safety alert symbol.





Indicates a potentially hazardous situation which, if not avoided, could result in death or serious (irreversible) injury

Indicates a potentially hazardous situation which, if not avoided, may result in moderate or minor (reversible) injury

NOTICE

Indicates a potentially hazardous situation which, if not avoided, may result in damage to the product and/or its functions or of a property in its proximity.

1.4.2 Safety alert symbol



Use of the safety alert symbol indicates a **risk of injury**. Observe all measures that are marked with the safety alert symbol in order to avoid death or injury.

1.4.3 Pictograms

Warning signs			
Electrical hazard	Hot surface	Explosive atmosphere	Stability hazard
Lifting hazard	Risk of corrosion and / or chemical burns	Suffocation hazard	Harmful substances
Biohazard	Pollution Hazard		
Mandatory action signs	1		
Mandatory regulation	Read operating instructions	Disconnect the power plug	Lift with several persons
Lift with mechanical	Environment protection	Wear protective gloves	Wear safety goggles
assistance			Wear safety goggies
Prohibition signs			
Do NOT touch	Do NOT spray with water		
Information to	be observed in order to ens	ure optimum function of the	e product.

1.4.4 Word message panel structure

Type / cause of hazard.

Possible consequences.


- $\ensuremath{\varnothing}$ Instruction how to avoid the hazard: prohibition.
- > Instruction how to avoid the hazard: mandatory action.


Observe all other notes and information not necessarily emphasized in the same way, in order to avoid disruptions that could result in direct or indirect injury or property damage.


1.5 Localization / position of safety labels on the chamber

The following labels are located on the chamber:

Pictograms	(Warning signs)	Service label
	 Hot surface ED, FD, FED: outer chamber door BD, BF: on the glass door handle On chamber rear next to the exhaust duct 	Service - Hotline International: +49 (0) 7462 / 2005-555 USA Toll Free: +1 866 885 9794 от: +1 631 224 4340 Россия и СНГ: +7 495 98815 17 service@binder-world.com www.binder-world.com
	 Read operating manual UL chamber: outer chamber door BD, BF with optional interior socket: below the interior socket 	







|--|

Replace safety labels that are no longer legible. Contact BINDER Service for these replacements.

1.6 Type plate

The type plate is located on the left-hand side of the chamber, bottom right-hand.

Nominal temp. IP protection Safety device Class	572 °F 20 DIN 12880 2.0	1,30 kW / 5,7 A 230 V / 50 Hz 230 V / 60 Hz 1 N ~			
Art. No. Project No.	9010-0305				
Built	2021	Drying and hear BINDE Im Mittl	R GmbH leren Ösch 5	FED 115	Serial No. 00000000000000
			Tuttlingen / Germany inder-world.com	E3.1	Made in Germany





Figure 3: Type plate (example BD 115-230V optional chamber)

Indications of the type plate (example)

Indication		Information	
BINDER		Manufacturer: BINDER GmbH	
BD 115		Model designation	
Incubator		Chamber name: Incubator	
Drying and heating over	en	Chamber name: Drying and heating oven	
Serial No.	000000000000	Serial No of the chamber	
Built	2021	Year of construction	
Nominal temperature	100 °C 212 °F	Nominal temperature	
IP protection	20	IP type of protection acc. to EN 60529	
Temp. safety device	DIN 12880	Temperature safety device acc. to standard DIN 12880	
Class	3.1	Class of temperature safety device	
Art. No.	9110-0305	Art. no. of the chamber	
Project No.		Optional: Special application acc. to project no.	
1,30 kW		Nominal power	
5,7 A		Nominal current	
230 V / 50 Hz		Nominal voltage ± 10%	
230 V / 60 Hz		at the indicated power frequency	
1 N ~		Current type	
With option internal socket: Nominal power: 0,85 kW		With option internal socket: increased total nominal power	

Symbols on the type plate


Symbol	Valid for	Information
CE	All chambers	CE conformity marking
	All chambers	Electrical and electronic equipment manufactured / placed on the market in the EU after 13 August 2005 and be dis- posed of in separate collection according to Directive 2012/19/EU on waste electrical and electronic equipment (WEEE).
	All chambers	GS mark of conformity of the "VDE Prüf- und Zertifizier- ungsinstitut" (Testing and Certification Institute of the Associ- ation for Electrical, Electronic and Information Technologies
EAC	Not for UL cham- bers	The chamber is certified according to Customs Union Tech- nical Regulation (CU TR) for the Eurasian Economic Union (Russia, Belarus, Armenia, Kazakhstan Kyrgyzstan).
LISTED LABORATORY EQUIPMENT	UL chambers only	 The chamber is certified by Underwriters Laboratories Inc.[®] according to the following standards: UL 61010-1, 3rd Edition, 2012-05, rev. 2015-07 CAN/CSA-C22.2 No. 61010-1, 3rd Edition, 2012-05, rev. 2015-07


1.7 General safety instructions on installing and operating the chambers

With regard to operating the chambers and to the installation location, please observe the local and national regulations relevant for your country (for Germany: DGUV guidelines 213-850 on safe working in laboratories, issued by the employers' liability insurance association).


BINDER GmbH is only responsible for the safety features of the chamber provided skilled electricians or qualified personnel authorized by BINDER perform all maintenance and repair, and if components relating to chamber safety are replaced in the event of failure with original spare parts.

To operate the chamber, use only original BINDER accessories or accessories from third-party suppliers authorized by BINDER. The user is responsible for any risk caused by using unauthorized accessories.





The chambers do not dispose of any measures of explosion protection.



Any solvent contained in the charging material must not be explosive or inflammable. I.e., irrespective of the solvent concentration in the steam room, NO explosive mixture with air must form. The temperature inside the chamber must lie below the flash point or below the sublimation point of the charging material. Familiarize yourself with the physical and chemical properties of the charging material, as well as the contained moisture constituent and its behavior with the addition of heat energy.

Familiarize yourself with any potential health risks caused by the charging material, the contained moisture constituent or by reaction products that may arise during the temperature process. Take adequate measures to exclude such risks prior to putting the chamber into operation.



The chambers were produced in accordance with VDE regulations and were routinely tested in accordance to VDE 0411-1 (IEC 61010-1).

During and shortly after operation, the temperature of the inner surfaces almost equals the set-point. The glass doors and glass door handles (BD, BF), inner chamber, exhaust duct, door window (option), and the door gaskets will become hot during operation.

Danger of burning by touching hot chamber parts during or after operation.
Burns.
Ø Do NOT touch the glass doors, inner surfaces, exhaust duct, door window, access ports, door gaskets, or the charging material during operation.
\oslash BF, FD, FED: Do NOT place the power cable over the door gap when the chamber is hot after operation.

1.8 Intended use

$\sim \circ$

Observing the instructions in this operating manual and conducting regular maintenance work (chap. 14.2) is part of the intended use.

Any use of the chambers that does not comply with the requirements specified in this Operating Manual shall be considered improper use.

Other applications than those described in this chapter are not approved.

Use

The chambers are suitable for exact tempering of harmless materials and for drying and heat treatment of solid or pulverized charging material, as well as bulk material, using the supply of heat. They can be used to dry e.g. glassware, and for warm storage of liquids in containers.

DO NOT use the device for drying processes that release so large amounts of water vapor that condensation occurs.


Because of their precise temperature accuracy the incubators BD and BF are especially useful for incubation of cultures at a standard temperature of 37 $^{\circ}$ C / 98.6 $^{\circ}$ F.

Requirements for the chamber load

Any solvent must not be explosive or flammable. A mixture of any component of the charging material with air must NOT be explosive. The operating temperature must lie below the flash point or below the sublimation point of the charging material. Any component of the charging material must NOT be able to release toxic gases.

The loading material shall not contain any corrosive ingredients that may damage the machine components made of stainless steel, aluminum, and copper. Such ingredients include in particular acids and halides. Any corrosive damage caused by such ingredients is excluded from liability by BINDER GmbH.

The chamber does not dispose of any measures of explosion protection.



 \varnothing Do NOT introduce any substance which could lead to release of toxic gases.

Contamination of the chamber by toxic, infectious or radioactive substances must be prevented

Danger of intoxication and infection through contamination of the chamber with toxic, infectious or radioactive substances. Damages to health.
 Protect the interior of the chamber from contamination by toxic, infectious or radioactive substances.
Take suitable protective measures when introducing and removing toxic, infectious or radioactive material

In case of foreseeable use of the chamber there is no risk for the user through the integration of the chamber into systems or by special environmental or operating conditions in the sense of EN 61010-1:2010. For this, the intended use of the chamber and all its connections must be observed.

Connect only external devices to the chamber interfaces Ethernet (regular with FED, optional with BD, BF, ED, FD) and USB which are compliant with the standards EN 61010-1:2010 or EN 60950-1:2006 mod.

Medical devices

The chambers are not classified as medical devices as defined by the Medical Device Directive 93/42/EEC and Regulation (EU) No 2017/745.



Due to the special demands of the Medical Device Directive (MDD), these chambers are not qualified to perform sterilization of medical devices as defined by the directive 93/42/EEC.

Personnel Requirements

Only trained personnel with knowledge of the Operating Manual can set up and install the chamber, start it up, operate, clean, and take it out of operation. Service and repairs call for further technical requirements (e.g. electrical know-how), as well as knowledge of the service manual.

Installation site requirements

The chambers are designed for setting up inside a building (indoor use).

The requirements described in the Operating Manual for installation site and ambient conditions (chap. 3.4) must be met.



WARNING: If customer should use a BINDER chamber running in non-supervised continuous operation, we strongly recommend in case of inclusion of irrecoverable specimen or samples to split such specimen or samples and store them in at least two chambers, if this is feasible.

1.9 Foreseeable Misuse

Other applications than those described in chap. 1.8 are not approved.

This expressly includes the following misuses (the list is not exhaustive), which pose risks despite the inherently safe construction and existing technical safety equipment:

- Non-observance of Operating Manual
- Non-observance of information and warnings on the chamber (e.g. control unit messages, safety identifiers, warning signals)
- Installation, startup, operation, maintenance and repair by untrained, insufficiently qualified, or unauthorized personnel
- Missed or delayed maintenance and testing
- Non-observance of traces of wear and tear
- Insertion of materials excluded or not permitted by this Operating Manual.
- Non-compliance with the admissible parameters for processing the respective material.
- Installation, testing, service or repair in the presence of solvents
- Installation of replacement parts and use of accessories and operating resources not specified and authorized by the manufacturer
- Bypassing or changing protective systems, operation of the chamber without the designated protective systems
- Non-observance of messages regarding cleaning and disinfection of the chamber.
- Spilling water or cleaning agent on the chamber, water penetrating into the chamber during operation, cleaning or maintenance
- · Cleaning activity while the chamber is turned on
- Operation of the chamber with a damaged housing or damaged power cord
- · Continued operation of the chamber during an obvious malfunction



- Insertion of objects, particularly metallic objects, in louvers or other openings or slots on the chamber
- Human error (e.g. insufficient experience, qualification, stress, exhaustion, laziness)

To prevent these and other risks from incorrect operation, it is recommended the operator issue operating instructions and standard operating procedures (SOPs).

1.10 Residual Risks

The unavoidable design features of a chamber, as well as its proper field of application, can also pose risks, even during correct operation. These residual risks include hazards which, despite the inherently safe design, existing technical protective equipment, safety precautions and supplementary protective measures, cannot be ruled out.

Messages on the chamber and in the Operating Manual warn of residual risks. The consequences of these residual risks and the measures required to prevent them are listed in the Operating Manual. Moreover, the operator must take measures to minimize hazards from unavoidable residual risks. This includes, in particular, issuing operating instructions.

The following list summarizes the hazards against which this Operating Manual and the Service Manual warn, and specifies protective measures at the appropriate spots:

Unpacking, Transport, Installation

- Sliding or tilting the chamber
- Setup of the chamber in unauthorized areas
- Installation of a damaged chamber
- Installation of a chamber with damaged power cord
- Inappropriate site of installation
- Missing protective conductor connection

Normal operation

- Assembly errors
- Contact with hot surfaces on the housing
- Contact with hot surfaces in the interior and inside of doors
- Emission of non-ionizing radiation from electrical operating resources
- Contact with live parts in normal state

Cleaning and Decontamination

- Penetration of water into the chamber
- Inappropriate cleaning and decontamination agents
- Enclosure of persons in the interior

Malfunction and Damage

- Continued operation of the chamber during an obvious malfunction or outage of the heating system
- Contact with live parts during error status
- Operation of a unit with damaged power cord

Maintenance

- Maintenance work on live parts.
- Execution of maintenance work by untrained/insufficiently qualified personnel
- Electrical safety analysis during annual maintenance not performed

Trouble-shooting and Repairs

- Non-observance of warning messages in the Service Manual
- Trouble-shooting of live parts without specified safety measures
- Absence of a plausibility check to rule out erroneous inscription of electrical components
- Performance of repair work by untrained/insufficiently qualified personnel
- Inappropriate repairs which do not meet the quality standard specified by BINDER
- Use of replacement parts other than BINDER original replacement parts
- Electrical safety analysis not performed after repairs

2. Chamber description

BINDER incubators BD and BF and drying and heating ovens ED, FD and FED are equipped with an electronic PID-controller with digital display.

The incubators BD and BF indicate the temperature with an accuracy of a tenth of a degree. The drying and heating ovens ED, FD and FED indicate the temperature with an accuracy of one degree.

All chambers are heated electrically. Incubators BD and drying and heating ovens ED are ventilated naturally. Incubators BF and drying and heating ovens FD and FED are ventilated by fan-assisted, forced-air circulation.

The concept of air conduction guarantees high level of spatial and time-based temperature precision, thanks to the direct and distributed air circulation into the interior. With BF, FD and FED, the fan supports exact attainment and maintenance of the desired temperature accuracy.

The chambers are regularly equipped with an overtemperature safety device class 1 acc. to DIN12880:2007 and with an overtemperature safety controller (overtemperature temperature safety device class 2 or class 3.1 acc. to DIN12880:2007), see chap. 7).

Material: The inner chamber and the inside of the doors are made of stainless steel V2A (German material no. 1.4301, US equivalent AISI 304 and material no. 1.4016, US equivalent AISI 430). Drying and heating ovens ED, FD and FED: When operating the chambers at temperatures above 150 °C / 302 °F, the impact of the oxygen in the air may cause discoloration of the metallic surfaces (yellowish-brown or blue) by natural oxidation processes. These colorations are harmless and will in no way impair the function or quality of the chamber. The housing is RAL 7035 powder-coated. All corners and edges are also completely coated.

All chamber functions are easy and comfortable to use thanks to their clear arrangement. Major features are easy cleaning of all chamber parts and avoidance of undesired contamination.

The chambers are regularly (FED) or optionally equipped with an Ethernet interface for computer communication, e.g. via the APT-COM[™] 4 Multi Management Software (option, chap. 12.1) and with a USB interface to read out the measured values in real time.

The models size 720 are equipped with four castors. Both front castors can be locked by brakes.

Temperature ranges see technical data (chap. 16.4 - 16.8).

2.1 Chamber overview




Figure 4: Overview, closed chamber (chamber with single door)




Figure 5: Overview, open chamber with glass door (chamber with single door) (BD, BF)

- (1) Triangular instrument panel with controller R4 and USB interface
- (2) Door handle
- (3) Outer door
- (4) Glass door handle (BD and BF)
- (5) Glass door (BD and BF)
- (6) Rack

2.2 Triangular instrument panel




Figure 6: Triangular instrument panel

- (7) Controller display
- (8) Functional controller buttons
- (9) USB interface
- (10) Switch for Interior lighting (option)

2.3 Main power switch ED, FD, FED 720

The chambers ED, FD, FED size 720 are equipped with a main power switch located on the chamber rear.



Off



On

Figure 7: Main power switch on the chamber rear


3. Completeness of delivery, transportation, storage, and installation

3.1 Unpacking, and checking equipment and completeness of delivery

After unpacking, please check the chamber and its optional accessories, if any, based on the delivery receipt for completeness and for transportation damage. Inform the carrier immediately if transportation damage has occurred.

The final tests of the manufacturer may have caused traces of the racks on the inner surfaces. This has no impact on the function and performance of the chamber.

Please remove any transportation protection devices and adhesives in/on the chamber and on the doors and take out the operating manuals and accessory equipment.



If you need to return the chamber, please use the original packing and observe the guidelines for safe lifting and transportation (chap. 3.2).

For disposal of the transport packing, see chap. 15.1.

Note on second-hand chambers (Ex-Demo-Units):


Second-hand chambers are chambers that have been used for a short time for tests or exhibitions. They are thoroughly tested before resale. BINDER ensures that the chamber is technically sound and will work flawlessly.

Second-hand chambers are marked with a sticker on the chamber door. Please remove the sticker before commissioning the chamber.

3.2 Guidelines for safe lifting and transportation

The front castors of chambers size 720 can be blocked by brakes. Please move the chambers with castors only when empty and on an even surface, otherwise the castors may be damaged. After operation, please observe the guidelines for temporarily decommissioning the chamber (chap. 15.2).

BINDER



• Permissible ambient temperature range during transport: -10 °C to +60 °C / 14 °F to 140 °F.

You can order transport packing and pallets for transportation purposes from BINDER Service.

3.3 Storage


Intermediate storage of the chamber is possible in a closed and dry room. Observe the guidelines for temporary decommissioning (chap. 15.2).

- Permissible ambient temperature range during storage: -10 °C to +60 °C / 14 °F to 140 °F.
- Permissible ambient humidity: max. 70 % r.H., non-condensing


When after storage in a cold location you transfer the chamber to its warmer installation site, condensation may form. Before start-up, wait at least one hour until the chamber has attained ambient temperature and is completely dry.

3.4 Location of installation and ambient conditions

Set up the chamber on an even and non-flammable surface, free from vibration and in a well-ventilated, dry location and align it using a spirit level. The site of installation must be capable of supporting the chamber's weight (see technical data, chap. 16.4 to 16.7). The chambers are designed for setting up inside a building (indoor use).



Do not install or operate the chamber in potentially explosive areas.



Ambient conditions

• Permissible ambient temperature range during operation: +18 °C up to +40 °C / 64.4 °F to 104 °F. At elevated ambient temperature values, fluctuations in temperature can occur.



The ambient temperature should not be substantially higher than the indicated ambient temperature of +25 °C / 77 °*F* to which the specified technical data relate. For other ambient conditions, deviations from the indicated data are possible.

- Permissible ambient humidity: 70 % r.H. max., non-condensing.
- Installation height: max. 2000 m / 6562 ft. above sea level.

Minimum distances

When placing several chambers of the same size side by side, maintain a minimum distance of 250 mm / *9.84 in* between each chamber. Wall distances: rear 160 mm / *6.30 in*, sides 100 mm / *3.94 in*. Spacing above the chamber of at least 100 mm / *3.94 in* must also be accounted for.

- When placing several chambers of the same size side by side, maintain a minimum distance of 250 mm / 9.84 *in* between each chamber.
- Wall distances: rear 160 mm / 6.30 in, sides 100 mm / 3.94 in.
- Spacing above the chamber of at least 100 mm / 3.9 in must also be accounted for.

Two devices up to size 115 can be stacked on top of each other. For this purpose place rubber pads under all feet of the upper chamber to prevent the device from slipping.



Chambers sizes 260 and 720 must NOT be stacked.



Other requirements

To completely separate the chamber from the power supply, you must disconnect the power plug. Install the chamber in a way that the power plug is easily accessible and can be easily pulled in case of danger. Do not conduct the power cable above the exhaust duct.

4. Installation

4.1 Mounting the tilt protection holders (chambers with window)

For chambers equipped with the option "door with window" it is recommended to install the supplied tilt protection.

Scope of delivery of tilt protection kit (Art.no. 8009-0870):

- 2 screws
- 2 tilt protection holders

Preparing the tilt protection holders

• Depending on the desired wall distance, you can bend the tilt protection holders accordingly.




Figure 8: Variable length of the tilt protection holder depending on the bend

Installation on the chamber

- Plug the two tilt protection holders each with the tab on the provided spot on the edge of the rear panel. The screw holes in the rear wall and in the tilt protection holders must align.
- Fix the tilt protection holders each with one of the supplied screws on the chamber rear wall.

Wall mounting

• Then fix both tilt protection holders on the wall, each with 2 screws Ø 6mm suitable for the wall (B)

4.2 Electrical connection

The chambers are supplied ready for connection and.

The chambers ED720 / FD720 / FED720 / FED720UL come with a fixed power connection cable of at least 1800 mm / 70.87 *in* in length. The other chambers come with an IEC connector plug.

Model	Power plug / power cable	Nominal voltage +/- 10% at the indicated power frequency	Current type	Chamber fuse
BD056-230V BF056-230V	Grounded plug	230 V at 50 Hz 230 V at 60 Hz	1N~	6,3 A
ED056-230V FD056-230V FED056-230V	Grounded plug	230 V at 50 Hz 230 V at 60 Hz	1N~	6,3 A
BD115-230V BF115-230V	Grounded plug	230 V at 50 Hz 230 V at 60 Hz	1N~	6,3 A
ED115-230V FD115-230V FED115-230V	Grounded plug	230 V at 50 Hz 230 V at 60 Hz	1N~	6,3 A
BD260-230V BF260-230V	Grounded plug	230 V at 50 Hz 230 V at 60 Hz	1N~	8,0 A
ED260-230V FD260-230V FED260-230V	Grounded plug	230 V at 50 Hz 230 V at 60 Hz	1N~	12,5 A
BD720-230V BF720-230V	Grounded plug	230 V at 50 Hz 230 V at 60 Hz	1N~	12,5 A
ED720-400V FD720-400V FED720-400V	Grounded plug	400 V at 50 Hz 400 V at 60 Hz	3N~	
BD056UL-120V BF056UL-120V	NEMA 5-15P SJT 14AWG*3C	120 V at 50 Hz 120 V at 60 Hz	1N~	12,5 A
ED056UL-120V FD056UL-120V FED056UL-120V	NEMA 5-15P SJT 14AWG*3C	120 V at 50 Hz 120 V at 60 Hz	1N~	12,5 A
BD115UL-120V BF115UL-120V	NEMA 5-15P SJT 14AWG*3C	120 V at 50 Hz 120 V at 60 Hz	1N~	12,5 A
ED115UL-120V FD115UL-120V FD115UL-120V	NEMA 5-15P SJT 14AWG*3C	120 V at 50 Hz 120 V at 60 Hz	1N~	12,5 A
BD260UL-120V BF260UL-120V	NEMA 5-15P SJT 12AWG*3C	120 V at 50 Hz 120 V at 60 Hz	1N~	12,5 A
ED260UL-240V FD260UL-240V FED260UL-240V	NEMA 6-20P SJT 12AWG*3C	240 V at 50 Hz 240 V at 60 Hz	2~	
BD720UL-240V BF720UL-240V	NEMA 6-20P SJT 12AWG*3C	240 V at 50 Hz 240 V at 60 Hz	2~	
FED720UL-208V	NEMA L21-20P	208 V at 50 Hz 208 V at 60 Hz	3N~	



• The domestic socket must also provide a protective conductor. Make sure that the connection of the protective conductor of the domestic installations to the chamber's protective conductor meets the latest technology. The protective conductors of the socket and plug must be compatible!



Electrical hazard due to missing protective conductor connection. Deadly electric shock.

- Make sure that the chamber's power plug and the power socket match and securely connect the electrical protective conductors of the chamber and the house installation.
- Only use original connection cables from BINDER according to the above specification.

UL chambers: Use only a UL Listed Power supply cord (UL category ELBZ) according to the above specification. For outside USA use a certified power supply cord according to national requirements.

• Prior to connection and start-up, check the power supply voltage. Compare the values to the specified data located on the chamber's type plate (left-hand side of the chamber, chap. 1.6).



NOTICE

Danger of incorrect power supply voltage due to improper connection. Damage to the chamber.

- > Check the power supply voltage before connection and start-up.
- Compare the power supply voltage with the data indicated on the type plate.
- When connecting, please observe the regulations specified by the local electricity supply company as well as the local or national electrical regulations (VDE directives for Germany).
- Observe a sufficient current protection according to the number of devices that you want to operate. We recommend the use of a residual current circuit breaker.
- BF, FD, FED: Do not place the power cable over the door gap when the chamber is hot after operation.
- Pollution degree (acc. to IEC 61010-1): 2
- Over-voltage category (acc. to IEC 61010-1): II

See also electrical data (chap. 16.4 to 16.7).



To completely separate the chamber from the power supply, you must disconnect the power plug. Install the chamber in a way that the power plug is easily accessible and can be easily pulled in case of danger.

4.3 Connection to an exhaust/ventilation system (optional)


Active suction from the chamber must only be effected together with external air. Therefore, the chamber's exhaust air duct shall not be immediately connected to an active exhaust system.

When connecting to an active exhaust system, proceed as follows:

• Perforate the connecting piece between the exhaust air duct and the exhaust system.

Or

• Use an exhaust air funnel placed in a distance of 3-5 cm / 1 to 2 in from the exhaust air duct. The funnel's opening must be at least twice as large as the diameter of the exhaust air duct.



The exhaust duct on the chamber rear will become hot during operation



Danger of burning by touching hot chamber parts during operation. Burns.

CAUTION


 \varnothing Do NOT touch the exhaust duct during operation.

4.4 Inserting the racks


Observe the correct orientation of the racks:

Standard rack: The lateral brackets must be above the rack surface when inserting the rack.

Optional heavy load rack: The lateral brackets must be below the rack surface when inserting the rack.




Standard rack



Optional heavy load rack

Figure 9: Correct orientation when inserting the racks



5. Start up

Insert the plug into a suitable socket (chap. 4.2).

BF, FD, FED size 720: Turn the chamber at the main power switch (chap. 2.3).



If there is no other indication on the controller than the standby symbol, press the standby button until the display lights up.

The controller now shows normal display (chap. 6.2). If a timer function was active prior to turning off the chamber, it is shown in the controller display.



Warming chambers may release odors in the first few days after commissioning. This is not a quality defect. To reduce odors quickly we recommend heating up the chamber to its nominal temperature for one day and in a well-ventilated location.

5.1 Behavior when opening the door

BD, ED: Depending on the temperature, heating performance may be adapted when opening the door.

BF, FD, FED: When opening the door, heating and fan turn off as long as the door remains open.

5.2 Performance during and after power failures

During a power failure, all controller functions are shut down.

After the power returns, all functions return to the same status the chamber had before power failure. The controller continues to function in the original operating mode it was in previously before the power failure occurred.

• Performance after power failure in Standby mode

Control is deactivated

• Performance after power failure that occurred while equilibrating to the set-point:

The set-points are immediately resumed.


• Performance after power failure during timer operation

The program is resumed at the point where the interruption occurred with the latest set-points reached during the program run.

5.3 Loading

When loading the chamber, observe the maximum permissible load per rack and the maximum permissible total load (see tecnical data, chap. 16.4 to 16.8).

Observe the correct orientation of the racks (chap. 4.4).




6. Overview and general settings on the R4 controller

6.1 Controller overview



Buttons on t	he controller
$\triangle \nabla$	The arrow buttons serve to navigate and to enter the values
ок	The OK button serves to select the parameters and to confirm the entered values
IJ	The Back button serves to reach the preceding level
С	If the Standby button is pressed down for approx. 3 seconds, the display changes to standby mode. To activate the display, press down the standby button again for approx. 3 seconds
	Display in standby mode with standby symbol



6.2 Normal display

Setp 25%	Normal display with chambers without fan (BD, ED) or with fixed fan speed (FD)
Fan 40%	Normal display with chambers with adjustable fan speed (BF, FED)

	From Normal display
without fan 5x with fan 6x	with the arrow-up button to the user menu
ок	Confirm with OK.
Password	Enter the password (factory setting: 00 00) and confirm each entry with OK.
4 x 🛆	with the arrow-up button to the language setting menu.
Lansuase	The current menu language is shown.
ок	Press OK to select the menu language.
	The setting flashes.
	Select the setting with the arrow buttons
ок	and confirm with OK.
2x 🖻	Back to Normal display.

6.3 Setting the menu language

There are the following options:



6.4 Setting date and time

	From Normal display
without fan 5x with fan 6x	with the arrow-up button to the user menu
ок	Confirm with OK.
Password	Enter the password (factory setting: 00 00) and confirm each entry with OK.
14 Jul2015	The current date is shown.

ок	Press OK to set the year .
	The setting flashes.
	Enter the year with the arrow buttons (any setting)
ок	and confirm with OK.

ок	Press OK to set the month .
Month	The setting flashes.
	Enter the month with arrow buttons (1 to 12)
ок	and confirm with OK.



Without the optional real time clock, these settings must be repeated when the power supply is interrupted..

ок	Press OK to set the day .
Day	The setting flashes.
	Enter the day with arrow buttons (1 to 31)
ок	and confirm with OK.

ок	Press OK to set the hour .
Hours	The setting flashes.
	Enter the hour with arrow buttons (0 to 23)
ок	and confirm with OK.

ок	Press OK to set the minute .
Minute	The setting flashes.
	Enter the minute with arrow buttons (0 to 59)
ок	and confirm with OK.

2x 🗾 Back to Normal display.

6.5 Selecting the temperature unit

You can chose between degrees Celsius °C and degrees Fahrenheit °F.

If the unit is changed, the temperature set-point and limits are converted accordingly.

Also when specifying the ramp function (see chap. 9) this setting is accordingly taken as the basis.

}	C = degrees Celsius F= degrees Fahrenheit	0 °C = 31°F	Conversion:
29	F= degrees Fahrenheit	100 °C = 212°F	[Value in °F] = [Value in °C] * 1.8 + 32

	From Normal display
without fan 5x with fan 6x	with the arrow-up button to the user menu
ок	Confirm with OK.
Password	Enter the password (factory setting: 00 00) and confirm each entry with OK.
Δ	With the arrow-up button to the temperature unit selection menu.
Unit	The current temperature unit is shown.
ок	Press OK to select the temperature unit.
Unit	The setting flashes.
	Select the setting with arrow buttons
ок	and confirm with OK.
2x 🔁	Back to Normal display.

From Normal display with the arrow-up button to the Set-point entry menu. Image: Set Point Set Po			
Image: Set Point is displayed. Image: Set Point is display. Image: Set Point is displayed. Image: Set Point flashes		From Normal display	
Sete Press OK to enter the temperature set-point. Image: Sete The temperature set-point flashes. Image: Sete Enter the temperature set-point with arrow buttons with an accuracy of a tenth of a degree (BD, BF) or of one degree (ED, FD, FED) Image: Sete Enter the temperature set-point with arrow buttons with an accuracy of a tenth of a degree (BD, BF) or of one degree (ED, FD, FED) Image: Sete and confirm with OK. Image: Sete Back to Normal display. Image: Sete with chambers with adjustable fan speed (BF, FED): Image: Sete go on to enter the fan speed. Image: Sete The fan speed set-point is displayed. Image: Sete Press OK to enter the fan speed Image: Sete The fan speed set-point flashes Image: Sete The fan speed set-point flashes Image: Sete The fan speed in steps of 10 % with arrow buttons		with the arrow-up button to the Set-point entry menu.	
Image: Set	Set _P	The current temperature set-point is displayed.	
Set Image: Se	ок	Press OK to enter the temperature set-point.	
of a degree (BD, BF) or of one degree (ED, FD, FED) and confirm with OK. Back to Normal display. or with chambers with adjustable fan speed (BF, FED): go on to enter the fan speed. Image: Second	SetP	The temperature set-point flashes.	
Back to Normal display. or with chambers with adjustable fan speed (BF, FED): go on to enter the fan speed. Image: Family or the fan speed set-point is displayed. Image: Press OK to enter the fan speed Image: Press OK to enter the fan speed set-point flashes Image: Press OK to enter the fan speed in steps of 10 % with arrow buttons			
or with chambers with adjustable fan speed (BF, FED): go on to enter the fan speed. Image: Second	ок	and confirm with OK.	
go on to enter the fan speed. Image: Fan	D	Back to Normal display.	
Image: Constraint of the fan speed set-point is displayed. Image: Constraint of the fan speed Image: Constraint of the fan speed Image: Constraint of the fan speed set-point flashes Image: Constraint of the fan speed in steps of 10 % with arrow buttons	or	with chambers with adjustable fan speed (BF, FED):	
Fan Press OK to enter the fan speed Image: Constant of the fan speed set-point flashes The fan speed set-point flashes Image: Constant of the fan speed in steps of 10 % with arrow buttons Adjust the fan speed in steps of 10 % with arrow buttons		go on to enter the fan speed.	
Image: Constraint of the fan speed set-point flashes Image: Constraint of the fan speed in steps of 10 % with arrow buttons	Fan	The fan speed set-point is displayed.	
Adjust the fan speed in steps of 10 % with arrow buttons	ок	Press OK to enter the fan speed	
	Fark	The fan speed set-point flashes	
and confirm with OK.		Adjust the fan speed in steps of 10 % with arrow buttons	
	ок	and confirm with OK.	
Back to Normal display.		Back to Normal display.	


Check and/or adjust the safety controller following any changes of the set-point (chap. 7).

6.6 Set-point entry for temperature and fan speed

6.6.1 Set-point entry for temperature in two-door chambers (ED, FD, FED 720)

With two-door chambers (ED, FD, FED 720) the maximum adjustable temperature set-point depends on the ambient temperature:

- Range from 18 °C up to 26 °C: maximum temperature set-point: 300 °C
- Range from above 26 °C up to 40 °C: maximum temperature set-point decreasing with increasing ambient temperature





Ambient temperature	Maximum temperature set-point
18 °C up to 26 °C	300 °C
28 °C	298 °C
30 °C	296 °C
32 °C	294 °C
34 °C	292 °C
36 °C	290 °C
38 °C	288 °C
40 °C	286 °C

This ensures the maximum lifetime of the controller.

6.7 Adjusting the air flap position

Opening the air flap in the exhaust duct serves to adjust the air change.

The position of the air flap in the exhaust duct serves to adjust the fresh air entry. With the open the air flap, fresh air can enter through the fresh air tube. For chambers with fan, fan operation will increase fresh air entry.

If the air flap is completely open, the spatial temperature accuracy can be negatively influenced.

	From Normal display	
without fan 2x with fan 3x	with the arrow-up button to the Adjusting the air flap position menu.	
Flap	The current air flap position is shown.	
ок	Press OK to select the air flap position.	
Flap	The setting flashes.	
	Select the position with arrow buttons	
ок	and confirm with OK.	
Ŋ	Back to Normal display.	

There are the following options:

Air flap closed	Air flap slightly opened	Air flap half opened	Air flap almost open	Air flap open
Flap	Image: constraint of the second se	₩ Flap	₩ °F ₩ € Flap	Flap

The setting can be done in steps of 15°.

As part of an automatic function test, the exhaust air flap opens briefly every 24 hours.

6.8 Changing the passwords for user level and general controller functions

In this menu you can change the passwords for access to the user menu and to all controller functions.

You can set two passworts for different access levels:

L1 (level 1): The password enables access control to the user level

L2 (level 2): The password enables access control to all controller functions

Factory setting for both passwords is 00 00 (no password assigned).

As soon as a password has been assigned, access to the respective functions is blocked and only available after enering the correct password.

	From Normal display		
without fan 5x with fan 6x	with the arrow-up button to the user menu		
ок	Confirm with OK.		
Password	Enter the password (factory setting: 00 00) and confirm each entry with OK.		
2 x 🔼	With the arrow-up button to the password setting menu .		
Password	Password L1 for access to the user level. Confirm with OK.		
ELAN S	The current password level L1 is shown. The setting flashes. You can change beween L1 and L2 with arrow buttons.		
	Select the setting with arrow buttons (if desired)		
ок	and confirm with OK.		
Assword	The current password for the selected password level is shown. The left two digits are flashing.		
	Enter the desired numbers with arrow buttons,		
ок	confirm with OK and go on.		

BINDER

	The right two digits of the password are flashing.	
	Enter the desired numbers with arrow buttons	
ок	and confirm with OK.	
Password	The modified password (L1 or L2 depending on the selection) is shown (example: L1).	
	If you want to change beween L1 and L2 , confirm with OK. Thereafter you can change to the other password level and also modify the password.	
2x 🔁	Back to Normal display.	



Keep in mind any modification of the password. There is no access to the user menu without the correct password L1. Without the correct password L2 access control to all controller functions is blocked.

7. Overtemperature protection

7.1 Overtemperature protective device (class 1)

The chambers are regularly equipped with an overtemperature protective device (safety device class 1 acc. to DIN 12880:2007). It serves to protect the chamber, its environment and the contents against exceeding the maximum permissible temperature. When a defined temperature is reached, which is approx. by 20 °C to 30 °C above the chamber's nominal temperature, the overtemperature protective device turns off the heating.

Cut-off temperature values:

BD, BF: 120 °C

ED 260, ED 720: 320 °C

ED 56, ED 115, FD 56, FED 56: 330 °C

FD 115, FD 260, FD 720, FED 115, FED 260, FED 720: 350 °C

The message "Overtemperature" is displayed on the controller.



If the overtemperature protective device class 1 has turned off the heating, proceed as follows:

- Disconnect the chamber from the power supply for at least 10 seconds (pull the power plug).
- If appropriate, have an expert examine and rectify the cause of the fault.
- Let the chamber cool down
- Restart the chamber.

As soon as the inner chamber temperature after restart is below the defined cut-off temperature of the overtemperature protective device class 1, the alarm message is deleted automatically.

Reset temperature values:

BD, BF: 90 °C ED 260, ED 720: 220 °C ED 56, ED 115, FD 56, FED 56: 230 °C FD 115, FD 260, FD 720, FED 115, FED 260, FED 720: 250 °C

7.2 Safety controller

The chambers are regularly equipped with an adjustable electronic safety controller. It serves to protect the chamber, its environment and the contents against exceeding the maximum permissible temperature. Please observe the regulations applicable to your country (for Germany: DGUV guidelines 213-850 on safe working in laboratories, issued by the employers' liability insurance association).

Depending on the chamber type the safety controller acts as an over temperature safety device class 2 ("temperature limiter") or class 3.1 ("temperature protection") acc. to DIN 12880:2007.

Check the setting regularly and adjust it following any changes of the set-point.

• Safety controller class 2 ("temperature limiter") with ED, FD and FED

The safety controller class 2 limits the temperature inside the chamber to the entered safety controller set-point. In the event of a fault (if this maximum temperature is exceeded) the safety controller completely turns off the heating until manual reset. This status is reported visually by an alarm message and, in case of the option audible alarm with activated buzzer (chap. 7.6) additionally by the buzzer sounding.

If the safety controller class 2 has turned off the heating, we recommend proceeding as follows:

- Disconnect the chamber from the power supply.
- Have an expert examine and rectify the cause of the fault.
- Restart the chamber
- Reset the alarm message

• Safety controller class 3.1 ("temperature protection") with BD and BF

The safety controller class 3.1 limits the temperature inside the chamber to the entered safety controller set-point. In the event of a fault (if this maximum temperature is exceeded), it takes over the control to this value. This status is reported visually by an alarm message and, in case of the option audible alarm with activated buzzer (chap. 7.6) additionally by the buzzer sounding.

The safety controller keeps control of the chamber until the chamber temperature cools down below the safety controller set-point value.

If the safety controller class 3.1 has taken over control, we recommend proceeding as follows:

- Disconnect the chamber from the power supply.
- Have an expert examine and rectify the cause of the fault.
- Restart the chamber
- Reset the alarm message

Function check:

Check the safety controller at appropriate intervals for its functionality. It is recommended that the authorized operating personnel should perform such a check, e.g., before starting a longer work procedure.

7.3 Setting the safety controller set-point

A limit temperature is entered as the safety controller set-point , i.e. the absolute maximum permitted temperature value.

Example: Temperature set-point 45 °C, safety controller set-point 50 °C.

Regularly check the safety controller setting relating to the entered temperature set-point
 Set the safety controller set-point by approx. 2 °C to 5 °C above the desired temperature set-point.

without fan 4x with fan 5x	From Normal display with the arrow-up button to the Safety controller set-point setting menu.
TLim	The current safety controller set-point is shown (class 2 "temperature limiter" or class 3.1 "temperature protection" depend- ing on the chamber type).
ок	Press OK to enter the safety controller set-point.
	The safety controller set-point flashes.
	Enter the safety controller set-point with arrow buttons: 10 °C up to 100 °C (with an accuracy of a tenth of a degree) with BD, BF 10 °C up to 300 °C (with an accuracy of one degree) with ED, FD, FED
ок	and confirm with OK.
n	Back to Normal display.

7.4 Alarm message and proceeding in case of an alarm

The alarm icon flashes in the display. In case of the option audible alarm with activated buzzer (chap. 7.6) the buzzer sounds.

• Safety controller class 2 ("temperature limiter")



The alarm icon flashes in the display. In case of the option audible alarm with activated buzzer (chap. 7.6) the buzzer sounds.

The heating turns off.

Resetting the alarm:

With option audible alarm with activated buzzer: Mute the buzzer pressing the OK button.

As soon as the inner chamber temperature has cooled down below the safety controller set-point, the alarm icon is lit permanently. You can reset the alarm message on the controller. To do this, reset the alarm message in the safety controller menu with the OK button. The heating is released and temperature control is resumed by the controller.

• Safety controller class 3.1 ("temperature protection")



The alarm icon flashes in the display. In case of the option audible alarm with activated buzzer (chap. 7.6) the buzzer sounds.

The heating turns off.

Resetting the alarm:

With option audible alarm with activated buzzer: Mute the buzzer pressing the OK button.

As soon as the inner chamber temperature has cooled down below the safety controller set-point, You can reset the alarm message in the safety controller menu with the OK button, The heating is released and temperature control is resumed by the controller.

Note:

When the safety controller class 2 or class 3.1 had been activated, you should disconnect the chamber from the power supply and have an expert examine and rectify the cause of the fault.

7.5 Function check

Check the temperature safety device class 2 or class 3.1 at appropriate intervals for its functionality. It is recommended that the authorized operating personnel should perform such a check, e.g., before starting a longer work procedure

7.6 Disconnectable audible over-temperature alarm (option)

This option permits activating an audible signal:

If the buzzer is activated, an audible signal sounds when the limit temperature set at the safety controller is exceeded. This happens in addition to the alarm message on the controller display.

Turning off the audible alarm does not influence the safety controller's function.

	From Normal display
without fan 5x with fan 6x	with the arrow-up button to the user menu
ок	Confirm with OK.
Password	Enter the password (factory setting: 00 00) and confirm each entry with OK.
5 x 🔼	With the arrow-up button to the alarm buzzer setting menu
Buzzer	The current setting is shown.
ок	Press OK to select the alarm buzzer setting
BUZI	The setting flashes.
	Select the setting with arrow buttons
ок	and confirm with OK.
n	Back to Normal display.

There are the following options

Alarm buzzer On	Alarm buzzer Off
Buzzer	Buzzer

8. Timer functions

8.1 Selecting the timer function

There are up to three 3 timer functions:

DelayOff	Timer function "Delayed Off"The selected timer run-time immediately starts running down.When the timer expires the heating turns off. Chambers with fan: The fan may continue working according to the selections made.	
DelayTemp.	Timer function "Temperature dependent Delayed Off" The selected timer run-time only starts running down, when the actual value reached or exceeds the selected set-point. When the timer expires the heating turns off. Chambers with fan: The fan may continue working according to the selections made.	
DelayOn	Timer function "Delayed On" The selected timer run-time immediately starts running down, the heating turns off. Chambers with fan: The fan may be working according to the selections made. When the timer expires the heating turns on and remains in continuous opera- tion.	

The chambers BD, ED and FD offer the timer function "Delayed Off"

The chambers BF and FED offer all the three timer functions.

Stage	Heating	Fan (Chamber with fixed fan speed: FD)	Fan (Chambers with adjustable fan speed: BF, FED)
Timer function "De	elayed Off"		
Timer running	Control to the temperature set-point	On (100 %)	Rotation speed according to fan speed set-point
After the timer ex- pired	Off	On (100 %) or Off (0 %) acc. to selection	Rotation speed according to setting of timer function
Timer function "Te	emperature dependent Del	ayed Off"	
Possibly heating- up phase until the temperature set- point is reached	Control to the temperature set-point		Rotation speed according to fan speed set-point
Timer running	Control to the temperature set-point		Rotation speed according to fan speed set-point
After the timer ex- pired	Off		Rotation speed according to setting of timer function
Timer function "Delayed On"			
Timer running	Off		Rotation speed according to setting of timer function
After the timer expired	Control to the temperature set-point		Rotation speed according to fan speed set-point

General information on the setting:

In the setting menus of the timer functions, it is always required to confirm **all** parameters with OK, otherwise all entries made will be lost.

Timer run-time is set in days, hours, and minutes. If days have been entered, they are shown in the controller display preceded by an underscore:



Setting: 0 days (not shown), 10 hours, 30 minutes




Setting: 2 days, 10 hours (minutes not shown)

8.2 Timer function "Delayed Off"

8.2.1 Entry and activation of the timer run-time and fan setting

	From Normal display
	with the arrow-down button to the Timer function "Delayed Off" menu (with connected USB device: press the arrow-down button twice)
DelayOff	Current Timer function "Delayed Off"
ок	Confirm with OK and go on to enter the days of the timer run-time .
Days	The current timer run-time (days) is shown. The day value flashes.
	Enter the days of the timer run-time (0 up to 9)
ок	confirm with OK and go on to enter the hours of the timer run-time
Hours	The hour value flashes.
	Enter the hours of the timer run-time (0 up to 23)
ок	confirm with OK and go on to enter the minutes of the timer run-time



Chambers without fan (BD, ED):	
Image: Constraint of the second secon	The timer function "Delayed Off" is activated
D	Back to Normal display.

Chamber with fixed fan speed (FD): Select whether the fan shall operate or not after the timer has expired	
	The current setting of the fan operation is shown
	Select fan operation: On (100 %) or Off (0 %)
ок	and confirm with OK.
Image: Constraint of the second se	The timer function "Delayed Off" is activated
D	Back to Normal display.

Chambers with adjustable fan speed (BF, FED): Enter the fan speed set-point valid for the time after the timer has expired.	
Fand	The current fan speed set-point is shown
	Enter the fan speed in steps of 10 % 0 % and 40 % up to 100 %
ок	and confirm with OK.
@	The timer function "Delayed Off" is activated
n	Back to Normal display.

Normal display during timer operation with timer function "Delayed Off"



The timer run-time until turning off the heating is running down.

Normal display after the timer has expired:

© DelayOff	Chambers without fan (BD, ED): Timer function "Delayed Off". The timer has expired. The heating is off.
Off FOn	Chamber with fixed fan speed (FD): Timer function "Delayed Off". The timer has expired. The heating is off. The fan is operating (On) or not (Off), as selected
Coff 100%	Chambers with adjustable fan speed (BF, FED): Timer function "Delayed Off". The timer has expired. The heating is off. The fan operates with the selected fan speed.

When the timer has expired, the heating is off. The chamber cools down to ambient temperature.

To restart the chamber you need to turn off the timer function (chap. 8.2.2)

8.2.2 Turning off the timer function or changing the settings

To turn off the timer function "Delayed Off" during the time when it is still running, set all time values (days, hours, minutes) to zero in the according menu (terminating the timer function). As long as the timer is running, the timer function settings can be subsequently modified in this menu.

After the timer expired you can turn off the timer function by pressing the OK button. Alternatively you can also set the time to zero in the according menu.

BINDER

Turning off the timer function when the timer has expired

In Normal display press the OK button.

Alternatively you can set the time to zero:

	From Normal display
	with the arrow-down button to the Timer function "Delayed Off" menu (with connected USB device: press the arrow-down button twice)
DelayOff Or	Chambers without fan (BD, ED): Timer function "Delayed Off", the timer has expired. The heating is off.
Image: Constraint of the second secon	Chamber with fixed fan speed (FD): Timer function "Delayed Off", the timer has expired. The heating is off. The fan is operating (On) or not (Off), as selected
or	
	Chambers with adjustable fan speed (BF, FED):
00000000000000000000000000000000000000	Timer function "Delayed Off", the timer has expired. The heating is off. The fan operates with the selected fan speed.
ок	Confirm the timer function "Delayed Off" with OK
see chap. 8.2.1	Set all time values to zero
2x 🔁	Back to Normal display.

Changing or terminating the timer function when the timer is running:

	From Normal display
\bigtriangledown	with the arrow-down button to the Timer function "Delayed Off" menu (with connected USB device: press the arrow-down button twice)
Fan	Chambers without fan (BD, ED): Timer function "Delayed Off", the timer is running.
or Off Fon or	Chamber with fixed fan speed (FD): Timer function "Delayed Off", the timer is running.
@	Chambers with adjustable fan speed (BF, FED): Timer function "Delayed Off", the timer is running.
ок	Confirm the timer function "Delayed Off" with OK
see chap. 8.2.1	Modify the time values or set them to zero to terminate the timer function
2x 🔁	Back to Normal display.

8.3 Timer function "Temperature dependent Delayed Off" (BF, FED)

8.3.1 Entry and activation of the timer run-time, fan setting and set-point entry

	From Normal display
3 x 🔽	with the arrow-down button to the Timer function "Temperature dependent Delayed Off " menu (with connected USB device: press the arrow-down button 4 times)
DelayTemp.	Current Timer function "Temperature dependent Delayed Off"
ок	Confirm with OK and go on to enter the days of the timer run-time .
Days	The current timer run-time (days) is shown. The day value flashes.
	Enter the days of the timer run-time (0 up to 9)
ок	confirm with OK and go on to enter the hours of the timer run-time
Hours	The hour value flashes.
\triangle	Enter the hours of the timer run-time (0 up to 23)
ок	confirm with OK and go on to enter the minutes of the timer run-time
	The minute value flashes.
	Enter the minutes of the timer run-time (0 up to 59)
ок	confirm with OK and go on to enter the temperature set-point

BINDER

Set P	The temperature set-point flashes. When the set-point is reached, the timer starts running down.
	Enter the temperature set-point,
ок	confirm with OK and go on to enter the fan speed set-point valid for the time af- ter the timer has expired
Fand	The current fan speed is shown
	Enter the fan speed in steps of 10 % 0 % and 40 % up to 100 %
ок	and confirm with OK.
Image: Constraint of the second sec	The timer function "Temperature dependent Delayed Off" is activated
D	Back to Normal display.

Normal display during timer operation with timer function "Temperature dependent Delayed Off"

₩	Normal display while the timer is running
CG	If the entered set-point exceeds the current actual value, the chamber heats
Toff 40	up.
↓ ↓ ↓ ↓ ↓ ↓ ↓ ↓ ↓ ↓ ↓ ↓	The defined time only begins to run when the current value is by 1 °C below the set point. After the defined time has expired, the heating is turned off. Normal display during the timer run-time. The timer run-time until turning off the heating (and possibly the fan) is running.

Normal display after the timer has expired:

G Off 100%	Timer function "Temperature dependent Delayed Off". The timer has expired. The heating is off. The fan operates with the selected fan speed.
------------	-------------------------------------------------------------------------------------------------------------------------------------------------

When the timer has expired, the heating is off. The chamber cools down to ambient temperature.

To restart the chamber you need to turn off the timer function (chap. 8.3.2).

8.3.2 Turning off the timer function or changing the settings

To turn off the timer function "Temperature dependent Delayed Off" during the time when it is still running, set all time values (days, hours, minutes) to zero in the according menu (terminating the timer function). As long as the timer is running, the timer function settings can be subsequently modified in this menu.

After the timer expired you can turn off the timer function by pressing the OK button. Alternatively you can also set the time to zero in the according menu.

Turning off the timer function when the timer has expired

In Normal display press the OK button.

Alternatively you can set the time to zero:

	From Normal display
3 x 🔽	with the arrow-down button to the Timer function "Temperature dependent Delayed Off" menu (with connected USB device: press the arrow-down button 4 times)
Coff tOFF	Timer function "Temperature dependent Delayed Off", the timer has expired
ок	Confirm the timer function "Temperature dependent Delayed Off" with OK
see chap. 8.3.1	Set all time values to zero
Ŋ	Back to Normal display

Changing or terminating the timer function when the timer is running:

	From Normal display	
3 x 🔽	with the arrow-down button to the Timer function "Temperature dependent Delayed Off" menu (with connected USB device: press the arrow-down button 4 times)	
B Toff 0:08	Timer function "Temperature dependent Delayed Off", the timer is running	
ок	Confirm the timer function "Temperature dependent Delayed Off" with OK	
see chap. 8.3.1	Modify the time values or set all time values to zero to terminate the timer func- tion	
Ŋ	Back to Normal display	

8.4 Timer function "Delayed On" (BF, FED)

8.4.1 Entry and activation of the timer run-time and fan setting

	From Normal display	
2 x 🔽	with the arrow-down button to the Timer function "Delayed On" menu (with connected USB device: press the arrow-down button 3 times)	
DelayOn	Current timer function "Delayed On"	
ок	Confirm with OK and go on to enter the fan speed set-point valid during the time the timer is running	
FanOper	The current fan speed set-point is shown	
	Enter the fan speed in steps of 10 % 0 % and 40 % up to 100 %	
ок	Confirm with OK and go on to enter the days of the timer run-time .	
Days	The current timer run-time (days) is shown. The day value flashes.	
	Enter the days of the timer run-time (0 up to 9)	
ок	confirm with OK and go on to enter the hours of the timer run-time	
Hours	The hour value flashes.	
	Enter the hours of the timer run-time (0 up to 23),	
ок	confirm with OK and go on to enter the minutes of the timer run-time	

BINDER

	The minute value flashes.
	Enter the minutes of the timer run-time (0 up to 59)
ок	confirm with OK and go on to enter the temperature set-point
Set#	The temperature set-point flashes.
	Enter the temperature set-point,
ок	confirm with OK and go on to enter the fan speed set-point valid for the time af- ter the timer has expired
Fart	The current fan speed set-point is shown
	Enter the fan speed in steps of 10 % - 0 % and 40 % up to 100 %
ок	and confirm with OK.
© 0n 100%	The timer function "Delayed On" is activated
2x 🔁	Back to Normal display.

Normal display during timer operation with timer function "Delayed On"

122		
6		
	0n	00:10

The Timer run-time until turning on the heating is running. Timer function "Delayed On". The heating is off, temperature approximates ambient temperature.

Normal display after the timer has expired:



The timer has expired. The time function is off. The heating is active to equilibrate the temperature set-point. The fan operates with the selected fan speed.

8.4.2 Changing the settings

After the timer expired, the timer function "Delayed On" deactivates, therefore turning the function off is not required.

As long as the timer is still running, the timer function settings can be subsequently modified in this menu.

To terminate the timer function, all time values (days, hours, minutes) must be set to zero in the according menu.

Changing or terminating the timer function when the timer is running:

© 0n 00:04	From Normal display during timer operation
2 x 🔽	with the arrow-down button to the Timer function "Delayed On" menu. (with connected USB device: press the arrow-down button 3 times)
© 100%	Timer function "Delayed On", the timer is running. The fan operates with the selected fan speed.
ок	Confirm the timer function "Delayed On" with OK
see chap. 8.4.1	Modify the time values or set all time values to zero to terminate the timer function
2x 🔁	Back to Normal display.

8.5 Temperature programming example (BF, FED)

The chamber shall heat up to a temperature of 50 $^\circ$ C, maintain this temperature for three hours and then turn off.

Proceeding: Select timer function "Temperature dependent Delayed Off" (chap. 8.3) and perform the following settings:

- Enter a timer run-time of 3 hours
- Enter the set point 50 °C
- Specify the fan speed after the timer expires

9. Ramp function

9.1 General information

You can program temperature ramps in order to extend heating up times. This may be necessary in some cases to prevent temperature stress in the material during the heating up phase. Temperature ramps should only be used if required. Using them may result in considerably slowing down the heating up times. When the ramp function is turned off, the chamber will heat up with its maximum heating capacity.

The entry means the nominal value gradient and limits the maximum temperature increase to this value. Due to the heat and evaporation energy assumed by the drying material, smaller temperature gradients may also result.

A ramp proceeds from a previously entered set-point to a new, higher one. The temperature must be equilibrated to the start set-point. Perform the setting in the following 3 steps:

- 1. Enter the temperature set-point as start ramp set-point and let the temperature equilibrate to this value
- 2. Define the temperature increase (**ramp gradient**) in °C/min or in °F/min in the setting menu "Ramp function"

You can select a gradient from "0.0" up to "1.0" or from "1" up to "10" according to the chamber type.

When setting the gradient to "0.0" or "0", ramp function is turned off. The chamber will then heat up with its maximum heating capacity.

The chamber will try to heat up according to the entered gradient, i.e. with a speed of xx degree per minute. A heating-up rate of 0.4 °C/min for the incubators BD and BF resp. 4 °C/min for the heating and drying ovens ED, FD and FED can be regarded as a realistic maximum.

3. Enter the target ramp set-point in the "ramp function" setting menu.

As soon as the entries have been adopted, the ramp function is activated. The chamber heats up with the entered gradient, if the set ramp target value lies above the actual temperature value.


During ramp operation the **effective ramp set-point** continually rises in accordance to the entered gradient from the previously entered set-point to the new one. The actual value follows this continually changing effective ramp set-point. As soon as the ramp target value is reached, this temperature is maintained constant.

The actual temperature value, the selected gradient, and the target value are shown in Normal display. The effective ramp set-point can be seen through the temperature set-point function.

	From Normal display
without fan 3x with fan 4x	with the arrow-up button to the Ramp function menu.
	Ramp function (not programmed)
ок	Press OK to enter the gradient in degree per minute.
mire	The gradient flashes.
	Enter the gradient (0 up to 9)
ок	confirm with OK and go on toenter the target ramp set-point
G SetP	The target ramp set-point value flashes.
	Enter the target ramp set-point
ок	and confirm with OK.
°C ~ 1 40	Selected ramp with gradient 1 and target ramp set-point 40 °C (example)
D	Back to Normal display.
<pre></pre>	Normal display showing the actual temperature, the gradient and the target ramp set-point (example: 40 °C)

9.2 Setting and displaying the ramp function

Temperature course with ramp function



9.3 Displaying the effective ramp set-point and changing the target ramp setpoint

- 1 40	From Normal display showing the actual temperature (example: 27 $^{\circ}$ C), the gradient (example: 1) and the target ramp set-point (example: 40 $^{\circ}$ C)
Δ	to the display of the effective ramp set-point
SetP	The effective ramp set-point is displayed
ок	Press OK to set the target ramp set-point
Setr	The target ramp set-point value flashes.
	Enter the target ramp set-point
ок	and confirm with OK.
SetP	The effective ramp set-point is displayed
2x 🗾	Back to Normal display.
	Normal display showing the actual temperature, the gradient and the modi- fied target ramp set-point (example: 50 °C)

9.4 Turning off the ramp function

To turn off the ramp function, the gradient must be set to zero in the according menu. The set-point can be entered as desired.

	From Normal display	
without fan 3x with fan 4x	with the arrow-up button to the Ramp function menu	
¢ • • 1 50	Ramp function with programmed gradient and target ramp set-point	
ок	Press OK to enter the gradient in degree per minute.	
/min	The gradient flashes.	
	Enter the gradient zero (turning off the ramp function)	
ок	confirm with OK and go on to the ramp target value	
Set P	The target ramp set-point value flashes.	
	Enter the target ramp set-point	
ок	and confirm with OK.	
	Ramp function (not programmed)	
D	Back to Normal display.	
- 1 40	Normal display showing the actual temperature, the gradient and the target ramp set-point (example: 40 $^\circ\text{C})$	
2x 🔁	Back to Normal display.	

Instead of turning off the ramp function (gradient = 0), you can also modify the settings of the gradient and target ramp set-point in this menu.

10. Data recording via USB interface

The USB interface located in the triangular instrument panel serves to read out the measured values, which are put out in real time. Following data is recorded: Timer, actual temperature value, temperature set-point, Object temperature sensor (chambers with option Object temperature display), Fan (chambers with fan), air flap position, safety controller, analog output (option), heating regulation ratio.



Connect only USB sticks to the USB interface.


Data are stored directly on the USB stick. They are issued in the selected language as a spreadsheet with the file extension ".csv" and can be further processed in the desired program.

10.1 Starting data recording

Connect the USB stick to the interface located in the triangular instrument panel.

	From Normal display	
	with the arrow-down button to the USB menu	
	Current state of data recording: stopped	
ок	Confirm with OK.	

Chambers without the real time clock option do not save date and time after shut-down. To make sure that recorded data is provided with the correct date and time, with these chambers it is required to enter date and time first :



After entry of the minute and confirmation with OK, further setting in the USB menu continues.

Inter	The data recording interval flashes.
	Enter the interval in minutes
ок	and confirm with OK.
	The current state of data recording (stopped) flashes.
	Change to start data recording.

	The current state of data recording (started) flashes	
ок	Confirm with OK.	
	Data recording is running.	
	Back to Normal display.	

Data recording continues also during stand-by mode of the chamber. Disconnecting the power supply interrupts data recording, which continues after the power returns. To terminate data recording, stop it via the menu (chap. 10.2).

10.2 Terminating data recording

	From Normal display
\bigtriangledown	with the arrow-down button to the USB menu
	Current state of data recording: running
ок	Confirm with OK.
	The current state of data recording (running) flashes
	Change to stop data recording
	The current state of data recording (stopped) flashes
ок	Confirm with OK.
	Data recording is stopped.
D	Back to Normal display.



11. Network configuration for chambers with Ethernet interface

The settings of this submenu are used for networking chambers with an Ethernet interface, e.g. to connect them with BINDER's APT-COM[™] 4 Multi Management Software (option, chap. 12.1).

This menu subsequently offers the following settings:

- Display of the chamber's MAC address (no setting)
- Switching on and off the DHCP state
- Entering the IP address
- Entering the **subnet** mask number
- Entering the gateway number

	From Normal display			
without fan 5x with fan 6x	with the arrow-up button to the user menu			
ок	Confirm with OK.			
Password	Enter the password (factory setting: 00 00) and confirm each entry with OK.			
5 x 🔼	with the arrow-up button to the Ethernet menu.			
ок	Confirm with OK.			
MRC_Hish	The first digits of the MAC address are shown.			
ок	Continue with OK. The middle digits of the MAC address are shown.			
ок	Continue with OK. The last digits of the MAC address are shown.			


Displaying the chamber's MAC address serves to identify the chamber in the Ethernet network.

Example: 00-0F-67-0F-42-40



ок	Continue with OK to switching on and off the DHCP state	
	The current DHCP state is displayed (on). The setting flashes.	
	Select the setting with the arrow buttons.	


There are the following options:



To configure the network settings, the DHCP state must be switched off. Otherwise, the DHCP-server would assign the network configuration..

ок	Confirm with OK and go on to enter the IP address .
IPAdr_A	The first digits of the current IP address are shown. The setting flashes.
	Enter the desired values with the arrow buttons
	and confirm with OK.
ок	The next digits of the IP address are shown.
_	Enter the desired values accordingly.

Example value: 0.0.0.0



BINDER

ок	Confirm with OK and go on to enter the subnet mask .
	The first digits of the current subnet mask are shown. The setting flashes.
	Enter the desired values with the arrow buttons
ок	and confirm with OK. The next digits of the subnet mask are shown. Enter the desired values accordingly.

Example value: 255.255.255.0

255	255	255	0
Subnet_H	Subnet_H	Subnet_L	Subnet_LL

ок	Confirm with OK and go on to enter the gateway .
GatWay	The first digits of the current gateway are shown. The setting flashes.
	Enter the desired values with the arrow buttons
	and confirm with OK.
ок	The next digits of the gateway are shown.
	Enter the desired values accordingly.

Example value: 0.0.0.0



ок	Confirm with OK.
2x 🔁	Back to Normal display.

12. Options

12.1 APT-COM[™] 4 Multi Management Software (option)

The chamber is regularly (FED) or optionally equipped with an Ethernet interface on the chamber rear that can connect the BINDER APT-COM[™] 4 Multi Management Software. The actual temperature value is given at adjustable intervals. Programming can be performed graphically via PC. Up to 100 chambers can be cross-linked. The MAC Address is indicated in the **Ethernet** menu (chap. 11). For further information, please refer to the APT-COM[™] 4 operating manual.

To establish a connection via the chamber's Ethernet interface, the chamber must be turned off.

12.2 Object temperature display with additional Pt100 temperature sensor (option)

With this option an additional flexible temperature sensor Pt100 measures the chamber temperature or the temperature of the charging material which is shown on the controller. The sensor-top protective tube of the flexible Pt100 can be immersed into liquid substances.

The object temperature display enables the determination of the actual temperature of the charging material during the whole process. The object temperature is displayed in the controller in Normal Display.



Chamber with option object temperature display: Actual temperature value and actual object temperature value

Technical data of thePt100 sensor:

- Three-wire technique
- Class B (DIN EN 60751)
- Temperature range up to 320 °C / 608 °F
- Stainless steel protective tube 45 mm length, material no. 1.4501

12.3 Analog output for temperature (option)

With this option the chamber is equipped with an analog output 4-20 mA for temperature. This output permits transmitting data to external data registration systems or devices.

The connection is carried out as a DIN socket at the rear of the chamber as follows:



ANALOG OUTPUT 4-20 mA DC

PIN 1: Temperature – PIN 2: Temperature + **Temperature range:** BD, BF: 0 °C / 32 °F to +100 °C / 212 °F ED, FD, FED: 0 °C / 32 °F to +300 °C / 572 °F A suitable DIN plug is enclosed.

Figure 11: Pin allocation of DIN socket for option analog outputs

12.4 HEPA fresh air filter (option for FD, FED)

With this option, the introduced fresh air is cleaned by means of a high efficiency submicron particulate air filter type HEPA class H 14 (acc. to DIN EN 1822:2009). Replace the filter insert, if necessary, by removing the metal cover of the filter at the left side of the chamber (Art. No. 6014-0003).

12.5 Water protected disconnectable internal socket (option BD, BF)


You can turn on or off the disconnectable water protected internal socket (CA3GD) for line voltage inside the chamber by the switch in the triangular instrument panel, independent of the incubator operating or not. Thus, devices operated inside the incubator can be started or stopped without any need to open the chamber doors.

The internal socket provides a protective conductor, therefore it is possible to connect appliances with protection class I.


The internal socket is splash proof.

IP system of protection 67 230 V 1N ~ 50-60 Hz. Charge max. 500 W

The maximum permitted operating temperature with this option is 90 °C / 194 °F.



Heat emission of electrical devices connected inside the chamber may modify the temperature range.



Controller shutdown by the standby button has no effect on the internal socket. The internal socket can remain switched on even though controller is shut down.



12.6 Mostly gas-tight version (option for BF, FD, FED)

With this option the chamber is additionally sealed, so the loss when introducing gases is decreased. The chamber is not completely gas-tight, so it is impossible to establish overpressure. The sealing diminishes the release of vapors via the housing that may be set free from the charging material when heated. Carry-ing-off via the regular evacuation duct, e.g. into a waste air installation, is likely to further reduce emissions.

The chamber is not completely gas-tight. Gases from inside the chamber can escape into the surrounding atmosphere.
 Observe the occupational exposure limit OEL for the released substance set by the national authorities (formerly maximum permitted workplace concentration). Respect the relevant regulations.
 Any harmful gas that might escape has to be led out via good room ventilation or a suitable exhaust system. Place the chamber, if necessary, below a gas vent.

The air flap does not close the exhaust duct completely. The delivered plug serves to avoid emerging of vapors or loss of introduced inert gas, if any, via the exhaust duct. Due to special demands of heat resistance, use the delivered plug only.



≻

NOTICE

BINDER

Danger of inflammation when using an inappropriate plug. Damage to the chamber and ist surroundings.

Use ONLY the supplied plug to close the exhaust duct.

For drying purpose, please remove the plug in order to permit dissipation of the generated vapor, which would lead to condensation in the inner chamber.

12.7 Inert gas connection with mostly gas-tight version (option for BF, FD, FED)

With this option the chamber is additionally sealed, so the loss when introducing inert gases is decreased. For details on the mostly gas-tight version please refer to chap. 12.7.

The chamber is equipped with two ports for inert gas (nitrogen or noble gases).

The ports are located **on the top panel in the middle** and **on the chamber rear at the bottom right**. Each of these ports can be used as inlet or outlet, depending on the nature of the inert gas:

- lighter gas (nitrogen, helium): lower port as inlet
- heavy gas (e.g. argon): upper port as inlet


This distinction is important when operating with a reduced fan speed.

Connection

Observe the legal requirements and relevant standards and regulations for the safe handling of gas cylinders and inert gases.

r	1	
3	G	eneral information for safe handling of gas cylinders:
3	•	Store and use gas cylinders only in well ventilated areas.
	•	Open the gas cylinder valve slowly to avoid pressure surges
	•	Secure gas cylinders during storage and use against falling (chaining).
	•	Transport gas cylinders with a cylinder cart, do not carry, roll, or throw them
	•	Always close the valve even with apparently empty cylinders; screw on the cap when not in use. Return gas cylinders with the valve closed
	•	Do not open gas cylinders by force. Mark them when damaged
	•	Observe relevant regulations for dealing with gas cylinders.

Connect a flexible gas tube to the gas hose connection adapter (diameter 10mm), which is used for gas inlet, and secure it with hose clamps (hose and hose clamps are not enclosed). There is a constant gas flow after establishing the connection.




Use a pressure reducer and make sure to avoid any excessive outlet pressure when connecting the gas hose to the chamber.



The chamber is not entirely gas-tight. Inert gases from inside the chamber can escape into the surrounding atmosphere.

Inert gases in high concentrations is hazardous to health. They are colorless and almost odorless and therefore practically imperceptible. Inhalation of inert gases can cause drowsiness up to respiratory arrest. When the O_2 content of the air decreases below 18%, there is risk of death from lack of oxygen. Any gas that might escape has to be led out via good room ventilation or a suitable exhaust system.

BINDER



Inert gases, which are heavier than air, may accumulate in low-lying areas of the installation site.

The "Mostly gas-tight version" reduces the loss of gas.

Setting (sample values):

If you want to flush the chamber with an air exchange rate of 1 per hour, set the flow rate on the pressure reducer according to the interior volume.

Chamber with 56 I internal volume: The flow rate corresponding to 56 I / h is 0.9 I / min.

Chamber with 115 I internal volume: The flow rate corresponding to 115 I / h is 1.9 I / min.

Chamber with 260 I internal volume: The flow rate corresponding to 260 I / h is 4.3 I / min.

The air flap does not close the exhaust duct completely. The delivered plug serves to avoid loss of introduced inert gas via the exhaust duct. Due to special demands of heat resistance, use the delivered plug only.



NOTICE


Danger of inflammation when using an inappropriate plug.
Damage to the chamber and ist surroundings.
> Use ONLY the supplied plug to close the exhaust duct.

For drying purpose, please remove the plug in order to permit dissipation of the generated vapor, which would lead to condensation in the inner chamber.

13. Cleaning and decontamination

Clean the chamber after each use in order to prevent potential corrosion damage by ingredients of the loading material.

Prior to renewed startup, allow the chamber to completely dry after all cleaning and decontamination measures.



13.1 Cleaning

Disconnect the chamber from the power supply before cleaning. Disconnect the power plug.

The interior of the chamber must be kept clean. Thoroughly remove any residues of the charging material

Wipe the surfaces with a moistened towel. In addition, you can use the following cleaning agents:

Exterior surfaces inner chamber racks door gaskets	Standard commercial cleaning detergents free from acid or halides. Alcohol-based solutions. We recommend using the neutral cleaning agent Art. No. 1002-0016.
Instrument panel	Standard commercial cleaning detergents free from acid or halides. We recommend using the neutral cleaning agent Art. No. 1002-0016.
Zinc coated hinge parts rear chamber wall	Standard commercial cleaning detergents free from acid or halides. Do NOT use a neutral cleaning agent on zinc coated surfaces.

Do not use cleaning agents that may cause a hazard due to reaction with components of the device or the charging material. If there is doubt regarding the suitability of cleaning products, please contact BINDER service.

We recommend using the neutral cleaning agent Art. No. 1002-0016 for a thorough cleaning. Any corrosive damage that may arise following use of other cleaning agents is excluded from liability by BINDER GmbH.


Any corrosive damage caused by a lack of cleaning, is excluded from liability by BINDER GmbH.



NOTICE

Danger of corrosion by using unsuitable cleaners. Damage to the chamber.

- Ø Do NOT use acidic or chlorine cleaning detergents.
- Ø Do NOT use a neutral cleaning agent on other kind of surfaces e.g., the zinc coated hinge parts or the rear chamber wall.



With every decontamination method, always use adequate personal safety controls.

Following cleaning, leave the chamber door open or remove the plugs of the optional access port.



The neutral cleaning agent may cause health problems in contact with skin and if ingested. Follow the operating instructions and safety hints labeled on the bottle of the neutral cleaning agent.

Recommended precautions: To protect the eyes use sealed protective goggles. Wear gloves. Suitable protective gloves in full contact with media: butyl or nitrile rubber, penetration time >480 minutes.

Danger of chemical burns through contact with skin or ingestion of the neutral cleaning agent.
Skin and eye damage. Environmental damage.
\varnothing Do not ingest the neutral cleaning agent. Keep it away from food and beverages.
arnothing Do NOT empty the neutral cleaning agent into drains.
Wear protective gloves and goggles.
Avoid skin contact with the neutral cleaning agent.

13.2 Decontamination / chemical disinfection


The operator must ensure that proper decontamination is performed in case a contamination of the chamber by hazardous substances has occurred.

Disconnect the chamber from the power supply prior to decontamination. Pull the power plug.

Do not use decontamination agents that may cause a hazard due to reaction with components of the device or the charging material. If there is doubt regarding the suitability of cleaning products, please contact BINDER service.

You can use the following disinfectants:

Inner chamber	Standard commercial surface disinfectants free from acid or halides.	
	Alcohol-based solutions.	
	We recommend using the disinfectant spray Art. No. 1002-0022.	



With every decontamination / disinfection method, always use adequate personal safety controls.

In case of impurity of the interior with biological or chemical hazardous material, there are three possible procedures depending on the type of contamination and of the charging material.

- The drying and heating ovens ED, FD and FED can be hot air sterilized at 190 °C / 374 °F for at least 30 minutes. All inflammable goods must be removed from the interior before. With the incubators BD and BF it is possible to perform a hot-air disinfection at 100 °C / 212 °F.
- 2. Spray the inner chamber with an appropriate disinfectant.

Before start-up, the chamber must be absolute dry and ventilated, because explosive gases may form during the decontamination process.

3. BD, ED: If necessary, have strongly contaminated inner chamber parts removed by an engineer for cleaning, or have them exchanged. Sterilize the inner chamber parts in a sterilizer or autoclave. You can also remove and sterilize the racks.



In case of eye contact, the disinfectant spray may cause eye damage due to chemical burns. Follow the operating instructions and safety hints labeled on the bottle of the disinfectant spray.

Recommended precautions: To protect the eyes use sealed protective goggles.



Danger of chemical burns through eye contact with the disinfectant spray.

Eye damage. Environmental damage

- \varnothing Do NOT empty the disinfectant into drains.
- Wear protective goggles.



After using the disinfectant spray, allow the chamber to dry thoroughly, and aerate it sufficiently.

14. Maintenance and service, troubleshooting, repair, testing

14.1 General information, personnel qualification

• Maintenance

See chap. 14.2.

• Simple troubleshooting

Chap. 14.3 describes troubleshooting by operating personnel. It does not require technical intervention into the chamber, nor disassembly of chamber parts.

For personnel requirements please refer to chap. 1.1.

Detailed troubleshooting

If errors cannot be identified with simple troubleshooting, further troubleshooting must be performed by BINDER Service or by BINDER qualified service partners or technicians, in accordance with the description in the Service Manual.

For personnel requirements please refer to the Service Manual.

Repair

Repair of the chamber can be performed by BINDER Service or by BINDER qualified service partners or technicians, in accordance with the description in the Service Manual.

After maintenance, the chamber must be tested prior to resuming operation.

• Electrical testing

To prevent the risk of electrical shock from the electrical equipment of the chamber, an annual repeat inspection as well as a test prior to initial startup and prior to resuming operation after maintenance or repair, are required. This test must meet the requirements of the competent public authorities. We recommend testing under DIN VDE 0701-0702:2008 in accordance with the details in the Service Manual.

For personnel requirements please refer to the Service Manual.

14.2 Maintenance intervals, service



Electrical hazard during live maintenance work.

Deadly electric shock.

- \oslash The chamber must NOT become wet during operation or maintenance works.
- \varnothing Do NOT remove the rear panel of the chamber.
- Disconnect the chamber before conducting maintenance work. Turn off the main power switch (ED, FD, FED size 720) and pull the power plug.
- Make sure that all maintenance work will be conducted by licensed electricians or experts authorized by BINDER.

Ensure regular maintenance work is performed at least once a year.

The warranty becomes void if maintenance work is conducted by non-authorized personnel.



Replace the door gasket only when cold. Otherwise, the door gasket may become damaged.

We recommend taking out a maintenance agreement. Please consult BINDER Service.

BINDER telephone hotline: +49 (0) 7462 2005 555 **BINDER** fax hotline: +49 (0) 7462 2005 93555 BINDER e-mail hotline: customerservice@binder-world.com BINDER service hotline USA: +1 866 885 9794 or +1 631 224 4340 x3 (toll-free in the USA) BINDER service hotline Asia Pacific: +852 390 705 04 or +852 390 705 03 BINDER service hotline Russia and CIS +7 495 988 15 16 **BINDER** Internet website http://www.binder-world.com **BINDER** address BINDER GmbH, post office box 102, D-78502 Tuttlingen

International customers, please contact your local BINDER distributor.

14.3 Simple troubleshooting

Defects and shortcomings can compromise the operational safety of the chamber and can lead to risks and damage to equipment and persons. If there are is a technical fault or shortcoming, take the chamber out of operation and inform BINDER Service. If you are not sure whether there is a technical fault, proceed according to the following list. If you cannot clearly identify an error or there is a technical fault, please contact BINDER Service.



Only qualified service personnel authorized by BINDER must perform repair. Repaired chambers must comply with the BINDER quality standards.

Fault description	Possible cause	Required measures			
General	General				
		Check connection to power supply.			
	No power supply.	ED, FD, FED size 720: Check whether the chamber is turned on at the main power switch.			
	Wrong voltage.	Check power supply for correct voltage (chap. 4.2).			
Chamber without function.	Chamber fuse has responded.	Check chamber fuse and replace it if appropriate (chap. 16.3). If it re- sponds again, contact BINDER ser- vice.			
	Defective controller.	Contact BINDER service.			
	Overtemperature protective device class 1 has turned off the chamber.	Disconnect the chamber from the power supply for at least 10 sec- onds and let it cool down. If the de- vice responds again, contact BINDER service			
Chamber without function, the standby icon is displayed.	Chamber in standby mode.	Press down the standby button until the display lights up.			
Temperature					
	Chamber door not properly closed.	Completely close chamber door.			
Set-point temperature is not	Defective door gasket.	Replace door gasket,			
reached after specified time.	Controller not adjusted.	Calibrate and adjust controller.			
	Wrong voltage.	Check the power supply for correct voltage (chap. 4.2).			
	Defective controller.				
Chamber heating permanently,	Defective Pt 100 sensor.	Contact BINDER service.			
set-point not held.	Defective semiconductor relay.				
	Controller not adjusted.	Calibrate and adjust controller.			



Fault description	Possible cause	Required measures				
Temperature (continued)						
Chamber doesn't heat up. Heat-	Defective heating element.	Contact BINDER service.				
ing icon is displayed.	Defective semiconductor relay.					
Chamber doesn't heat up. No	Timer run off.	Re-program the timer or turn it off.				
heating icon in the display. Con-	Defective semiconductor relay.	Contact BINDER service.				
troller display working.	Defective controller.	Contact BINDER Service.				
<i>BD, BF:</i> Alarm message "!TProt" is displayed	Safety device class 3.1 has re- sponded.	Check the settings of the tempera- ture set-point and of the safety de- vice class 3.1 (chap. 7.3).				
<i>ED, FD, FED:</i> Chamber without function. Alarm message "!TLim" is displayed	Safety device class 2 has turned off the chamber.	Let cool down the chamber. Check the settings of the temperature set- point and of the safety device class 2 (chap. 7.3). If appropriate, select suitable limit value.				
	Defective safety device class 2.	Contact BINDER service.				
Deviations from the indicated heating-up times.	Chamber fully loaded.	Load the chamber less or consider longer heating-up times.				
Controller						
Message "1999" in the controller display	Sensor rupture between sen- sor and controller.	Contact BINDER service.				
Fault description	Possible cause	Required measures				
Miscellaneous						
<i>BF, FD, FED:</i> The fan doesn't turn or turns too slowly.	<i>BF, FED:</i> Fan speed set too low	<i>BF, FED:</i> Set fan speed to 100%				
	Defective fan.	Contact BINDER service.				

14.4 Sending the chamber back to BINDER GmbH

If you return a BINDER product to us for repair or any other reason, we will only accept the product upon presentation of an **authorization number** (RMA number) that has previously been issued to you. An authorization number will be issued after receiving your complaint either in writing or by telephone **prior** to your sending the BINDER product back to us. The authorization number will be issued following receipt of the information below:

- BINDER product type and serial number
- Date of purchase
- Name and address of the dealer from which you bought the BINDER product
- Exact description of the defect or fault
- Complete address, contact person and availability of that person
- Exact location of the BINDER product in your facility
- A contamination clearance certificate (chap. 18) must be faxed in advance

The authorization number must be applied to the packaging in such a way that it can be easily recognized or be recorded clearly in the delivery documents.

For safety reasons we cannot accept a chamber delivery if it does not carry an authorization number.

Return address: BINDER GmbH, Abteilung Service Gänsäcker 16, 78502 Tuttlingen, Germany

15. Disposal

15.1 Disposal of the transport packing

Packing element	Material	Disposal
Straps to fix packing on pallet	Plastic	Plastic recycling
Wooden transport box (option)	Non-wood (compressed match- wood, IPPC standard)	Wood recycling
with metal screws	Metal	Metal recycling
Pallet (from size 115 on)	Solid wood (IPPC standard)	Wood recycling
Transport box	Cardboard	Paper recycling
with metal clamps	Metal	Metal recycling
Top cover (size 720)	Cardboard	Paper recycling
Edge protection	Styropor [®] or PE foam	Plastic recycling
Protection of doors and racks	PE foam	Plastic recycling
Bag for operating manual	PE foil	Plastic recycling
Insulating air cushion foil (packing of optional accessories)	PE foil	Plastic recycling

If recycling is not possible, all packing parts can also be disposed of with normal waste.

15.2 Decommissioning

- ED, FD, FED size 720: Turn off the chambers at the main power switch (chap. 2.3).
- Disconnect the chamber from the power supply (pull the power plug).
- With option inert gas connection (chap. 12.8): Close the inert gas supply and remove the gas connection.



- Temporal decommissioning: See indications for appropriate storage, chap. 3.3.
- Final decommissioning: Dispose of the chamber as described in chap. 15.3 to 15.5.

15.3 Disposal of the chamber in the Federal Republic of Germany

According to Annex I of Directive 2012/19/EU of the European Parliament and of the Council on waste electrical and electronic equipment (WEEE), BINDER devices are classified as "monitoring and control instruments" (category 9) only intended for professional use". They must not be disposed of at public collecting points.

The chambers bear the symbol for the marking of electrical and electronic equipment manufactured / placed on the market in the EC after 13 August 2005 and be disposed of in separate collection according to Directive 2012/19/EU on waste electrical and electronic equipment (WEEE) and German national law for electrical and electronic equipment (Elektro- und Elektronikgerätegesetz, ElektroG). WEEE marking: crossed-out wheeled bin with solid bar under. A significant part of the materials must be recycled in order to protect the environment.



BINDER

At the end of the device's service life, have the device disposed of according to the German national law for electrical and electronic equipment (Elektro- und Elektronikgerätegesetz, ElektroG from 20 October 2015, BGBI. I p. 1739) or contact BINDER service who will organize taking back and disposal of the chamber according to the German national law for electrical and electronic equipment (Elektro- und Elektronikgerätegesetz, Elektro- und Elektro- Elektro- und



NOTICE


Danger of violation against existing law if not disposed of properly. Failure to comply with applicable law.

- \varnothing Do NOT dispose of BINDER devices at public collecting points.
- Have the device disposed of professionally at a recycling company which is certified according to the German national law for electrical and electronic equipment (Elektround Elektronikgerätegesetz, ElektroG from 20 October 2015, BGBI. I p. 1739). or
- Instruct BINDER Service to dispose of the device. The general terms of payment and delivery of BINDER GmbH apply, which were valid at the time of purchasing the chamber.

Certified companies disassemble waste (used) BINDER equipment in primary substances for recycling according to Directive 2012/19/EU. The devices must be free from toxic, infectious or radioactive substances in order to eliminate any health hazards to the employees of the recycling companies.

Prior to handing the chamber over to a recycling company, it is the user's responsibility that it is free from toxic, infectious or radioactive substances.

- Prior to disposal, clean all introduced or residual toxic substances from the chamber.
- Prior to disposal, disinfect the chamber from all sources of infection. Be aware that sources of infection may also be located outside the inner chamber.
- If you cannot safely remove all toxic substances and sources of infection from the chamber, dispose of it as special waste according to national law.
- Fill out the contamination clearance certificate (chap. 18) and enclose it with the chamber.



A chamber from which all toxic substances or sources of infection cannot be safely removed must be considered as "special" waste according to national law. Dispose of it accordingly.


15.4 Disposal of the chamber in the member states of the EU except for the Federal Republic of Germany

According to Annex I of Directive 2012/19/EU of the European Parliament and of the Council on waste electrical and electronic equipment (WEEE), BINDER devices are classified as "monitoring and control instruments" (category 9) only intended for professional use". They must not be disposed of at public collecting points.

The chambers bear the symbol for the marking of electrical and electronic equipment manufactured / placed on the market in the EC after 13 August 2005 and be disposed of in separate collection according to the Directive 2012/19/EU on waste electrical and electronic equipment (WEEE). WEEE marking: crossed-out wheeled bin with solid bar under.



At the end of the device's service life, notify the distributor who sold you the device, who will take back and dispose of the chamber according to the Directive 2012/19/EU on waste electrical and electronic equipment (WEEE).




If your distributor is not able to take back and dispose of the chamber, please contact BINDER service.

Certified companies disassemble waste (used) BINDER equipment in primary substances for recycling according to Directive 2012/19/EU. The devices must be free from toxic, infectious or radioactive substances in order to eliminate any health hazards to the employees of the recycling companies.

Prior to handing the chamber over to a recycling company, it is the user's responsibility that it is free from toxic, infectious or radioactive substances.

- Prior to disposal, clean all introduced or residual toxic substances from the chamber.
- Prior to disposal, disinfect the chamber from all sources of infection. Be aware that sources of infection may also be located outside the inner chamber.
- If you cannot safely remove all sources of infection and toxic substances from the chamber, dispose of it as special waste according to national law.
- Fill out the contamination clearance certificate (chap. 18) and enclose it with the chamber.



WARNING

Danger of intoxication and infection through contamination of the chamber with toxic, infectious or radioactive substances.

Ø NEVER take a chamber contaminated with toxic substances or sources of infection for recycling according to Directive 2012/19/EU.

- Prior to disposal, remove all toxic substances and sources of infection from the chamber.
- A chamber from which all toxic substances or sources of infection cannot be safely removed must be considered as "special" waste according to national law. Dispose of it accordingly.

15.5 Disposal of the chamber in non-member states of the EU



NOTICE

Danger of violation against existing law if not disposed of properly. Failure to comply with applicable law. Alteration of the environment.

- For final decommissioning and disposal of the chamber, please contact BINDER service.
- > Follow the statutory regulations for appropriate, environmentally friendly disposal.

16. Technical description

16.1 Factory calibration and adjustment

This chamber was calibrated and adjusted in the factory. Calibration and adjustment were performed using standardized test instructions, according to the QM DIN EN ISO 9001 system applied by BINDER (certified since December 1996 by TÜV CERT). All test equipment used is subject to the administration of measurement and test equipment that is also constituent part of the BINDER QM DIN EN ISO 9001 systems. They are controlled and calibrated to a DKD-Standard at regular intervals.

16.2 Definition of usable volume

The usable volume illustrated below is calculated as follows:





Figure 12: Determination of the useable volume

The technical data refers to the defined usable volume.



16.3 Over current protection

Single-phase devices are protected by one (UL chambers) or two miniature fuses against over current, accessible from the outside. The miniature fuses are located at the rear of the chamber above the power cable connection. Each fuse holder is equipped with a fuse clip 5 mm x 20 mm (cUL-Version 6,3x32 mm). A fuse may be replaced only with a substitute of the same ratings. Refer to the technical data of the respective device type.

Two-phase devices are equipped with a resettable miniature circuit breaker (combination element).

Three-phase devices are equipped with an internal miniature circuit breaker.

16.4 BD technical data

Chamber size)		BD 56	BD 115	BD 260	BD 720
Exterior dime						-
Width, net		mm / inch	560 / 22.01	710 / 27.95	815 / 32.09	1165 / 45.87
Height, gross	(incl. feet)	mm / inch	625 / 24.60	735 / 28.94	965 / 37.99	1590 / 62.60
Depth, net	· · · · ·	mm / inch	565 / 22.24	605 / 23.82	760 / 29.92	816 / 32.13
	incl. door handle, ex-		0.40 / 05 00	000 / 00 77	045 / 00 00	070 / 0 / 05
haust duct)	,	mm / <i>inch</i>	640 / 25.20	680 / 26.77	815 / 32.09	870 / 34.25
Wall clearance	e rear (minimum)	mm / inch	160 / <i>6.30</i>	160 / <i>6.30</i>	160 / 6.30	160 / 6.30
Wall clearance	e side (minimum)	mm / <i>inch</i>	100 / 3.94	100 / 3.94	100 / 3.94	100 / 3.94
Exhaust duct,	outer diameter	mm / <i>inch</i>	52 / 2.05	52 / 2.05	52 / 2.05	52 / 2.05
Doors						
Number of doo	ors		1	1	1	2
Number of inn	er glass doors		1	1	1	2
Interior dime	nsions					
Width		mm / inch	360 / 14.17	510 / 20.08	610 / 24.02	960 / 37.80
Height		mm / inch	420 / 16.54	530 / 20.87	760 / 29.92	1280 / 50.39
Depth		mm / inch	380 / 14.96	420 / 16.53	550 / 21.65	605 / 23.81
Interior volume	9	l / cu.ft.	57 / 2.01	112 / 3.96	255 / 9.01	737 / 26.03
Steam space	volume	l / cu.ft.	63 / 2.22	127 / 4.49	279 / 9.85	791 / 27.93
Racks						
Quantity of rac	cks (regular)		2	2	2	2
Quantity of rac	cks (max.)		4	5	8	16
Max. load per	rack (standard rack)	Kg / Ibs	30 / 66	30 / 66	40 / 88	45 / 99
Max. load per (optional perfo		Kg / Ibs	35 / 77	35 / 77	35 / 77	35 / 77
Max. load per (optional heav		Kg / Ibs	50 / 110	70 / 154	80 / 176	70 / 154
Permitted tota	l load	Kg / Ibs	70 / 154	150 / 330	270 / 595	315 / 694
Weight						
Weight (empty	()	Kg / Ibs	38 / 84	54 / 119	85 / 187	170 / 375
Temperature	data					
Temperature range	from … degrees above ambient	°C / °F	5/9	5 / 9	5 / 9	5 / 9
Tange	up to	°C / °F	100 / 212	100 / 212	100 / 212	100 / 212
Temperature f at 37 °C / 98.6		± K	0.2	0.1	0.2	0.1
Temperature u at 37 °C / 98.6	uniformity (variation) S °F	±Κ	0.4	0.4	0.4	0.7
Heating up tim	ne to 37 °C / 98.6 °F	minutes	45	55	65	70
	e after door was sec at 37 °C / 98.6 °F	minutes	16	16	19	23
Electrical dat	a (model versions BD0)56-230V, E	BD115-230V, B	D260-230V, B	D720-230V)	
System of pro 60529	tection acc. to EN	IP	20	20	20	20
Nominal at 50) Hz power frequency	V	230	230	230	230
voltage) Hz power frequency	V	230	230	230	230
Current type			1N~	1N~	1N~	1N~
Nominal powe	r	kW	0.30	0.35	0.85	1.65
Nominal powe socket	r with option internal	kW	0.80	0.85	1.35	2.15



Chambe	r size		BD 56	BD 115	BD 260	BD 720
Electrica	I data (continued)		I		I	1
Power pl	ug of the power cable			Ground	ed plug	
Chamber fuse (external) 5x20 mm / A		6.3	6.3	8.0	12.5	
Overtem class 1	perature protective device	°C	120	120	120	120
Installatio	Installation category acc. to IEC 61010-1		II	II	II	II
Pollution	degree acc. to IEC 61010-1		2	2	2	2
	t electrical data for BD-UL ersions BD056UL-120V, BD				L-240V)	
Nominal	at 50 Hz power frequency	V	120	120	120	240
voltage (+/-10%)	at 60 Hz power frequency	V	120	120	120	240
Current type		1N~	1N~	1N~	2~	
Nominal	power	kW	0,30	0,35	0,95	1,75
Power pl	ug of the power cable	NEMA	5-15P	5-15P	5-15P	6-20P
	r fuse (external) i / 250V / time-lag (T)	А	12.5	12.5	12.5	
Miniature	e circuit breaker (internal)	А				16
Environ	nent-specific data					
Energy c 98.6 °F	onsumption at 37 °C /	Wh/h	25	25	40	78

All indications are average values, typical for chambers produced in series. We reserve the right to change technical specifications at any time.

F

If the chamber is fully loaded, the specified heating up times may vary according to the load.

With option interior socket: If electrical devices are connected and operating inside the chamber, the temperature range may be modified due to heat emission.

16.5 BF technical data

Chamber size		BF 56	BF 115	BF 260	BF 720
Exterior dimensions					
Width, net	mm / inch	560 / 22.01	710 / 27.95	815 / 32.09	1165 / <i>45</i> .87
Height, gross (incl. feet)	mm / inch	625 / 24.60	735 / 28.94	965 / 37.99	1590 / 62.60
Depth, net	mm / inch	565 / 22.24	605 / 23.82	760 / 29.92	816 / 32.13
Depth, gross (incl. door handle and exhaust duct)	mm / <i>inch</i>	640 / 25.20	680 / 26.77	815 / 32.09	870 / 34.25
Wall clearance rear (minimum)	mm / inch	160 / <i>6.30</i>	160 / <i>6.30</i>	160 / 6.30	160 / <i>6.30</i>
Wall clearance side (minimum)	mm / inch	100 / 3.94	100 / 3.94	100 / 3.94	100 / 3.94
Exhaust duct, outer diameter	mm / <i>inch</i>	52 / 2.05	52 / 2.05	52 / 2.05	52 / 2.05
Doors					
Number of doors		1	1	1	2
Number of inner glass doors		1	1	1	2
Interior dimensions					
Width	mm / <i>inch</i>	400 / 15.75	550 / 21.65	650 / 25.59	1000 / 39.37
Height	mm / inch	440 / 17.32	550 / 21.65	780 / 30.71	1300 / <i>51.18</i>
Depth	mm / inch	340 / 13.39	380 / 14.96	510 / 20.08	560 / 22.05



Chamber size)		BF 56	BF 115	BF 260	BF 720
	nsions (continued)			21		
Interior volume		l / cu.ft.	59 / 2.08	114 / 4.03	257 / 9.08	728 / 25.71
Steam space		/ cu.ft.	66 / 2.33	127 / 4.49	279 / 9.85	791 / 27.93
Racks		.,				
Quantity of rac	cks (regular)		2	2	2	2
Quantity of rac	, , ,		4	5	8	16
	rack (standard rack)	Kg / Ibs	30 / 66	30 / 66	40 / 88	45 / 99
Max. load per (optional perfo	rack	Kg / Ibs	35 / 77	35 / 77	35 / 77	35 / 77
Max. load per (optional heav		Kg / Ibs	50 / 110	70 / 154	80 / 176	70 / 154
Permitted tota		Kg / Ibs	70 / 154	150 / 330	270 / 595	315 / 694
Weight				L	L	
Weight (empty	()	Kg / Ibs	39 / 86	54 / 119	85 / 187	166 / 366
Temperature	data					
Temperature	from … degrees above ambient	°C / °F	7 / 12.6	8 / 14.4	7 / 12.6	10 / <i>18</i>
range	up to	°C / °F	100 / 212	100 / 212	100 / 212	100 / 212
Temperature f 98.6 °F	luctuation at 37 °C /	±Κ	0.1	0.1	0.1	0.1
Temperature u at 37 °C / 98.6	uniformity (variation) S °F	±Κ	0.3	0.3	0.3	0.3
Heating up tim	ne to 37 °C / 98.6 °F	minutes	8	8	8	15
	e after door was sec at 37 °C / 98.6 ° <i>F</i>	minutes	3	4	4	4
Electrical dat	a (model versions BF0	56-230V, B	F115-230V, B	F260-230V, BI	720-230V)	
System of pro 60529	tection acc. to EN	IP	20	20	20	20
voltage) Hz power frequency	V	230	230	230	230
(+/-10%) at 60) Hz power frequency	V	230	230	230	230
Current type			1N~	1N~	1N~	1N~
Nominal powe	r	kW	0.40	0.40	0.90	1.75
Nominal powe socket	r with option internal	kW	0.90	0.90	1.40	2.25
Power plug of	the power cable			Ground	ed plug	
Chamber fuse 5x20 mm / 250	(external) DV / time-lag (T)	А	6.3	6.3	8.0	12.5
Overtemperati class 1	ure protective device	°C	120	120	120	120
Installation cat	tegory acc. to IEC 610	10-1	I	II	II	II
Pollution degre	ee acc. to IEC 61010-1		2	2	2	2
	trical data for BF-UL ns BF056UL-120V, BF				-240V)	
Nominal at 50) Hz power frequency	V	120	120	120	240
voltage (+/-10%) at 60) Hz power frequency	V	120	120	120	240
Current type			1N~	1N~	1N~	2~
Power plug of	the power cable	NEMA	5-15P	5-15P	5-15P	6-20P
Nominal powe		kW	0.40	0.40	1.00	1.85
Chamber fuse 5x20 mm / 250	(external) DV / time-lag (T)	А	12.5	12.5	12.5	
Miniature circu	uit breaker (internal)	А				16



Chamber size		BF 56	BF 115	BF 260	BF 720
Environment-specific data					
Noise level (mean value)	dB (A)	43	43	43	43
Energy consumption at 37 °C / 98.6 °F	Wh/h	60	60	70	130

All indications are average values, typical for chambers produced in series. We reserve the right to change technical specifications at any time.

If the chamber is fully loaded, the specified heating up times may vary according to the load.

With option interior socket: If electrical devices are connected and operating inside the chamber, the temperature range may be modified due to heat emission.

16.6 ED technical data

Chamber size		ED 56	ED 115	ED 260	ED 720
Exterior dimensions					
Width, net	mm / <i>inch</i>	560 / 22.01	710 / 27.95	815 / 32.09	1165 / <i>45</i> .87
Height, gross (incl. feet)	mm / inch	625 / 24.60	735 / 28.94	965 / 37.99	1590 / 62.60
Depth, net	mm / inch	565 / 22.24	605 / 23.82	760 / 29.92	816 / 32.13
Depth, gross (incl. door handle and exhaust duct)	mm / <i>inch</i>	640 / 25.20	680 / 26.77	815 / 32.09	870 / 34.25
Wall clearance rear (minimum)	mm / <i>inch</i>	160 / <i>6.30</i>	160 / 6.30	160 / 6.30	160 / 6.30
Wall clearance side (minimum)	mm / <i>inch</i>	100 / 3.94	100 / 3.94	100 / 3.94	100 / 3.94
Exhaust duct, outer diameter	mm / inch	52 / 2.05	52 / 2.05	52 / 2.05	52 / 2.05
Doors					
Number of door(s)		1	1	1	2
Interior dimensions					
Width	mm / inch	360 / 14.17	510 / 20.08	610 / 24.02	960 / 37.80
Height	mm / inch	420 / 16.54	530 / 20.87	760 / 29.92	1280 / 50.39
Depth	mm / inch	380 / 14.96	425 / 16.73	550 / 21.65	610 / 24.02
Interior volume	I / cu.ft.	57 / 2.01	114 / <i>4</i> .03	255 / 9.01	749 / 26.45
Steam space volume	l / cu.ft.	63	127 / 4.49	273 / 9.64	791 / 27.93
Racks					
Quantity of racks (regular)		2	2	2	2
Quantity of racks (max.)		4	5	8	16
Max. load per rack (standard rack)	Kg / Ibs	30 / 66	30 / 66	44 / 88	45 / 99
Max. load per rack (optional perfo- rated shelf)	Kg / Ibs	35 / 77	35 / 77	35 / 77	35 / 77
Max. load per rack (optional heavy load rack)	Kg / Ibs	50 / 110	70 / 154	80 / 176	70 / 154
Permitted total load	Kg / Ibs	70 / 154	150 / 330	270 / 595	315 / 694
Weight					
Weight (empty)	Kg / Ibs	39 / 86	54 / 119	85 / 187	169 / 373



Chamber size			ED 56	ED 115	ED 260	ED 720
Temperature of	data					
Temperature	from degrees above ambient	°C / °F	5 / 9	5 / 9	5 / 9	5 / 9
range	up to	°C / °F	300 / 572	300 / 572	300 / 572	300 / 572
Temperature fl 302 °F	uctuation at 150 °C /	±Κ	0.5	0.4	1.0	0.8
Temperature u at 150 °C / 302	niformity (variation) ? ° <i>F</i>	±Κ	2.5	2.0	2.0	3,2
Heating up time	e to 150 °C / 302 °F	minutes	45	45	55	85
Recovery time opened for 30 : <i>302 °F</i>		minutes	35	16	20	25
Electrical data (model version	i s ED056-230V, ED11	5-230V, ED	260-230V, ED	0720-400V)		
IP system of pr 60529	otection acc. to EN	IP	20	20	20	20
Nominal at 50	Hz power frequency	V	230	230	230	400
voltage (+/-10%) at 60	Hz power frequency	V	230	230	230	400
Current type			1N~	1N~	1N~	3N~
Nominal power	-	kW	1.05	1.25	2.25	4.10
Power plug of t	the power cable			Ground	ed plug	
Chamber fuse 5x20 mm / 250	(external) V / time-lag (T)	А	6.3	6.3	12.5	
Miniature circu	it breaker (internal)	А				16
Overtemperatu class 1	re protective device	°C	330	330	330	330
Installation cate	egory acc. to IEC 610	10-1		II	II	
Pollution degre	e acc. to IEC 61010-1		2	2	2	2
	r ical data for ED-UL s ED056UL-120V, ED					
Nominal at 50	Hz power frequency	V	120	120	240	
voltage (+/-10%) at 60	Hz power frequency	V	120	120	240	
Current type			1N~	1N~	2~	
Power plug of t	the power cable	NEMA	5-15P	5-15P	6-20P	
Nominal power	-	kW	1.15	1.35	2.45	
Chamber fuse 5x20 mm / 250	(external) V / time-lag (T)	А	12.5	12.5		
Miniature circu (combination e	it breaker lement, external)	А			16	
Environment-	specific data					
Energy consun 302 °F	nption at 150 °C /	Wh/h	180	250	370	700

All indications are average values, typical for chambers produced in series. We reserve the right to change technical specifications at any time.

F

If the chamber is fully loaded, the specified heating up times may vary according to the load.

16.7 FD technical data

Chamber size		FD 56	FD 115	FD 260	FD 720
Exterior dimensions		1000	10110	10100	10120
Width, net	mm / inch	560 / 22.01	710 / 27.95	815 / 32.09	1165 / 45.87
Height, gross (incl. feet)	mm / inch	625 / 24.60	735 / 28.94	965 / 37.99	1590 / 62.60
Depth, net	mm / inch	565 / 22.24	605 / 23.82	760 / 29.92	816 / 32.13
Depth, gross (incl. door handle and					
exhaust duct)	mm / inch	640 / 25.20	680 / 26.77	815 / 32.09	870 / 34.25
Wall clearance rear (minimum)	mm / <i>inch</i>	160 / <i>6.30</i>	160 / <i>6.30</i>	160 / 6.30	160 / 6.30
Wall clearance side (minimum)	mm / <i>inch</i>	100 / 3.94	100 / 3.94	100 / 3.94	100 / 3.94
Exhaust duct, outer diameter	mm / <i>inch</i>	52 / 2.05	52 / 2.05	52 / 2.05	52 / 2.05
Doors					
Number of door(s)		1	1	1	2
Interior dimensions					
Width	mm / <i>inch</i>	400 / 15.75	550 / 21.65	650 / 25.59	1000 / 39.37
Height	mm / <i>inch</i>	440 / 17.32	550 / 21.65	780 / 30.71	1300 / <i>51.18</i>
Depth	mm / <i>inch</i>	345 / 13.58	385 / 15.16	510 / 20.08	570 / 22.44
Interior volume	/ cu.ft.	60 / 2.12	116 / <i>4.10</i>	257 / 9.08	741 / 26.17
Steam space volume	/ cu.ft.	67 / 2.37	127 / 4.49	279 / 9.85	791 / 27.93
Racks					
Quantity of racks (regular)		2	2	2	2
Quantity of racks (max.)		4	5	8	16
Max. load per rack (standard rack)	Kg / Ibs	30 / 66	30 / 66	44 / 88	45 / 99
Max. load per rack (optional perfo- rated shelf)	Kg / Ibs	35 / 77	35 / 77	35 / 77	35 / 77
Max. load per rack (optional heavy load rack)	Kg / Ibs	50 / 110	70 / 154	80 / 176	70 / 154
Permitted total load	Kg / Ibs	70 / 154	150 / 330	270 / 595	315 / 694
Weight					
Weight (empty)	Kg / Ibs	39 / 86	54 / 119	85 / 187	166 / 366
Temperature data					
Temperature from degrees above ambient	°C / °F	10 / <i>18</i>	10 / <i>18</i>	10 / <i>18</i>	12 / 21.6
range up to	°C / °F	300 / 572	300 / 572	300 / 572	300 / 572
Temperature fluctuation at 150 °C / 302 °F	± K	0.3	0.3	0.5	0.6
Temperature uniformity (variation) at 150 °C / <i>302 °F</i>	± K	1.7	1.7	1.9	2.5
Heating up time to 150 °C / 302 °F	minutes	15	19	20	25
Recovery time after door was opened for 30 sec at 150 °C / 302 °F	minutes	4	5	6	6
Ventilation data					
Air change at 100 °C / 212 °F	x/h	80	32	9	5
Electrical data (model versions FD056-230V, FD11	•	•			
,	5-230V, FD	260-230V, FD	720-400V)		
System of protection acc. to EN 60529	5-230V, FD	260-230V, FD 20	720-400V) 20	20	20
				20 230	20 400



Chambe	r size		FD 56	FD 115	FD 260	FD 720
	II data (continued) ersions FD056-230V, FD11	5-230V, FD	260-230V, FD	720-400V)		
Current t	уре		1N~	1N~	1N~	3N~
Nominal	power	kW	1.10	1.30	2.30	4.50
Power pl	ug of the power cable			Ground	ed plug	
	fuse (external) / 250V / time-lag (T)	А	6.3	6.3	12.5	
Miniature	circuit breaker (internal)	А				16
Overtem class 1	perature protective device	°C	330	350	330	330
Installatio	on category acc. to IEC 610	10-1		II	II	
Pollution	degree acc. to IEC 61010-1		2	2	2	2
	electrical data for FD-UL ersions FD056UL-120V, FD					
Nominal	at 50 Hz power frequency	V	120	120	240	
voltage (+/-10%)	at 60 Hz power frequency	V	120	120	240	
Current t	уре		1N~	1N~	2~	
Nominal	power	kW	1.20	1.40	2.50	
Power pl	ug of the power cable	NEMA	5-15P	5-15P	6-20P	
	⁻ fuse (external) n / 250V / time-lag (T)	А	12.5	12.5		
Miniature tion elem	circuit breaker (combina- ent)	А			16	
Environ	nent-specific data					
Noise lev	rel (mean value)	dB (A)	43	43	43	43
Energy c 302 °F	onsumption at 150 °C /	Wh/h	300	340	420	800

All indications are average values, typical for chambers produced in series. We reserve the right to change technical specifications at any time.

If the chamber is fully loaded, the specified heating up times may vary according to the load.

16.8 FED technical data

Chamber size		FED 56	FED 115	FED 260	FED 720
Exterior dimensions					
Width, net	mm / <i>inch</i>	560 / 22.01	710 / 27.95	815 / 32.09	1165 / <i>45</i> .87
Height, gross (incl. feet)	mm / inch	625 / 24.60	735 / 28.94	965 / 37.99	1590 / 62.60
Depth, net	mm / inch	565 / 22.24	605 / 23.82	760 / 29.92	816 / 32.13
Depth, gross (incl. door handle and exhaust duct)	mm / <i>inch</i>	640 / 25.20	680 / 26.77	815 / 32.09	870 / 34.25
Wall clearance rear (minimum)	mm / <i>inch</i>	160 / <i>6.30</i>	160 / <i>6.30</i>	160 / 6.30	160 / 6.30
Wall clearance side (minimum)	mm / inch	100 / 3.94	100 / 3.94	100 / 3.94	100 / 3.94
Exhaust duct, outer diameter	mm / <i>inch</i>	52 / 2.05	52 / 2.05	52 / 2.05	52 / 2.05
Doors					
Number of door(s)		1	1	1	2



Chamber size Interior dimensions Width Height Depth Interior volume Steam space volume Racks Quantity of racks (regular Quantity of racks (regular Quantity of racks (max.) Max. load per rack (stand Max. load per rack (option rated shelf) Max. load per rack (option load rack) Permitted total load Weight Weight (empty) Temperature data Temperature fluctuation a 302 °F from above a up to	lard rack) nal perfo- nal heavy	mm / inch mm / inch I / cu.ft. I / cu.ft. Kg / lbs Kg / lbs Kg / lbs Kg / lbs	FED 56 400 / 15.75 440 / 17.32 345 / 13.58 60 / 2.12 67 / 2.37 2 2 4 30 / 66 35 / 77 50 / 110	FED 115 550 / 21.65 550 / 21.65 385 / 15.16 116 / 4.10 127 / 4.49 2 5 30 / 66 35 / 77 70 / 154	FED 260 650 / 25.59 780 / 30.71 510 / 20.08 257 / 9.08 279 / 9.85 2 2 8 44 / 88 35 / 77 20 / 470	FED 720 1000 / 39.37 1300 / 51.18 570 / 22.44 741 / 26.17 791 / 27.93 2 16 45 / 99 35 / 77
Width Height Depth Interior volume Steam space volume Racks Quantity of racks (regular Quantity of racks (max.) Max. load per rack (stand Max. load per rack (option rated shelf) Max. load per rack (option load rack) Permitted total load Weight Weight (empty) Temperature data Temperature data Temperature fluctuation a 302 °F	lard rack) nal perfo- nal heavy	mm / inch mm / inch I / cu.ft. I / cu.ft. Kg / lbs Kg / lbs Kg / lbs	440 / 17.32 345 / 13.58 60 / 2.12 67 / 2.37 2 4 30 / 66 35 / 77 50 / 110	550 / 21.65 385 / 15.16 116 / 4.10 127 / 4.49 2 5 30 / 66 35 / 77	780 / 30.71 510 / 20.08 257 / 9.08 279 / 9.85 2 2 8 44 / 88 35 / 77	1300 / 51.18 570 / 22.44 741 / 26.17 791 / 27.93 2 16 45 / 99
Height Depth Interior volume Steam space volume Racks Quantity of racks (regular Quantity of racks (max.) Max. load per rack (stand Max. load per rack (option rated shelf) Max. load per rack (option load rack) Permitted total load Weight Weight (empty) Temperature data Temperature data Temperature fluctuation a 302 °F	lard rack) nal perfo- nal heavy	mm / inch mm / inch I / cu.ft. I / cu.ft. Kg / lbs Kg / lbs Kg / lbs	440 / 17.32 345 / 13.58 60 / 2.12 67 / 2.37 2 4 30 / 66 35 / 77 50 / 110	550 / 21.65 385 / 15.16 116 / 4.10 127 / 4.49 2 5 30 / 66 35 / 77	780 / 30.71 510 / 20.08 257 / 9.08 279 / 9.85 2 2 8 44 / 88 35 / 77	1300 / 51.18 570 / 22.44 741 / 26.17 791 / 27.93 2 16 45 / 99
Depth Interior volume Steam space volume Racks Quantity of racks (regular Quantity of racks (max.) Max. load per rack (stand Max. load per rack (option rated shelf) Max. load per rack (option load rack) Permitted total load Weight Weight (empty) Temperature data Temperature data Temperature fluctuation a 302 °F	lard rack) nal perfo- nal heavy	mm / inch / cu.ft. / cu.ft. Kg / lbs Kg / lbs Kg / lbs	345 / 13.58 60 / 2.12 67 / 2.37 2 4 30 / 66 35 / 77 50 / 110	385 / 15.16 116 / 4.10 127 / 4.49 2 5 30 / 66 35 / 77	510 / 20.08 257 / 9.08 279 / 9.85 279 / 9.85 2 8 44 / 88 35 / 77	570 / 22.44 741 / 26.17 791 / 27.93 2 16 45 / 99
Interior volume Steam space volume Racks Quantity of racks (regular Quantity of racks (max.) Max. load per rack (stand Max. load per rack (option rated shelf) Max. load per rack (option load rack) Permitted total load Weight Weight (empty) Temperature data Temperature data Temperature fluctuation a 302 °F	lard rack) nal perfo- nal heavy	/ cu.ft. / cu.ft. Kg / lbs Kg / lbs Kg / lbs	60 / 2.12 67 / 2.37 2 4 30 / 66 35 / 77 50 / 110	116 / 4.10 127 / 4.49 2 5 30 / 66 35 / 77	257 / 9.08 279 / 9.85 2 8 44 / 88 35 / 77	741 / 26.17 791 / 27.93 2 16 45 / 99
Steam space volume Racks Quantity of racks (regular Quantity of racks (max.) Max. load per rack (stand Max. load per rack (option rated shelf) Max. load per rack (option Max. load per rack (option rated shelf) Max. load per rack (option load rack) Permitted total load Weight Weight (empty) Temperature data range from above a up to Temperature fluctuation a 302 °F	lard rack) nal perfo- nal heavy	I / cu.ft. Kg / lbs Kg / lbs Kg / lbs	67 / 2.37 2 4 30 / 66 35 / 77 50 / 110	127 / 4.49 2 5 30 / 66 35 / 77	279 / 9.85 2 8 44 / 88 35 / 77	791 / 27.93 2 16 45 / 99
Racks Quantity of racks (regular Quantity of racks (max.) Max. load per rack (stand Max. load per rack (option rated shelf) Max. load per rack (option Ioad rack) Permitted total load Weight Weight (empty) Temperature data range from above a up to Temperature fluctuation a 302 °F	lard rack) nal perfo- nal heavy	Kg / Ibs Kg / Ibs Kg / Ibs	2 4 30 / 66 35 / 77 50 / 110	2 5 30 / 66 35 / 77	2 8 44 / 88 35 / 77	2 16 45 / 99
Quantity of racks (regular Quantity of racks (max.) Max. load per rack (stand Max. load per rack (option rated shelf) Max. load per rack (option load rack) Permitted total load Weight Weight (empty) Temperature data Temperature data Temperature from above a up to Temperature fluctuation a 302 °F	lard rack) nal perfo- nal heavy	Kg / Ibs Kg / Ibs	4 30 / 66 35 / 77 50 / 110	5 30 / 66 35 / 77	8 44 / 88 35 / 77	16 45 / 99
Quantity of racks (max.) Max. load per rack (stand Max. load per rack (option rated shelf) Max. load per rack (option load rack) Permitted total load Weight Weight (empty) Temperature data Temperature data Temperature form above a up to Temperature fluctuation a 302 °F	lard rack) nal perfo- nal heavy	Kg / Ibs Kg / Ibs	4 30 / 66 35 / 77 50 / 110	5 30 / 66 35 / 77	8 44 / 88 35 / 77	16 45 / 99
Max. load per rack (stand Max. load per rack (option rated shelf) Max. load per rack (option load rack) Permitted total load Weight Weight (empty) Temperature data Temperature data Temperature from above a up to Temperature fluctuation a 302 °F	nal perfo- nal heavy	Kg / Ibs Kg / Ibs	30 / 66 35 / 77 50 / 110	30 / 66 35 / 77	44 / 88 35 / 77	45 / 99
Max. load per rack (option rated shelf) Max. load per rack (option load rack) Permitted total load Weight Weight (empty) Temperature data Temperature data Temperature from above a up to Temperature fluctuation a 302 °F	nal perfo- nal heavy	Kg / Ibs Kg / Ibs	35 / 77 50 / 110	35 / 77	35 / 77	
rated shelf) Max. load per rack (option load rack) Permitted total load Weight Weight (empty) Temperature data Temperature data range Temperature fluctuation a 302 °F	nal heavy	Kg / Ibs	50 / 110			35 / 77
load rack) Permitted total load Weight Weight (empty) Temperature data Temperature range Temperature fluctuation a 302 °F		-		70 / 154	00/470	
Weight Weight (empty) Temperature data Temperature range Image Temperature fluctuation a 302 °F		Kg / Ibs		-	80 / 176	70 / 154
Weight (empty) Temperature data Temperature range Temperature fluctuation a 302 °F			70 / 154	150 / 330	270 / 595	315 / 694
Temperature dataTemperature rangefrom above a up toTemperature fluctuation a 302 °F						
Temperature range Temperature fluctuation a 302 °F	_	Kg / Ibs	39 / 86	54 / 119	85 / 187	162 / 357
Temperature above a range up to Temperature fluctuation a 302 °F	-					•
Temperature fluctuation a 302 °F	degrees	°C / °F	10 / 18	10 / <i>18</i>	10 / <i>18</i>	12 / 21.6
Temperature fluctuation a 302 °F	ambient	077	10778	10778	10778	12/21.0
302 [°] F		°C / °F	300 / 572	300 / 572	300 / 572	300 / 572
Tomporature	at 150 °C /	±Κ	0.3	0.3	0.5	0.6
Temperature uniformity (v at 150 °C / 302 °F	variation)	±Κ	1.4	1.2	1.6	2.0
Heating up time to 150 °C	C/302 °F	minutes	15	19	20	25
Recovery time after door opened for 30 sec at 150 ° <i>F</i>		minutes	4	5	6	6
Ventilation data						•
Air change at 100 °C	/ 212 °F	x/h	80	32	9	5
Electrical data (model versions FED056-	-230V, FED	115-230V,	FED260-230V	, FED720-400 ^v	√)	
System of protection acc. 60529	to EN	IP	20	20	20	20
Nominal at 50 Hz power	frequency	V	230	230	230	400
voltage (+/-10%) at 60 Hz power	frequency	V	230	230	230	400
Current type			1N~	1N~	1N~	3N~
Nominal power		kW	1.10	1.30	2.30	4.50
Power plug of the power	cable			Ground	ed plug	•
Chamber fuse (external) 5x20 mm / 250V / time-la	g (T)	А	6.3	6.3	12.5	
Miniature circuit breaker ((internal)	А				16
Overtemperature protecti class 1	,	°C	330	350	330	330
Installation category acc.	to IEC 610	10-1	II	II	II	11
Pollution degree acc. to II			2	2	2	2



Chambe	r size		FED 56	FED 115	FED 260	FED 720
	t electrical data for FED-UI ersions FED056UL-120V, F					
Nominal	at 50 Hz power frequency	V	120	120	240	208
voltage (+/-10%)	at 60 Hz power frequency	V	120	120	240	208
Current t	уре		1N~	1N~	2~	3N~
Nominal	power	kW	1.20	1.40	2.50	4.50
Power pl	ug of the power cable	NEMA	5-15P	5-15P	6-20P	L21-20P
	r fuse (external) n / 250V / time-lag (T)	А	12.5	12.5		
	e circuit breaker (combina- ient, external)	А			16	
Miniature	e circuit breaker (internal)	А				16
Environ	ment-specific data					
Noise lev	vel (mean value)	dB (A)	43	43	43	43
Energy c 302 °F	onsumption at 150 °C /	Wh/h	300	340	420	800

All indications are average values, typical for chambers produced in series. We reserve the right to change technical specifications at any time.

If the chamber is fully loaded, the specified heating up times may vary according to the load.

16.9 Equipment and options (extract)



To operate the chamber, use only original BINDER accessories or accessories / components from third-party suppliers authorized by BINDER. The user is responsible for any risk arising from using unauthorized accessories.

	BD	BF	ED	FD	FED
Standard equipment			20		
Microprocessor temperature controller	\checkmark		\checkmark	\checkmark	\checkmark
One timer function: Delayed Off	\checkmark		\checkmark	\checkmark	
Three timer functions: Delayed On, Delayed Off and Temperature dependent Delayed Off		\checkmark			\checkmark
Adjustable ramp function	\checkmark		\checkmark	\checkmark	\checkmark
Temperature safety controller class 3.1 acc. to DIN 12880:2007	\checkmark	\checkmark			
Temperature safety controller class 2 acc. to DIN 12880:2007			\checkmark	\checkmark	\checkmark
Inner glass door	\checkmark				
USB interface to read out the measured values	\checkmark		\checkmark	\checkmark	\checkmark
Communication interface Ethernet					\checkmark
Exhaust duct, internal diameter 50 mm / <i>1.97 inches</i> , with ad- justable ventilation slide	\checkmark	\checkmark	\checkmark	\checkmark	\checkmark
Adjustable air change by means of rear exhaust duct (50 mm)	\checkmark	\checkmark	\checkmark	\checkmark	\checkmark
2 racks, chrome-plated	\checkmark		\checkmark	\checkmark	\checkmark

B	IN	D	Ε	R

Options / accessories		l i			
Rack, chrome-plated or stainless steel					\checkmark
Perforated rack, stainless steel	\checkmark	\checkmark	\checkmark	\checkmark	\checkmark
Heavy load rack, stainless steel	\checkmark	\checkmark	\checkmark	\checkmark	\checkmark
Access ports with various diameters, with silicone plug	\checkmark	\checkmark		\checkmark	\checkmark
Door with window			\checkmark	\checkmark	\checkmark
Interior lightning			\checkmark	\checkmark	\checkmark
Communication interface Ethernet	\checkmark				
Battery backed real-time clock	\checkmark				\checkmark
Rubber pads for safe stacking (5 pieces)	\checkmark				\checkmark
Object temperature display with additional Pt100 temperature sensor	\checkmark	\checkmark	\checkmark	\checkmark	\checkmark
Disconnectable water protected internal socket (CA3GD) for line voltage inside the chamber, IP type of protection 67, 230 V 1N \sim 50-60 Hz. Max. load 500 W	\checkmark	\checkmark			
Analog output 4-20 mA for temperature with 6 pole DIN socket, DIN plug included	\checkmark	\checkmark	\checkmark	\checkmark	\checkmark
HEPA Fresh air filter, class H 14 (DIN EN 1822:2009)				\checkmark	\checkmark
Mostly gas-tight version				\checkmark	\checkmark
Inert gas connection (gas inlet and outlet), with mostly gas- tight version		\checkmark		\checkmark	\checkmark
Disconnectable audible over-temperature alarm	\checkmark			\checkmark	\checkmark
FKM door gasket (temperature resistant up to 200 °C / 392 ° <i>F</i> max.)			\checkmark	\checkmark	\checkmark
Measurement of air change rate acc. to ASTM D5374				\checkmark	\checkmark
Factory calibration certificate	\checkmark	\checkmark		\checkmark	\checkmark
Extension to factory calibration certificate (additional values)	\checkmark	\checkmark		\checkmark	\checkmark
Measuring protocol acc. to DIN 12880:2007	\checkmark	\checkmark		\checkmark	\checkmark
Qualification folder	\checkmark	\checkmark		\checkmark	\checkmark
Neutral cleaning agent (liquid concentrate)	\checkmark	\checkmark		\checkmark	\checkmark
Stable table on wheels with castors and locking brakes					\checkmark

16.10 Accessories and spare parts (extract)

BINDER GmbH is responsible for the safety features of the chamber only, provided skilled electricians or qualified personnel authorized by BINDER perform all maintenance and repair, and if components relating to chamber safety are replaced in the event of failure with original spare parts. The user is responsible for any risks arising from using unauthorized accessories/components.

Chamber size	56	115	260	720
Description		Art.	No.	
Rack, chrome-plated BD, ED	6004-0174	6004-0175	6004-0177	6004-0179
Rack, chrome-plated BF, FD, FED	6004-0166	6004-0167	6004-0169	6004-0171
Rack, stainless steel BD, ED	6004-0158	6004-0159	6004-0161	6004-0163
Rack, stainless steel BF, FD, FED	6004-0150	6004-0151	6004-0153	6004-0155
Perforated rack, stainless steel BD, ED	6004-0190	6004-0191	6004-0193	6004-0195
Perforated rack, stainless steel BF, FD, FED	6004-0182	6004-0183	6004-0185	6004-0187
Rack, heavy load, stainless steel BD, ED	6004-0201	6004-0202	6004-0203	6004-0205
Rack, heavy load, stainless steel BF, FD, FED	6004-0198	6004-0199	6004-0200	6004-0204



Chamber size	56	115	260	720
Description		Art.	No.	
Door gasket, silicone	6005-0254	6005-0255	6005-0258	6005-0260
Door gasket made of FKM (temperature re- sistant up to 200 °C / 392 °F max.)	6005-0265	6005-0266	6005-0268	6005-270
Stable table on wheels with castors and locking brakes	9051-0005	9051-0005	9051-0006	
Chamber fuse 5x20mm 250V 6,3 A time lag (T)	5006-0092	5006-0092		
Chamber fuse 5x20mm 250V 8,0 A time lag (T)			5006-0093	
Chamber fuse 5x20mm 250V 12,5 A time lag (T)	5006-0096	5006-0096	5006-0096	

Description	Art. No.
HEPA fresh air filter (replacement) class H 14 (DIN EN 1822:2009)	6014-0003
Rubber pads for safe stacking (5 pieces)	8012-1887
Neutral cleaning agent, 1 kg	1002-0016

For information on components not listed here, please contact BINDER Service.


Validation service	Art. no.
Qualification folder IQ-OQ (printed version)	7007-0001
Qualification folder IQ-OQ (digital version)	7057-0001
Qualification folder IQ-OQ-PQ (printed version)	7007-0005
Qualification folder IQ-OQ-PQ (digital version)	7057-0005
Execution of IQ-OQ	DL420300
Execution of IQ-OQ-PQ	DL440500

Chamber type	BD	BF	ED	FD	FED
Calibration service			ArtNo.		
Calibration of temperature including certifi- cate (1 measuring point)	DL300101	DL300101	DL300101	DL300101	DL300101
Spatial temperature measurement including certificate (9 measuring points)	DL300109	DL300109	DL300109	DL300109	DL300109
Spatial temperature measurement including certificate (18 measuring points)	DL300118	DL300118	DL300118	DL300118	DL300118
Spatial temperature measurement including certificate (27 measuring points)	DL300127	DL300127	DL300127	DL300127	DL300127
Measurement of air ventilation acc. to ASTM D 5374, including certificate			DL330000	DL330000	DL330000

38


56,5

16.11 Dimensions size 56




120
















38

57





16.13 Dimensions size 260







16.14 Dimensions size 720







BD / BF / ED / FD / FED (E3.1) 02/2022



17. Certificates and declarations of conformity

17.1 EU Declaration of Conformity for BD

Binder State Binder State Вилона EU-Konformitätserklärung / EU Declaration of Conformity / Déclaration de UE / Dichiarazione di conformità UE / Де Contract sun EU Hersteller / Manufacturer / Fabricant / Fabricante / Fabricante / Производитель BINDER GmbH Anschrift / Address / Adresse / Dirección / Indirizzo / Agpec Im Mittleren Ösch 5, 78532 Tuttlingen, Incubators with natural convection Incubators with natural convection Incubators with natural convection Incubators de convección natural Incubatoria a convection naturale Инкубаторы с естественной конвекци Typenbezeichnung / Type / Type / Tipo / Tipo / Tim BD 56, BD 115, BD 260, BD 720 (E3.1 9010-0325, 9010-0326, 91190-0325, 9010-0326, 911900-0325, 9110-0325, 9010-0326, 911900-0325, 9110-0325, 9010-0326, 911900-0331, 9110-0331, 9010-0332, 9110-0331, 9010-0332, 9110-0331, 9010-0332, 9110-0331, 9010-0332, 9110-0331, 9010-0332, 9110-0332, 9010-0332, 9110-0332, 9010-0332, 9110-0332, 9010-0332, 9110-0332, 9010-0332, 9100-0332, 9110-0332, 9010-0332, 9110-0332, 9010-0332, 9110-0332, 9010-0332, 9110-0332, 9010-0332, 9110-0332, 9010-0332, 9110-0332, 9010-0332, 9110-0332, 9010-0332, 9110-0332, 9010-0332, 9110-0332, 9010-0332, 9110-0332, 9010-0332, 9110-0332, 9010-0332, 9110-0332, 9010-0332, 9110-0332, 9010-0332, 9110-0332, 9010-0332, 9110-0332, 9010-0332, 9110-0332, 9010-0332, 9110-0332, 9010-0332, 9110-0332, 9010-0332, 9110-0332, 9010-0332, 9110-0332, 9010-0332, 9110-0332, 9010-0332, 9110-0332, 9010-0332, 9110-0332, 9010-0332, 9110-0332, 9010-0332, 9110-0332, 9110-0332, 9110-0332, 9110-0332, 9110-0332, 9110-0332, 9110-0332, 9110-0332, 9110-0332, 9110-0332, 9110-0332, 9110-0332, 9110-0332, 9110-0332, 9110-0332, 9110-0332, 9110-0332, 9110-0332, 9110-0332, 9110-0332, 9110-0332, 9110-0332, 9110-0332, 9110-0332, 9110-0332, 911	Germany
UE / Declaración de conformidad UE / Dichiarazione di conformità UE / Де cooтветствия EU Негsteller / Manufacturer / Fabricant / Fabricante / Производитель Anschrift / Address / Adresse / Dirección / Indirizzo / Aqpec BINDER GmbH Produkt / Product / Produit / Producto / Prodotto / Продукт Im Mittleren Ösch 5, 78532 Tuttlingen, Incubators with natural convection Incubateurs à convection naturelle Incubadoras de convección natural Incubatori a convection naturelle Incubadoras de convección naturale Инкубаторы с естественной конвекц Typenbezeichnung / Type / Type / Tipo / Tipo / Tun BD 56, BD 115, BD 260, BD 720 (E3.1) Art. No. / Art. no. / Réf. / Art. Nº / Art. n. / № арт. 9010-0323, 9110-0323, 9010-0326, 91 9010-0329, 9110-0329, 9010-0330, 91	Germany
UE / Declaración de conformidad UE / Dichiarazione di conformità UE / Де cooтветствия EU Негsteller / Manufacturer / Fabricant / Fabricante / Fabricante / Производитель Anschrift / Address / Adresse / Dirección / Indirizzo / Aqpec BINDER GmbH Produkt / Product / Produit / Producto / Prodotto / Продукт Im Mittleren Ösch 5, 78532 Tuttlingen, Incubators with natural convection Incubators with natural convection Incubators de convección naturale Incubatori a convección naturale Инкубаторы с естественной конвекц Typenbezeichnung / Type / Type / Tipo / Tipo / Tun BD 56, BD 115, BD 260, BD 720 (E3.1) Art. No. / Art. no. / Réf. / Art. Nº / Art. n. / № арт. 9010-0323, 9110-0323, 9010-0324, 94 9010-0329, 9110-0329, 9010-0330, 91 9010-0329, 9110-0329, 9010-0330, 91	Germany
UE / Declaración de conformidad UE / Dichiarazione di conformità UE / Де cooтветствия EU Негsteller / Manufacturer / Fabricant / Fabricante / Fabricante / Производитель Anschrift / Address / Adresse / Dirección / Indirizzo / Aqpec BINDER GmbH Produkt / Product / Produit / Producto / Prodotto / Продукт Im Mittleren Ösch 5, 78532 Tuttlingen, Incubators with natural convection Incubators with natural convection Incubators de convección naturale Incubatori a convección naturale Инкубаторы с естественной конвекц Typenbezeichnung / Type / Type / Tipo / Tipo / Tun BD 56, BD 115, BD 260, BD 720 (E3.1) Art. No. / Art. no. / Réf. / Art. Nº / Art. n. / № арт. 9010-0323, 9110-0323, 9010-0324, 94 9010-0329, 9110-0329, 9010-0330, 91 9010-0329, 9110-0329, 9010-0330, 91	Germany
UE / Declaración de conformidad UE / Dichiarazione di conformità UE / Де cooтветствия EU Негsteller / Manufacturer / Fabricant / Fabricante / Fabricante / Производитель Anschrift / Address / Adresse / Dirección / Indirizzo / Aqpec BINDER GmbH Produkt / Product / Produit / Producto / Prodotto / Продукт Im Mittleren Ösch 5, 78532 Tuttlingen, Incubators with natural convection Incubateurs à convection naturelle Incubadoras de convección natural Incubatori a convección naturale Инкубаторы с естественной конвекц Typenbezeichnung / Type / Type / Tipo / Tipo / Tun BD 56, BD 115, BD 260, BD 720 (E3.1) Art. No. / Art. no. / Réf. / Art. Nº / Art. n. / № арт. 9010-0323, 9110-0323, 9010-0324, 94 9010-0329, 9110-0329, 9010-0330, 91 9010-0329, 9110-0329, 9010-0330, 91	Germany
Fabbricante / Производитель Anschrift / Address / Adresse / Dirección / Indirizzo / Aqpec Im Mittleren Ösch 5, 78532 Tuttlingen, Produkt / Product / Produit / Producto / Prodotto / Продукт Inkubatoren mit freier Konvektion Incubators with natural convection Incubators de convección naturale Incubatori a convección naturale Инкубаторы с естественной конвекц Typenbezeichnung / Type / Type / Tipo / Tipo / Tип BD 56, BD 115, BD 260, BD 720 (E3.1 9010-0323, 9110-0323, 9010-0324, 91 9010-0329, 9110-0329, 9010-0330, 91	
Fabbricante / Производитель Anschrift / Address / Adresse / Dirección / Indirizzo / Aqpec Im Mittleren Ösch 5, 78532 Tuttlingen, Produkt / Product / Produit / Producto / Prodotto / Продукт Inkubatoren mit freier Konvektion Incubators with natural convection Incubators de convección natural Incubatori a convección naturale Инкубаторы с естественной конвекц Typenbezeichnung / Type / Type / Tipo / Tipo / Tип BD 56, BD 115, BD 260, BD 720 (E3.1 9010-0323, 9110-0323, 9010-0324, 91 9010-0329, 9110-0329, 9010-0330, 91	
Адрес Produkt / Product / Produit / Producto / Prodotto / Продукт Inkubatoren mit freier Konvektion Incubators with natural convection Incubators a convection naturale Incubatori a convection naturale Инкубаторы с естественной конвекц Typenbezeichnung / Type / Type / Tipo / Timo / Tim BD 56, BD 115, BD 260, BD 720 (E3.1 9010-0323, 9110-0323, 9010-0324, 91 9010-0325, 9110-0325, 9010-0326, 91 9010-0329, 9110-0329, 9010-0330, 91	
Продукт Incubators with natural convection Incubateurs à convection naturale Incubatori a convección natural Incubatori a convección naturale Инкубаторы с естественной конвекц Турельеzеichnung / Type / Type / Tipo / Тип BD 56, BD 115, BD 260, BD 720 (E3.1 Art. No. / Art. no. / Réf. / Art. № / Art. n. / № арт. 9010-0323, 9110-0323, 9010-0324, 91 9010-0325, 9110-0325, 9010-0326, 91 9010-0329, 9110-0329, 9010-0330, 91	ией
Incubateurs à convection naturelle Incubadoras de convección natural Incubatori a convezione naturale Инкубаторы с естественной конвекц Typenbezeichnung / Type / Type / Tipo / Тип BD 56, BD 115, BD 260, BD 720 (E3.1 BD 56, BD 115, BD 260, BD 720 (E3.1 Art. No. / Art. no. / Réf. / Art. N° / Art. n. / № арт. 9010-0323, 9110-0323, 9010-0324, 91 9010-0325, 9110-0325, 9010-0326, 91 9010-0329, 9110-0329, 9010-0330, 91	ией
Incubatori a convezione naturale Инкубаторы с естественной конвекц Турепьеzeichnung / Type / Type / Tipo / Тип BD 56, BD 115, BD 260, BD 720 (E3.1 Art. No. / Art. no. / Réf. / Art. N° / Art. n. / № арт. 9010-0323, 9110-0323, 9010-0324, 91 9010-0325, 9110-0325, 9010-0326, 91 9010-0329, 9110-0329, 9010-0330, 91	ией
Инкубаторы с естественной конвекц Турепbezeichnung / Туре / Туре / Тіро / Тіпо / Тип BD 56, BD 115, BD 260, BD 720 (E3.1 Art. No. / Art. no. / Réf. / Art. № / Art. n. / № арт. 9010-0323, 9110-0323, 9010-0324, 91 9010-0325, 9110-0325, 9010-0326, 91 9010-0329, 9110-0329, 9010-0330, 91	ией
Art. No. / Art. no. / Réf. / Art. N° / Art. n. / № арт. 9010-0323, 9110-0323, 9010-0324, 91 9010-0325, 9110-0325, 9010-0326, 91 9010-0325, 9110-0325, 9010-0326, 91 9010-0329, 9110-0329, 9010-0330, 91 9010-0329, 9110-0329, 9010-0330, 91	
9010-0325, 9110-0325, 9010-0326, 91 9010-0329, 9110-0329, 9010-0330, 91)
9010-0329, 9110-0329, 9010-0330, 91	10-0324
	10-0326
9010-0331, 9110-0331, 9010-0332, 91	
	10-0332
Die oben beschriebenen Produkte sind konform mit folgenden EU-Richtlinien:	
The products described above are in conformity with the following EU Directives:	
Les produits décrits ci-dessus sont conformes aux directives UE suivantes:	
Los productos descritos arriba cumplen con las siguientes directivas de la UE:	
I prodotti sopra descritti sono conformi alle seguenti direttive UE: Продукты, указанные выше, полностью соответствуют следующим EU руководствам:	
продукты, указанные выше, полностые соответствуют следующим со руководствам.	
• 2014/35/EU	
Niederspannungsrichtlinie 2014/35/EU / Low voltage directive 2014/35/EU / Directive b 2014/35/UE / Directiva sobre baja tensión 2014/35/UE / Direttiva Bassa tensione 2014/ Директива по низкому напряжению 2014/35/EU	
• 2014/30/EU	
EMV-Richtlinie 2014/30/EU / EMC Directive 2014/30/EU / Directive CEM 2014/30/UE / 2014/30/UE / Direttiva EMC 2014/30/UE / Директива ЭМС 2014/30/EU	Directiva CEI
• 2011/65/EU, (EU) 2015/863	
RoHS-Richtlinien 2011/65/EU und (EU) 2015/863 / RoHS Directives 2011/65/EU and (Directives RoHS 2011/65/UE et (UE) 2015/863 / Directivas RoHS 2011/65/UE y (UE) 2 artiting RoHS 2011/65/UE et (UE) 2015/863 / Directivas RoHS 2011/65/UE y (EU) 2015/863 / Directivas RoHS 2015/863 / Directivas ROHS 2015/863 / Directivas ROHS 2015/863 / Dir	015/863 / Di-
rettive RoHS 2011/65/UE et (UE) 2015/863 / Директивы RoHS 2011/65/EU и (EU) 207	5/803
1/2	
1 / 2 BINDER GmbH Postfach 102 D-78502 Tuttlingen Anschrift: BINDER GmbH Im Mittleren Ösch 5 D-78532 Tuttlingen	



Best conditions for your success

Die oben beschriebenen Produkte tragen entsprechend die Kennzeichnung CE. The products described above, corresponding to this, bear the CE-mark. Les produits décrits ci-dessus, en correspondance, portent l'indication CE. Los productos descritos arriba, en conformidad, llevan la indicación CE. I prodotti sopra descritti, conformi a quanto sopra, portano il marchio CE. Данные продукты в соответствии с изложенным выше маркированы знаком CE.

Die oben beschriebenen Produkte sind konform mit folgenden harmonisierten Normen: The products described above are in conformity with the following harmonized standards: Les produits décrits ci-dessus sont conformes aux normes harmonisées suivantes: Los productos descritos arriba cumplen con las siguientes normas: I prodotti sopra descritti sono conformi alle seguenti normative armonizzate: Продукты, указанные выше, полностью соответствуют следующим стандартам:

Sicherheit / Safety / Sécurité / Seguridad / Sicurezza / Нормативы по безопасности

- EN 61010-1:2010
- EN 61010-2-010:2014

EMV / EMC / CEM / CEM / EMC / ЭМС

• EN 61326-1:2013

RoHS

EN IEC 63000:2018

78532 Tuttlingen, 28.01.2022 BINDER GmbH

110.

P. Wimmer Vice President Vice President Vice président Vicepresidente vicepresidente Bице-президент



Глава департамента R&D

2/2

 BINDER GmbH
 Postfach 102
 D-78502 Tuttlingen
 Anschrift: BINDER
 GmbH
 Im Mittleren
 Ösch 5
 D-78532 Tuttlingen

 Kontakt:
 Telefon: +49 (0) 74 62 / 20 05 - 0
 |
 Telefax: +49 (0) 74 62 / 20 05 - 100
 info@binder-world.com
 |
 www.binder-world.com

 Geschäftsführung:
 Dipl-Ing. Peter M. Binder
 |
 Amtsgericht Stuttgart, HRB 727150
 |
 Stz der Gesellschaft: Tuttlingen

 Bankverbindung:
 Kreissparkasse Tuttlingen
 IBAN-Code: DE05643 500700 000002266
 |
 SWIFT-Code: SOLA DE S1TUT

 Deutsche Bank Tuttlingen
 IBAN-Code: DE56653 70075 0213870900
 |
 SWIFT-Code: DEUT DE \$\$653

 Altgeräte-Entsorgung gemäß WEEE-Reg.-Nr. DE 37004983

 SWIFT-Code: DEUT DE \$\$653

17.2 EU Declaration of Conformity for BF

	Best conditions for your s
СС EU-Konformitätserklärung / EU Decl. UE / Declaración de conformidad UE соответствия EU	aration of Conformity / Déclaration de conformité : / Dichiarazione di conformità UE / Декларация
Hersteller / Manufacturer / Fabricant / Fabricante / Fabbricante / Производитель	BINDER GmbH
Anschrift / Address / Adresse / Dirección / Indirizzo / Адрес	Im Mittleren Ösch 5, 78532 Tuttlingen, Germany
Produkt / Product / Produit / Producto / Prodotto / Продукт	Inkubatoren mit Umluft Incubators with forced convection Incubateurs à convection forcée Incubadoras de convección forzada Incubatori a convezione forzata Инкубаторы с принудительной циркуляцией воздуха
Typenbezeichnung / Type / Type / Tipo / Tipo / Тип	BF 56, BF 115, BF 260, BF 720 (E3.1)
Art. No. / Art. no. / Réf. / Art. № / Art. n. / № apτ.	9010-0313, 9110-0313, 9010-0314, 9110-0314 9010-0315, 9110-0315, 9010-0316, 9110-0316 9010-0319, 9110-0319, 9010-0320, 9110-0320
Die oben beschriebenen Produkte sind konform mit The products described above are in conformity with Les produits décrits ci-dessus sont conformes aux d Los productos descritos arriba cumplen con las sigu	n the following EU Directives: lirectives UE suivantes:
The products described above are in conformity with Les produits décrits ci-dessus sont conformes aux d Los productos descritos arriba cumplen con las sigu I prodotti sopra descritti sono conformi alle seguenti Продукты, указанные выше, полностью соответс • 2014/35/EU Niederspannungsrichtlinie 2014/35/EU / Low vo	folgenden EU-Richtlinien: n the following EU Directives: lirectives UE suivantes: lientes directivas de la UE: direttive UE: твуют следующим EU руководствам: Ditage directive 2014/35/EU / Directive basse tension /35/UE / Direttiva Bassa tensione 2014/35/UE /
The products described above are in conformity with Les produits décrits ci-dessus sont conformes aux d Los productos descritos arriba cumplen con las sigu l prodotti sopra descritti sono conformi alle seguenti Продукты, указанные выше, полностью соответс • 2014/35/EU Niederspannungsrichtlinie 2014/35/EU / Low vo 2014/35/UE / Directiva sobre baja tensión 2014 Директива по низкому напряжению 2014/35/EU • 2014/30/EU	folgenden EU-Richtlinien: n the following EU Directives: lirectives UE suivantes: lientes directivas de la UE: direttive UE: твуют следующим EU руководствам: oltage directive 2014/35/EU / Directive basse tension /35/UE / Direttiva Bassa tensione 2014/35/UE / EU 14/30/EU / Directive CEM 2014/30/UE / Directiva CEM
The products described above are in conformity with Les produits décrits ci-dessus sont conformes aux d Los productos descritos arriba cumplen con las sigu l prodotti sopra descritti sono conformi alle seguenti Продукты, указанные выше, полностью соответс • 2014/35/EU Niederspannungsrichtlinie 2014/35/EU / Low vo 2014/35/UE / Directiva sobre baja tensión 2014 Директива по низкому напряжению 2014/35/E • 2014/30/EU EMV-Richtlinie 2014/30/EU / EMC Directive 201 2014/30/UE / Direttiva EMC 2014/30/UE / Дире • 2011/65/EU, (EU) 2015/863 RoHS-Richtlinien 2011/65/EU und (EU) 2015/86	folgenden EU-Richtlinien: n the following EU Directives: lirectives UE suivantes: uientes directivas de la UE: direttive UE: твуют следующим EU руководствам: oltage directive 2014/35/EU / Directive basse tension /35/UE / Direttiva Bassa tensione 2014/35/UE / EU 14/30/EU / Directive CEM 2014/30/UE / Directiva CEM ектива ЭМС 2014/30/EU 63 / RoHS Directives 2011/65/EU and (EU) 2015/863 / Di-
The products described above are in conformity with Les produits décrits ci-dessus sont conformes aux d Los productos descritos arriba cumplen con las sigu l prodotti sopra descritti sono conformi alle seguenti Продукты, указанные выше, полностью соответс • 2014/35/EU Niederspannungsrichtlinie 2014/35/EU / Low vo 2014/35/UE / Directiva sobre baja tensión 2014 Директива по низкому напряжению 2014/35/E • 2014/30/EU EMV-Richtlinie 2014/30/EU / EMC Directive 2012 2014/30/UE / Direttiva EMC 2014/30/UE / Дире • 2011/65/EU, (EU) 2015/863 RoHS-Richtlinien 2011/65/EU und (EU) 2015/863 rettive RoHS 2011/65/UE et (UE) 2015/863 / Ди	folgenden EU-Richtlinien: n the following EU Directives: lirectives UE suivantes: uientes directivas de la UE: direttive UE: твуют следующим EU руководствам: oltage directive 2014/35/EU / Directive basse tension /35/UE / Direttiva Bassa tensione 2014/35/UE / EU 14/30/EU / Directive CEM 2014/30/UE / Directiva CEM 2014/30/EU / Directive CEM 2014/30/UE / Directiva CEM 2014/30/EU / Directives 2011/65/EU and (EU) 2015/863 / Directivas RoHS 2011/65/EU и (EU) 2015/863
The products described above are in conformity with Les produits décrits ci-dessus sont conformes aux d Los productos descritos arriba cumplen con las sigu l prodotti sopra descritti sono conformi alle seguenti Продукты, указанные выше, полностью соответс • 2014/35/EU Niederspannungsrichtlinie 2014/35/EU / Low vo 2014/35/UE / Directiva sobre baja tensión 2014 Директива по низкому напряжению 2014/35/E • 2014/30/EU EMV-Richtlinie 2014/30/EU / EMC Directive 2012 2014/30/UE / Direttiva EMC 2014/30/UE / Дире • 2011/65/EU, (EU) 2015/863 RoHS-Richtlinien 2011/65/EU und (EU) 2015/863 rettive RoHS 2011/65/UE et (UE) 2015/863 / Ди	folgenden EU-Richtlinien: n the following EU Directives: lirectives UE suivantes: uientes directivas de la UE: direttive UE: твуют следующим EU руководствам: oltage directive 2014/35/EU / Directive basse tension /35/UE / Direttiva Bassa tensione 2014/35/UE / EU 14/30/EU / Directive CEM 2014/30/UE / Directiva CEM ектива ЭМС 2014/30/EU 63 / RoHS Directives 2011/65/EU and (EU) 2015/863 / Di-



Best conditions for your success

Die oben beschriebenen Produkte tragen entsprechend die Kennzeichnung CE. The products described above, corresponding to this, bear the CE-mark. Les produits décrits ci-dessus, en correspondance, portent l'indication CE. Los productos descritos arriba, en conformidad, llevan la indicación CE. I prodotti sopra descritti, conformi a quanto sopra, portano il marchio CE. Данные продукты в соответствии с изложенным выше маркированы знаком CE.

Die oben beschriebenen Produkte sind konform mit folgenden harmonisierten Normen: The products described above are in conformity with the following harmonized standards: Les produits décrits ci-dessus sont conformes aux normes harmonisées suivantes: Los productos descritos arriba cumplen con las siguientes normas: I prodotti sopra descritti sono conformi alle seguenti normative armonizzate: Продукты, указанные выше, полностью соответствуют следующим стандартам:

Sicherheit / Safety / Sécurité / Seguridad / Sicurezza / Нормативы по безопасности

- EN 61010-1:2010
- EN 61010-2-010:2014

EMV/EMC/CEM/CEM/EMC/ЭMC

• EN 61326-1:2013

RoHS

• EN IEC 63000:2018

78532 Tuttlingen, 28.01.2022 BINDER GmbH

11ª. 18

P. Wimmer Vice President Vice President Vice président Vicepresidente исергезidente Вице-президент



Director R & D Chef de service R&D Responsable I & D Direttore R & D Глава департамента R&D

2/2

 BINDER GmbH
 Postfach 102
 D-78502 Tuttlingen
 Anschrift: BINDER
 GmbH
 Im Mittleren
 Ösch 5
 D-78532 Tuttlingen

 Kontakt:
 Telefon: +49 (0) 74 62 / 20 05 - 0
 |
 Telefax: +49 (0) 74 62 / 20 05 - 10
 info@binder-world.com
 |
 www.binder-world.com

 Geschäftsführung:
 Dijl-Ing. Peter M. Binder
 |
 Amtsgericht Stuttgart, HRB 727150
 |
 Siz der Gesellschaft: Tuttlingen

 Bankverbindung:
 Kreissparkasse Tuttlingen
 IBAN-Code: DE05643 500700 000002266
 |
 SWIFT-Code: SOLA DE S1TUT

 Deutsche Bank Tuttlingen
 IBAN-Code: DE56653 70075 0213870090
 |
 SWIFT-Code: DEUT DE \$\$653

 Altgeräte-Entsorgung gemäß WEEE-Reg.-Nr. DE 37004983
 |
 SWIFT-Code: DEUT DE \$\$653

17.3 EU Declaration of Conformity for ED

	Best conditions for your s
EU-Konformitätserklärung / EU UE / Declaración de conformida соответствия EU	Declaration of Conformity / Déclaration de conformité ad UE / Dichiarazione di conformità UE / Декларация
Hersteller / Manufacturer / Fabricant / Fabri- cante / Fabbricante / Производитель	BINDER GmbH
Anschrift / Address / Adresse / Dirección / In- dirizzo / Agpec	Im Mittleren Ösch 5, 78532 Tuttlingen, Germany
Produkt / Product / Produit / Producto / Pro- dotto / Продукт	Trocken- und Wärmeschränke mit freier Konvektion Drying and heating ovens with natural convection Etuves de chauffage et de séchage à convection naturelle Estufas de secado y calentamiento de convección natural Stufe per essiccazione e riscaldamento a convezione natural rale Сушильные и сухожаровые шкафы с естественной конвекцией
Туреnbezeichnung / Туре / Туре / Тіро / Тіро / Тип	ED 56, ED 115, ED 260, ED 720 (E3.1)
Art. No. / Art. no. / Réf. / Art. № / Art. n. / № арт.	9010-0333, 9110-0333, 9010-0334, 9110-0334 9010-0335, 9110-0335, 9010-0336, 9110-0336 9010-0339, 9110-0339, 9010-0340, 9110-0340 9010-0341, 9110-0341
Die oben beschriebenen Produkte sind konfor	m mit folgenden EU-Richtlinien:
The products described above are in conformi Les produits décrits ci-dessus sont conformes Los productos descritos arriba cumplen con la l prodotti sopra descritti sono conformi alle seg Продукты, указанные выше, полностью соо • 2014/35/EU Niederspannungsrichtlinie 2014/35/EU / L 2014/35/UE / Directiva sobre baja tensión Директива по низкому напряжению 201 • 2014/30/EU EMV-Richtlinie 2014/30/EU / EMC Directi 2014/30/UE / Direttiva EMC 2014/30/UE / • 2011/65/EU, (EU) 2015/863 RoHS-Richtlinien 2011/65/EU und (EU) 20 Directives RoHS 2011/65/UE et (UE) 2015	m mit folgenden EU-Richtlinien: ty with the following EU Directives: aux directives UE suivantes: s siguientes directivas de la UE: guenti direttive UE: тветствуют следующим EU руководствам: ow voltage directive 2014/35/EU / Directive basse tension 2014/35/UE / Direttiva Bassa tensione 2014/35/UE / 4/35/EU ve 2014/30/EU / Directive CEM 2014/30/UE / Directiva CEM
The products described above are in conformi Les produits décrits ci-dessus sont conformes Los productos descritos arriba cumplen con la l prodotti sopra descritti sono conformi alle seg Продукты, указанные выше, полностью соо • 2014/35/EU Niederspannungsrichtlinie 2014/35/EU / L 2014/35/UE / Directiva sobre baja tensión Директива по низкому напряжению 201 • 2014/30/EU EMV-Richtlinie 2014/30/EU / EMC Directi 2014/30/UE / Direttiva EMC 2014/30/UE / • 2011/65/EU, (EU) 2015/863 RoHS-Richtlinien 2011/65/EU und (EU) 20 Directives RoHS 2011/65/UE et (UE) 2015	m mit folgenden EU-Richtlinien: ty with the following EU Directives: aux directives UE suivantes: s siguientes directivas de la UE: guenti direttive UE: тветствуют следующим EU руководствам: ow voltage directive 2014/35/EU / Directive basse tension 2014/35/UE / Direttiva Bassa tensione 2014/35/UE / 4/35/EU ve 2014/30/EU / Directive CEM 2014/30/UE / Directiva CEM / директива ЭМС 2014/30/EU 015/863 / RoHS Directives 2011/65/EU and (EU) 2015/863 5/863 / Directivas RoHS 2011/65/UE y (UE) 2015/863 / Di-



Best conditions for your success

Die oben beschriebenen Produkte tragen entsprechend die Kennzeichnung CE. The products described above, corresponding to this, bear the CE-mark. Les produits décrits ci-dessus, en correspondance, portent l'indication CE. Los productos descritos arriba, en conformidad, llevan la indicación CE. I prodotti sopra descritti, conformi a quanto sopra, portano il marchio CE. Данные продукты в соответствии с изложенным выше маркированы знаком CE.

Die oben beschriebenen Produkte sind konform mit folgenden harmonisierten Normen: The products described above are in conformity with the following harmonized standards: Les produits décrits ci-dessus sont conformes aux normes harmonisées suivantes: Los productos descritos arriba cumplen con las siguientes normas: I prodotti sopra descritti sono conformi alle seguenti normative armonizzate: Продукты, указанные выше, полностью соответствуют следующим стандартам:

Sicherheit / Safety / Sécurité / Seguridad / Sicurezza / Нормативы по безопасности

- EN 61010-1:2010
- EN 61010-2-010:2014

EMV / EMC / CEM / CEM / EMC / ЭМС

• EN 61326-1:2013

RoHS

• EN IEC 63000:2018

78532 Tuttlingen, 28.01.2022 BINDER GmbH

110. 15

P. Wimmer Vice President Vice President Vice président Vicepresidente исергезidente Вице-президент



Letter F & E Director R & D Chef de service R&D Responsable I & D Direttore R & D Глава департамента R&D

2/2

 BINDER GmbH
 Postfach 102
 D-78502 Tuttlingen
 Anschrift: BINDER
 GmbH
 Im Mittleren
 Ösch 5
 D-78532 Tuttlingen

 Kontakt:
 Telefon: +49 (0) 74 62 / 20 05 - 0
 |
 Telefax: +49 (0) 74 62 / 20 05 - 10
 info@binder-world.com
 |
 www.binder-world.com

 Geschäftsführung:
 Dijl-Ing. Peter M. Binder
 |
 Amtsgericht Stuttgart, HRB 727150
 |
 Siz der Gesellschaft: Tuttlingen

 Bankverbindung:
 Kreissparkasse Tuttlingen
 IBAN-Code: DE05643 500700 000002266
 |
 SWIFT-Code: SOLA DE S1TUT

 Deutsche Bank Tuttlingen
 IBAN-Code: DE56653 70075 0213870090
 |
 SWIFT-Code: DEUT DE \$\$653

 Altgeräte-Entsorgung gemäß WEEE-Reg.-Nr. DE 37004983
 |
 SWIFT-Code: DEUT DE \$\$653

17.4 EU Declaration of Conformity for FD

UE / Declaración de conformida	Declaration of Conformity / Déclaration de conformité ad UE / Dichiarazione di conformità UE / Декларация
соответствия EU	
Hersteller / Manufacturer / Fabricant / Fabri- cante / Fabbricante / Производитель	BINDER GmbH
Anschrift / Address / Adresse / Dirección / In- dirizzo / Адрес	Im Mittleren Ösch 5, 78532 Tuttlingen, Germany
Produkt / Product / Produit / Producto / Pro- dotto / Продукт	Trocken- und Wärmeschränke mit Umluft Drying and heating ovens with forced convection Etuves de chauffage et de séchage à convection forcée Estufas de secado y calentamiento de convección forzada Stufe per essiccazione e riscaldamento a convezione for- zata Сушильные и сухожаровые шкафы с принудительной конвекцией
Typenbezeichnung / Type / Type / Tipo / Tipo / Тип	FD 56, FD 115, FD 260, FD 720 (E3.1)
Art. No. / Art. no. / Réf. / Art. № / Art. n. / № арт.	9010-0303, 9110-0303, 9010-0304, 9110-0304 9010-0305, 9110-0305, 9010-0306, 9110-0306 9010-0309, 9110-0309, 9010-0310, 9110-0310 9010-0311, 9110-0311
2014/35/UE / Directiva sobre baja tensión Директива по низкому напряжению 201	ty with the following EU Directives: aux directives UE suivantes: s siguientes directivas de la UE: guenti direttive UE: тветствуют следующим EU руководствам: ow voltage directive 2014/35/EU / Directive basse tension 2014/35/UE / Direttiva Bassa tensione 2014/35/UE /
The products described above are in conformi Les produits décrits ci-dessus sont conformes Los productos descritos arriba cumplen con la l prodotti sopra descritti sono conformi alle seg Продукты, указанные выше, полностью соог • 2014/35/EU Niederspannungsrichtlinie 2014/35/EU / L 2014/35/UE / Directiva sobre baja tensión Директива по низкому напряжению 2014 • 2014/30/EU	ty with the following EU Directives: aux directives UE suivantes: s siguientes directivas de la UE: guenti direttive UE: тветствуют следующим EU руководствам: ow voltage directive 2014/35/EU / Directive basse tension 2014/35/UE / Direttiva Bassa tensione 2014/35/UE / 4/35/EU ve 2014/30/EU / Directive CEM 2014/30/UE / Directiva CEM
The products described above are in conformi Les produits décrits ci-dessus sont conformes Los productos descritos arriba cumplen con la l prodotti sopra descritti sono conformi alle seg Продукты, указанные выше, полностью соог • 2014/35/EU Niederspannungsrichtlinie 2014/35/EU / L 2014/35/UE / Directiva sobre baja tensión Директива по низкому напряжению 201- • 2014/30/EU EMV-Richtlinie 2014/30/EU / EMC Directir 2014/30/UE / Direttiva EMC 2014/30/UE / • 2011/65/EU, (EU) 2015/863 RoHS-Richtlinien 2011/65/EU und (EU) 2015/ Directives RoHS 2011/65/UE et (UE) 2015/	ty with the following EU Directives: aux directives UE suivantes: s siguientes directivas de la UE: guenti direttive UE: тветствуют следующим EU руководствам: ow voltage directive 2014/35/EU / Directive basse tension 2014/35/UE / Direttiva Bassa tensione 2014/35/UE / 4/35/EU ve 2014/30/EU / Directive CEM 2014/30/UE / Directiva CEM
The products described above are in conformi Les produits décrits ci-dessus sont conformes Los productos descritos arriba cumplen con la l prodotti sopra descritti sono conformi alle seg Продукты, указанные выше, полностью соог • 2014/35/EU Niederspannungsrichtlinie 2014/35/EU / L 2014/35/UE / Directiva sobre baja tensión Директива по низкому напряжению 201- • 2014/30/EU EMV-Richtlinie 2014/30/EU / EMC Directir 2014/30/UE / Direttiva EMC 2014/30/UE / • 2011/65/EU, (EU) 2015/863 RoHS-Richtlinien 2011/65/EU und (EU) 2015/ Directives RoHS 2011/65/UE et (UE) 2015/	ty with the following EU Directives: aux directives UE suivantes: s siguientes directivas de la UE: guenti direttive UE: тветствуют следующим EU руководствам: ow voltage directive 2014/35/EU / Directive basse tension 2014/35/UE / Direttiva Bassa tensione 2014/35/UE / 4/35/EU ve 2014/30/EU / Directive CEM 2014/30/UE / Directiva CEI директива ЭМС 2014/30/EU 015/863 / RoHS Directives 2011/65/EU and (EU) 2015/863 5/863 / Directivas RoHS 2011/65/UE y (UE) 2015/863 / Di-



Best conditions for your success

Die oben beschriebenen Produkte tragen entsprechend die Kennzeichnung CE. The products described above, corresponding to this, bear the CE-mark. Les produits décrits ci-dessus, en correspondance, portent l'indication CE. Los productos descritos arriba, en conformidad, llevan la indicación CE. I prodotti sopra descritti, conformi a quanto sopra, portano il marchio CE. Данные продукты в соответствии с изложенным выше маркированы знаком CE.

Die oben beschriebenen Produkte sind konform mit folgenden harmonisierten Normen: The products described above are in conformity with the following harmonized standards: Les produits décrits ci-dessus sont conformes aux normes harmonisées suivantes: Los productos descritos arriba cumplen con las siguientes normas: I prodotti sopra descritti sono conformi alle seguenti normative armonizzate: Продукты, указанные выше, полностью соответствуют следующим стандартам:

Sicherheit / Safety / Sécurité / Seguridad / Sicurezza / Нормативы по безопасности

- EN 61010-1:2010
- EN 61010-2-010:2014

EMV / EMC / CEM / CEM / EMC / ЭМС

• EN 61326-1:2013

RoHS

• EN IEC 63000:2018

78532 Tuttlingen, 28.01.2022 BINDER GmbH

110. 60

P. Wimmer Vice President Vice President Vice président Vicepresidente исергезidente Вице-президент



Глава департамента R&D

2/2

 BINDER GmbH
 Postfach 102
 D-78502 Tuttlingen
 Anschrift: BINDER
 GmbH
 Im Mittleren
 Ösch 5
 D-78532 Tuttlingen

 Kontakt:
 Telefon: +49 (0) 74 62 / 20 05 - 0
 |
 Telefax: +49 (0) 74 62 / 20 05 - 100
 info@binder-world.com
 |
 www.binder-world.com

 Geschäftsführung:
 Dijl-Ing. Peter M. Binder
 |
 Amtsgericht Stuttgart, HRB 727150
 |
 Sitz der Gesellschaft: Tuttlingen

 Bankverbindung:
 Kreissparkasse Tuttlingen
 IBAN-Code: DE05643 500700 000002266
 |
 SWIFT-Code: SOLA DE S1TUT

 Deutsche Bank Tuttlingen
 IBAN-Code: DE56633 70075 0213870000
 |
 SWIFT-Code: DEUT DE \$\$653

 Altgeräte-Entsorgung gemäß WEEE-Reg.-Nr. DE 37004983
 |
 SWIFT-Code: DEUT DE \$\$653

17.5 EU Declaration of Conformity for FED

	BINDE Best conditions for your s
	EU Declaration of Conformity / Déclaration de conformité nidad UE / Dichiarazione di conformità UE / Декларация
Hersteller / Manufacturer / Fabricant / Fabricante / Fabbricante / Производитель	BINDER GmbH
Anschrift / Address / Adresse / Dirección / Indirizzo / Адрес	Im Mittleren Ösch 5, 78532 Tuttlingen, Germany
Produkt / Product / Produit / Producto / Prodotto / Продукт	Trocken- und Wärmeschränke mit Umluft Drying and heating ovens with forced convection Etuves de chauffage et de séchage à convection forcée Estufas de secado y calentamiento de convección forzada Stufe per essiccazione e riscaldamento a convezione forzata Сушильные и сухожаровые шкафы с принудительной конвекцией
Туреnbezeichnung / Туре / Туре / Тіро / Тіро / Тип	FED 56, FED 115, FED 260, FED 720 (E3.1)
Art. No. / Art. no. / Réf. / Art. № / Art. n. / № арт.	9010-0293, 9110-0293, 9010-0294, 9110-0294 9010-0295, 9110-0295, 9010-0296, 9110-0296 9010-0299, 9110-0299, 9010-0300, 9110-0300 9010-0301, 9110-0301, 9010-0302, 9110-0302
	form mit folgonden ELL Dichtlinion
Die oben beschriebenen Produkte sind kon The products described above are in confo Les produits décrits ci-dessus sont conform Los productos descritos arriba cumplen cor I prodotti sopra descritti sono conformi alle Продукты, указанные выше, полностью с	rmity with the following EU Directives: nes aux directives UE suivantes: n las siguientes directivas de la UE:
The products described above are in confo Les produits décrits ci-dessus sont conform Los productos descritos arriba cumplen cor l prodotti sopra descritti sono conformi alle Продукты, указанные выше, полностью с • 2014/35/EU	rmity with the following EU Directives: nes aux directives UE suivantes: n las siguientes directivas de la UE: seguenti direttive UE:
The products described above are in confo Les produits décrits ci-dessus sont conform Los productos descritos arriba cumplen cor l prodotti sopra descritti sono conformi alle Продукты, указанные выше, полностью с • 2014/35/EU Niederspannungsrichtlinie 2014/35/EU	rmity with the following EU Directives: nes aux directives UE suivantes: n las siguientes directivas de la UE: seguenti direttive UE: соответствуют следующим EU руководствам: / Low voltage directive 2014/35/EU / Directive basse tension ión 2014/35/UE / Direttiva Bassa tensione 2014/35/UE /
The products described above are in confo Les produits décrits ci-dessus sont conform Los productos descritos arriba cumplen cor l prodotti sopra descritti sono conformi alle Продукты, указанные выше, полностью с • 2014/35/EU Niederspannungsrichtlinie 2014/35/EU 2014/35/UE / Directiva sobre baja tens Директива по низкому напряжению 2 • 2014/30/EU EMV-Richtlinie 2014/30/EU / EMC Directiva	rmity with the following EU Directives: hes aux directives UE suivantes: h las siguientes directivas de la UE: seguenti direttive UE: соответствуют следующим EU руководствам: / Low voltage directive 2014/35/EU / Directive basse tension ión 2014/35/UE / Direttiva Bassa tensione 2014/35/UE / 2014/35/EU
The products described above are in confo Les produits décrits ci-dessus sont conform Los productos descritos arriba cumplen cor l prodotti sopra descritti sono conformi alle Продукты, указанные выше, полностью с • 2014/35/EU Niederspannungsrichtlinie 2014/35/EU 2014/35/UE / Directiva sobre baja tens Директива по низкому напряжению 2 • 2014/30/EU	rmity with the following EU Directives: hes aux directives UE suivantes: h las siguientes directivas de la UE: seguenti direttive UE: соответствуют следующим EU руководствам: / Low voltage directive 2014/35/EU / Directive basse tension ión 2014/35/UE / Direttiva Bassa tensione 2014/35/UE / 2014/35/EU
 The products described above are in confo Les produits décrits ci-dessus sont conform Los productos descritos arriba cumplen cor l prodotti sopra descritti sono conformi alle Продукты, указанные выше, полностью с 2014/35/EU Niederspannungsrichtlinie 2014/35/EU 2014/35/UE / Directiva sobre baja tens Директива по низкому напряжению 2 2014/30/EU EMV-Richtlinie 2014/30/EU / EMC Dire 2014/30/UE / Direttiva EMC 2014/30/U 2011/65/EU, (EU) 2015/863 RoHS-Richtlinien 2011/65/EU und (EU Directives RoHS 2011/65/UE et (UE) 2 	rmity with the following EU Directives: hes aux directives UE suivantes: h las siguientes directivas de la UE: seguenti direttive UE: соответствуют следующим EU руководствам: / Low voltage directive 2014/35/EU / Directive basse tension ión 2014/35/UE / Direttiva Bassa tensione 2014/35/UE / 2014/35/EU
 The products described above are in confo Les produits décrits ci-dessus sont conform Los productos descritos arriba cumplen cor l prodotti sopra descritti sono conformi alle Продукты, указанные выше, полностью с 2014/35/EU Niederspannungsrichtlinie 2014/35/EU 2014/35/UE / Directiva sobre baja tens Директива по низкому напряжению 2 2014/30/EU EMV-Richtlinie 2014/30/EU / EMC Dire 2014/30/UE / Direttiva EMC 2014/30/U 2011/65/EU, (EU) 2015/863 RoHS-Richtlinien 2011/65/EU und (EU Directives RoHS 2011/65/UE et (UE) 2 	mity with the following EU Directives: hes aux directives UE suivantes: h las siguientes directivas de la UE: seguenti direttive UE: соответствуют следующим EU руководствам: / Low voltage directive 2014/35/EU / Directive basse tension ión 2014/35/UE / Direttiva Bassa tensione 2014/35/UE / 2014/35/EU ective 2014/30/EU / Directive CEM 2014/30/UE / Directiva CEM ective 2014/30/EU / Directive CEM 2014/30/UE / Directiva CEM E / Директива ЭМС 2014/30/EU) 2015/863 / RoHS Directives 2011/65/EU and (EU) 2015/863 / Directivas RoHS 2011/65/UE y (UE) 2015/863 / Di-





P. Wimmer Vice President Vice President Vice président Vicepresidente исергезидент J. Bollaender Leiter F & E Director R & D Chef de service R&D Responsable I & D Direttore R & D Глава департамента R&D

2/2

 BINDER GmbH
 Postfach 102
 D-78502 Tuttlingen
 Anschrift: BINDER
 GmbH
 Im Mittleren
 Ösch 5
 D-78532 Tuttlingen

 Kontakt:
 Telefon: +49 (0) 74 62 / 20 05 - 0
 |
 Telefax: +49 (0) 74 62 / 20 05 - 100
 |
 info@binder-world.com
 |
 www.binder-world.com

 Geschäftsführung:
 Dipl-Ing. Peter M. Binder
 |
 Antsgericht Stuttgart, HRB 727150
 |
 Sitz der Gesellschaft: Tuttlingen

 Bankverbindung:
 Kreissparkasse Tuttlingen
 IBAN-Code: DE56633 70075 0213870900
 |
 SWIFT-Code: SOLA DE S1TUT

 Deutsche Bank
 Tuttlingen
 IBAN-Code: DE56633 70075 0213870900
 |
 SWIFT-Code: DEUT DE S563

 Altgeräte-Entsorgung gemäß WEEE-Reg.-Nr. DE 37004983
 |
 SWIFT-Code: DEUT DE S5653



17.6 Certificate for the GS mark of conformity of the "VDE Prüf- und Zertifizierungsinstitut" (Testing and Certification Institute of the Association for Electrical, Electronic and Information Technologies)



1792300-2945-0010 / 241177 / TL4 / ZIE 2019-01-30 2016-09-2 Dieses Blatt gilt nur in Verbindung mit Blatt 1 des Zeichengenehmigungsausweises Nr. 40045043. This supplement is only valid in conjunction with page 1 of the Certificate No. 40045043. Wärmeschrank, Labor Heating cabinet, laboratory Typ(en) / Type(s) A) BF056 B) BF115 C) BF260 D) BD056 E) BD115 F) BD280 G) ED056 K) FD115 D) FD056 K) FD115 D) FD056 K) FD115 D) FD056 K) FD115 D) FD280 M) FED056 N) FED15 O) FD720 R) BF720 R) BT720 S) FD720 Bererkung BF - Brutschrank mit Lüfter zur Lufturwälzung / Incubator with fan for air circulation Remarkur BF - Brutschrank mit Lüfter zur Lufturwälzung / Incubator with fan for air circulation FD - mit Lüfter zur Lufturwälzung / Incubator with fan for air circulation FD - mit Lüfter zur Lufturwälzung / Incubator with fan for air circulation FD - mit Lüfter zur Lufturwälzung / Incubator with	VD Zei	E Prüf- und chengenehr	Zertifizierungsinstitut nigung	Ausweis <i>Certifica</i> 40045	ate No. Page
1792300-2945-0010 / 241177 / TL4 / ZIE 2019-01-30 2016-09-2 Dieses Blatt gilt nur in Verbindung mit Blatt 1 des Zeichengenehmigungsausweises Nr. 40045043. This supplement is only valid in conjunction with page 1 of the Certificate No. 40045043. Wärmeschrank, Labor Heating cabinet, laboratory Typ(en) / Type(s) A) BF056 B) BF115 C) BF260 D) BD056 E) BD115 F) BD280 G) ED056 K) FD115 D) FD056 K) FD115 D) FD280 R) FED15 O) FED260 R) BF720 R) BF720 S) FD720 Bererkung BF - Brutschrank mit Lüfter zur Lufturmwälzung / Incubator with fan for air circulation Remarkung BF - Brutschrank mit Lüfter zur Lufturmwälzung / Incubator with fan for air circulation FD - mit Lüfter zur Lufturmwälzung / unbator with fan for air circulation PD - BT20 Bererkung BF - Brutschrank mit Lüfter zur Lufturmwälzung / Incubator with fan for air circulation FD - mit Lüfter zur Lufturmwälzung / Incubator with fan for air circulation FD - mit Lüfter zur Lufturmwälzung / Incubator with FD - mit Lüfter zur Luftu				holder	
This supplement is only valid in conjunction with page 1 of the Certificate No. 40045043. Wärmeschrank, Labor Heating cabinet, laboratory Typ(en) / Type(s) A) BF056 B) BF115 C) BF280 D) BD056 E) BD115 F) BD260 G) ED056 H) ED115 I) ED260 J) FD056 K) FD115 L) FD280 M) FED115 O) FED280 M] FED210 BF720 C) BD720 R] ED720 R] ED720 R] ED720 BF-Srutschrank mit Lüfter zur Luftunwälzung / Incubator with fan for air circulation FO- mit Lüfter zur Luftunwälzung / Incubator with FD- mit Lüfter zur Luftunwälzung / Incubator			7 / TL4 / ZIE		Datum <i>/ Date</i> 2016-09-28
Heating cabinet, laboratory Typ(en) / Type(s) A) BF056 B) BF115 C) BF260 D) BD056 E) BD115 F) BD260 G) ED056 H) ED115 I) ED260 J) FD260 M) FED260 N) FD260 N) FED260 P) BF720 Q) BD720 R) ED720 S) FD720 Bemerkung Remark BF - Brutschrank mit Lüfter zur Luftumwälzung / Incubator with fan for air circulation B) - Brutschrank / Incubator E) - Universal Warmeschrank, natürliche Luftumwälzung / Universal Warmeschrank, natürliche Luftumwälzung / Universal Warmeschrank, natürliche Luftumwälzung / Universal Warmeschrank, natürliche zur Luftumwälzung / Universal Warmeschrank mit Düfter zur Luftumwälzung / Incubator with fan for air circulation F) - mit Lüfter zur Luftumwälzung / Incubator with fan for air circulation F) - mit Lüfter zur Luftumwälzung / Incubator with					
A) BF056 B) BF115 C) BF260 D) BD056 E) BD115 F) BD260 G) ED056 H) ED115 I) ED260 J) FD056 K) FD115 L) FD260 M) FED156 N) FED156 O) FED260 P) BF720 Q) BD720 R) ED720 S) FD720 S) FD720 S) FD720 S) FD720 S) FD720 S) FD720 Bernerkung Remark BF - Brutschrank mit Lüfter zur Luftumwälzung / Incubator with fan for air circulation BD - Brutschrank / Incubator ED - Universal heating cabinet, natural air circulation FD - mit Lüfter zur Luftumwälzung / Unubator with fan for air circulation BD - Brutschrank / Incubator ED - Universal heating cabinet, natural air circulation FD - mit Lüfter zur Luftumwälzung / Unubator with fan for air circulation FD - mit Lüfter zur Luftumwälzung / Unubator with FD - mit Lüfter zur Luftumwälzung / Unubator with FD - mit Lüfter zur Luftumwälzung / Unubator with FD - mit Lüfter zur Luftumwälzung / Incubator with					
 B BF115 BF260 BD056 BD115 BD260 G) ED056 H) ED115 I) ED260 J) FD056 K) FD115 L) FD260 M) FED056 N) FED115 O) FED260 P) BF720 Q) BD720 R) ED720 S) FD720 T) FD720 BEmerkurg Remark BF - Brutschrank mit Lüfter zur Luftumwälzung / Incubator with fan for air circulation FD - Brutschrank / Incubator BD - Brutschrank mit Lüfter zur Luftumwälzung / Incubator with fan for air circulation FD - Brutschrank mit Lüfter zur Luftumwälzung / Universal heating cabinet, natural air circulation FD - mit Lufter zur Luftumwälzung / Universal heating cabinet, natural air circulation FD - mit Lufter zur Luftumwälzung / Universal heating cabinet, natural air circulation FD - mit Lufter zur Luftumwälzung / Universal heating cabinet, natural air circulation FD - mit Lufter zur Luftumwälzung / Universal heating cabinet, natural air circulation FD - mit Lufter zur Luftumwälzung / Universal heating cabinet, natural air circulation FD - mit Lufter zur Luftumwälzung / Universal heating cabinet, natural air circulation FD - mit Lufter zur Luftumwälzung / Universal heating cabinet, natural air circulation FD - mit Lufter zur Luftumwälzung / Universal heating cabinet, natural air circulation FD - mit Lufter zur Luftumwälzung / Universal heating cabinet, natural air circulation FD - mit Lufter zur Luftumwälzung / Universal heating cabinet, natural air circulation FD - mit Lufter zur Luftumwälzung / Universal heating cabinet, natural air circulation FD - mit Lufter zur Luftumwälzung / Universal heating cabinet, natural air circulation FD - mit Lufter zur Luftumwälzung / Universal heating cabinet, natural air circulation FD - mit Lufter zur Luftumwälzung / Universal heating cabinet, natural air circulation FD - mit Lufter zur Luftumwälzung / Universal heating cabinet, natural air circulation FD - mit Lufter zur Luftumwälzung / Universal heating	Тур(ег	n) / <i>Type(s)</i>			
Bemerkung BF - Brutschrank mit Lüfter zur Luftumwälzung / Incubator with Remark BD - Brutschrank / Incubator ED - Universalwärmeschrank, natürliche Luftumwälzung / Universal heating cabinet, natürliche icrculation FD - mit Lüfter zur Luftumwälzung / with fan for air circulation FED - mit Lüfter zur Luftumwälzung und Drehzahlregelung / with fan for air circulation BF - Brutschrank mit Lüfter zur Luftumwälzung / Incubator with	B)C)D)E)F)G)H)J)KL)MN)OP)Q)R)	BF115 BF260 BD056 BD115 BD260 ED056 ED115 ED260 FD056 FD115 FD260 FED156 FED115 FED260 BF720 BD720 ED720			
Fortsetzung siehe Blatt 3 / Continued on page 3	Fortse	rkung rk etzung siehe Blatt 3 /	fan for air circulation BD - Brutschrank / Incubator ED - Universalwärmeschrank, natürli Universal heating cabinet, natural air FD - mit Lüfter zur Luftumwälzung / with fan for air circulation and speed BF - Brutschrank mit Lüfter zur Luftur fan for air circulation Die zwei / drei folgenden Ziffern beze Innenraumvolumen / The two / three following digits are sig volume	che Luftumwälzung / circulation vith fan for air circulation ind Drehzahlregelung / regulation mwälzung / Incubator with ichnen das	

Zeichengenehm	Zertifizierungsinstitut ligung		Certificate No. 40045043	Page 3
Name und Sitz des Genehmigungs-In Binder GmbH, Im Mittleren Ös	habers / Name and registered seat of the Certificate ch 5, 78532 Tuttlingen	e holder		
Aktenzeichen <i>/ File ref.</i> 1792300-2945-0010 / 241177	/ TL4 / ZIE	letzte Änderung / <i>u</i> 2019-01-30	• CARRON (1997) (1997) (1997) (1997)	um / <i>Date</i> 16-09-28
	mit Blatt 1 des Zeichengenehmigungsauswe onjunction with page 1 of the Certificate No. 40			
PAK AfPS GS 2014:01PAK	Das Produkt entspricht den Anforde PAK-Dokument AfPS GS 2014:01P.			
PAH AfPS GS 2014:01PAH	The product is in accordance with th PAH-document AfPS GS 2014:01P/	ne requirements of AH.		
Weitere Angaben Further information	Siehe Anlage Nr. 1 von 2019-01-30 See Appendix No. 1 dated 2019-01-			
Section TL4				
	VDE Prùf- und Zertifizierungsinstitut GmbH * Testing and Certif			

VDE Prüf- un Zeichengene	nd Zertifizierungsinstitut Phmigung	Cer	weis-Nr. / <i>tificate No</i> .)45043	Beiblatt / <i>Suppleme</i> i
	ungs-Inhabers / Name and registered seat of the Certificate en Ösch 5, 78532 Tuttlingen	holder		
Aktenzeichen / <i>File ref.</i> 1792300-2945-0010 / 24	11177 / TL4 / ZIE	letzte Änderung / <i>updat</i> 2019-01-30		m / <i>Date</i> 6-09-28
	teil des Zeichengenehmigungsausweises Nr. 400450 <i>he Certificate No. 40045043.</i>	43.		
Wärmeschrank, La Heating cabinet, la				
Fertigungsstätte(n) Place(s) of manufactur	e			
Referenz/ <i>Reference</i> 30007949	Binder GmbH Gänsäcker 16 78532 Tuttlingen			
VDE Prüf- und Zertifizier VDE Testing and Certific Fachgebiet TL4 Section TL4				
				<u>^</u>
		cation Institute		

VDE Prüf- und Zertifizierungsinstitut Zeichengenehmigung	Ausweis-Nr. / Infoblatt / Certificate No. Info shee 40045043
Name und Sitz des Genehmigungs-Inhabers / <i>Name and registered seat of the Certificate</i> Binder GmbH, Im Mittleren Ösch 5, 78532 Tuttlingen	holder
Aktenzeichen / <i>File ref.</i> 1792300-2945-0010 / 241177 / TL4 / ZIE	letzte Änderung / updatedDatum / Date2019-01-302016-09-28
Dieses Blatt gilt nur in Verbindung mit Blatt 1 des Zeichengenehmigungsauswei This supplement is only valid in conjunction with page 1 of the Certificate No. 40	
Genehmigung zum Benutzen des auf Seite 1 abgebildeten markenr VDE:	echtlich geschützten Zeichens des
Grundlage für die Benutzung sind die Allgemeinen Geschäftsber Zertifizierungsinstitut GmbH (www.vde.com/AGB-Institut). Das Recht z bezeichnete Firma mit den genannten Fertigungsstätten und die zugeordneten Bezeichnungen. Die Fertigungsstätte muss so eing Herstellung der geprüften und zertifizierten Ausführung gewährleistet ist Die Genehmigung ist so lange gültig wie die VDE-Bestimmungen gelten haben, sofern sie nicht auf Grund anderer Bedingungen aus der VDE F zurückgezogen werden muss. Der Gültigkeitszeitraum einer VDE-GS-Zeichengenehmigung kann gesetzlichen und / oder normativen Änderungen kann die VDE-GS- einem früheren als dem angegebenen Datum verlieren. Produkte, die das Biozid Dimethylfumarat (DMF) enthalten, dürfen 2009/251/EG nicht mehr in den Verkehr gebracht oder auf dem Markt be Der VDE-Zeichengenehmigungsausweis wird ausschließlich auf der ers	zur Benutzung erstreckt sich nur auf die oben aufgeführten Produkte mit der erichtet sein, dass eine gleichmäßig , die der Zertifizierung zugrunde geleger Prüf- und Zertifizierungsordnung (PM102 auf Antrag verlängert werden. Be Zeichengenehmigung ihre Gültigkeit zu gemäß der Kommissionsentscheidung ereitgestellt werden. ten Seite unterzeichnet.
Approval to use the legally protected Mark of the VDE as shown on Basis for the use are the general terms and conditions of the V (www.vde.com/terms-institute). The right to use the mark is granted of named places of manufacture and the listed products with the related ty shall be equipped in a way that a constant manufacturing of the certified The approval is valid as long as the VDE specifications are in force unless it is withdrawn according to the VDE Testing and Certification Pro The validity period of a VDE-GS-Mark Approval may be prolonged on ro or normative requirements, the validity period of a VDE-GS-Mark Appro Products containing the biocide dimethylfumarate (DMF) may not be market according to the Commission Decision 2009/251/EC. The approval is solely signed on the first page.	VDE Testing and Certification Institute only to the mentioned company with the pereferences. The place of manufacture I construction is assured. on which the certification is based on pocedure (PM102E). equest. In case of changes in legal and val may be shortened.

17.7 Certificate of Compliance for the UL Certification Mark from Underwriters Laboratories

Certificate Number	2019-2-26-E200795
Report Reference	E200795-D1002-1/A0/C3-UL
Issue Date	2019-2-26
Issued to:	Binder GmbH
Applicant Company:	Im Mittleren Oesch 5 Tuttlingen, 78532 DE
Listed Company:	Same as Applicant
This is to certify that	Heating cabinet, laboratory
representative samples of	BF056-UL, BF115-UL, BF260-UL, BF720-UL, BD056-UL,
	BD115-UL, BD260-UL, BD720-UL, ED056-UL, ED115-UL,
	ED260-UL, FD056-UL, FD115-UL,FD260-UL, FED056-UL, FED115-UL, FED260-UL, FED720-UL
	Have been investigated by UL in accordance with the Standard(s) indicated on this Certificate.
Standard(s) for Safety:	UL 61010-1, 3rd Edition, May 11, 2012, Revised July 15 201
)(4)(4)(4)(4)(4)(4)(4)(4)(4)(4)(4)(4)(4)	CAN/CSA-C22.2 No. 61010-1-12, 3rd Edition, Revision date July 2015
Additional Standards:	IEC 61010 2-010: 2014 (Third Edition)
Additional Information:	See the UL Online Certifications Directory at
Only those products bearing the UL Certific Certification and Follow-Up Service.	www.ul.com/database for additional information. ation Mark should be considered as being covered by UL's
Look for the UL Certification Mark on the pr	roduct
eYuleYuleYuleYu	
This is to certify that representative samples according to the current UL requirements.	s of the product as specified on this certificate were tested
3 mars of the	
Bruce Mahrenholz, Assistant Chief Engineer, Global Inspection and Field Service	
Joseph Hosey, General Manager, Director of Sales – Canada, UNDERWRITERS L	ABORATORIES OF CANADA INC.

18. Contamination clearance certificate

18.1 For chambers located outside the USA and Canada

Declaration with regard to safety and health

Erklärung zur Sicherheit und gesundheitlichen Unbedenklichkeit

The German Ordinance on Hazardous Substances (GefStofV), and the regulations regarding safety at the workplace, require that this form be filled out for all products that are returned to us, so that the safety and health of our employees can be warranted.

Die Sicherheit und Gesundheit unserer Mitarbeiter, die Gefahrstoffverordnung GefStofV und die Vorschriften zur Sicherheit am Arbeitsplatz machen es erforderlich, dass dieses Formblatt für alle Produkte, die an uns zurückgeschickt wird.



In the absence of a completely filled out form, a repair is not possible. Ohne Vorliegen des vollständig ausgefüllten Formblattes ist eine Reparatur nicht möglich.

 A completely filled out form should be transmitted by Fax (+49 (0) 7462 2005 93555) or by letter in advance to us, so that this information is available before the equipment/component part arrives. A second copy of this form should accompany the equipment/component part. Eventually the carrier should be informed.

Eine vollständig ausgefüllte Kopie dieses Formblattes soll per Telefax (Nr. +49 (0) 7462 2005 93555) oder Brief vorab an uns gesandt werden, so dass die Information vorliegt, bevor das Gerät/Bauteil eintrifft. Eine weitere Kopie soll dem Gerät/Bauteil beigefügt sein. Ggf. ist auch die Spedition zu informieren.

Incomplete information or non-conformity with this procedure will inevitably lead to substantial delays in
processing. We hope you will have understanding for this measure, which lies outside of our area of
influence, and that you will help us to speed up this procedure.

Unvollständige Angaben oder Nichteinhalten dieses Ablaufs führen zwangsläufig zu beträchtlichen Verzögerungen in der Abwicklung. Bitte haben Sie Verständnis für Maßnahmen, die außerhalb unserer Einflussmöglichkeiten liegen und helfen Sie mit, den Ablauf beschleunigen.

• Please fill out this form completely.

Bitte unbedingt vollständig ausfüllen!

1.	Chamber/ component part / type: / Gerät / Bauteil / Typ:
2.	Serial No./ Serien-Nr.:
3.	Details about utilized substances / biological substances / Einzelheiten über die eingesetzten Substanzen/biologische Materialien:
3.1	Designations / Bezeichnungen:
a)	
b)	
c)	
3.2	Safety measures required for handling these substances / Vorsichtsmaßnahmen beim Umgang mit diesen Stoffen:
a)	
b)	
c)	

3.3	Measures to be taken in case of skin contact or release into the atmosphere / Maßnahmen bei Personenkontakt oder Freisetzung:
a)	-
b)	
c)	
d)	
3.4	Other important information that must be taken into account / Weitere zu beachtende und wichtige Informationen:
a)	
b)	
c)	
4.	Declaration on the risk of these substances (please checkmark the applicable items) / Erklärung zur Gefährlichkeit der Stoffe (bitte Zutreffendes ankreuzen) :
□ 4.1	For non toxic, non radioactive, biologically harmless materials / für nicht giftige, nicht radio- aktive, biologisch ungefährliche Stoffe:
	rewith guarantee that the above-mentioned chamber / component part… / Wir versichern, dass rät/Bauteil
	not been exposed to or contains any toxic or otherwise hazardous substances / weder giftige noch stige gefährliche Stoffe enthält oder solche anhaften.
	t eventually generated reaction products are non-toxic and also do not represent a hazard / auch entstandene Reaktionsprodukte weder giftig sind noch sonst eine Gefährdung darstellen.
	ntual residues of hazardous substances have been removed / evtl. Rückstände von Gefahrstoffen ent- t wurden.
4.2	For toxic, radioactive, biologically harmful or hazardous substances, or any other ha- zardous materials / für giftige, radioaktive, biologisch bedenkliche bzw. gefährliche Stoffe oder an- derweitig gefährliche Stoffe.
We he	rewith guarantee that … / Wir versichern, dass …
pon plet	e hazardous substances, which have come into contact with the above-mentioned equipment/com- ent part, have been completely listed under item 3.1 and that all information in this regard is com- e / die gefährlichen Stoffe, die mit dem o.g. Gerät/Bauteil in Kontakt kamen, in 3.1 aufgelistet sind und alle aben vollständig sind.
	t the chamber /component part has not been in contact with radioactivity / das Gerät/Bauteil nicht mit ioaktivität in Berührung kam
5. I	Kind of transport / transporter / Transportweg/Spediteur:
Transp	ort by (means and name of transport company, etc.) Versendung durch (Name Spediteur o.ä.)
Date of	f dispatch to BINDER GmbH / Tag der Absendung an BINDER GmbH:

We herewith declare that the following measures have been taken / Wir erklären, dass folgende Maßnahmen getroffen wurden:
Hazardous substances were removed from the chamber / component part, so that no hazard exists for corresponding persons in the handling or repair of these items / das Gerät/Bauteil wurde von Gefahrstoffen befreit, so dass bei Handhabung/Reparaturen für die betreffenden Person keinerlei Gefährdung besteht
The chamber was securely packaged and properly identified / das Gerät wurde sicher verpackt und vollständig gekennzeichnet.
 Information about the hazardousness of the shipment (if required) has been provided to the transporter / der Spediteur wurde (falls vorgeschrieben) über die Gefährlichkeit der Sendung informiert.
We herewith commit ourselves and guarantee that we will indemnify BINDER GmbH for all damages that are a consequence of incomplete or incorrect information provided by us, and that we will exempt BINDER GmbH from eventual damage claims by third parties./ Wir versichern, dass wir gegenüber BINDER für jeden Schaden, der durch unvollständige und unrichtige Angaben entsteht, haften und BINDER gegen eventuell entstehende Schadenansprüche Dritter freistellen.
We are aware that, in accordance with Article 823 of the German Civil Code (BGB), we are directly liable with regard to third parties, in this instance especially the employees of BINDER GmbH, who have been entrusted with the handling / repair of the chamber / component. / Es ist uns bekannt, dass wir gegenüber Dritten – hier insbesondere mit der Handhabung/Reparatur des Geräts/des Bauteils betraute Mitarbeiter der Firma BINDER - gemäß §823 BGB direkt haften
Name:
Position:
Date / Datum:
Signature / Unterschrift:
Company stamp / Firmenstempel:

A 39	Equipment that is returned to the factory for repair must be accompanied by a completely filled out contamination clearance certificate. For service and maintenance works on site, such a contamination clearance certificate must be submitted to the service technician before the start of the works. No repair or maintenance of the equipment is possible, without a properly filled out contamination clearance certificate.
-------------	-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

18.2 For chambers located in the USA and Canada

Product Return Authorization Request

Please complete this form and the Customer Decontamination Declaration (next 2 pages) and attach the required pictures. E-mail to: IDL_SalesOrderProcessing_USA@binder-world.com

After we have received and reviewed the complete information we will decide on the issue of a RMA number. Please be aware that size specifications, voltage specifications as well as performance specifications are available on the internet at <u>www.binder-world.us</u> at any time.

	Please fill:		
Reason for return request	O Duplicate order		
	O Duplicate	e shipment	
	O Demo		Page one completed by sales
	O Power P	lug / Voltage	115V / 230 V / 208 V / 240V
	O Size doe	s not fit space	
	O Transpo	rt Damage	Shock watch tripped? (pictures)
	O Other (sp	pecify below)	
Is there a replacement PO?	O Yes	O No	
If yes -> PO #			
If yes -> Date PO placed			
Purchase order number			
BINDER model number			
BINDER serial number			
Date chamber was received			
Was the chamber unboxed?	O Yes	O No	
Was the chamber plugged in?	O Yes	O No	
Was the chamber in opera- tion?	O Yes	O No	
Pictures of chamber at-	O Yes	O No	Pictures have to be attached!
tached? Pictures of Packaging at- tached?	O Yes	O No	

Take notice of shipping laws and regulations.

	Customer Contact Information	Distributor Contact Information
Name		
Company		
Address		
Phone		
E-mail		

Customer (End User) Decontamination Declaration

Health and Hazard Safety declaration

To protect the health of our employees and the safety at the workplace, we require that this form is completed by the user for all products and parts that are returned to us. (Distributors or Service Organizations cannot sign this form)

NO RMA number will be issued without a completed form. Products or parts returned to our NY warehouse without a RMA number will be refused at the dock.

A second copy of the completed form must be attached to the outside of the shipping box.

1.	Chamber/ component part / type:
2.	Serial No.
3.	List any exposure to hazardous liquids, gasses or substances and radioactive material
3.1	List with MSDS sheets attached where available or needed
(if ther	e is not enough space available below, please attach a page):
a)	
b)	
c)	
3.2	Safety measures required for handling the list under 3.1
a)	
b)	
c)	
3.3	Measures to be taken in case of skin contact or release into the atmosphere:
a)	
b)	
c)	
d)	
3.4	Other important information that must be considered:
a)	
b)	
c)	

4.	Declaration of Decontamination		
	kic, radioactive, biologically and chemically harmful or hazardous substances, or any other		
	lous materials.		
We he	reby guarantee that		
4.1 Any hazardous substances, which have come into contact with the above-mentioned equipment / component part, have been completely listed under item 3.1 and that all information in this regard is complete.			
	e chamber /component part has not been in contact with radioactivity		
4.3 Ar	ny Hazardous substances were removed from the chamber / component part, so that no hazard kists for a persons in the shipping, handling or repair of these returned chamber		
or	ne chamber was securely packaged in the original undamaged packaging and properly identified in the outside of the packaging material with the chamber designation, the RMA number and a opy of this declaration.		
4.5 Sł	nipping laws and regulations have not been violated.		
l hereb	y commit and guarantee that we will indemnify BINDER Inc. for all damages that are a con-		
	nce of incomplete or incorrect information provided by us, and that we will indemnify and		
	armless BINDER Inc. from eventual damage claims by third parties.		
Name:			
Positio	n:		
Compa	iny:		
Addres	s:		
Phone	#•		
FIIONE	<i>π</i>		
Email:			
Date:			
Signatu	ıre:		



Equipment returned to the NY warehouse for repair must be accompanied by a completed customer decontamination declaration. For service and maintenance works on site, such a customer decontamination declaration must be submitted to the service technician before the start of work. No repair or maintenance of the equipment is possible without a completed form.