DOC022.98.80205

**LANGE** 

**2100AN**

11/2014, Edition 6

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## Specifications

Specifications are subject to change without notice.

Specification	Details
Measurement method	Nephelometric
Regulatory	Meets EPA Method 180.1 ASTM D7315 - Standard Test Method for Determination of Turbidity Above 1 Turbidity Unit (TU) in Static Mode ASTM D6655 - Standard Test Method for Determination of Turbidity Below 5 NTU in Static Mode
Light source	Tungsten filament lamp
Measurement modes	NTU, NEP (Nephelo), EBC, Abs (absorbance), %T (% transmittance), CU (color units) and two user-defined units

Specification	Details
Range	NTU (Ratio on, manual range): 0–0.999, 0–9.99, 0–99.9, 0–10,000 NTU (Ratio on, auto range): 0–10,000 auto decimal NTU (Ratio off): 0–40 Nephelo (Ratio on, manual range): 0–9.99, 0–99.9, 0–67,000 Nephelo (Ratio on, auto range): 0–67,000 auto decimal Nephelo (Ratio off): 0–268 EBC (Ratio on, manual range): 0–0.999, 0–9.99, 0–99.9, 0–2450 EBC (Ratio on, auto range): 0–2450 auto decimal EBC (Ratio off): 0–9.8 Absorbance (manual range): 0–0.999, 0–2.00 Absorbance (auto range): 0–2.00 Transmittance (%): 1.0–100 Color units (at 455 nm): 0–500
Accuracy <sup>1, 2, 3, 4</sup>	Ratio on: $\pm 2\%$ of reading plus 0.01 NTU from 0–1000 NTU, $\pm 5\%$ of reading from 1000–4000 NTU, $\pm 10\%$ of reading from 4000–10,000 NTU Ratio off: $\pm 2\%$ of reading plus 0.01 NTU from 0–40 NTU Absorbance: $\pm 0.005$ Abs from 0–1 Abs at 455 nm Transmittance: 0.12% T from 10–100% T at 455 nm Color: $\pm 2$ CU from 0–30 (calibrated at 15 CU), $\pm 5$ CU from 0–500 CU (calibrated at 500 CU)

- <sup>1</sup> Turbidity specifications identified using USEPA filter assembly, recently prepared formazin standard and matched 25-mm sample cells.
- <sup>2</sup> Use of a flow cell is necessary to get color measurement specifications.
- <sup>3</sup> Intermittent electromagnetic radiation of 3 volts/meter or greater may cause slight accuracy shifts.
- <sup>4</sup> Reference conditions: 23  $\pm$  2 °C, 50%  $\pm$  10% RH noncondensing, 115/230 VAC, 50/60 Hz

Specification	Details
Resolution	Turbidity: 0.001 NTU/EBC Nephelo: 0.01 Absorbance: 0.001 Abs Transmittance: 0.1% T Color: 1 CU
Repeatability	±1% of reading or 0.01 NTU, whichever is greater (under reference conditions)
Response time	Signal averaging off: 6.8 seconds Signal averaging on: 14 seconds (when 10 measurements are used to calculate the average)
Stabilization time	Ratio on: 30 minutes after start-up Ratio off: 60 minutes after start-up
Reading modes	Manual or auto range, signal averaging on and adjustable or off, Ratio on or off
Power requirement	115–230 VAC, 50/60 Hz (automatic power selection) 28 W maximum
Pollution degree/installation category	2; II
Protection Class	1
Operating conditions	Temperature: 0 to 40 °C (32 to 104 °F) Relative humidity: 0–90% at 25 °C, 0–75% at 40 °C, noncondensing Altitude: 2000 m (6560 ft) maximum Indoor use only
Storage conditions	–40 to 60 °C (–40 to 140 °F), instrument only
Printer	Built-in (thermal, 58-mm, up to 28 column)
Interface	RS232C serial interface by way of DB9 subminiature D-shell connector for data output to computer or printer, and data input (command). No handshaking.

Specification	Details
Air purge	Dry nitrogen or instrument grade air (ANSI MC 11.1, 1975) 0.1 scfm at 69 kPa (10 psig); 138 kPa (20 psig) maximum Hose barb connection for 1/8-inch tubing
Sample cells	Round cells 95 x 25 mm (3.74 x 1 in.) borosilicate glass with rubber-lined screw caps <b>Note:</b> <i>Smaller sample cells (less than 25 mm) can be used when a cell adapter is used.</i>
Sample requirements	25 mm sample cell: 20 mL minimum 0 to 95 °C (32 to 203 °F)
Enclosure	High-impact polycarbonate plastic
Dimensions	30.5 x 40 x 15.6 cm (12.0 x 15.7 x 6.1 in.)
Weight	3.8 kg (8.5 lb)
Certification	CE, cETLus

## General information

In no event will the manufacturer be liable for direct, indirect, special, incidental or consequential damages resulting from any defect or omission in this manual. The manufacturer reserves the right to make changes in this manual and the products it describes at any time, without notice or obligation. Revised editions are found on the manufacturer's website.

## Additional information

Additional information is available on the manufacturer's website.

## Safety information

### NOTICE

The manufacturer is not responsible for any damages due to misapplication or misuse of this product including, without limitation, direct, incidental and consequential damages, and disclaims such damages to the full extent permitted under applicable law. The user is solely responsible to identify critical application risks and install appropriate mechanisms to protect processes during a possible equipment malfunction.

Please read this entire manual before unpacking, setting up or operating this equipment. Pay attention to all danger and caution statements. Failure to do so could result in serious injury to the operator or damage to the equipment.

Make sure that the protection provided by this equipment is not impaired. Do not use or install this equipment in any manner other than that specified in this manual.

## Use of hazard information

### ▲ DANGER

Indicates a potentially or imminently hazardous situation which, if not avoided, will result in death or serious injury.

### ▲ WARNING

Indicates a potentially or imminently hazardous situation which, if not avoided, could result in death or serious injury.

### ▲ CAUTION

Indicates a potentially hazardous situation that may result in minor or moderate injury.

### NOTICE

Indicates a situation which, if not avoided, may cause damage to the instrument. Information that requires special emphasis.

## Precautionary labels

Read all labels and tags attached to the instrument. Personal injury or damage to the instrument could occur if not observed. A symbol, if noted

on the instrument, will be included with a danger or caution statement in the manual.

	This symbol, if noted on the instrument, references the instruction manual for operation and/or safety information.
	Electrical equipment marked with this symbol may not be disposed of in European public disposal systems after 12 August of 2005. In conformity with European local and national regulations (EU Directive 2002/96/EC), European electrical equipment users must now return old or end-of-life equipment to the Producer for disposal at no charge to the user. <b>Note:</b> For return for recycling, please contact the equipment producer or supplier for instructions on how to return end-of-life equipment, producer-supplied electrical accessories, and all auxiliary items for proper disposal.

## Certification

### Canadian Radio Interference-Causing Equipment Regulation, IECs-003, Class A:

Supporting test records reside with the manufacturer.

This Class A digital apparatus meets all requirements of the Canadian Interference-Causing Equipment Regulations.

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### FCC Part 15, Class "A" Limits

Supporting test records reside with the manufacturer. The device complies with Part 15 of the FCC Rules. Operation is subject to the following conditions:

1. The equipment may not cause harmful interference.
2. The equipment must accept any interference received, including interference that may cause undesired operation.

Changes or modifications to this equipment not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment. This equipment has been tested and found to comply with the limits for a Class A digital device, pursuant to Part 15 of the FCC rules. These limits are designed to provide reasonable

protection against harmful interference when the equipment is operated in a commercial environment. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instruction manual, may cause harmful interference to radio communications. Operation of this equipment in a residential area is likely to cause harmful interference, in which case the user will be required to correct the interference at their expense. The following techniques can be used to reduce interference problems:

1. Disconnect the equipment from its power source to verify that it is or is not the source of the interference.
2. If the equipment is connected to the same outlet as the device experiencing interference, connect the equipment to a different outlet.
3. Move the equipment away from the device receiving the interference.
4. Reposition the receiving antenna for the device receiving the interference.
5. Try combinations of the above.

## Product overview


The 2100AN laboratory turbidimeter measures turbidity in NTUs (nephelometric turbidity units), NEPs (nephelos) and EBCs (European Brewing Convention units). NEPs and EBCs are calculated using the conversion factors of 6.7 nephelos per 1.0 NTU and 0.245 EBCs per 1.0 NTU. The 2100AN turbidimeter also measures absorbance, transmittance and color (APHA Pt-Co Method).

In addition, two user-defined measurement units can be specified. The application specific mode of operation uses the nephelometric optical system and the NTU measurement mode.

The turbidimeter has a built-in printer and an RS232 output for connection to a printer, data logger or computer and a recorder output.


The turbidimeter contains a real-time clock with battery. The clock provides a time-date stamp on all data transmitted to the built-in printer or to external devices by way of the RS232 interface (i.e., measurements and calibration records).

**Figure 1 Front overview**

1 Mode display: shows the calibration standard number, setup number or sample number	5 Cover for the sample cell compartment
2 Keypad	6 Printer cover
3 Sample cell holder	7 Eight-digit LED display
4 Light shield	

**Figure 2 Back overview**


1 Power cord connector	5 Air purge fitting
2 Fuse holder	6 Recorder output jack for a chart recorder (0 to 1 V output)
3 Power switch	7 Remote cable jack for flow valve module connection to the automatic flow cell (low pressure)
4 DB9 connector for RS232 cable	8 Lamp access cover

## Product components

Refer to [Figure 3](#) to make sure that all components have been received. If any of these items are missing or damaged, contact the manufacturer or a sales representative immediately.

**Figure 3 Instrument components**

1 2100AN turbidimeter	7 StablCal® Calibration kit
2 USEPA filter assembly	8 Gelex® secondary turbidity standardization kit <sup>5</sup>
3 455 nm filter assembly	9 Dust cover
4 Oiling cloth	10 Printer paper roll (2x) <sup>6</sup>
5 Six 1" sample cells (30 mL) with caps	11 Power cord
6 Silicone oil	

<sup>5</sup> Supplied with 4700100 only.

<sup>6</sup> Do not remove the plastic wrapper from the paper rolls until the paper is installed.

## Installation

### ⚠ DANGER

Multiple hazards. Only qualified personnel must conduct the tasks described in this section of the document.

## Put paper in the printer

### NOTICE


Use only the provided thermal paper. Use of other thermal paper may cause poor print quality and decrease the life of the print-head.

#### Notes:

- Do not rub the thermal paper with a hard object.
  - Do not use chemical paste on thermal paper.
  - A red line on the edge of the thermal paper shows when the paper supply is low.
1. Cut the end of the paper with scissors to make an arrow shape.
  2. Open the printer cover.
  3. Put the point of the thermal paper in the paper entrance slot.
  4. Push the paper through until the point of the paper comes out the exit slot.
  5. Pull the paper out of the exit slot until the full width of the paper is past the exit slot.
  6. Put the paper roll in the printer.
  7. Put the thermal paper through the slot in the printer cover, then close the printer cover.

## User interface

Figure 4 Keypad

1 ENTER key	7 FLOW key
2 EDIT (arrow) keys	8 SIGNAL AVG key
3 SAMPLE key	9 PRINT key
4 RANGE key	10 LINE FEED key
5 UNITS/Exit key	11 SETUP key
6 RATIO key	12 CAL/Zero key

Table 1 Key descriptions


Key	Description
	Enters the value on the display. Starts the measurement of a calibration standard. Clears data from the buffer.
	Changes the numbers and/or letters on the display. Steps through the calibration standards. The right arrow key moves the cursor to the previous or next digit.

**Table 1 Key descriptions (continued)**

Key	Description
SAMPLE	Starts the changing of the sample number shown on the mode display.
RANGE	Selects automatic or manual ranging.
UNITS Exit	Selects the unit of measure. Exits Calibration or Setup mode without saving changes.
RATIO	Turns Ratio on or off.
FLOW	Turns on or off the Flow mode of operation. Used only with the automated flow cell.
SIGNAL AVG	Turns signal averaging on or off.
PRINT	Sends the data that is on the display to a printer or computer. Sends a calibration data report to a printer or computer when in Calibration mode. Sends diagnostic results to a printer or computer if held down when the instrument is turned on. Provides a print of the setup commands when in Setup mode. Turns the print interval feature on or off if the instrument has been configured with a printer interval.
LINE FEED	Moves the printer paper forward one line.
SETUP	Turns on Setup mode and starts the selection of the setup number on the mode display.
CAL Zero	Starts a calibration when in NTU, NEP or EBC mode. Starts analytical zeroing when in CU, %T or Abs mode.

**Figure 5 Indicator lights**

1 Lamp icon light	7 PRINT light
2 "CAL?" light	8 SIGNAL AVG light
3 "Manual" light	9 CAL/Zero light
4 "Auto" light	10 SETUP light
5 RATIO light	11 SAMPLE light
6 FLOW light	

**Table 2 Light descriptions**

Light	Description
	Illuminated when the instrument light source is on. Flashes when there is not sufficient light for measurement.
CAL?	"CAL?" is shown during a calibration if the calibration data is not within the acceptable range. Flashes when the instrument should be calibrated. <i>Note: The CAL? light applies when the USEPA filter and a 25-mm sample cell are used. Ignore the CAL? light if illuminated during calibration when a different filter or a smaller sample cell is used. Push UNITS/Exit to start measurements.</i>
Manual	"MANUAL" is shown above the Range Mode label when the instrument is in manual ranging mode.
Auto	"AUTO" is shown below the Range Mode label when the instrument is in auto ranging mode.
RATIO	Illuminated when Ratio is on.
FLOW	Illuminated when the Flow mode of operation is selected. Flashes when the flow cycle is done.
PRINT	Illuminated when the printer interval feature is selected. Flashes when a print interval has been selected but is not active.
SIGNAL AVG	Illuminated when signal averaging is on.
CAL Zero	Illuminated when Calibration or Zeroing mode is selected.
SETUP	Illuminated when Setup mode is selected.
SAMPLE	Illuminated when Sample mode is selected.

## Startup

### Turn the instrument on

1. Put the instrument on a stable, level surface that is free of vibration. Do not put in direct sunlight.
2. Make sure that there is air circulation around the instrument. Keep the back and area below the instrument free of material that could decrease air flow through the vents.
3. Connect the power cord to the power plug on the back of the instrument.
4. Connect the power cord to a power socket with ground contact.
5. Push the power switch on the back of the instrument to turn the instrument on.

### Turn the keypad sound off (optional)

By default, the instrument makes an audible sound when a key is pushed. To turn the keypad sound off:

1. Push **SETUP**. The SETUP light turns on.
2. Use the arrow keys to select 00.
3. Push **ENTER**.
4. Use the arrow keys to select the sound option:

Option	Description
<b>BEEP ON</b>	An audible sound is made when a key is pushed.
<b>BEEP OFF</b>	No sound is made when a key is pushed.

5. Push **ENTER**.
6. Push **SETUP**.

### Set the date and time

1. Push **SETUP**. The SETUP light turns on.
2. Use the arrow keys to select an option:

Option	Description
<b>05</b>	Sets the hours and minutes (HH-MM).

Option	Description
06	Sets the month and day (MM-DD).
07	Sets the year (YY).

3. Push **ENTER**.
4. Use the arrow keys to change the value.
5. Push **ENTER**.
6. Push **SETUP**.

### Show the current time (optional)

1. Push **SETUP**. The SETUP light turns on.
2. Use the arrow keys to select 08.
3. Push **ENTER**. The current time is shown on the display (HH-MM-SS).
4. Push **SETUP**.

## Standard operation

### Calibrate the turbidimeter with StablCal® Standards

Calibrate the turbidimeter before it is used for the first time using the StablCal® sealed vial standards provided. As an alternative, calibration can be done with recently prepared formazin standards.

Calibrate the turbidimeter at least every 3 months or as specified by the regulating authority when data is used for USEPA reporting.

The instrument is ready for calibration 60 minutes after start-up. Keep the instrument on 24 hours a day if the instrument is used regularly.

**Note:** *Unknown results may occur if standards other than the recommended calibration points are used. The recommended calibration points (< 0.1, 20, 200, 1000, 4000 and 7500 NTU) provide the best calibration accuracy. Use of standards other than StablCal, or user-prepared formazin, may result in less accurate calibrations. The manufacturer cannot guarantee the performance of the instrument if calibrated with co-polymer styrenedivinylbenzene beads or other suspensions.*

### Prepare the StablCal standards

When received and at intervals:

1. Clean the exterior surface of the StablCal vials with laboratory glass cleaning detergent.
2. Rinse the vials with distilled or deionized water.
3. Dry the vials with a lint-free cloth.

**Note:** *Never shake or invert the < 0.1 NTU standard. If the standard has been mixed or shaken, do not move the vial for 15 minutes or more before using.*

**Note:** *Do not remove the caps from the sealed vials.*

Make sure that the StablCal standards are at ambient instrument temperature before use (and no greater than 40 °C (104 °F)).

Mix the standards before use:

1. Open the case lid. Remove the < 0.1 NTU standard from the plastic case.
2. Leave the other standards in the case. Close the case lid.
3. Shake the case vigorously for at least 10 seconds.
4. Let the standards stand with no movement for 3–5 minutes before use.


### Calibration notes

- Make sure that the instrument is in the same ambient conditions as where it is used.
- Make sure that the standards are at the same ambient temperature as the instrument before use.
- Use only the provided silicone oil. This silicone oil has the same refractive index as the vial glass and masks minor glass differences and scratches.
- Store the oiling cloth in a plastic storage bag to keep the cloth clean.
- If power is lost during calibration, the new calibration data is lost and the last calibration data is used. To exit a calibration and not save the new values, push **UNITS/Exit**.
- In Calibration mode, automatic range and signal averaging on are selected. When calibration is completed, all operational modes go back to the last settings.
- All nephelometric (turbidity units of measure) calibrations are done at the same time.


- Ratio-on and Ratio-off calibration data is measured and recorded at the same time.
- The 7500-NTU standard does not have to be measured during calibration if turbidity less than 4000 NTU will be measured. Push **CAL/Zero** after the 4000 NTU standard is measured to complete the calibration procedure.

- Clean the USEPA filter assembly before doing a primary calibration, or at least every 3 months (which is the USEPA-recommended primary calibration interval).


### StabCal calibration procedure


1. Remove the filter assembly. Refer to [Change the filter assembly](#) on page 16.

2. Clean the lens of the USEPA filter assembly. Refer to [Clean the filter assembly](#) on page 17.


3. Hold the tab of the USEPA filter assembly so that the arrows point toward the front of the instrument. Push the filter assembly fully in the housing.

4. Push **CAL/Zero**. The CAL/Zero light turns on, and the mode display shows "00". The NTU value of the dilution water that was used in the previous calibration is shown on the display.


5. Get the < 0.1 NTU vial. Clean the vial with a soft, lint-free cloth to remove water spots and fingerprints. Do not invert the vial.

6. Apply a small bead of silicone oil from the top to the bottom of the vial.


7. Use the oiling cloth to apply the oil equally to the surface of the vial. Remove the excess oil. Make sure that the vial is almost dry.


8. Put the vial in the sample cell holder with the triangle on the vial aligned with the reference mark on the sample cell holder. Close the cover.


9. Push **ENTER**.  
The instrument display counts down, then measures the standard.  
The next expected standard (e.g., 20.00) is shown. The mode display shows "01".

10. Remove the vial from the sample cell holder.

11. Do steps 5–10 for the other StablCal vials (from lowest to highest NTU standard).  
The mode display shows "00" after the last vial is measured.

12. Push **CAL/Zero**.  
The instrument saves the new calibration data and goes back to Measurement mode.

## StablCal standards storage

- Do not move a StablCal standard to a different container for storage. Keep StablCal standards in the plastic case provided with the cover closed.
- Store at 5 to 25 °C (41 to 77 °F).
- For long-term storage (more than one month between use), keep at 5 °C (41 °F).

## Turbidity measurement

### ▲ WARNING

Potential explosion and fire hazard. This instrument is for measuring water based samples. Do not measure solvent or combustible based samples.

For accurate turbidity readings use clean sample cells and remove air bubbles.

## Measurement notes

Proper measurement techniques are important in minimizing the effects of instrument variation, stray light and air bubbles. For accurate and repeatable measurements:

### Instrument

- Make sure that the instrument is on a level, stationary surface that is free of vibration during the measurement.
- The USEPA filter assembly is required for turbidity measurements reported for United States Environmental Protection Agency (USEPA), National Primary Drinking Water Regulations (NPDWR) or National Pollutant Discharge Elimination System (NPDES) permits.
- Turn the instrument on 30 minutes (Ratio on) or 60 minutes (Ratio off) before measurement. Keep the instrument on 24 hours a day if the instrument is used regularly.
- Always close the sample compartment lid during measurement, calibration and storage.
- Remove the sample cell from the instrument and turn off the instrument if the instrument is stored for an extended time period (more than a month).

- Keep the sample compartment lid closed to keep dust and dirt out.

### Sample cells

- Always cap the sample cell to prevent spillage of the sample into the instrument.
- Always use clean sample cells in good condition. Dirty, scratched or damaged cells can result in readings that are not accurate.
- Make sure that cold samples do not “fog” the sample cell.
- Store sample cells filled with distilled or deionized water and cap tightly.
- For the best accuracy, use a single sample cell for every measurement or a flow cell.

### Turbidity measurement procedure

1. Rinse a clean, empty sample cell two times with the solution to be measured and drain to waste. Fill to the line (about 30 mL) with sample and immediately put the cap on the sample cell.


2. Clean the sample cells with a soft, lint-free cloth to remove water spots and fingerprints.


3. Apply a small bead of silicone oil from the top to the bottom of the sample cells.

4. Use the oiling cloth provided to apply the oil equally to the surface of the sample cells. Remove the excess oil. Make sure that the sample cells are almost dry.

5. Gently and slowly invert the sample cell to fully mix the sample. Be careful not to add air bubbles.


6. Put the sample cell in the sample cell holder with the triangle on the sample cell aligned with the reference mark on the sample cell holder. Close the cover.

**Note:** As an alternative, matched sample cells may be used for measurements but do not provide as good of accuracy or precision as a single indexed sample cell or flow cell. When using matched sample cells, align the orientation mark on the sample cell with the reference mark on the sample cell holder.

### Measurement

- Measure samples immediately to prevent temperature changes and settling. Before a measurement is taken, always make sure that the sample is homogeneous throughout.
- Avoid sample dilution when possible.
- Avoid instrument operation in direct sunlight.

7. Read and record the value when stable.

**Note:** To print or send (via RS232) a measurement record, push **PRINT**.

## Measurement techniques

Measurements may be made with different operation mode settings and optional accessories.

Calibrate the instrument whenever the sample cell pathlength is changed.

### Manual or automatic ranging

The manufacturer recommends that ranging be set to automatic for most measurements.

The setting can be changed at any time during sample measurement.

Push **RANGE** repeatedly to step the instrument from automatic ranging to manual ranging and then scroll through the manual range settings.

"MANUAL" is shown above the Range Mode label on the instrument when manual ranging is selected. "AUTO" is shown below the Range Mode label on the instrument when automatic ranging is selected.

Notes:

- When manual ranging is selected, the display flashes all 9s when the sample being measured is greater than the selected range. The display flashes all 0s when the sample measured is less than the selected range.

- When automatic ranging is selected, the display flashes 9s when the sample is greater than the maximum range of the instrument. The display flashes 9s when Ratio is off and the measurement is greater than 40 NTUs (268 nephelos or 9.8 EBCs). Turn Ratio on to increase the range.
- When automatic ranging is selected, the display flashes all 0s when the measurement is less than the range of the instrument or a negative value. Calibrate the instrument. When measuring color, absorbance or transmittance, set the zero reference point again.

### Signal averaging on or off

Signal averaging corrects for reading fluctuations that are caused by random drifting particles in the sample. When signal averaging is on, an average reading is calculated every 3 seconds and shown on the display.

The manufacturer recommends that signal averaging be on for most measurements.

Push **SIGNAL AVG** to turn signal averaging on or off. The SIGNAL AVG light turns on when signal averaging is on.

Push **ENTER** when signal averaging is on to erase data in the signal averaging buffer and provide an immediate update on the display as necessary. This is especially useful when measuring samples with large differences in turbidity.

To change the number of measurements that are used to calculate the average reading (default=10):

1. Push **SETUP**. The SETUP light turns on.
2. Select 09 using the arrow keys.
3. Push **ENTER**.
4. Use the arrow keys to select the number of measurements—1 to 15.  
*Note: If a number greater than 15 is selected, 15 measurements will be used.*
5. Push **ENTER**.
6. Push **SETUP**.

### Ratio on or off

Ratio on provides very good linearity, calibration stability and a wide measurement range. Ratio on helps correct for interference when color is present in the sample that absorbs at the wavelength of incident light.

The manufacturer recommends that Ratio on be used for most measurements. Ratio must be on to measure samples greater than 40 NTUs (268 nephelos or 9.8 EBCs).

Ratio can be turned on for NTU, EBC, NEP and ASC -1- and -2- measurements.

Push **RATIO** to turn Ratio on or off. The Ratio light is on when Ratio is on.

Notes:

- If the sample being measured is greater than 40 NTU (or equivalent) and Ratio is off, the display will show 9s and the RATIO light will flash. Push **RATIO** to turn Ratio on and remove the over-range condition.
- Measurements with Ratio on and measurements with Ratio off are almost the same for turbidity measurements that are less than 40 NTU if interferences caused by color or light absorbing particles are not present.

## Maintenance

### ⚠ DANGER

Multiple hazards. Only qualified personnel must conduct the tasks described in this section of the document.

### Clean the instrument

Keep the instrument clean to get continuous and accurate operation.

### NOTICE

Never use cleaning agents such as turpentine, acetone or similar products to clean the instrument including the keypad.

1. Turn the instrument off and disconnect the power cord.
2. Clean the surface of the instrument with a soft, moist cloth and a weak soap solution.
3. Dry the surface of the instrument with a lint-free cloth.

### Change the filter assembly

### NOTICE

The filter assembly is fragile and must be handled with care to prevent damage.

1. Hold the tab of the filter assembly and pull straight up and out of the instrument.
2. Store the filter assembly in a clean container.
3. Before installation, clean the lens of the filter assembly. Refer to [Clean the filter assembly](#) on page 17.
4. Hold the tab of the filter assembly with the arrows pointing toward the front of the instrument.
5. Push the filter assembly fully into the housing.

## Clean the filter assembly

**Note:** Be careful not to push the lens out of the filter assembly.

1. Clean both sides of the lens of the filter assembly with glass cleaner, lens cleaner or isopropyl alcohol, and a cotton-tipped swab or lens tissue.
2. Inspect the filter glass for scratches or other damage.
3. If a cloudy circle is seen around the edge of the filter, the filter material is delaminating. Replace the filter assembly.

## Replace the lamp

### ⚠ CAUTION

Wear protective eye wear when the lamp is turned on and the lamp cover is removed.

### ⚠ CAUTION





Burn hazard. The lamp must be cool before removal from the instrument.

#### Notes:

- Replace the lamp with the same size, style and electrical rating (4708900).
- Do not touch the lamp as oil from skin will damage the lamp. Clean the lamp with alcohol as necessary.
- Either lamp lead can be put in either terminal block position.
- Turn the instrument on 30 minutes (Ratio on) or 60 minutes (Ratio off) before measurement or calibration.
- Calibrate the instrument after the lamp is replaced.

To remove the lamp, refer to the illustrated steps.

To install the lamp, do the illustrated steps in the opposite direction.

## Replace a fuse


⚠ DANGER	
	<p>Fire hazard. Use the same type and current rating to replace fuses.</p>

Replacement parts:

- Fuse for 115 V operation, time-delay, 250 V, 1.6 A (3030700), or
- Fuse for 230 V operation, time-delay, 250 V, 1.6 A (3030600)

To replace a fuse, refer to the illustrated steps in [Figure 6](#).

**Figure 6** Replace a fuse

## Troubleshooting

Refer to the tables in this section for error codes, diagnostic codes, common problem messages or symptoms, possible causes and corrective actions.

### Error codes

[Table 3](#) lists the error codes shown for different conditions. Error codes identify instrument malfunction or operator error.

The instrument continues operation in an error condition.

Push **ENTER** to clear an error code from the display.

**Note:** Any calibration being calculated when an error occurs, is discarded. The old calibration is kept.

**Table 3** Error codes

Error	Description	Solution
ERR01	The turbidity of the dilution water is greater than 0.5 NTU.	Start the calibration again with lower turbidity dilution water. <b>Note:</b> Ignore ERR01 when the sample cell diameter is less than 25 mm. Push <b>UNITS/Exit</b> to go back to measurement mode.
ERR02	<ul style="list-style-type: none"> <li>• Two calibration standards have the same value.</li> <li>• The difference between two calibration standards is less than 60.0 NTU.</li> <li>• The turbidity of Standard 1 is too low (less than 10 NTU).</li> </ul>	<ol style="list-style-type: none"> <li>1. Inspect the preparation of standards.</li> <li>2. Do the calibration again.</li> </ol> <p><b>Note:</b> Ignore ERR02 when the sample cell diameter is less than 25 mm. Push <b>UNITS/Exit</b> to go back to measurement mode.</p>

**Table 3 Error codes (continued)**

Error	Description	Solution
<b>ERR03</b>	Low light error	<ol style="list-style-type: none"> <li>1. Put the sample in the instrument again.</li> <li>2. Make sure that the lamp icon light is on.</li> <li>3. Make sure that an object is not in the light path.</li> <li>4. Do sample dilution if necessary.</li> </ol> <p><i>Note: If this error occurs when a filter assembly other than the USEPA filter assembly is installed, the filter assembly should not be used for turbidity measurements.</i></p>
<b>ERR04</b>	Memory malfunction	<ol style="list-style-type: none"> <li>1. Turn the instrument off and then back on.</li> <li>2. Contact Technical Support if the error occurs again.</li> </ol>
<b>ERR05</b>	A/D is over the range	<ol style="list-style-type: none"> <li>1. Make sure that the light shield is closed.</li> <li>2. Contact Customer Service if necessary.</li> </ol>
<b>ERR06</b>	A/D is under the range	<ol style="list-style-type: none"> <li>1. Make sure that no object is in the light path.</li> <li>2. Contact Customer Service if necessary.</li> </ol>
<b>ERR07</b>	Light leak	<ol style="list-style-type: none"> <li>1. Make sure that the cover for the sample cell compartment is closed.</li> <li>2. Turn the instrument off and then back on.</li> </ol>
<b>ERR09</b>	Printer time out error or paper in the internal printer can not move	<ol style="list-style-type: none"> <li>1. Gently pull up on the paper in the internal printer to remove the obstruction.</li> <li>2. Make sure that the external printer is connected correctly.</li> <li>3. Make sure that the external printer is selected (online).</li> </ol>

**Table 3 Error codes (continued)**

Error	Description	Solution
<b>ERR10</b>	System voltage out of range	<ol style="list-style-type: none"> <li>1. Turn the instrument off and then back on.</li> <li>2. Contact Customer Service if the error occurs again.</li> </ol>
<b>ERR11</b>	System loop test error	<ol style="list-style-type: none"> <li>1. Turn the instrument off and then back on.</li> <li>2. Contact Customer Service if the error occurs again.</li> </ol>
<b>ERR12</b>	ASC units name error	Enter an application specific calibration (ASC) unit name that is not one of the default units (i.e., NTU or EBC).
<b>ERR14</b>	Invalid time error	The time must be between 00-00 and 23-59.
<b>ERR15</b>	Invalid date error	The date must be between 01-00 and 12-31.

## Diagnostic codes

[Table 4](#) lists the diagnostic codes that are used to get information about instrument operation when instrument operation is in doubt.

To do a diagnostic test:

1. Push **SETUP**.
2. Use the arrow keys to enter a diagnostic code.
3. Push **ENTER** to show the diagnostic value.
4. Push **UNITS/Exit** to go back to Measurement mode.

*Note: To print a diagnostic report, hold down **PRINT**, then turn the instrument on.*

**Table 4 Diagnostic codes**

Code	Display	Description
21	"PRINT TST"	Printer test
22	Test results are shown.	Display test

**Table 4 Diagnostic codes (continued)**

Code	Display	Description
23	Test results are shown.	Keyboard test
24	Test results are shown.	Memory test

sample with the least absorbance. Record the reading as negative absorbance.

## Delete calibration data

To delete any calibration data entered by the user:

1. Turn off the instrument.
2. Push and hold **CAL/Zero**.
3. Turn on the instrument.

The CAL? light flashes. The instrument starts in Calibration mode.

4. Calibrate the instrument before use.

## Flashing 9s

When manual ranging is selected, the display will flash all 9s when the sample being measured is greater than the selected range.

When automatic ranging is selected, the display will flash 9s when the sample is greater than the maximum range of the instrument. The display will also flash 9s if Ratio is off and the measurement is greater than 40 NTUs (268 nephelos or 9.8 EBCs). Turn Ratio on.

## Flashing 0s

When manual ranging is selected, the display will flash all 0s when the sample measured is less than the selected range.

When automatic ranging is selected, the display will flash all 0s when the measurement is less than the range of the instrument or a negative value. Calibrate the instrument.

- When measuring color, absorbance or transmittance, set the zero reference point again.
- When measuring absorbance, make sure that the reading is positive. To measure samples with negative absorbance, set the analytical zero using the sample with the greatest absorbance and read the

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P.O. Box 389, Loveland, CO 80539-0389 U.S.A.
(800) 227-4224 (U.S.A. only)