Oxford A & B Series Instructions

Semimikrowaage SBC 22 SBC 22 Laboratory Balances

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Before reading these instructions, fold out the first page showing a general view of the balance.

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Warranty

Do not miss out on the benefits of our extended warranty. Please complete the warranty registration card and return it to the address listed below (see supplement).

In the unlikely event that your balance should require servicing within the 3-year warranty period, please contact the office listed below for your country, and state the model number, serial number and the problem:

Federal Republic of Germany: SCALTEC INSTRUMENTS GmbH Rudolf-Diesel-Strasse 1 37308 Heiligenstadt, Germany Phone: (0)3606-6701-0 Fax: (0)3606-6701-60 North America: SCALTEC INSTRUMENTS INC P.O. Box 9 Brentwood, New York 11717 Phone: 1-800-7893332 Fax: 516-254-0918

All other countries: SCALTEC INSTRUMENTS GmbH Rudolf-Diesel-Strasse 1 37308 Heiligenstadt, Germany Phone: (0)3606-6701-0 Fax: (0)3606-6701-60

Storage and Shipping Conditions

Allowable storage temperature: +5 °C \ldots +40 °C $_{+41}$ °F \ldots +104 °F

After unpacking the balance, please check it immediately for any visible damage as a result of rough handling during shipment. If this is the case, proceed as directed in the section entitled "Safety Inspection."

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Ambient Conditions

The SCALTEC SBA/SBC Balance is designed to provide reliable weighing results under normal ambient conditions in the laboratory and in industry.

Therefore, please observe the following when choosing a suitable place of installation:

- Set up the balance on a stable, even surface
- Protect the balance from aggressive chemical vapors
- Protect the balance from drafts that come from open windows and doors
- Avoid exposing the balance to extreme vibrations
- Avoid exposing the balance to extreme heat or to direct sunlight

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Information on Radio Frequency Interference

Warning!

This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instruction manual, may cause interference to radio communications. It has been tested and found to comply with the limits for a Class A computing device pursuant to Subpart J of Part 15 of FCC rules, which are designed to provide reasonable protection against such interference, when operated in a commercial environment. Operation of this equipment in a residential area is likely to cause interference, in which case the user, at his own expense, will be required to take whatever measures may be required to correct the interference. Einsetzen der geeichten Waage im eichpflichtigen Verkehr Vor dem Einsatz im eichpflichtigen Verkehr ist die Waage am Aufstellort zu justieren (s. Kapitel »Justieren« ab Seite 17).

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This balance is not allowed to be used for weighing goods intended for direct sale to the public. The type-approval certificate for verification applies only to non-automatic weighing instruments; for automatic operation with or without auxiliary measuring devices, you must comply with the regulations of your country applicable to the place of installation of your balance. A suitable thermometer and barometer are recommended for monitoring the ambient conditions.

For balances of accuracy class ①, a thermometer and barometer are recommended for monitoring the ambient conditions. The temperature range indicated on the verification ID label must not be exceeded during operation. The balance must warm up for at least 24 hours after initial connection to AC power or after a relatively long power outage.

The manufacturer complies with EC Directive No. 90/384/ EEC for non-automatic weighing instruments, which has been in effect since January 1, 1993, within the Single European Market as well as the accreditation of the quality management system of the manufacturer by Lower Saxony's Regional Administrative Department of Legal Metrology (Niedersächsische Landesverwaltungsamt-Eichwesen) from February 15, 1993.

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Important Note Concerning All Verified Balances Approved for Use as Legal Measuring Instruments:

Provided that an official lead seal is required for the verified balance, a control seal is affixed to the balance. Unauthorized attempts to remove this seal will irreversibly damage it. If you break the seal, the validity of the verification will become void, and you must have your balance re-verified.

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Connecting the SBC 22 Balance to the Electronic Box:

- Plug the cable into the socket of the electronics box.

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- Schirmring
- Unterschale
- Waagschale

Preparing the SBA 31/32/33 and the SBC 22/31/32/33 with a Draft Shield Chamber

Place the following components inside the chamber in the order given:

- Shield disk
- Pan support disk
- Weighing pan

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Preparing the SBA 41/42 and the SBC 41/42 with a Round Glass Draft Shield

 Place the shield disk on the balance. Turn the disk counterclockwise until it is securely fastened.

Place components listed below on the balance in the order given:

- Pan support
- Weighing pan
- Glass draft shield cylinder
- Draft shield cover

Preparing the SBA 53 and the SBC 53 with a Round Weighing Pan

 Place the shield disk on the balance. Turn the disk counterclockwise until it is securely fastened.

Place components listed below on the balance in the order given:

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Preparing the SBA 51/52/61/62 with a Rectangular Weighing Pan

- Place the weighing pan (1) on the balance

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Connecting the Balance to AC Power

The balance is powered by an AC adapter. Make sure that the voltage rating printed on this unit is identical to your local line voltage.

Plug the cord of the AC adapter into the AC jack (12); then plug the AC adapter into a wall outlet (mains).

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Leveling the Balance Using the Level Indicator

At the place of installation, level the balance using the leveling feet (5) so that the air bubble is centered within the circle of the level indicator (16).

Retract the two auxiliary feet located at the back of the balance.

Alle Waagen nach der Libelle einstellen: Anheben der Waage: Stellfüße rechtsherum drehen. Senken der Waage: Stellfüße linksherum drehen. Use the level indicator as a guide to level the balance: To lift the balance, extend the leveling feet (turn clockwise). To lower the balance, retract the leveling feet (turn counterclockwise).

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After retracting the rear auxiliary feet, extend them until they touch the surface on which the balance rests.

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Operating the Balance

Warmup Time

To deliver exact results, the balance must warm up for at least 30 minutes after initial connection to AC power or after a relatively long power outage. Only after this time will the balance have reached the required operating temperature.

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Turning the Display On or Off (Standby Mode)

Press the ON/OFF key (9) to turn the display on or off.

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Self-Test

After the balance has been turned on, an automatic self-test of the balance's electronic circuitry is performed. At the end of the self-test, a zero readout is displayed. This means that the balance is ready to operate.

Important Note Concerning Verified Balances Approved for Use as Legal Measuring Instruments:

For verified balances that have a verification scale interval "e" which is greater than the scale interval "d," the last digit on the display is bordered.

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Taring

A weight can be determined accurately only from a defined zero point. Press TARE **(8)** to zero the weight display. You can tare within the entire weighing range of the balance.

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Important Note Concerning Verified Balances Approved for Use as Legal Measuring Instruments:

The small circle in the weight display (on the left) shows that the balance is exactly tared to "0" (\pm 0.25 of a scale interval).

Simple Weighing

Place your sample on the weighing pan **(1)** to determine the weight. Read off the weight indicated on the display only after the weight unit "g" or a different unit selected appears as the stability symbol.

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During calibration, the sensitivity of the balance is adjusted to changes in the ambient conditions.

Calibrate your new balance at the place of installation after each warmup period and before the first weight measurement. In verified models the balance should be calibrated at least once a day.

The balance offers you various calibration functions. You can select these functions by setting the appropriate menu code.

All SBC models have a built-in calibration weight. Menu code selection: 1 9 3 *

The calibration weight is internally applied by servomotor and removed at the end of calibration.

Using Verified SBC Balances as Legal Measuring Instruments

Before using your balance as a legal measuring instrument, you must carry out an "Internal Calibration" operation at the place of installation after the warmup period.

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If any interference affects calibration, you may obtain a brief display of the error message "Err 02." In this case, tare and press the CAL key again.

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Important Note

During calibration, the weighing pan must be unloaded.

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External Calibration**

Menu code selection: 1 9 1*

Use only calibration weights with an accuracy equal to or better than the readability of your balance. Refer to the list of calibration weights under "Accessories" on page 73.

When a zero readout is displayed, press the CAL key **(10)**. This starts calibration. The calibration weight value in grams is then displayed.

Errors or interference at the start of the calibration routine are indicated by the error message "Err 02." If this is the case, tare and press the CAL key again when a zero readout appears.

Center the calibration weight on the weighing pan. The balance then calibrates automatically. When calibration is completed, the value of the calibration weight and the stability symbol "g" are displayed.

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- * = factory setting for SBC balances with a readability of ≤0.1 mg
- ** = not applicable to verified balances of accuracy class I

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Calibration Test for SBC Models with a Built-in Calibration Weight

Menu code selection: 1 9 4

Rather substantial changes in barometric pressure and temperature may affect the display response of the balance. To ensure that you obtain the full accuracy of the balance, even when you use the entire weighing range, we have implemented a calibration test function. It will help you quickly to decide whether or not you need to recalibrate your balance (for example, to maintain the same accuracy during long-term weighing series).

Unload the balance and tare.

When a zero readout is displayed, press the CAL **(10)** key. The built-in calibration weight is now internally applied by servomotor. At the same time, **CAL** is displayed. After the display has stabilized, the deviation of the momentary readout from the target weight (in grams only) is indicated.

If external interference affects the calibration test, you may obtain a brief display of the error message "Err 02." In this case, tare and then press the CAL key again.

Important Note

The balance should be calibrated if the deviation of the readout from zero is more than the reproducibility specified for your balance in the "Technical Specifications" starting on page 61.

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To calibrate the balance with the built-in weight: press the CAL key briefly

or

to quit the calibration test: hold down the CAL key for more than 2 seconds

Blocking the Calibration Functions

You can block the calibration functions by setting code 1 9 7 (when the menu access switch **(4)** is locked).

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Data Interface Port

If you wish to record weighing data using a SCALTEC printer, plug the printer connector into the interface port **(15)** of the balance.

Important Note Concerning Verified Balances Approved for Use as Legal Measuring Instruments:

When using the balance as a legal measuring instrument, you may connect it only to auxiliary measuring devices that are legally permitted. However, if you are not using the balance as a legal measuring instrument, you may connect it to any peripheral devices that you wish, e.g., additional remote displays, personal computer, etc. Such peripheral devices must be marked to indicate that they may not be connected when the balance is being used as a legal measuring instrument.

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Important Note

Make sure to unplug the balance from AC power before you connect or disconnect a peripheral device (printer or PC) to or from the interface port.

To print data on hard copy or output them on a computer screen, press the PRINT key (7).

For information about the data output parameters and data ID codes, see page 36. For details on the data interface (such as the data output or input formats, pin assignment, etc.), contact SCALTEC.

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Below-Balance Weighing

A port for a below-balance weighing hanger is located on the bottom of the balance.

For Balances with an Analytical Draft Shield Chamber:

To open the below-balance port, turn the cover plate.

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For Balances of the SBA/SBC Series (Except for SBA 61/62):

To open the below-balance port, remove the cover plate from the bottom of the balance.

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Now you can attach a sample using a suspension wire, for example. Common applications for below-balance weighing include density determination and immersing a sample in a special atmosphere (medium for reaction).

Important Note Concerning Verified Balances Approved for Use as Legal Measuring Instruments:

The below-balance weighing port may not be opened or used when an approved balance is being operated as a legal measuring instrument.

Important Note

When you use a below-balance weighing hanger, you must install a shield to protect against drafts.

Fastening an Antitheft Locking Device*

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To fasten an antitheft locking device, use the lug **(17)** on the rear panel of the balance.

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Troubleshooting Guide

Problem	Causes	Solution
No segments appear on the weight display (11)	 No AC power is available The AC adapter is not plugged in 	 Check the AC power supply Plug in the AC adapter
No segments appear on the weight display after adjustment/calibration?	 The surface on which the balance rests is not stable Internal stability has not been reached 	 Make sure that the ambient conditions are stable Prevent vibrations from affecting the surface on which the scale rests Close the draft shield
The weight display shows "H"	- The load exceeds the capacity of the balance	– Unload the balance
The weight display shows "L" or "Err 54"	 The weighing pan (1) and/or the pan support disk is not in place 	 Position the pan and/or the pan support disk (depending on the balance model)
The weight display briefly shows "Err 02"	 The display did not show a zero readout when the CAL key (10) was pressed to calibrate The balance is loaded 	 Press the TARE key (8); then press the CAL key again Unload the balance
The special code "�" does not go out on the weight display	 None of the keys has been pressed since the balance was turned on 	– Press a key
The weight readout changes constantly	 Unstable ambient conditions Too much vibration or the balance is exposed to a draft 	 Set up the balance in another area Access the menu to select the correct code for the weighing environment
The weight readout is obviously wrong	 The balance is not calibrated The balance was not tared before weighing The air bubble of the level indicator (16) is not centered within the circle 	 Calibrate the balance (see page 17) Tare before weighing Level the balance (see page 12)

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Servicing

Regular servicing by a SCALTEC service technician will extend the service life of your balance and increase the measuring accuracy. SCALTEC can offer you service contracts, with your choice of regular maintenance intervals ranging from 1 month to 2 years.

Safety Inspection

If there is any indication that safe operation of the balance with the AC adapter is no longer warranted, turn off the power and disconnect the equipment from AC power immediately. Lock the equipment in a secure place to ensure that it cannot be used for the time being.

In this case, notify your SCALTEC customer service office.

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The **factory settings** of the menu codes are identified by an "*" in this instruction manual. You can select the functions not identified by an "*" by setting the respective menu codes. You will find these codes in the section entitled "Balance Operating Parameters," which starts on page 31.

Important Note Concerning Verified Balances Approved for Use as Legal Measuring Instruments:

The balance operating menu can also be changed when the balance is being used as a legal measuring instrument. Codes that are not permitted for operation of the balance as a legal measuring instrument are blocked and cannot be selected. The balance operating menu on verified balances cannot be locked with the menu access switch ("-" not displayed).

Changing Menu Code Settings

There are three steps to changing a code:

- Accessing the menu
- Setting a code
- Confirming and storing this code

The keys have special functions for setting menu codes:

- CAL (10) = Increases a number by one with each press (the numbers change in cycles)
- TARE (8) = Confirms a code setting; stores a code setting; and leaves the menu
- PRINT (7) = Moves to the right stepwise, and to the left from the third number to the first one

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Accessing the Menu

- Press ON/OFF to turn off the balance.
 Turn the balance back on. While all segments are displayed, briefly press the TARE key (8).
- Release the TARE key when "I" is displayed.
- If "-" is displayed in front of the left-hand number, remove the protective cap from the rear panel of the balance, and move the menu access switch (4) in the direction of the arrow.
- Press CAL until "8" appears
- Press PRINT until the 2nd number of the code appears
- Press CAL until "5" appears
- Press PRINT until the 3rd number appears (when you move to this number, the previously set menu code will appear)
- Press CAL to select "4"
- Confirming a Menu Code Setting Press TARE to confirm the code you have just set (this is indicated by the "o" after the code).
- To **store** the new menu code setting, press TARE for more than 2 seconds.

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Important Note

The current code setting in the balance operating menu is identified by a small, superscript "o" after the last number. When you access the operating menu, the previously set code will be displayed after you have selected the righthand number.

If you would like to change several menu code settings, you do not have to press the TARE key after each change to leave the balance operating menu.

Please do not forget to relock the balance operating menu. The symbol "-" indicates that the menu is locked.

Undoing All Menu Code Changes - Reset Function

The reset function lets you undo all menu code changes, which means that you will obtain the original factory-set menu codes identified by an "*." To use this function, select code $9 - 1^{\circ}$.

Reset Function	Code
On	9 1
Off	9 2*

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* = factory setting; depends on the balance model in some cases

** = setting not applicable to verified balances of accuracy class \square

 not applicable to the verified balances approved for use as legal measuring instruments

Auto-Zero-Funktion	Code	Auto Zero Function	Code
Auto-Zero ein	1 6 1*	Auto Zero on	1 6 1*
Auto-Zero aus	1 6 2	Auto Zero off	1 6 2
Justier-/Kalibrier- und Linearisierfunktionen mit CAL	Code	Calibration and Linearization Functions Using CAL	Code
Externes Justieren	1 9 1*/**	External calibration	1 9 1*
Internes Justieren für SBC-Modelle	1 9 3*	Internal calibration for SBC models	1 9 3*
Kalibrieren für SBC-Modelle	1 9 4	Calibration test for SBC models	1 9 4
Internes Linearisieren bei SBC 22	196	Internal linearization for SBC 22	196
Justierfunktion gesperrt	1 9 7	Calibration functions blocked	1 9 7

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 $^{*}\,$ = factory setting; depends on the balance model in some cases ** = setting not applicable to verified balances of accuracy class 1

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The FUNCTION key (6) lets you switch back and forth (toggle) between two weight units. The 2nd unit is identified by the display symbol "R1."

	Code		
FUNCTION key blocked	2	1	1*
Mass unit conversion by toggling	2	1	2*

Weight Units

The "basic weight unit" is the unit in which your balance will weigh the moment you turn it on. This unit is selected for the first range. You can select a different unit for each of the two ranges by setting the appropriate menu codes.

	Symbol	Conversion Factor	Code						
	-	1 g =	1s	t Ra	inge		2n	d Ra	ange R1
Grams (o)	0	1.	1	7	1		3	1	1
Grams	g	1.	1	7	2*		3	1	2 *
Kilograms ¹)	kg	0.001	1	7	3		3	1	3
Carats**	ct	5.	1	7	4		3	1	4
Pounds**	lb	0.0022046226	1	7	5		3	1	5
Ounces**	OZ	0.035273962	1	7	6		3	1	6
Troy ounces**	ozt	0.032150747	1	7	7		3	1	7
Hong Kong taels*	*tl	0.02671725	1	7	8		3	1	8
Singapore taels*	* tl	0.02646063	1	7	9		3	1	9
Taiwanese taels**	* tl	0.02666666	1	7	10		3	1	10
Grains**	GN	15.43235835	1	7	11		3	1	11
Pennyweights**	dwt	0.643014931	1	7	12		3	1	12

	Symbol	Conversion Factor		Coc	de	е			
	-	1 g =	1s	t Rá	ange		2n	d R	ange R1
Milligrams ²)	mg	1000.	1	7	13		3	1	13 *
Parts per pound*	*0	1.1287667712	1	7	14		3	1	14
Chinese taels**	tl	0.02645547175	1	7	15		3	1	15
Mommes**	m	0.2667	1	7	16		3	1	16
Austrian carats**	k	5.	1	7	17		3	1	17
Tola**	t	0.0857333381	1	7	18		3	1	18
Baht**	b	0.06578947436	1	7	19		3	1	19
Mesghal**	m	0.217	1	7	20		3	1	20
Display Mode						Сос	de		
			1s	t Rá	ange		2n	d R	ange R1
Highest possible accuracy ³)		1	8	1*		3	2	1*	
Rounding factor 2	2**		1	8	3		3	2	3

1 8 4

1 8 5

3 2 4

3 2 5

= factory setting; depends on the balance model in some cases

- ** = not applicable to verified balances approved for use as legal measuring instruments
- ¹) = not applicable to verified balances of accuracy class \square

Rounding factor 5**

Rounding factor 10**

²) = not applicable to verified balances of accuracy class \square

³) = applicable only to verified balances approved for use as legal measuring instruments

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Baud Rate Co				
150 baud	5	1	1	
300 baud	5	1	2	
600 baud	5	1	3	
1,200 baud	5	1	4 *	
2,400 baud	5	1	5	
4,800 baud	5	1	6	
9,600 baud	5	1	7	
19,200 baud	5	1	8	
Parity	Code			
Mark parity	5	2	1	
Space parity	5	2	2	
Odd parity	5	2	3 *	
Even parity	5	2	4	
Number of Stop Bits	Со			
1 stop bit	5	3	1*	
2 stop bits	5	3	2	
Handshake Mode		de		
Software handshake	5	4	1	
Hardware handshake with 2 characters after CTS	5	4	2 *	
Hardware handshake with 1 character after CTS	5	4	3	

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Data Output Parameter	Сс	ode	
Print on request regardless of stability	6	1	1
Print on request after stability, with storage of the function	6	1	2 *
Print on request at stability, without storage of the function	6	1	3
Auto print regardless of stability	6	1	4
Auto print at stability	6	1	5
Auto Print	Сс	ode	
Start/stop auto print using the PRINT key	6	2	1
Auto print not stoppable	6	2	2 *
Data Output at Defined Intervals	Сс	ode	
1 display update	6	3	1*
2 display updates	6	3	2
Automatic Taring after Data Output	Сс	ode	
Data output without automatic taring	6	4	1 *
Data output with automatic taring	6	4	2
Automatic Output of the Application Parameters			
Off	7	1	1*
On	7	1	2

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Additional Functions

Menu Access Function

You can define the function of the menu access switch by setting the code for the balance operating menu to "accessible." The balance operating menu will then be accessible at all times, which means that you can change the menu codes at any time regardless of the setting of the menu access switch.

Access to Balance Operating Menu			
Accessible (codes can be changed)	8	1	1
Locked (codes can be read only)	8	1	2 *
Key Functions			
Accessible	8	3	1*
Blocked (except for ON/OFF)	8	3	2

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Menu Access Switch Function on Verified Balances of Accuracy Class ${\rm I\!I}$

The balance operating menu can also be changed when the balance is being used as a legal measuring instrument. Codes that are not permitted for operation of the balance as a legal measuring instrument are blocked and cannot be selected. The balance operating menu on verified balances approved for use as legal measuring instruments cannot be locked with the menu access switch ("--" not displayed). For balances of accuracy class (II), the menu access switch is adjusted to lock the "External Calibration" function after verification. To do this, the switch must be moved to the left. A control seal is affixed to the balance to lock the menu access switch. Unauthorized attempts to remove this seal will irreversibly damage it. If you break the seal, the validity of the verification will become void, and you must have your balance re-verified.

Power-on Mode			
(Power) off -> on <-> standby	8	5	1*
on <> standby	8	5	3
Automatic power-on	8	5	4

Application Programs

In addition to the functions implemented for weighing, SCALTEC balances offer you a variety of application programs.

To select an application program or a specific function, set the appropriate code in the menu.

Additional settings for running a particular program and displaying or printing the data on hard copy are listed in a table of codes.

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Tare Memory

Code 2 1 6

Press the FUNCTION key (6) to store the tare weight. The balance is now automatically tared so you can weigh starting with a zero readout. If you have stored a value in the tare memory, a "1" will be printed or output after the ID code of this value to identify it as such on the printout or on the computer screen.

"Tare - Net - Gross" Weights

Menu codes used in the example:	Code			
Tare memory Automatic output of all parameters	2 1 6 7 1 2			
With data ID code	722			
Print tare memory	7 3 2			
As an alternative, you can select this output parameter:				
Print net value	7 3 1*			

Step/Key		Re	adout		Data Out	out/Printo	out	
CF ¹), TARE	-		0.00 g					
Place conta on pan; press FUN	ainer CTION	+	22.65 g 0.00 g r	NET	т 1	+	22.65	g
Fill contair with sampl press PRIN	ner e; IT	+	250.24 g	NET	N 1	+	250.24	g
CF ¹), PRIN	IT	+	272.89 g		Ν	+	272.89	g
The data ID) codes				mean:			
т1	+	2	22.65	g	tare weig (weight v	ht stored alue)	in the memory	
N 1	+	25	50.24	g	net weigh	nt when ta	re weight is store	ed
N	+	27	2.89	g	gross we	ight = tare	e + net	

Application: Obtaining printouts of tare, net and gross weights

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Practical Example: Net Total

Menu codes used in the example:	Сс	Code			
Tare memory	2	1	6		
Automatic data output	7	1	2		
With data ID code	7	2	2		
Printout of individual components/tare weights	7	3	1*		
As an alternative, you can select this output parameter:					
Printout of the net total weight/total tare weight	7	3	2		
Application: Simple compounding and formulating of covers	lcompo	non	łc		

Application: Simple compounding and formulating of several components with additive storage, automatic taring and outputting of component weight

Step/Key	Re	adout	Data Output/Printout			
Place container on pan	+	22.65 g				
CF ¹), TARE		0.00 g				
Weigh in 1st component; store: FUNCTION	+	4.61 g 0.00 g NET	N 1	+	4.61	g
Weigh in 2nd component; store: FUNCTION	+	60.33 g 0.00 g NET	N 1	+	60.33	g
Weigh in and store additional components: FUNCTION		xxx.xx g Net	N 1	+	xxx.xx	g
Finish weighing in components; output total weight: CF ¹), PRINT	+	272.89 g	N	+	272.89	g
The data ID codes			mean:			
N1 + N +	6 27	0.33 g 2.89 g	compor total cor (net tota	nent weight mponents w al weight)	(net individual v veighed-in	weight)

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Weighing in Percent

Code 2 1 5

This application program enables you to obtain weight readouts in percent which are in proportion to a reference weight. The reference weight readout is stored as a menu-defined percentage (factory setting: 100%).

Changing the Reference Percentage

You can change the reference percentage in cycles. Choose from the following settings: 5, 10, 20, 50 and 100.

To activate the	Hold down FUNCTION (6) for more than 2 seconds until
change function:	"rEF 100%" appears on the display
To change the setting:	Briefly press FUNCTION
To store the percentage:	Hold down FUNCTION for more than 2 seconds to
	store it permanently (in the non-volatile memory) after
	you turn off the power
This setting is not cance	eled by the reset code $91^{\circ}!$

5			
Storage Parameter for the Reference Weight/Value	Со		
With full accuracy according to the internal resolution	3	5	1*
According to the display accuracy	3	5	2
Display Parameter for Readouts in Percent	Со	de	
Without a decimal place	3	6	1
With one decimal place	3	6	2 *
With two decimal places	3	6	3 *
With three decimal places	3	6	4

If the weight stored is too light to be displayed, the number of decimal places is automatically decreased.

* = factory setting; depends on the balance model in some cases

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Practical Example:	Determination	of the Residual	Weight in Percent

Menu codes used in the example:	Code				
Weighing in percent	2 1 5				
Reference % and weight	7 1 2				
With data ID code	722				

Application: Quick determination of the dry weight of a sample

Step/Key	Rea	adout	Data Outpu	t/Printout		
Place container on pan	+	22.65 g				
CF ¹), TARE		0.00 g				
Add prepared sample to container; FUNCTION) + +	4.61 g 100.0 %	pRef Wxx%	+ +	100 4.61	% g
For moisture analysis, press TARE at this point		0.0 %				
Remove container; dry sample		XX.X %				
Reweigh sample in container	+	72.5 %				
PRINT	+	72.5 %	Prc	+	72.5	%

Step/Key		Rea	dout		Data Outpu	ut/Printo	ut
Press CF to residual wei readout and reference va	clear ight Ilue	+	3.34 g				
PRINT		+	3.34 g		Ν	+	3.34 g
The data ID	codes				mean:		
pRef Wxx% Prc N	+ + + +		100 4.61 72.5 3.34	% g % g	reference p net referen calculated net residua	percentaq ce weigh percenta al weight	ge it xx% ige

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Counting

The counting program allows automatic conversion of weights into piece counts based on a reference sample weight. A weight readout is stored as a reference sample quantity (factory setting: 10 pcs = pieces).

Changing the Reference Sample Quantity

You can change the reference sample quantity in cycles. Choose from the following settings for the reference sample quantity: 5, 10, 20, 50 and 100.

To activate the	Hold down FUNCTION (6) for more than 2 seconds until
change function:	"rEF pcs" appears on the display
To change the setting:	Briefly press FUNCTION again
To store the quantity:	Hold down FUNCTION for more than 2 seconds to store
	it permanently (in the non-volatile memory) after you turn
	off the power
This setting is not can	celed by the reset code $9 = -1^{\circ}$

Storage Parameter for the Peference Sample Weight

Storage Parameter for the Reference Sample weight	Code			
With full accuracy according to the internal resolution	3	5	1*	
According to the display accuracy	3	5	2	

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Practical Example: Parts Counting		Press					
Menu codes used in the example:	Code	FUNCTION	rEF 20 pcs				
Counting Reference quantity and ref. weight With data ID code	2 1 4 7 1 2 7 2 2	To store: hold down FUNCTION	I				
Application: Counting bulk quantities of items with the same we	ight	2 seconds					
Step/Key Readout Data Output/Printout		Press FUNCTION					
Place container on pan + 22.65 g		to confirm reference piece count	+ 10 pcs	nRef wRef	+ + 5	10 .6546	pcs g
CF ¹), TARE 0.00 g		Fill container with					
Add 10 counted parts rEF 10 pcs		desired quantity					
Change the ref-		counting them)	+ 500 pcs				
erence sample atv. if desired:		PRINT	+ 500 pcs	Qnt	+	500	pcs
To use 20 as the reference qty., add 10 more parts to container.		Press CF to clear the weight readout and the reference value PRINT	+ 2827.35 g + 2827.35 g				
ref. qty. to 20,		The data ID codes		mean:			
hold down FUNCTION for > 2 seconds rEF 10 pcs		nRef + wRef + 5 Qnt + N + 28	10 p .6546 g 500 p 27.35 g	c s reference p reference w c s calculated net weight	iece count /eight piece count		
) = press CF for more than two seconds 							

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7 Animal Weighing/Averaging

Code: 217

Use this program to determine the weights of live animals or weights under unstable ambient conditions.

In this program, the balance calculates the weight as the average of a selectable number of individual weighing operations by the balance. These weighing operations are also known as "subweighing operations."

During averaging, the number of individual subweighing operations left is shown in the application display in a "countdown" mode. Once all subweighing operations have been performed, the calculated mean value is indicated as a stable readout on the weight display.

Manual or Automatic Start Mode

Depending on the menu code you select, animal weighing will be started automatically or manually by pressing a key. If you select the "automatic mode," you will have to press FUNCTION to start the program for averaging the first weight. You can press CF to interrupt a weight measurement in progress at any time. In the "automatic mode," the symbol "**AUTO**" will be displayed during weight measurement. The result is locked into the display. The "mouse" symbol or "**AUTO**" flashes during this time. The readout will stop flashing after you have unloaded the balance, and the next weight will be displayed.

Animal weighing/averaging will start by:			
Manual mode	3	8	1
Automatic mode	3	8	2*

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Delayed Start Mode

A rule of thumb to go by for selecting the right setting to weigh animals is: the more active an animal is, the greater the difference must be between two successive subweights measured.

Depending on individual requirements, starting the averaging operation can be delayed either in the automatic or manual mode until the animal you are weighing has calmed down to a certain degree.

In this case, the start criterion is defined by the difference between two successive subweights measured.

If the animal moves, the start criterion is not met; therefore, averaging will not start.

Once the animal has calmed down, the program checks whether two measured subweights are within the previously selected range. If so, the actual averaging operation will be started.

Delay start	Сс	ode		
Slight	3	7	1	
Average	3	7	2 *	
Considerable	3	7	3	

During averaging, the number of subweighing operations left to perform is shown on the weight display (countdown mode).

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Changing the Number of Subweighing Operations

You can change the number of subweighing operations used to average a weight. Change this number in cycles. You can choose from the following settings: 5, 10, 20, 50 and 100.

To activate the change function:	Hold down FUNCTION (6) for more than
0	2 seconds until "rEF 10" appears on the display
	If you enter a wrong number, press CF to clear it.
To change the setting:	Briefly press FUNCTION
To store the number:	Hold down FUNCTION for more than 2 seconds
	to store this number permanently (in the
	non-volatile memory) after you turn off the power

This setting is **not** canceled by the reset code $9 - -1^{\circ}!$

Storage Threshold for the Automatic Start Mode

To obtain an added measure of reliability in the automatic mode in order to avoid an "incorrect start," a weight value must correspond to a minimum load of 100 display increments.

Once averaging has been completed, the program will stop until the balance is unloaded to half the value (50 display increments) of the storage threshold.

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Practical Example: Animal Weighing with Automatic	Start Mode	Step/Key	Readout	Data Output/Printo	out
Menu codes used in the example:	Place animal				
Animal weighing Delay start until dif. is average	2 1 7 3 7 2*	weighing bowl on balance			
Automatic start mode	3 8 2*	CF ¹), TARE	0.00 g		
Automatic output of all parameters With data ID code	7 1 2 7 2 2	Change number - of subweighing	Ŭ		
As an alternative you can select this code: Code		operations; e.g.,			
Other delay start	37 x				
Manual start mode	3 8 1	- FUNCTION for more			
You can change the number of subweighing operations.		than 2 seconds	ref 10		
If you do not enter a number, averaging will be done based on the shown on the display (factory setting: 10 subweighing operatio	ne number Ins).	Press FUNCTION	ref 20		
Once you press FUNCTION the first time, all further weighing o will be performed automatically after they have met the start crit	To store: hold down FUNCTION				
Application: Automatic weighing of animals based on 20 subv	for more than 2 seconds	ref 20			
		Place 1st animal in bowl			
		FUNCTION	20 465.20 g	mdef +	20
				x-Net +	465.20 g

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Application:

Use of the balance in quality assurance systems in compliance with ISO, GLP, GMP, EN, etc.

The balance can record all completed calibration and adjustment operations and print out data in compliance with the requirements of Good Laboratory Practice (GLP). The balance, interfaced with a data printer or a computer, creates a document that records the date, time, serial number and model number, making it possible to clearly trace data to the balance that generated it and the time at which it was generated.

Select the GLP/GMP-compliant printout or record mode by setting the respective code in the balance operating menu:

GLP/GMP-compliant printout/record mode	Code	
Off	8 10	1*
Only for the calibration and linearization functions	8 10	2
Always on	8 10	3

The following menu code setting must be selected in order to obtain a GLP/GMP-compliant printout/record:

With data ID code	7	2	2	
-------------------	---	---	---	--

7 2 2

Important Note

GLP/GMP printouts/records will not be generated if the factory setting, code 7 2 1, is selected. In addition, do not select the "Auto print" data output parameter (code 6 1 4 or 6 1 5).

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Operating the Balance with a Device that Has a GLP/GMP Printing or Recording Capability:

You can connect a special SCALTEC printer (order no. SPR 01) to SCALTEC balances for generating GLP/GMP-compliant printouts.

The GLP/GMP-compliant SCALTEC printer provides you with the following functions:

- Turns the GLP/GMP function on or off
- Date/time
- ID no. for identifying a workstation/operator
- Printouts with balance-specific data

Printout/Record for Calibration and Linearization Functions

A printout or record is generated at the end of the following functions:

- All calibration and linearization functions

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The printout can have the following lines:

SCA Model	ALTEC SBC22 70500146
I d	
Date : Start: Cal. :	21-Jun-97 10:05 Extern
Set. :	200.00000 g
End : Name :	10:05

Balance family and manufacturer
Balance model
Serial number of the balance
Space for entering a workstation or operator ID
Current date
Time at which the application started
Calibration mode (in this case, "external calibration")
Calibration weight value (only for "external calibration")
End of application
Space for signature of the operator responsible

Additional information for other calibration modes:

Cal. :	Intern		: Interna
Lin. :	Intern		: Interna
Lin. :	Extern		: Extern
Stat.:	Complete		: Status
			or line
Cal. :	Test		: Calibra
Diff.:-	0.00004	g	: Data m
			the cal

ternal calibration ternal linearization kternal linearization tatus comment for calibration linearization functions alibration test ata measured during e calibration test

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Data Printout/Record

To have a data record printed out, perform the following:

- Press PRINT (after you have turned on the balance or cleared a function by pressing CF) to output the printout heading
- Press PRINT to output additional data
- To end printout generation and recording of data, press CF

Generation of a GLP/GMP-compliant printout or record is also ended when a calibration operation is started.

For application programs, the reference data (parameters) can be included in the printout/record.

Automatic Output of the Reference Data (Parameters)	Code		
Off	7	1	1*
Reference %/qty. and reference weight	7	1	2

To have a data record printed out, perform the following:

 Press FUNCTION to output the printout heading and reference data (the reference data will be stored at the same time)

If you input and store new reference values while a GLP/GMP-compliant record is being printed out, the new reference data will be output.

- Press PRINT to output weighing data
- Press CF to end printout generation (this also ends once a calibration or linearization operation has been started)
- Then press CF to clear the reference data for the application programs

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The printout of the data record can have the following lines:

SCALTEC Model SBC22 S/N 70500146 Id	: Balance family and manufacturer : Balance model : Serial number of the balance : Space for entering a workstation or operator ID
Date: 10-Jun-97 Start: 11:15 Ser.: nRef + 10	: Current date : Time at which the application started : Space for entering a project number p c s : Reference data (in this case, "counting" – see also page 46)
Qnt + 500	p c s : Measured values (in this case, "calculated piece count")
End : 12:15 Name :	: End of application : Space for signature of the operator responsible

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Technical Specifications

Model		SBA 31	SBA 32	SBA 33	SBA 41	SBA 42
Weighing capacity	g	220	120	65	320	120
Readability	g	0.0001	0.0001	0.0001	0.001	0.001
Tare range (by subtraction)	g	- 220	-120	-65	-320	-120
Reproducibility (standard deviation)*	g	≤±0.0001	≤±0.0001	≤±0.0001	≤±0.003	≤±0.003
Linearity	g	≤±0.0002	≤±0.0002	≤±0.0002	≤±0.002	≤±0.002
Response time (average)	S	≤2	≤2	≤2	≤1.5	≤1.5
Adaptation to ambient conditions and application requirements		By selection of 1 of 4 optim	mized filter levels			
Display update (depends on the filter level selected)	S	0.1 - 0.4				
Allowable ambient operating temperature	°C	+5+40	+5+40	+5+40	0+40	0+40
Operating temperature range	°C	+10+30				
Sensitivity drift within +10 +30 °C	/°C	$\leq \pm 2 \cdot 10^{-6}$				
Pan size	mm	Ø 80	Ø 80	Ø 80	Ø 115	Ø 115
Weighing chamber height (effect. dimens.)	mm	225	225	225	-	-
Net weight, approx.	kg	5.4	5.4	5.4	2.8	2.8
Selectable weight units		g, kg, ct, lb, oz, ozt, tlh, tls	s, tlt, GN, dwt, mg, /lb, tlc, m	nom, K, tol, bat, MS		
Selectable application programs		Mass unit conversion by t	toggling; tare memory; net t	total; weighing in percent; a	nimal weighing	
Automatic zero tracking function (can be turned off by menu code)		Standard feature				
External calibration weight value (of at least accuracy class)	g	200 (E2)	100 (E2)	50 (E2)	200 (F1)	100 (F1)
Power requirements	V~	Via AC adapter 230 VAC o	or 115 VAC, -20% +15%	6		
Frequency	Hz	48-63				
Power consumption (average)	VA	16: maximum; 8: average	е			
Hours of operation with fully charged SBP 01 external battery pack, approx.:	h	32				
Built-in interface port		RS-232 C-S/V24-V28; 7- software/hardware hands	bit; -even, -mark, -odd, -sp hake	bace; Transmission rates: 15	50 19,200 baud; 1 or 2 sto	op bits;
Standard features/equipment supplied:						
Analytical draft shield chamber / Dust cover		Х	Х	Х	-	-
Round glass draft shield		-	-	-	Х	Х
AC adapter, varies acc. to country		Х	Х	Х	Х	Х
Level indicator		Х	Х	Х	Х	Х
Hanger for below-balance weighing		Х	х	Х	Х	X

* = standard deviation of the reproducibility acc. to DIN 8120, Part 3

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Technical Specifications

Model		SBA 51	SBA 52	SBA 53	SBA 61	SBA 62
Weighing capacity	g	3200	2220	620	7200	4200
Readability	g	0.01	0.01	0.01	0.1	0.1
Tare range (by subtraction)	g	-3200	-2220	-620	-7200	-4200
Reproducibility (standard deviation)*	g	≤±0.03	≤±0.03	≤±0.03	≤±0.2	≤±0.2
Linearity	g	≤±0.02	≤±0.02	≤±0.01	≤±0.2	≤±0.1
Response time (average)	S	≤1.5	≤1.5	≤1.5	≤1	≤1
Adaptation to ambient conditions and application requirements		By selection of 1 of 4 optim	mized filter levels			
Display update (depends on the filter level selected)	S	0.1-0.4				
Allowable ambient operating temperature	°C	0+40				
Operating temperature range	°C	+10+30				
Sensitivity drift within +10 +30 °C	/°C	$\leq \pm 2 \cdot 10^{-6}$	$\leq \pm 2 \cdot 10^{-6}$	≤±2·10 ⁻⁶	≤±5 · 10 ⁻⁶	$\leq \pm 5 \cdot 10^{-6}$
Pan size	mm	180 x 180	180 x 180	Ø 115	210 x 180	210x180
Net weight, approx.	kg	3.5	3.5	1.9	2.8	2.8
Selectable weight units		g, kg, ct, lb, oz, ozt, tlh, tls	s, tlt, GN, dwt, mg, /lb, tlc, m	nom, K, tol, bat, MS		
Selectable application programs		Mass unit conversion by I	oggling; tare memory; net t	otal; weighing in percent; a	nimal weighing	
Automatic zero tracking function (can be turned off by menu code)		Standard feature				
External calibration weight value (of at least accuracy class)	g	2000 (F1)	2000 (F1)	500 (F2)	5000 (F2)	2000 (F2)
Power requirements	V~	Via AC adapter 230 VAC o	or 115 VAC; –20% +15%	ò		
Frequency	Hz	48-63				
Power consumption (average)	VA	16: maximum; 8: average	5			
Hours of operation with fully charged SBP 01 external battery pack, approx.	h	25	25	25	45	45
Built-in interface port		RS-232 C-S/V24-V28; 7- software/hardware hands	bit; -even, -mark, -odd, -sp hake	ace; Transmission rates: 15	019,200 baud; 1 or 2 sto	op bits;
Standard features/equipment supplied:						
AC adapter, varies acc. to country		Х	Х	Х	Х	Х
Level indicator		Х	Х	Х	Х	Х
Hanger for below-balance weighing		Х	Х	Х	-	-

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Technical Specifications Verified Balances with EC Type Approval

Model		SBC 22
Accuracy class		
Weighing capacity Max.*	g	60/210
Scale interval d*	mg	0.01/0.1
Verification scale interval e*	g	0.001
Minimum capacity Min.*	g	0.001
Tare range (by subtraction)	g	-210
Range of use according to Directive No. 90/384/EEC*	g	0.001–210
Response time (average)	S	≤12/3
Adaptation to ambient conditions and application requirements		By selection of 1 of 4 optimized filters
Display update (depends on the filter level selected)	S	0.2-0.4
Allowable ambient operating temperature	°C	+5+40
Operating temperature range	°C	+15+25
Pan size	mm	Ø 80
Weighing chamber height (effect. dimens.)	mm	225
Net weight approx.	kg	6,6
Selectable weight units		g, mg
Selectable application programs		Mass unit conversion by toggling; tare memory; net total; weighing in percent; animal weighing
Automatic zero tracking function (can be turned off by menu code)		Standard feature
Power requirements	V~	Via AC adapter 230 VAC or 115 VAC, -20% +15%
Frequency	Hz	48-63
Power consumption (average)	VA	16: maximum; 8: average
Hours of operation with fully charged SBP 01 external battery pack, approx.:	h	25
Built-in interface port		RS-232 C-S/V24-V28; 7-bit; -even, -odd; Transmission rates: 150 19200 baud; 1 or 2 stop bits;
		software/hardware handshake
Standard features/equipment supplied:		
Dust cover		X
AC adapter, varies acc. to country		X
Built-in calibration weight		X
Level indicator		x

* = EC Directive No. 90/384/EEC for Non-Automatic Weighing Instruments applicable in the E.U. and within the European Economic Area

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Technical Specifications

Model		SBC 31	SBC 32	SBC 33	SBC 41	SBC 42		
Accuracy class								
Weighing capacity Max.*	g	220	120	65	320	120		
Scale interval d*	mg	0.1	0.1	0.1	1	1		
Verification scale interval e*	g	0.001	0.001	0.001	0.01	0.01		
Minimum capacity Min.*	g	0.01	0.01	0.01	0.02	0.02		
Tare range (by subtraction)	g	-220	-120	-65	-320	-120		
Range of use according to Directive No. 90/384/EEC*	g	0.01-220	0.01-120	0.01–65	0.02-320	0.02-120		
Response time (average)	S	≤2	≤2	≤2	≤1.5	≤1.5		
Adaptation to ambient conditions and application requirements		By selection of 1 of 4 optimized filters						
Display update (depends on the filter level selected)	S	0.1-0.4						
Allowable ambient operating temperature	°C	+5+40	+5+40	+5+40	0+40	0+40		
Operating temperature range	°C	+15+25	+15+25	+15+25	+10+30	+10+30		
Pan size	mm	Ø 80	Ø 80	Ø 80	Ø 115	Ø 115		
Weighing chamber height	mm	225	225	225	-	-		
Net weight approx.	kg	5.3	5.3	5.3	4.7	4.7		
Selectable weight units		g, mg	g, mg	g, mg	g, kg	g, kg		
Selectable application programs		Mass unit conversion by toggling; tare memory; net total; weighing in percent; animal weighing						
Automatic zero tracking function (can be turned off by menu code)		Standard feature						
Power requirements	V~	Via AC adapter 230 VAC or 115 VAC, –20% +15%						
Frequency	Hz	48-63						
Power consumption (average)	VA	16: maximum; 8: average						
Hours of operation with fully charged SBP 01 external battery pack, approx.: h		32						
Built-in interface port		RS-232 C-S/V24-V28; 7-bit; -even, -odd; Transmission rates: 15019200 baud; 1 or 2 stop bits; software/hardware handshake						
Standard features/equipment supplied:								
Analytical draft shield		Х	Х	Х	-	-		
Dust cover		Х	Х	Х	-	-		
Round glass draft shield		-	-	-	Х	Х		
Built-in calibration weight		х	х	Х	х	х		
AC adapter, varies acc. to country		Х	Х	Х	Х	Х		
Level indicator		х	Х	Х	Х	Х		

* = EC Directive No. 90/384/EEC for Non-Automatic Weighing Instruments applicable in the E.U. and within the European Economic Area

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Technical Specifications

Model		SBC 51	SBC 52	SBC 53	SBC 61	
Accuracy class					I	
Weighing capacity Max.*	g	3200	2200	620	6200	
Scale interval d*	g	0.01	0.01	0.01	0.1	
Verification scale interval e*	g	0.1	0.1	0.1	1	
Minimum capacity Min.*	g	0.5	0.5	0.5	5	
Tare range (by subtraction)	g	-3200	-2200	-620	-6200	
Range of use according to Directive No. 90/384/EEC*	g	0.5-3200	0.5-2200	0.5–620	5–6200	
Response time (average)	S	≤1.5	≤1.5	≤1.5	≤1	
Adaptation to ambient conditions and application requirements		By selection of 1 of 4 optir	mized filters			
Display update (depends on the filter level selected)	S	0.1-0.4				
Allowable ambient operating temperature	°C	0+40				
Operating temperature range	°C	+10+30				
Pan size	mm	180 x 180	180 x 180	Ø 115	180 x 180	
Net weight approx.	kg	4.4	4.2	4.7	4.4	
Selectable weight units		g, kg				
Selectable application programs		Mass unit conversion by t	oggling; tare memory; net t	otal; weighing in percent; ar	nimal weighing	
Automatic zero tracking function (can be turned off by menu code)		Standard feature				
Power requirements	V~	Via AC adapter 230 VAC or 115 VAC, –20% +15%				
Frequency	Hz	48-63				
Power consumption (average)	VA	16: maximum; 8: average				
Hours of operation with fully charged SBP 01 external battery pack, approx.:	h	25				
Built-in interface port		RS-232 C-S/V24-V28; 7-bit; -even, -odd; Transmission rates: 15019200 baud; 1 or 2 stop bits; software/hardware handshake				
Standard features/equipment supplied:						
AC adapter, varies acc. to country		х	х	х	х	
Built-in calibration weight		Х	Х	Х	х	
Level indicator		х	Х	Х	x	

* = EC Directive No. 90/384/EEC for Non-Automatic Weighing Instruments applicable in the E.U. and within the European Economic Area

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Council Directive 89/336/EU "Electomagnetic compatibility (EMC)"

Important Note:

The operator shall be responsible for any modifications to SCALTEC equipment and for any connections of cables or equipment not supplied by SCALTEC and must check and, if necessary, correct these modifications and connections. On request, SCALTEC will provide information on the minimum operating specifications (in accordance with the Standards listed above for defined immunity to interference).

73/23/EU "Electrical equipment designed for use within certain voltage limits"

If you use electrical equipment in installations and under ambient conditions requiring higher safety standards, you must comply with the provisions as specified in the applicable regulations for installations in your country.

For information on the current legal requirements of your country, please contact your local SCALTEC customer service office.

Additional Directive for Weighing Instruments Used in Legal Metrology: Directive 90/384/EU "Non-Automatic Weighing Instruments"

This directive also regulates the performance of the EC verification by the manufacturer, provided that an EC Type-Approval Certificate has been issued (see page 72) and the manufacturer has been accredited by an officer of a Notified Body registered at the Commission of the European Community for performing such verification. On February 15, 1993, the manufacturers received accreditation for performing verification.

Subsequent Verifications within the European Countries

The validity of the verification will become void in accordance with the national regulations of the country in which the balance is used. For information on verification and legal regulations currently applicable in your country, and to obtain the names of the persons to contact, please contact your local SCALTEC office, dealer or service center.

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The weighing instruments were initially verified under number 0111 by the manufacturer at the place of manufacture.

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