

UltiMate 3000

System Installation and Application



Revision: 2.1

Date: September 2005

Now sold under the Thermo Scientific brand



© 2005 Dionex



Table of Contents

| 1 | Intr | oduction | 1 |
|---|-------|---|----|
| | 1.1 | General | 1 |
| | 1.2 | Safety Precautions | |
| | 1.3 | Unpacking | 2 |
| | 1.4 | Supported System Configurations | |
| | 1.5 | Intended Use | |
| 2 | Inst | allation | 7 |
| | 2.1 | Facility Requirements | 7 |
| | 2.1 | Rear Panel Connectors | |
| | 2.2.1 | | |
| | 2.2.2 | | |
| | 2.3 | Opening the Front Panel | |
| | 2.4 | Preparing the System | |
| 3 | Ann | lication Examples and Fluid Connections | 13 |
| _ | 3.1 | Manual Injection (Basic Configuration) | |
| | 3.1 | Direct Injection | |
| | 3.3 | Preconcentration | |
| | 3.4 | 2D LC Salt Plug | |
| | 3.5 | Comprehensive 2D LC | |
| | 3.6 | Parallel LC | |
| 4 | Auto | omated Control by Chromeleon | 45 |
| | 4.1 | General | |
| | 4.2 | Assigning a Flow Splitter to a Pump | |
| | 4.3 | Recording the Pump Pressure as a Separate Channel | |
| | 4.4 | Synchronizing the Autosampler with a Pump | |
| | 4.5 | Displaying the Degasser-Related Parameters | |
| _ | | | |
| 5 | Soft | Key Menus | 55 |
| 6 | Tro | ıbleshooting | 57 |
| 7 | , D | A. B. C. A. | (1 |
| 7 | Kou | tine Maintenance | 61 |
| 8 | Acce | essories and Spare Parts | 63 |
| | 8.1 | Basic Configuration Kits | 64 |
| | 8.1.1 | | |
| | 8.1.2 | Biocompatible Devices | 66 |
| | 8.2 | Direct Injection Kits | |
| | 8.2.1 | Standard Devices | 68 |
| | 8.2.2 | Biocompatible Devices | 70 |

| 10 | Index | 80 |
|----|-----------------------------|----|
| 9 | Technical Information | 87 |
| | 8.6.2 Biocompatible Devices | 83 |
| | 8.6.1 Standard Devices | |
| 8 | 3.6 Parallel LC Kits | |
| | 8.5.2 Biocompatible Devices | 80 |
| | 8.5.1 Standard Devices | |
| 8 | 3.5 Comprehensive 2D LC Kit | 79 |
| | 8.4.2 Biocompatible Devices | |
| | 8.4.1 Standard Devices | 77 |
| 8 | 3.4 2D LC Salt Plug Kit | 77 |
| | 8.3.2 Biocompatible Devices | 74 |
| | 8.3.1 Standard Devices. | 72 |
| 8 | 3.3 Preconcentration Kits | 72 |

1 Introduction

1.1 General

When working with analytical instrumentation, you should know the potential hazards of using chemical solvents.

This document is intended as a supplementary guide to the *Operating Instructions* for the individual modules of an UltiMate 3000 system (also called "the system" throughout the manual), providing information about the installation and application of the entire system. The descriptions in the manual apply to the standard (stainless steel) and biocompatible versions of the system modules. If some detail applies to only one version, the version is identified by name.

In order to obtain a full understanding of the system, Dionex recommends that you review this document and the *Operating Instructions* for the system modules thoroughly before beginning installation and operation of the system.

At various points throughout the document, messages of particular importance are indicated by certain symbols:

Please note: Indicates general information intended to optimize the performance of

the instrument.

Important: Indicates that failure to take note of the accompanying information may

result in damage to the instrument.

Important: Indique que ne pas tenir compte de l'information jointe peut

endommager l'instrument.

Warning: Indicates that failure to take note of the accompanying information may

result in personal injury.

Avertissement: Indique que ne pas tenir compte de l'information jointe peut entraîner

des blessures corporelles.

The information contained in this document is subject to change without notice and should not be construed as a commitment by Dionex. Dionex assumes no responsibility for any errors that may appear in this document. This document is believed to be complete and accurate at the time of publication. In no event shall Dionex be liable for incidental or consequential damages in connection with or arising from the use of this document.

CHROMELEON® is registered trademark and UltiMateTM is a trademark of Dionex. All other trade or company names mentioned are subject to the copyright and the property and trademark rights of the respective companies.

All rights reserved, including those for photomechanical reproduction and storage on electronic media. No part of this publication may be copied or distributed, transmitted, transcribed, stored in a retrieval system, or transmitted into any human or computer language, in any form or by any means, electronic, mechanical, magnetic, manual, or otherwise, or disclosed to third parties without the express written permission of Dionex.

1.2 Safety Precautions

To avoid the possibility of personal injury and/or damage to the UltiMate system, observe the safety precautions for the individual system modules (\rightarrow Operating Instructions for the instruments).

1 Important: When connecting the capillaries, make sure that the connectors are free from contaminants. Even minute particles may cause damage to the flow splitter, flow control valve, and column.

Important: Lorsque vous connectez les capillaires, assurez-vous que les raccords sont exempts de tout contaminant. Même d'infimes particules peuvent causer des dommages au système (ex. diviseur de débit, vanne de régulation de débit et colonne).

Please note: Use only standard HPLC solvents (HPLC-grade: 0.2 μm, filtered) and buffers that are compatible with components in the flow path of the flow manager. Note the special properties of the solvents such as viscosity, boiling point, UV absorption (UV/VIS detector), and refractive index (refractive index detector).

Buffer concentration: Typically up to 1 mol/L.

1.3 Unpacking

All electrical and mechanical components of the system modules are carefully tested before the instruments are shipped from the factory. After unpacking, please inspect the modules for any signs of mechanical damage that may have occurred during transit.

Please note: Immediately report any shipping damage to both the incoming carrier and Dionex. Shipping insurance will compensate for the damage only if

reported immediately.

Please note: Keep the original shipping container and the packing material. They provide excellent protection for the instrument in case of future transit. Shipping the instrument in any other packaging automatically voids the product warranty.

Unpack the instruments as follows:

- Place the shipping container on the floor and remove the accessories kit and the power cord.
- Grasp the instrument by the sides. Slowly and carefully, pull unit out of the shipping container and place it on a stable surface.

To prevent the instrument from falling, always lift the unit by the sides. Do not lift the unit by the packaging material or the front panel.

⚠ Important: Afin d'empêcher l'instrument de tomber, saisissez-la par les côtés. Ne soulevez l'instrument à l'aide du matériau d'emballage ou par la porte du panneau avant.

- Remove the foam inserts.
- Check off the contents of the accessories kit pack against the accessories list.

⚠ Important: To prevent damage to the instrument when lifting or moving, always hold

it by the bottom or sides.

Important: Afin d'éviter d'endommager l'instrument lorsque que vous le soulevez ou

le déplacez, saisissez-la toujours par le bas ou les côtés de l'unité.

1.4 Supported System Configurations

For an overview of the currently available system modules, refer to the table below. If you have any questions, do not hesitate to contact your Dionex sales representative or distributor.

Please note: For more information about the configuration of the individual modules, refer to the related *Operating Instructions*.

| Application | UltiMate 3000 System Module | Part No. |
|-------------|---|------------------------|
| | SRD-3600—Solvent rack with 6 channel vacuum degasser (intended, e.g., for use with an LPG-3600) —or— | 5035.9230 |
| | SRD-3300 Solvent rack with 3 channel vacuum degasser (intended, e.g., for use with an LPG-3300) | 5035.9240 |
| | LPG-3600—Ternary dual low-pressure micro gradient pump (with loading pump and micro pump in one enclosure) Standard pump (LPG-3600M) Biocompatible pump (LPG-3600MB) —or— | 5035.0035 5037.0035 |
| | LPG-3300—Ternary low-pressure micro gradient pump) Standard pump (LPG-3300M) Biocompatible pump (LPG-3300MB) | 5035.0040 5037.0040 |
| | UVD-3000—UV detector with nano flow cell | 5073.0010 |
| Nano HPLC | FLM-3100(B) with two 10-port micro switching valves) Standard flow manager (FLM-3100) Biocompatible flow manager (FLM-3100B) —or— | 5720.0010 5721.0010 |
| | FLM-3200(B) with 10-port micro and 6-port nano switching valves Standard flow manager (FLM-3200) Biocompatible flow manager (FLM-3200B) —or— | 5720.0020 5721.0020 |
| | FLM-3300(B) with one 10-port micro switching valve Standard flow manager (FLM-3300) Biocompatible flow manager (FLM-3300B) | 5720.0030 5721.0030 |
| | Note: These flow managers are fitted with a splitter cartridge for a split ratio of 1:1000. | |

| Application | UltiMate 3000 System Module | Part No. |
|--------------------|---|------------------------|
| | UVD-3000—UV detector with capillary flow cell | 5073.0020 |
| Capillary HPLC | FLM-3100(B) with two 10-port micro switching valves) Standard flow manager (FLM-3100) Biocompatible flow manager (FLM-3100B) —or— | 5720.0015 5721.0015 |
| | FLM-3200(B) with 10-port micro and 6-port nano switching valves Standard flow manager (FLM-3200) Biocompatible flow manager (FLM-3200B) —or— FLM-3300(B) with one 10-port micro switching valve | 5720.0025 5721.0025 |
| | Standard flow manager (FLM-3300) Biocompatible flow manager (FLM-3300B) | 5720.0035 5721.0035 |
| | Note: These flow managers are fitted with a splitter cartridge for a split ratio of 1:100. | |
| | UVD-3000—UV detector with micro flow cell | 5073.0030 |
| Micro HPLC | FLM-3100(B) with two 10-port micro switching valves) Standard flow manager (FLM-3100) Biocompatible flow manager (FLM-3100B) —or— | 5720.0018 5721.0018 |
| | FLM-3200(B) with 10-port micro and 6-port nano switching valves Standard flow manager (FLM-3200) Biocompatible flow manager (FLM-3200B) —or— | 5720.0028 5721.0028 |
| | FLM-3300(B) with one 10-port micro switching valve Standard flow manager (FLM-3300) Biocompatible flow manager (FLM-3300B) | 5720.0038 5721.0038 |
| | Note: These flow managers are fitted with a splitter cartridge for a split ratio of 1:6. | |
| | UVD-3000—UV detector with nano flow cell | 5073.0010 |
| Monolithic HPLC | FLM-3100(B) with two 10-port micro switching valves) Standard flow manager (FLM-3100) Biocompatible flow manager (FLM-3100B) —or— | 5720.0015 5721.0015 |
| | FLM-3200(B) with 10-port micro and 6-port nano switching valves Standard flow manager (FLM-3200) Biocompatible flow manager (FLM-3200B) —or— | 5720.0025 5721.0025 |
| | FLM-3300(B) with one 10-port micro switching valve Standard flow manager (FLM-3300) Biocompatible flow manager (FLM-3300B) | 5720.0035 5721.0035 |
| | Note: These flow managers are fitted with a splitter cartridge for a split ratio of 1:100. | |
| | WPS-3000(B)—Well-plate autosampler with temperature control Standard autosampler (WPS-3000) Biocompatible autosampler (WPS-3000B) —or— | 5820.0010 5821.0010 |
| | WPS-3000T(B)—Well-plate autosampler without temperature control Standard autosampler (WPS-3000) Biocompatible autosampler (WPS-3000B) | 5820.0020 5821.0020 |

1.5 Intended Use

The UltiMate 3000 system is designed to perform equally well as system for routine analyses or as a sophisticated research instrument for use in capillary, nano, and micro HPLC (high performance liquid chromatography) applications.

The system is controlled by the **Chromeleon** Chromatography Management System. However, it can also be operated with other data systems, such as Analyst[®] (Applied Biosystems/MDS Sciex), HyStarTM (Bruker Daltonics), or Xcalibur[®] (Thermo Electron Corporation). To do so, the appropriate add-on module to Chromeleon is required.

Please note that the system may be operated only using the accessories originally supplied with the system modules and within its technical specification (\rightarrow page 87).

If there is any question regarding appropriate usage, contact Dionex before proceeding.

Dionex cannot be held liable for any damage, material or otherwise, resulting from inappropriate or improper use of the system.

| UltiMate 3000 - System Installation and Apple | ication | |
|---|---------|--|
| | | |
| | | |
| | | |
| | | |
| | | |
| | | |
| | | |
| | | |
| | | |
| | | |
| | | |
| | | |
| | | |
| | | |
| | | |
| | | |
| | | |
| | | |
| | | |
| | | |
| | | |
| | | |
| | | |
| | | |
| | | |
| | | |
| | | |
| | | |
| | | |
| | | |
| | | |
| | | |
| | | |
| | | |
| | | |
| | | |
| | | |
| | | |
| | | |
| | | |
| | | |
| | | |
| | | |
| | | |
| | | |

2 Installation

2.1 Facility Requirements

After unpacking the modules of the UltiMate 3000 system, allow the units to warm up for approximately 4 hours before connecting them to the power supply. This delay allows any condensation that might have occurred during shipping to evaporate. After 4 hours, check the modules; if the condensation is still there, allow them to continue to warm up (without connecting them to the mains) until the condensation is completely gone.

Install the system in the laboratory on a stable surface that is free of vibration. Make sure that the surface is resistant to solvents. Avoid locations with extreme changes in temperature (such as direct sunlight or drafts) and high humidity. Allow sufficient clearance behind the system for power connections and ventilation.

Dionex recommends that you stack the individual modules of an UltiMate 3000 system, for example, as shown below. Note that the arrangement of the system modules depends on the application (→ Application Examples and Fluid Connections, page 13).



Fig. 1: Example for an UltiMate 3000 system

2.2 Rear Panel Connectors




Fig. 2: Rear panel view of an UltiMate 3000 system

Please note: Fig. 2 illustrates how the system modules are interconnected. This applies accordingly to for any other system structure.

2.2.1 Power Connection

Use the power cords from in the accessories kits to connect the modules to the mains. The instruments are equipped with a standard power supply unit with an automatic voltage selector. Thus, no adjustment is required to adapt the line voltage to local voltage requirements.

The solvent rack is usually connected to the mains via the pump $(\rightarrow \text{Fig. 2})$. However, an external power supply unit (part no. 6510.0004) is available as an option.

2.2.2 Interfaces for Device Control

Apart from the solvent rack, all modules of an UltiMate 3000 system can be connected to the server PC separately. However, Dionex recommends interconnecting the modules, and then connecting the system to the Chromeleon server PC via only one USB connection, e.g., from the autosampler (\rightarrow Fig. 2).

Please note: The detector is equipped with an RS-232 interface. Connect the detector to the autosampler of the UltiMate 3000 system (→ Fig. 2), using the RS-232/USB adaptor cable from the detector's accessories kit. For information about the USB installation of the adaptor cable, refer to the *Operating Instructions* for the detector.

Please note: The solvent rack is operated with Chromeleon via the pump.

For information about the connection ports on the system modules, as well as the connector pinout, refer the *Operating Instructions* for the respective instrument.

2.3 Opening the Front Panel

To provide easy access to the components in the instrument compartments, the front panels of the autosampler, detector, pump, and solvent rack tilt upward. The open cover locks in its topmost position. The front panel of the flow manager tilts downward.

⚠ **Important:** When lifting or moving an instrument, always lift by the bottom or sides. Lifting the instrument by the front panel may damage the front panel door.

⚠ Important: Lorsque vous soulevez ou déplacez un instrument, saisissez la toujours par le dessous ou les côtés l'unité. Soulever l'unité par le panneau avant risque d'endommager la porte du panneau avant.

2.4 Preparing the System

Important: The system modules are filled with 2-propanol when being shipped from the factory. During initial operation of the modules, make sure that the solvents used are miscible. Otherwise, use an appropriate intermediate solvent.

⚠ Important: Les instruments dans un système UltiMate 3000 sont stockés sous 2propanol. Lors du démarrage initial des instruments, assurez-vous que les solvants utilisés soient miscibles avec le 2-propanol. Sinon, suivez les étapes intermédiaires appropriées.

- Observe the information about the facility requirements when installing the system.
- Connect drainage tubing from the devices accessories kits to the outlet of the modules to direct liquids that might have gathered in the enclosure to the waste. The outlets are on the bottom right of the devices. (Please note: The detector has no outlet.)
 - Please note: You may connect the tubing as shown in Fig. 3. To do so, use the connecting pieces from the devices' accessories kits.




Fig. 3: Connecting drainage tubing

For the system structure shown in Fig. 3, Dionex recommends cutting the tubing to the following lengths:

| Description | Length (mm) |
|---|-------------|
| Tubing between solvent rack and pump | 215 |
| Tubing between pump and flow manager | 312 |
| Tubing between flow manager and autosampler | 332 |

Important: To prevent damage to the instrument, make sure that no part of the tubing is placed higher than the connection port.

⚠ Important: Afin d'éviter d'endommager l'instrument, assurez-vous qu'aucune pièce des tubes n'est placée plus haute que le port de raccordement.

Please note: In addition to the outlets described above, the flow manager and the solvent rack are equipped with another outlet on the rear panel. The WPS-3000T autosampler with temperature control has an additional outlet at the bottom right (→ Fig. 3). Connect the required tubing as described in the *Operating Instructions* for the devices.

- Prepare the individual modules as described in the *Operating Instructions* for the respective instrument.
- Connect the modules on the rear panel as described in Rear Panel Connectors (\rightarrow page 8).
- For examples of the fluid connections between the system modules, refer to Application Examples and Fluid Connections (→ page 13).
- Use the power cords shipped with the instruments to connect the units to the mains. Turn on the instruments by pressing the power switch on the rear panel.

| UltiMate 3000 - System Installation and Application | |
|---|--|
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |

3 Application Examples and Fluid Connections

Different system configurations allow you to use the UltiMate 3000 system for various applications. For information about the related system setup and the fluid connections in the system, refer to:

- Manual Injection (\rightarrow page 16)
- **Direct Injection** (\rightarrow page 20)
- Preconcentration (\rightarrow page 24)
- **2D** LC Salt Plugs (\rightarrow page 28)
- Comprehensive 2D LC (→ page 33)
- Parallel LC (\rightarrow page 39)

When installing the fluid connections in an UltiMate 3000 system, keep the following in mind:

⚠ Important: Different fitting systems are used in an UltiMate 3000 system. Therefore, install the capillaries and fittings only at the positions for which they are intended.

Important: Différents types de raccords sont utilisés dans le système UltiMate 3000. Par conséquent, installez les capillaires et les raccords uniquement dans aux endroits où ils sont prévus.

⚠ **Important:** When you operate biocompatible system modules, observe the instructions for connecting the capillaries and other components. Refer to the related sections of the Operating Instructions for the module(s).

Avant de accorder les capillaires où des autres éléments à un module biocompatible, référez-vous aux sections correspondantes dans le mode d'emploi pour le(s) module(s).

⚠ **Important:** Dionex recommends starting the fluid connections for the system modules at the pump. Connect the capillary to the pump first and rinse thoroughly before connecting the other end of the capillary to another module.

⚠ Important: Lorsque vous connectez les capillaires, respectez la séquence suivante: Raccordez le capillaire à la pompe. Rincez le capillaire avant de raccorder l'autre fin du capillaire à l'autre unité.

⚠ Important: When connecting the capillaries, make sure that the connectors are free from contaminants. Even minute particles may cause damage to the flow splitter, flow control valve, and column.

⚠ Important: Lorsque vous connectez les capillaires, assurez-vous que les raccords sont exempts de tout contaminant. Même d'infimes particules peuvent causer des dommages au système (ex. diviseur de débit, vanne de régulation de débit et colonne). **⚠** Important: Use only the capillaries shipped with the modules and original Dionex spare capillaries. Important: Utilisez uniquement les capillaires fournis avec les unités et les capillaires de rechange d'origine Dionex. ⚠ Important: To connect the capillaries to an injection valve or selector valve, install only the ferrules and fittings shipped with the capillaries. ✓!\ Important: Pour brancher les capillaires à une vanne d'injection ou de sélection, installez uniquement les ferrules et les raccords livrés avec la vanne.

⚠ Important: Reuse used fittings and ferrules only for the same capillary connection. This is to avoid increased dead volume.

Important: La réutilisation des raccords et ferrules n'est possible que pour la connexion capillaire d'origine, afin d'éviter l'apparition de volumes morts

Please note: The outlet block inlets of the pump are fitted with frit holders and inline filters. The pump is shipped with stainless steel frits (porosity: 0.5 µm). Check the permeability of the filter frits at regular intervals. When the pump delivers water at a flow rate of 2 ml/min and when the outlet is open, the pressure should not exceed 10 bar. If necessary, replace the filter frits (part no. 6000.0045 = 10 filter frits with SS $0.5 \mu m$) as described in the *Operating Instructions* for the pump.

Please note: Keep in mind that the flow splitter in the flow manager, the flow cell in the detector and the associated capillaries depend on the specific application and that it might be necessary to replace them as necessary.

The following flow splitters and flow cells are available:

| Description | Part No. |
|--|------------|
| Splitter cartridge (standard) for nano HPLC applications (split ratio 1:1000) | 6720.3150A |
| Splitter cartridge (standard) for capillary HPLC applications (split ratio 1:100) | 6720.3160A |
| Splitter cartridge (standard) for micro HPLC applications (split ratio 1:6) | 6720.3170A |
| Splitter cartridge (biocompatible) for nano HPLC applications (split ratio 1:1000) | 6721.3150A |
| Splitter cartridge (biocompatible) for capillary HPLC applications (split ratio 1:100) | 6721.3160A |
| Splitter cartridge (biocompatible) for micro HPLC applications (split ratio 1:6) | 6721.3170A |
| Flow cell for nano HPLC applications | 6073.0002 |
| Flow cell for capillary HPLC applications | 6073.0003 |
| Flow cell for micro HPLC applications | 6073.0004 |

For nano HPLC applications, the inner diameter of the capillaries is 20 $\mu m.$ For capillary HPLC applications, the inner diameter is 50 $\mu m.$

3.1 Manual Injection (Basic Configuration)

For this application, the sample is injected manually and directed to the separation column via the switching valve in the flow manager.

Dionex recommends that you stack the system modules as shown below:



Fig. 4: Module arrangement for manual injection

The system includes the modules listed in the table below. In addition, the appropriate Basic Configuration Kit (\rightarrow section 8.1, page 64) is required.

| System Module | | Par | t No. | |
|--|------------------------|------------------------|------------------------|------------------------|
| | Nano HPLC | Capillary HPLC | Micro HPLC | Monolith. HPLC |
| SRD-3300 | | 5035 | 5.9240 | |
| LPG-3300 | 503 | 5.0040 (standard) or | 5037.0040 (biocom | patible) |
| FLM-3300 Standard Biocompatible | 5720.0030 5021.0030 | 5720.0035 5021.0035 | 5720.0038 5021.0038 | 5720.0035 5021.0035 |
| Basic Configuration Kit Standard Biocompatible | 6720.0065 6721.0065 | 6720.0082 6721.0082 | 6720.0083 6721.0083 | 6720.0084 6721.0084 |

Fig. 5 and the related table provide an overview of the fluid connections for manual injection. For an example of the fluid connections for nano HPLC, also refer to Fig. 6 (\rightarrow page 19).




Fig. 5: Fluid connections for manual injection (overview)

| No. | Description | | Part | Part No.1 | |
|-----|--|-----------------------------------|--|---|-----------------------------------|
| | | Nano HPLC | Capillary HPLC | Micro HPLC | Monolith. HPLC |
| 1 | Solvent lines (analytical) | | 0030 | 6030.2548 | |
| 2 | Solvent lines from degasser to proportioning valve (set of 3 including the appropriate fittings) | | 6030 | 6030.2547 | |
| 3 | Connection proportioning valve to working pump head (including fittings) | | 9803 | 6035.2514 | |
| 4 | Capillary from working pump head to equilibration pump head (including fittings) | | 06030 | 6030.2515 | |
| 5 | Capillary from equilibration head to outlet block (including fittings) | | 9803 | 6035.3010 | |
| 9 | Inline filter (SS 0.5 µm) in frit holder on outlet block | | 0009 | 6000.0045 | |
| 7 | Capillary from pump to flow splitter (short connection) | 922 | 6035.2553 ² (standard) or 6037.2553 ² (biocompatible) | 037.2553 ² (biocomp | atible) |
| 8 | Capillary from switching valve to flow splitter (both in flow manager) Standard Biocompatible | $\frac{6720.0033^2}{6721.0033^2}$ | $6720.0034^2 \\ 6721.0034^2$ | $\frac{6720.0035^2}{6721.0035^2}$ | 6720.0034^2 6721.0034^2 |
| 6 | Capillary from port 6 to port 8 of the switching valve Standard Biocompatible | $\frac{6820.0015^2}{6821.0015^2}$ | 6820.0016^2 6821.0016^2 | $\frac{6820.0018^2}{6821.0018^2}$ | $\frac{6820.0015^2}{6821.0015^2}$ |
| 10 | Connection from switching valve (flow manager) to separation column | Connect the coluconnection | Connect the column to the switching valve. If necessary, use the appropriate connection parts from the respective Basic Configuration kit. | valve. If necessary, u ective Basic Configu | se the appropriate uration kit. |
| 111 | Separation column | 160321^2 | 160295^2 | 160282^{2} | 161409^2 |
| 12 | Separation column outlet | | All columns have | All columns have an inbuilt outlet. | |
| 13 | Capillary for manual injection | part no. 6 | Included in Manual Injection Port, part no. 6720.9007^2 (standard) or 6721.9007^2 (biocompatible) | Included in Manual Injection Port, 9007 ² (standard) or 6721.9007 ² (bioc | ompatible) |

¹ The part number refers to the packing unit. For more information, contact your Dionex Sales Representative. ² These parts are included in the Basic Configuration Kit for the related application (→ section 8.1, page 64).

The part numbers in Fig. 6 refer to nano HPLC (standard devices):




Fig. 6: Fluid connections for manual injection (nano HPLC, standard devices)

3.2 Direct Injection

For this application, the sample is drawn through the needle into the sample loop while the autosampler's injection valve is in **Load** position. The injection valve is switched to **Inject** position, directing the sample flow from the sample loop onto the separation column.

Dionex recommends that you stack the system modules as shown below:



Fig. 7: Module arrangement for direct injection

The system includes the modules listed in the table below. In addition, the appropriate Direction Injection Kit (\rightarrow section 8.2, page 68) is required.

| System Module | | Part | No. | |
|---|------------------------|--|------------------------|------------------------|
| | Nano HPLC | Capillary HPLC | Micro HPLC | Monolith. HPLC |
| SRD-3300 | | 5035.9 | 9240 | |
| LPG-3300 | 5035. | .0040 (standard) or 50 | 037.0040 (biocom | patible) |
| Option: UVD-3000 with flow cell | 5073.0010 | 5073.0020 | 5073.0030 | 5073.0010 |
| FLM-3300 Standard Biocompatible | 5720.0030 5721.0030 | 5720.0035 5721.0035 | 5720.0038 5721.0038 | 5720.0035 5721.0035 |
| WPS-3000 <i>or</i> WPS-3000 T | | .0010 (standard) or 58 .0020 (standard) or 58 | , | |
| Direct Injection Kit Standard Biocompatible | 6720.0042 6721.0042 | 6720.0043 6721.0043 | 6720.0046 6721.0046 | 6720.0044 6721.0044 |

Fig. 8 and the related table provide an overview of the fluid connections for direct injection. For an example of the fluid connections for nano HPLC, also refer to Fig. 9 (\rightarrow page 23).




Fig. 8: Fluid connections for direct injection (overview)

| No. | Description | | Part No. ¹ | No.1 | |
|-----|--|-----------------------------|---|--|--|
| | | Nano HPLC | Capillary HPLC | Micro HPLC | Monolith. HPLC |
| 1 | Solvent lines (analytical) | | 6030 | 6030.2548 | |
| 2 | Solvent lines from degasser to proportioning valve (set of 3 including the appropriate fittings) | | 6030.2547 | 2547 | |
| ε | Connection proportioning valve to working pump head (including fittings) | | 6035. | 6035.2514 | |
| 4 | Capillary from working pump head to equilibration pump head (including fittings) | | 0209 | 6030.2515 | |
| 5 | Capillary from equilibration head to outlet block (including fittings) | | 6035. | 6035.3010 | |
| 9 | Inline filter (SS 0.5 μm) in frit holder on outlet block | | 0009 | 6000.0045 | |
| 7 | Capillary from pump to flow splitter Long connection (system with detector) short connection (system without detector) | 6035. 6035. | 6035.2550^2 (standard) or 6037.2550^2 (biocompatible) 6035.2553^2 (standard) or 6037.2553^2 (biocompatible) | 037.2550 ² (biocomp 037.2553 ² (biocomp | atible) atible) |
| 8 | Capillary from separation column to detector flow cell | The f | The flow cell is shipped with appropriate capillaries. | ith appropriate capil | laries. |
| 6 | Separation column | 160321^2 | 160295^2 | 160282^2 | 161409^2 |
| 10 | Capillary from autosampler switching valve to separation column (PEEKsil) Standard Biocompatible | 6720.0024^2 6721.0024^2 | 6720.0025^2 6721.0025^2 | $6720.0026^2 \\ 6721.0026^2$ | 6720.0025^2 6721.0025^2 |
| 11 | Capillary from autosampler switching valve to flow splitter (PEEKsil) Standard Biocompatible | 6720.0027^2 6721.0027^2 | 6720.0028^2 6721.0028^2 | $6720.0029^2 \\ 6721.0029^2$ | 6720.0028^2 6721.0028^2 |
| 12 | Standard needle (2.4 µl, fused silica) | 6850 | 6820.3010 (standard) or 6821.3010 (biocompatible) | 821.3010 (biocompa | ıtible) |
| 13 | Sample loop Standard Biocompatible | 6820 6821 | 6820.0015 ² 6821.0015 ² | $6820.0018^{2} \\ 6821.0018^{2}$ | 6820.0015 ² 6821.0015 ² |
| | | | | | |

¹ The part number refers to the packing unit. For more information, contact your Dionex Sales Representative. ² These parts are included in the Direct Injection Kit for the related application $(\rightarrow$ section 8.2, page 68).

The part numbers in Fig. 9 refer to nano HPLC (standard devices):




Fig. 9: Fluid connections for direct injection (nano HPLC, standard devices)

3.3 Preconcentration

For this application, the sample is drawn through the needle into the sample loop while the autosampler's injection valve is in **Load** position. The injection valve is switched to the **Inject** position, allowing the loading pump to direct the sample flow from the sample loop onto the precolumn for sample enrichment. The switching valve in the flow manager is switched from **10-1** to **1-2**. The micro pump and the flow splitter direct the sample onto the separation column. Dionex recommends that you stack the system modules as shown below:



Fig. 10: Module arrangement for preconcentration

The system includes the modules listed in the table below. In addition, the appropriate Preconcentration Kit (\rightarrow section 8.3, page 72) is required:

| System Module | | Part | No. | |
|---|------------------------|--|------------------------|------------------------|
| | Nano HPLC | Capillary HPLC | Micro HPLC | Monolith. HPLC |
| SRD-3600 | | 5035.9 | 9230 | |
| LPG-3600 | 5035. | 0035 (standard) or 50 | 037.0035 (biocom | patible) |
| Option: UVD-3000 with flow cell | 5073.0010 | 5073.0020 | 5073.0030 | 5073.0010 |
| FLM-3300 Standard Biocompatible | 5720.0030 5721.0030 | 5720.0035 5721.0035 | 5720.0038 5721.0038 | 5720.0035 5721.0035 |
| WPS-3000 <i>or</i> WPS-3000 T | | 0010 (standard) or 58 0020 (standard) or 58 | | |
| Preconcentration Kit Standard Biocompatible | 6720.0047 6721.0047 | 6720.0048 6721.0048 | 6720.0050 6721.0050 | 6720.0049 6721.0049 |

Fig. 11 and the related table provide an overview of the fluid connections for preconcentration. For an example of the fluid connections for nano HPLC, also refer to Fig. 12 (\rightarrow page 27).




Fig. 11: Fluid connections for preconcentration (overview)

| No. | Description | | Part No.1 | |
|-----|---|--------------------------------|---|---------------------------------|
| | | Nano HPLC | Capillary HPLC | Monolith. HPLC |
| 1 | Solvent lines (analytical) | | 6030.2548 | |
| 2 | Solvent lines from degasser to proportioning valve (set of 3 including the appropriate fittings) | | 6030.2547 | |
| 3 | Connection proportioning valve to working pump head (including fittings) | | 6035.2514 | |
| 4 | Capillary from working pump head to equilibration pump head (including fittings) | | 6030.2515 | |
| 5 | Capillary from equilibration head to outlet block (including fittings) | | 6035.3010 | |
| 9 | Inline filter (SS 0.5 µm) in frit holder on outlet block | | 6000.0045 | |
| 7 | Capillary from pump to flow splitter Long connection (system with detector) Short connection (system without detector) | 6035.2550 ² (| 6035.2550^2 (standard) or 6037.2550^2 (biocompatible) 6035.2553^2 (standard) or 6037.2553^2 (biocompatible) | ocompatible) ocompatible) |
| ∞ | Capillary from flow splitter to flow manager switching valve Standard Biocompatible | 6720.0033^2 6721.0033^2 | $\frac{6720.0034^2}{6721.0034^2}$ | $6720.0034^2 \\ 6721.0034^2$ |
| 6 | Precolumn holder (with the appropriate capillaries) and Precolumn | $\frac{6720.0012^2}{160454^2}$ | $6720.0012^2 \\ 160454^2$ | $\frac{163972^{2}}{163972^{2}}$ |
| 10 | Separation column (The column is connected directly to the switching valve.) | 160321^2 | 160295^2 | 161409^2 |
| 11 | Capillary from the separation column to the detector flow cell | The flow ce | The flow cell is shipped with appropriate capillaries. | e capillaries. |
| 12 | Capillary from loading pump outlet to autosampler Long connection (system with detector) Short connection (system without detector) | 6720.0032 ² (| 6720.0032^2 (standard) or 6721.0032^2 (biocompatible) 6720.0031 (standard) or 6721.0031 (biocompatible) | ocompatible) compatible) |
| 13 | Capillary from autosampler to flow manager switching valve | 6720.0038 ² (| 6720.0038^2 (standard) or 6721.0038^2 (biocompatible) | ocompatible) |
| 14 | Sample loop | 6820.0018 ² (| 6820.0018^2 (standard) or 6821.0018^2 (biocompatible) | ocompatible) |
| 15 | Standard needle (2.4 µl, fused silica) | 6820.3010 | 6820.3010 (standard) or 6821.3010 (biocompatible) | compatible) |

¹ The part number refers to the packing unit. For more information, contact your Dionex Sales Representative. ² These parts are included in the Preconcentration Kit for the related application (\rightarrow section 8.3, page 72).

The part numbers in Fig. 12 refer to nano HPLC (standard devices):




Fig. 12: Fluid connections for preconcentration (nano HPLC, standard devices)

3.4 2D LC Salt Plug

This application supports online two-dimensional separation of peptides: While the autosampler's injection valve is in **Load** position, the sample is drawn through the needle into the sample loop. The injection valve is switched to the **Inject** position, allowing the loading pump to direct the sample onto the ion exchange column (SCX). The right valve in the flow manager is in position **10-1**. The part of the sample that is not retained is trapped on the reverse phase pre-column; the left valve in the flow manager is in position **10-1**.

When loading is finished, the ion exchange column is removed from the fluidic path by switching the right valve to position 1-2. This allows desalting the sample collected on the trap column with the loading pump.

After the desalting step, the left valve is switched to position 1-2 to run the analytical gradient through the trap column and the reverse phase separation column.

When the first fraction is separated, both valves are switched to position 10-1 and a salt solution is injected. A new fraction of peptides elutes from the ion exchange column to the trap column. The desalting and analytical steps described above are carried out.

The process is iterated with increasing concentrations of salt injected until all the peptides from the ion exchange column are eluted and separated on the reverse phase column.

Dionex recommends that you stack the system modules as shown below:



Fig. 13: Module arrangement for 2D LC salt plug

The system includes the modules listed in the table below. In addition, the 2D LC Salt Plug Kit (\rightarrow section 8.4, page 77) is required:

| System Module | Part No. |
|--|---|
| | Nano HPLC |
| SRD-3600 | 5035.9230 |
| LPG-3600 | 5035.0035 (standard) or 5037.0035 (biocompatible) |
| Option: UVD-3000 with flow cell | 5073.0010 |
| FLM-3100 | 5720.0010 (standard) or 5721.0010 (biocompatible) |
| WPS-3000 or WPS-3000 T | 5820.0010 (standard) or 5821.0010 (biocompatible) 5820.0020 (standard) or 5821.0020 (biocompatible) |
| 2D LC Salt Plug Kit Standard Biocompatible | 6720.0051 6721.0051 |

Fig. 14 and the related table provide an overview of the fluid connections for 2D LC salt plug. Also, refer to Fig. 15 (\rightarrow page 32).




Fig. 14: Fluid connections for 2D LC salt plug (overview)

| No. | Description | Part No. ¹ |
|-----|---|---|
| | | Nano HPLC |
| 1 | Solvent lines (analytical) | 6030.2548 |
| 2 | Solvent lines from degasser to proportioning valve (set of 3 including the appropriate fittings) | 6030.2547 |
| 8 | Connection proportioning valve to working pump head (including fittings) | 6035.2514 |
| 7 | Capillary from working pump head to equilibration pump head (including fittings) | 6030.2515 |
| 5 | Capillary from equilibration head to outlet block (including fittings) | 6035.3010 |
| 9 | Inline filter (SS 0.5 µm) in frit holder on outlet block | 6000.0045 |
| L | Capillary from pump to flow splitter Long connection (system with detector) Short connection (system without detector) | 6035.2550 ² (standard) or 6037.2550 ² (biocompatible) 6035.2553 ² (standard) or 6037.2553 ² (biocompatible) |
| 8 | Capillary from flow splitter to left switching valve in flow manager | 6720.0033 ² (standard) or 6021.0033 (biocompatible) |
| 6 | Precolumn holder (with appropriate capillaries) and Trap column | $6720.0012^2 \\ 160454^2$ |
| 10 | Separation column (The column is connected directly to the switching valve.) | 160321^2 |
| 11 | Capillary from separation column to detector flow cell | The flow cell is shipped with appropriate capillaries. |
| 12 | PEEKsil bridge from left to right switching valve in the flow manager | 6720.0060^2 (standard) or 6721.0060^2 (biokompatibel) |
| 13 | Ion exchange column (SCX) | 162152 ² |
| 14 | Capillary from loading pump outlet block to autosampler Long connection (system with detector) Short connection (system without detector) | 6720.0032 ² (standard) or 6721.0032 ² (biocompatible) 6720.0031 (standard) or 6721.0031 (biocompatible) |
| 15 | Sample loop (20 µl, PEEKsil) | 6820.0018 ² (standard) or 6821.0018 ² (biocompatible) |
| 16 | Capillary from autosampler to right switching valve in flow manager | 6720.0038 ² (standard) or 6721.0038 ² (biocompatible) |
| 17 | Standard needle (2.4 µl, fused silica) | 6820.3010 (standard) or 6821.3010 (biocompatible) |

¹ The part number refers to the packing unit. For more information, contact your Dionex Sales Representative. ² These parts are included in the 2D LC Salt Plug Kit (\rightarrow section 8.4, page 77).

Fig. 15, too, illustrates the fluid connections for 2D LC salt plug. The part numbers in the picture refer to the standard devices.




Fig. 15: Fluid connections for 2D LC salt plug (nano HPLC, standard devices)

3.5 Comprehensive 2D LC

This application supports online comprehensive two-dimensional separation of peptides: While the autosampler's injection valve is in **Load** position, the sample is drawn through the needle into the sample loop. The injection valve is switched to the **Inject** position, allowing the capillary pump to direct the sample to the ion exchange column (SCX).

The peptides that are not retained are collected on trap column 1 (\rightarrow Fig. 17, no. 12) that is inline with the ion exchange column. Valve 2 is in position **1-2**. In parallel, valve 1 is in position **1-2** to flush the trap column 2 (\rightarrow Fig. 17, no. 15) via the loading pump.

At the end of the loading step, valve 2 is switched to position **10-1**, setting trap column 2 in line with the ion exchange column. Trap column 1 is now in line with the loading pump for a desalting step. When desalting is finished, valve 1 is switched to position **10-1**. This connects trap column 1 and the reverse phase column, thus allowing the separation of the peptides collected on trap column 1. During that time, a salt gradient is started. The peptides eluting from the ion exchange column are trapped on trap column 2, collecting a new peptides fraction.

At the end of the separation on the reverse phase column, valve 1 is switched to position **1-2** to condition trap column 1 via the loading pump (acetonitrile is removed). When trap column 1 is ready again, valve 2 is switched to position **1-2**. Trap column 1 is flushed with the salt buffer while trap column 2 is equilibrated with the loading solvent.

At this stage, a new salt gradient starting from the highest salt percentage reached in the previous one can be started to load peptides on trap column 1. Simultaneously, the peptides previously collected on the trap column 2 are being separated, with valve 1 in position 10-1.

The process can be iterated, running an increasing salt concentration on the ion exchange column. In this way, peptides are collected on one trap column while a reverse phase separation is performed for the peptides collected on the other trap column during the previous iteration.

- Please note: The right choice of the time allowed for desalting and reconditioning the trap columns is crucial for this application.
- Please note: The application comprises 2 pumps (LPG-3300 and LPG-3600) and two FLM-3300 flow managers (→ Fig. 16, page 34). Therefore, adapt the names under which the devices are identified in the installation environment and in the Chromeleon client in the Chromeleon Server Configuration program. Keep in mind that you may have to adapt the links to the associated control panels as well.

SRD-3600 UltiMate 3000 Solvent Rack LPG-3600 Pump UltiMate 3000 SRD-3300 Solvent Rack UVD-3000 UV Detector LPG-3300 Pump -FLM-3300 Flow Manager WPS-3000 or WPS-3000 T Autosampler FLM-3300

Dionex recommends that you arrange the system modules as shown below:

Fig. 16: Module arrangement for comprehensive 2D LC

Flow Manager

The system includes the modules listed in the table below. In addition, the Comprehensive 2D LC Kit (\rightarrow section 8.5, page 79) is required:

| System Module | Part No. |
|---------------------------------|--|
| | Capillary/Nano HPLC |
| SRD-3300 | 5035.9240 |
| LPG-3300 | 5035.0040 (standard) or 5037.0040 (biocompatible) |
| WPS-3000 or WPS-3000 T | 5820.0010 (standard) or 5821.0010 (biocompatible) 5820.0020 (standard) or 5821.0020 (biocompatible) |
| SRD-3600 | 5035.9230 |
| LPG-3600 | 5035.0035 (standard) or 5037.0035 (biocompatible) |
| Option: UVD-3000 with flow cell | 5073.0010 |
| FLM-3300 (Nano) | 5720.0030 (standard) or 5721.0030 (biocompatible) |
| FLM-3300 (Capillary) | 5720.0035 (standard) or 5721.0035 (biocompatible) |
| Comprehensive 2D LC Kit | 6720.0053 (standard) or 6721.0053 (biocompatible) |

Fig. 17 and the related table provide an overview of the fluid connections for 2D LC comprehensive for capillary (1st dimension) nano (2nd dimension) HPLC. Also, refer to Fig. 18 (\rightarrow page 38).




Fig. 17: Fluid connections for comprehensive 2D LC (overview)

| No. | Description | Part No. ¹ |
|-----|---|---|
| | | Capillary/Nano HPLC |
| 1 | Solvent lines (analytical) | 6030.2548 |
| 2 | Solvent lines from degasser to proportioning valve (set of 3 including the appropriate fittings) | 6030.2547 |
| 3 | Connection proportioning valve to working pump head (including fittings) | 6035.2514 |
| 4 | Capillary from working pump head to equilibration pump head (including fittings) | 6030.2515 |
| 5 | Capillary from equilibration head to outlet block (including fittings) | 6035.3010 |
| 9 | Inline filter (SS 0.5 μm) in frit holder on outlet block | 6000.0045 |
| 7 | Capillary from capillary pump (LPG-3600) to flow splitter Long connection (system with detector) Short connection (system without detector) | 6035.2556 ² (standard) or 6037.2556 (biocompatible) 6035.2554 ² (standard) or 6037.2554 (biocompatible) |
| 8 | Capillary from autosampler to flow splitter (cap or micro) | 6720.0028 ² (standard) or 6721.0028 (biocompatible) |
| 6 | Capillary from autosampler to ion exchange column (SCX) | 6720.0058 ² (standard) or 6721.0058 ² (biocompatible) |
| 10 | Ion exchange column (SCX) | 162122^{2} |
| 11 | Capillary from ion exchange column (SCX) to switching valve | The ion exchange column has an inbuilt outlet. |
| 12 | Precolumn holder (with appropriate capillaries) and Trap column | $6720.0012^2 \\ 160454^2$ |
| 13 | Capillary from pump loading pump (LPG-3300) to switching valve (nano FLM) | 6720.0057^2 (standard) or 6721.0057^2 (biocompatible) |
| 14 | PEEKsil bridges from switching valve (nano FLM) to switching valve (cap/micro FLM) | $(2x) 6720.0062^2$ (standard) or 6721.0062^2 (biocompatible) |
| 15 | Precolumn holder (with appropriate capillaries) and Trap column | $\frac{6720.0012^2}{160454^2}$ |
| 16 | PEEKsil bridge: Nano flow manager: from port 4 to port 7 Cap/micro flow manager: from port 3 to port 8 | 6720.0061^{2} (standard) or 6721.0061^{2} (biocompatible) 6720.0061^{2} (standard) or 6721.0061^{2} (biocompatible) |
| 17 | Separation column (The column is connected directly to the switching valve.) | 160321^2 |
| 18 | Capillary from separation column to detector flow cell | The flow cell is shipped with appropriate capillaries. |

| | No. Description | Part No. |
|----------------------|---|---|
| | | Capillary/Nano HPLC |
| 19 Ca _l | Capillary from analytical pump (LPG-3600) to flow splitter | |
| Lo: Sh | Long connection (system with detector) Short connection (system without detector) | 6035.2536 (standard) or 6037.2536 (biocompatible) 6035.2554 ² (standard) or 6037.2554 (hiocompatible) |
| | | (conditions) |
| 20 Ca _l | Capillary from nano flow splitter to switching valve | 6720.0033 ² (standard) or 6021.0033 (biocompatible) |
| 21 Sar | 21 Sample loop (20 µl, PEEKsil) | 6820.0018^2 (standard) or 6821.0018^2 (biocompatible) |
| 22 Sta | Standard needle (2.4 µl, fused silica) | 6820.3010 (standard) or 6821.3010 (biocompatible) |

¹ The part number refers to the packing unit. For more information, contact your Dionex Sales Representative.
² These parts are included in the Comprehensive 2D LC Kit (\rightarrow section 8.5, page 79).

Fig. 18, too, illustrates the fluid connections for comprehensive 2D LC for capillary (1st dimension) and nano (2nd dimension) HLPC. The part numbers in the picture refer to the standard devices.




Fig. 18: Fluid connections for comprehensive 2D LC (standard devices)

3.6 Parallel LC

This application allows you to increase the throughput for peptide separations. While the autosampler's injection valve is in **Load** position, the sample is drawn through the needle into the sample loop. The injection valve is switched to the **Inject** position, allowing the loading pump to direct the sample from the sample loop onto trap column 1 (in the FLM-3200; \rightarrow Fig. 20, no. 9) for sample preconcentration.

Valve 1 (left switching valve in the FLM-3200) is in position **10-1**. In parallel, valve 2 (in the FLM-3300) is in position **1-2** to allow the wash step to elute remaining hydrophobic compounds from trap column 2 (in the FLM-3300; \rightarrow Fig. 20, no. 11) and separation column 2 (in the FLM-3300; \rightarrow Fig. 20, no. 19). The nano valve (right switching valve in the FLM-3200) is in position **6-1** to direct the eluent from the wash-step to the detector.

When the sample is loaded on trap column 1, valve 1 is switched to position **1-2**. The analytical pump 1 then directs the sample from trap column 1 to separation column 1 (in the FLM-3200; \rightarrow Fig. 20, no. 14). At the time when the first peaks are expected, the nano valve is switched to position **1-2** to direct the eluent to the detector. In the meantime, separation column 2 is reconditioned. Shortly before the separation ends, valve 2 is switched to position **10-1**. Trap column 2 is conditioned with the loading solvent.

As soon as the gradient over separation column 1 is finished, a new sample is injected. Valve 1 is in position **1-2.** Valve 2 is in position **10-1** to load the sample on trap column 2. At the same time, separation column 1 is washed and its eluent is collected, with the nano valve being in position **1-2**.

When the trap column 2 is loaded, valve 2 is switched to position **1-2**. This connects the separation column 2 with the trap column 2 (analytical pump). At the end of the wash step for separation column 1, the nano valve is switched to position **6-1**. The separation column 1 is reconditioned while the separation runs on separation column 2. Before the experiment ends, valve 1 is switched to position **10-1** to prepare for the next injection.

This sequence, which can be repeated as often as needed, allows running the analysis on one column and simultaneously washing and equilibrating the other column.

Please note: The application comprises two pumps (LPG-3300 and LPG-3600) and two flow managers (FLM-3200 and FLM-3300; → Fig. 19, page 40). Therefore, adapt the names under which the devices are identified in the installation environment and in the Chromeleon client in the Chromeleon Server Configuration program. Keep in mind that you may have to adapt the links to the associated control panels as well.

Dionex recommends that you arrange the system modules as shown below:




Fig. 19: Module arrangement for parallel LC

The system includes the modules listed in the table below. In addition, the appropriate Parallel LC Kit (\rightarrow section 8.6, page 81) is required:

| System Module | | Part No. | | |
|--|--|------------------------------|-------------|--|
| | Nano/Nano Capillary/Capillary Monolith./Monolith. HPLC HPLC HPLC | | | |
| SRD-3300 | 5035.9240 | | | |
| LPG-3300 | 5035.0040 (standard) or 5037.0040 (biocompatible) | | | |
| WPS-3000 <i>or</i> WPS-3000 T | 5820.0010 (standard) or 5821.0010 (biocompatible) 5820.0020 (standard) or 5821.0020 (biocompatible) | | | |
| SRD-3600 | | 5035.9230 | | |
| LPG-3600 | 5035.0035 (s | standard) or 5037.0035 (biod | compatible) | |
| Option: UVD-3000 with flow cell | 5073.0010 5073.0020 5073.0010 | | | |
| FLM-3200 Standard Biocompatible | 5720.0020 5720.0025 5720.0025 5721.0020 5721.0025 5721.0025 | | | |
| FLM-3300 Standard Biocompatible | 5720.0030 5720.0035 5720.0035 5721.0030 5721.0035 5721.0035 | | | |
| Parallel LC Kit Standard Biocompatible | 6720.0054 6720.0055 6721.0056 6721.0054 6721.0055 6721.0056 | | | |

Fig. 20 and the related table provide an overview of the fluid connections for parallel LC. For an example of the fluid connections for nano HPLC, also refer to Fig. 21 (\rightarrow page 44).




Fig. 20: Fluid connections for parallel LC (overview)

| • | Description | | Part No.1 | |
|----------------|---|--|---|--|
| | | Nano/Nano HPLC | Capillary/Capillary HPLC | Monolith./Monolith. HPLC |
| 1 | Solvent lines (analytical) | | 6030.2548 | |
| 2 S | Solvent lines from degasser to proportioning valve (set of 3 including the appropriate fittings) | | 6030.2547 | |
| 3 C | Connection proportioning valve to working pump head (including fittings) | | 6035.2514 | |
| 4 C | Capillary from working pump head to equilibration pump head (including fittings) | | 6030.2515 | |
| 5 C | Capillary from equilibration head to outlet block (including fittings) | | 6035.3010 | |
| П 9 | Inline filter (SS 0.5 µm) in frit holder on outlet block | | 6000.0045 | |
| 7 C | Capillary from loading pump (LPG-3300) to autosampler | 6720.0057 ² (s | 6720.0057 ² (standard) or 6721.0057 ² (biocompatible) | ocompatible) |
| 8 | Capillary from autosampler to FLM-3200 left switching valve | 6720.0038 ² (s | 6720.0038 ² (standard) or 6721.0038 ² (biocompatible) | ocompatible) |
| Б Б | Precolumn holder (with appropriate capillaries) and Trap column | $6720.0012^2 \\ 160454^2$ | $\frac{6720.0012^2}{160454^2}$ | 164972 ² |
| 10 P | PEEKsil bridge from left switching valve (FLM-3200) to FLM-3300 switching valve | 6720.0059^{2} (| 6720.0059^2 (standard) or 6721.0059^2 (biocompatible) | ocompatible) |
| 11 P | Precolumn holder (with appropriate capillaries) and Trap column | $6720.0012^2 \\ 160454^2$ | $\frac{6720.0012^2}{160454^2}$ | 164972 ² |
| 12 C | Capillary from analytical pump 1 (LPG-3600-left pump) to FLM-3200 flow splitter Long connection (system with detector) Short connection (system without detector) | 6035.2556 ² (6035.2554 ² (| 6035.2556 ² (standard) or 6037.2556 ² (biocompatible) 6035.2554 ² (standard) or 6037.2554 ² (biocompatible) | ocompatible) ocompatible) |
| 13 F S B | FLM-3200: Capillary from flow splitter to left switching valve Standard Biocompatible | 6720.0033^{2} 6721.0033^{2} | 6720.0034^2 6721.0034^2 | $6720.0034^2 \\ 6721.0034^2$ |
| 14 S | Separation column (The column is connected directly to the switching valve.) | 160321^2 | 160295^2 | 16409^{2} |
| 15 F | FLM-3200: Capillary from the separation column to the right switching valve | Connect the column to t connection p | Connect the column to the switching valve. If necessary, use the appropriate connection parts from the respective Parallel LC kit. | sary, use the appropriate rallel LC kit. |
| 16 C | Capillary from right switching valve (FLM-3200) to detector flow cell | The flow ce | The flow cell is shipped with appropriate capillaries. | e capillaries. |

| No. | No. Description | | Part No. 1 | |
|-----|---|--|---|------------------------------|
| | | Nano/Nano HPLC | Capillary/Capillary HPLC | Monolith./Monolith. HPLC |
| 17 | Capillary from analytical pump 2 (LPG-3600, right pump) to FLM-3300 flow splitter Long connection (system with detector) Short connection (system without detector) | 6035.2556 ² 6035.2554 ² | 6035.2556 ² (standard) or 6037.2556 (biocompatible) 6035.2554 ² (standard) or 6037.2554 (biocompatible) | ocompatible) |
| 18 | FLM-3300: Capillary from flow splitter to switching valve Standard Biocompatible | 6720.0033^2 6721.0033^2 | $6720.0034^2 \\ 6721.0034^2$ | $6720.0034^2 \\ 6721.0034^2$ |
| 19 | Separation column (The column is connected directly to the switching valve.) | 160321^2 | 160295^2 | 16409^2 |
| 20 | Connection of separation column (FLM-3300) to right switching valve (FLM-3200) | The separation colu | The separation columns have an inbuilt fused silica outlet (for direct connection to the switching valve). | lica outlet (for direct ve). |
| 21 | Sample loop (20 µl, PEEKsil) | 6820.0018^{2} | 6820.0018^2 (standard) or 6821.0018^2 (biocompatible) | ocompatible) |
| 22 | Standard needle (2.4 µl, fused silica) | 6820.3010 (| 6820.3010 (standard) or 6821.3010 (biocompatible) | compatible) |

¹ The part number refers to the packing unit. For more information, contact your Dionex Sales Representative. ² These parts are included in the Parallel LC Kit for the related application (\rightarrow section 8.6, page 81).

The part numbers in Fig. 21 refer to nano HPLC (standard devices):




Fig. 21: Fluid connections for parallel LC (nano HPLC; standard devices)

4 Automated Control by Chromeleon

1 Please note: The descriptions below refer to Chromeleon 6.70.

4.1 General

The modules of the UltiMate 3000 system can be controlled by the Chromeleon Chromatography Management System. To control the modules, a Chromeleon version 6.60 Service Pack 3 or higher and a Timebase Class 1 Chromeleon license are required. (If you have any questions, please contact your Dionex sales representative.)

Dionex recommends that you follow the installation steps below:

- Install the Chromeleon software before connecting the system modules to the USB port on the Chromeleon server PC. This will ensure that most of the necessary USB drivers are copied to appropriate locations (except the drivers for the RS-232/USB interface cable, see below).
- Connect the solvent rack to the 15-pin D-sub port on the pump.
- Connect the modules that are controlled via a USB connection (LPG-3600, FLM-3100, and WPS-3000) to each other. Use the USB cables that are shipped with the modules. Dionex recommends connecting all modules to one "central module" (e.g. the autosampler) and connecting this module to the CM Server PC (→ Fig. 2, page 8).
- Power up the "central module." The Plug&Play Manager (Hardware Wizard) guides you
 through the remaining USB installation for this module. Continue with powering up and
 installing the other USB modules. For more information about the USB installation steps,
 refer to the Operating Instructions for the modules.
- Connect the UVD-3000's RS-232 port to the USB hub on the "central module," using the RS-232/USB interface cable that is shipped with the detector. The Plug&Play Manager (Hardware Wizard) guides you through the remaining USB installation for the interface cable. The driver is available on the Chromeleon software CD in the Drivers\W&T RS232-USB Interface Cable directory. During installation, a new COM port is created in the Windows Device Manager. This port is required when you configure the UVD-3000 in Chromeleon. For more information, refer to the *Operating Instructions* for the detector. An alternative is to connect the detector directly to a COM port on the Chromeleon Server PC, using a 1:1 modem cable (Dionex part no.: 8914.0143).
- Install and configure the system modules in the Chromeleon Server Configuration program as described in the *Operating Instructions* for the module. You can add the modules in any order to the timebase. However, some settings affect more than one module. Therefore, the preferred order of installation is as follows (→ Fig. 22, page 46):
 - Flow Manager (FLM-3x00)
 - Pump (LPG-3x00): Assign a flow splitter in the flow manager to a pump as described in section 4.2 (→ page 47)
 - Autosampler (WPS-3000): Synchronize the injection with the pump as described in section $4.4 (\rightarrow page 51)$

- Detector (UVD-3000): Use the COM port created during installation of the RS-232/ USB interface cable.
- Solvent Rack (SRD-3x00): No configuration needed

If you install the modules in a different order, reopen the pump's and sampler's **Properties** pages and make the necessary configuration changes.




Fig. 22: Installation of an UltiMate 3000 system in the Server Configuration program (here: for Direct Injection, Preconcentration, or 2D LC Salt Plug)

Please note: The *Operating Instructions* for the modules also provide an overview of the different commands and properties supported by Chromeleon for the respective instrument. For more information, refer to the *Chromeleon online Help* and the *Chromeleon manuals*.

Please note: Operational and/or Performance Qualification allows you to check and document the quality of your HPLC system. All required materials and detailed instructions are available from Dionex on request.

4.2 Assigning a Flow Splitter to a Pump

If your system includes an FLM-3000 series flow manager, you can use the **Devices** tab page of the pump's properties to specify whether a flow splitter is connected to the pump.




Fig. 23: Devices page (here: for an LPG-3600)

Please note: Flow splitter assignment is possible only if both the pump and the flow manager are installed in Chromeleon.

The default entry is <None>. To change the setting, click the arrow next to the uses FLM-3x00 Flow Splitter input field and select the splitter from the drop-down.

- Please note: When you assign a flow splitter to a pump, the pump flow settings apply directly to the flow on the splitter outlet, i.e., the split ratio is considered automatically.
- Please note: For the standard configuration of the UltiMate 3000 system with an LPG-3600 pump, Dionex recommends assigning the flow splitter to the right pump (MicroPump).

Flow splitter assignment also affects the properties that are available for the pump in the **Commands** dialog box:

| Chromeleon | Description |
|-----------------|--|
| Flow and/or | Flow rate in μl/min |
| MasterFlow | Pumps without flow splitter assignment: |
| | Flow indicates the total flow through the pump. |
| | MasterFlow is not available. |
| | Pumps with flow splitter assignment: |
| | Flow indicates the flow through the column. |
| | MasterFlow indicates the flow through the master pump, i.e., the total |
| | flow which is the flow before the flow splitter. |
| Pressure and/or | Indicates the pressure |
| MasterPressure | Pumps without flow splitter assignment: |
| | Pressure indicates the pump pressure. |
| | MasterPressure is not available. |
| | Pumps with flow splitter assignment: |
| | Pressure indicates the current column pressure. |
| | MasterPressure indicates the pressure of the master pump, i.e., the |
| | pressure before the flow splitter. |

For more information, refer to the *Operating Instructions* for the pump and/or flow manager or to the Chromeleon *online Help*.

4.3 Recording the Pump Pressure as a Separate Channel

If your system includes a flow manager, the **Pressure Signal** check box is selected by default on the **Configuration** tab page in the flow manager's properties.

Accept this setting if you want to record the column pressure as a separate channel.




Fig. 24: Configuration tab page

Please note: When your UltiMate 3000 system includes one flow manager: Accept the default name (ColumnPressure), under which the column pressure signal is identified in the installation environment.

When your system includes more than one flow manger: Adapt the name for all additional flow managers according to your application. Keep in mind that you have to re-link the corresponding controls on the related control panel(s).

If the column pressure is recorded as a separate channel, **ColumnPressure** appears in the **Commands** dialog box, also.

| Chromeleon | Description |
|----------------|---|
| ColumnPressure | Click the "+"sign beside the name to display the items underneath: |
| | Delta —reports the signal's slope, i.e., the difference between the current value and the value one second ago. This is useful for triggers. |
| | Signal —has the following entries: Value—reports the current signal value (read-only), UpperLimit—sets the upper signal limit, and LowerLimit—sets the lower signal limit. |
| | AcqOn—starts data acquisition. AcqOff—terminates data acquisition. |
| | Retention —reports the retention time of the signal (read-only). |
| | MaxAutoStep —sets the maximum step rate for Auto Step Mode (range: 0.15.1 s; default: 5.1 s). |
| | Step —sets the step for data acquisition (range: 0.014.80 s; Auto selects the best step dynamically.) |
| | Average—averages all measured values over the step interval. (The default setting is: On. Off records only the last point of each interval.) |

For more information, refer to the *Operating Instructions* for the flow manager and/or to the Chromeleon *online* Help.

4.4 Synchronizing the Autosampler with a Pump

If your system includes a WPS-3000 or WPS-3000T autosampler, you can synchronize the injection command of the autosampler with the strokes of an LPG-3600 or LPG-3300 pump. Synchronization ensures that all injections are performed at the same phase of the pump cycle.

On the **Segments** tab page of the autosampler's properties, determine whether synchronization shall be performed.

Please note: Dionex recommends always enabling synchronization as synchronization helps enhancing the analytical results considerably.

Click the arrow of the **Sync Inject with pump** field and select an option from the drop-down list:




Fig. 25: Segments tab page (here: when the UltiMate 3000 system includes an LPG-3600 pump)

- Please note: The options are available only if both the autosampler and the pump are installed in Chromeleon.
- Select **None** to disable synchronization.
- When the UltiMate 3000 system includes an LPG-3600 pump: Select **LPG3x00_L_STRK** to enable synchronization with the left pump. To enable synchronization with the right pump, select **LPG3x00_R_STRK**.
- When the UltiMate 3000 system includes an LPG-3300 pump: Select **LPG3x00 STRK** to enable synchronization.
- Please note: You can change the default synchronization assignment for a specific application, using the PumpDevice and SyncWithPump properties in the Commands dialog box under Sampler on the control panel or in the program file (PGM).

To disable synchronization for a specific application, set **SyncWithPump** to **Off**. Use the **PumpDevice** property to select the pump for which synchronization shall be performed. (The drop-down list contains names of the pumps installed in the Server Configuration program.)

Changing the settings for these properties from the **Commands** dialog box does not change the standard synchronization setting on the **Segments** tab page.

For more information, refer to the *Operating Instructions* for the autosampler and/or pump or to the Chromeleon *online Help*.

4.5 Displaying the Degasser-Related Parameters

If your system includes an SRD-3000 series solvent rack, the degasser can be controlled with Chromeleon via the pump. In this case, the following degasser-related parameters are listed in the **Commands** dialog box under the pump to which the solvent rack is connected:

| Chromeleon | Description |
|-----------------|---|
| Degasser | Turns the solvent rack degasser on or off. For pumps operated with an SRD-3000 solvent rack, the setting should always be On . For pumps that are operated without degasser or with an external third-party degasser, select Off . |
| DegasserVacuum | Indicates the status of the degasser vacuum in the solvent rack (OK or NotOK ; read-only). |
| SolventRackLeak | Indicates the status of the leak sensor in the solvent rack (Leak or NoLeak; read-only). |

Please note: Dionex recommends always leaving the degasser on.

Turning off the pump, automatically turns off the degasser, too.

| UltiMate 3000 - System Installation and Application | on |
|---|----|
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |

5 Soft Key Menus

For all UltiMate 3000 system modules except the solvent rack and UV detector, four magnetic buttons (soft keys) under the front panel display provide access to various menus, allowing you to change certain settings and/or perform certain commands.

A magnetic Dionex menu pen (part no. 6300.0100) is included in each module's accessories kit. In an UltiMate 3000 system, another menu pen is installed in the solvent rack (\rightarrow Fig. 26).




Fig. 26: Menu pen in the solvent rack

Touch the front panel just under the display with the Dionex menu pen to show the soft key menus (\rightarrow Fig. 27):




Fig. 27: Accessing the soft key menus (here: on the autosampler)

| UltiMate 3000 - System Installation and Ap | plication | |
|--|-----------|--|
| | | |
| | | |
| | | |
| | | |
| | | |
| | | |
| | | |
| | | |
| | | |
| | | |
| | | |
| | | |
| | | |
| | | |
| | | |
| | | |
| | | |
| | | |
| | | |
| | | |
| | | |
| | | |
| | | |
| | | |
| | | |
| | | |
| | | |
| | | |
| | | |
| | | |
| | | |
| | | |
| | | |
| | | |
| | | |
| | | |
| | | |
| | | |
| | | |
| | | |
| | | |
| | | |
| | | |
| | | |
| | | |

6 Troubleshooting

The following table provides a summary of possible operating problems, lists probable causes, and suggests remedial actions:

| Problem | Probable Cause | Remedial Action |
|---|--|---|
| No display on the LCD. | The instrument is not connected to the mains. | Connect the power cord. |
| | The power is turned off. | Turn on the power. |
| | The instrument is in standby mode. | Press the Standby key on the instrument's front panel. |
| | The brightness and/or contrast is not adjusted correctly. | Adjust the brightness and/or contract (→ Contrast and Brightness properties in Chromeleon). |
| | The fuses blow. | Replace the fuses (→ Operating Instructions for the instrument) |
| | Replacement fuse blows immediately. | Contact Dionex Service. |
| | Error in the electronic system. | Contact Dionex Service. |
| The instrument does not work correctly when controlled by Chromeleon. | There is no connection between the instrument and the Chromeleon server PC. | Check the connection to the server PC. |
| | The USB port on the computer is not ready for operation. | Check the USB port on the computer. |
| The system has very high backpressure. | One or more capillaries in the system are blocked or damaged by bending. | Check the capillaries in the system, remove the blockage, or replace the capillaries. |
| The system has very high backpressure at the column and at the pump. | The column is contaminated or blocked. | Rinse or replace the column. |
| The system has very high backpressure at the pump (the pressure at the column is normal). | The filter frit (porosity: 0.5µm) on the pump's high-pressure side is dirty. | Replace the filter frit. |
| High baseline drift | The column is contaminated. | Clean or replace the column. |
| | The system is not sufficiently equilibrated. | Flush the system until equilibration. |
| | The solvent is degraded. | Use fresh solvent and check the eluent filter frits. In aqueous solvents, growth of microorganisms is possible. |
| | The environmental conditions are unstable. | Make sure that the temperature and the humidity are constant. Use fresh solvent the solvent |
| | Solvent is degraded. | Use Iresh solvent the solvent |
| | The detector is not warmed up. | Allow the full detector warm-up time. |
| | Detector problem. | \rightarrow Detector manual. |

| Problem | Probable Cause | Remedial Action |
|--|---|--|
| High noise level, pulsation | There are pressure fluctuations from the pump. | Prime the pump; check general function (→ <i>Operating Instructions</i> for the pump). |
| | There are gas bubbles in the system. | Prime the system (\rightarrow <i>Operating Instructions</i> for the pump). |
| | The solvent is degraded. | Use fresh solvent. |
| | The detector is defective. | Contact Dionex Service. |
| | The lamp is too old. | Replace the detector lamp. |
| Peak Broadening, increased dead time | The capillary's inner diameter is too large. | Change the capillary. |
| | The filter frits on the solvent lines are clogged. | Check the filter for permeability. Replace the filter frit if necessary (→ Operating Instructions for the pump). |
| | The capillaries are clogged. | Replace the capillaries. |
| | The sample loop is clogged. | Replace the sample loop (→ Operating Instructions for the autosampler). |
| | The proportioning valve is defective. | Contact Dionex Service. |
| | The column is overloaded or contaminated. | Clean or replace the column. |
| | The solvent is degraded. | Use fresh solvent. |
| | The flow splitter is clogged or defective. | Replace the splitter cartridge (|
| Triangular peaks | The column is overloaded (the protein concentration is too high). | Dilute the sample. |
| Additional peaks appear in the injection peak. | With gradients, the equilibration time after the flush cycle is too short or the dead volume is too high. | Increase the equilibration time and/or eliminate possibly existing dead volume. |
| Poor reproducibility | The autosampler draws air from the vial. | There is not enough sample in the vial, the settings for the SampleHeight and/or TransLiquidHeight autosampler parameter are incorrect, or there are too many replicates (\rightarrow Operating Instructions for the autosampler). |
| | There are gas bubbles in the syringe. | Flush the syringe (\rightarrow <i>Operating Instructions</i> for the autosampler). |
| | There is an gas bubble in the flow path. | Perform a wash cycle (→ Operating Instructions for the autosampler). |
| | The autosampler, the injection valve, or the syringe valve is not tight. | → Operating Instructions for the autosampler. |

| Problem | Probable Cause | Remedial Action |
|-------------------------------------|--|---|
| Poor reproducibility (cont'd) | Carry-over occurs in the system. | Flush the needle using an appropriate solvent. |
| | The capillary connections are not installed properly or they are not | Check and tighten the capillary connections. |
| | tight. | Exchange the needle if necessary |
| | | (→ Operating Instructions for the autosampler). |
| | There are dead volumes in the capillary connections. | Replace the fittings. |
| | The pump head seals are not tight. | Replace the seals (\rightarrow Operating Instructions for the pump). |
| | The gradient is irreproducible. | Change the gradient. |
| | The sample is unstable and decomposes. | Use new sample or change the conditions. |
| | The gradient is not reproducible. | Check the pump function and degassing. |
| | | Check the suction frits for obstruction and exchange as necessary. |
| | The environmental conditions are unstable. | Make sure that the temperature and air humidity are constant. (Use a column thermostat.) |
| | Contamination occurs somewhere in the system. | Flush the system using an appropriate solvent. |
| No flow | The system is leaking. | Find and eliminate the leak. |
| | The check valves are not | Correctly install or replace the |
| | installed properly in the pump | check valves $(\rightarrow Operating$ |
| | (not in direction of flow) or defective. | <i>Instructions</i> for the pump). |
| | There is air in the solvent or in | Purge the pump (\rightarrow <i>Operating</i> |
| | the pump head(s). | Instructions for the pump). |
| | The pump's inline filter is clogged. | Replace the filter (\rightarrow Operating Instructions for the pump). |
| Flow fluctuation | The inlet path is blocked. | Check the pump's inlet lines, filter, proportioning valve etc. for signs of blockage. |
| | There is air in the inlet path. | Purge the pump (→ Operating Instructions for the pump) and check the degasser. |
| | The pump's check valves are dirty or defective. | Clean or replace the check valves if necessary (→ <i>Operating Instructions</i> for the pump). |
| | The pump head seals are not tight. | Replace the seals (\rightarrow Operating Instructions for the pump). |
| Poor degassing | There is a leak in the capillaries or solvent lines or there are loose connections. | Inspect the capillary and solvent tube connections for leakage; tighten loose fitting connections. Reduce the flow rate. |
| LIV cional drift during the | The flow rate is too high. | Increase or reduce the TFA |
| UV signal drift during the gradient | In applications with TFA, the TFA concentration is incorrect. If the signal drift upward, solvent B contains too much TFA and vice | concentration as needed. |

| Problem | Probable Cause | Remedial Action |
|---|--|--|
| After switching from the loading pump flow to the nano flow, there is a drop in the UV signal, and the signal does not return to the initial level fast enough. | There is dead volume in the connection to the precolumn. | Slide the capillary into the column as far as it will go. Check and replace the precolumn if necessary. |
| | The precolumn is defective. | Replace the precolumn. |
| No peaks or only few, poorly resolved peaks in the chromatogram | The precolumn is too short, the flow rate of the loading pump is too high and/or the loading time is too long. (The three items interact: If the precolumn is too short, peak may be lost because they are not retained. If the loading time is too long or if the loading flow is too high, hydrophilic peaks may be lost. This also affects the reproducibility because the proteins are transported a long way into the precolumn.) | Consider the three items. Use a longer precolumn, reduce the flow rate from the loading pump and/or reduce the loading time. |
| Reproducible ghost peaks in the chromatogram. | The degasser channels are contaminated. | Rinse the degasser channels (→ Operating Instructions for the solvent rack). |
| | The solvents are degraded or dirty or their purity is insufficient. | Use fresh solvents. |
| | Contamination occurs somewhere in the system. | Flush the system using an appropriate solvent. |
| The temperature in the flow manager does not change for some time although the temperature set point has not been reached. | The ambient temperature is too high. | Reduce the ambient temperature (e.g., by ventilating the room). |
| | The ventilation slots on the sides and/or bottom of the instrument are obstructed. | Make sure that the ventilation slots are not obstructed in any way. |
| | The front panel door is not completely closed. | When closing the door, make sure that the locking mechanism locks the door. |
| | The capillaries are placed in such a manner that ambient air can enter the column chamber. | Make sure that the capillaries rest flat on the edge of the housing. |
| | The door seal is damaged. | Contact Dionex Service. |

Please note: For information about the error messages that may appear in the Chromeleon Audit Trail if an error occurs or if communication between the instrument and Chromeleon cannot be established, refer to the associated sections in the *Operating Instructions* for the instruments.

7 Routine Maintenance

The system modules of the UltiMate 3000 system are made of high-quality components and materials to minimize maintenance requirements. The painted surfaces, as well as the display, are relatively resistant to weak acids, alkali, and organic solvents. Nevertheless, immediately wipe up all liquids spilled onto the modules' surface, using lint-free cloth or paper. If surfaces are exposed for longer periods, these liquids can cause damage.

The *Operating Instructions* for the individual instruments provide instructions for shutting down the instrument and describe all maintenance procedures that can be carried out by the user. Dionex personnel should perform any additional servicing, as well as annular inspections to detect contamination, wear, etc. If unexpected problems occur, please contact Dionex Service.

Please note: Before you return any instrument to Dionex for repair, contact Dionex Service or your local distributor. An RMA (Return Material Authorization) number is required in order to track your instrument. Always use the original packaging when shipping the module. Shipping the instrument in anything other than the original packaging will void the warranty. Refer to the warranty statement in the terms of sale for more information.

| UltiMate 3000 - System Installation and Application | | | |
|---|--|--|--|
| | | | |
| | | | |
| | | | |
| | | | |
| | | | |
| | | | |
| | | | |
| | | | |
| | | | |
| | | | |
| | | | |
| | | | |
| | | | |
| | | | |
| | | | |
| | | | |
| | | | |
| | | | |
| | | | |
| | | | |
| | | | |
| | | | |
| | | | |
| | | | |
| | | | |
| | | | |
| | | | |
| | | | |
| | | | |
| | | | |
| | | | |
| | | | |
| | | | |
| | | | |
| | | | |
| | | | |
| | | | |
| | | | |
| | | | |
| | | | |
| | | | |
| | | | |
| | | | |
| | | | |
| | | | |

8 Accessories and Spare Parts

Dionex accessories and spare parts are always maintained at the latest technical standard. Therefore, part numbers are subject to alteration. However, updated parts will always be compatible with the parts they replace.

- Please note: The part numbers always refer to the packing unit. Unless stated otherwise, the packing unit is one unit.
- Please note: For more information about which accessories and/or spare parts are available for the individual modules, refer to the *Operating Instructions* for the instruments.
- Please note: Use original Dionex spare parts only. Substituting non-Dionex parts or using non-Dionex accessories may impair the performance of the system.

8.1 Basic Configuration Kits

8.1.1 Standard Devices

Nano HPLC:

| Description | Part No. | Qty in Kit |
|---|-----------|------------|
| Basic Configuration Kit for Nano HPLC, including: | 6720.0065 | |
| Capillary from pump to flow manager: | | |
| Long connection (system with detector) | 6035.2550 | 1 |
| Short connection (system without detector) | 6035.2553 | 1 |
| Manual Injection Port | 6720.9007 | 1 |
| Capillary from flow splitter to switching valve (PEEKsil, 20 µm ID, 30 cm long) | 6720.0033 | 1 |
| Teflon [®] tubing (500 μm ID, 100 cm long) for use as waste tubing | 6720.0077 | 2 |
| Universal fitting (1/16", finger-tight, extra long thread; 4 fittings) | 6720.0015 | 1 |
| PEEK sleeves for connection with fused silica capillaries (280 µm O.D., set of 5) | 6720.0064 | 1 |
| Ferrule (1/16"; stainless steel, set of 10) | 6720.0017 | 5 ferrules |
| Nut (1/16", 19 mm long, stainless steel, set of 10) | 6720.0019 | 5 nuts |
| Teflon [®] tubing (250 μm ID, low pressure connection of 280 μm OD fused silica capillaries, set of 5) | 6720.0030 | 1 |
| Separation column (C18 PepMap100, 75µm ID, 15 cm long) | 160321 | 1 |
| Syringe (50 µl, with special needle) | 713.80565 | 1 |
| Cytochrome C Digest test sample | 161089 | 1 |
| Polypropylene vials with glass insert (set of 25 vials) | 6820.0027 | 1 |
| Caps for polypropylene vials (set of 25 caps) | 6820.0028 | 1 |
| Sample loop (PEEKsil, 1µl) | 6820.0015 | 1 |
| End cap (set of 5 caps) | 164140 | 1 |

Capillary HPLC:

| Description | Part No. | Qty in Kit |
|---|-----------|------------|
| Basic Configuration Kit for Capillary HPLC, including: | 6720.0082 | |
| Capillary from pump to flow manager: | | |
| Long connection (system with detector) | 6035.2550 | 1 |
| Short connection (system without detector) | 6035.2553 | 1 |
| Capillary from flow splitter to switching valve | 6720.0034 | 1 |
| (PEEKsil, 50 μm ID, 30 cm long) | | |
| Manual Injection Port | 6720.9007 | 1 |
| Teflon [®] tubing (500 μm ID, 100 cm long) for use as waste tubing | 6720.0077 | 2 |
| Universal fitting (1/16", finger-tight, extra long thread; 4 fittings) | 6720.0015 | 1 |
| Universal fitting (1/16", finger-tight, extra long body; set of 2) | 6720.0072 | 1 |
| Teflon® tubing (250 μm ID, low pressure connection of 280 μm OD fused | 6720.0030 | 1 |
| silica capillaries, set of 5) | | |
| Column (C18 PepMap100, 300μm ID, 15 cm long) | 160295 | 1 |
| Syringe (50 μl, with special needle) | 713.80565 | 1 |
| Cytochrome C Digest test sample | 161089 | 1 |
| Adaptor for capillary column connection | 6720.0040 | 1 |
| Polypropylene vials with glass insert (set of 25 vials) | 6820.0027 | 1 |

| Description | Part No. | Qty in Kit |
|---|-----------|------------|
| Caps for polypropylene vials (set of 25 caps) | 6820.0028 | 1 |
| Sample loop (PEEKsil, 5µl) | 6820.0016 | 1 |
| End cap (set of 5 caps) | 164140 | 1 |

Micro HPLC:

| Description | Part No. | Qty in Kit |
|---|-----------|------------|
| Basic Configuration Kit for Micro HPLC, including: | 6720.0083 | |
| Capillary from pump to flow manager: | | |
| Long connection (system with detector) | 6035.2550 | 1 |
| Short connection (system without detector) | 6035.2553 | 1 |
| Capillary from flow splitter to switching valve | 6720.0035 | 1 |
| (PEEKsil, 75 μm ID, 30 cm long) | | |
| Manual Injection Port | 6720.9007 | 1 |
| Teflon [®] tubing (500 μm ID, 100 cm long) for use as waste tubing | 6720.0077 | 2 |
| Universal fitting (1/16", finger-tight, extra long thread; 4 fittings) | 6720.0015 | 1 |
| Separation column (C18 PepMap100, 1.0 mm ID, 15 cm long) | 160282 | 1 |
| Syringe (50 µl, with special needle) | 713.80565 | 1 |
| Cytochrome C Digest test sample | 161089 | 1 |
| Polypropylene vials with glass insert (set of 25 vials) | 6820.0027 | 1 |
| Caps for polypropylene vials (set of 25 caps) | 6820.0028 | 1 |
| Sample loop (PEEKsil, 20 µl) | 6820.0018 | 1 |
| End cap (set of 5 caps) | 164140 | 1 |

Monolithic HPLC:

| Description | Part No. | Qty in Kit |
|--|-----------|------------|
| Basic Configuration Kit for monolithic HPLC, including: | 6720.0084 | |
| Capillary from pump to flow manager: | | |
| Long connection (system with detector) | 6035.2550 | 1 |
| Short connection (system without detector) | 6035.2553 | 1 |
| Capillary from flow splitter to switching valve (PEEKsil, 50 µm ID, 30 cm long) | 6720.0034 | 1 |
| Manual Injection Port | 6720.9007 | 1 |
| Teflon [®] tubing (500 μm ID, 100 cm long) for use as waste tubing | 6720.0077 | 2 |
| Universal fitting (1/16", finger-tight, extra long thread; 4 fittings) | 6720.0015 | 1 |
| Peek sleeves for connections with fused silica capillaries (360 µm OD, set of 5) | 6720.0078 | 1 |
| Ferrule (1/16"; stainless steel, set of 10) | 6720.0017 | 5 ferrules |
| Nut (1/16", 19 mm long, stainless steel, set of 10) | 6720.0019 | 5 nuts |
| Union (micro-tight), including 2 fittings and 1 gauge plug | 6720.0074 | 1 |
| PEEK sleeves for connections with micro-tight union (280 μm OD; set of 10) | 6720.0075 | 1 |
| PEEK sleeves for connections with micro-tight union (380 μm OD; set of 10) | 6720.0076 | 1 |
| Monolithic capillary column (200μm, ID 5 cm long) | 160409 | 1 |
| Syringe (50 µl, with special needle) | 713.80565 | 1 |
| Cytochrome C Digest test sample | 161089 | 1 |
| Polypropylene vials with glass insert (set of 25 vials) | 6820.0027 | 1 |
| Caps for polypropylene vials (set of 25 caps) | 6820.0028 | 1 |
| Sample loop (PEEKsil, 1µl) | 6820.0015 | 1 |
| End cap (set of 5 caps) | 164140 | 1 |

8.1.2 Biocompatible Devices

Nano HPLC:

| Description | Part No. | Qty in Kit |
|---|-----------|------------|
| Basic Configuration Kit (biocompatible) for Nano HPLC, including: | 6721.0065 | |
| Tool for preassembly of biocompatible fittings | 6000.0065 | 1 |
| Capillary from pump to flow manager: | | |
| Long connection (system with detector) | 6037.2550 | 1 |
| Short connection (system without detector) | 6037.2553 | 1 |
| Manual Injection Port | 6721.9007 | 1 |
| Capillary from flow splitter to switching valve (PEEKsil, 20 µm ID, 30 cm long) | 6721.0033 | 1 |
| Teflon [®] tubing (500 μm ID, 100 cm long) for use as waste tubing | 6720.0077 | 2 |
| Universal fitting (1/16", finger-tight, extra long thread; 4 fittings) | 6720.0015 | 1 |
| PEEK sleeves for connection with fused silica capillaries (280 µm O.D., set of 5) | 6720.0064 | 1 |
| Universal fitting (PEEK, 1/16"), long hex nut with ferrule (set of 2 each) | 6721.0017 | 3 |
| Teflon [®] tubing (250 μm ID, low pressure connection of 280 μm OD fused silica capillaries, set of 5) | 6720.0030 | 1 |
| Separation column (C18 PepMap100, 75μm ID, 15 cm long) | 160321 | 1 |
| Syringe (50 µl, with special needle) | 713.80565 | 1 |
| Cytochrome C Digest test sample | 161089 | 1 |
| Polypropylene vials with glass insert (set of 25 vials) | 6820.0027 | 1 |
| Caps for polypropylene vials (set of 25 caps) | 6820.0028 | 1 |
| Sample loop (PEEKsil, 1µl) | 6821.0015 | 1 |
| End cap (set of 5 caps) | 164140 | 1 |

Capillary HPLC:

| Description | Part No. | Qty in Kit |
|---|-----------|------------|
| Basic Configuration Kit (biocompatible) for Capillary HPLC, including: | 6721.0082 | |
| Capillary from pump to flow manager: | | |
| Long connection (system with detector) | 6037.2550 | 1 |
| Short connection (system without detector) | 6037.2553 | 1 |
| Capillary from flow splitter to switching valve | 6721.0034 | 1 |
| (PEEKsil, 50 μm ID, 30 cm long) | | |
| Manual Injection Port | 6721.9007 | 1 |
| Teflon [®] tubing (500 μm ID, 100 cm long) for use as waste tubing | 6720.0077 | 2 |
| Universal fitting (1/16", finger-tight, extra long thread; 4 fittings) | 6720.0015 | 1 |
| Universal fitting (1/16", finger-tight, extra long body; set of 2) | 6720.0072 | 1 |
| Teflon [®] tubing (250 μm ID, low pressure connection of 280 μm OD fused | 6720.0030 | 1 |
| silica capillaries, set of 5) | | |
| Column (C18 PepMap100, 300μm ID, 15 cm long) | 160295 | 1 |
| Syringe (50 μl, with special needle) | 713.80565 | 1 |
| Cytochrome C Digest test sample | 161089 | 1 |
| Adaptor for capillary column connection | 6721.0040 | 1 |
| Polypropylene vials with glass insert (set of 25 vials) | 6820.0027 | 1 |
| Caps for polypropylene vials (set of 25 caps) | 6820.0028 | 1 |
| Sample loop (PEEKsil, 5µl) | 6821.0016 | 1 |
| End cap (set of 5 caps) | 164140 | 1 |

Micro HPLC:

| Description | Part No. | Qty in Kit |
|---|-----------|------------|
| Basic Configuration Kit (biocompatible) for Micro HPLC, including: | 6721.0083 | |
| Tool for preassembly of biocompatible fittings | 6000.0065 | 1 |
| Capillary from pump to flow manager: | | |
| Long connection (system with detector) | 6037.2550 | 1 |
| Short connection (system without detector) | 6037.2553 | 1 |
| Capillary from flow splitter to switching valve | 6721.0035 | 1 |
| (PEEKsil, 75 μm ID, 30 cm long) | | |
| Manual Injection Port | 6721.9007 | 1 |
| Teflon [®] tubing (500 μm ID, 100 cm long) for use as waste tubing | 6720.0077 | 2 |
| Universal fitting (1/16", finger-tight, extra long thread; 4 fittings) | 6720.0015 | 1 |
| Separation column (C18 PepMap100, 1.0 mm ID, 15 cm long) | 160282 | 1 |
| Syringe (50 μl, with special needle) | 713.80565 | 1 |
| Cytochrome C Digest test sample | 161089 | 1 |
| Polypropylene vials with glass insert (set of 25 vials) | 6820.0027 | 1 |
| Caps for polypropylene vials (set of 25 caps) | 6820.0028 | 1 |
| Sample loop (PEEKsil, 20 μl) | 6821.0018 | 1 |
| End cap (set of 5 caps) | 164140 | 1 |

Monolithic HPLC:

| Description | Part No. | Qty in Kit |
|--|------------------------|------------|
| Basic Configuration Kit (biocompatible) for monolithic HPLC, including: | 6721.0084 | |
| Tool for preassembly of biocompatible fittings | 6000.0065 | 1 |
| Capillary from pump to flow manager: Long connection (system with detector) Short connection (system without detector) | 6037.2550 6037.2553 | 1 1 |
| Capillary from flow splitter to switching valve (PEEKsil, 50 µm ID, 30 cm long) | 6721.0034 | 1 |
| Manual Injection Port | 6721.9007 | 1 |
| Teflon [®] tubing (500 μm ID, 100 cm long) for use as waste tubing | 6720.0077 | 2 |
| Universal fitting (1/16", finger-tight, extra long thread; 4 fittings) | 6720.0015 | 1 |
| Peek sleeves for connections with fused silica capillaries (360 µm OD, set of 5) | 6720.0078 | 1 |
| Universal fitting (PEEK, 1/16"), long hex nut with ferrule (set of 2 each) | 6721.0017 | 3 |
| Union (micro-tight), including 2 fittings and 1 gauge plug | 6720.0074 | 1 |
| PEEK sleeves for connections with micro-tight union (280 μm OD; set of 10) | 6720.0075 | 1 |
| PEEK sleeves for connections with micro-tight union (380 μm OD; set of 10) | 6720.0076 | 1 |
| Monolithic capillary column (200μm, ID 5 cm long) | 160409 | 1 |
| Syringe (50 µl, with special needle) | 713.80565 | 1 |
| Cytochrome C Digest test sample | 161089 | 1 |
| Polypropylene vials with glass insert (set of 25 vials) | 6820.0027 | 1 |
| Caps for polypropylene vials (set of 25 caps) | 6820.0028 | 1 |
| Sample loop (PEEKsil, 1µl) | 6821.0015 | 1 |
| End cap (set of 5 caps) | 164140 | 1 |

8.2 Direct Injection Kits

8.2.1 Standard Devices

Nano HPLC:

| Description | Part No. | Qty in Kit |
|---|-----------|------------|
| Direct Injection Kit for Nano HPLC, including: | 6720.0042 | |
| Capillary from pump to flow manager: | | |
| Long connection (system with detector) | 6035.2550 | 1 |
| Short connection (system with detector) | 6035.2553 | 1 |
| Capillary from autosampler to flow manager (PEEKsil, 20µm, 50 cm long) | 6720.0027 | 1 |
| Capillary from autosampler to flow cell (PEEKsil, 20 µm ID, 50 cm long), including 1/16" zero dead volume union | 6720.0024 | 1 |
| PEEK sleeves for connection with fused silica capillaries (280 µm O.D., set of 5) | 6720.0064 | 1 |
| Ferrule (1/16"; stainless steel, set of 10) | 6720.0017 | 5 ferrules |
| Nut (1/16", 19 mm long, stainless steel, set of 10) | 6720.0019 | 5 nuts |
| Teflon [®] tubing (250 μm ID, low pressure connection of 280 μm OD fused silica capillaries, set of 5) | 6720.0030 | 1 |
| Separation column (C18 PepMap100, 75µm ID, 15 cm long) | 160321 | 1 |
| Cytochrome C Digest test sample | 161089 | 1 |
| Polypropylene vials with glass insert (set of 25 vials) | 6820.0027 | 1 |
| Caps for polypropylene vials (set of 25 caps) | 6820.0028 | 1 |
| End cap (set of 5 caps) | 164140 | 1 |

Capillary HPLC:

| Description | Part No. | Qty in Kit |
|---|-----------|------------|
| Direct Injection Kit for Capillary HPLC, including: | 6720.0043 | |
| Capillary from pump to flow manager: | | |
| Long connection (system with detector) | 6035.2550 | 1 |
| Short connection (system with detector) | 6035.2553 | 1 |
| Capillary from autosampler to flow manager (PEEKsil, 50µm ID, 50 cm long) | 6720.0028 | 1 |
| Capillary from autosampler to separation column (PEEKsil, 50 µm Id, 50 cm | 6720.0025 | 1 |
| long), including 1/16" zero dead volume union | | |
| Universal fitting (1/16", finger-tight, extra long body; set of 2) | 6720.0072 | 1 |
| Teflon [®] tubing (250 μm ID, low pressure connection of 280 μm OD fused | 6720.0030 | 1 |
| silica capillaries, set of 5) | | |
| Column (C18 PepMap100, 300µm ID, 15 cm long) | 160295 | 1 |
| Cytochrome C Digest test sample | 161089 | 1 |
| Polypropylene vials with glass insert (set of 25 vials) | 6820.0027 | 1 |
| Caps for polypropylene vials (set of 25 caps) | 6820.0028 | 1 |
| End cap (set of 5 caps) | 164140 | 1 |

Micro HPLC:

| Description | Part No. | Qty in Kit |
|---|-----------|------------|
| Direct Injection Kit for Micro HPLC, including: | 6720.0046 | |
| Capillary from pump to flow manager: | | |
| Long connection (system with detector) | 6035.2550 | 1 |
| Short connection (system with detector) | 6035.2553 | 1 |
| Capillary from autosampler to flow manager (PEEKsil, 75µm ID, 50 cm long) | 6720.0029 | 1 |
| Capillary from autosampler to separation column (PEEKsil, 75 µm ID, 50 cm | 6720.0026 | 1 |
| long), including 1/16" zero dead volume union | | |
| Universal fitting (1/16", finger-tight, extra long thread; 4 fittings) | 6720.0015 | 1 |
| Separation column (C18 PepMap100, 1.0 mm ID, 15 cm long) | 160282 | 1 |
| Cytochrome C Digest test sample | 161089 | 1 |
| Sample loop (PEEKsil, 20 μl) | 6820.0018 | 1 |
| Polypropylene vials with glass insert (set of 25 vials) | 6820.0027 | 1 |
| Caps for polypropylene vials (set of 25 caps) | 6820.0028 | 1 |
| End cap (set of 5 caps) | 164140 | 1 |

| Description | Part No. | Qty in Kit |
|---|-----------|------------|
| Direct Injection Kit for monolithic HPLC, including: | 6720.0044 | |
| Capillary from pump to flow manager: | | |
| Long connection (system with detector) | 6035.2550 | 1 |
| Short connection (system without detector) | 6035.2553 | 1 |
| Capillary from autosampler to flow manager (PEEKsil, 50µm ID, 50 cm long) | 6720.0028 | 1 |
| Capillary from autosampler to separation column (PEEKsil, 50 µm ID, 50 cm long), including 1/16" zero dead volume union | 6720.0025 | 1 |
| Peek sleeves for connections with fused silica capillaries (360 µm OD, set of 5) | 6720.0078 | 1 |
| Ferrule (1/16"; stainless steel, set of 10) | 6720.0017 | 5 ferrules |
| Nut (1/16", 19 mm long, stainless steel, set of 10) | 6720.0019 | 5 nuts |
| Union (micro-tight), including 2 fittings and 1 gauge plug | 6720.0074 | 1 |
| PEEK sleeves for connections with micro-tight union (280 μm OD; set of 10) | 6720.0075 | 1 |
| PEEK sleeves for connections with micro-tight union (380 μm OD; set of 10) | 6720.0076 | 1 |
| Monolithic capillary column (200μm, ID 5 cm long) | 160409 | 1 |
| Cytochrome C Digest test sample | 161089 | 1 |
| Polypropylene vials with glass insert (set of 25 vials) | 6820.0027 | 1 |
| Caps for polypropylene vials (set of 25 caps) | 6820.0028 | 1 |
| End cap (set of 5 caps) | 164140 | 1 |

8.2.2 Biocompatible Devices

Nano HPLC:

| Description | Part No. | Qty in Kit |
|---|-----------|------------|
| Direct Injection Kit (biocompatible) for Nano HPLC, including: | 6721.0042 | |
| Tool for preassembly of biocompatible fittings | 6000.0065 | 1 |
| Capillary from pump to flow manager: | | |
| Long connection (system with detector) | 6037.2550 | 1 |
| Short connection (system with detector) | 6037.2553 | 1 |
| Capillary from autosampler to flow manager (PEEKsil, 20µm, 50 cm long) | 6721.0027 | 1 |
| Capillary from autosampler to flow cell (PEEKsil, 20 μm ID, 50 cm long), including 1/16" PEEK union | 6721.0024 | 1 |
| PEEK sleeves for connection with fused silica capillaries (280 µm O.D., set of 5) | 6720.0064 | 1 |
| Universal fitting (PEEK, 1/16"), long hex nut with ferrule (set of 2 each) | 6721.0017 | 3 |
| Teflon [®] tubing (250 μm ID, low pressure connection of 280 μm OD fused silica capillaries, set of 5) | 6720.0030 | 1 |
| Separation column (C18 PepMap100, 75μm ID, 15 cm long) | 160321 | 1 |
| Cytochrome C Digest test sample | 161089 | 1 |
| Polypropylene vials with glass insert (set of 25 vials) | 6820.0027 | 1 |
| Caps for polypropylene vials (set of 25 caps) | 6820.0028 | 1 |
| End cap (set of 5 caps) | 164140 | 1 |

Capillary HPLC:

| Description | Part No. | Qty in Kit |
|---|------------------------|------------|
| Direct Injection Kit (biocompatible) for Capillary HPLC, including: | 6721.0043 | |
| Tool for preassembly of biocompatible fittings | 6000.0065 | 1 |
| Capillary from pump to flow manager: Long connection (system with detector) Short connection (system with detector) | 6037.2550 6037.2553 | 1 |
| Capillary from autosampler to flow manager (PEEKsil, 50µm ID, 50 cm long) | 6721.0028 | 1 |
| Capillary from autosampler to separation column (PEEKsil, 50 µm Id, 50 cm long), including 1/16" PEEK union | 6721.0025 | 1 |
| Universal fitting (1/16", finger-tight, extra long body; set of 2) | 6720.0072 | 1 |
| Teflon [®] tubing (250 μm ID, low pressure connection of 280 μm OD fused silica capillaries, set of 5) | 6720.0030 | 1 |
| Column (C18 PepMap100, 300μm ID, 15 cm long) | 160295 | 1 |
| Cytochrome C Digest test sample | 161089 | 1 |
| Polypropylene vials with glass insert (set of 25 vials) | 6820.0027 | 1 |
| Caps for polypropylene vials (set of 25 caps) | 6820.0028 | 1 |
| End cap (set of 5 caps) | 164140 | 1 |

Micro HPLC:

| Description | Part No. | Qty in Kit |
|---|------------------------|------------|
| Direct Injection Kit (biocompatible) for Micro HPLC, including: | 6721.0046 | |
| Tool for preassembly of biocompatible fittings | 6000.0065 | 1 |
| Capillary from pump to flow manager: Long connection (system with detector) Short connection (system with detector) | 6037.2550 6037.2553 | 1 |
| Capillary from autosampler to flow manager (PEEKsil, 75µm ID, 50 cm long) | 6721.0029 | 1 |
| Capillary from autosampler to separation column (PEEKsil, 75 μ m ID, 50 cm long), including 1/16" PEEK union | 6721.0026 | 1 |
| Universal fitting (1/16", finger-tight, extra long thread; 4 fittings) | 6720.0015 | 1 |
| Separation column (C18 PepMap100, 1.0 mm ID, 15 cm long) | 160282 | 1 |
| Cytochrome C Digest test sample | 161089 | 1 |
| Sample loop (PEEKsil, 20 µl) | 6821.0018 | 1 |
| Polypropylene vials with glass insert (set of 25 vials) | 6820.0027 | 1 |
| Caps for polypropylene vials (set of 25 caps) | 6820.0028 | 1 |
| End cap (set of 5 caps) | 164140 | 1 |

| Description | Part No. | Qty in Kit |
|---|-----------|------------|
| Direct Injection Kit (biocompatible) for monolithic HPLC, including: | 6721.0044 | |
| Tool for preassembly of biocompatible fittings | 6000.0065 | 1 |
| Capillary from pump to flow manager: | | |
| Long connection (system with detector) | 6037.2550 | 1 |
| Short connection (system without detector) | 6037.2553 | 1 |
| Capillary from autosampler to flow manager (PEEKsil, 50µm ID, 50 cm long) | 6721.0028 | 1 |
| Capillary from autosampler to separation column (PEEKsil, 50 µm ID, 50 cm long), including 1/16" PEEK union | 6721.0025 | 1 |
| Peek sleeves for connections with fused silica capillaries (360 µm OD, set of 5) | 6720.0078 | 1 |
| Universal fitting (PEEK, 1/16"), long hex nut with ferrule (set of 2 each) | 6721.0017 | 3 |
| Union (micro-tight), including 2 fittings and 1 gauge plug | 6720.0074 | 1 |
| PEEK sleeves for connections with micro-tight union (280 μm OD; set of 10) | 6720.0075 | 1 |
| PEEK sleeves for connections with micro-tight union (380 μm OD; set of 10) | 6720.0076 | 1 |
| Monolithic capillary column (200μm, ID 5 cm long) | 160409 | 1 |
| Cytochrome C Digest test sample | 161089 | 1 |
| Polypropylene vials with glass insert (set of 25 vials) | 6820.0027 | 1 |
| Caps for polypropylene vials (set of 25 caps) | 6820.0028 | 1 |

8.3 Preconcentration Kits

8.3.1 Standard Devices

Nano HPLC:

| Description | Part No. | Qty in Kit |
|--|-----------|------------|
| Preconcentration Kit for Nano HPLC, including: | 6720.0047 | |
| Capillary from pump to flow manager: | | |
| Long connection (system with detector) | 6035.2550 | 1 |
| Short connection (system without detector) | 6035.2553 | 1 |
| Capillary from loading pump outlet block to autosampler switching valve (PEEK, 130 µm ID, 75 cm long) | 6720.0032 | 1 |
| Capillary from autosampler to flow manager switching valve (PEEKsil, 75 µm ID, 50 cm long) | 6720.0038 | 1 |
| Capillary from flow splitter to switching valve (PEEKsil, 20 µm ID, 30 cm long) | 6720.0033 | 1 |
| Capillary to nano column in flow manager (PEEKsil, 20 µm ID, 50 cm long, including 1/16" zero dead volume fitting) | 6720.0024 | 1 |
| Capillary to nano column in flow manager (PEEKsil, 20 µm ID, 50 cm long) | 6720.0027 | 1 |
| PEEK sleeves for connection with fused silica capillaries (280 µm O.D., set of 5) | 6720.0064 | 1 |
| Ferrule (1/16"; stainless steel, set of 10) | 6720.0017 | 5 ferrules |
| Nut (1/16", 19 mm long, stainless steel, set of 10) | 6720.0019 | 5 nuts |
| Universal fitting (1/16", finger-tight, extra long thread; 4 fittings) | 6720.0015 | 1 |
| Teflon [®] tubing (250 μm ID, low pressure connection of 280 μm OD fused silica capillaries, set of 5) | 6720.0030 | 1 |
| Teflon [®] tubing (500 μm ID, 100 cm long) for use as waste tubing | 6720.0077 | 2 |
| μ-Precolumn [™] holder (5 mm), including connecting tubing (30 μm ID) | 6720.0012 | 1 |
| Precolumn (C18 PepMap100, 300 μm ID, 5 mm long, set of 5) | 160454 | 1 |
| Separation column (C18 PepMap100, 75µm ID, 15 cm long) | 160321 | 1 |
| Sample loop (20 µl, PEEK) | 6820.0018 | 1 |
| Cytochrome C Digest test sample | 161089 | 1 |
| Polypropylene vials with glass insert (set of 25 vials) | 6820.0027 | 1 |
| Caps for polypropylene vials (set of 25 caps) | 6820.0028 | 1 |
| End cap (set of 5 caps) | 164140 | 1 |

Capillary HPLC:

| Description | Part No. | Qty in Kit |
|--|------------------------|------------|
| Preconcentration Kit for Capillary HPLC, including: | 6720.0048 | |
| Capillary from pump to flow manager: Long connection (system with detector) Short connection (system without detector) | 6035.2550 6035.2553 | 1 |
| Capillary from loading pump outlet block to autosampler switching valve (PEEK, 130 µm ID, 75 cm long) | 6720.0032 | 1 |
| Capillary from autosampler to flow manager switching valve (PEEKsil, 75 µm ID, 50 cm long) | 6720.0038 | 1 |
| Capillary from flow splitter to switching valve (PEEKsil, 50 µm ID, 30 cm long) | 6720.0034 | 1 |
| Adaptor from capillary column to flow manager switching valve (PEEKsil, 50 µm ID, 10 cm long, 1/16") | 6720.0040 | 1 |

| Description | Part No. | Qty in Kit |
|---|-----------|------------|
| Capillary to capillary column in flow manager (PEEKsil, 20 µm ID, 50 cm long, including 1/16" zero dead volume fitting) | 6720.0025 | 1 |
| Capillary to capillary column in flow manager (PEEKsil, 20 µm ID, 50 cm long) | 6720.0028 | 1 |
| Universal fitting (1/16", finger-tight, extra long body; set of 2) | 6720.0072 | 1 |
| Teflon® tubing (250 μm ID, low pressure connection of 280 μm OD fused silica capillaries, set of 5) | 6720.0030 | 1 |
| μ-Precolumn [™] holder (5 mm), including connecting tubing (30 μm ID) | 6720.0012 | 1 |
| Precolumn (C18 PepMap100, 300 μm ID, 5 mm long, set of 5) | 160454 | 1 |
| Teflon® tubing (500 μm ID, 100 cm long) for use as waste tubing | 6720.0077 | 2 |
| Universal fitting (1/16", finger-tight, extra long thread; 4 fittings) | 6720.0015 | 1 |
| Ferrule (1/16"; stainless steel, set of 10) | 6720.0017 | 5 ferrules |
| Nut (1/16", 19 mm long, stainless steel, set of 10) | 6720.0019 | 5 nuts |
| Column (C18 PepMap100, 300μm ID, 15 cm long) | 160295 | 1 |
| Sample loop (20 µl, PEEK) | 6820.0018 | 1 |
| Cytochrome C Digest test sample | 161089 | 1 |
| Polypropylene vials with glass insert (set of 25 vials) | 6820.0027 | 1 |
| Caps for polypropylene vials (set of 25 caps) | 6820.0028 | 1 |
| End cap (set of 5 caps) | 164140 | 1 |

| Description | Part No. | Qty in Kit |
|---|-----------|------------|
| Preconcentration Kit for Monolithic HPLC, including: | 6720.0049 | |
| Capillary from pump to flow manager: | | |
| Long connection (system with detector) | 6035.2550 | 1 |
| Short connection (system without detector) | 6035.2553 | 1 |
| Capillary from loading pump outlet block to autosampler switching valve (PEEK, 130 µm ID, 75 cm long) | 6720.0032 | 1 |
| Capillary from autosampler to flow manager switching valve (PEEKsil, 75 µm ID, 50 cm long) | 6720.0038 | 1 |
| Capillary from flow splitter to switching valve (PEEKsil, 50 µm ID, 30 cm long) | 6720.0034 | 1 |
| Peek sleeves for connections with fused silica capillaries (360 µm OD, set of 5) | 6720.0078 | 1 |
| Capillary to capillary column in flow manager (PEEKsil, 20 µm ID, 50 cm long, including 1/16" zero dead volume fitting) | 6720.0025 | 1 |
| Capillary to capillary column in flow manager (PEEKsil, 20 µm ID, 50 cm long) | 6720.0028 | 1 |
| Ferrule (1/16"; stainless steel, set of 10) | 6720.0017 | 5 ferrules |
| Nut (1/16", 19 mm long, stainless steel, set of 10) | 6720.0019 | 5 nuts |
| Union (micro-tight), including 2 fittings and 1 gauge plug | 6720.0074 | 1 |
| PEEK sleeves for connections with micro-tight union (280 μm OD; set of 10) | 6720.0075 | 1 |
| PEEK sleeves for connections with micro-tight union (380 μm OD; set of 10) | 6720.0076 | 1 |
| Teflon® tubing (500 μm ID, 100 cm long) for use as waste tubing | 6720.0077 | 2 |
| Universal fitting (1/16", finger-tight, extra long thread; 4 fittings) | 6720.0015 | 1 |
| Monolithic column (200 μm ID, 5 cm long) | 161409 | 1 |
| Monolithic trap column (200 μm ID, 5 mm long) | 163972 | 1 |
| Sample loop (20 µl, PEEK) | 6820.0018 | 1 |

| Description | Part No. | Qty in Kit |
|---|-----------|------------|
| Cytochrome C Digest test sample | 161089 | 1 |
| Polypropylene vials with glass insert (set of 25 vials) | 6820.0027 | 1 |
| Caps for polypropylene vials (set of 25 caps) | 6820.0028 | 1 |
| End cap (set of 5 caps) | 164140 | 1 |

8.3.2 Biocompatible Devices

Nano HPLC:

| Description | Part No. | Qty in Kit |
|---|-----------|------------|
| Preconcentration Kit (biocompatible) for Nano HPLC, including: | 6721.0047 | |
| Tool for preassembly of biocompatible fittings | 6000.0065 | 1 |
| Capillary from pump to flow manager: | | |
| Long connection (system with detector) | 6037.2550 | 1 |
| Short connection (system without detector) | 6037.2553 | 1 |
| Capillary from loading pump outlet block to autosampler switching valve (PEEK, 130 µm ID, 75 cm long) | 6721.0032 | 1 |
| Capillary from autosampler to flow manager switching valve (PEEKsil, 75 µm ID, 50 cm long) | 6721.0038 | 1 |
| Capillary from flow splitter to switching valve (PEEKsil, 20 µm ID, 30 cm long) | 6721.0033 | 1 |
| Capillary to nano column in flow manager (PEEKsil, 20 µm ID, 50 cm long, including 1/16" PEEK fitting) | 6721.0024 | 1 |
| Capillary to nano column in flow manager (PEEKsil, 20 µm ID, 50 cm long) | 6721.0027 | 1 |
| PEEK sleeves for connection with fused silica capillaries (280 µm O.D., set of 5) | 6720.0064 | 1 |
| Universal fitting (PEEK, 1/16"), long hex nut with ferrule (set of 2 each) | 6721.0017 | 3 |
| Universal fitting (1/16", finger-tight, extra long thread; 4 fittings) | 6720.0015 | 1 |
| Teflon [®] tubing (250 μm ID, low pressure connection of 280 μm OD fused silica capillaries, set of 5) | 6720.0030 | 1 |
| Teflon [®] tubing (500 μm ID, 100 cm long) for use as waste tubing | 6720.0077 | 2 |
| μ-Precolumn [™] holder (5 mm), including connecting tubing (30 μm ID) | 6720.0012 | 1 |
| Precolumn (C18 PepMap100, 300 μm ID, 5 mm long, set of 5) | 160454 | 1 |
| Separation column (C18 PepMap100, 75µm ID, 15 cm long) | 160321 | 1 |
| Sample loop (20 µl, PEEK) | 6821.0018 | 1 |
| Cytochrome C Digest test sample | 161089 | 1 |
| Polypropylene vials with glass insert (set of 25 vials) | 6820.0027 | 1 |
| Caps for polypropylene vials (set of 25 caps) | 6820.0028 | 1 |
| End cap (set of 5 caps) | 164140 | 1 |

Capillary HPLC:

| Description | Part No. | Qty in Kit |
|--|------------------------|------------|
| Preconcentration Kit (biocompatible) for Capillary HPLC, including: | 6721.0048 | |
| Tool for preassembly of biocompatible fittings | 6000.0065 | 1 |
| Capillary from pump to flow manager: Long connection (system with detector) Short connection (system without detector) | 6037.2550 6037.2553 | 1 |
| Capillary from loading pump outlet block to autosampler switching valve (PEEK, 130 µm ID, 75 cm long) | 6721.0032 | 1 |
| Capillary from autosampler to flow manager switching valve (PEEKsil, 75 µm ID, 50 cm long) | 6721.0038 | 1 |
| Capillary from flow splitter to switching valve (PEEKsil, 50 µm ID, 30 cm long) | 6721.0034 | 1 |
| Adaptor from capillary column to flow manager switching valve (PEEKsil, 50 μm ID, 10 cm long, 1/16") | 6721.0040 | 1 |
| Capillary to capillary column in flow manager (PEEKsil, 20 µm ID, 50 cm long, including 1/16" PEEK fitting) | 6721.0025 | 1 |
| Capillary to capillary column in flow manager (PEEKsil, 20 µm ID, 50 cm long) | 6721.0028 | 1 |
| Universal fitting (1/16", finger-tight, extra long body; set of 2) | 6720.0072 | 1 |
| Teflon [®] tubing (250 μm ID, low pressure connection of 280 μm OD fused silica capillaries, set of 5) | 6720.0030 | 1 |
| μ-Precolumn [™] holder (5 mm), including connecting tubing (30 μm ID) | 6720.0012 | 1 |
| Precolumn (C18 PepMap100, 300 μm ID, 5 mm long, set of 5) | 160454 | 1 |
| Teflon [®] tubing (500 μm ID, 100 cm long) for use as waste tubing | 6720.0077 | 2 |
| Universal fitting (1/16", finger-tight, extra long thread; 4 fittings) | 6720.0015 | 1 |
| Universal fitting (PEEK, 1/16"), long hex nut with ferrule (set of 2 each) | 6721.0017 | 3 |
| Column (C18 PepMap100, 300μm ID, 15 cm long) | 160295 | 1 |
| Sample loop (20 μl, PEEK) | 6821.0018 | 1 |
| Cytochrome C Digest test sample | 161089 | 1 |
| Polypropylene vials with glass insert (set of 25 vials) | 6820.0027 | 1 |
| Caps for polypropylene vials (set of 25 caps) | 6820.0028 | 1 |
| End cap (set of 5 caps) | 164140 | 1 |

| Description | Part No. | Qty in Kit |
|--|------------------------|------------|
| Preconcentration Kit (biocompatible) for Monolithic HPLC, including: | 6721.0049 | |
| Tool for preassembly of biocompatible fittings | 6000.0065 | 1 |
| Capillary from pump to flow manager: Long connection (system with detector) Short connection (system without detector) | 6037.2550 6037.2553 | 1 1 |
| Capillary from loading pump outlet block to autosampler switching valve (PEEK, 130 µm ID, 75 cm long) | 6721.0032 | 1 |
| Capillary from autosampler to flow manager switching valve (PEEKsil, 75 µm ID, 50 cm long) | 6721.0038 | 1 |
| Capillary from flow splitter to switching valve (PEEKsil, 50 µm ID, 30 cm long) | 6721.0034 | 1 |
| Peek sleeves for connections with fused silica capillaries (360 µm OD, set of 5) | 6720.0078 | 1 |

| Description | Part No. | Qty in Kit |
|---|-----------|------------|
| Capillary to capillary column in flow manager (PEEKsil, 20 µm ID, 50 cm long, including 1/16" PEEK fitting) | 6721.0025 | 1 |
| Capillary to capillary column in flow manager (PEEKsil, 20 µm ID, 50 cm long) | 6721.0028 | 1 |
| Universal fitting (PEEK, 1/16"), long hex nut with ferrule (set of 2 each) | 6721.0017 | 3 |
| Union (micro-tight), including 2 fittings and 1 gauge plug | 6720.0074 | 1 |
| PEEK sleeves for connections with micro-tight union (280 μm OD; set of 10) | 6720.0075 | 1 |
| PEEK sleeves for connections with micro-tight union (380 μm OD; set of 10) | 6720.0076 | 1 |
| Teflon [®] tubing (500 μm ID, 100 cm long) for use as waste tubing | 6720.0077 | 2 |
| Universal fitting (1/16", finger-tight, extra long thread; 4 fittings) | 6720.0015 | 1 |
| Monolithic column (200 μm ID, 5 cm long) | 161409 | 1 |
| Monolithic trap column (200 μm ID, 5 mm long) | 163972 | 1 |
| Sample loop (20 µl, PEEK) | 6821.0018 | 1 |
| Cytochrome C Digest test sample | 161089 | 1 |
| Polypropylene vials with glass insert (set of 25 vials) | 6820.0027 | 1 |
| Caps for polypropylene vials (set of 25 caps) | 6820.0028 | 1 |
| End cap (set of 5 caps) | 164140 | 1 |

8.4 2D LC Salt Plug Kit

8.4.1 Standard Devices

Nano HPLC:

| Description | Part No. | Qty in Kit |
|---|-----------|------------|
| 2D LC Salt Plug Kit for Nano HPLC, including: | 6720.0051 | |
| Capillary from pump to flow manager: | | |
| Long connection (system with detector) | 6035.2550 | 1 |
| Short connection (system without detector) | 6035.2553 | 1 |
| Capillary from loading pump outlet block to autosampler switching valve (PEEK, 130 µm ID, 75 cm long) | 6720.0032 | 1 |
| Capillary from autosampler to flow manager switching valve (PEEKsil, 75 µm ID, 50 cm long) | 6720.0038 | 1 |
| Capillary from flow splitter to switching valve (PEEKsil, 20 µm ID, 30 cm long) | 6720.0033 | 1 |
| Capillary to nano column in flow manager (PEEKsil, 20 µm ID, 50 cm long, including 1/16" zero dead volume fitting) | 6720.0024 | 1 |
| Capillary to nano column in flow manager (PEEKsil, 20 µm ID, 50 cm long) | 6720.0027 | 1 |
| PEEK sleeves for connection with fused silica capillaries (280 µm O.D., set of 5) | 6720.0064 | 1 |
| Ferrule (1/16"; stainless steel, set of 10) | 6720.0017 | 5 ferrules |
| Nut (1/16", 19 mm long, stainless steel, set of 10) | 6720.0019 | 5 nuts |
| Universal fitting (1/16", finger-tight, extra long thread; 4 fittings) | 6720.0015 | 1 |
| Teflon [®] tubing (250 μm ID, low pressure connection of 280 μm OD fused silica capillaries, set of 5) | 6720.0030 | 1 |
| Teflon [®] tubing (500 μm ID, 100 cm long) for use as waste tubing | 6720.0077 | 2 |
| PEEKsil bridge between left and right flow manager switching valves | 6720.0060 | 1 |
| (75 µm ID, 50 cm long), including appropriate fittings (stainless steel) | | |
| µ-Precolumn [™] holder (5 mm), including connecting tubing (30 μm ID) | 6720.0012 | 1 |
| Precolumn (trap column, C18 PepMap100, 300 µm ID, 5 mm long, set of 5) | 160454 | 1 |
| Separation column (C18 PepMap100, 75µm ID, 15 cm long) | 160321 | 1 |
| Ion exchange column (300 μm ID x 10 cm) packed with Poros 10 S with connections, 130 μm ID PEEK inlet (30cm) and outlet (10 cm) | 162152 | 1 |
| Sample loop (20 μl, PEEK) | 6820.0018 | 1 |
| Proteine Mixture Digest test sample | 160188 | 1 |
| Cytochrome C Digest test sample | 161089 | 1 |
| Polypropylene vials with glass insert (set of 25 vials) | 6820.0027 | 1 |
| Caps for polypropylene vials (set of 25 caps) | 6820.0028 | 1 |
| End cap (set of 5 caps) | 164140 | 1 |

8.4.2 Biocompatible Devices

Nano HPLC:

| Description | Part No. | Qty in Kit |
|---|-----------|------------|
| 2D LC Salt Plug Kit (biocompatible) for Nano HPLC, including: | 6721.0051 | |
| Tool for preassembly of biocompatible fittings | 6000.0065 | 1 |
| Capillary from pump to flow manager: | | |
| Long connection (system with detector) | 6037.2550 | 1 |
| Short connection (system without detector) | 6037.2553 | 1 |
| Capillary from loading pump outlet block to autosampler switching valve (PEEK, 130 µm ID, 75 cm long) | 6721.0032 | 1 |
| Capillary from autosampler to flow manager switching valve (PEEKsil, 75 µm ID, 50 cm long) | 6721.0038 | 1 |
| Capillary from flow splitter to switching valve (PEEKsil, 20 µm ID, 30 cm long) | 6721.0033 | 1 |
| Capillary to nano column in flow manager (PEEKsil, 20 µm ID, 50 cm long, including 1/16" PEEK fitting) | 6721.0024 | 1 |
| Capillary to nano column in flow manager (PEEKsil, 20 µm ID, 50 cm long) | 6721.0027 | 1 |
| PEEK sleeves for connection with fused silica capillaries (280 µm O.D., set of 5) | 6720.0064 | 1 |
| Universal fitting (PEEK, 1/16"), long hex nut with ferrule (set of 2 each) | 6721.0017 | 3 |
| Universal fitting (1/16", finger-tight, extra long thread; 4 fittings) | 6720.0015 | 1 |
| Teflon [®] tubing (250 μm ID, low pressure connection of 280 μm OD fused silica capillaries, set of 5) | 6720.0030 | 1 |
| Teflon [®] tubing (500 μm ID, 100 cm long) for use as waste tubing | 6720.0077 | 2 |
| PEEKsil bridge between left and right flow manager switching valves | 6721.0060 | 1 |
| (75 μm ID, 50 cm long), including appropriate fittings (stainless steel) | | |
| μ-Precolumn [™] holder (5 mm), including connecting tubing (30 μm ID) | 6720.0012 | 1 |
| Precolumn (trap column, C18 PepMap100, 300 µm ID, 5 mm long, set of 5) | 160454 | 1 |
| Separation column (C18 PepMap100, 75µm ID, 15 cm long) | 160321 | 1 |
| Ion exchange column (300 μm ID x 10 cm) packed with Poros 10 S with connections, 130 μm ID PEEK inlet (30cm) and outlet (10 cm) | 162152 | 1 |
| Sample loop (20 µl, PEEK) | 6821.0018 | 1 |
| Proteine Mixture Digest test sample | 160188 | 1 |
| Cytochrome C Digest test sample | 161089 | 1 |
| Polypropylene vials with glass insert (set of 25 vials) | 6820.0027 | 1 |
| Caps for polypropylene vials (set of 25 caps) | 6820.0028 | 1 |
| End cap (set of 5 caps) | 164140 | 1 |

8.5 Comprehensive 2D LC Kit

8.5.1 Standard Devices

Capillary/Nano HPLC:

| Description | Part No. | Qty in Kit |
|---|------------------------|------------|
| 2D LC Comprehensive Kit for Capillary/Nano HPLC, including: | 6720.0053 | |
| Capillary from loading pump to flow manager switching valve (PEEK, 130 µm ID, 45 cm long), including appropriate fitting connections | 6720.0057 | 1 |
| Capillary from flow manager to autosampler (Cap; PEEKsil, 50 µm ID, 50 cm long) | 6720.0028 | 1 |
| Capillary from autosampler to capillary column in flow manager (50 µm ID, 70 cm long), including zero dead volume union (1/16") | 6720.0058 | 1 |
| Capillary from pump to flow manager Long connection (system with detector) Short connection (system without detector) | 6035.2556 6035.2554 | 2 2 |
| Capillary from nano splitter outlet to switching valve (PEEKsil, 20 µm ID, 30 cm long) | 6720.0033 | 1 |
| Capillary to nano column in flow manager (PEEKsil, 20 µm ID, 50 cm long, including 1/16" zero dead volume fitting) | 6720.0024 | 1 |
| Capillary to nano column in flow manager (PEEKsil, 20 µm ID, 50 cm long) | 6720.0027 | 1 |
| PEEKsil bridge between left and right switching valves in flow manager (30 µm ID, 15 cm long), including appropriate fitting connection (stainless steel) | 6720.0061 | 2 |
| PEEKsil bridge between left and right switching valves in flow manager (30 µm ID, 30 cm long), including appropriate fitting connection (stainless steel) | 6720.0062 | 2 |
| PEEK sleeves for connection with fused silica capillaries (280 µm O.D., set of 5) | 6720.0064 | 1 |
| Ferrule (1/16"; stainless steel, set of 10) | 6720.0017 | 1 |
| Nut (1/16", 19 mm long, stainless steel, set of 10) | 6720.0019 | 1 |
| Universal fitting (1/16", finger-tight, extra long thread; 4 fittings) | 6720.0015 | 1 |
| Teflon [®] tubing (500 μm ID, 100 cm long) for use as waste tubing | 6720.0077 | 2 |
| Ion exchange column (Poros 10S, ID 300 μm x 15 cm) | 162122 | 1 |
| μ-Precolumn [™] holder (5 mm), including connecting tubing (30 μm ID) | 6720.0012 | 2 |
| Precolumn (trap column, C18 PepMap100, 300 μm ID, 5 mm long, set of 5) | 160454 | 2 |
| Separation column (C18 PepMap100, 75µm ID, 15 cm long) | 160321 | 1 |
| Sample loop (20 μl, PEEK) | 6820.0018 | 1 |
| Proteine Mixture Digest test sample | 160188 | 1 |
| Cytochrome C Digest test sample | 161089 | 1 |
| Polypropylene vials with glass insert (set of 25 vials) | 6820.0027 | 1 |
| Caps for polypropylene vials (set of 25 caps) | 6820.0028 | 1 |
| End cap (set of 5 caps) | 164140 | 1 |

8.5.2 Biocompatible Devices

Capillary/Nano HPLC:

| Description | Part No. | Qty in Kit |
|---|------------------------|------------|
| 2D LC Comprehensive Kit (biocompatible) for Capillary/Nano HPLC, | 6721.0053 | |
| including: | | |
| Tool for preassembly of biocompatible fittings | 6000.0065 | 1 |
| Capillary from loading pump to flow manager switching valve (PEEK, 130 µm ID, 45 cm long), including appropriate fitting connections | 6721.0057 | 1 |
| Capillary from flow manager to autosampler (Cap; PEEKsil, 50 μm ID, 50 cm long) | 6721.0028 | 1 |
| Capillary from autosampler to capillary column in flow manager (50 µm ID, 70 cm long), including PEEK union (1/16") | 6721.0058 | 1 |
| Capillary from pump to flow manager Long connection (system with detector) Short connection (system without detector) | 6037.2556 6037.2554 | 2 2 |
| Capillary from nano splitter outlet to switching valve (PEEKsil, 20 µm ID, 30 cm long) | 6721.0033 | 1 |
| Capillary to nano column in flow manager (PEEKsil, 20 µm ID, 50 cm long, including 1/16" PEEK fitting) | 6721.0024 | 1 |
| Capillary to nano column in flow manager (PEEKsil, 20 µm ID, 50 cm long) | 6721.0027 | 1 |
| PEEKsil bridge between left and right switching valves in flow manager (30 µm ID, 15 cm long), including appropriate fitting connection (stainless steel) | 6721.0061 | 1 |
| PEEKsil bridge between left and right switching valves in flow manager (30 µm ID, 30 cm long), including appropriate fitting connection (stainless steel) | 6721.0062 | 1 |
| PEEK sleeves for connection with fused silica capillaries (280 µm O.D., set of 5) | 6720.0064 | 1 |
| Universal fitting (PEEK, 1/16"), long hex nut with ferrule (set of 2 each) | 6721.0017 | 3 |
| Universal fitting (1/16", finger-tight, extra long thread; 4 fittings) | 6720.0015 | 1 |
| Teflon [®] tubing (500 μm ID, 100 cm long) for use as waste tubing | 6720.0077 | 2 |
| Ion exchange column (Poros 10S, ID 300 μm x 15 cm) | 162122 | 1 |
| µ-Precolumn [™] holder (5 mm), including connecting tubing (30 μm ID) | 6720.0012 | 1 |
| Precolumn (trap column, C18 PepMap100, 300 μm ID, 5 mm long, set of 5) | 160454 | 2 |
| Separation column (C18 PepMap100, 75μm ID, 15 cm long) | 160321 | 1 |
| Sample loop (20 μl, PEEK) | 6821.0018 | 1 |
| Proteine Mixture Digest test sample | 160188 | 1 |
| Cytochrome C Digest test sample | 161089 | 1 |
| Polypropylene vials with glass insert (set of 25 vials) | 6820.0027 | 1 |
| Caps for polypropylene vials (set of 25 caps) | 6820.0028 | 1 |
| End cap (set of 5 caps) | 164140 | 1 |

8.6 Parallel LC Kits

8.6.1 Standard Devices

Nano/Nano HPLC:

| Description | Part No. | Qty in Kit |
|--|-----------|------------|
| Parallel LC Kit, Nano/Nano, including: | 6720.0054 | |
| Capillary from loading pump to autosampler | 6720.0057 | 1 |
| (PEEK, 130 μm ID, 45 cm long), including appropriate fitting connections | | |
| Capillary from autosampler to FLM-3200 flow manager (PEEKsil, ID 75 μm, 50 cm long) | 6720.0038 | 1 |
| Capillary from pump to flow manager: | | |
| Long connection (system with detector) | 6035.2556 | 2 |
| Short connection (system without detector) | 6035.2554 | 2 |
| Capillary from nano flow splitter to flow manager switching valve (PEEKsil, 20 µm ID, 30 cm long) | 6720.0033 | 2 |
| Capillary to nano column in flow manager (PEEKsil, 20 µm ID, 50 cm long, including 1/16" zero dead volume fitting) | 6720.0024 | 1 |
| Capillary to nano column in flow manager (PEEKsil, 20 µm ID, 50 cm long) | 6720.0027 | 1 |
| PEEKsil bridge from flow manager to flow manager (75 μm ID, 30 cm long), including appropriate fitting connections (stainless steel) | 6720.0059 | 1 |
| PEEK sleeves for connection with fused silica capillaries (280 µm O.D., set of 5) | 6720.0064 | 2 |
| Ferrule (1/16"; stainless steel, set of 10) | 6720.0017 | 1 |
| Nut (1/16", 19 mm long, stainless steel, set of 10) | 6720.0019 | 1 |
| PEEK sleeve (1/32", 300 μm ID, 3 cm; set of 6) | 6720.0079 | 1 |
| Fitting (1/32") for nano switching valve (set of 6) | 6720.0080 | 1 |
| Teflon [®] tubing (500 μm ID, 100 cm long) for use as waste tubing | 6720.0077 | 3 |
| Fused silica tubing, 5m (75 μ m ID \pm 3 μ m; 280 μ m OD \pm 10 μ m) | 6720.0081 | 1 |
| Universal fitting (1/16", finger-tight, extra long thread; 4 fittings) | 6720.0015 | 1 |
| µ-Precolumn [™] holder (5 mm), including connecting tubing (30 μm ID) | 6720.0012 | 2 |
| Precolumn (trap column, C18 PepMap100, 300 μm ID, 5 mm long, set of 5) | 160454 | 2 |
| Separation column (C18 PepMap100, 75µm ID, 15 cm long) | 160321 | 2 |
| Sample loop (20 µl, PEEK) | 6820.0018 | 1 |
| Cytochrome C Digest test sample | 161089 | 1 |
| Polypropylene vials with glass insert (set of 25 vials) | 6820.0027 | 1 |
| Caps for polypropylene vials (set of 25 caps) | 6820.0028 | 1 |
| End cap (set of 5 caps) | 164140 | 1 |

Capillary/Capillary HPLC:

| Description | Part No. | Qty in Kit |
|--|-----------|------------|
| Parallel LC Kit, Cap/Cap, including: | 6720.0055 | |
| Capillary from pump to flow manager: | | |
| Long connection (system with detector) | 6035.2556 | 2 |
| Short connection (system without detector) | 6035.2554 | 2 |
| Capillary from loading pump to autosampler (PEEK, 130 µm ID, 45 cm long), including appropriate fitting connections | 6720.0057 | 1 |
| Capillary from autosampler to FLM-3200 flow manager (PEEKsil, 75 µm ID, 50 cm long) | 6720.0038 | 1 |
| PEEKsil bridge between the flow managers (75 μm ID, 30 cm long), including appropriate fitting connections (stainless steel) | 6720.0059 | 1 |
| Capillary from capillary flow splitter to switching valve (PEEKsil, 50 µm ID, 30 cm long) | 6720.0034 | 2 |
| Adaptor from capillary column to flow manager switching valve (PeekSil, 50 µm ID, 10 cm long, 1/16") | 6720.0040 | 2 |
| Capillary to capillary column in flow manager (PEEKsil, 20 µm ID, 50 cm long, including 1/16" zero dead volume fitting) | 6720.0025 | 1 |
| Capillary to capillary column in flow manager (PEEKsil, 20 µm ID, 50 cm long) | 6720.0028 | 1 |
| Universal fitting (1/16", finger-tight, extra long body; set of 2) | 6720.0072 | 2 |
| Teflon® tubing (250 μm ID, low pressure connection of 280 μm OD fused silica capillaries, set of 5) | 6720.0030 | 1 |
| Teflon® tubing (500 μm ID, 100 cm long) for use as waste tubing | 6720.0077 | 3 |
| Universal fitting (1/16", finger-tight, extra long thread; 4 fittings) | 6720.0015 | 5 ferrules |
| Ferrule (1/16"; stainless steel, set of 10) | 6720.0017 | 5 nuts |
| Nut (1/16", 19 mm long, stainless steel, set of 10) | 6720.0019 | 1 |
| μ-Precolumn [™] holder (5 mm), including connecting tubing (30 μm ID) | 6720.0012 | 2 |
| Precolumn (trap column, C18 PepMap100, 300 μm ID, 5 mm long, set of 5) | 160454 | 2 |
| Column (C18 PepMap100, 300μm ID, 15 cm long) | 160295 | 2 |
| Sample loop (20 μl, PEEK) | 6820.0018 | 1 |
| Cytochrome C Digest test sample | 161089 | 1 |
| Polypropylene vials with glass insert (set of 25 vials) | 6820.0027 | 1 |
| Caps for polypropylene vials (set of 25 caps) | 6820.0028 | 1 |
| End cap (set of 5 caps) | 164140 | 1 |

| Description | Part No. | Qty in Kit |
|--|-----------|------------|
| Parallel LC Kit, 2x Monolithic, including: | 6720.0056 | |
| Capillary from loading pump to autosampler (PEEK, 130 µm ID, 45 cm long), including appropriate fitting connections | 6720.0057 | 1 |
| Capillary from autosampler to FLM-3200 flow manager (PEEKsil, 75 µm ID, 50 cm long) | 6720.0038 | 1 |
| Capillary from pump to flow manager | | |
| Long connection (system with detector) | 6035.2556 | 2 |
| Short connection (system without detector) | 6035.2554 | 2 |
| Capillary from capillary flow splitter to switching valve (PEEKsil, 50 µm ID, 30 cm long) | 6720.0034 | 2 |
| PEEKsil bridge between the flow managers (ID 75 μm, 30 cm long), including appropriate fitting connections (stainless steel) | 6720.0059 | 1 |

| Description | Part No. | Qty in Kit |
|---|-----------|------------|
| PEEK sleeves for connection with fused silica capillaries (280 µm O.D., set of 5) | 6720.0064 | 1 |
| Capillary to capillary column in flow manager (PEEKsil, 20 µm ID, 50 cm long, including 1/16" zero dead volume fitting) | 6720.0025 | 1 |
| Capillary to capillary column in flow manager (PEEKsil, 20 µm ID, 50 cm long) | 6720.0028 | 1 |
| Ferrule (1/16"; stainless steel, set of 10) | 6720.0017 | 1 |
| Nut (1/16", 19 mm long, stainless steel, set of 10) | 6720.0019 | 1 |
| PEEK sleeve (1/32", ID 300 μm, 3 cm; set of 6) | 6720.0079 | 1 |
| PEEK sleeve (1/32", ID 400 μm, 3 cm; set of 6) | 6720.0041 | 1 |
| Fitting (1/32") for nano switching valve (set of 6) | 6720.0080 | 1 |
| Teflon [®] tubing (500 μm ID, 100 cm long) for use as waste tubing | 6720.0077 | 3 |
| Fused silica tubing, 5m (75 μ m ID \pm 3 μ m; 280 μ m OD \pm 10 μ m) | 6720.0081 | 1 |
| Universal fitting (1/16", finger-tight, extra long thread; 4 fittings) | 6720.0015 | 1 |
| Monolithic column (200 μm ID 5 cm long) | 161409 | 2 |
| Monolithic trap column (ID 200 μm, 5 mm long, PS-DVB) | 163972 | 2 |
| Sample loop (20 µl, PEEK) | 6820.0018 | 1 |
| Cytochrome C Digest test sample | 161089 | 1 |
| Polypropylene vials with glass insert (set of 25 vials) | 6820.0027 | 1 |
| Caps for polypropylene vials (set of 25 caps) | 6820.0028 | 1 |
| End cap (set of 5 caps) | 164140 | 1 |

8.6.2 Biocompatible Devices

Nano/Nano HPLC:

| Description | Part No. | Qty in Kit |
|--|-----------|------------|
| Parallel LC Kit (biocompatible), Nano/Nano, including: | 6721.0054 | |
| Tool for preassembly of biocompatible fittings | 6000.0065 | 1 |
| Capillary from loading pump to autosampler (PEEK, 130 µm ID, 45 cm long), including appropriate fitting connections | 6721.0057 | 1 |
| Capillary from autosampler to FLM-3200 flow manager (PEEKsil, ID 75 µm, 50 cm long) | 6721.0038 | 1 |
| Capillary from pump to flow manager: | | |
| Long connection (system with detector) | 6037.2556 | 2 |
| Short connection (system without detector) | 6037.2554 | 2 |
| Capillary from nano flow splitter to flow manager switching valve (PEEKsil, 20 µm ID, 30 cm long) | 6721.0033 | 2 |
| Capillary to nano column in flow manager (PEEKsil, 20 µm ID, 50 cm long, including 1/16" PEEK fitting) | 6721.0024 | 1 |
| Capillary to nano column in flow manager (PEEKsil, 20 µm ID, 50 cm long) | 6721.0027 | 1 |
| PEEKsil bridge from flow manager to flow manager (75 μm ID, 30 cm long), including appropriate fitting connections (stainless steel) | 6721.0059 | 1 |
| PEEK sleeves for connection with fused silica capillaries (280 µm O.D., set of 5) | 6720.0064 | 2 |
| Universal fitting (PEEK, 1/16"), long hex nut with ferrule (set of 2 each) | 6721.0017 | 5 |
| PEEK sleeve (1/32", 300 μm ID, 3 cm; set of 6) | 6720.0079 | 1 |
| Fitting (1/32") for nano switching valve (set of 6) | 6720.0080 | 1 |

| Description | Part No. | Qty in Kit |
|--|-----------|------------|
| Teflon [®] tubing (500 μm ID, 100 cm long) for use as waste tubing | 6720.0077 | 3 |
| Fused silica tubing, 5m (75 μ m ID \pm 3 μ m; 280 μ m OD \pm 10 μ m) | 6720.0081 | 1 |
| Universal fitting (1/16", finger-tight, extra long thread; 4 fittings) | 6720.0015 | 1 |
| µ-Precolumn [™] holder (5 mm), including connecting tubing (30 μm ID) | 6720.0012 | 2 |
| Precolumn (trap column, C18 PepMap100, 300 µm ID, 5 mm long, set of 5) | 160454 | 2 |
| Separation column (C18 PepMap100, 75µm ID, 15 cm long) | 160321 | 2 |
| Sample loop (20 μl, PEEK) | 6821.0018 | 1 |
| Cytochrome C Digest test sample | 161089 | 1 |
| Polypropylene vials with glass insert (set of 25 vials) | 6820.0027 | 1 |
| Caps for polypropylene vials (set of 25 caps) | 6820.0028 | 1 |
| End cap (set of 5 caps) | 164140 | 1 |

Capillary/Capillary HPLC:

| Description | Part No. | Qty in Kit |
|--|-----------|------------|
| Parallel LC Kit (biocompatible), Cap/Cap, including: | 6721.0055 | |
| Tool for preassembly of biocompatible fittings | 6000.0065 | 1 |
| Capillary from pump to flow manager: | | |
| Long connection (system with detector) | 6037.2556 | 2 |
| Short connection (system without detector) | 6037.2554 | 2 |
| Capillary from loading pump to autosampler (PEEK, 130 µm ID, 45 cm long), including appropriate fitting connections | 6721.0057 | 1 |
| Capillary from autosampler to FLM-3200 flow manager (PEEKsil, 75 µm ID, 50 cm long) | 6721.0038 | 1 |
| PEEKsil bridge between the flow managers (75 μm ID, 30 cm long), including appropriate fitting connections (stainless steel) | 6721.0059 | 1 |
| Capillary from capillary flow splitter to switching valve (PEEKsil, 50 µm ID, 30 cm long) | 6721.0034 | 2 |
| Adaptor from capillary column to flow manager switching valve (PeekSil, 50 µm ID, 10 cm long, 1/16") | 6721.0040 | 2 |
| Capillary to capillary column in flow manager (PEEKsil, 20 µm ID, 50 cm long, including 1/16" PEEK fitting) | 6721.0025 | 1 |
| Capillary to capillary column in flow manager (PEEKsil, 20 µm ID, 50 cm long) | 6721.0028 | 1 |
| Universal fitting (1/16", finger-tight, extra long body; set of 2) | 6720.0072 | 2 |
| Teflon [®] tubing (250 μm ID, low pressure connection of 280 μm OD fused silica capillaries, set of 5) | 6720.0030 | 1 |
| Teflon [®] tubing (500 μm ID, 100 cm long) for use as waste tubing | 6720.0077 | 3 |
| Universal fitting (1/16", finger-tight, extra long thread; 4 fittings) | 6720.0015 | 5 ferrules |
| Universal fitting (PEEK, 1/16"), long hex nut with ferrule (set of 2 each) | 6721.0017 | 5 |
| μ-Precolumn [™] holder (5 mm), including connecting tubing (30 μm ID) | 6720.0012 | 2 |
| Precolumn (trap column, C18 PepMap100, 300 μm ID, 5 mm long, set of 5) | 160454 | 2 |
| Column (C18 PepMap100, 300μm ID, 15 cm long) | 160295 | 2 |
| Sample loop (20 µl, PEEK) | 6821.0018 | 1 |
| Cytochrome C Digest test sample | 161089 | 1 |
| Polypropylene vials with glass insert (set of 25 vials) | 6820.0027 | 1 |
| Caps for polypropylene vials (set of 25 caps) | 6820.0028 | 1 |
| End cap (set of 5 caps) | 164140 | 1 |

| Description | Part No. | Qty in Kit |
|--|-----------|------------|
| Parallel LC Kit (biocompatible), 2x Monolithic, including: | 6721.0056 | |
| Tool for preassembly of biocompatible fittings | 6000.0065 | 1 |
| Capillary from loading pump to autosampler (PEEK, 130 µm ID, 45 cm long), including appropriate fitting connections | 6721.0057 | 1 |
| Capillary from autosampler to FLM-3200 flow manager (PEEKsil, 75 µm ID, 50 cm long) | 6721.0038 | 1 |
| Capillary from pump to flow manager | | |
| Long connection (system with detector) | 6037.2556 | 2 |
| Short connection (system without detector) | 6037.2554 | 2 |
| Capillary from capillary flow splitter to switching valve (PEEKsil, 50 µm ID, 30 cm long) | 6721.0034 | 2 |
| PEEKsil bridge between the flow managers (ID 75 μm, 30 cm long), including appropriate fitting connections (stainless steel) | 6721.0059 | 1 |
| PEEK sleeves for connection with fused silica capillaries (280 µm O.D., set of 5) | 6720.0064 | 1 |
| Capillary to capillary column in flow manager (PEEKsil, 20 µm ID, 50 cm long, including 1/16" PEEK fitting) | 6721.0025 | 1 |
| Capillary to capillary column in flow manager (PEEKsil, 20 µm ID, 50 cm long) | 6721.0028 | 1 |
| Universal fitting (PEEK, 1/16"), long hex nut with ferrule (set of 2 each) | 6721.0017 | 5 |
| PEEK sleeve (1/32", ID 300 μm, 3 cm; set of 6) | 6720.0079 | 1 |
| PEEK sleeve (1/32", ID 400 μm, 3 cm; set of 6) | 6720.0041 | 1 |
| Fitting (1/32") for nano switching valve (set of 6) | 6720.0080 | 1 |
| Teflon [®] tubing (500 μm ID, 100 cm long) for use as waste tubing | 6720.0077 | 3 |
| Fused silica tubing, 5m (75 μ m ID \pm 3 μ m; 280 μ m OD \pm 10 μ m) | 6720.0081 | 1 |
| Universal fitting (1/16", finger-tight, extra long thread; 4 fittings) | 6720.0015 | 1 |
| Monolithic column (200 μm ID 5 cm long) | 161409 | 2 |
| Monolithic trap column (ID 200 μm, 5 mm long, PS-DVB) | 163972 | 2 |
| Sample loop (20 µl, PEEK) | 6821.0018 | 1 |
| Cytochrome C Digest test sample | 161089 | 1 |
| Polypropylene vials with glass insert (set of 25 vials) | 6820.0027 | 1 |
| Caps for polypropylene vials (set of 25 caps) | 6820.0028 | 1 |
| End cap (set of 5 caps) | 164140 | 1 |

| UltiMate 3000 - System Installation and Application | | | |
|---|--|--|--|
| | | | |
| | | | |
| | | | |
| | | | |
| | | | |
| | | | |
| | | | |
| | | | |
| | | | |
| | | | |
| | | | |
| | | | |
| | | | |
| | | | |
| | | | |
| | | | |
| | | | |
| | | | |
| | | | |
| | | | |
| | | | |
| | | | |
| | | | |
| | | | |
| | | | |
| | | | |
| | | | |
| | | | |
| | | | |
| | | | |
| | | | |
| | | | |
| | | | |
| | | | |
| | | | |
| | | | |
| | | | |
| | | | |
| | | | |
| | | | |
| | | | |
| | | | |
| | | | |
| | | | |
| | | | |

9 Technical Information

| Column flow range: | 50 nl/min to 2.5 mL/min using predefined splitter cartridges and active flow control |
|--|---|
| Flow control: | Electronically controlled and actively balanced flow splitter, independent from solvent composition and column backpressure |
| Gradient delay time: | Typically less than 1.5 min to splitter outlet at 200 nl/min |
| Sample thermostatting: | 4 to 45°C, max. 22 °C below ambient |
| Thermostatted column compartment and switching valves: | 5–70 °C (max. 15 °C below ambient) |
| UV detection: | Dedicated flow cells for nano-, capillary- and micro HPLC with 10 mm path lengths |
| User input and display: | Large LCDs, programmable to show important system parameters Soft keys for operation during installation and maintenance |
| System control: | All functions and parameters are software controlled via USB 1.1 connection |
| GLP features: | In Chromeleon: Automatic equipment qualification (AutoQ), System Wellness monitoring and recording of all system parameters in audit trail by Chromeleon software |
| MS software interfacing: | Via the appropriate add-ons to Chromeleon, single-point control for Analyst [®] (Applied Biosystems/ MDS Sciex) HyStar [™] (Bruker Daltonics) Xcalibur [®] (Thermo Electron Corporation) |

Technical information: May 2005.

All technical specifications are subject to change without notice.

For information about the technical specification on the individual system modules, refer to the *Operating Instructions* for the instrument.

| UltiMate 3000 - System Installation and Application | on |
|---|----|
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |

10 Index

| 2D LC Salt Plug28, 77 | $^{\prime}$ |
|-----------------------------|------------------------------|
| \mathbf{A} | Injection Command |
| Accessories 63 | |
| Applications 13 | |
| 2D LC Salt Plug28 | |
| Comprehensive 2D LC33 | |
| Direct Injection20 |) |
| Manual Injection16 | Leak |
| Parallel LC39 | Location |
| Preconcentration | |
| В | Maintenance |
| | Manual Injection 16.6/ |
| Basic Configuration | Menus 55 |
| C | O |
| Chromeleon | Operational Qualification 46 |
| Control45, 47, 49, 51, 53 | Operational Qualification |
| Error Messages60 |) P |
| Column Pressure |) |
| Comprehensive 2D LC | Parallel LC |
| Configurations | |
| Connections | Power Connection |
| Device Control. | |
| Fluidics 13 | Preparation10 |
| Power | |
| Rear Panel | R |
| D | Rear Panel (Connections) |
| Degasser53 | \mathbf{S} |
| Degasser Parameters | Safety Precautions2 |
| DegasserVacuum | Soft Keys |
| Direct Injection | Solvent Rack 53 |
| _ | Spare Parts 63 |
| ${f E}$ | Specification 87 |
| Error Messages (Chromeleon) | |
| F | T |
| _ | |
| Facility Requirements | |
| Flow Splitter | |
| Fluid Connections | |
| Front Panel | U |
| | Unpacking2 |
| | |

| UltiMate 3000: System Installation and Application | |
|--|--|
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |