

AAnalyst 700

Atomic Absorption Spectrometer

SYSTEM DESIGN

The AAnalyst™ 700 is a fully-integrated benchtop design atomic absorption spectrometer, incorporating all spectrometer and atomizer components in a single instrument.

OPTICAL SYSTEM

Real-time double-beam optical system with protectively-coated, front-surfaced, reflecting optics. Optical system sealed within protective cover.

Monochromator	Littrow design with motorized drive for automatic wavelength selection and peaking.
Wavelength range:	190 - 900 nm
Diffraction grating:	1800 lines/mm blazed at 236 nm and 597 nm
Grating area:	64 x 72 mm
Reciprocal linear dispersion:	1.6 nm/mm (nominal)
Focal length:	267 mm
Spectral bandwidths:	0.2, 0.7 and 2.0 nm, dual height; motorized slit drive for automatic slit selection

Detector	Wide-range segmented solid-state detector, including a built-in low-noise CMOS charge amplifier array.
-----------------	--

Automatic Lamp Selection	8-lamp holder with built-in power supplies for hollow cathode and electrodeless discharge lamps. Computer-controlled lamp selection and alignment via WinLab32™ software. Operating parameters are automatically set when using PerkinElmer Lumina™ hollow cathode lamps.
---------------------------------	---

FLAME SYSTEM

Gas Controls	Fully computer-controlled with oxidant and fuel monitoring. Keyboard-actuated remote ignition system with air-acetylene. Acetylene flow is automatically adjusted prior to the oxidant change when switching to or from nitrous oxide-acetylene operation. TotalFlow™ control of the oxidant and fuel gases for constant fuel:oxidant ratio.
---------------------	--

Safety Functions	Interlocks prevent ignition if the proper burner head, the nebulizer/end cap, or the burner drain system are not correctly installed; the liquid level in the drain vessel is incorrect; or gas pressures are too low. Interlocks automatically shut down gases if a flame is not detected, or if any of the other interlock functions are activated. Safe shutdown in the event of a power failure.
-------------------------	--

Burner System	Premix design that can be moved automatically into the sample compartment via software control and a motorized carriage. Alignment of the flame in the light beam is fully automatic and optimized, using a motorized burner mount for vertical and horizontal burner adjustments. The burner is equipped with a high-strength inert mixing chamber, angled to ensure proper drainage. Includes adjustable high-precision nebulizer and an all-titanium, 10-cm, single-slot burner head for air-acetylene operation.
----------------------	--

BACKGROUND CORRECTION

Built-in continuum source double-beam background correction using a high-intensity deuterium arc lamp.

GRAPHITE FURNACE

Built-in computer-controlled Heated Graphite Atomizer (HGA™) which can be moved automatically into the sample compartment and positioned via software control and a motorized carriage. An external protective gas stream around the graphite tube prevents the entrance of outside air to maximize tube life. An internal purge gas goes through the graphite tube to remove the volatilized matrix vapors during drying and thermal pretreatment. The two gas streams are independently computer-controlled. Pneumatic opening and closing of the furnace for easy tube change.

FURNACE FEATURES

Program Flexibility	Analytical programs with up to 12 steps can be set up. Each step can be programmed with the following parameters:
Temperature:	Ambient up to 3000 °C in steps of 10 °C
Ramp and Hold Time:	0 to 99 sec in steps of 1 sec
Internal Gas Flow:	0, 50, or 150 mL/min; can be switched over to another type of gas
Required Inert Gas:	Argon. Inlet pressure 300 kPa (3 bar) minimum. Maximum gas consumption is 1220 mL/min.
Water Coolant:	An optional circulatory cooling system is available for use with the AAnalyst 700. When operating the AAnalyst 700 without the circulatory cooling system, cooling water meeting the following specifications should be used: Sediment-free drinking water; 20-40 °C; flow rate not less than 2 L/min; pressure between 2.5 and 4.5 bar; pH between 6.5 and 7.5; hardness not greater than 14°dH or 100 ppm.

FURNACE AUTOSAMPLER

Sampler Table	Removable sample tray with 88 and 146 sampling positions and 1 overflow container for pipet washing. Minimum sample requirement: ca. 0.1 mL.
Dispensable Volume	Sample and Reagent: 1 – 99 µL, selectable in increments of 1 µL.
Electronics	The autosampler is powered from the spectrometer and is software-controlled.

DATA CONTROL SYSTEM

Complete PC control using WinLab32 software. Provides complete control of the instrument and its major accessories plus data handling and storage.

Data Handling	Instrument readings linear in absorbance (-0.500 A to +2.000 A), concentration or emission intensity. Integration times operator-selectable from 0.1 to 60 sec. in increments of 0.1 sec. Reading modes include time-averaged integration, non-averaged integration (peak area), and peak height measurement. Includes built-in statistics. Up to thirty (30) standards and a choice of calibration equations. Reslope of the analytical curve using a single operator-selected calibration standard. Built-in IEEE-488 interface for computer connection and use of optional accessories.
----------------------	--

DIMENSIONS

110 cm wide x 65 cm high x 70 cm deep (104 cm deep with furnace autosampler)

Weight	147 kg (without controller and cooling system)
---------------	--

POWER REQUIREMENTS

230 V (+/-10%), 50 Hz (+/-1%), 3600 VA (maximum); 240V (+/-10%), 60 Hz., 4000 VA (maximum)

Electrical Protection	As defined in EN 61010-1; Insulation Class I; Installation Category II; Pollution Degree 2.
------------------------------	---

TECHNICAL STANDARDS

Certification	Designed and tested to be in compliance with the legal requirements for laboratory instruments. The instrument was developed and produced in compliance with ISO 9001. The WinLab32 software provides required control parameters for GLP and instrument performance verification.
Safety Standards	EN 61010-1, EN 61010-2-061, CSA C22.2 No. 1010.1, CSA C22.2 No. 1010.2.061. The instruments bear the CE Mark and the CSA/NRTL Certification Mark.
EMC Standards	EN 61326, EN 55011, EN 61000-3-2, EN 61000-3-3

ENVIRONMENTAL REQUIREMENTS

Ambient Temperature	+15 °C to +35 °C
Relative Humidity	20% to 80%, non-condensing

COOLING SYSTEM

Optional self-priming recirculatory system with fan-assisted heat exchanger for constant cooling of the graphite furnace. Water temperature during operation approx. 36 °C; water flow 2.5 L/min.

Power Requirements	230 V (+/- 10%), 50/60 Hz (+/- 1%); approx. 140 VA
Dimensions	20 cm wide x 375 cm high x 50 cm deep
Weight	18 kg with coolant

PerkinElmer Life and
Analytical Sciences
710 Bridgeport Avenue
Shelton, CT 06484-4794 USA
Phone: (800) 762-4000 or
(+1) 203-925-4602
www.perkinelmer.com

