
PERKINELMER AANALYST 700 USER MANUAL

AAnalyst 700

Atomic Absorption Spectrometer

User's Guide

AAAnalyst 700

Atomic Absorption Spectrometer

User's Guide

PerkinElmer Bodenseewerk is registered for the design and manufacture of laboratory analytical equipment under the quality requirements of BS EN ISO 9001.

Certificate No. FM 22178

Notice

The information contained in this document is subject to change without notice.

Release history

Part Number	Release	Publication Date
0993-5255	1	January 1998
	2	September 1998
	3	June 1999
	4	March 2000
	5	September 2000

Trademarks

PerkinElmer is a trademark of PerkinElmer, Inc.

AAAnalyst, AA WinLab, Lumina, and HGA are trademarks of affiliates of PerkinElmer Instruments LLC.

Windows is a trademark and Microsoft is a registered trademark of Microsoft Corporation.

Registered names, trademarks, etc. used in this document, even when not specifically marked as such, are not to be considered unprotected by law.

Copyright information

This document contains proprietary information that is protected by copyright.

All rights are reserved. No part of this document may be reproduced in any form whatsoever or translated into any language without the prior written permission of PerkinElmer Instruments LLC or one of its affiliates.

Copyright ©1998–2000 by affiliates of PerkinElmer Instruments LLC

Printed in the Federal Republic of Germany

Technical Documentation

PerkinElmer Bodenseewerk

D-88647 Ueberlingen, Federal Republic of Germany

Contents

page

Chapter 1 Safety and Regulatory Information

Safety practices and conventions	1-3
Correct use of the instrument	1-4
Electrical safety	1-5
Operating conditions	1-7
Laboratory hygiene	1-8
Laboratory ventilation	1-8
Safe handling of chemicals	1-9
Waste disposal	1-10
UV radiation	1-11
High temperatures	1-12
Handling compressed gases	1-13
Safety practices for flame atomization	1-16
Safety practices for electrothermal atomization	1-22
Safety practices for FIAS and mercury hydride systems	1-24
References for laboratory safety practice	1-26
Regulatory information	1-27
IEC 1010 compliance	1-27
Electrical protection	1-27
Electromagnetic compatibility (EMC)	1-28
Symbols used on the instrument	1-29
Warning Markings on the instrument	1-30
Warning markings at the front of the spectrometer	1-30
Warning markings on the rear panel	1-32
Warning markings on the deuterium lamp cover	1-35

	<i>page</i>
Chapter 2 Spectrometer Installation	
Before you install the spectrometer	2-3
Moving and reinstalling the spectrometer	2-4
Preparing the spectrometer to be moved	2-4
Moving the spectrometer	2-5
Reinstalling the spectrometer	2-6
Connecting the burner gases	2-7
Connecting the furnace gases	2-10
To connect the drain system to the burner	2-15
To install the furnace autosampler	2-15
Connecting cooling sources	2-16
Electrical connections	2-22
Chapter 3 Preparing the System for Analyses	
Switching on the system	3-3
Installing lamps	3-7
Installing sample trays	3-11
Shutting down the system	3-12
Emergency shutdown	3-14
Chapter 4 Spectrometer: System Maintenance	
PerkinElmer Service	4-3
Maintenance checklist	4-4
Spectrometer maintenance procedures	4-5
Cleaning the windows of the optical system	4-5
Maintaining the electrical components	4-7
Changing the air filters	4-10
Replacing the deuterium lamp	4-11
Maintaining the acetylene supply	4-14
Chapter 5 Parts and Supplies	
Parts provided	5-3
Replacement parts	5-4
Obtaining supplies, replacement parts, and accessories	5-4

	<i>page</i>
Chapter 6 Laboratory Requirements	
Introduction	6-3
Operating conditions	6-3
Space requirements for the spectrometer system	6-5
Electrical requirements	6-8
Gas supply requirements	6-9
Burner gases	6-10
Furnace gases	6-14
Cooling the spectrometer system	6-16
Lamps	6-18
Computer and printer	6-18
Chapter 7 Spectrometer: System Description	
The AAnalyst 700 Spectrometer System	7-3
The optical system	7-6
Switches and connectors on the spectrometer	7-8
Technical data	7-14
Chapter 8 Translations of Warnings	
Translations of Warnings	8-3
Symbols used on the instrument	8-45
Index	
Customer Service	

Contents

Safety and Regulatory Information

1

Contents	page
Safety practices and conventions	1-3
Correct use of the instrument	1-4
Electrical safety	1-5
Operating conditions	1-7
Laboratory hygiene	1-8
Laboratory ventilation	1-8
Safe handling of chemicals	1-9
Waste disposal	1-10
UV radiation	1-11
High temperatures	1-12
Handling compressed gases	1-13
Safety practices for flame atomization	1-16
Safety practices for electrothermal atomization	1-22
Safety practices for FIAS and mercury hydride systems	1-24
References for laboratory safety practice	1-26
Regulatory information	1-27
IEC 1010 compliance	1-27
Electrical protection	1-27
Electromagnetic compatibility (EMC)	1-28
Symbols used on the instrument	1-29
Warning Markings on the instrument	1-30
Warning markings at the front of the spectrometer	1-30
Warning markings on the rear panel	1-32
Warning markings on the deuterium lamp cover	1-35

Safety practices and conventions

This instrument has been designed to comply with a wide variety of international standards governing the safety of laboratory equipment. In routine use, the instrument poses virtually no risk to you. If you take some simple, common-sense precautions, you can maintain the continued safe operation of the instrument.

The guide provided with the instrument contains information and warnings that must be followed by the user to ensure safe operation and to maintain the instrument in a safe condition. This advice is intended to supplement, not supersede, the normal safety code of behavior prevailing in the country of operation.

The information provided does not cover every safety procedure that should be practiced. Ultimately, maintenance of a safe laboratory environment is the responsibility of the user and the user's organization.

Possible hazards that could harm the user or result in damage to the instrument are clearly stated at appropriate places throughout this guide.

Any of the following safety conventions can be used throughout this guide:

*This symbol alerts you to situations that could result in **personal injury** to yourself or other persons.*

Details about these circumstances are in a box like this one.

*This symbol alerts you to the risk of **electric shock** that could result in **personal injury** to yourself or other persons.*

Details about these circumstances are in a box like this one.

*This symbol alerts you to the risk of **hot surfaces** that could cause **personal injury** to yourself or other persons.*

Details about these circumstances are in a box like this one.

*This symbol alerts you to the risk of **ultraviolet radiation** that could cause **eye damage** to yourself or other persons. Details about these circumstances are in a box like this one.*

Caution

*The term **CAUTION** alerts you to situations that could result in **serious damage to the instrument** or other equipment. Details about these circumstances are described in a message similar to this one.*

Correct use of the instrument

Before you install or use your instrument, and in order to get the best results, you should be familiar with all of the instruments in the system and know how to operate them. You should also be aware of the safety procedures in force in your laboratory, especially those concerning atomic spectrometry instruments. Consult the guides supplied with the instruments before you start.

If you use the instrument in a manner not specified in the guide, or if you use it for a purpose other than that intended, you may damage the instrument, or compromise your own, or someone else's, safety.

This instrument should only be operated by persons who are suitably qualified and have received adequate training.

Electrical safety

Warning: Unauthorized Adjustments and Servicing

W1.2

Only a PerkinElmer service engineer or similarly trained and authorized person should be permitted to service the instrument.

- *Do not attempt to make adjustments, replacements, repairs, or modifications to this instrument except as described in the documentation supplied with the instrument.*

Warning: Electrical Hazard

W2.1

Any interruption of the protective conductor inside or outside the instrument or disconnection of the protective conductor (earth/ground) terminal is likely to make the instrument dangerous.

Intentional interruption is prohibited.

Lethal voltages are present in the instrument

- *Even with the power switch off, line power voltages can still be present within the instrument.*
- *When the instrument is connected to line power, terminals may be live, and opening covers or removing parts (except those to which access can be gained without use of a tool) is likely to expose live parts.*
- *Capacitors inside the instrument may still be charged even if the instrument has been disconnected from all voltage sources.*

- To ensure satisfactory and safe operation of the instrument, it is essential that the line power cord is connected to a properly installed line power outlet that has a protective conductor (earth/ground).
- Do not attempt to make internal adjustments or replacements except as directed in the guide provided with the instrument.
- Do not operate the instrument with any covers or parts removed.
- Disconnect the instrument from all voltage sources before opening it for any adjustment, replacement, maintenance, or repair.
If, afterwards, the opened instrument must be operated for further adjustment, maintenance, or repair, this must only be done by a skilled person who is aware of the hazard involved.

- Use only fuses with the required current rating and of the specified type for replacement.

Do not use makeshift fuses or short-circuit the fuse holders.

- Whenever it is likely that the instrument is no longer electrically safe for use, make the instrument inoperative and secure it against any unauthorized or unintentional operation.

The instrument is likely to be electrically unsafe when it:

- Shows visible damage.
- Fails to perform the intended measurement.
- Has been subjected to prolonged storage under unfavorable conditions.
- Has been subjected to severe transport stresses.

Note: *The spectrometer is provided with an operational on/off switch and with a circuit breaker.*

The operational on/off switch turns the instrument on and off, but it does not disconnect the instrument from the line power supply.

The circuit breaker functions as a 'disconnecting device' and disconnects the instrument entirely from the line power supply.

Further information on these two switches is provided in 'Switches and connectors on the spectrometer' starting at page 7-6.

Operating conditions

**Warning: Explosive Atmosphere**

W1.3

This instrument is not designed for operation in an explosive atmosphere.

The instrument will operate correctly under the following conditions:

- Indoors.
- Ambient temperature +15 °C to +35 °C (+59 °F to +95 °F), with a maximum change not exceeding 2.8 °C (5 °F) per hour.
- Ambient relative humidity 20% to 80%, without condensation.
- Altitude in the range 0 m to 2000 m (sea level to 6500 feet).
- The location must be free of dust, smoke, and corrosive fumes.
- You can store the instrument safely under the following conditions:
 - Ambient temperature –20 °C to +60 °C (–4 °F to +140 °F).
 - Ambient relative humidity 20% to 80%, without condensation.
 - Altitude in the range 0 m to 12000 m (sea level to 40 000 feet).

Note: *When you remove the instrument from storage and before you put it into operation, allow it to stand for at least a day under the approved operating conditions.*

Laboratory hygiene

- Keep the work area scrupulously clean to avoid contaminating your samples and to maintain a safe working environment. Clean up spilled chemicals immediately and dispose of them properly.
- Do not allow smoking in the work area. Smoking is a source of significant contamination and also a potential route for ingesting harmful chemicals.
- Do not store, handle, or consume food in the work area.

Laboratory ventilation

Toxic combustion products, metal vapor, ozone, etc., can be generated by the system, depending on the type of analyses being performed.

- You must provide an efficient laboratory ventilation system to remove toxic products generated during instrument operation.

Safe handling of chemicals

Some chemicals used with the instrument may be hazardous or may become hazardous after completion of an analysis.

The responsible body¹ must take the necessary precautions to ensure that the surrounding workplace is safe and that instrument operators are not exposed to hazardous levels of toxic substances (chemical or biological) as defined in the applicable national, state, and local health and safety regulations and laws. Venting for fumes and disposal of wastes must be in accordance with all national, state, and local health and safety regulations and laws.

- Use, store, and dispose of chemicals in accordance with the manufacturer's recommendations and the applicable national, state, and/or local regulations.
- Wear appropriate eye protection at all times while handling chemicals. Depending on the types of chemicals you are handling, wear safety glasses with side shields, or goggles, or a full-face shield.
- Wear suitable protective clothing, including gloves if necessary, resistant to the chemicals you are handling.
- When preparing chemical solutions, always work in a fume hood that is suitable for the chemicals you are using.
- Perform sample preparation away from the instrument to minimize corrosion and contamination.
- Clean up spills immediately using the appropriate equipment and supplies, such as spill cleanup kits.
- Do not put open containers of solvent near the instrument.
- Store solvents in an approved cabinet (with the appropriate ventilation) away from the instrument.

1. Definitions from IEC 1010-1:

Responsible body: 'individual or group responsible for the use and maintenance of equipment, and for ensuring that operators are adequately trained.'

Operator: 'person operating equipment for its intended purpose.'

Waste disposal

Contents of waste containers

The materials that you collect in waste containers may include small amounts of the substances that were analyzed and other chemicals used in the analyses. If these materials are toxic, corrosive, or contain organics you may have to treat the collected effluent as hazardous waste. Refer to your local safety regulations for proper disposal procedures.

Disposing of defective lamps

Hollow cathode lamps and electrodeless discharge lamps contain small quantities of the lamp element in a very pure form. If a lamp containing toxic elements becomes defective or otherwise unusable, you must treat it as hazardous waste and dispose of it accordingly. A licensed company in the chemical waste disposal business can provide lamp disposal in accordance with environmental regulatory requirements. Please note that PerkinElmer does not take back defective lamps for disposal.

In addition, hollow cathode lamps and deuterium lamps are maintained under reduced pressure. Handle and dispose of them correctly to minimize the implosion risk.

For information, the quantities of analyte material contained in the cathodes of hollow cathode lamps are listed in the table below:

Quantity of material in the cathode	Single-element hollow cathode lamps	Multi-element hollow cathode lamps
less than 5 g	As, Au, B, Ba, Be, Ca, Dy, Er, Eu, Ga, Gd, Ge, Hf, Ho, In, Ir, K, La, Li, Mg, Na, Nd, Pd, Pr, Pt, Re, Rh, Ru, Sc, Se, Sm, Sn, Ta, Tb, Tm, Yb, Y	Na-K, Pt-Ru
5–10 g	Ag, Al, Bi, Cd, Co, Cr, Mn, Mo, Ni, Sb, Si, Sr, Te, V, W, Zn, Zr	Ca-Mg, Ca-Zn, Ag-Au, Sn-Te, Ca-Mg-Zn, Ca-Mg-Al
10–15 g	Cu, Fe, Hg, Nb, P, Pb, Ti, Tl	all other multi-element HCLs

For EDLs, the quantity of analyte material used is much less than the quantities used in HCLs. Typically, only several milligrams of material are placed in the EDL bulb. No PerkinElmer EDLs (System 1 or 2) contain more than 0.05 g of analyte material.

UV radiation

You should be aware of the health hazard presented by UV radiation.

- When the instrument is on, do not remove any covers unless specifically instructed to do so in the guide since otherwise you may be exposed to potentially hazardous UV radiation.
- Always wear UV-absorbing eye protection when viewing any of these sources:
 - The flame, especially the nitrous oxide/acetylene flame.
 - The graphite furnace when it is heated to incandescence.
 - Hollow cathode or electrodeless discharge lamps.
 - The deuterium background correction lamp.

High temperatures

Burner system

The flame can generate temperatures of up to 2800 °C.

- Do not touch the burner head until it has cooled to room temperature.

Graphite furnace

The HGA graphite furnace can generate temperatures of up to 3000 °C.

- Do not touch any part of the graphite furnace until it has cooled to room temperature.
- Do not attempt to inject a sample into the graphite tube when the furnace is hot. You may be subjected to fumes from the sudden vaporization of the sample, and may damage the pipet.

Quartz tube atomizer cell

The quartz tube atomizer cell (QTA-cell) used for the hydride-generation technique (option) can reach temperatures of up to 1000 °C.

- Do not touch any part of the heating mantle (electrically heated systems), or the cell holder (flame-heated systems), or the QTA-cell until they have cooled to room temperature.

Handling compressed gases

Note: The permanent installation of gas supplies is the responsibility of the user and should conform to local safety and building regulations.

Summary of gas hazards

Hazards associated with the different gases used in atomic absorption spectrometry are presented in the table below:

Gas	Suffocation	Explosion	Spontaneous Decomposition or Combustion
Air	—	—	—
Argon	x	—	—
Nitrogen	x	—	—
Acetylene (ethyne, C ₂ H ₂)	x	x	x
Nitrous Oxide (dinitrogen monoxide, N ₂ O)	x	—	x

Contact the gas supplier for a safety data sheet containing detailed information on the potential hazards associated with the gas.

Warning: Compressed Gases

High pressure gas cylinders can be dangerous if mishandled or misused.

- *Always handle gas cylinders with caution and observe your local regulations for the safe handling of gas cylinders.*

Identifying cylinders

- Legibly mark cylinders to clearly identify the contents and status (full, empty, etc.). Use the chemical name or commercially accepted name for the gas.

Storing cylinders

- Store cylinders in accordance with the applicable national, state, and/or local regulations and standards.
- When gas cylinders are stored in confined areas, such as a storage room, make sure that ventilation is adequate to prevent toxic or explosive accumulations of gas.
The storage room should be well ventilated and dry. This is particularly important in confined areas.
- Do not store cylinders near elevators, gangways, or in locations where heavy moving objects may strike or fall against them.
- Use and store cylinders away from exits and exit routes.
- Locate gas cylinders away from heat sources, including heat lamps.
Compressed gas cylinders should not be subjected to temperatures above 52 °C (125 °F).
- Do not allow ignition sources in the storage area and keep cylinders away from readily ignitable substances such as gasoline or waste, or combustibles in bulk, including oil.
- Store all gas cylinders only in a vertical position, with the valve cap in place, and fastened securely to an immovable bulkhead or a permanent wall.
- If you are storing cylinders outdoors, store them above ground on a suitable floor where they are protected against temperature extremes (including the direct rays of the sun).

Handling cylinders

- Move cylinders with a suitable hand truck after making sure that the valve cap is securely in place and that the cylinder is properly fastened to the hand truck.
- Use only approved regulators, tubing, and hose connectors.
- Arrange gas hoses where they will not be damaged or stepped on and where things will not be dropped on them.
- Do not attempt to refill gas cylinders.
- Check the condition of pipes, hoses, and connectors regularly, and replace any damaged parts.
- Perform periodic gas leak tests at all joints and seals of the gas system by applying an approved gas leak detection solution.

Safety practices for flame atomization

Safety interlocks

PerkinElmer provides a number of safety interlocks on the burner system to monitor gas pressure and check for the proper setup of the burner head, nebulizer, and drain system. In addition a flame sensor checks that the flame is burning.

- Do not attempt to defeat these interlocks; you may compromise your own, or someone else's, safety.

Safe use of burner gases

The seepage of fuel gas or fuel gas mixture from the burner system, the drain system, the gas control system, or the gas connections constitutes a serious fire hazard.

- Make sure that there are no breaks or leaks in any of these systems and that all the seals are correctly installed and in good condition.
- Regularly check for leaks at joints and seals using an approved leak test solution.
- When you perform any maintenance or installation procedures, follow the instructions in the guide exactly.
- Do not attempt to service the gas control system yourself.
A PerkinElmer service engineer, or similarly authorized and trained person, must perform the work.
- When you shut down the instrument, for example at the end of the working day, shut all the gas lines at the cylinder or regulator valves. Bleed the lines between the regulator and instrument to atmosphere before switching off the ventilation system.

Air supply

An air compressor is recommended for the air supply to the burner system. Cylinders of compressed air should only be used as an emergency or short-term solution.

- If you use cylinder air, make sure that it is compressed air and not oxygen-enriched air (e.g. breathing air); the use of oxygen-enriched air can cause a flashback of the flame.
- **Never** use oxygen as the oxidant since this can cause an explosion.

Nitrous oxide

(dinitrogen monoxide, N_2O)

- Tubing and fittings carrying nitrous oxide must be free of all oil, grease, and organic materials.
Spontaneous combustion may occur if nitrous oxide comes into contact with these materials.
- Use a double-stage or heated regulator for nitrous oxide to prevent freezing of the diaphragm and a loss of pressure regulation.
- Do not store nitrous oxide in close proximity to flammable gases, such as acetylene or hydrogen.

Acetylene

- Always use ‘Atomic Absorption Grade’ acetylene (ethyne) dissolved in acetone (propan-2-one). This is satisfactory for all PerkinElmer AA spectrometers.
- Do not store acetylene in close proximity to oxidizing gases, such as nitrous oxide.
- Use approved tubing and fittings for acetylene. Never use copper fittings since acetylene forms an explosive compound with copper.
- Prevent acetylene from coming into contact with copper, silver, mercury or gaseous chlorine.
- Periodically check for the presence of acetylene in the laboratory atmosphere, especially near the ceiling.
- For acetylene, set an outlet gauge pressure of between 90 kPa and 100 kPa (0.9–1.0 bar, 13–14.5 psig). Never allow the outlet gauge pressure to exceed 103 kPa (1.03 bar, 15 psig); acetylene can explode spontaneously above this pressure.
- Solvent Carryover.

If the acetylene cylinder pressure falls to below 600 kPa / 6.0 bar / 87 psig (at 20 °C), some of the solvent used to stabilize the acetylene (such as acetone) may be carried over into the burner. This could influence the characteristics of the burner with a resulting influence on the analytical results.

Change to a new acetylene cylinder when the cylinder pressure falls to below this value.

Safe operation of the flame

- Before you ignite the flame make sure that:
 - The laboratory fume ventilation system is operating;
 - The burner head is installed correctly;
 - The burner end cap is secured firmly;
 - The fuel and oxidant tubing fittings are properly connected;
 - The burner unit contains the necessary solvent-resistant components if you intend to analyze samples containing organic solvents;
 - The atomizer compartment door is closed;
 - The pH of the liquid in the drain vessel is greater than pH 10 if you intend to aspirate cyanide solutions.
- Do not leave the flame unattended. Always make sure that there is a fire extinguisher near the instrument.
- Never change the gas pressure or shut a gas valve while the flame is burning.
- Do not allow the burner head slot to become blocked. This can cause a flashback of the flame.
- Do not place open containers of flammable liquids and solvents near to the flame. Be especially careful with samples that contain highly volatile solvents.
- When you have analyzed samples containing organic solvents, flush all traces of solvent from the burner system.

The drain system

- Always place the drain vessel in a well ventilated place underneath the spectrometer, **in full view** while you are working with the spectrometer. This prevents the build-up of potentially hazardous gases, and allows you to see the liquid level.
- Never use a glass drain vessel.
- Do not kink, fold, or loop the drain tube.
- Never place the drain tube directly into a laboratory sink.
- Empty the drain vessel frequently; especially do not allow the drain vessel to fill with organic solvent.
- Thoroughly flush the drain system, as described in the guide, when you shut down the system.
- Check the condition of the drain tube regularly, especially if you use organic solvents. Replace the drain tube when it first shows signs of cracking or discoloration.

Flashbacks

A flashback is an explosion of the premixed gases in the spray chamber. The safety interlocks normally prevent conditions that could cause a flashback, but if a flashback should ever occur:

- Check that the burner head slot is clean.
- Make sure that you are using compressed air and not oxygen-enriched air.
- Check the flow spoiler or impact bead for damage, and replace any damaged parts.
- Carefully check and tighten all connections to, and seals on, the burner unit, and check the gas flow settings before you attempt to ignite the flame. Seepage of fuel and oxidant from the spray chamber may cause a flashback.
- Check that the burner head is seated correctly.

Hazards with Flame Atomization

Toxic combustion products:

Toxic combustion products can be generated by the system, depending on the type of analyses being performed.

- You must provide an efficient laboratory ventilation system to remove toxic products generated during instrument operation.

High temperatures:

The flame can generate temperatures of up to 2800 °C.

- To avoid serious burns, never touch the burner head until it has cooled to room temperature.

UV radiation:

The flame, especially the nitrous oxide/acetylene flame, emits ultraviolet radiation.

- Keep the atomizer compartment door closed when the flame is burning and never directly view the flame unless you are wearing UV-absorbing glasses.

Acetylides:

If you have aspirated high concentrations of copper, silver, or mercury solutions into an acetylene flame, unstable acetylides may have formed in the spray chamber. If permitted to dry, these compounds may explode.

- Aspirate solution continuously to prevent residues from drying.
- Thoroughly flush the spray chamber and drain system with water immediately after such an analysis.

Cyanide solutions:

If you intend to aspirate cyanide solutions, make sure that the pH of the liquid in the drain vessel is greater than pH 10; toxic hydrogen cyanide gas is formed when cyanides contact acidic solutions.

- Never allow solutions containing cyanides to mix with acidic solutions.

Safety practices for electrothermal atomization

Safety interlocks

PerkinElmer provides a number of safety interlocks on the graphite furnace (electrothermal atomization) system that monitor the gas pressure and the temperature of certain system components.

- Do not attempt to defeat these interlocks; you may compromise your own, or someone else's, safety.

Inert gas

The graphite furnace requires a supply of inert gas. PerkinElmer recommends **argon**.

Although nitrogen can be used, it is not recommended for the inert gas supply. At furnace temperatures above 2300 °C, nitrogen reacts with graphite to form toxic cyanogen gas (CN)₂.

Safety checks

- Before you start an analysis, make sure that:
 - The laboratory fume ventilation system is operating;
 - The inert gas supply is connected and set to the correct pressure;
 - The cooling water supply is turned on;
 - The autosampler waste bottle is empty.

Hazards with electrothermal atomization

Toxic products:

Toxic metal vapors, ozone, etc., can be generated by the system, depending on the type of analyses being performed.

- You must provide an efficient laboratory ventilation system to remove toxic products generated during instrument operation.
- If you use nitrogen as the inert gas, toxic cyanogen gas (CN)₂ will be produced when the furnace temperature is above 2300 °C.

High temperatures:

The HGA graphite furnace can generate temperatures of up to 3000 °C

- To avoid serious burns, never touch the graphite components until the furnace has cooled to room temperature.
- Do not attempt to inject a sample into the graphite tube when the furnace is hot. You may be subjected to fumes from the sudden vaporization of the sample, and may damage the pipet.

UV radiation:

The graphite furnace emits ultraviolet radiation when it is heated to incandescence.

- Never directly view the graphite furnace unless you are wearing UV-absorbing glasses.

Safety practices for FIAS and mercury hydride systems

These systems are offered as options for the spectrometer.

Safety checks

- Before you start an analysis, make sure that:
 - The laboratory fume ventilation system is operating;
 - The tubing in the fluid system is not constricted since this could lead to a pressure build-up in the system;
 - The drain tubing is inserted in the waste container and that this container is empty.
- If you are using a flame to heat the QTA-cell or if you intend to perform flame-FIAS analyses, refer to the information under *Safety Practices for Flame Atomization* before you ignite the flame.

Hazards with FIAS and mercury hydride systems

Toxic products:

Toxic metal vapors, etc., can be generated by the system, depending on the type of analyses being performed.

- You must provide an efficient laboratory ventilation system to remove toxic products generated during instrument operation.

High temperatures:

The QTA-cell used for the hydride-generation technique can reach temperatures of up to 1000 °C.

- Do not touch any part of the heating mantle (electrically heated systems), or the cell holder (flame-heated systems), or the QTA-cell until they have cooled to room temperature.

Hazardous chemicals:

Some of the chemicals required for analyses are corrosive and/or toxic.

- Refer to the information under *Safe Handling of Chemicals* and to the specific warnings in the documentation provided with the system before you start an analysis.

References for laboratory safety practice

- Bretherik, L., *Bretherik's Handbook of Reactive Chemical Hazards*, 4th ed., Butterworth & Co. Ltd., London, UK, 1990.
- *Safe Practices in Chemical Laboratories*, Royal Society of Chemistry, London, UK, 1989.
- *Hazards in the Chemical Laboratory*, 5th edition, Luxon, S.G., ed., Royal Society of Chemistry, London, UK, 1992.
- *CRC Handbook of Laboratory Safety*, 3rd edition, Furr, K., ed., The Chemical Rubber Co. Press, Florida, USA, 1990.
- *Prudent Practices for Handling Hazardous Chemicals in Laboratories*, National Research Council, National Academic Press, Washington D.C., USA, 1981.
- *Sax's Dangerous Properties of Industrial Materials*, 9th edition, Sax, N.I. and Lewis, R.J., eds., Wiley, New York, USA, 1998.
- Pohanish, R.P. and Greene, S.A., *Hazardous Materials Handbook*, Wiley, New York, USA, 1996.
- Compressed Gas Association, Inc., Arlington, VA 22202, USA, various publications.
- Data Sheets provided by chemical manufacturers, e.g.:
 - USA, Material Safety Data Sheets (MSDS),
 - FRG, DIN-Sicherheitsblätter,
 - UK, Hazard Data Sheets.

Regulatory information

IEC 1010 compliance

This instrument has been designed and tested in accordance with EN 61010-1:1993 (IEC 1010-1): *Safety requirements for electrical equipment for measurement, control, and laboratory use*, including Amendments 1 and 2 to this standard, and EN 61010-2-061:1996 (IEC 1010-2-061:1995): *Particular requirements for laboratory atomic spectrometers with thermal atomization and ionization*.

Electrical protection

Insulation

Class I as defined in IEC 1010-1.

Installation category

This instrument is able to withstand transient overvoltage according to Installation Category II as defined in IEC 1010-1 and IEC 664.

Pollution degree

This equipment will operate safely in environments that contain nonconductive foreign matter and condensation up to Pollution Degree 2 as defined in IEC 1010-1 and IEC 664.

Electromagnetic compatibility (EMC)

European Union (EMC Directives)

This instrument has been designed and tested to meet the requirements of the EU Directives 89/336/EEC and 92/31/EEC. It complies with the generic EMC standard EN 61326-1:1997 *Electrical equipment for measurement, control, and laboratory use – EMC requirements – General requirements*.

This instrument has passed the following EMC tests:

Emission of conducted and radiated noise	Electromagnetic compatibility
EN 55011:96	EN61 000-4-2:95
EN 61 000-3-2:95	EN61 000-4-3:95
EN 61 000-3-3:95	EN61 000-4-4:95
	EN61 000-4-5:95
	EN61 000-4-6:96

United States (FCC)

This instrument is classified as a digital device used exclusively as industrial, commercial, or medical equipment. It is exempt from the technical standards specified in Part 15 of FCC Rules and Regulations, based on Section 15.103 (c).

Japan (FCC)

This instrument has been tested and found to comply with the limits of a Class A digital device, pursuant to Part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference when the equipment is operated in a commercial environment. This equipment generates, uses, and can radiate radio frequency energy and, if not installed and used in accordance with the instruction manual, may cause harmful interference to radio communications. Operation of this equipment in a residential area is likely to cause harmful interference in which case the user will be required to correct the interference at his own expense. Changes or modifications not expressly approved by the manufacturer could void the user's authority to operate the equipment.

Symbols used on the instrument

Any of the following symbols can be used on the instrument:

This symbol indicates that there is a potential hazard and that the user must look in the guide for an explanation of the hazard and how to avoid it.

This symbol indicates the risk of electric shock.

This symbol indicates the risk of hot surfaces.

This symbol indicates the off position of the main power switch or circuit breaker. With the switch in this position, the instrument is disconnected entirely from the line power supply.

This symbol indicates the on position of the main power switch or circuit breaker. With the switch in this position, the instrument is connected to the line power supply but is not necessarily switched on and operational.

This symbol indicates the operational on/off switch.

This symbol indicates alternating current.

Warning Markings on the instrument

Warning markings at the front of the spectrometer

A

Risk of hot surfaces

Risque de surfaces chaudes

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Oberflächen

B

Warning: Moving Parts – Risk of Injury

Some moving parts of the instrument are accessible in normal operation.

Keep hands, clothing and other objects away from the moving parts of the instrument.

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
C

Nebulizer clamp must be lowered and locked over the nebulizer flange.

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B

Warning: Moving Parts – Risk of Injury


Some moving parts of the instrument are accessible in normal operation.
Keep hands, clothing and other objects away from the moving parts of the instrument.

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
D

Power outlets

Power outlets for a printer, computer, or monitor.

Only connect a printer, a computer, or a monitor to the respective outlet when the instrument requires 230 V, 50/60 Hz.


E

Cooling system outlet

Connect the cooling system to this power outlet only.

Do not connect the cooling system to any other supply.

Do not connect other equipment to this power outlet.

F

Burner Gases

Never set the outlet pressure to a value exceeding the maximum pressure (P_{max}).

Warning: Acetylene – Explosion Hazard

Acetylene can decompose explosively at pressures higher than 103 kPa / 1.03 bar / 15 psig.

Always make sure that the acetylene outlet gauge pressure is below this value.

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Warning: Flashback Hazard

PerkinElmer burner systems are designed for use with compressed air. The use of oxygen can cause an explosion in the burner system, and oxygen-enriched air can cause a flashback of the flame.

Never use oxygen or oxygen-enriched air with PerkinElmer burner systems. Use only compressed air.

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
Warning

Disconnect supply before servicing

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Warning


All field wiring must have insulation
suitable for at least 300 V

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Warning markings on the deuterium lamp cover

The warning markings below are visible when you open the top cover of the spectrometer.

High Temperatures – Risk of Burns

The lamp is hot.

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UV Radiation – Risk of Eye Damage

The lamp emits intense UV radiation which can damage your eyes.

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Electrical Hazard

The deuterium lamp operates at high voltage.

To avoid electrical shock, disconnect the line power cord before opening the cover.

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Spectrometer Installation

2

Contents	page
Before you install the spectrometer	2-3
Moving and reinstalling the spectrometer	2-4
Preparing the spectrometer to be moved	2-4
Moving the spectrometer	2-5
Reinstalling the spectrometer	2-6
Connecting the burner gases	2-7
Connecting the furnace gases	2-10
To connect the drain system to the burner	2-15
To install the furnace autosampler	2-15
Connecting cooling sources	2-16
Electrical connections	2-22

Before you install the spectrometer

When the spectrometer arrives

- Check the exterior of the shipping crate for obvious signs of damage. If the crate has been damaged, file a claim with the authorized carrier immediately and inform your nearest PerkinElmer office.
- Contact your local PerkinElmer office to arrange for the installation.

Before the service engineer arrives

- You must prepare a suitable site for the spectrometer and any accessories. See *Chapter 6, Laboratory Requirements*, and the guides for any accessories.
- You must make sure that the fume ventilation system and the necessary gas supply lines and regulators are installed.
- Although you should not perform the initial installation yourself, we recommend that you read the installation instructions to learn about the spectrometer and its requirements before you start using it. These instructions are also provided in case you need to alter any of the connections or move the instrument.

When the service engineer arrives

- The service engineer will remove the spectrometer from the crate and check that all parts have been delivered. If anything is missing or damaged, the service engineer will make sure that the respective parts are supplied.
- The service engineer will remove the shipping locks, perform pre-installation checks, install the spectrometer, and make sure that everything functions correctly.

Moving and reinstalling the spectrometer

We recommend that only a PerkinElmer service engineer or a PerkinElmer trained and authorized person moves and reinstalls the spectrometer. Should you attempt to do this yourself, PerkinElmer cannot be held responsible for any damage or malfunctioning of the spectrometer.

The instructions here are primarily intended for occasions where you may want to move the spectrometer **short** distances to a new location within the laboratory. Read the instructions carefully and be sure you understand them thoroughly before you attempt to move the spectrometer. Read the detailed information given in the relevant sections of this chapter before you reinstall the spectrometer. If you intend to transport the spectrometer over a longer distance, contact your local PerkinElmer office since internal shipping locks must be set via the service software, inside the spectrometer.

Preparing the spectrometer to be moved

1. Close the AA WinLab application, then shut down the Windows operating system.
2. Switch off the spectrometer, computer and all the accessories.
3. Disconnect all the electrical cables.
4. Remove the furnace autosampler from the spectrometer (*see 'Removing the autosampler from the spectrometer' in the Graphite Furnace User's Guide*).
5. Disconnect the burner gas supplies.
6. Disconnect the drain system from the burner.
7. Disconnect the furnace gas supplies.
8. Disconnect the cooling water supply.
9. Prepare the computer and other system components for transport as described in their own guides.
10. Attach the original handles securely. Only use the original handles and screws.

Note: *Make sure that the handles are properly secure before you move the spectrometer.*

Moving the spectrometer

Warning: Risk of injury through lifting

W1.4

The instrument weighs approximately 150 kg.

Improper lifting can cause painful and sometimes permanent back injury.

- *Use proper lifting techniques and always use the handles to lift or move the instrument.*

Caution: Risk of damage to the instrument

The spectrometer is a precision optical instrument. Move the spectrometer very carefully and protect it from vibrations and jolts. Do not tilt the spectrometer.

If the spectrometer must be tilted or mechanically lifted, contact your PerkinElmer service representative before moving the spectrometer. The shipping locks must be properly secured to prevent any damage occurring if the spectrometer is tilted or jolted.

Carefully lift the spectrometer onto a movable table to transport it to the new location.

Reinstalling the spectrometer

Reinstall the spectrometer by following exactly the procedures described in this chapter, as summarized in the steps below:

1. Prepare the new site for the spectrometer; see *Chapter 6, Laboratory Requirements*.

Make sure that there will be enough space at the rear and sides of the spectrometer system for air to circulate freely.

2. Remove the handles. Save the handles and screws.
Insert the plastic plugs provided (B015-3803) into the screw holes.
3. Connect the burner and furnace gas supplies.
4. Connect the drain system to the burner.
5. Install the furnace autosampler.
6. Connect the cooling water supply.
7. Make sure that the main power switches of all the instruments are in the off position, then connect the electrical cables.

Connecting the burner gases

The connectors for the burner gases are located on the rear panel of the spectrometer.

What you need to provide:

- Suitable supplies of acetylene, air, and nitrous oxide. The gases must conform with the requirements listed in *Chapter 6, Laboratory Requirements*.
- Suitable pressure regulators and on/off valves near to the spectrometer for the gas supplies.

The gas supply systems must comply with your local safety regulations. Refer also to the information given in *Chapter 1, Safety and Regulatory Information*.

Warning: Flashback Hazard

W5.9

PerkinElmer burner systems are designed for use with compressed air. The use of oxygen can cause an explosion in the burner system, and oxygen-enriched air can cause a flashback of the flame.

- *Never use oxygen or oxygen-enriched air with PerkinElmer burner systems. Use only compressed air.*

What is provided with the spectrometer:

- Gas hose assemblies for acetylene, air, and nitrous oxide.
The fittings on the gas hose assemblies are suitable for connection to most commercially available pressure regulators. These are R 3/8" standard threaded fittings. The fuel (acetylene) fitting has a left-hand thread, while the oxidant (air, nitrous oxide) fittings have right-hand threads.
If you require an adapter for your regulator, contact the gas supplier or the regulator manufacturer, or your local PerkinElmer office.

To connect the gases

Refer to Figure 2-1 when you connect the burner gases.

Air

1. Make sure that the air supply is turned off.
2. Use the black hose assembly to connect a source of clean, dry, and filtered air to the fitting labelled **Air/Luft** on the rear panel of the spectrometer.

Nitrous oxide

1. Make sure that the nitrous oxide supply is turned off.
2. Use the blue N_2O hose assembly to connect a source of clean, dry nitrous oxide to the fitting labelled **N₂O** on the rear panel of the spectrometer.

Acetylene

1. Make sure that the acetylene supply is turned off.
2. Use the red C_2H_2 hose assembly to connect a source of clean, dry acetylene to the fitting labelled **C₂H₂** on the rear panel of the spectrometer.


Figure 2-1. Connecting the Burner Gases

Connecting the furnace gases

The connectors for the furnace gases are located on the rear panel of the spectrometer.

Note: *The furnace will operate correctly only if there is gas pressure at both the Normal Gas and Special Gas inlets. If you are using an inert gas only you must connect the gas to both inlets using the Y-connector.*

What you need to provide:

- A suitable supply of argon (or nitrogen) and if necessary a suitable supply of a special gas. The gases must conform with the requirements listed in *Chapter 6, Laboratory Requirements*.
- Suitable pressure regulators and on/off valves near to the spectrometer for the inert and special gas supplies.
- A tube cutter for the special gas tube. This is an optional accessory; Part No. 0992-3183.

The gas supply systems must comply with your local safety regulations. Refer also to the information given under *Chapter 1, Safety and Regulatory Information*.

What is provided with the spectrometer:

- About 10 m of special gas tube, B050-8264.
- A Y-connector with automatic push-fit connectors, 0992-0179.
- Two screw connectors, R ¼" for the gas regulator. These are not supplied with all instrument versions.

To connect argon (or nitrogen) only

Refer to Figure 2-2 on the following page when you connect the furnace gas.

1. Make sure that the argon (or nitrogen) supply is turned off.
2. Use the special tube cutter to cut two, 10 cm lengths of tube from the special gas tube. The ends of the tubes must be cut squarely.
3. Connect the long gas tube to the argon (or nitrogen) pressure regulator. Use the R ¼" screw connector provided or use a connector or adapter suitable for the gas regulator. Make sure that all the connections are secure.
4. Pass the gas tube to the rear of the spectrometer. Make sure that it cannot be damaged or strained.
5. Push the two short tubes and the long tube into the connectors on the Y-connector as shown on the following page.
6. Push the free ends of the two short gas tubes into the **Normal Gas** and **Special Gas** connectors on the manifold as shown.
Push the tubes into the connectors until the tubes butt up to the end of the connector.

Note: *To disconnect the gas tube from a connector:
Carefully push the black locking ring of the connector into the connector and pull the tube out of the connector.*


Figure 2-2. Connecting Argon (or Nitrogen) only

To connect both argon and a special gas

Refer to Figure 2-3 on the following page when you connect the furnace gases.

1. Make sure that the gas supplies are turned off.
2. Use the special tube cutter to cut the special gas tube into two pieces for the two gas supplies. The ends of the tubes must be cut squarely.
3. Connect the gas tubes to the argon (or nitrogen) and special gas pressure regulators. Use the R ¼" screw connectors provided or use a connector or adapter suitable for the gas regulator. Make sure that all the connections are secure.
4. Pass the gas tubes to the rear of the spectrometer. Make sure that they cannot be damaged or strained.
5. Push the free ends of the gas supply tubes into the **Normal Gas** and **Special Gas** connectors on the manifold as shown.
Push the tubes into the connectors until the tubes butt up to the end of the connector.

Note: *To disconnect the gas tube from a connector:
Carefully push the black locking ring of the connector into the connector and pull the tube out of the connector.*


Figure 2-3. Connecting Argon (or Nitrogen) and a Special Gas

To connect the drain system to the burner

Procedures for connecting the drain system to the burner are described in *Maintaining the drain system* in the *Burner System User's Guide*.

To install the furnace autosampler

Procedures for demounting and reinstalling the furnace autosampler are described in *Installing the furnace autosampler* in the *Graphite Furnace User's Guide*.

Connecting cooling sources

Installing the cooling system

PerkinElmer offers a recirculating cooling system as an option for the spectrometer system. This cooling system includes a pressure-sensitive safety device that switches off the pump if the internal pressure in the circulatory system increases to above a safe value (see 'Cooling the furnace' in the *Graphite Furnace User's Guide*).

Connecting and disconnecting the water hoses

The inlet and outlet hoses are permanently attached to the cooling system. They have automatic connectors which are marked with the direction of the water flow.


1. Place the cooling system in a convenient position **on the floor** underneath the spectrometer. Do not place the cooling system on the same bench as the spectrometer since vibrations can adversely influence the spectrometer.

Make sure that:

- You can reach the on/off switch at the rear.
 - There is space at the rear and sides of the cooling system for air to circulate freely.
 - There is space to remove the stopper and fill the system.
 - You can see the water level indicator.
2. Pass the two hoses from the cooling system to the rear of the spectrometer. Make sure they are not strained.
 3. Push the hose connectors onto the correct fittings on the spectrometer as shown on the following page.
The connectors are spring loaded. Push them straight onto the fittings until they click into place.

Note: To disconnect the water hoses:

Carefully pull the knurled locking ring down slightly, to release the lock, and pull the connector off the fitting on the spectrometer.

Connecting the hoses



Figure 2-4. Connecting the Cooling System

Filling the cooling system

The cooling system is self-priming; air in the system is removed automatically as the cooling system operates. The coolant is a 1:10 solution of glycerol in water (1 part glycerol plus 9 parts water).

Note: *The cooling system contains 0.5 L glycerol when it leaves the factory. (See 'Cooling system maintenance procedures' in the Graphite Furnace User's Guide for the procedure to follow when you empty and refill the cooling system. Do not use antifreeze.)*

1. Remove the stopper from the top of the cooling system – twist and pull.
2. Fill the system to the **Max** mark with sediment-free drinking water – about 4.5 L.
Do not use deionized water since this can cause corrosion.
3. Replace the stopper securely.

Water level indicator

Connecting external cooling sources

If you intend to run the spectrometer system without the PerkinElmer cooling system, you must connect another suitable source of cooling, such as the main water supply or a laboratory thermostat. If you connect the spectrometer system to the main water supply, make sure that the water pressure is as given in Table 2-1.


Figure 2-5. Connecting the Main Water Supply

A water hose kit with pressure relief valve is provided with the spectrometer system which you must install in the water supply line to the spectrometer. This valve opens if the water pressure to the furnace system exceeds 4 bar. This feature eliminates the risk of damage to the furnace system through excess water pressure.

The cooling water must conform with the requirements listed in '*External cooling*' on page 6-16.

Caution**Water Pressure – Leaking hose fittings or burst water hoses**

If the water inlet pressure is too high, hose fittings may leak or hoses may even burst.

- *For reliable and safe operation, you must install the pressure relief valve kit B313-0124 in the water supply line.*

The hoses provided in the kit are fitted with automatic connectors that incorporate stop valves; these valves prevent water from running out when the hoses are disconnected. The hoses are marked with the direction of the water flow.

Table 2-1. Pressure for the Cooling Water Supply

Outlet Pressure		
$P_{\min}$	$P_{\max}$	
250	350	kPa
2.5	3.5	bar
36	50	psig

- Check the water pressure to your laboratory and make sure that the water inlet pressure is set to a value between the minimum value ($P_{\min}$) and the maximum value ($P_{\max}$).
 - Make sure that the outlet pressure is **never** set to a value exceeding the maximum pressure ($P_{\max}$).
1. The PerkinElmer service engineer will connect the pressure relief valve to the water tap (faucet).
 2. Make sure that the drain hose from the pressure relief valve hangs down into a suitable laboratory sink or drain.
 3. Connect the inlet hose (marked **WATER IN**) to the fitting at the bottom of the pressure relief valve.
 4. Connect the other end of the inlet hose to the **Furnace Cooling** fitting marked **In** on the spectrometer.
The connectors are spring loaded. Push them straight onto the fittings until they click into place.

5. Connect the outlet hose (marked **WATER OUT**) to the **Furnace Cooling** fitting marked **Out** on the spectrometer.
- Note:** To disconnect the water hoses:
Carefully pull the knurled locking ring down slightly, to release the lock, and pull the connector off the fitting on the spectrometer.*
6. We recommend that you fit a proprietary flow indicator in the outlet hose so that you can easily monitor that cooling water is flowing.
 7. Place the free end of the outlet hose in a suitable laboratory sink or drain. Make sure that the hose cannot slip out.


Figure 2-6. Connecting Cooling Water Hoses

Electrical connections

Warning: Electrical Hazard

W2.3

- To prevent potential injury to yourself and damage to the instrument, switch OFF all instruments in the system and disconnect them from the line power supply before you perform the procedures described below.

Read the following instructions very carefully before connecting any cables.

See *Chapter 6, Laboratory Requirements* for the specification of the line power supply.

Summary

1. Make sure that the spectrometer and the cooling system (where applicable) are switched off.
2. Switch the circuit breaker on the rear panel of the spectrometer to the off position (marked ○).
3. Switch off the computer and all other instruments that you will connect to the spectrometer or computer. See the guides for the individual instruments for the procedure to use.
4. Disconnect the line power cords of the spectrometer and all other instruments from the line power supply.
5. Connect the instrument control cables and power cords as shown on the following page.
6. When you have completed all connections, switch on the system (see *'Switching on the system'*, page 3-3).


Figure 2-7. Electrical Connections

Preparing the System for Analyses

3

Contents	page
Switching on the system	3-3
Installing lamps	3-7
Installing sample trays	3-11
Shutting down the system	3-12
Emergency shutdown	3-14

Switching on the system

1. Read the safety information at the front of this guide before you operate the system.
2. Make sure that the spectrometer system and other ancillary instruments are correctly installed and connected to their power supplies.
Make sure that the coolant hoses are connected correctly to the spectrometer. Where applicable, make sure that the cooling system is filled to the **Max** mark (see 'Filling the cooling system', page 2-18).
3. Switch on the fume ventilation system for the work area.
4. Turn on the burner gases and adjust the outlet gauge pressures to the recommended values (see Table 3-1 on page 3-5).
Always set the same outlet gauge pressure for air and nitrous oxide.
5. Turn on the inert (and special) gas for the furnace and adjust the outlet gauge pressures to the recommended values (see Table 3-2 on page 3-6).

Note: You must always turn on the air supply since air pressure is required for the internal pneumatic system. If you intend to work with only one atomizer, you only need to additionally turn on the gas supplies for that atomizer.

6. Turn on the cooling water to the system:
 1. If you are using the circulatory cooling system, switch it on by pressing the power switch located on the rear panel.
You can leave this switch permanently in the on position since power is supplied by the spectrometer.
 2. If you are using the main water supply, set a flowrate of approximately 2.5 L/min.
Note: If the water is too cold or flowing too fast, moisture may condense on the graphite contacts or on the furnace windows. Moisture on the graphite contacts leads to shorter lifetimes of the graphite tube and contacts, while moisture on the windows leads to a reduction in the radiant energy.
7. Switch the circuit breaker on the rear panel of the spectrometer to the on position (marked I).
You can leave this switch permanently in the on position during routine operation; you switch the spectrometer system on and off via the operational on/off switch located on the front panel.

8. Switch on the computer.
If your computer has a power management system, switch this off.
9. Install a Lumina hollow cathode lamp in any holder in the lamp compartment (*procedure: page 3-7*).
A hollow cathode lamp must be installed in at least one holder before you switch on the spectrometer.
10. Switch on the spectrometer system with the operational on/off switch located on the front panel.
11. If you need the diluter, autosampler, or flow injection system, switch on the respective systems.
12. If you have changed the configuration of the system, such as installed a different autosampler or added or removed any accessories, use the **Reconfigure** utility to reconfigure the software (*see the Online Help in AA WinLab*).
13. Start AA WinLab: In the AA WinLab group, double click on the AA WinLab icon.
14. Set up the various parts of the system:
 1. Install any further lamps that you require for the analyses (*procedure: page 3-7*).
 2. Set up the burner system
(*see 'Setting up the system for analyses using the flame technique' in the Burner System User's Guide*).
 3. Set up the graphite furnace system
(*see 'Setting up the system for analyses using the furnace technique' in the Graphite Furnace User's Guide*).

Procedures for shutting down the system are described on page 3-12.

In an emergency, you can shut down the system as described on page 3-14.

Table 3-1. Required Pressures for the Burner Gases

Gas	Outlet Gauge Pressure					
	kPa		bar		psig	
	$P_{\min}$	$P_{\max}$	$P_{\min}$	$P_{\max}$	$P_{\min}$	$P_{\max}$
Acetylene (C ₂ H ₂)	90	100	0.9	1.0	13.0	14.5
Air	450	500	4.5	5.0	65.0	72.5
Nitrous oxide (N ₂ O)	450	500	4.5	5.0	65.0	72.5

- Set the outlet gauge pressure for each gas to a value just below the maximum value ($P_{\max}$).
- Always set the same outlet gauge pressure for air and nitrous oxide.
- **Never** set the outlet gauge pressure to a value exceeding the maximum pressure ($P_{\max}$).
- Always turn on the air supply since air pressure is required for the spectrometer's pneumatic system.

**Warning: Acetylene – Explosion Hazard**

W5.7

Acetylene can decompose explosively at pressures higher than 103 kPa/
1.03 bar/15 psig.

- Always make sure that the acetylene outlet gauge pressure is below this value.

Table 3-2. Required Pressures for the Furnace Gases

Gas	Outlet Gauge Pressure					
	kPa		bar		psig	
	$P_{\min}$	$P_{\max}$	$P_{\min}$	$P_{\max}$	$P_{\min}$	$P_{\max}$
Argon	350	400	3.5	4.0	50.0	58.0
Nitrogen	350	400	3.5	4.0	50.0	58.0
Special	350	400	3.5	4.0	50.0	58.0

- Set the outlet gauge pressure for each gas to a value just below the maximum value ($P_{\max}$).
- **Never** set the outlet gauge pressure to a value exceeding the maximum pressure ($P_{\max}$).

Warning: Toxic Gases

W6.1

If you use nitrogen as the inert gas, toxic cyanogen gas (CN)₂ will be produced when the furnace temperature is above 2300 °C.


- *Make sure that the main laboratory fume ventilation system removes all gases and fumes emitted by the furnace.*

Installing lamps

Lamp compartment

The lamp compartment is located on the left-hand side of the spectrometer, beneath a hinged cover. The compartment is easily accessible from the front when the cover is open.

Lamp compartment cover

The spectrometer is designed to operate with PerkinElmer Lumina hollow cathode lamps (HCLs) and coded electrodeless discharge lamps (EDLs).

The lamp compartment has holders for up to eight lamps.

You can install Lumina HCLs in all eight holders. The holders incorporate the receptacles for connection of the Lumina lamps.

You can install EDLs in holders 1–4 only. The sockets for the power supplies to the EDLs are located on the left-hand side of the lamp compartment.



Figure 3-1. Lamp Holders in the Lamp Compartment

Installing and connecting lamps

HCLs

Note: A hollow cathode lamp must be installed in at least one holder before you switch on the spectrometer.


1. Open the lamp compartment cover by lifting it at the front and swinging it back.
2. Slide the Lumina lamp for the required element into the selected lamp holder and make sure that the integrated socket slides fully home over the receptacle on the lamp holder.
3. Close the lamp compartment cover.


EDLs

You must use the lamp holders numbered 1 through 4 for EDLs.

1. Open the lamp compartment cover by lifting it at the front and swinging it back.
2. If necessary, assemble the lamp by sliding the lamp driver into the lamp sleeve for the required element. Depress the locking pin as the sleeve slides into place and make sure that the locking pin extends through the locating hole in the sleeve.

3. Slide the lamp carefully into the selected lamp holder until it butts up against the end-stop inside the holder.
4. Insert the coding plug over the receptacle on the lamp holder.
5. Insert the lamp connector into the correct EDL socket at the left-hand side of the lamp compartment. Use the receptacle that has the same number as the lamp holder.
6. Close the lamp compartment cover.

Removing lamps

Do not unplug and remove lamps while they are switched on.

1. On the **Toolbar**, click on **Lamps**. The Align Lamps window appears.
2. Make sure that the lamps you want to remove are switched off:
The **On** buttons must *not* be green. If necessary, click on the respective **On** buttons.
3. Close the Align Lamps window.
4. Unplug and remove the lamps. For EDLs, remove the coding plugs as well.

Tips for installing lamps

- Use only PerkinElmer lamps. You can use both hollow cathode and electrodeless discharge lamps.
- Do not touch the front window of the lamp; perspiration or other contamination can reduce the intensity of the radiation.
- If you install more than one lamp containing a particular element, for example, a single-element copper lamp and a multi-element lamp containing copper, the system will use the lamp in the holder with the lower number when you perform determinations of that element. Make sure that you place the preferred lamp in the holder with the lower number.
- With EDLs, connect each lamp plug to the socket that has the same number as the holder of the lamp.
- With EDLs, make sure that you use the correct coding plug and connect it to the receptacle on the lamp holder containing the EDL.
- If you use the wrong coding plug or connect the coding plug to the wrong receptacle, you will cause the spectrometer to malfunction and may cause irreparable damage to the lamp.

Installing sample trays

A sample tray can be inserted into one of three locations in front of the spectrometer – on the left-hand side, in front of the atomizer compartment, and on the right-hand side. You can leave a tray permanently installed on the left-hand side if you wish. You cannot install a tray in front of the atomizer compartment or on the right-hand side as long as the furnace autosampler is installed. When the autosampler is swivelled to the standby position at the right, you can install a tray in front of the atomizer compartment.

The sample trays are robust and can support an accessory, such as a diluter, or a flame autosampler, or a flow injection system.

To install the tray

1. Hold the tray with both hands and position it in front of the location into which you wish to insert it.
2. Insert the tray support posts into the guides at the base of the instrument and slide in the tray until it hits the stop. Then gently lift the front of the tray and slide it finally into place. You can feel it lock into place.

To remove the tray

1. Gently lift the tray at the front and then pull it back to slide it out of the guides.

Shutting down the system

Note: In an emergency, switch off the spectrometer using the operational on/off switch on the front panel.

The flame is automatically extinguished and the gases shut down in a safe sequence.

See page 3-14 for the procedure to follow after an emergency shutdown.

You must shut down each atomizer correctly before you shut down the spectrometer.

To shut down the flame system

1. With the flame still burning, aspirate the correct rinsing solutions to rinse the burner system (see 'Extinguishing the flame' in the Burner System User's Guide for detailed information).

Either aspirate the solutions manually or use the flame autosampler.

2. Extinguish the flame and bleed the gas lines:
 1. In the **Flame Control** window, click on the **Flame** icon.
The system uses a predetermined shutdown sequence to extinguish the flame safely.
 2. Shut down the gases to the spectrometer at source.
 3. In the **Flame Control** window, click on **Bleed Gases** to depressurize the gas lines.

To shut down the furnace system

1. Flush the furnace autosampler rinsing system (see 'Autosampler maintenance procedures' in the Graphite Furnace User's Guide for detailed information).
2. Shut down the furnace gases to the spectrometer at source.

To shut down the spectrometer

1. Exit AA WinLab:
In the **File** menu, click on **Exit**.
2. Switch off the spectrometer at the operational on/off switch on the front panel.
3. Switch off any accessories.
4. If you do not want to use the computer for other tasks, shut down the computer and printer as described in their user's guides.
5. Empty the burner drain vessel (*see 'Emptying the drain vessel' in the Burner System User's Guide*).
6. Empty the furnace autosampler waste bottle.

Note: *Dispose of hazardous or corrosive solutions properly and refer to your local safety regulations for proper disposal procedures.*

Emergency shutdown

If a situation arises where you have to shut down the system quickly and you do not have time to carry out the correct procedure, you can use the procedure described below.

Note: *This procedure is for use in an emergency only. Do not use it as a routine shutdown procedure.*

1. Switch off the spectrometer at the operational on/off switch on the front panel.
2. Shut down the gases to the spectrometer at source.

When you return to the spectrometer, proceed as follows to shut down the system properly:

3. If you were previously running concentrated salt solutions, organic solvents, or solutions containing cyanide, copper, silver, or mercury:
 1. Remove the burner head (*see the Burner System User's Guide*).
 2. Pour a large quantity of deionized water slowly through the neck of the spray chamber to thoroughly flush the chamber and the drain system.
4. Check the burner head.

If there are deposits along the slot, clean the burner head (*see the Burner System User's Guide*).
5. Reinstall the burner head.
6. Switch on and set up the spectrometer system.
7. Ignite the flame.
8. Extinguish the flame properly.

Spectrometer: System Maintenance

4

Contents	page
PerkinElmer Service	4-3
Maintenance checklist	4-4
Spectrometer maintenance procedures	4-5
Cleaning the windows of the optical system	4-5
Maintaining the electrical components	4-7
Changing the air filters	4-10
Replacing the deuterium lamp	4-11
Maintaining the acetylene supply	4-14

PerkinElmer Service

**Warning: Unauthorized Adjustments and Servicing**

W1.2

Only a PerkinElmer service engineer or similarly trained and authorized person should be permitted to service the instrument.

- *Do not attempt to make adjustments, replacements, repairs, or modifications to this instrument except as described in the documentation supplied with the instrument.*

If the spectrometer system does not function correctly...

- Switch off the spectrometer system and make sure that it cannot be operated.
- Contact your local PerkinElmer office.

You should only perform the maintenance procedures described in this guide.

For any other maintenance or service contact your local PerkinElmer office to arrange for a service engineer to visit.

In particular you should only allow a PerkinElmer service engineer or similarly trained and authorized person to perform any work on:

- electrical components inside the instrument;
- the gas connections and gas system pneumatics;
- the monochromator and all other optical system components.

Before the service engineer arrives:

1. Make sure that the instrument and work area are clean.
2. Empty corrosive solutions from the burner drain vessel.
Then add water to the drain trap to activate the function.
3. Empty the furnace autosampler waste bottle.

Troubleshooting

See the *Troubleshooting Guide* supplied with the AA WinLab operating system. This guide includes information on problems that can arise with both the software and the spectrometer system.

Maintenance checklist

Maintenance checklists for the spectrometer system are presented below. Maintenance checklists for the burner system are presented in the *Burner System User's Guide* and for the graphite furnace, cooling system, and the furnace autosampler in the *Graphite Furnace User's Guide*.

General preventative maintenance procedures	
Instrument and work area	For safety reasons, and to avoid contaminating new samples, you should make sure that the instrument and the work area are always absolutely clean. This is especially important when you are working with toxic substances and when you are measuring trace amounts of any elements. Wipe up spills immediately before they can cause further contamination or damage. Before you leave the instrument standing idle for any period of time, make sure that all the components are clean.
Fume ventilation	Make sure that the main fume ventilation system functions correctly. Always switch it on when you are using the instrument.
Environmental conditions	Make sure that the environmental conditions are as described in <i>Chapter 6, Laboratory Requirements</i> .
Acetylene cylinder	Always change to a new acetylene cylinder when the cylinder pressure falls to below 600 kPa (6 bar; 87 psig) (<i>more information: page 4-14</i>).

Spectrometer, maintenance checklist	
External surfaces	Wipe over the surfaces with a lint-free cloth moistened with a dilute solution of laboratory detergent. The surfaces are resistant to dilute acids and alkalis, and to a lesser extent to strong acids and alkalis and organic solvents.
Air filters	Fit new filters yearly or more often in a dusty environment (<i>procedure: page 4-10</i>).
Atomizer compartment windows, lamp end windows	These windows do not normally require cleaning but if they become marked, e.g., with fingerprints, you can clean them (<i>procedure: page 4-5</i>).

Spectrometer maintenance procedures

Procedures described

- Cleaning the windows.
- Changing the air filters.
- Resetting the circuit breaker and changing the fuses.
- Replacing the deuterium lamp.

Cleaning the windows of the optical system

Caution

Risk of damage to optical components

- *Never touch or clean the surfaces of the grating or mirrors. If mirror surfaces are scratched they must be resurfaced. Do not allow cleaning fluids to come into contact with any of the reflecting surfaces.*
-

Most of the optical components are under the main cover of the spectrometer and are protected by the black optics covers.

Usually these components require no maintenance, provided the laboratory fume ventilation system is working correctly and the environmental conditions are as described in *Chapter 6, Laboratory Requirements*.

The only optical parts that are exposed are:

- The atomizer compartment windows.
- The windows of the lamps.
- The furnace windows. Cleaning the furnace windows is described under *System Maintenance* in the *Graphite Furnace User's Guide*.

These windows are important parts of the whole optical system. Handle them carefully, as you would any sensitive optical component. Do not directly touch the windows.

If mirror surfaces need cleaning, this must be done by a PerkinElmer service engineer or a similarly authorized and trained person who knows the special techniques required.

To clean the windows:

This must be done only by a skilled person who knows the special techniques required.

Note: *The more you rub the surface of the windows, the more chance there is of scratching them. You must change scratched windows.*

W3.1

Warning: UV Radiation – Risk of Eye Damage

The lamps may emit UV radiation which can damage your eyes.

- *Do not gaze into a lighted lamp.*
- *Always wear UV-absorbing safety glasses when looking at the radiation from the lamps.*

1. Blow dust off optical surfaces with clean dry air before you attempt to clean them. Do not rub the surfaces.
2. Clean the window surfaces with a soft lint-free cloth moistened with spectroscopic grade alcohol. The window must dry quickly to avoid leaving drying marks.

Keep the cleaning solutions and water away from other parts of the instrument.

Maintaining the electrical components

Warning: Electrical Hazard

W2.3

- *To prevent potential injury to yourself and damage to the instrument, switch OFF all instruments in the system and disconnect them from the line power supply before you perform the procedures described below.*

Parts that may be serviced by the user

The only parts that you should attempt to maintain are those described below. The electronic components in the instrument do not contain any parts which you can maintain yourself. If there is a problem, contact your local PerkinElmer office.

If any fuses blow frequently or if the automatic circuit breaker switch is activated frequently, contact your local PerkinElmer office for advice.

Resetting the automatic circuit breaker

The spectrometer has an automatic circuit breaker which prevents the system from drawing too much current.

When the circuit breaker is activated, the switch flips to the off position (marked O). This can happen while an analysis is running.

To reset the switch, wait a few seconds then put the switch to the on position (marked I).

If this switch is activated frequently, contact your local PerkinElmer office for advice.

Changing the fuses

Warning: Electrical Hazard

W2.9

- *Use only fuses with the correct current rating and of the specified type for replacement.*
- *Do not use makeshift fuses and do not short-circuit the fuse holders.*

1. Switch off the spectrometer and all other instruments in the system as described in the relevant guides.
2. Switch the circuit breaker on the rear panel of the spectrometer to the off position (marked O).
3. Disconnect the spectrometer's line power cord from the line power supply.
4. Disconnect all the line power cords from the sockets at the rear of the spectrometer.
5. To change the fuses numbered F1, F2, F3, F4:
 1. Insert the tip of a suitable screwdriver in the slot of the respective fuse holder and unscrew the fuse holder.
 2. Remove the blown fuse, then insert a new one of the correct type and rating. The correct rating is shown on the label of the fuse on the rear panel of the spectrometer.
 3. Screw the fuse holder into the correct socket.
6. Reconnect the line power cords correctly.


Figure 4-1. Automatic Circuit Breaker and Fuses on the rear panel


Changing the air filters

To make sure that the electrical components in the spectrometer are properly cooled, the air filters must allow sufficient air into the instrument. Check the filters regularly and change them if they are dirty.

You should not need to change the filters more than once every year if the environmental conditions are as described under *Chapter 6, Laboratory Requirements*.

Replacement filters are available; see *Chapter 5, Parts and Supplies* for the part number.

1. Carefully pull off the snap-on plastic grid that holds the air filter in place.
2. Replace the filter and push the grid back into place.

Replacing the deuterium lamp

The deuterium lamp for the background corrector has a limited lifetime. After a prolonged period of operation the lamp will exhibit decreasing radiant intensity and will ultimately fail.

The lifetime of the deuterium lamp will be prolonged if you do not switch it on and off frequently.

You can replace the deuterium lamp assembly yourself. The replacement lamp is fitted in a prealigned mount.


Figure 4-2. Location of Deuterium Lamp inside the Spectrometer

Warning: Electrical Hazard

W2.3

- To prevent potential injury to yourself and damage to the instrument, switch OFF all instruments in the system and disconnect them from the line power supply before you perform the procedures described below.

1. Switch off the spectrometer at the operational on/off switch on the front panel.
2. Switch the circuit breaker on the rear panel of the spectrometer to the off position (marked O).
3. Disconnect the spectrometer from the line power supply.


Figure 4-3. Deuterium Lamp Cover and Bracket Assembly

4. Open the top cover of the spectrometer:
 1. Remove the two securing screws; one on the left and one on the right.
 2. Swing back the cover to the open position.
5. Remove the deuterium lamp assembly:
 1. Undo the two captive screws on the deuterium lamp cover.
 2. Carefully lift up the deuterium lamp cover.
The deuterium lamp assembly is attached to a bracket underneath the cover.
 3. Carefully place the complete assembly on top of the black optics cover.
 4. Pull out the deuterium lamp plug from the socket in the lamp compartment.
Do not pull on the deuterium lamp leads.
6. Slacken the thumbscrew on the bracket and remove the old lamp assembly.
Save the thumbscrew for use with the new lamp assembly.

7. Unpack the new lamp assembly, taking care to hold it only by the metal mount to prevent fingermarks on the lamp envelope.
8. Slip the slot at the base of the lamp mount over the stud on the lamp bracket and secure the lamp assembly with the thumbscrew saved in step 6.
9. Wipe the lamp envelope with a soft cloth moistened with alcohol to remove dirt, since this would otherwise be burned in when the lamp is hot.
10. Insert the deuterium lamp plug into the socket in the lamp compartment. The plug is asymmetric and can only be inserted in one direction.
11. Reinstall the lamp assembly:
 1. Carefully lower the deuterium lamp assembly into the lamp compartment until the lamp cover fits in place.
 2. Screw in the two captive screws on the lamp cover.
12. Close the top cover:
 1. Swing down the cover to the closed position.
 2. Screw in the two securing screws; one on the left and one on the right.

Since the deuterium lamp is fitted with as prealigned mount, you do not normally need to align the lamp.

Note: Do not attempt to make adjustments to the lamp bracket assembly. If there is indication that the alignment is not correct (e.g., low energy for the BG bar graph in the Align Lamps or Continuous Graphics windows), please contact your local PerkinElmer office.

Maintaining the acetylene supply

Warning: Acetylene – Explosion Hazard

W5.7

Acetylene can decompose explosively at pressures higher than 103 kPa / 1.03 bar / 15 psig.

- *Always make sure that the acetylene outlet gauge pressure is below this value.*
- You should depressurize both the fuel and oxidant gas hoses at the end of the working day or if you are not going to use the instrument for an extended period. Close the cylinder valves to avoid the possibility of gas hoses being subjected to the full cylinder pressure in the event that a pressure regulator should fail.
- Acetylene is supplied dissolved in a solvent such as acetone, and a small amount of solvent carryover with the acetylene is normal. However, as tank pressure falls, the relative amount of solvent entering the gas stream increases and can give erratic results, particularly for elements such as calcium, tin, chromium, molybdenum and others whose sensitivity is highly dependent on the fuel:oxidant ratio. The pressure drop is not linear with gas removal since the acetylene is dissolved in a solvent, and a pressure of 600 kPa (6 bar; 85 psig) at ambient temperature indicates that the cylinder is nearly empty. For this reason, replace acetylene cylinders when the cylinder pressure drops to about 600 kPa (6 bar; 85 psig).
- Always store and use acetylene cylinders vertically, rather than horizontally, to prevent liquid acetone from reaching the cylinder valve.
- The practice of ‘cracking the valve’ of an acetylene tank (that is, opening the valve slightly for a very short period previous to attaching the regulator) is not recommended. While such an action will clear the valve opening of dust or dirt particles and may remove acetone from the cylinder valve, it is a potentially hazardous practice and one which should never be attempted in the presence of an open flame, sparks or other possible sources of ignition.

Parts and Supplies

5

Contents	page
Parts provided	5-3
Replacement parts	5-4
Obtaining supplies, replacement parts, and accessories	5-4

Parts provided

PerkinElmer reserves the right to alter the schedule of parts provided without prior notice.

Qty.	Item	Part No.
1	AAAnalyst 700 atomic absorption spectrometer.	B312-0082
1	Read cable.	B312-1302
3	Special line cords	B310-0003
1 pack	Replacement air filters	B050-2706
4	Plugs, for handle screw holes	B015-3803
1 pack	Replacement fuses	B312-1303
1	AS-800 furnace autosampler	B300-1500
1	Swing-arm mount for autosampler.	B312-0922
	<i>Further parts provided with the furnace autosampler are listed under 'Parts and Supplies' in the Graphite Furnace User's Guide.</i>	
1	Fuel hose assembly (red)	B015-4984
1	Air hose assembly (black)	B015-4987
1	Nitrous oxide hose assembly (blue)	B015-4990
10 m	Special gas tube (FLEXO)	B050-8264
1	Y-connector.	0992-0179
1	Hose kit.	B311-0036
1	Drain vessel	B312-1306
1	Sample tube	B015-6272
5	Corkpreen gaskets for burner end cap	0047-2014
1	O-ring for burner end cap	0990-2147
1	O-ring for the spray chamber	0990-2219
1	Burner head cleaning tool.	0303-1573
1 pack	Graphite components	B300-0408
5	Graphite tubes with integrated platforms	B300-0407
1	Tube alignment tool	B008-0112
1 set	Tools for changing graphite contacts	B312-1301
1	Documentation package	B312-0299

Replacement parts

Obtaining supplies, replacement parts, and accessories

Supplies, replacement parts, and accessories can be ordered directly from PerkinElmer, using the part numbers quoted in the guides provided with the instrument.

e-ssentials

PerkinElmer Instruments on-line catalog, your single source for quality supplies and accessories.

<http://essentials.perkinelmer.com/>

PerkinElmer's catalog service offers a full selection of high-quality supplies. To place an order for supplies and many replacement parts, request a free catalog, or ask for information:

- If you are located *within the U.S.*, call toll free **1-800-762-4002**, 8 a.m. to 8 p.m. EST. Your order will be shipped promptly, usually within 24 hours.
- If you are located *outside of the U.S.*, call your local PerkinElmer sales or service office.

What is listed in this section

This section lists the most important replacement parts for your spectrometer system. Parts listed under *Parts provided* are not repeated in this section. Information about other parts and accessories for the system is given in the current PerkinElmer brochures and price lists.

Background corrector

Description	Part No.
Deuterium lamp	B014-8615

Laboratory Requirements

6

Contents	page
Introduction	6-3
Operating conditions	6-3
Space requirements for the spectrometer system	6-5
Electrical requirements	6-8
Gas supply requirements	6-9
Burner gases	6-10
Furnace gases	6-14
Cooling the spectrometer system	6-16
Lamps	6-18
Computer and printer	6-18

Introduction

The AAnalyst 700 spectrometer has some special requirements that you must consider when setting up the laboratory. The laboratory must fulfil these requirements, described in this chapter, before you can operate the spectrometer correctly.

Brochure no. L-301 – *Preparing Your Laboratory for PerkinElmer Atomic Spectroscopy Instruments* – also contains useful information. This brochure is available from PerkinElmer.

Operating conditions

**Warning: Explosive Atmosphere**

W1.3

This instrument is not designed for operation in an explosive atmosphere.

The instrument will operate correctly under the following conditions:

- Indoors.
- Ambient temperature +15 °C to +35 °C (+59 °F to +95 °F), with a maximum change not exceeding 2.8 °C (5 °F) per hour.
- Ambient relative humidity 20% to 80%, without condensation.
- Altitude in the range 0 m to 2000 m (sea level to 6500 feet).
- The location must be free of dust, smoke, and corrosive fumes.

Note: *When you remove the instrument from storage and before you put it into operation, allow it to stand for at least a day under the approved operating conditions.*

Bench and location requirements

The location you choose for the spectrometer must fulfil the criteria listed below:

- Place the spectrometer on a sturdy bench or trolley (cart) capable of sustaining the weight. The work surface must be level, flat, clean, dry, and free from vibration. It should be strong and stiff enough to bear the weight of the spectrometer system over a long period without warping or bending. A movable table is very convenient, but make sure that the wheels can be securely locked.
- Position the bench or trolley out of direct sunlight and away from radiators and heaters.
- Position the system near to the electricity and gas supply points, and close to a source of cooling water.
- Position the system so that the atomizer compartment is located beneath a suitable fume ventilation system.
- Make sure that there is space at the rear and sides of the system for air to circulate freely.
- Make sure that there is sufficient room for you to work comfortably with the instruments, and that you can reach the connections at the rear of the instruments and the electricity and gas supply points.

Fume ventilation

You must provide an efficient fume ventilation system to remove the fumes generated by the atomizer.

The specification for a suitable system is described in brochure no. L-301 – *Preparing Your Laboratory for PerkinElmer Atomic Spectroscopy Instruments*. This brochure is available from PerkinElmer.

Warning: Toxic Gases – Fume Ventilation System

Without adequate ventilation, potentially toxic vapors can build up in the laboratory.

- *Your laboratory must have a reliable fume ventilation system before you use this instrument.*

Exhaust venting is important for the following reasons:

- It will protect laboratory personnel and the instrument from toxic vapors, which may be produced by some samples.
- It will minimize the effects of room drafts and the laboratory atmosphere on flame stability.
- It will remove heat produced by the flame.

We strongly recommend that you do not place the spectrometer in a chemical hood. If a chemical hood must be used, make arrangements to prevent contamination due to corrosive vapors and backdrafts from other hoods. **Do not** perform sample preparation in the same hood as the spectrometer.

Space requirements for the spectrometer system

The illustrations on the following pages show the dimensions of the spectrometer and the space that you require to accommodate the spectrometer, computer, and ancillary equipment.

Computer and printer

Place the computer near to the spectrometer in a convenient position for working with both instruments. A movable table, specially designed for a computer and printer, may be most convenient.

The printer needs to be in a convenient place with enough room for the paper supply to stand and feed smoothly into the printer.

Accessories

Make sure that there is sufficient room for any accessories near to the spectrometer.

For example, we recommend that you place the (optional) FIAS flow injection system on a movable table in front of the atomizer compartment to minimize the length of the sample transfer tube. Make sure that the wheels of the table can be securely locked.


Figure 6-1. Space Requirements for the Spectrometer System (1)


Figure 6-2. Space Requirements for the Spectrometer System (2)

Electrical requirements

The AAnalyst 700 spectrometer requires a stable electrical supply of 230 V, 50 Hz or 60 Hz, capable of delivering 20 A, that has a correctly wired protective earthing system (ground connection).

The line power cable has a 3-pole connector (line + neutral + protective earth), conforming to IEC 309-2, rated for 20 A. The cable is 3 meters long.

The line power supply should be free of line transients in excess of 50 V peak. If the electrical supply voltage produces large AC line voltage fluctuations, a qualified electrician should install a voltage regulator between the electrical outlet and the instrument.

The power demand of the graphite furnace is approximately 3600 VA, the maximum graphite tube temperature is 3000 °C. Monitoring of the energy required for each user-defined furnace program reliably prevents the line power supply from being overloaded.

The line power supply point must be within 3 meters of the rear of the spectrometer.

We recommend that you connect the spectrometer, computer, printer, and any accessories to the same phase of the line power supply and the same protective earth.

The line power supply must conform with local safety regulations and be checked by a qualified electrician before you connect the instrument to line power.

See the guides for the computer and printer and any other instruments for information about their power requirements.

Gas supply requirements

The permanent installation of gas supplies is the responsibility of the user and should conform to local safety and building regulations. PerkinElmer recommends that cylinders be stored outside the laboratory and that gases be piped to the spectrometer system in approved gas lines.

- The user must provide the gas supplies, regulators, connectors, and valves. A double-stage or heated regulator is recommended for nitrous oxide. Either the gas supplier or PerkinElmer can supply the correct regulators and valves.
- You must be able to reach the on/off valves easily and see the pressure indicators.
- Consult *Chapter 1, Safety and Regulatory Information* for recommended safety procedures.

**Warning: Compressed Gases**

High pressure gas cylinders can be dangerous if mishandled or misused.

- *Always handle gas cylinders with caution and observe your local regulations for the safe handling of gas cylinders.*

Burner gases

Two types of flame are commonly used for atomic absorption analyses with flame atomization: the air/acetylene flame and the nitrous oxide/acetylene flame. The specifications and requirements for the burner gases are given below:

Table 6-1. Specifications and Requirements of the Burner Gases

Gas	Purity	Outlet Gauge Pressure					
		kPa		bar		psig	
		$P_{\min}$	$P_{\max}$	$P_{\min}$	$P_{\max}$	$P_{\min}$	$P_{\max}$
Acetylene (ethyne, C ₂ H ₂)	99.996%, atomic absorption grade in acetone	90	100	0.9	1.0	13.0	14.5
Air	oxygen ≤20%; particle free; water free	450	500	4.5	5.0	65.0	72.5
Nitrous oxide (dinitrogen monoxide, N ₂ O)	99.996%	450	500	4.5	5.0	65.0	72.5

- Set the outlet gauge pressure for each gas to a value just below the maximum value ($P_{\max}$).
- Always set the same outlet gauge pressure for air and nitrous oxide.
- **Never** set the outlet pressure to a value exceeding the maximum pressure ($P_{\max}$).

Warning: Acetylene – Explosion Hazard

WS.7

Acetylene can decompose explosively at pressures higher than 103 kPa/ 1.03 bar/15 psig.

- Always make sure that the acetylene outlet gauge pressure is below this value.

Air supply

For the burner system, an air supply that is clean, dry, and oil-free is required. The air supply should provide a minimum of 28 liters/min. (1 cubic foot/min.) at a minimum pressure of 450 kPa (65 psig). The air supply should be suitably filtered to remove oil, water, and solid particles. Since compressors are generally noisy, we recommend that you locate the compressor outside the laboratory.

Cylinders of compressed air should only be used as an emergency or short-term solution.

Warning: Flashback Hazard

W5.9

PerkinElmer burner systems are designed for use with compressed air. The use of oxygen can cause an explosion in the burner system, and oxygen-enriched air can cause a flashback of the flame.

- *Never use oxygen or oxygen-enriched air with PerkinElmer burner systems. Use only compressed air.*

Air Compressor: A small oil-less air compressor meeting the air supply requirements is recommended. The use of an air dryer filter is strongly recommended to remove entrained water, oil, water aerosols, and solid particles from compressed air lines. Water in the air supply can cause severe damage to the pneumatics.

Caution

Risk of damage to the instrument

Water in the air supply can damage to the pneumatics.

- **Use only dry air.** *Bleed compressor tanks to remove water at least once a week.*
-

Acetylene supply

The acetylene supply should be of ‘atomic absorption’ grade.

Warning: Acetylene – Explosion Hazard

Acetylene can decompose explosively at pressures higher than 103 kPa/ 1.03 bar/15 psig.

- *Always make sure that the acetylene outlet gauge pressure is below this value.*

Grade: Suitable acetylene typically has a minimum purity specification of 99.6% with the actual assay being about 99.8%. Acetylene is supplied dissolved in acetone to prevent explosive decomposition. You should not use acetylene that is dissolved in any other solvent. Experience indicates that solvents can be carried through the regulator with the acetylene and can weaken, to the point of rupture, the walls of the tubing carrying the acetylene to the burner, causing a potential hazard.

Cylinder Specifications: A size 1A acetylene cylinder contains about 8500 liters (300 cubic feet) of acetylene and usually lasts about 30 hours of burning time with an air-acetylene flame. With a nitrous oxide-acetylene flame, the same size cylinder will last about 10 hours.

Flash Arrestor: PerkinElmer recommends the use of the Flash Arrestor (P/N N930-0068). The flash arrestor provides a positive shutoff of gas by checking the reverse flow under all conditions. It also stops the gas supply when a flashback occurs. The flash arrestor can be directly connected to the acetylene fuel line.

Piping/Tubing/Fittings: Use galvanized iron, steel, or wrought iron tubing. Never use copper tubing with acetylene. Joints may be welded or made up of threaded or flanged fittings, typically stainless steel or aluminum, or brass composed of less than 65% copper. Rolled, forged, or cast steel or malleable iron fittings may also be used. Cast iron fittings cannot be used safely for acetylene lines.

Warning: Risk of explosion

Acetylene can react with copper to form a readily explosive compound.

- *Never use copper tubing or fittings for acetylene gas lines.*

Nitrous oxide supply

The nitrous oxide supply should be at least 99% pure. A heated regulator is recommended.

Cylinder Specifications: Cylinders of nitrous oxide (99.0% minimum purity) are available from local suppliers.

Nitrous oxide is supplied in the liquid state, initially at a pressure of about 5000 kPa (750 psig). A size 1A cylinder of nitrous oxide contains about 14800 liters (520 cubic feet) and will typically last for 10–12 hours of burning time.

Warning: Nitrous oxide is a strong oxidant and an asphyxiant

Nitrous oxide is an asphyxiant; it can be considered a potential hazard whenever it is stored or used in a confined space.

Spontaneous combustion may occur if nitrous oxide comes into contact with grease, oil, or other organic material.

- *Store nitrous oxide cylinders in a well-ventilated area outside the laboratory.*
- *Do not store nitrous oxide cylinders in close proximity to flammable gases (e.g., acetylene).*
- *Make sure tubing and fittings carrying nitrous oxide are free of grease, oil, or other organic material.*

Heated Regulator: Since the nitrous oxide is in liquid form, the pressure gauge does not give a true indication of how much nitrous oxide remains in the cylinder until the pressure starts to fall rapidly as the residual gas is drawn off. When nitrous oxide is rapidly removed from the cylinder, the expanding gas causes cooling of the cylinder pressure regulator and the regulator diaphragm can freeze. This can create erratic flame conditions or, in the most extreme case, a flashback.

It is therefore advisable to use a regulator with a built-in heater. A heated nitrous oxide pressure regulator is available from PerkinElmer.

Furnace gases

Inert gas

The graphite furnace requires a supply of inert gas to prevent the tube and the analyte atoms from being oxidized when the tube is heated.

Argon is the recommended inert gas.

If you cannot obtain argon, you can use purified nitrogen, but nitrogen does not provide such efficient protection as argon at temperatures above 2000 °C.

Nitrogen can also react with the samples and the graphite tube at typical atomization temperatures.

Warning: Toxic Gases

W6.1

If you use nitrogen as the inert gas, toxic cyanogen gas (CN)₂ will be produced when the furnace temperature is above 2300 °C.

- *Make sure that the main laboratory fume ventilation system removes all gases and fumes emitted by the furnace.*

Special gas

For some applications you may want to use a special gas. A typical example is the use of air during one of the pyrolysis steps to better ash organic materials.

Table 6-2. Specifications and Requirements of the Furnace Gases

Gas	Purity	Outlet Gauge Pressure						Max. Flowrate
		kPa		bar		psig		L/min
		$P_{\min}$	$P_{\max}$	$P_{\min}$	$P_{\max}$	$P_{\min}$	$P_{\max}$	1.2
Argon	99.996%	350	400	3.5	4.0	50.0	58.0	1.2
Nitrogen	99.996%	350	400	3.5	4.0	50.0	58.0	1.2
Special	particle free; water free	350	400	3.5	4.0	50.0	58.0	1.2

- Set the outlet gauge pressure for each gas to a value just below the maximum value ($P_{\max}$).
- **Never** set the outlet pressure to a value exceeding the maximum pressure ($P_{\max}$).

Cooling the spectrometer system

The graphite furnace uses water cooling to allow the graphite tube to return to its equilibrium temperature (approximately ambient) quickly between measurements (see *'Cooling the furnace' in the Graphite Furnace User's Guide*). Cooling coils are also incorporated in the spectrometer's power supply to maintain the electronic components at a comfortable working temperature.

Cooling system

PerkinElmer offers a recirculating cooling system as an option. It provides coolant at a constant temperature and flow rate, to allow the most reproducible analytical conditions. It also makes the system independent of external water supplies and saves water (see *'Cooling the furnace' in the Graphite Furnace User's Guide*). This cooling system includes a pressure-sensitive safety device that opens a bypass if the internal pressure in the circulatory system increases to above a safe value.

Power requirements: 230 V; 50/60 Hz (electrical power is supplied by the spectrometer)

Part Number: B019-8072

External cooling

If you intend to run the spectrometer system without the PerkinElmer cooling system, you must connect another suitable source of cooling, such as the main water supply or a laboratory thermostat. If you connect the spectrometer system to the main water supply, make sure that the water pressure is as given in Table 6-3.

Table 6-3. Requirements for the Cooling Water Supply

Water quality	sediment-free drinking water
Temperature	20–40 °C
Flow rate	not less than 2 L/min
Pressure	2.5–3.5 bar (36–50 psig)
pH	6.5–7.5
Hardness	not greater than 14 °dH or 100 ppm

Unfiltered industrial water should **not** be used.

A water hose kit with pressure relief valve is provided with the spectrometer system which you must permanently install in the water supply line to the system. This valve opens if the water pressure to the system exceeds 4 bar. This feature eliminates the risk of damage to the water hoses and hose fittings in the system through excess water pressure.

The water pressure to the system should be at least 2.5 bar (36 psig). If the pressure is lower, the flow rate of water through the system is reduced and heat will not be removed quickly enough.

The maximum water pressure to the system is 3.5 bar (50 psig); the pressure relief valve opens automatically if the water inlet pressure exceeds 4 bar – the slightly higher pressure of 4 bar allows for minor daily fluctuations in the pressure of the main water supply.

Lamps

These are often called primary source lamps, atomic absorption source lamps, spectroscopic lamps, or line sources.


You need to provide the lamps for your particular applications.

Brochure no. L-656 – *Hollow Cathode Lamps and Electrodeless Discharge Lamps for Atomic Absorption*, contains descriptions of all the lamps that are available. This brochure is available from PerkinElmer.

Computer and printer

PerkinElmer offers a suitable computer for controlling the spectrometer and graphite furnace, and a suitable printer. Contact your local PerkinElmer office for more information.

PerkinElmer cannot accept responsibility for incorrect functioning of the spectrometer system if you are not using the computer offered with the spectrometer by PerkinElmer.

Spectrometer: System Description

7

Contents	page
The AAnalyst 700 Spectrometer System	7-3
The optical system	7-6
Switches and connectors on the spectrometer	7-8
Technical data	7-14

The AAnalyst 700 Spectrometer System

Concept

The AAnalyst 700 spectrometer from PerkinElmer is a high performance atomic absorption spectrometer incorporating a burner system for flame atomization and a graphite furnace for electrothermal atomization. The spectrometer is capable of fully automatic single or multi-element analyses. The complete system is operated from the associated computer using the AA WinLab software.

Features of the AAnalyst 700 System

- Double-beam optical system optimized for maximum radiation throughput.
- Two fully integrated atomizers: a burner system for flame atomization and a graphite furnace for electrothermal atomization. Atomizer changeover is automatic at the click of a mouse button.
- Longitudinally-heated graphite furnace featuring full Stabilized Temperature Platform Furnace (STPF) technology.
- Deuterium background correction system.
- Automated self-optimization of the burner system.
- High capacity furnace autosampler.
- Easy-to-use AA WinLab software with which you control the spectrometer system. AA WinLab also incorporates a range of powerful quality control functions.

System components

The AAnalyst 700 is a compact benchtop spectrometer. The major components of the system are:

- The spectrometer, containing the optical system, the electronics, and the atomizer compartment.
- The burner system, mounted on a carriage that can be driven into and out of the atomizer compartment by software command (*see 'Burner design' in the Burner System User's Guide for more information*).
- The graphite furnace, mounted on a carriage that can be driven into and out of the atomizer compartment by software command (*see 'Furnace design' in the Graphite Furnace User's Guide for more information*).
- The furnace autosampler, mounted on a swing arm in front of the atomizer compartment (*see 'Autosampler design' in the Graphite Furnace User's Guide for more information*).
- The (optional) furnace cooling system (*see 'Cooling the furnace' in the Graphite Furnace User's Guide for more information*).
- The computer and the associated AA WinLab software.
- The printer; this is optional but strongly recommended.

System control

To control the spectrometer and perform your analyses, you use the AA WinLab application program, which runs under the Microsoft Windows operating system. The computer keyboard, the screen, and the mouse form the user interface. When you switch on the spectrometer system it performs an initialization procedure to check the instrument. While the spectrometer is operating, the control system continually monitors many instrument functions. It generates appropriate messages when a problem occurs and prevents potentially hazardous situations from arising.

You can obtain more information about the application program in the AA WinLab online help.

Accessories

PerkinElmer offers autosamplers for operation with the burner system, flow injection systems, diluters, and much more. Please contact your local PerkinElmer office to obtain more information.

The optical system

The AAnalyst 700 uses a double-beam optical system as depicted in Figure 7-1. Front-surfaced, reflecting optical components are used throughout. The mirror surfaces are coated with silica for durability.

The lamp compartment has positions for up to eight spectral sources. The radiation from the requested source lamp is reflected by a movable mirror into the optical path.

A beam splitter splits the radiation into a sample beam and a reference beam. The sample beam passes through the atomizer while the reference beam passes around the atomizer and is conducted by a fiber optic. The sample and reference beams are measured simultaneously.

The Littrow-configuration monochromator has a large area reflecting grating blazed at two wavelengths; 236 nm for operation in the UV spectral range and 597 nm for operation in the visible spectral range. The monochromator offers a choice of three spectral slit widths. Two slit heights are available for each slit width. The high slits are used with the flame technique while the low slits are used with the graphite furnace technique. The detector is a segmented array of discrete photodiodes.

For continuum source background correction, the radiation from the deuterium lamp is merged with the source radiation at the beam splitter.


Figure 7-1. Optical System

Switches and connectors on the spectrometer

Operational on/off switches


Figure 7-2. Switches and Ports on the Front Panel

Spectrometer

You normally switch the system on and off with the operational on/off switch on the front panel of the spectrometer.

Note: The spectrometer is provided with an operational on/off switch and with a circuit breaker.

The operational on/off switch turns the instrument on and off, but it does not disconnect the instrument from the line power supply.


The circuit breaker functions as a 'disconnecting device' and disconnects the instrument entirely from the line power supply.

Cooling system

The (optional) cooling system requires a 230 V line power supply, which is supplied and controlled from the spectrometer. A special line power cord is provided for this purpose. Do not connect the cooling system to an external power source.

You can leave the on/off switch on the cooling system in the 'on' position since the system is controlled from the spectrometer.

Cooling System

Accessory ports

Three accessory ports are located on the front panel of the spectrometer behind slide-up covers. You use these ports to connect certain accessories to the spectrometer, such as the furnace autosampler. The ports provide the electrical supply and control signals for operation of the accessories.

Electrical connectors and controls on the rear panel

Item	Designation	Function
1	Instrument	Socket for a read cable to an external instrument to activate a read cycle.
2	Communication	Socket for the communication cable connecting the spectrometer to the computer.
3	Cooling System ~230 V 50/60 Hz	Power outlet for the (optional) cooling system. Connect the cooling system to this socket only. Do not connect other equipment to this power outlet.
4	 ~230 V 50/60 Hz	Power outlet for a printer.
5	 ~230 V 50/60 Hz	Power outlet for a computer.
6	 ~230 V 50/60 Hz	Power outlet for a monitor.
7		Circuit breaker. An automatic switch that prevents the system from drawing too much current. It can also be used to isolate the spectrometer system from the line power supply.
8		Captive line power cord for the complete system.
9	F1	Fuse for the cooling system power outlet; 4 A slow-blow (T).
10	F2	Fuse for the printer power outlet; 4 A slow-blow (T).
11	F3	Fuse for the computer power outlet; 4 A slow-blow (T).
12	F4	Fuse for the monitor power outlet; 4 A slow-blow (T).


Figure 7-3. Electrical Connectors and Switches on the Rear Panel

The gas and cooling connectors

Item	Designation	Function
Furnace Cooling		
1	In	Fitting for connection of the cooling water inlet hose to the spectrometer.
2	Out	Fitting for connection of the cooling water outlet hose from the spectrometer.
Furnace Gases		
3	Ar/N₂ Normal Gas	Connector for the connection of the normal gas supply to the furnace (argon or nitrogen).
4	Special Gas	Connector for the connection of the special gas supply to the furnace. If you are not using a special gas, you must also connect the normal gas to this connector.
Burner Gases		
5	Air/Luft	Fitting for connection of the air supply to the spectrometer. The air supply is required for the burner and also for the pneumatic system in the spectrometer.
6	N₂O	Fitting for connection of the nitrous oxide supply to the burner.
7	C₂H₂	Fitting for connection of the acetylene supply to the burner.


Figure 7-4. Gas and Cooling Connectors on the Rear Panel


Technical data

Principle	Atomic absorption spectrometer. Computer controlled, using a special application program running under a graphical user interface. Real-time double-beam, all-reflecting optical system. Front-surfaced optics protected by a silica coating. AC modulated spectral line sources.
Monochromator	Littrow configuration grating monochromator. Grating: ruled area: 64 mm x 72 mm, 1800 lines/mm blaze wavelengths: 236 nm and 597 nm Focal length: 267 mm Slit widths: 0.2, 0.7 and 2.0 mm Wavelength range: 190–870 nm Reciprocal linear dispersion: 1.6 nm/mm (nominal)
Detector	Segmented discrete solid-state detector with CMOS charge amplifier array; software-controlled automatic gain control (AGC).
Lamps	Line radiation sources. Positions for up to 8 lamps in lamp compartment. Single- and multi-element hollow cathode lamps (HCLs) and up to 4 electrodeless discharge lamps (EDLs) can be used.
Background correction	Continuum source (deuterium lamp).
Atomizers	Flame atomization: Premix burner system with pneumatic nebulizer and interchangeable burner heads; automatic changeover between air/acetylene and nitrous oxide/acetylene flames. Electrothermal atomization: Longitudinally-heated graphite furnace.

Signal handling	Time average, peak average, peak height, peak area. Measurement ranges: - absorbance: -0.5 to 2.000 - concentration: 0.00001 to 9999 concentration units - expansion factor: 0.01 to 100 - read time: 1 to 60 s - read delay: 0 to 60 s
Calibration	Automatic zeroing, resloping, and recalibration. Calibration/recalibration using up to 15 calibration solutions ('standards'). Calibration coefficients determined by a least squares procedure. Analytical curve technique: - Linear with zero or non-zero intercept. - Non-linear, 2-coefficient, or 3-coefficient, zero or non-zero intercept. Analyte addition technique: - Zero or non-zero intercept. - Additions calibrate (zero intercept)
Application program	Special application program: AA WinLab. Online help function. Hard disk or floppy storage of measurement data, signals, and completed analytical methods. Special Reporter, Reformatter, and Library Manager applications for generating data reports and reformatting data for export to other data handling programs.
Graphics	High resolution display of the signals. Simultaneous display of the corrected analyte signal and the background signal. Calibration curve display. Storage and recall of transient absorbance signals. Hard copy printout of screen graphics on an associated printer.
Power requirements	230 V single-phase, $50/60$ Hz, power consumption 5000 VA.

Spectrometer System Description

Electrical protection	As defined in EN 61010-1-1993 (IEC 1010-1): Insulation: Class I Insulation category (overvoltage category): II Pollution degree: 2		
Safety standards	As defined in EN 61010-1-1993 (IEC 1010-1: 1990 + A1 + A2) and EN 61010-2-061:1996. CSA standards: CAN/CSA-C22.2 No. 1010.1-92, Amendment 1 (1992), 2 (1995), and CAN/CSA 61010-2-061 (1995). UL standards: UL 3101-1 1993 and UL 3111-1 1994. Compliant to the requirements of Canadian Standards Association (CSA) and Underwriters Laboratory (UL) as listed in File No. LR 45111. Certificate of Approval and GS mark from the LGA Bavaria, FRG, Institute of Equipment Safety Testing and Certification, and CSA-NRTL/C mark from CSA, Canada, File No. LR 45111.		
EMC standards	As defined in EN 61326-1:1997.		
Environmental requirements	Ambient temperature: +15 °C to +35 °C 20–80% relative humidity; non-condensing		
Dimensions	Width: 1095 mm	with handles: 1220 mm	
	Height: 655 mm	with door open: 865 mm	
	Depth: 700 mm	with handles: 770 mm	
		with autosampler: 1040 mm	
Mass (weight)	147 kg Spectrometer alone, not including the computer, autosampler, and cooling system.		

Translations of Warnings

8

Contents	page
Translations of Warnings	8-3
Symbols used on the instrument	8-45

Translations of Warnings

This section contains translations of the safety conventions and warnings used in this guide.

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This symbol alerts you to situations that could result in **personal injury** to yourself or other persons.
Details about these circumstances are in a box like this one.

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This symbol alerts you to the risk of **electric shock** that could result in **personal injury** to yourself or other persons.

Details about these circumstances are in a box like this one.

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This symbol alerts you to the risk of **hot surfaces** that could cause **personal injury** to yourself or other persons.
Details about these circumstances are in a box like this one.

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This symbol alerts you to the risk of **ultraviolet radiation** that could cause **eye damage** to yourself or other persons.

Details about these circumstances are in a box like this one.

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Caution:

The term CAUTION alerts you to situations that could result in **serious damage to the instrument** or other equipment.

Details about these circumstances are described in a message similar to this one.

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
Index

A

AAAnalyst 700
 concept 7-3
 features 7-3
 system components 7-4
 system control 7-4
 technical data 7-14
 accessory ports 7-9
 acetylene
 connecting 2-8
 cylinder specifications 6-12
 explosion hazard 6-12
 fittings 6-12
 flash arrestor 6-12
 grade 6-12
 safe use 1-18
 tubing 6-12
 warning, copper 6-12
 acetylene supply
 maintaining 4-14
 acetylides 1-21
 air filters
 changing 4-10
 air supply
 compressor requirements 6-11
 connecting 2-8
 warning 6-11
 argon
 connecting 2-11
 automatic circuit breaker
 description 7-8, 7-10
 resetting 4-8

B

burner drain system
 connecting 2-15
 burner gases
 connecting 2-7
 specifications and requirements 6-10
 burner system
 drain system 1-20
 safe operation 1-19

C

changing air filters 4-10
 changing fuses 4-8
 checklist
 general maintenance 4-4
 spectrometer maintenance 4-4
 circuit breaker *see* automatic circuit breaker
 cleaning windows 4-5
 compressed gases 1-13
 connecting
 acetylene 2-8
 air 2-8
 argon 2-11
 burner drain system 2-15
 burner gases 2-7
 external cooling sources 2-19
 furnaces gases 2-10
 nitrogen 2-11
 nitrous oxide 2-8
 special gas 2-13
 connecting furnace gas tubes 2-11
 connecting water hoses 2-16
 connections, electrical 2-22

connectors on the spectrometer 7-8
 cooling 7-12
 electrical 7-10
 gas 7-12
 cooling
 graphite furnace 6-16
 spectrometer system 6-16
 cooling connectors 7-12
 cooling system (recirculating) 6-16
 connecting and disconnecting water hoses 2-16
 filling 2-18
 furnace cooling 6-16
 installing 2-16
 spectrometer cooling 6-16
 correct use of the instrument 1-4
 customer service – *see end page*
 cyanide solutions 1-21

D

deuterium lamp
 replacing 4-11
 dimensions of the system 6-5
 dimensions, spectrometer 7-16
 disconnecting furnace gas tubes 2-11
 disconnecting water hoses 2-16
 disposing of contents of waste containers 1-10
 disposing of defective lamps 1-10
 drain system 1-20

E

- electrical components
 - maintaining 4-7
- electrical connections 2-22
- electrical connectors 7-10
- electrical protection 1-27, 7-16
- electrical requirements 6-8, 7-15
- electrical safety 1-5
- electrodeless discharge lamps
 - installing and connecting 3-8
- electromagnetic compatibility 1-28, 7-16
- electrothermal atomization
 - hazards 1-23
 - safety practices 1-22
- emergency shutdown 3-14
- environment
 - operating conditions 1-7
 - storage conditions 1-7
- environmental requirements 7-16
- external cooling 6-16
 - connecting 2-19
 - furnace 6-16
 - spectrometer 6-16

F

- FIAS and mercury/hydride systems
 - hazards 1-25
 - safety practices 1-24
- filling the recirculating cooling system 2-18
- flame
 - safe operation 1-19
- flame atomization
 - safe operation of the flame 1-19
 - safe use of burner gases 1-16
 - safety practices 1-16

- flashbacks 1-20
- fume ventilation 6-4
- furnace autosampler
 - installing 2-15
- furnace gas tubes
 - connecting 2-11
 - disconnecting 2-11
- furnace gases
 - connecting 2-10
 - specifications and requirements 6-15
- fuses
 - changing 4-8

G

- gas connectors 7-12
- gas hazards 1-13
- gas supply requirements 6-9
 - acetylene 6-12
 - air 6-11
 - burner gases 6-10
 - furnace gases 6-14
 - nitrous oxide 6-13

H

- hazards with electrothermal atomization
 - high temperatures 1-23
 - toxic products 1-23
 - UV radiation 1-23
- hazards with FIAS and mercury/hydride systems
 - hazardous chemicals 1-25
 - high temperatures 1-25
 - toxic products 1-25
- hazards with flame atomization
 - acetylides 1-21
 - cyanide solutions 1-21
 - high temperatures 1-21
 - toxic combustion products 1-21
 - UV radiation 1-21

- heated regulator for nitrous oxide 6-13
- high temperatures 1-12
- hollow cathode lamps
 - installing 3-8
- hygiene 1-8

I

- IEC 1010 compliance 1-27
- installation procedures 2-3
- installing
 - recirculating cooling system 2-16
 - sample trays 3-11
- installing and connecting lamps 3-8
- installing furnace autosampler 2-15

L

- laboratory hygiene 1-8
- laboratory requirements
 - cooling system (recirculating) 6-16
 - cooling the spectrometer system 6-16
 - fume ventilation 6-4
 - gas supplies 6-9
 - operating conditions 6-3
 - water cooling 6-16
- laboratory ventilation 1-8
- lamp compartment 3-7
- lamps
 - installing and connecting 3-8
 - removing 3-10

M

maintaining electrical components 4-7
 maintaining the acetylene supply 4-14
 maintenance checklists 4-4
 maintenance procedures 4-5
 markings on the instrument 1-30
 mass 7-16
 moving the spectrometer 2-4

N

nitrogen
 connecting 2-11
 nitrous oxide
 connecting 2-8
 heated regulator 6-13
 warnings 6-13

O

operating conditions 1-7, 6-3
 operational on/off switches 7-8
 optical system 7-6

P

parts and supplies 5-1
 PerkinElmer service 4-3
 power requirements 7-15

R

Reconfigure utility 3-4
 references for laboratory safety practice 1-26
 regulatory information
 electromagnetic compatibility 1-28, 7-16
 removing lamps 3-10
 replacement parts 5-4
 catalog service 5-4

replacing air filters 4-10
 replacing fuses 4-8
 replacing the deuterium lamp 4-11
 required pressures
 burner gases 6-10
 furnace gases 6-15
 resetting the circuit breaker 4-8

S

safe handling of chemicals 1-9
 safe use of burner gases
 acetylene 1-18, 6-12
 air 1-17, 6-11
 nitrous oxide 1-17, 6-13
 safety
 compressed gases 1-13
 correct use of the instrument 1-4
 cyanide solutions 1-21
 electrical safety 1-5
 handling chemicals 1-9
 high temperatures 1-12
 operation of the flame 1-19
 UV radiation 1-11
 safety and EMC standards 1-28, 7-16
 safety conventions 1-3
 translations 8-3
 see also symbols
 safety practices
 electrothermal atomization 1-22
 FIAS and mercury/hydride systems 1-24
 flame atomization 1-16
 references 1-26
 sample trays
 installing 3-11

shutting down the system 3-12
 space requirements 6-5
 special gas
 connecting 2-13
 spectrometer
 dimensions 7-16
 mass (weight) 7-16
 moving 2-4
 spectrometer maintenance 4-5
 storage conditions 1-7
 switches on the spectrometer 7-8
 operational on/off switch 7-8
 switching on the system 3-3
 symbols used in the guide 1-3
 symbols used on the instrument 1-29
 translations 8-45
 system
 shutting down 3-12

T

technical data
 spectrometer 7-14
 translations
 safety conventions 8-3
 symbols used on instrument 8-45
 warnings in guide 8-3
 troubleshooting 4-3
 tubing
 acetylene supply 6-12

U

UV radiation 1-11

V

ventilation, laboratory 1-8

W

warning markings on the instrument 1-30

warnings

acetylene and copper 6-12

acetylene explosion hazard
6-12

acetylene pressure 6-12

air supply 6-11

nitrous oxide 6-13

warnings in guide

translations 8-3

waste disposal

contents of waste containers
1-10

defective lamps 1-10

water cooling 6-16

connecting 2-19

water hoses

connecting 2-16

disconnecting 2-16

weight 7-16

windows

cleaning 4-5

Customer Service

Company Name and Address

PerkinElmer Analytical Instruments – Headquarters

761 Main Avenue, Norwalk
Connecticut 06859-0001, U.S.A.

Tel.: 203-762-4000, Fax.: 203-762-4050

PerkinElmer on the Internet

<http://www.perkinelmer.com>

E-mail: info@perkin-elmer.com

Supplies, replacement parts, and accessories

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761 Main Avenue, Norwalk, CT 06859-0010, U.S.A.,