

## Thermo Scientific FOCUS GC

Full-sized GC performance

Specifically targeted for routine quality control and high productivity laboratories, the Thermo Scientific FOCUS GC is a compact, single-channel gas chromatograph. Conceived and built around the "core" hardware and unique features of the Thermo Scientific TRACE GC Ultra, the FOCUS GC™ incorporates new, best of breed electronics. The FOCUS GC is the reference for precision and accuracy of results, as well as a standard for durability in QA/QC and high throughput environments.



### Thermo Scientific AI/AS 3000 II Autosampler

Combined with the 8-sample AI 3000 II Autoinjector or the 105-sample AS 3000 II Autosampler, the FOCUS GC is positioned in line with your laboratory's throughput requirements at any moment in time, thus preserving scalability.

### Thermo Scientific TriPlus Autosampler

For additional productivity and flexibility, the Thermo Scientific TriPlus Autosampler is the ideal solution, able to simultaneously automate two adjacent FOCUS GC gas chromatographs. Scalability is preserved, with capacities for up to 300 samples and multiple injection modes.



### Compact Size

Small footprint for bench space optimization: 430 x 510 x 350 mm (HxDxW)

### Optimized User Interface

Allows simple control of the main functions of the GC while preventing access to the more critical parameters of the method from unauthorized personnel. In line with the requirements of QA/QC environments.

### Automatic Column Characterization

Automatic measurement of the effective column's resistance for exceptional column-to-column and instrument-to-instrument retention time reproducibility at no additional cost.

### Wide Oven

Wide oven with front opening allows easy accessibility for straightforward installation of capillary columns mounted on standard diameter cages. Oven dimensions are 225 x 225 x 140 mm (HxWxD)

### Additional Instrument Control

Instrument start-up and real-time diagnosis, on-line help describing each parameter, logbook with deviation records, and local keyboard lock provide complete and safe system handling.

### Leak Check

Automatic leak check function for rapid routine system control.

## Thermo Scientific FOCUS GC Gas Chromatograph

### Features and Technical Specifications

#### Specifications

##### Power Supply

110 V and 230 V power supply

##### Base Unit

Height 430 mm, Depth 510 mm, and Width 350 mm

Weight: 30 kg

##### Column Oven

Compatible with standard capillary column cages

Temperature range: from few degrees above ambient to 350 °C

Program rates: 0.1 to 120 °C/min

Ramps: 7

Plateaus: 8

Typical cool-down time: 350 °C to 50 °C in 270 seconds

Typical heat-up time: 50 °C to 350 °C in 270 seconds

##### Injectors

Split-Splitless Temperature range: 50 °C to 375 °C

Purged-Packed Temperature range: 50 °C to 400 °C. Suitable for both wide bore and packed columns

##### Digital Pressure/Flow Controller (1000 kPa)

High pressure pneumatic module

Ambient temperature compensation: Yes

Ambient pressure compensation: Yes

Gas saver: Yes

Automatic column characterization: Yes

Automatic leak check: Yes

Pressure surge: Yes

##### Base Body with Gas Control for FID

Includes a complete base body with all relevant pneumatic controls for three gas lines

The temperature of the base body is finely controlled in the range of 50 °C to 375 °C in 1 °C increments

The pneumatic module includes three pressure controllers, calibrated restrictors, and on/off valves for H<sub>2</sub>, Air, and make-up gases

##### Detectors

##### MDA

##### Linearity

##### Selectivity or additional features

Flame Ionization Detector

2 x 10<sup>-12</sup> gC/sec

Better than 10<sup>7</sup>

Flameout detection and timed programming capability. Acquisition rate: 300 Hz

Thermal Conductivity Detector

600 pg Ethane/mL He

10<sup>6</sup>

Automated software switch function

##### Gas Sampling Valve

6-port sampling/switching valve available as an optional accessory for petrochemical and industrial applications. Equipped with selectable sampling loops.

##### Autosamplers

FOCUS GC is compatible with AI 3000 II Autoinjector, AS 3000 II Autosampler, and TriPlus Autosampler



[www.thermoscientific.com](http://www.thermoscientific.com)

©2011 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries. Specifications, terms and pricing are subject to change. Not all products are available in all countries. Please consult your local sales representative for details.

**Africa-Other** +27 11 570 1840  
**Australia** +61 3 9757 4300  
**Austria** +43 1 333 50 34 0  
**Belgium** +32 53 73 42 41  
**Canada** +1 800 530 8447  
**China** +86 10 8419 3588  
**Denmark** +45 70 23 62 60

**Europe-Other** +43 1 333 50 34 0  
**Finland/Norway/Sweden** +46 8 556 468 00  
**France** +33 1 60 92 48 00  
**Germany** +49 6103 408 1014  
**India** +91 22 6742 9434  
**Italy** +39 02 950 591

**Japan** +81 45 453 9100  
**Latin America** +1 561 688 8700  
**Middle East** +43 1 333 50 34 0  
**Netherlands** +31 76 579 55 55  
**New Zealand** +64 9 980 6700  
**Russia/CIS** +43 1 333 50 34 0  
**South Africa** +27 11 570 1840

**Spain** +34 914 845 965  
**Switzerland** +41 61 716 77 00  
**UK** +44 1442 233555  
**USA** +1 800 532 4752

COMPANY WITH QUALITY MANAGEMENT  
 SYSTEM CERTIFIED BY DNV  
 = ISO 9001:2008 =

Thermo Fisher Scientific S.p.A.  
 Milan, Italy is ISO Certified.



Thermo Fisher Scientific (Shanghai)  
 Instruments Co., Ltd. is ISO certified.

PS10056\_E 05/11C

**Thermo**  
 SCIENTIFIC