

Cole-Parmer®

Cole-Parmer® Programmable Vacuum Ovens

- Easily manage up to eight programs with eight steps each
- PID temperature controller provides accurate and reliable temperature control
- Unit stops heater automatically if temperature is out of control
- Full-view window allows continuous monitoring
- Rounded inner chamber corners offer easy cleaning
- Rubber O-ring and door lock handle create a secure seal to prevent air leakage or unauthorized access
- Stainless steel inner chamber and powder-coated steel plate exterior resist corrosion for longer lifetime



coleparmer.com

USA:
+1.800.323.4340
+1.847.549.7600

Canada: +1.800.363.5900
China: +86.21.5109.9909
France: +33 (0) 1486 37800
Germany: +49 (0) 9377 92030

India: +1.800.266.1244
Italy: +39 (0) 1313 89513
UK: +44 (0) 1480.272279
All other countries: +1.847.549.7600

Cole-Parmer is an
Antylia Scientific company
Find out more at antylia.com





Cole-Parmer® Programmable Vacuum Ovens

Easily set up 8 programs for parameters such as temperature, time, and fan speed. Our ovens are designed with independent vacuum interface and vents, an air inlet for easy access to inert gases such as nitrogen, and a solenoid valve that controls vacuum degree and vacuum pump. Multiple safety features include temperature deviation alarm, overtemperature set point alarm, overcurrent protection, and auto-start after power loss/return.

Each oven includes stainless steel shelves with mounting brackets, connecting vacuum pipes, and wrench. Additional shelves can be ordered separately (see table).

Specifications

Temperature range: 10°C above ambient to 200°C (50 to 392°F)

Temperature resolution: 0.1°C (0.1°F)

Temperature uniformity: ±10%

Temperature stability: ±0.1°C (±1.8°F)

Temperature display: LCD

Timer range: 0 to 99 hrs, 59 mins

Vacuum range: 0 to 29.9" Hg (133 Pa)

Vacuum display: dial gauge



Capacity cu ft (L)	Rise time to 100°C	No. of shelves included	No. of shelf positions	Dimensions (W x H x D)		Power			Oven	Additional shelf
				Chamber	Overall	VAC	Hz	Watts	Item number	Item number
0.85 (24)	90 mins	2	3	11.8" x 10.8" x 11.8"	17.5" x 24.4" x 22.8"	120	60	700	ML-52411-30	ML-52411-53
				(30 x 27.5 x 30 cm)	(44.5 x 62 x 58 cm)	220	50		ML-52411-31	
1.8 (53)	80 mins	2	5	16.3" x 13.5" x 14.5"	22.8" x 27.75" x 26"	120	60	1400	ML-52411-32	ML-52411-54
				(41.5 x 34.5 x 37 cm)	(58 x 70.5 x 67.5 cm)	220	50		ML-52411-33	
3.2 (91)	80 mins	3	6	17.7" x 17.7" x 17.7"	24" x 31.9" x 31.7"	120	60	2000	ML-52411-34	ML-52411-55
				(45 x 45 x 45 cm)	(61 x 81 x 80.5 cm)	220	50		ML-52411-35	

coleparmer.com

USA:
+1.800.323.4340
+1.847.549.7600

4751

Canada: +1.800.363.5900
China: +86.21.5109.9909
France: +33 (0) 1486 37800
Germany: +49 (0) 9377 92030

India: +1.800.266.1244
Italy: +39 (0) 1313 89513
UK: +44 (0) 1480.272279
All other countries: +1.847.549.7600

Cole-Parmer is an
Antylia Scientific company
Find out more at antylia.com

