


Smart vs. Smart PRO comparison



| Controller | Smart | Smart PRO |
|------------------------|--|--|
| Display | 4.3" + touch screen | 7" + touch screen |
| Network | LAN | LAN and Wi-Fi |
| USB | YES saving registration data saving events | YES saving registration data saving events uploading programs |
| Keypad | Numeric | Alphanumeric |
| Languages | PL, EN, RU, CZ, IT, PT, UA, FR, ES | PL, EN, RU, CZ, IT, PT, UA, FR, ES |
| Main Screen | Dashboard (all relevant data visible from one main window) | Dashboard (all relevant data visible from one main window) |
| Users | - | 5 |
| Users account types | - | User / Admin / Super Admin |
| Programs | 5 | 40 |
| Program name | Free number assigned | Any |
| Priority | Parameters | Parameters, time |
| Segments | 6 | 100 |
| Light control | Only ON/OFF (FOT) | YES (FIT) |
| Schedule | - | 10 schedules |
| Data registry | max. 10,000 measurement data stored for a maximum of 6 months | max. 10,000 measurement data stored for a maximum of 12 months |
| Events registry | YES | YES |
| Statistics | YES - only the current cycle | YES - from every segment and program cycle |
| Temp. protection class | 1.0 or 2.0 (3.1, 3.2, 3.3 - option) | 3.1 or 3.2 or 3.3 |
| Quick Note | - | Ability to enter user text notes |
| Graph | - | YES |
| Mail notifications | - | Alarm notifications |
| Unit name | Fixed (serial number) | Editable |
| Alarm Bar | YES | YES |
| Quick Change | YES | YES |
| Software for PC | LabDesk (option) | LabDesk |