

## CB 210 - CO2 incubator with hot air sterilization

Precision equipment for cell cultivation of the future. It has a drift-free infrared CO2 measuring system to ensure stable pH values, condensation-free Permadry<sup>™</sup> interior chamber with integrated shelf supports, easy to clean, where the surface area is kept to a minimum in order to reduce potential sites for germs to establish; an effective sterilization program that operates at 180 °C in compliance with standards, for cultivation without any contamination. All of these features along with the absolutely precise temperature accuracy, which you have come to expect from us.



## Performance features and equipment:

- Electronically controlled APT.line<sup>™</sup> preheating chamber assuring temperature accuracy and reproducible results
- Temperature range 7 °C (13 °F) above ambient temperature up to 60 °C (140 °F)
- MCS controller for temperature and CO2 concentration
- Features:
  - User friendly LCD screen
  - · Easy-to-read menu guide
  - Integrated electronic chart recorder
  - · Variety of options for the graphic display of process parameters
  - Real-time clock
- Standard compliant 180 °C (356 °F) hot air sterilization using an overnight cycletime (acc to. DIN 58947, Pharmacopoeias, ANSI)
- VENTAIR™ Jacket System
- Drift free infrared CO2 measurement system
- Gas mixing head
- Permadry™ system condensation-free double-pan humidification system, maintains dry interior walls
- Weldless deep-drawn inner chamber made of stainless steel with integrated shelf support system
- Electronic self-diagnostic system for errors with optical and audible alarm, as well as relay contact for central monitoring
- Independent adjustable temperature safety device class 3.1, providing full protection against chamber over-temperature, with visual and audible temperature alarm
- Tightly-fitting inner glass door
- RS 422 interface for use with optional GMP/GLP and FDA guideline 21 CFR Part 11 compliant APT-COM™ DataControlSystem software
- 3 perforated shelves made of stainless steel for CB 150, CB 210 with Standard equipment and O2-Control 2 perforated shelves made of stainless steel for CB 53 with Standard equipment and O2-Control
- Units are stackable with stacking adapter
- Door Lock
- BINDER test certificate

## Technical specification CB 210

	CB 210
Exterior dimensions	CB 210
Width (mm/inch)	740 / 29.1
Height (inclusive feet) (mm/inch)	1069 / 42.1
Depth (mm/inch)	715 / 28.2
plus connections (mm/inch)	60 / 2.36
plus door handle and I-Triangle (mm/inch)	54 / 2.13
Wall clearance rear (mm/inch)	100 / 3.94
Wall clearance side (mm/inch)	50 / 2.0
Wan cicaranee side (inimitien)	0072.0
Interior dimensions	
Width (mm/inch)	560 / 22.1
Height (mm/inch)	750 / 29.5
Depth (mm/inch)	500 / 19.7
Interior volume (I/cu.ft.)	210 / 7.5
Perforated shelves, stainless steel (number standard/max.)	3 / 8
Dimensions of perforated shelves, Width (mm/inch)	555,5 / 21.9
Dimensions of perforated shelves, Depth (mm/inch)	444 / 17.5
Weight of the unit (empty) (kg/lbs.)	121 / 267
Temperature / CO2 data	
Temperature range 7 °C (13 °F) above ambient up to °C/°F	60 / 140
Temperature variation at 37 °C (98.6 °F) (± °C)	0.4
Temperature fluctuation (± °C)	0.1
Recovery time after door was opened for 30 sec. 1)	
at 37 °C (98.6 °F) (Min.)	3
CO2-range (Vol% CO2)	0-20
Setting accuracy (Vol% CO2)	0.1
Recovery time after door was opened for 30 sec. 1)	
up to 5 Vol% (min)	3
CO2-measurement	IR
Connection hose nozzle DN6 for CO2 for hose with internal diameter (mm/inch)	6 / 0.24
Humidity average value (% RH)	95
O2-range (Vol% O2)	0.2-95
Setting accuracy (Vol% O2)	0.1
Recovery time 1)	100
from 20 vol % up to 0.2 vol % O2 (Min.)	120
from 20 vol % up to 5 vol % O2 (Min.)	64
from 20 vol % up to 10 vol % O2 (Min.)	31
from 20 vol % up to 15 vol % O2 (Min.)	14
from 20 vol % up to 30 vol % O2 (Min.)	
from 20 vol % up to 50 vol % O2 (Min.)	25
from 20 vol % up to 80 vol % O2 (Min.) O2-measurement	75 ZrO2
Connection hose nozzle DN6 for O2 / N2 for hose with internal diameter (mm/inch)	6 / 0.24
Electrical data	
Housing protection acc. to EN 60529	IP 20
Nominal voltage (±10 %) 50/60 Hz (V)	230 / 115
Nominal power (W)	1500
Energy consumption at 37 °C (98.6 °F) (W)	140





1) up to 98 % of the set value All technical data are specified for units with standard equipment at an ambient temperature of 25°C and a voltage fluctuation of ±10 %. The temperature data are determinated in accordance to factory standard following DIN 12880 respecting the recommended wall clearances of 10 % of the height, width and depth of the inner chamber. All indications are average values, typical for units produced in series. We reserve the right to alter technical specifications at all times.





Gas-tight, divided inner glass door

Maintains consistent climatic conditions within the incubator. Minimal loss of incubator atmosphere, heat and, CO2 during loading, as well as fast recovery times.



Silicon access ports

Silicon access ports are provided for inserting external measuring devices. They have a diameter of 30 mm and can be sealed with a silicone plug on both sides, for positioning either on the left or the right rear side.



Stacking frame

Used for stacking two BINDER CO2 Incubators which are thermally isolated from each other.

## Options and accessories CB series



	CB 210
Silicone access port closable with two silicone lids 30 mm (1.18 inch), left side	0
Silicone access port closable with two silicone lids 30 mm (1.18 inch), right side	0
Silicone access port on both sides closable with two silicone lids 30 mm (1.18 inch), back	0
Keyboard lock	0
Gastank connection kit for CO2, consisting of a gastank pressure regulator, max. pressure 10 part, with connection parts, and a (5 m / 16.4 ft.) hose	0
Gastank connection kit for O2, consisting of a gastank pressure regulator, max. pressure 10 part, with connection parts, and a (5 m / 16.4 ft.) hose	0
Gastank connection kit for N2, consisting of a gastank pressure regulator, max. pressure 10 part, with connection parts, and a (5 m / 16.4 ft.) hose	0
BINDER Gas Supply Service: External gastank changer for connection of 2 gastanks, either CO2 or N2, with audible and visual alarm and zero-voltage alarm output	0
Ethernet interface instead of the RS 422 interface	0
4 - 20 mA analog output for temperature and CO2 measurements (e.g. chart recorder connection), with 6-pin DIN socket (output not adjustable)	0
nterior power socket with ON / OFF switch and LEMO connector (maximum power rating 230 V, 1N, 50 / 60 Hz)	0
ndependent electronic safety system Intelligent Fail Safe. Unique safety plus for continous monitoring of the CO2 control, preventing any unnoticed deviations of the CO2 concentration from the set point. Note: Not possible with options access port or divided glass door	0
Factory calibration certificate for temperature and CO2. Temperature measurement in center / CO2 measurement performed using analyzed test gas at 37 °C (98.6 °F) and 5 % CO2	0
Factory calibration certificate for the O2 control option. O2 measurement performed using analyzed test gas at 1% O2	0
Temperature measurement acc. to DIN 12880 (27 measuring points) at 37 °C (98.6 °F) or at specified temperature with measuring protocol and certificate.	0
Manual for Primary Human Cell Culture	0
Current feedthrough (8-pin) for low voltage with LEMO socket (coverable) and LEMO plug (loading capacity max. 24 V AC/DC - 2A)	
nternal CO2 gastank changer, for connecting two gas gastanks. Offers precise control over the current condition of CO2 supply to incubator by means of alarm messaging and incident reporting	0
nternal CO2 gastank changer, for connecting two gas gastanks, with external connection for up to one additional CO2 incubator. Offers precise control over the current condition of CO2 supply to incubator by means of alarm messaging and incident reporting	0
nternal O2 and N2 gastank changer for connecting two gastanks of either gas	0
Shelf, perforated, stainless steel	0
Divided shelf, stainless steel for divided gas - tight inner glass door	0
Base on castors	0
Stacking frame vibration - free, on castors with stop brake, for direct and safe stacking of 2 CB ncubators with wedge equipment.	0
Stacking adapter for direct thermal decoupled stacking of two CB 210 CO2 incubators	0
CELLROLL. Modular and expandable roller bottle system for cell cultivation, for 6 roller bottles	0



Connect With Us

