

Choose the right Hydrolab Water Purifier for your lab

Is your feed water prefiltered or straight from the tap?

I'm using pre-filtered water

Tap water

What level of purification do you need to produce?

1 (Ultrapure)

2 (Deionised)

3 (Filtered)

What level of reporting / control do you need?

What level of reporting / control do you need?

What level of reporting / control do you need?

Advanced

Basic

HLP*

Advanced

Basic

Advanced

Ultra

- Feed with Type 2 or 3
- Produces Type 1
- Type C controller for full monitoring and reporting

Spring

- Feed with tap water
- Produces Type 2 and 1
- Type C controller for full monitoring and reporting

HLP

- Feed with tap water
- Produces Type 3 and 2 as standard. Optional UV upgrade produces Type 1
- Basic Type A controller

Technical

- Feed with tap water
- Produces Type 3
- Choose basic Type A controller, or Type C for full monitoring and reporting