


Choose the right Hydrolab Water Purifier for your lab

Is your feed water prefiltered or straight from the tap?



What level of purification do you need to produce?



What level of reporting / control do you need?

What level of reporting / control do you need?

What level of reporting / control do you need?

Advanced

Basic

HLP*

Advanced

Basic

Advanced

Ultra

- Feed with Type 2 or 3
- Produces Type 1
- Type C controller for full monitoring and reporting

Spring

- Feed with tap water
- Produces Type 2 and 1
- Type C controller for full monitoring and reporting

HLP

- Feed with tap water
- Produces Type 3 and 2 as standard. Optional UV upgrade produces Type 1
- Basic Type A controller

Technical

- Feed with tap water
- Produces Type 3
- Choose basic Type A controller, or Type C for full monitoring and reporting