

Cryolab Series

High Thermal Efficiency

2 Litre to 50 Litre Capacity

Ideal for Lab Applications

5 Year Vacuum Warranty

Strong and Lightweight

Unique Pouring Trolley

LN2 Dispenser With Hose or Fixed

Dispensing Systems

Super Insulation



Cryolab Series

All Cryolabs are designed and manufactured to the very latest standards and offer strength and durability within a lightweight high-efficiency package. State of the art super-insulation and ultra-high vacuum techniques ensure unrivalled thermal efficiency and vacuum life.

Product Features



Cryolab 35 fitted with LN2 Dispenser.



Cryolab 35 fitted with Dispenser and Hose



Dispenser shelf & Mini Pouring Spout .

(Part Number 9701051 CL35/50)
(Part Number 9701050 CL25)

To be used in conjunction with a withdrawal device.



Transport trolley for the Cryolab 25.

The pouring Trolley is 1033mm high x 472mm wide.



Roller Base to suit Cryolab 20, 35, & 50 models.

Roller Base 110mm high.



Cryolab 10 Withdrawal Device. Including a swan neck dispenser.

With a range of storage capacities from 2 to 50 litres and an array of accessories designed to facilitate safe and reliable liquid nitrogen handling the Cryolab Series are ideal for laboratory applications.



Cryolab Series Technical Data

MODEL	Cryolab 2	Cryolab 5	Cryolab 10	Cryolab 20	Cryolab 25	Cryolab 35	Cryolab 50
PART NUMBER	9901015	9901030	9901040	9901050	9901062	9901080	9901090
LN2 Capacity (litres)	2	5	10	20	25	35	50
Evaporation Rate (litres/day)	0.08	0.18	0.2	0.18	0.25	0.27	0.28
Static Holding Time (days)	25 ¹	28 ¹	50 ¹	111 ¹	100 ¹	130 ¹	178 ¹
Weight Empty (Kg)	2.8	4.1	5.6	9.1	7.4	12.2	15.5
Weight Full (kg)	4.4	8.1	13.7	25.3	27.6	40.5	55.9
Neck diameter (mm)	32	51	51	51	51	63	63
Overall Height (mm)	385	495	560	570	585	650	810
Overall Diameter (mm)	185	216	305	407	395	407	407
Roller base Part Number Height 110mm	N/A	N/A	N/A	9902082	9901063	9902082	9902082
LN2 Dispenser Part Number	N/A	N/A	9701306	9701056	9701056	9701057	9701057
LN2 Dispenser with flexible hose Part Number	N/A	N/A	N/A	9701157	9701157	9701152	9701152
Pouring Trolley Part Number	N/A	N/A	N/A	N/A	9901064	N/A	N/A

Add 110mm to the height if you include a roller base.
The pouring Trolley is 1033mm high x 472mm wide.
Note: ¹Working times will vary according to container usage.

In the interest of our continued product improvement, Statebourne reserve the right to amend the specification of our products.

Dispensing Systems

We have put together a range of dispensing systems for the Cryolab Range to help you work out the best possible package to suit your application.

Please request a full set of Dispensing System pdfs



Statebourne Cryogenics
18 Parsons Road,
Parsons Industrial Estate,
Washington,
Tyne and Wear
NE37 1EZ
England.



Web: www.statebourne.com

Tel: +44 (0) 191 416 4104

Fax: +44 (0) 191 415 0369

Email: info@statebourne.com

