

COOLSAFE 55-4

SPECIFICATIONS



CoolSafe 55-4	Unit	Part no 7001000060
Digital read-out of temperature		Yes
Digital read-out of pressure		No
Ultimate temperature at 20°C room temperature	°C	-55
Power requirements	V/Hz	230/50-60
Power consumption	W	500
Insulation	mm	90
Condenser capacity	kg	3
Condenser capacity per 24 hours	kg	2,5
Total volume of Condenser	L	4
Weight net	kg	40
Cabinet dimensions (DxWxH)	mm	500x400x520
Cabinet material		Polyester coated steel
Condenser material		Stainless steel AISI 316
Cooling media		R 507
Vacuum Valve		Solenoid
Drain tap		Yes
Keypad for control options		No
R 485/USB port for Labogene software		No