



Upright Laboratory Control Plus Refrigerator 288L

GDPA code: **LSFSR288UK** Lab SF SD UP Ref

SKU code: **444410661**

- ✓ 2°C to 10°C temperature range
- ✓ Spark free interior
- ✓ External digital controller and temperature display
- ✓ Fitted lock and 2 keys
- ✓ Audible and visual temperature alarms



The Experts in
Medical Refrigeration

www.lec-medical.co.uk
☎ 0344 815 3742

1 of 3



Specification: Upright Laboratory Control Plus Refrigerator 288L

Gross Capacity (litres)	288
Net Capacity (litres)	284
Unit External Width (mm)	595mm
Unit External Depth (mm)	634mm
Unit External Height (mm)	1570mm
Product Weight (kg)	58kg
Packed Width (mm)	657mm
Packed Depth (mm)	734mm
Packed Height (mm)	1585mm
Packed Weight (kg)	63kg
Country of Origin	UK
Internal Dimensions (mm) (WxDxH)	480mm W x 1382mm H x 445mm D (step 255mm H x 175mm D)
Temperature Range	2°C to 10°C
Colour	White/green
External Finish	Painted Steel
Internal Finish	Plastic
Type of plastic	SRL800 (HIPS/ABS mix
Door Finish	Painted Steel
Door Type	Solid
Door Style	Hinged, reversible
Door Open Depth (mm)	1185mm
Replacement Door Seal	No
Open Door Sensor	Yes
Lock	Yes, 2 Keys
Adjustable Feet (number & position)	Yes - 2 front, 2 rear rollers
Porthole	No
Porthole Size	N/A
Porthole Location	N/A
Temperature Display on Panel	Yes
Temperature Alarm	Yes
Remote Alarm Contacts	Yes
Audible Alarms	Yes
Visual Alarms	Yes
Number of Shelves	5
Number of Shelf Positions	12
Shelf Size (mm) (WxDxH)	540mm W x 6mm H x 497mm D
Shelf Material	Wire
Number of Baskets	0
Basket Material	N/A
Basket Size (mm) (WxDxH)	N/A
Interior Light	No
Light Type	—
Light Activation	—
Fan	No
Spark Free	Yes
Auto or Manual Defrost	Auto
Controller Type	Digital, external
Number of Probes	1
Probe Type	Air
Climate Class	SN
Refrigerant Type	R600a
CFC Free	Yes
HFC Free	Yes

Specification: Upright Laboratory Control Plus Refrigerator 288L

HCFC Free	Yes
Ammonia Free	Yes
UK 3 Pin Plug	Yes
Other Plug Types Available	Yes
Noise Level (dB)	<45dB
Energy Consumption (kwh)	—
Voltage Supply	220 -240
Aux Output	No
Insulation Thickness (mm)	60mm
Insulation Type	Cyclopentane
Blowing Agent	Yes
Anti Microbial Handle	Yes
Optional Drawer System Available	No
Wall Brackets	No
Digital Lock	No