

Thermo Scientific RT Touch Series Magnetic Stirrers

Stay in control with our low profile digital display stirrer, ideal for routine stirring applications where visible stir speed is required.

FEATURES

- Digital control
- Gradual start-up stirring to prevent splashing
- Designed for stir rod support
- Low profile lightweight design with small footprint; two sizes available
- Strong magnetic coupling to ensure consistent stirring during experimental procedures
- Speed Control: 30 to 2000rpm
- Two non-slip silicone plate covers provided with all units (1 black, 1 white)



4

Ordering Information and Specifications

Cat. No.	88880013	88880014	88880015	88880016
Plug Configuration				
Speed Control	Digital		Digital	
Stirring Capacity (L)	4L		5L	
Speed Range	30 to 2000rpm		30 to 2000rpm	
Maximum Load (vessel)	20kg / 44.09lbs		25kg / 55.12lbs	
Body Material	Polypropylene		Polypropylene	
Top Plate Dimensions	170mm / 6.69in diameter		220mm / 8.69in diameter	
Overall Dimensions (WxDxH)	210 x 280 x 72mm 8.27 x 11.02 x 2.82in		260 x 360 x 72mm 10.24 x 14.17 x 2.82in	
Power Consumption	6W		6W	
Power	100– 240V, 50/60Hz		100 – 240V, 50/60Hz	
Certifications	CE		CE	
Warranty	2 Years		2 Years	

Accessories for RT Touch Series Magnetic

Accessories for RT Touch Series Magnetic	Cat. No.
Stirrers	88880145
Non-slip silicone plate cover for items 88880013 and 88880014 (White color, Silicone Ø170cm)	88880146

Non-slip silicone plate cover for items

88880015 and 88880016 (White color, Silicone Ø220cm)



Non-slip Silicone Plate Covers