

QUICK  **CLEAR**™

Suction Machines



GB  **UK**
HEALTHCARE

quick-clear.com



GBUK Healthcare is the UK's trusted name for medical suction equipment and consumables.

As a leading supplier of electric suction machines, we have extensive experience in providing bespoke solutions for hospital, community, homecare and commercial customers.

In addition, our comprehensive suction portfolio features market leading products such as the Serres disposable suction liner system, the Tendertip range of open and closed suction catheters and the Caretip range of open suction catheters.

QUICK **CLEAR**™

QuickClear is the new range of suction machines from GBUK Healthcare which has been designed to meet the needs of today's healthcare professionals.

Featuring the latest technology, this modern range of static and portable units has been specially constructed to combine outstanding performance, reliability and affordability.

Ideal for hospital, community, emergency, or commercial use, the QuickClear range is fully supported by GBUK Healthcare's own service department. Fully trained engineers provide a flexible on site repair and maintenance service to ensure that customer needs are speedily met.

Outstanding Performance
Reliability
Affordability



Dynamic

Compact and versatile with exceptional battery performance

Ideal for hospital, community and homecare

- ✓ Two and a half hour battery run time – reliable suctioning
- ✓ Robust, compact construction – knock proof and easy to transport
- ✓ High quality, double piston pump – powerful suction performance
- ✓ One touch on/off button with vacuum pre-set option – easy to use
- ✓ Battery, power and charging indicators – simple status check



Status indicators

Technical Specification

Vacuum Range	0 – 605 mmHg 0 – 80 kPa
Max. Flow Rate	20 L/min
Dimensions	310x240x190 mm
Weight	4.9 kg
Voltage	230V (50–60Hz)
Run Time (Battery)	up to 2.5 hours
Charge Time	4 hours
Max. Noise Level	60 dB

Order Code

Description

GBUK-11-1215(R)	QuickClear Dynamic (with 1 litre reusable canister)
GBUK-11-1215(S)	QuickClear Dynamic (with 1 litre Serres canister)
GBUK-11-1215(V)	QuickClear Dynamic (with 1 litre V canister fitting)

Vortex

**Strong, reliable,
maintenance free running**

**Ideal for multiple applications in hospital,
community and commercial settings**

- ✓ Efficient, fan cooled suction unit – reliable continuous performance
- ✓ Vibration free running at all vacuum levels – unobtrusive operation
- ✓ Heavy duty, oil free double piston pump – maintenance free
- ✓ Simple wipe clean controls – reduced contamination risk

Overflow protection



Technical Specification

Vacuum Range	0 – 680 mmHg 0 – 90 kPa
Max. Flow Rate	38 L/min
Dimensions	440x300x310 mm
Weight	12 kg
Voltage	230V (50–60Hz)
Max. Noise Level	39.4 dB

Order Code

Description

GBUK-11-1112(R)	QuickClear Vortex (with 2 litre reusable canister)
GBUK-11-1112(S)	QuickClear Vortex (with 2 litre Serres canister)
GBUK-11-1112(V)	QuickClear Vortex (with 2 litre V canister fitting)

Surgical

Ultra quiet, powerful performance

Ideal for all surgical procedures including liposuction

- ✓ High powered, fan cooled suction unit – for continuous theatre use
- ✓ Ultra quiet performance – less noise than a PC
- ✓ Stainless steel eurorail – fits any size or combination of canisters
- ✓ Integral vacuum changeover valve – limitless collection capacity

Vacuum shift



Optional shelf



Optional storage basket



Technical Specification

Vacuum Range	0 – 680 mmHg 0 – 90 kPa
Max. Flow Rate	38 L/min
Dimensions	480x900x460 mm
Voltage	230V (50–60Hz)
Max. Noise Level	39.4 dB

Order Code	Description
GBUK-11-1122(R)	QuickClear Surgical (with 2 x 2 litre reusable canisters)
GBUK-11-1122(S)	QuickClear Surgical (with 2 x 2 litre Serres canisters)
GBUK-11-1122(V)	QuickClear Surgical (with 2 x 2 litre V canister fittings)
GBUK-11-5116	Optional shelf
GBUK-11-5105	Optional storage basket

Rescue

**Slim, lightweight design
for maximum portability**

**Ideal for community, homecare
use and first line response**

- ✓ **Slim and lightweight –
easy to carry, transport and store**
- ✓ **Efficient battery performance –
fast, effective airway clearance**
- ✓ **Rugged, modern design – shock
resistant and easy to clean**
- ✓ **Includes mains charger, carry bag
and strap – ready to go in an instant**



Shoulder strap



Carry bag



Technical Specification

Vacuum Range	0 – 610 mmHg 0 – 81 kPa
Max. Flow Rate	25 L/min
Dimensions	240x125x216 mm
Weight	2.8 kg
Voltage	230V (50–60Hz)/12V DC
Battery	Ni-Mh 3.3A (10 unit pk)
Run Time (Battery)	up to 45 mins
Charge Time	4 hours

Order Code	Description
GBUK-11-90818(R)	QuickClear Rescue (with 1200ml reusable canister)
GBUK-11-90818(S)	QuickClear Rescue (with 1 litre Serres canister)
GBUK-11-90818(V)	QuickClear Rescue (with 1 litre V canister fitting)

High Suction Controller (0 – 760 mmHg)

- ✓ Lightweight, plastic construction, with high quality 'WIKA' vacuum gauge
- ✓ Mechanical overflow safety trap
- ✓ ISO 10079-3 compliant. 7 year warranty

Order Code	Description
PSR102	QuickClear High Suction Controller (with direct, BS probe)
PSR103	QuickClear High Suction Controller (rail-mounted)
PSROT	Mechanical Overflow Trap



Oxygen Flow Meter (0.15 L/Min)

- ✓ Aluminium-bodied, with durable polycarbonate inner and outer tubes
- ✓ Calibrated against 60 psi inlet pressure and 100% tested to ensure accuracy
- ✓ ISO 15002 compliant. 7 year warranty

Order Code	Description
POF101	QuickClear Oxygen Flow Meter, Single (with direct, BS probe)
POF102	QuickClear Oxygen Flow Meter, Single (rail-mounted)
POF103	QuickClear Oxygen Flow Meter, Twin (with direct, BS probe)
POF104	QuickClear Oxygen Flow Meter, Twin (rail-mounted)



Protection Filters

- ✓ In-line protection for pipeline and suction equipment
- ✓ Fit common electric suction filter mounts
- ✓ Bacteria and hydrophobic versions

Order Code	Description
GBUK-1002	QuickClear Bacteria Filter
GBUK-1003	QuickClear Hydrophobic Filter
GBUK-91418	QuickClear Rescue Hydrophobic Filter



GBUK-1002



GBUK-1003



GBUK-91418