



www.ietLtd.com Proudly serving laboratories worldwide since 1979

CALL 001.847.913.0777 for Refurbished & Certified Lab Equipment

Shimadzu SCL-10Avp

SCL-10AVP

Liquid Chromatographs



Shimadzu's SCL-10AVP System Controller directly controls all VP Series components from its simplified display menu for easy, user-friendly operation. All Shimadzu HPLC system modules are conveniently linked to the controller by fiber optic cables for swift 'Plug-and-Play' set-up. Fiber optic interfacing enables low noise, high-speed data transmission which enhances system reliability and sensitivity. An optional SCSI interface allows multiple VP Series HPLC systems to be linked.

Features

New User-Friendly Interface • Large, Easy-To-Read Characters and Icons • Customizable Display Menu and Start-Up Screen.	The new controller features many improvements including a larger, more legible display and customizable menus which simplify system operation for novice and advanced users. Menus can now be tailored to see only the parameters that you care about. Function keys, detailed help information (with valid parameter ranges), and a graphical user interface streamline parameter setup and instrument operation.
--	--

• On-Screen "Help" • Graphical User Interface	
GLP/GMP Compliance	Built-in validation support functions help you comply with GLP/GMP regulations and enhance productivity by reducing the labor required to perform validation tasks. A record of each module's operation is preserved and can be reviewed on-screen.
Transfer and Store Methods	Save method parameters to the built-in 3.5" disk drive for reference and for secure transfer to other VP Series HPLC systems. ROM version updates are downloaded from disk for easy updating.
Simple Mode for Routine Analysis	Operate your isocratic system with a simplified parameter setup. Enter pump flow rate, detector wavelength, column oven temperature and autoinjector batch schedule parameters, and begin your run.
Advanced Mode for Full Control	The SCL-10AVP's advanced mode accesses full VP Series capability, including the autoinjector's automated advanced sample pretreatment functionality. All operational parameters are controlled by the SCL-10AVP.



www.IetLtd.com Proudly serving laboratories worldwide since 1979

CALL 001.847.913.0777 for Refurbished & Certified Lab Equipment